

Item No.	Product Name	Gene Names	Uniport ID	Test Methods	Standard	Sensitivity	Assay Range
FY-EU1001	EPA(Eicosapentaenoic Acid) ELISA Kit	EPA		Competitive	20 ng/mL	0.13 ng/mL	0.31-20 ng/mL
FY-EU1002	PGE1(Prostaglandin E1) ELISA Kit	PGE1		Competitive	2000 pg/mL	9.55 pg/mL	31.25-2000 pg/mL
FY-EU1003	FLAG-Tag(DYKDDDDK-Tag Protein) ELISA Kit	FLAG-Tag		Competitive	1000 ng/mL	4.51 ng/mL	15.63-1000 ng/mL
FY-EU1004	F-E3(Free Estriol) ELISA Kit	F-E3		Competitive	40 ng/mL	0.218 ng/mL	0.63-40 ng/mL
FY-EU1005	F-TESTO(FreeTestosterone) ELISA Kit	F-TESTO		Competitive	2000 pg/mL	18.75 pg/mL	31.25-2000 pg/mL
FY-EU1006	VA(Vitamin A) ELISA KitKit	VA		Competitive	1000 ng/mL	9.38 ng/mL	15.63-1000 ng/mL
FY-EU1007	4-HNE(4-Hydroxynonenal) ELISA Kit	4-HNE		Competitive	40 ng/mL	0.38 ng/mL	0.63-40 ng/mL
FY-EU1008	GnRH(GonadotropinReleasing Hormone) ELISA Kit	GnRH		Competitive	1000 pg/mL	9.38 pg/mL	15.63-1000 pg/mL
FY-EU1009	SP(Substance P) ELISA KitKit	SP		Competitive	5000 pg/mL	46.88 pg/mL	78.13-5000 pg/mL
FY-EU1010	MEK(Methionine-Enkephalin) ELISA Kit	MEK		Competitive	8000 pg/mL	75.00 pg/mL	125.00-8000 pg/mL
FY-EU1011	LTE4(Leukotriene E4) ELISA Kit	LTE4		Competitive	2000 pg/mL	9.41 pg/mL	31.25-2000 pg/mL
FY-EU1012	10-HETE(10-Hydroxyeicosatetraenoic Acid) ELISA Kit	10-HETE		Competitive	800 ng/mL	3.95 ng/mL	12.5-800 ng/mL
FY-EU1013	TGF-β1(Transforming Growth Factor Beta 1) ELISA Kit	TGF-β1		Sandwich	2000 pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EU1014	AcCoA(Acetyl CoA) ELISA Kit	AcCoA	P14611	Sandwich	20 nmol/mL	0.119 nmol/mL	0.32-20 nmol/mL
FY-EU1015	25-OH-D(25 Hydroxy Vitamin D) ELISA Kit	25-OH-D		Competitive	200 ng/mL	0.91 ng/mL	3.13-200 ng/mL
FY-EU1016	FSH(Follicle Stimulating Hormone) ELISA Kit	FSH		Competitive	100mIU/mL	0.67 mIU/mL	1.57-100mIU/mL
FY-EH1017	Human ANG2(Angiopoietin 2) ELISA Kit	ANG2	O15123	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-EM1018	Mouse ANG2(Angiopoietin 2) ELISA Kit	ANG2	O35608	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1019	Rat ANG2(Angiopoietin 2) ELISA Kit	ANG2	O35462	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EU1020	THB(Tetrahydrobiopterin) ELISA Kit	THB		Competitive	100 ng/mL	0.51 ng/mL	1.57-100 ng/mL
FY-EU1021	Cit(Citrulline) ELISA Kit	Cit		Competitive	100 ng/mL	0.51 ng/mL	1.57-100 ng/mL
FY-EU1022	GSSG(Oxidized Glutathione) ELISA Kit	GSSG		Competitive	20000 ng/mL	91.3 ng/mL	312.5-20000 ng/mL
FY-EU1023	FK506(Tacrolimus) ELISA Kit	FK506		Competitive	10 ng/mL	0.053ng/mL	0.16-10 ng/mL
FY-EU1024	CTN(Creatinine) ELISA Kit	CTN		Competitive	200 μg/mL	0.89 μg/mL	3.13-200 μg/mL
FY-EU1025	11-KT(11-keto Testosterone) ELISA Kit	11-KT		Competitive	10000 pg/mL	15.29 pg/mL	156.25-10000 pg/mL
FY-EU1026	ACTb(Actin Beta) ELISA Kit	ACTb		Competitive	10 ng/mL	0.057 ng/mL	0.16-10 ng/mL
FY-EU1027	Lipoteichoic Acids(LTA) ELISA Kit	Lipoteichoic Acids		Competitive	20 ng/mL	0.13 ng/mL	0.31-20 ng/mL
FY-EU1028	Biotin ELISA Kit	Biotin		Competitive	200 pg/mL	4.7 pg/mL	3.13-200 pg/mL
FY-EH1029	Human BAFF/CD257(B-Cell Activating Factor) ELISA Kit	BAFF/CD257	Q9Y275	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM1030	Mouse BAFF/CD257(B-cell Activating Factor) ELISA Kit	BAFF/CD257	Q9WU72	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER1031	Rat BAFF/CD257(B-Cell Activating Factor) ELISA Kit	BAFF/CD257	Q9WU72	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EU1032	MGO(Methylglyoxal) ELISA Kit	MGO		Competitive	100 ng/mL	0.51 ng/mL	1.57-100 ng/mL
FY-EU1033	HS(Heparan Sulfate) ELISA Kit	HS;HS; HS3ST1; 3OST; 3OS		Competitive	20 ng/mL	0.19 ng/mL	0.31-20 ng/mL
FY-EU1034	LTD4(Leukotriene D4) ELISA Kit	LTD4		Sandwich	20 ng/mL	0.082 ng/mL	0.31-20 ng/mL
FY-EU1035	LTD4(Leukotriene D4) ELISA Kit	LTD4		Competitive	20 ng/mL	0.082 ng/mL	0.31-20 ng/mL
FY-EU1036	DOPAC (3,4-dihydroxyphenyl acetic acid) Elisa Kit.	DOPAC		Competitive	100 ng/mL	0.59 ng/mL	1.57-100 ng/mL
FY-EU1037	DHEA-S(Dehydroepiandrosterone sulfate) ELISA Kit	DHEA-S		Competitive	400 ng/mL	4.25 ng/mL	6.25-400 ng/mL
FY-EU1038	SB(Stercobilin) ELISA Kit	SB		Competitive	100 ng/mL	0.56 ng/mL	1.57-100 ng/mL

FY-EU1039	CH(Cholesterol) ELISA Kit	CH		Competitive	500 µg/mL	2.48 µg/mL	7.82-500 µg/mL
FY-EU1040	SMD(Spermidine) ELISA Kit	SMD		Competitive	200 ng/mL	0.87 ng/mL	3.13-200 ng/mL
FY-EU1041	VB7(Vitamin B7) ELISA Kit	VB7		Competitive	20000 pg/mL	90.2 pg/mL	312.5-20000 pg/mL
FY-EU1042	Gly(Glycine) ELISA Kit	Gly		Competitive	100 µg/mL	0.55 µg/mL	1.57-100 µg/mL
FY-EH1043	Human BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	BMP-2	P12643	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM1044	Mouse BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	BMP-2		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER1045	Rat BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	BMP-2	P49001	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EU1046	CS(Chondroitin Sulfate) ELISA Kit	CS		Competitive	10000 pg/mL	15.17 pg/mL	156.25-10000 pg/mL
FY-EU1047	LRT(Liraglutide) ELISA Kit	LRT		Competitive	1000 pg/mL	4.64 pg/mL	15.63-1000 pg/mL
FY-EU1048	S1P(Sphingosine-1-Phosphate) ELISA Kit	S1P		Competitive	1000 ng/mL	4.62 ng/mL	15.63-1000 ng/mL
FY-EU1049	CBG(Cinobufagin) ELISA Kit	CBG		Competitive	20000 ng/mL	91.3 ng/mL	312.5-20000 ng/mL
FY-EU1050	Tyr(Tyrosol) ELISA Kit	Tyr		Competitive	30000 ng/mL	142.7 ng/mL	468.75-30000 ng/mL
FY-EU1051	GTD(Gastrodin) ELISA Kit	GTD		Competitive	20000 ng/mL	87.5 ng/mL	312.5-20000 ng/mL
FY-EU1052	DC(Deoxycholate) ELISA Kit	DC		Competitive	1000 ng/mL	4.99 ng/mL	15.63-1000 ng/mL
FY-EU1053	VB2(Vitamin B2) ELISA Kit	VB2		Competitive	600 ng/mL	2.55 ng/mL	9.38-600 ng/mL
FY-EU1054	RVT(Resveratrol) ELISA Kit	RVT		Competitive	20000 ng/mL	89.3 ng/mL	312.5-20000 ng/mL
FY-EH1055	Human BMP-4(Bone Morphogenetic Protein 4) ELISA Kit	BMP-4	P12644	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1056	Mouse BMP-4(Bone Morphogenetic Protein 4) ELISA Kit	BMP-4		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1057	Rat BMP-4(Bone Morphogenetic Protein 4) ELISA Kit	BMP-4	Q06826	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EU1058	PEG(Polyethylene Glycol) ELISA Kit	PEG		Competitive	1000 ng/mL	4.72 ng/mL	15.63-1000 ng/mL
FY-EU1059	Trp(Tryptophan) ELISA Kit	Trp		Competitive	100 µg/mL	0.62 µg/mL	1.57-100 µg/mL
FY-EU1060	PCr(Phosphocreatine) ELISA Kit	PCr		Competitive	10000 ng/mL	49.3 ng/mL	156.25-10000 ng/mL
FY-EU1061	PIP3(Phosphatidylinositol Trisphosphate) ELISA Kit	PIP3		Competitive	500 ng/mL	2.83 ng/mL	7.82-500 ng/mL
FY-EU1062	GS(Guanosine) ELISA Kit	GS		Competitive	20000 pg/mL	82.7 pg/mL	312.5-20000 pg/mL
FY-EU1063	NTM(Natamycin) ELISA Kit	NTM		Competitive	300 ng/mL	1.45 ng/mL	4.69-300 ng/mL
FY-EU1064	VMA(Vanillylmandelic Acid) ELISA Kit	VMA		Competitive	10000 ng/mL	43.7 ng/mL	156.25-10000 ng/mL
FY-EU1065	PP(Dipalmitin) ELISA Kit	PP		Competitive	1000 pg/mL	4.94 pg/mL	15.63-1000 pg/mL
FY-EU1066	PI(Phosphatidylinositol) ELISA Kit	PI		Competitive	10000 pg/mL	48.6 pg/mL	156.25-10000 pg/mL
FY-EU1067	CP(Cloprostenol) ELISA Kit	CP		Competitive	10000 pg/mL	44.3 pg/mL	156.25-10000 pg/mL
FY-EH1068	Human BMP-7(Bone Morphogenetic Protein 7) ELISA Kit	BMP-7	P18075	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1069	Mouse BMP-7(Bone Morphogenetic Protein 7) ELISA Kit	BMP-7	P23359	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1070	Rat BMP-7(Bone Morphogenetic Protein 7) ELISA Kit	BMP-7		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EU1071	Glu(Glutamic Acid) ELISA Kit	Glu		Competitive	10000 ng/mL	48.2 ng/mL	156.25-10000 ng/mL
FY-EU1072	Tau(Taurine) ELISA Kit	Tau		Competitive	1000 ng/mL	4.39 ng/mL	15.63-1000 ng/mL
FY-EU1073	CNT(Carnitine) ELISA Kit	CNT		Competitive	1000 ng/mL	4.84 ng/mL	15.63-1000 ng/mL
FY-EU1074	Hep(Heparin) ELISA Kit	Hep		Competitive	1000 ng/mL	4.62 ng/mL	15.63-1000 ng/mL
FY-EU1075	LA(Lactic Acid) ELISA Kit	LA		Competitive	600 µg/mL	2.92 µg/mL	9.38-600 µg/mL
FY-EU1076	ALCAR(Acetylcarnitine) ELISA Kit	ALCAR		Competitive	5000 ng/mL	25.2 ng/mL	78.13-5000 ng/mL
FY-EU1077	TMAO(Trimethylamine-N-oxide) ELISA Kit	TMAO		Sandwich	1000 ng/mL	4.93 ng/mL	15.63-1000 ng/mL
FY-EU1078	ADP(Adenosine Diphosphate) ELISA Kit	ADP		Competitive	1500 ng/mL	7.17 ng/mL	23.44-1500 ng/mL
FY-EU1079	BA(Butyric Acid) ELISA Kit	BA		Competitive	10000 pg/mL	51.5 pg/mL	156.25-10000 pg/mL

FY-EU1080	TCA(Taurocholic Acid) ELISA Kit	TCA		Competitive	300 ng/mL	1.41 ng/mL	4.69-300 ng/mL
FY-EH1081	Human E-Cad(E-Cadherin) ELISA Kit	E-Cad	P12830	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1082	Mouse E-Cad(E-Cadherin) ELISA Kit	E-Cad	P09803	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER1083	Rat E-Cad(E-Cadherin) ELISA Kit	E-Cad	Q9ROT4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EU1084	SPH(Sphingomyelin) ELISA Kit	SPH		Competitive	400 µg/mL	1.85 µg/mL	6.25-400 µg/mL
FY-EU1085	PGH2(Prostaglandin H2) ELISA Kit	PGH2		Competitive	200 pg/mL	0.97 pg/mL	3.13-200 pg/mL
FY-EU1086	TG(Triglyceride) ELISA Kit	TG		Competitive	1000 µg/mL	1.8 µg/mL	15.63-1000 µg/mL
FY-EU1087	2-AG(2-Arachidonoylglycerol) ELISA Kit	2-AG		Competitive	300 ng/mL	1.41 ng/mL	4.69-300 ng/mL
FY-EU1088	LNG(Levonorgestrel) ELISA Kit	LNG		Competitive	100 ng/mL	0.53 ng/mL	1.57-100 ng/mL
FY-EU1089	12-HETE(12-Hydroxyeicosatetraenoic Acid) ELISA Kit	12-HETE		Competitive	1000 ng/mL	5.32 ng/mL	15.63-1000 ng/mL
FY-EU1090	5-HETE(5-Hydroxyeicosatetraenoic Acid) ELISA Kit	5-HETE		Competitive	800 ng/mL	3.95 ng/mL	12.5-800 ng/mL
FY-EU1091	PS(Phosphatidylserine) ELISA Kit	PS		Competitive	1200 ng/mL	5.27 ng/mL	18.75-1200 ng/mL
FY-EU1092	IP3(Inositol Triphosphate) ELISA Kit	IP3		Competitive	10000 pg/mL	49.3 pg/mL	156.25-10000 pg/mL
FY-EH1093	Human CA9(Carbonic Anhydrase IX) ELISA Kit	CA9	Q16790	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1094	Mouse CA9(Carbonic Anhydrase IX) ELISA Kit	CA9	Q8VHBS	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1095	Rat CA9(Carbonic Anhydrase IX) ELISA Kit	CA9		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EU1096	DAG(Diacylglycerol) ELISA Kit	DAG		Competitive	100000 pg/mL	141 pg/mL	1562.5-100000 pg/mL
FY-EU1097	LPC(Lysophosphatidylcholine) ELISA Kit	LPC		Competitive	20000 ng/mL	92.4 ng/mL	312.5-20000 ng/mL
FY-EU1098	cGMP(Cyclic Guanosine Monophosphate) ELISA Kit	cGMP		Competitive	50 pmol/mL	0.47 pmol/mL	0.78-50 pmol/mL
FY-EU1099	QA(Quinolinic Acid) ELISA Kit	QA		Competitive	100 ng/mL	0.55 ng/mL	1.57-100 ng/mL
FY-EU1100	DHA(Docosahexaenoic Acid) ELISA Kit	DHA		Competitive	1000 pg/mL	4.93 pg/mL	15.63-1000 pg/mL
FY-EU1101	KYNA(Kynurenic Acid) ELISA Kit	KYNA		Competitive	2000 ng/mL	9.24 ng/mL	31.25-2000 ng/mL
FY-EU1102	VK1(Vitamin K1) ELISA Kit	VK1		Competitive	10000 pg/mL	49.3 pg/mL	156.25-10000 pg/mL
FY-EU1103	UTP(Uridine Triphosphate) ELISA Kit	UTP		Competitive	5000 ng/mL	23.2 ng/mL	78.13-5000 ng/mL
FY-EU1104	VC(Vitamin C) ELISA Kit	VC		Competitive	40000 ng/mL	192.1 ng/mL	625-40000 ng/mL
FY-EU1105	5-HTP(5-Hydroxytryptophan) ELISA Kit	5-HTP		Competitive	400 ng/mL	1.77 ng/mL	6.25-400 ng/mL
FY-EH1106	Human CTSV(Cathepsin V) ELISA Kit	CTSV	O60911	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EU1107	HCy(Homocysteine) ELISA Kit	HCy		Competitive	8000 ng/mL	39.55 ng/mL	125-8000 ng/mL
FY-EU1108	LTC4(Leukotriene C4) ELISA Kit	LTC4		Competitive	30000 pg/mL	142.5 pg/mL	468.75-30000 pg/mL
FY-EU1109	cAMP(Cyclic Adenosine Monophosphate) ELISA Kit	cAMP		Competitive	20000 pg/mL	92.5 pg/mL	312.5-20000 pg/mL
FY-EU1110	LPA(Lysophosphatidic Acid) ELISA Kit	LPA		Competitive	10000 ng/mL	52.7 ng/mL	156.25-10000 ng/mL
FY-EU1111	SAM(S-Adenosyl Methionine) ELISA Kit	SAM		Competitive	1000 ng/mL	4.84 ng/mL	15.63-1000 ng/mL
FY-EU1112	PGI2(Prostacyclin) ELISA Kit	PGI2		Competitive	200 pg/mL	0.93 pg/mL	3.13-200 pg/mL
FY-EU1113	Bb(Bilirubin) ELISA Kit	Bb		Competitive	100 µg/mL	0.45 µg/mL	1.57-100 µg/mL
FY-EU1114	VD3(Vitamin D3) ELISA Kit	VD3		Competitive	400 ng/mL	2.14 ng/mL	6.25-400 ng/mL
FY-EU1115	HA(Hyaluronic Acid) ELISA Kit	HA		Competitive	400 ng/mL	1.81 ng/mL	6.25-400 ng/mL
FY-EU1116	DHVD3(1,25-Dihydroxyvitamin D3) ELISA Kit	DHVD3		Competitive	10000 pg/mL	49.5 pg/mL	156.25-10000 pg/mL
FY-EH1117	Human MIP-1α(Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	MIP-1α	P10147	Sandwich	1500pg/mL	14.06 pg/mL	23.44-1500pg/mL
FY-EM1118	Mouse MIP-1α(Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	MIP-1α	P10855	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1119	Rat MIP-1α(Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	MIP-1α	P50229	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EU1120	PGD2(Prostaglandin D2) ELISA Kit	PGD2		Competitive	1000 pg/mL	4.34 pg/mL	15.63-1000 pg/mL

FY-EU1121	Car(Carnosine) ELISA Kit	Car		Competitive	900 ng/mL	3.53 ng/mL	14.07-900 ng/mL
FY-EU1122	AED(Androstenediol) ELISA Kit	AED		Competitive	10000 pg/mL	53.4 pg/mL	156.25-10000 pg/mL
FY-EU1123	LTA4(Leukotriene A4) ELISA Kit	LTA4		Competitive	10000 pg/mL	46.3 pg/mL	156.25-10000 pg/mL
FY-EU1124	PAF(Platelet Activating Factor) ELISA Kit	PAF		Competitive	3000 pg/mL	14.4 pg/mL	46.88-3000 pg/mL
FY-EU1125	5-Me-THF(5-Methyltetrahydrofolate) ELISA Kit	5-Me-THF		Competitive	90 ng/mL	0.43 ng/mL	1.41-90 ng/mL
FY-EU1126	EE(Ethinylestradiol) ELISA Kit	EE		Competitive	2000 pg/mL	9.34 pg/mL	31.25-2000 pg/mL
FY-EU1127	8-epi-PGF2a(8-Epi Prostaglandin F2 Alpha) ELISA Kit	8-epi-PGF2a		Competitive	2000 pg/mL	9.55 pg/mL	31.25-2000 pg/mL
FY-EU1128	THFA(Tetrahydrofolic Acid) ELISA Kit	THFA		Competitive	100 ng/mL	0.57 ng/mL	1.57-100 ng/mL
FY-EU1129	Cor(Cortisol) ELISA Kit	Cor		Competitive	20 ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-EH1130	Human MIP-1β(Macrophage Inflammatory Protein 1 Beta) ELISA Kit	MIP-1β	P13236	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1131	Mouse MIP-1β(Macrophage Inflammatory Protein 1 Beta) ELISA Kit	MIP-1β	P14097	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1132	Rat MIP-1β(Macrophage Inflammatory Protein 1 Beta) ELISA Kit	MIP-1β	P50230	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EU1133	NAADP(Nicotinic Acid Adenine Dinucleotide Phosphate) ELISA Kit	NAADP		Competitive	10000 ng/mL	48.7 ng/mL	156.25-10000 ng/mL
FY-EU1134	VD(Vitamin D) ELISA Kit	VD		Competitive	400 ng/mL	2.19 ng/mL	6.25-400 ng/mL
FY-EU1135	DHT(Dihydrotestosterone) ELISA Kit	DHT		Competitive	2500 pg/mL	13.3 pg/mL	39.07-2500 pg/mL
FY-EU1136	Cort(Corticosterone) ELISA Kit	Cort		Competitive	100 ng/mL	0.51 ng/mL	1.57-100 ng/mL
FY-EU1137	E1S(Estrone Sulfate) ELISA Kit	E1S		Competitive	100 ng/mL	0.51 ng/mL	1.57-100 ng/mL
FY-EU1138	PTD(Pentosidine) ELISA Kit	PTD		Competitive	400 ng/mL	2.13 ng/mL	6.25-400 ng/mL
FY-EU1139	ADT(Androsterone) ELISA Kit	ADT		Competitive	90000 pg/mL	144 pg/mL	1406.25-90000 pg/mL
FY-EU1140	Arg(Arginine) ELISA Kit	Arg		Competitive	100 μg/mL	0.51 μg/mL	1.57-100 μg/mL
FY-EU1141	Pg(Progesterone) ELISA Kit	Pg		Competitive	5000 pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EU1142	HVA(Homovanillic Acid) ELISA Kit	HVA		Competitive	200 ng/mL	0.97 ng/mL	3.13-200 ng/mL
FY-EH1143	Human MCP-3(Monocyte Chemotactic Protein 3) ELISA Kit	MCP-3	P80098	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1144	Mouse MCP-3(Monocyte Chemotactic Protein 3) ELISA Kit	MCP-3	Q03366	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1145	Rat MCP-3(Monocyte Chemotactic Protein 3) ELISA Kit	MCP-3	Q9QXY8	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EU1146	VB1(Vitamin B1) ELISA Kit	VB1		Competitive	300 ng/mL	1.63 ng/mL	4.69-300 ng/mL
FY-EU1147	Ret(Retinol) ELISA Kit	Ret		Competitive	10000 ng/mL	45.1 ng/mL	156.25-10000 ng/mL
FY-EU1148	ATP(Adenosine Triphosphate) ELISA Kit	ATP		Competitive	1000 ng/mL	5.27 ng/mL	15.63-1000 ng/mL
FY-EU1149	SCT(Salmon Calcitonin) ELISA Kit	SCT		Competitive	20000 pg/mL	31.4 pg/mL	312.5-20000 pg/mL
FY-EU1150	rT3(Reverse Triiodothyronine) ELISA Kit	rT3		Competitive	6000 pg/mL	27.55 pg/mL	93.75-6000 pg/mL
FY-EU1151	CML(Carboxymethyl Lysine) ELISA Kit	CML		Competitive	5000 ng/mL	25.9 ng/mL	78.13-5000 ng/mL
FY-EU1152	ALLO(Allopregnanolone) ELISA Kit	ALLO		Competitive	100 ng/mL	0.55 ng/mL	1.57-100 ng/mL
FY-EU1153	HVD3(25-Hydroxyvitamin D3) ELISA Kit	HVD3		Competitive	200 ng/mL	0.91 ng/mL	3.13-200 ng/mL
FY-EU1154	VB6(Vitamin B6) ELISA Kit	VB6		Competitive	300 ng/mL	1.42 ng/mL	4.69-300 ng/mL
FY-EU1155	ACH(Acetylcholine) ELISA Kit	ACH		Competitive	1000 pg/mL	5.31 pg/mL	15.63-1000 pg/mL
FY-EH1156	Human ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit	ECF/CCL11	P51671	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1157	Mouse ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit	ECF/CCL11	P48298	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-ER1158	Rat ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit	ECF/CCL11	P97545	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EU1159	LTB4(Leukotriene B4) ELISA Kit	LTB4		Competitive	20000 pg/mL	89.3 pg/mL	312.5-20000 pg/mL
FY-EU1160	PGE2(Prostaglandin E2) ELISA Kit	PGE2		Competitive	200 pg/mL	0.97 pg/mL	3.13-200 pg/mL
FY-EU1161	NE(Noradrenaline) ELISA Kit	NE		Competitive	5000 pg/mL	26.1 pg/mL	78.13-5000 pg/mL

FY-EU1162	ALD(Aldosterone) ELISA Kit	ALD		Competitive	2000 pg/mL	9.41 pg/mL	31.25-2000 pg/mL
FY-EU1163	EPI(Epinephrine) ELISA Kit	EPI		Competitive	2000 pg/mL	9.74 pg/mL	31.25-2000 pg/mL
FY-EU1164	DA(Dopamine) ELISA Kit	DA		Competitive	1000 pg/mL	4.71 pg/mL	15.63-1000 pg/mL
FY-EU1165	ET(Epitestosterone) ELISA Kit	ET		Competitive	10000 pg/mL	15.29 pg/mL	156.25-10000 pg/mL
FY-EU1166	ASD(Androstenedione) ELISA Kit	ASD		Competitive	10000 pg/mL	48.8 pg/mL	156.25-10000 pg/mL
FY-EU1167	E3(Estriol) ELISA Kit	E3		Competitive	100 ng/mL	0.48 ng/mL	1.57-100 ng/mL
FY-EU1168	17-OHP(17-Hydroxyprogesterone) ELISA Kit	17-OHP		Competitive	20000 pg/mL	87.3 pg/mL	312.5-20000 pg/mL
FY-EH1169	Human TARC(Thymus Activation Regulated Chemokine) ELISA Kit	TARC	Q92583	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1170	Mouse TARC(Thymus Activation Regulated Chemokine) ELISA Kit	TARC	Q9WUZ6	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER1171	Rat TARC(Thymus Activation Regulated Chemokine) ELISA Kit	TARC	Q9ERE0	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EU1172	CNCbl(Cyanocobalamin) ELISA Kit	CNCbl		Competitive	10000 pg/mL	49.2 pg/mL	156.25-10000 pg/mL
FY-EU1173	VD2(Vitamin D2) ELISA Kit	VD2		Competitive	100 ng/mL	0.53 ng/mL	1.57-100 ng/mL
FY-EU1174	HA(Histamine) ELISA Kit	HA		Competitive	100 ng/mL	0.51 ng/mL	1.57-100 ng/mL
FY-EU1175	3-NT(3-Nitrotyrosine) ELISA Kit	3-NT		Competitive	300 ng/mL	1.35 ng/mL	4.69-300 ng/mL
FY-EU1176	TCPa(Alpha-Tocopherol) ELISA Kit	TCPa		Competitive	100 µg/mL	0.46 µg/mL	1.57-100 µg/mL
FY-EU1177	DDAVP(1-Desamino 8D Arginine Vasopressin) ELISA Kit	DDAVP		Competitive	10000 pg/mL	39.5 pg/mL	156.25-10000 pg/mL
FY-EU1178	MT(Melatonin) ELISA Kit	MT		Competitive	1000 pg/mL	4.94 pg/mL	15.63-1000 pg/mL
FY-EU1179	TXB2(Thromboxane B2) ELISA Kit	TXB2		Competitive	1500 pg/mL	7.2 pg/mL	23.44-1500 pg/mL
FY-EU1180	GA(Gibberellic Acid) ELISA Kit	GA		Competitive	10000 ng/mL	52.1 ng/mL	156.25-10000 ng/mL
FY-EH1181	Human SLC(Secondary Lymphoid Tissue Chemokine) ELISA Kit	SLC	O00585	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM1182	Mouse SLC(Secondary Lymphoid Tissue Chemokine) ELISA Kit	SLC	P84444(Ccl2	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1183	Rat SLC(Secondary Lymphoid Tissue Chemokine) ELISA Kit	SLC	A0A8I5ZU2'	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EU1184	IAA(Indole 3 Acetic Acid) ELISA Kit	IAA		Competitive	200 ng/mL	0.89 ng/mL	3.13-200 ng/mL
FY-EU1185	LXA4(Lipoxin A4) ELISA Kit	LXA4		Competitive	50 ng/mL	0.47 ng/mL	0.78-50 ng/mL
FY-EU1186	PGF2a(Prostaglandin F2 Alpha) ELISA Kit	PGF2a		Competitive	800 pg/mL	3.88 pg/mL	12.5-800 pg/mL
FY-EU1187	TXA2(Thromboxane A2) ELISA Kit	TXA2		Competitive	1000 pg/mL	4.39 pg/mL	15.63-1000 pg/mL
FY-EU1188	GSH(Glutathione) ELISA Kit	GSH		Competitive	100 ng/mL	0.48 ng/mL	1.57-100 ng/mL
FY-EU1189	ADMA(Asymmetrical Dimethylarginine) ELISA Kit	ADMA		Competitive	1000 ng/mL	4.51 ng/mL	15.63-1000 ng/mL
FY-EU1190	ABA(Abscisic Acid) ELISA Kit	ABA		Competitive	100 ng/mL	0.54 ng/mL	1.57-100 ng/mL
FY-EU1191	AA(Arachidonic Acid) ELISA Kit	AA		Competitive	200 µg/mL	0.94 µg/mL	3.13-200 µg/mL
FY-EU1192	1,5-AG(1,5-Anhydroglucitol) ELISA Kit	1,5-AG		Competitive	150 µg/mL	0.6 µg/mL	2.35-150 µg/mL
FY-EH1193	Human MPIF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	MPIF2	O00175	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM1194	Mouse MPIF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	MPIF2	Q9JKC0	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1195	Rat MPIF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	MPIF2		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EU1196	bHB(Beta-Hydroxybutyric Acid) ELISA Kit	bHB		Competitive	100 µg/mL	0.48 µg/mL	1.57-100 µg/mL
FY-EU1197	5-HIAA(5-Hydroxyindoleacetic Acid) ELISA Kit	5-HIAA		Competitive	400 ng/mL	1.93 ng/mL	6.25-400 ng/mL
FY-EU1198	E1(Estrone) ELISA Kit	E1		Competitive	2000 pg/mL	8.86 pg/mL	31.25-2000 pg/mL
FY-EU1199	Hyp(Hydroxyproline) ELISA Kit	Hyp		Competitive	5000 ng/mL	26.3 ng/mL	78.13-5000 ng/mL
FY-EU1200	FA(Folic Acid) ELISA Kit	FA		Competitive	10000 pg/mL	14.7 pg/mL	156.25-10000 pg/mL
FY-EU1201	MDA(Malondialdehyde) ELISA Kit	MDA		Competitive	2000 ng/mL	9.15 ng/mL	31.25-2000 ng/mL
FY-EU1202	5-HT(5-Hydroxytryptamine) ELISA Kit	5-HT		Competitive	900 ng/mL	3.97 ng/mL	14.07-900 ng/mL

FY-EU1203	DHEA(Dehydroepiandrosterone) ELISA Kit	DHEA		Competitive	10000 pg/mL	256 pg/mL	156.25-10000 pg/mL
FY-EU1204	ft4(Free Thyroxine) ELISA Kit	ft4		Competitive	90 pg/mL	0.48 pg/mL	1.41-90 pg/mL
FY-EU1205	ft3(Free Triiodothyronine) ELISA Kit	ft3		Competitive	100 pg/mL	0.94 pg/mL	1.56-100 pg/mL
FY-EH1206	Human MIP4α(Macrophage Inflammatory Protein 4 Alpha) ELISA Kit	MIP4α	P55774	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EU1207	AGE(Advanced Glycation End Product) ELISA Kit	AGE		Competitive	8000 ng/mL	38.2 ng/mL	125-8000 ng/mL
FY-EU1208	OT(Oxytocin) ELISA Kit	OT		Competitive	1000 pg/mL	5.51 pg/mL	15.63-1000 pg/mL
FY-EU1209	gABA(Gamma-Aminobutyric Acid) ELISA Kit	gABA		Competitive	2000 pg/mL	9.27 pg/mL	31.25-2000 pg/mL
FY-EU1210	T3(Triiodothyronine) ELISA Kit	T3		Competitive	10000 pg/mL	51.7 pg/mL	156.25-10000 pg/mL
FY-EU1211	Testo(Testosterone) ELISA Kit	Testo		Competitive	10 ng/mL	0.09 ng/mL	0.16-10ng/mL
FY-EU1212	OVA(Ovalbumin) ELISA Kit	OVA		Competitive	10000 pg/mL	7 pg/mL	156.25-10000 pg/mL
FY-EU1213	E2(Estradiol) ELISA Kit	E2		Competitive	1500 pg/mL	14.06 pg/mL	23.44-1500 pg/mL
FY-EU1214	T4(Thyroxine) ELISA Kit	T4		Competitive	300 ng/mL	1.42 ng/mL	4.69-300 ng/mL
FY-EU1215	BSA(Bovine Serum Albumin) ELISA Kit	BSA		Competitive	30000 ng/mL	161.2 ng/mL	468.75-30000 ng/mL
FY-EH1216	Human CTACK(Cutaneous T-cell Attracting Chemokine) ELISA Kit	CTACK	Q9Y4X3	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1217	Mouse CTACK(Cutaneous T-Cell Attracting Chemokine) ELISA Kit	CTACK	Q9Z1X0	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER1218	Rat CTACK(Cutaneous T-Cell Attracting Chemokine) ELISA Kit	CTACK		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EU1219	IS(serum indoxyl sulfate) ELISA Kit	IS		Competitive	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EU1220	PROA (Protein A) ELISA Kit	PROA		Competitive	3000 pg/mL	18.6 pg/mL	46.88-3000 pg/mL
FY-EU1221	PREG(Pregnenolone) ELISA Kit	PREG		Competitive	1000 pg/mL	6.12 pg/mL	15.63-1000 pg/mL
FY-EU1222	8-OHdG(8-Hydroxydeoxyguanosine) ELISA Kit	8-OHdG		Competitive	100 ng/mL	0.94 ng/mL	1.56~100 ng/mL
FY-EP1223	Porcine GSH(Glutathione) ELISA Kit	Porcine GSH		Competitive	100 ng/mL	0.48 ng/mL	1.57-100 ng/mL
FY-EP1224	Porcine MDA(Malondialdehyde) ELISA Kit	Porcine MDA		Competitive	2000 ng/mL	9.15 ng/mL	31.25-2000 ng/mL
FY-EH1225	Human FcεR II /CD23(Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit	FcεR II /CD23	P06734	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1226	Mouse FcεR II /CD23(Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit	FcεR II /CD23	P20693	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1227	Rat FcεR II /CD23(Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit	FcεR II /CD23		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EP1228	Porcine CGRP2(Calcitonin Gene Related Peptide 2) ELISA Kit	Porcine CGRP2		Sandwich	1000 pg/mL	4.94 pg/mL	15.63-1000 pg/mL
FY-EH1229	Human sCD40L(Soluble Cluster of Differentiation 40 Ligand) ELISA Kit	sCD40L	P29965	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1230	Human sCD163(Soluble Cluster of Differentiation 163) ELISA Kit	sCD163	Q86VB7	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ESH1231	Sheep NEFA(Non-Esterified Fatty Acid) ELISA Kit	NEFA		Competitive	20 ng/mL	0.152 ng/mL	0.32-20 ng/mL
FY-EH1232	Human GP39(Glycoprotein 39, Cartilage) ELISA Kit	GP39	P36222	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1233	Human CNTF(Ciliary Neurotrophic Factor) ELISA Kit	CNTF	P26441	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1234	Mouse CNTF(Ciliary Neurotrophic Factor) ELISA Kit	CNTF	P51642	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ER1235	Rat CNTF(Ciliary Neurotrophic Factor) ELISA Kit	CNTF	P20294	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ESH1236	Sheep NKB(Neurokinin B) ELISA Kit	NKB		Competitive	800 pg/mL	3.17 pg/mL	12.5-800 pg/mL
FY-EH1237	Human F3/TF(Tissue Factor) ELISA Kit	Human F3/TF		Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH1238	Human CRN(Corin) ELISA Kit	CRN	Q9Y5Q5	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1239	Human CX3CL1(CheMokine C-X3-C-Motif Ligand 1) ELISA Kit	CX3CL1	P78423	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1240	Mouse CX3CL1(CheMokine C-X3-C-Motif Ligand 1) ELISA Kit	CX3CL1	O35188	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1241	Rat CX3CL1(CheMokine C-X3-C-Motif Ligand 1) ELISA Kit	CX3CL1	O55145	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1242	Human GROα/CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit	GROα/CXCL1	P09341	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1243	Mouse GROα/CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit	GROα/CXCL1	P12850	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-ER1244	Rat GRO α /CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit	GRO α /CXCL1	P14095	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1245	Human GCP-2(Granulocyte Chemotactic Protein 2) ELISA Kit	GCP-2	P80162	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1246	Mouse GCP-2(Granulocyte Chemotactic Protein 2) ELISA Kit	GCP-2	P50228	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1247	Human IP-10/CXCL10(Interferon Gamma Induced Protein 10kDa) ELISA Kit	IP-10/CXCL10	P02778	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM1248	Mouse IP-10/CXCL10(Interferon Gamma Induced Protein 10kDa) ELISA Kit	IP-10/CXCL10	P17515	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1249	Rat IP-10/CXCL10(Interferon Gamma Induced Protein 10kDa) ELISA Kit	IP-10/CXCL10	P48973	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1250	Human I-TAC/CXCL11(Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit	I-TAC/CXCL11	O14625	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM1251	Mouse I-TAC/CXCL11(Interferon Inducible T-Cell Alpha Chemoattractant) ELISA Kit	I-TAC/CXCL11	Q9JHH5	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER1252	Rat I-TAC/CXCL11(Interferon Inducible T-Cell Alpha Chemoattractant) ELISA Kit	I-TAC/CXCL11		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EGP1253	Guinea Porcine ADMA(Asymmetrical Dimethylarginine) ELISA Kit	Guinea Porcine ADMA		Competitive	1000 ng/mL	4.51 ng/mL	15.63-1000 ng/mL
FY-EG1254	Goat vWF(Von Willebrand Factor) ELISA Kit	Goat vWF		Competitive	800 ng/mL	4.11 ng/mL	12.5-800 ng/mL
FY-EG1255	Goat ALB(Albumin) ELISA Kit	Goat ALB		Competitive	100 μ g/mL	0.59 μ g/mL	1.57-100 μ g/mL
FY-EH1256	Human BLC(B-Lymphocyte Chemoattractant) ELISA Kit	BLC	O43927	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1257	Mouse BLC(B-Lymphocyte Chemoattractant) ELISA Kit	BLC	O55038	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1258	Rat BLC(B-Lymphocyte Chemoattractant) ELISA Kit	BLC		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EZe1259	Zebrafish 5-HT(5-Hydroxytryptamine) ELISA Kit	Zebrafish 5-HT		Competitive	900 ng/mL	3.97 ng/mL	14.07-900 ng/mL
FY-EG1260	Goat NEFA(Non-Esterified Fatty Acid) ELISA Kit	Goat NEFA		Competitive	20 ng/mL	0.152 ng/mL	0.32-20 ng/mL
FY-EZe1261	Zebrafish GH(Growth Hormone) ELISA Kit	Zebrafish GH		Competitive	10000 pg/mL	55.7 pg/mL	156.25-10000 pg/mL
FY-EP1262	Porcine GnRH(Gonadotropin Releasing Hormone) ELISA Kit	Porcine GnRH		Competitive	1000 pg/mL	5.25 pg/mL	15.63-1000 pg/mL
FY-EH1263	Human Cyt-C(Cytochrome C) ELISA Kit	Cyt-C	P99999	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM1264	Mouse Cyt-C(Cytochrome C) ELISA Kit	Cyt-C	P62897	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER1265	Rat Cyt-C(Cytochrome C) ELISA Kit	Cyt-C		Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-EH1266	Human DKK1(Dickkopf Related Protein 1) ELISA Kit	DKK1	O94907	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1267	Mouse DKK1(Dickkopf Related Protein 1) ELISA Kit	DKK1	O54908	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER1268	Rat DKK1(Dickkopf Related Protein 1) ELISA Kit	DKK1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ESh1269	Sheep vWF(Von Willebrand Factor) ELISA Kit	vWF		Competitive	800 ng/mL	4.11 ng/mL	12.5-800 ng/mL
FY-EH1270	Human EGF(Epidermal Growth Factor) ELISA Kit	EGF	P01133	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-EM1271	Mouse EGF(Epidermal Growth Factor) ELISA Kit	EGF	P01132	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1272	Rat EGF(Epidermal Growth Factor) ELISA Kit	EGF	P01133	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EZe1273	Zebrafish DA(Dopamine) ELISA Kit	Zebrafish DA		Competitive	1000 pg/mL	4.71 pg/mL	15.63-1000 pg/mL
FY-ECH1274	Chicken 8-iso-PGF2 α (8-isoprostane) ELISA Kit	Chicken 8-iso-PGF2 α		Competitive	1000 pg/mL	9.375 pg/mL	15.625-1000 pg/mL
FY-EH1275	Human EGFR(Epidermal Growth Factor Receptor) ELISA Kit	EGFR	P00533	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1276	Mouse EGFR(Epidermal Growth Factor Receptor) ELISA Kit	EGFR	Q01279	Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-EH1277	Human EG-VEGF(Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit	EG-VEGF	P58294	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM1278	Mouse EG-VEGF(Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit	EG-VEGF	Q14A28	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1279	Rat EG-VEGF(Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit	EG-VEGF	Q8R414	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EP1280	Porcine CTXII(Cross Linked C-Telopeptide Of Type II Collagen) ELISA Kit	Porcine CTXII		Competitive	10000 pg/mL	52.3 pg/mL	156.25-10000 pg/mL
FY-EH1281	Human ET-1(Endothelin 1) ELISA Kit	ET-1	P05305	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-EM1282	Mouse ET-1(Endothelin 1) ELISA Kit	ET-1	P22387	Sandwich	50pg/mL	0.47 pg/mL	0.79-50pg/mL
FY-ER1283	Rat ET-1(Endothelin 1) ELISA Kit	ET-1	P22388	Sandwich	50pg/mL	0.47 pg/mL	0.79-50pg/mL
FY-ECH1284	Chicken AOPP(Advanced Oxidation Protein Products) ELISA Kit	Chicken AOPP		Competitive	10000 pg/mL	49.6 pg/mL	156.25-10000 pg/mL

FY-EH1285	Human EPCR(Endothelial Protein C Receptor) ELISA Kit	EPCR	Q9UNN8	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER1286	Rat EPCR(Endothelial Protein C Receptor) ELISA Kit	EPCR		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1287	Human FAS/CD95(Factor Related Apoptosis) ELISA Kit	FAS/CD95	P25445	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1288	Mouse FAS/CD95(Factor Related Apoptosis) ELISA Kit	FAS/CD95	P25446	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1289	Rat FAS/CD95(Factor Related Apoptosis) ELISA Kit	FAS/CD95	Q63199	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EP1290	Porcine CTXI(Cross Linked C-Telopeptide Of Type I Collagen) ELISA Kit	Porcine CTXI		Competitive	10000 pg/mL	52.9 pg/mL	156.25-10000 pg/mL
FY-ECH1291	Chicken IgY(Immunoglobulin of Yolk) ELISA Kit	Chicken IgY		Competitive	400 µg/mL	1.85 µg/mL	6.25-400 µg/mL
FY-EH1292	Human FASL/TNFSF6(Factor Related Apoptosis Ligand) ELISA Kit	FASL/TNFSF6	P48023	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1293	Mouse FASL/TNFSF6(Factor Related Apoptosis Ligand) ELISA Kit	FASL/TNFSF6	P41047	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1294	Rat FASL/TNFSF6(Factor Related Apoptosis Ligand) ELISA Kit	FASL/TNFSF6	P36940	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1295	Human AFGF/FGF1(Acidic Fibroblast Growth Factor 1) ELISA Kit	AFGF/FGF1	P05230	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1296	Mouse AFGF/FGF1(Acidic Fibroblast Growth Factor 1) ELISA Kit	AFGF/FGF1	P61148	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1297	Rat AFGF/FGF1(Acidic Fibroblast Growth Factor 1) ELISA Kit	AFGF/FGF1		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EZe1298	Zebrafish FSH(Follicle Stimulating Hormone) ELISA Kit	Zebrafish FSH		Competitive	180 ng/mL	0.87 ng/mL	2.82-180 ng/mL
FY-EZe1299	Zebrafish GnRH(Gonadotropin Releasing Hormone) ELISA Kit	Zebrafish GnRH		Competitive	16 ng/mL	0.1 ng/mL	0.25-16 ng/mL
FY-ECH1300	Chicken Neopterin ELISA Kit	Chicken Neopterin		Competitive	10 ng/mL	0.053ng/mL	0.16-10 ng/mL
FY-EG1301	Goat Testo(Testosterone) ELISA Kit	Goat Testo		Competitive	4000 pg/mL	17.7 pg/mL	62.5-4000 pg/mL
FY-EG1302	Goat DHEA(Dehydroepiandrosterone) ELISA Kit	Goat DHEA		Competitive	10000 pg/mL	256 pg/mL	156.25-10000 pg/mL
FY-EH1303	Human FGF19(Fibroblast Growth Factor 19) ELISA Kit	FGF19	O95750	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1304	Rat FGF19(Fibroblast Growth Factor 19) ELISA Kit	FGF19	Q8V181	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EP1305	Porcine Cor(Cortisol) ELISA Kit	Porcine Cor		Competitive	200 ng/mL	0.93 ng/mL	3.13-200 ng/mL
FY-EP1306	Porcine GAGs(Glycosaminoglycan) ELISA Kit	Porcine GAGs		Competitive	100 ng/mL	0.51 ng/mL	1.57-100 ng/mL
FY-ECH1307	Chicken CTXI(Cross Linked C-Telopeptide Of Type I Collagen) ELISA Kit	Chicken CTXI		Competitive	10000 pg/mL	52.9 pg/mL	156.25-10000 pg/mL
FY-ECH1308	Chicken T4(Thyroxine) ELISA Kit	Chicken T4		Competitive	300 ng/mL	1.42 ng/mL	4.69-300 ng/mL
FY-ECH1309	Chicken NE(Noradrenaline) ELISA Kit	Chicken NE		Competitive	5000 pg/mL	26.1 pg/mL	78.13-5000 pg/mL
FY-EH1310	Human FGF21(Fibroblast Growth Factor 21) ELISA Kit	FGF21	Q9NSA1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1311	Mouse FGF21(Fibroblast Growth Factor 21) ELISA Kit	FGF21	Q9JJN1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1312	Rat FGF21(Fibroblast Growth Factor 21) ELISA Kit	FGF21		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ECH1313	Chicken HA(Histamine) ELISA Kit	Chicken HA		Competitive	100 ng/mL	0.51 ng/mL	1.57-100 ng/mL
FY-ECH1314	Chicken 5-HT(5-Hydroxytryptamine) ELISA Kit	Chicken 5-HT		Competitive	900 ng/mL	3.97 ng/mL	14.07-900 ng/mL
FY-ECH1315	Chicken DA(Dopamine) ELISA Kit	Chicken DA		Competitive	1000 pg/mL	4.71 pg/mL	15.63-1000 pg/mL
FY-ECH1316	Chicken T3(Triiodothyronine) ELISA Kit	Chicken T3		Competitive	10000 pg/mL	51.7 pg/mL	156.25-10000 pg/mL
FY-EP1317	Porcine NEFA(Non-Esterified Fatty Acid) ELISA Kit	Porcine NEFA		Competitive	20 ng/mL	0.152 ng/mL	0.32-20 ng/mL
FY-ECH1318	Chicken Cort(Corticosterone) ELISA Kit	Chicken Cort		Competitive	100 ng/mL	0.51 ng/mL	1.57-100 ng/mL
FY-EZe1319	Zebrafish LH(Luteinizing Hormone) ELISA Kit	Zebrafish LH		Competitive	30000 pg/mL	145.2 pg/mL	468.75-30000 pg/mL
FY-EZe1320	Zebrafish Cor(Cortisol) ELISA Kit	Zebrafish Cor		Competitive	200 ng/mL	0.93 ng/mL	3.13-200 ng/mL
FY-ESh1321	Sheep PG(Progesterone) ELISA Kit	PG		Competitive	100 ng/mL	0.47 ng/mL	1.57-100 ng/mL
FY-EP1322	Porcine IP3(Inositol Triphosphate) ELISA Kit	Porcine IP3		Competitive	10000 pg/mL	49.7 pg/mL	156.25-10000 pg/mL
FY-EH1323	Human FST(Follistatin) ELISA Kit	FST	P19883	Sandwich	16ng/mL	0.15 ng/mL	0.25-16ng/mL
FY-EM1324	Mouse FST(Follistatin) ELISA Kit	FST	P47931	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-ER1325	Rat FST(Follistatin) ELISA Kit	FST	P21674	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EZe1326	Zebrafish MDA(Malondialdehyde) ELISA Kit	Zebrafish MDA	Competitive	1000 pg/mL	5.39 pg/mL	15.63-1000 pg/mL	
FY-EZe1327	Zebrafish GSH(Glutathione) ELISA Kit	Zebrafish GSH	Competitive	100 ng/mL	0.48 ng/mL	1.57-100 ng/mL	
FY-ESh1328	Sheep FSH(Follicle stimulating hormone) ELISA Kit	FSH	Competitive	200 ng/mL	1.032 ng/mL	3.13-200 ng/mL	
FY-ESh1329	Sheep LH(Luteinizing Hormone) ELISA Kit	LH	Competitive	30000 pg/mL	163.2 pg/mL	468.75-30000 pg/mL	
FY-ESh1330	Sheep E2(Estradiol) ELISA Kit	E2	Competitive	1000 pg/mL	4.47 pg/mL	15.63-1000 pg/mL	
FY-ECH1331	Chicken MDA(Malondialdehyde) ELISA Kit	Chicken MDA	Competitive	1000 pg/mL	5.39 pg/mL	15.63-1000 pg/mL	
FY-ESh1332	Sheep Cor(Cortisol) ELISA Kit	Cor	Competitive	200 ng/mL	0.93 ng/mL	3.13-200 ng/mL	
FY-EH1333	Human GDF15(Growth Differentiation Factor 15) ELISA Kit	GDF15	Q99988	Sandwich	1500pg/mL	14.06 pg/mL	23.44-1500pg/mL
FY-EM1334	Mouse GDF15(Growth Differentiation Factor 15) ELISA Kit	GDF15	Q9Z0J7	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1335	Rat GDF15(Growth Differentiation Factor 15) ELISA Kit	GDF15	Q9Z0J6	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EG1336	Goat Total cholesterol (TC) ELISA Kit	Goat Total cholesterol	Competitive	500 µg/mL	2.63 µg/mL	7.82-500 µg/mL	
FY-EZe1337	Zebrafish Testo(Testosterone) ELISA Kit	Zebrafish Testo	Competitive	4000 pg/mL	17.7 pg/mL	62.5-4000 pg/mL	
FY-EZe1338	Zebrafish E2(Estradiol) ELISA Kit	Zebrafish E2	Competitive	1000 pg/mL	4.45 pg/mL	15.63-1000 pg/mL	
FY-EZe1339	Zebrafish VTG(Vitellogenin) ELISA Kit	Zebrafish VTG	Competitive	2000 pg/mL	9.34 pg/mL	31.25-2000 pg/mL	
FY-ECH1340	Chicken HB(Hemoglobin) ELISA Kit	Chicken HB	Competitive	100 µg/mL	0.53 µg/mL	1.57-100 µg/mL	
FY-EH1341	Human HGF(Hepatocyte Growth Factor) ELISA Kit	HGF	P14210	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EM1342	Mouse HGF(Hepatocyte Growth Factor) ELISA Kit	HGF	Q08048	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1343	Rat HGF(Hepatocyte Growth Factor) ELISA Kit	HGF	P17945	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ECH1344	Chicken Bb(Bilirubin) ELISA Kit	Chicken Bb	Competitive	100 µg/mL	0.45 µg/mL	1.57-100 µg/mL	
FY-ECH1345	Chicken ALB(Albumin) ELISA Kit	Chicken ALB	Competitive	100 µg/mL	0.59 µg/mL	1.57-100 µg/mL	
FY-ECH1346	Chicken Cor(Cortisol) ELISA Kit	Chicken Cor	Competitive	200 ng/mL	0.93 ng/mL	3.13-200 ng/mL	
FY-ESh1347	Sheep IgG(Immunoglobulin G) ELISA Kit	IgG	Competitive	100 µg/mL	0.55 µg/mL	1.57-100 µg/mL	
FY-EG1348	Goat MDA(Malondialdehyde) ELISA Kit	Goat MDA	Competitive	1000 pg/mL	5.39 pg/mL	15.63-1000 pg/mL	
FY-EP1349	Porcine ASD(Androstenedione) ELISA Kit	Porcine ASD	Competitive	10000 pg/mL	48.8 pg/mL	156.25-10000 pg/mL	
FY-EP1350	Porcine AED(Androstenediol) ELISA Kit	Porcine AED	Competitive	10000 pg/mL	53.4 pg/mL	156.25-10000 pg/mL	
FY-EP1351	Porcine T4(Thyroxine) ELISA Kit	Porcine T4	Competitive	300 ng/mL	1.42 ng/mL	4.69-300 ng/mL	
FY-EP1352	Porcine T3(Triiodothyronine) ELISA Kit	Porcine T3	Competitive	10000 pg/mL	51.7 pg/mL	156.25-10000 pg/mL	
FY-EH1353	Human IFN-β(Interferon Beta) ELISA Kit	IFN-β	P01574	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1354	Mouse IFN-β(Interferon Beta) ELISA Kit	IFN-β	P01575	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1355	Rat IFN-β(Interferon Beta) ELISA Kit	IFN-β	P70499	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1356	Human IL-1α(Interleukin 1 Alpha) ELISA Kit	IL-1α	P01583	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-EM1357	Mouse IL-1α(Interleukin 1 Alpha) ELISA Kit	IL-1α	P01582	Sandwich	50pg/mL	0.47 pg/mL	0.79-50pg/mL
FY-ER1358	Rat IL-1α(Interleukin 1 Alpha) ELISA Kit	IL-1α	P16598	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ECH1359	Chicken HA(Hyaluronic Acid) ELISA Kit	Chicken HA	Competitive	400 ng/mL	1.81 ng/mL	6.25-400 ng/mL	
FY-EP1360	Porcine ADST(Androstenone) ELISA Kit	Porcine ADST	Competitive	2000nmol/ml	10.19nmol/ml	31.25-2000 nmol/ml	
FY-EP1361	Porcine ADT(Androsterone) ELISA Kit	Porcine ADT	Competitive	90000 pg/mL	144 pg/mL	1406.25-90000 pg/mL	
FY-EH1362	Human IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit	IL-1ra/IL-1F3	P18510	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1363	Mouse IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit	IL-1ra/IL-1F3	P25085	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1364	Rat IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit	IL-1ra/IL-1F3	P25086	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EP1365	Porcine Cort(Corticosterone) ELISA Kit	Porcine Cort	Competitive	100 ng/mL	0.51 ng/mL	1.57-100 ng/mL	
FY-ECH1366	Chicken GnRH(Gonadotropin Releasing Hormone) ELISA Kit	Chicken GnRH	Competitive	16 ng/mL	0.1 ng/mL	0.25-16 ng/mL	

FY-EH1367	Human IL-1R2(Interleukin 1 Receptor Type II) ELISA Kit	IL-1R2	P27930	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1368	Mouse IL-1R2(Interleukin 1 Receptor Type II) ELISA Kit	IL-1R2	P27931	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER1369	Rat IL-1R2(Interleukin 1 Receptor Type II) ELISA Kit	IL-1R2	P43303	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ECH1370	Chicken PGE2(Prostaglandin E2) ELISA Kit	Chicken PGE2		Competitive	25 ng/mL	0.1 ng/mL	0.4-25 ng/mL
FY-ECH1371	Chicken fT4(Free Thyroxine) ELISA Kit	Chicken fT4		Competitive	8 ng/mL	0.1 ng/mL	0.13-8 ng/mL
FY-EH1372	Human LIF(Leukemia Inhibitory Factor) ELISA Kit	LIF	P15018	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1373	Mouse LIF(Leukemia Inhibitory Factor) ELISA Kit	LIF	P09056	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1374	Rat LIF(Leukemia Inhibitory Factor) ELISA Kit	LIF		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ECH1375	Chicken fT3(Free Triiodothyronine) ELISA Kit	Chicken fT3		Competitive	8 ng/mL	0.1 ng/mL	0.13-8 ng/mL
FY-ECH1376	Chicken PG(Progesterone) ELISA Kit	Chicken PG		Competitive	100 ng/mL	0.47 ng/mL	1.57-100 ng/mL
FY-ECH1377	Chicken LH(Luteinizing Hormone) ELISA Kit	Chicken LH		Competitive	20 ng/mL	0.1 ng/mL	0.32-20 ng/mL
FY-ECH1378	Chicken FSH(Follicle Stimulating Hormone) ELISA Kit	Chicken FSH		Competitive	25 ng/mL	0.1 ng/mL	0.4-25 ng/mL
FY-EP1379	Porcine OB(Obestatin) ELISA Kit	Porcine OB		Competitive	5000 pg/mL	25.1 pg/mL	78.13-5000 pg/mL
FY-EP1380	Porcine AP17(Apelin 17) ELISA Kit	Porcine AP17		Competitive	2000 pg/mL	9.14 pg/mL	31.25-2000 pg/mL
FY-EH1381	Human TNFSF14(Tumor Necrosis Factor Ligand Superfamily, Member 14) ELISA Kit	TNFSF14	O43557	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1382	Mouse TNFSF14(Tumor Necrosis Factor Ligand Superfamily, Member 14) ELISA Kit	TNFSF14	Q9QYH9	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1383	Human MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit	MCSF	P09603	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1384	Mouse MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit	MCSF	P07141	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1385	Rat MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit	MCSF	Q8JZQ0	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1386	Human Pro-MMP-1(pro-Matrix Metalloproteinase 1) ELISA Kit	Pro-MMP-1	P03956	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1387	Human IL-3(Interleukin 3) ELISA Kit	IL-3	P08700	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1388	Mouse IL-3(Interleukin 3) ELISA Kit	IL-3	P01586	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1389	Rat IL-3(Interleukin 3) ELISA Kit	IL-3	P04823	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1390	Human IL-4(Interleukin 4) ELISA Kit	IL-4	P05112	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1391	Mouse IL-4(Interleukin 4) ELISA Kit	IL-4	P07750	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1392	Rat IL-4(Interleukin 4) ELISA Kit	IL-4	P20096	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1393	Human IL-13(Interleukin 13) ELISA Kit	IL-13	P35225	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1394	Mouse IL-13(Interleukin 13) ELISA Kit	IL-13	P20109	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1395	Rat IL-13(Interleukin 13) ELISA Kit	IL-13	P42203	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1396	Human IL-22(Interleukin 22) ELISA Kit	IL-22	Q9GZX6	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1397	Mouse IL-22(Interleukin 22) ELISA Kit	IL-22	Q9JJY9	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1398	Rat IL-22(Interleukin 22) ELISA Kit	IL-22	Q198B4	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1399	Human IL-23(Interleukin 23) ELISA Kit	IL-23	Q9GZX6	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-EM1400	Mouse IL-23(Interleukin 23) ELISA Kit	IL-23	Q9EQ14	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1401	Rat IL-23(Interleukin 23) ELISA Kit	IL-23	Q91Z84	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1402	Human IFN- γ (Interferon Gamma) ELISA Kit	IFN- γ	P01579	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1403	Mouse IFN- γ (Interferon Gamma) ELISA Kit	IFN- γ	P01580	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1404	Rat IFN- γ (Interferon Gamma) ELISA Kit	IFN- γ	P01581	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1405	Human TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	TNF- α	P01375	Sandwich	500pg/mL	3.91 pg/mL	7.81-500pg/mL
FY-EM1406	Mouse TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	Mouse TNF- α	P06804	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER1407	Rat TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	TNF- α	P16599	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-EH1408	Human VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit	VEGF-A	P15692	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1409	Mouse VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit	VEGF-A	Q00731	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1410	Rat VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit	VEGF-A	P16612	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1411	Human EREG(Epiregulin) ELISA Kit	EREG	O14944	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1412	Mouse EREG(Epiregulin) ELISA Kit	EREG	Q61521	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1413	Rat EREG(Epiregulin) ELISA Kit	EREG	A0A8I5ZT8:	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1414	Human Bcl-2(B-cell Leukemia/Lymphoma 2) ELISA Kit	Bcl-2	P10415	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1415	Mouse Bcl-2(B-cell Lymphoma/Leukemia 2) ELISA Kit	Bcl-2	P10417	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1416	Rat Bcl-2(B-cell Lymphoma/Leukemia 2) ELISA Kit	Bcl-2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1417	Human REN(Renin) ELISA Kit	REN	P00797	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1418	Mouse REN(Renin) ELISA Kit	REN	P06281	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1419	Rat REN(Renin) ELISA Kit	REN	P08424	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EP1420	Porcine BNP(Brain Natriuretic Peptide) ELISA Kit	Porcine BNP		Competitive	2000 pg/mL	9.27 pg/mL	31.25-2000 pg/mL
FY-EH1421	Human CTSC(Cathepsin C) ELISA Kit	CTSC	P53634	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1422	Human ApoA1(Apolipoprotein A1) ELISA Kit	ApoA1	P02647	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM1423	Mouse ApoA1(Apolipoprotein A1) ELISA Kit	ApoA1	Q00623	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER1424	Rat ApoA1(Apolipoprotein A1) ELISA Kit	ApoA1	P04639	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1425	Human MMP-25(Matrix Metalloproteinase 25) ELISA Kit	MMP-25	Q9NPA2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1426	Human MMP-12(Matrix Metalloproteinase 12) ELISA Kit	MMP-12	P39900	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1427	Mouse MMP-12(Matrix Metalloproteinase 12) ELISA Kit	MMP-12	P34960	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1428	Rat MMP-12(Matrix Metalloproteinase 12) ELISA Kit	MMP-12	Q63341	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1429	Human CT(Calcitonin) ELISA Kit	CT	P01258	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1430	Mouse CT(Calcitonin) ELISA Kit	CT	P70160	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1431	Rat CT(Calcitonin) ELISA Kit	CT	P01257	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EP1432	Porcine Anti-SST(Anti-Somatostatin Antibody) ELISA Kit	Porcine Anti-SST		Competitive	200 ng/mL	1.24 ng/mL	3.13-200 ng/mL
FY-EH1433	Human ACTH(Adrenocorticotrophic Hormone) ELISA Kit	ACTH	P63267	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1434	Mouse ACTH(Adrenocorticotrophic Hormone) ELISA Kit	ACTH		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1435	Rat ACTH(Adrenocorticotrophic Hormone) ELISA Kit	ACTH		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1436	Human PRL(Prolactin) ELISA Kit	PRL	P01236	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1437	Mouse PRL(Prolactin) ELISA Kit	PRL	P06879	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER1438	Rat PRL(Prolactin) ELISA Kit	PRL	P01237	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH1439	Human PAPP(A) (Pregnancy Associated Plasma Protein A) ELISA Kit	PAPP(A)	Q13219	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM1440	Mouse PAPP(A) (Pregnancy Associated Plasma Protein A) ELISA Kit	PAPP(A)	Q8R4K8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1441	Rat PAPP(A) (Pregnancy Associated Plasma Protein A) ELISA Kit	PAPP(A)		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1442	Human CSTA(Cystatin A) ELISA Kit	CSTA	P01040	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1443	Mouse CSTA(Cystatin A) ELISA Kit	CSTA	P56567	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1444	Human IL-1 β (Interleukin 1 Beta) ELISA Kit	IL-1 β	P01584	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM1445	Mouse IL-1 β (Interleukin 1 Beta) ELISA Kit	IL-1 β	P10749	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1446	Rat IL-1 β (Interleukin 1 Beta) ELISA Kit	IL-1 β	Q63264	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1447	Human IL-12(Interleukin 12) ELISA Kit	IL-12	P29459(IL12	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1448	Mouse IL-12(Interleukin 12) ELISA Kit	IL-12	P43432	Sandwich	50pg/mL	0.47 pg/mL	0.79-50pg/mL

FY-ER1449	Rat IL-12(Interleukin 12) ELISA Kit	IL-12	Q9R103, E1	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1450	Human IL-12 p40(Interleukin 12 p40) ELISA Kit	IL-12 p40	P29460	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1451	Mouse IL-12 p40(Interleukin 12 p40) ELISA Kit	IL-12 p40	P43432	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1452	Rat IL-12 p40(Interleukin 12 p40) ELISA Kit	IL-12 p40	Q62712	Sandwich	800pg/mL	7.5 pg/mL	12.5-800pg/mL
FY-EH1453	Human α -CTx(Alpha Crosslaps) ELISA Kit	α -CTx	P50983	Competitive	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EM1454	Mouse α -CTx(Alpha Crosslaps) ELISA Kit	α -CTx		Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER1455	Rat α -CTx(Alpha Crosslaps) ELISA Kit	α -CTx		Competitive	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EH1456	Human CD42(Cluster of Differentiation 42) ELISA Kit	CD42	P14770(GP9)	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH1457	Human HLA-B27(Human Leukocyte Antigen B27) ELISA Kit	HLA-B27	P01889	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1458	Human FGFR2(Fibroblast Growth Factor Receptor 2) ELISA Kit	FGFR2	P21802	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1459	Human MBP(Myelin Basic Protein) ELISA Kit	MBP	P02686	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1460	Mouse MBP(Myelin Basic Protein) ELISA Kit	MBP	P04370	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1461	Rat MBP(Myelin Basic Protein) ELISA Kit	MBP	P02688	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1462	Human TE(Telomerase) ELISA Kit	TE	O14746	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1463	Mouse TE(telomerase) ELISA Kit	TE		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-ER1464	Rat TE(Telomerase) ELISA Kit	TE		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1465	Human TFPI(Tissue Factor Pathway Inhibitor) ELISA Kit	TFPI	P10646	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1466	Mouse TFPI(Tissue Factor Pathway Inhibitor) ELISA Kit	TFPI		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER1467	Rat TFPI(Tissue Factor Pathway Inhibitor) ELISA Kit	TFPI	Q02445	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1468	Human TM(Thrombin) ELISA Kit	TM	P00734	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER1469	Rat TM(Thrombin) ELISA Kit	TM		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1470	Human FSTL1(Follistatin Like Protein 1) ELISA Kit	FSTL1	Q12841	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1471	Rat FSTL1(Follistatin Like Protein 1) ELISA Kit	FSTL1	Q62632	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1472	Human LOX(Lysyl Oxidase) ELISA Kit	LOX	P28300	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM1473	Mouse LOX(Lysyl Oxidase) ELISA Kit	LOX	P28301	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1474	Rat LOX(Lysyl Oxidase) ELISA Kit	LOX	P16636	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EP1475	Porcine PYY(Peptide YY) ELISA Kit	Porcine PYY		Competitive	1000 pg/mL	4.51 pg/mL	15.63-1000 pg/mL
FY-EH1476	Human GH(Growth Hormone) ELISA Kit	GH	P01241	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM1477	Mouse GH(Growth Hormone) ELISA Kit	GH	P06880	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1478	Rat GH(Growth Hormone) ELISA Kit	GH	P01244	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1479	Human COL4(Collagen Type IV) ELISA Kit	COL4		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM1480	Mouse COL4(Collagen Type IV) ELISA Kit	COL4		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER1481	Rat COL4(Collagen Type IV) ELISA Kit	COL4		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH1482	Human FN(Fibronectin) ELISA Kit	FN	P02751	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM1483	Mouse FN(Fibronectin) ELISA Kit	FN	P11276	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER1484	Rat FN(Fibronectin) ELISA Kit	FN	P04937	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH1485	Human IL-9(Interleukin 9) ELISA Kit	IL-9	P15248	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1486	Mouse IL-9(Interleukin 9) ELISA Kit	IL-9	P15247	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1487	Human PC I (Procollagen I) ELISA Kit	PC I		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1488	Human PCIII(Procollagen III) ELISA Kit	PCIII		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1489	Rat PCIII(Procollagen III) ELISA Kit	PCIII		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH1490	Human PIIIINP(N-Terminal Procollagen III Propeptide) ELISA Kit	PIIINP		Sandwich	1500pg/mL	14.06 pg/mL	23.44-1500pg/mL
FY-EM1491	Mouse PIIIINP(N-Terminal Procollagen III Propeptide) ELISA Kit	PIIINP		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1492	Rat PIIIINP(N-Terminal Procollagen III Propeptide) ELISA Kit	PIIINP		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EP1493	Porcine GLP2(Glucagon Like Peptide 2) ELISA Kit	Porcine GLP2		Competitive	10000 pg/mL	47.2 pg/mL	156.25-10000 pg/mL
FY-EH1494	Human TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	TIMP-1	P01033	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1495	Mouse TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	TIMP-1	P12032	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER1496	Rat TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	TIMP-1	P30120	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1497	Human ADCY2(Adenylate Cyclase 2) ELISA Kit	ADCY2	Q08462	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1498	Human AGGF1(Angiogenic Factor with G Patch and FHA Domains 1) ELISA Kit	AGGF1	Q8N302	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1499	Human C4a(Complement Component 4a) ELISA Kit	C4a	P0C0L4	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM1500	Mouse C4a(Complement Component 4a) ELISA Kit	C4a	P01029(678-	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER1501	Rat C4a(Complement Component 4a) ELISA Kit	C4a	P08649	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ESh1502	Sheep GHRL(Ghrelin) ELISA Kit	GHRL		Competitive	10000 pg/mL	48.3 pg/mL	156.25-10000 pg/mL
FY-EH1503	Human C5a(Complement Component 5a) ELISA Kit	C5a	P01031	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM1504	Mouse C5a(Complement Component 5a) ELISA Kit	C5a	P06684(679-	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER1505	Rat C5a(Complement Component 5a) ELISA Kit	C5a		Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH1506	Human IL-5(Interleukin 5) ELISA Kit	IL-5	P05113	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1507	Mouse IL-5(Interleukin 5) ELISA Kit	IL-5	P04401	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1508	Rat IL-5(Interleukin 5) ELISA Kit	IL-5	Q08125	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1509	Human IL-6R(Interleukin 6 Receptor) ELISA Kit	IL-6R	P08887	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1510	Mouse IL-6R(Interleukin 6 Receptor) ELISA Kit	IL-6R	P22272	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1511	Rat IL-6R(Interleukin 6 Receptor) ELISA Kit	IL-6R	P22273	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1512	Human FcγR I (Receptor I for the Fc Region of Immunoglobulin G) ELISA Kit	FcγR I	P12314	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH1513	Human NCAD(Neural Cadherin) ELISA Kit	NCAD	P19022	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1514	Mouse NCAD(Neural Cadherin) ELISA Kit	NCAD		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1515	Rat NCAD(Neural Cadherin) ELISA Kit	NCAD		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EG1516	Goat FSH(Follicle Stimulating Hormone) ELISA Kit	Goat FSH		Competitive	180 ng/mL	0.83 ng/mL	2.82-180 ng/mL
FY-EH1517	Human PDGFRL(Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit	PDGFRL	Q15198	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1518	Mouse PDGFRL(Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit	PDGFRL	Q6PE55	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER1519	Rat PDGFRL(Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit	PDGFRL	Q5RJP7	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EP1520	Porcine NT-ProBNP(N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit	Porcine NT-ProBNP		Competitive	10000 pg/mL	43.8 pg/mL	156.25-10000 pg/mL
FY-EP1521	Porcine GH(Growth Hormone) ELISA Kit	Porcine GH		Competitive	20000 pg/mL	91.4 pg/mL	312.5-20000 pg/mL
FY-ECH1522	Chicken IL17(Interleukin 17) ELISA Kit	Chicken IL17		Competitive	1000 pg/mL	25.4 pg/mL	15.63-1000 pg/mL
FY-EH1523	Human LUM(Lumican) ELISA Kit	LUM	P51884	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1524	Mouse LUM(Lumican) ELISA Kit	LUM	P51885	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1525	Human PRR4(Proline Rich Protein 4, Lacrimal) ELISA Kit	PRR4	Q16378	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1526	Human PSAP(Prosaposin) ELISA Kit	PSAP	P07602	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1527	Human PTGES(Prostaglandin E Synthase, Microsomal) ELISA Kit	PTGES	O14684	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM1528	Mouse PTGES(Prostaglandin E Synthase, Microsomal) ELISA Kit	PTGES	Q9JM51	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1529	Rat PTGES(Prostaglandin E Synthase, Microsomal) ELISA Kit	PTGES	Q9JHF3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1530	Human ISG15(Interferon Stimulated Gene 15) ELISA Kit	ISG15	P05161	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL

FY-EH1531	Human AIMP1(Aminoacyl tRNA Synthetase Complex Interacting Multifunctional Pr AIMP1		Q12904	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH1532	Human 11-DH-TXB2(11-Dehydro-thromboxane B2) ELISA Kit	11-DH-TXB2		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1533	Rat 11-DH-TXB2(11-Dehydrothromboxane B2) ELISA Kit	11-DH-TXB2		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EP1534	Porcine INHA(Inhibin A) ELISA Kit	Porcine INHA		Competitive	1000 pg/mL	5.16 pg/mL	15.63-1000 pg/mL
FY-EH1535	Human 15-LO(Arachidonate 15-Lipoxygenase) ELISA Kit	15-LO	P16050	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM1536	Mouse 15-LO(Arachidonate 15-Lipoxygenase) ELISA Kit	15-LO	P39654	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1537	Rat 15-LO(Arachidonate 15-Lipoxygenase) ELISA Kit	15-LO		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1538	Human BACE1(Beta-Site APP Cleaving Enzyme 1) ELISA Kit	BACE1	P56817	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM1539	Mouse BACE1(Beta-Site APP Cleaving Enzyme 1) ELISA Kit	BACE1	P56818	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1540	Rat BACE1(Beta-Site APP Cleaving Enzyme 1) ELISA Kit	BACE1	P56819	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1541	Human IL-32(Interleukin 32) ELISA Kit	IL-32	P24001	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1542	Human TNFRSF1A(Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit	TNFRSF1A	P19438	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM1543	Mouse TNFRSF1A(Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit	TNFRSF1A	P25118	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER1544	Rat TNFRSF1A(Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit	TNFRSF1A	P22934	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1545	Human IL-15(Interleukin 15) ELISA Kit	IL-15	P40933	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1546	Mouse IL-15(Interleukin 15) ELISA Kit	IL-15	P48346	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1547	Rat IL-15(Interleukin 15) ELISA Kit	IL-15	P97604	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1548	Human 5-LO(Arachidonate 5-Lipoxygenase) ELISA Kit	5-LO	P09917	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1549	Mouse 5-LO(Arachidonate 5-Lipoxygenase) ELISA Kit	5-LO	P48999	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1550	Rat 5-LO(Arachidonate 5-Lipoxygenase) ELISA Kit	5-LO		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1551	Human 5-NT(5-Nucleotidase) ELISA Kit	5-NT	P21589	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1552	Human α 1ACT(Alpha 1 Antichymotrypsin) ELISA Kit	α 1ACT	P01011	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH1553	Human IL-16(Interleukin 16) ELISA Kit	IL-16	Q14005	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1554	Mouse IL-16(Interleukin 16) ELISA Kit	IL-16	O54824	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1555	Rat IL-16(Interleukin 16) ELISA Kit	IL-16	D4A419	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1556	Human AR(Androgen Receptor) ELISA Kit	AR	P10275	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1557	Human ARO(Aromatase) ELISA Kit	ARO	P11511	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1558	Human BLVRB(Biliverdin Reductase B) ELISA Kit	BLVRB	P30043	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1559	Human IL-18(Interleukin 18) ELISA Kit	IL-18	Q14116	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1560	Mouse IL-18(Interleukin 18) ELISA Kit	IL-18	P70380	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1561	Rat IL-18(Interleukin 18) ELISA Kit	IL-18	P97636	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1562	Human IL-19(Interleukin 19) ELISA Kit	IL-19	Q9UHD0	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1563	Mouse IL-19(Interleukin 19) ELISA Kit	IL-19	Q8CJ70	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER1564	Rat IL-19(Interleukin 19) ELISA Kit	IL-19	D3ZF17	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1565	Human ACTG1(Actin Gamma 1) ELISA Kit	ACTG1	P63261	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1566	Human ACTN2(Actinin Alpha 2) ELISA Kit	ACTN2	P35609	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1567	Mouse ACTN2(Actinin Alpha 2) ELISA Kit	ACTN2	Q9JI91	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1568	Rat ACTN2(Actinin Alpha 2) ELISA Kit	ACTN2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1569	Human ACTN3(Actinin Alpha 3) ELISA Kit	ACTN3	Q08043	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1570	Human ADA(Adenosine Deaminase) ELISA Kit	ADA	P00813	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1571	Human ADAM10(A Disintegrin And Metalloprotease 10) ELISA Kit	ADAM10	O14672	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-EH1572	Human ADAM8(A Disintegrin And Metalloprotease 8) ELISA Kit	ADAM8	P78325	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1573	Human ADAM9(A Disintegrin And Metalloprotease 9) ELISA Kit	ADAM9	Q13443	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1574	Human ADAMTS4(ADAM with Thrombospondin Type 1 Motif 4) ELISA Kit	ADAMTS4	O75173	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1575	Human ADCY1(Adenylate Cyclase 1, Brain) ELISA Kit	ADCY1	Q08828	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EP1576	Porcine sIgA(Secretory Immunoglobulin A) ELISA Kit	Porcine sIgA		Competitive	100 ng/mL	0.53 ng/mL	1.57-100 ng/mL
FY-EP1577	Porcine AMH(Anti-Mullerian Hormone) ELISA Kit	Porcine AMH		Competitive	30000 pg/mL	133.8 pg/mL	468.75-30000 pg/mL
FY-EH1578	Human ADH1(Alcohol Dehydrogenase 1) ELISA Kit	ADH1	P07327	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1579	Human ADH(Antidiuretic Hormone) ELISA Kit	ADH	P01185(20~)	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1580	Mouse ADH(Antidiuretic Hormone) ELISA Kit	ADH	P35455	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1581	Rat ADH(Antidiuretic Hormone) ELISA Kit	ADH		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EG1582	Goat ACTH(Adrenocorticotrophic Hormone) ELISA Kit	Goat ACTH		Competitive	1000 pg/mL	4.92 pg/mL	15.63-1000 pg/mL
FY-EH1583	Human ADIPOR1(Adiponectin Receptor 1) ELISA Kit	ADIPOR1	Q96A54	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM1584	Mouse ADIPOR1(Adiponectin Receptor 1) ELISA Kit	ADIPOR1	Q91VH1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EGP1585	Guinea Porcine GT(Gastrin) ELISA Kit	Guinea Porcine GT		Competitive	1000 pg/mL	4.35 pg/mL	15.63-1000 pg/mL
FY-EGP1586	Guinea Porcine CCK8(Cholecystokinin 8, Octapeptide) ELISA Kit	Guinea Porcine CCK8		Competitive	1000 pg/mL	5.28 pg/mL	15.63-1000 pg/mL
FY-EH1587	Human ADIPOR2(Adiponectin Receptor 2) ELISA Kit	ADIPOR2	Q86V24	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EP1588	Porcine I-PTH(IntactParathormone) ELISA Kit	Porcine I-PTH		Competitive	1000 pg/mL	5.33 pg/mL	15.63-1000 pg/mL
FY-EP1589	Porcine GHRL(Ghrelin) ELISA Kit	Porcine GHRL		Competitive	10000 pg/mL	52.9 pg/mL	156.25-10000 pg/mL
FY-EG1590	Goat LH(Luteinizing Hormone) ELISA Kit	Goat LH		Competitive	30000 pg/mL	152.8 pg/mL	468.75-30000 pg/mL
FY-EH1591	Human ADM(Adrenomedullin) ELISA Kit	ADM	P35318	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1592	Mouse ADM(Adrenomedullin) ELISA Kit	ADM	P97297(95~)	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1593	Rat ADM(Adrenomedullin) ELISA Kit	ADM	P43145	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EGP1594	Guinea Porcine IgG2(Immunoglobulin G2) ELISA Kit	Guinea Porcine IgG2		Competitive	25000 ng/mL	124.5 ng/mL	390.63-25000 ng/mL
FY-EG1595	Goat GHRH(Growth Hormone Releasing Hormone) ELISA Kit	Goat GHRH		Competitive	1000 pg/mL	4.49 pg/mL	15.63-1000 pg/mL
FY-EP1596	Porcine AP(Aprotinin) ELISA Kit	Porcine AP		Competitive	100 ng/mL	0.47 ng/mL	1.57-100 ng/mL
FY-EP1597	Porcine ADM(Adrenomedullin) ELISA Kit	Porcine ADM		Competitive	1000 pg/mL	4.85 pg/mL	15.63-1000 pg/mL
FY-ECH1598	Chicken GC(Glucagon) ELISA Kit	Chicken GC		Competitive	1500 pg/mL	7.59 pg/mL	23.44-1500 pg/mL
FY-EP1599	Porcine SST(Somatostatin) ELISA Kit	Porcine SST		Competitive	500 pg/mL	2.74 pg/mL	7.82-500 pg/mL
FY-EH1600	Human ADRA1A(Adrenergic Receptor alpha-1A) ELISA Kit	ADRA1A	P35348	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EP1601	Porcine GHRH(Growth Hormone Releasing Hormone) ELISA Kit	Porcine GHRH		Competitive	1000 pg/mL	5.22 pg/mL	15.63-1000 pg/mL
FY-EP1602	Porcine CCK(Cholecystokinin) ELISA Kit	Porcine CCK		Competitive	1000 pg/mL	4.74 pg/mL	15.63-1000 pg/mL
FY-EP1603	Porcine D2D(D-Dimer) ELISA Kit	Porcine D2D		Competitive	5000 ng/mL	25.4 ng/mL	78.13-5000 ng/mL
FY-EH1604	Human ADRP(Adipose Differentiation Related Protein) ELISA Kit	ADRP	Q99541	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EP1605	Porcine vWF(Von Willebrand Factor) ELISA Kit	Porcine vWF		Competitive	1000 ng/mL	4.83 ng/mL	15.63-1000 ng/mL
FY-ECH1606	Chicken GH(Growth Hormone) ELISA Kit	Chicken GH		Competitive	10000 pg/mL	48.3 pg/mL	156.25-10000 pg/mL
FY-ECH1607	Chicken IgA(Immunoglobulin A) ELISA Kit	Chicken IgA		Competitive	100 ng/mL	0.48 ng/mL	1.57-100 ng/mL
FY-EH1608	Human IL-20(Interleukin 20) ELISA Kit	IL-20	Q9NYY1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM1609	Mouse IL-20(Interleukin 20) ELISA Kit	IL-20	Q9JKV9	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1610	Human DLK-1(Delta Like 1 Homolog) ELISA Kit	DLK-1	P80370	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1611	Human ACE2(Angiotensin I Converting Enzyme 2) ELISA Kit	ACE2	Q9BYF1	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER1612	Rat ACE2(Angiotensin I Converting Enzyme 2) ELISA Kit	ACE2	Q5EGZ1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-ECH1613	Chicken sIgA(Secretory Immunoglobulin A) ELISA Kit	Chicken sIgA		Competitive	50000 pg/mL	248.5 pg/mL	781.25-50000 pg/mL
FY-EH1614	Human BMP-3(Bone Morphogenetic Protein 3) ELISA Kit	BMP-3	P12645	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1615	Rat BMP-3(Bone Morphogenetic Protein 3) ELISA Kit	BMP-3		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1616	Human FLG-Ab(Filaggrin Antibody) ELISA Kit	FLG-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ECH1617	Chicken GHRL(Ghrelin) ELISA Kit	Chicken GHRL		Competitive	10000 pg/mL	46.4 pg/mL	156.25-10000 pg/mL
FY-EP1618	Porcine HB(Hemoglobin) ELISA Kit	Porcine HB		Competitive	300 µg/mL	1.74 µg/mL	4.69-300 µg/mL
FY-EH1619	Human CAPS(Calcyphosine) ELISA Kit	CAPS	Q13938	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1620	Human FABP4(Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit	FABP4	P15090	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-EM1621	Mouse FABP4(Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit	FABP4	P04117	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER1622	Rat FABP4(Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit	FABP4	P70623	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH1623	Human AFP-L3(Alpha-Fetoprotein Lens Culinaris Agglutinin 3) ELISA Kit	AFP-L3		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER1624	Rat AFP-L3(Alpha-Fetoprotein Lens Culinaris Agglutinin 3) ELISA Kit	AFP-L3		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1625	Human FUCA(Alpha-L-Fucosidase, Tissue) ELISA Kit	FUCA	P04066	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1626	Mouse FUCA(Alpha-L-Fucosidase, Tissue) ELISA Kit	FUCA	Q99LJ1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1627	Rat FUCA(Alpha-L-Fucosidase, Tissue) ELISA Kit	FUCA		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1628	Human CCNB1(Cyclin-B1) ELISA Kit	CCNB1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1629	Human AGC(Aggregan) ELISA Kit	AGC	P16112	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EP1630	Porcine IgG(Immunoglobulin G) ELISA Kit	Porcine IgG		Competitive	100 µg/mL	0.57 µg/mL	1.57-100 µg/mL
FY-EH1631	Human AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit	AGER	Q15109	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM1632	Mouse AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit	AGER	Q62151	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1633	Rat AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit	AGER	Q63495	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1634	Human CD3d(Cluster of Differentiation 3d) ELISA Kit	CD3d	P04234	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1635	Human AGR2(Anterior Gradient Protein 2) ELISA Kit	AGR2	O95994	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EGP1636	Guinea Porcine ALB(Albumin) ELISA Kit	Guinea Porcine ALB		Competitive	100 µg/mL	0.43 µg/mL	1.57-100 µg/mL
FY-EH1637	Human AGRN(Agrin) ELISA Kit	AGRN	O00468	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1638	Mouse AGRN(Agrin) ELISA Kit	AGRN	A2ASQ1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1639	Rat AGRN(Agrin) ELISA Kit	AGRN		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ECH1640	Chicken IgG(Immunoglobulin G) ELISA Kit	Chicken IgG		Competitive	100 µg/mL	0.57 µg/mL	1.57-100 µg/mL
FY-EGP1641	Guinea Porcine IgG(Immunoglobulin G) ELISA Kit	Guinea Porcine IgG		Competitive	100 µg/mL	0.48 µg/mL	1.57-100 µg/mL
FY-EH1642	Human AGT(Angiotensinogen) ELISA Kit	AGT	P01019	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM1643	Mouse AGT(Angiotensinogen) ELISA Kit	AGT	P11859	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER1644	Rat AGT(Angiotensinogen) ELISA Kit	AGT	P01015	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-EZel1645	Zebrafish EPI(Epinephrine) ELISA Kit	Zebrafish EPI		Competitive	2000 pg/mL	9.74 pg/mL	31.25-2000 pg/mL
FY-EZel1646	Zebrafish NE(Nor-Epinephrine) ELISA Kit	Zebrafish NE		Competitive	5000 pg/mL	26.1 pg/mL	78.13-5000 pg/mL
FY-EH1647	Human AIF1(Allograft Inflammatory Factor 1) ELISA Kit	AIF1	P55008	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1648	Rat AIF1(Allograft Inflammatory Factor 1) ELISA Kit	AIF1	P55009	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1649	Human ALDM(Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit	ALDM	P05091	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EP1650	Porcine GC(Glucagon) ELISA Kit	Porcine GC		Competitive	2000 pg/mL	11.25 pg/mL	31.25-2000 pg/mL
FY-EH1651	Human Alk-Smase(Alkaline Sphingomyelinase) ELISA Kit	Alk-Smase	Q6UWV6	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EP1652	Porcine ANG1(Angiopoietin1) ELISA Kit	Porcine ANG1		Competitive	200 ng/mL	1.17 ng/mL	3.13-200 ng/mL
FY-EH1653	Human INHB(Inhibin B) ELISA Kit	INHB	P05111(INH	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL

FY-EM1654	Mouse INHB(Inhibin B) ELISA Kit	INHB	Q04997(INH Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1655	Rat INHB(Inhibin B) ELISA Kit	INHB	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EG1656	Goat GHRL(Ghrelin) ELISA Kit	Goat GHRL	Competitive	10 ng/mL	0.071 ng/mL	0.16-10 ng/mL
FY-EG1657	Goat GC(Glucagon) ELISA Kit	Goat GC	Competitive	2000 pg/mL	12.3 pg/mL	31.25-2000 pg/mL
FY-EG1658	Goat Cor(Cortisol) ELISA Kit	Goat Cor	Competitive	200 ng/mL	0.93 ng/mL	3.13-200 ng/mL
FY-EH1659	Human CDK6(Cyclin Dependent Kinase 6) ELISA Kit	CDK6	Q00534 Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EP1660	Porcine FSH(Follicle Stimulating Hormone) ELISA Kit	Porcine FSH	Competitive	400 ng/mL	2.34 ng/mL	6.25-400 ng/mL
FY-EP1661	Porcine LH(Luteinizing Hormone) ELISA Kit	Porcine LH	Competitive	50 ng/mL	0.31 ng/mL	0.79-50 ng/mL
FY-EP1662	Porcine Testo(Testosterone) ELISA Kit	Porcine Testo	Competitive	20 ng/mL	0.11 ng/mL	0.32-20 ng/mL
FY-EP1663	Porcine CRH(Corticotropin Releasing Hormone) ELISA Kit	Porcine CRH	Competitive	10 ng/mL	0.061 ng/mL	0.16-10 ng/mL
FY-EP1664	Porcine ACTH(Adrenocorticotrophic Hormone) ELISA Kit	Porcine ACTH	Competitive	1000 pg/mL	5.81 pg/mL	15.63-1000 pg/mL
FY-EH1665	Human AMBP(Alpha-1-Microglobulin/Bikunin Precursor) ELISA Kit	AMBP	P02760 Sandwich	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-EP1666	Porcine TBARS(Thiobarbituric Acid Reactive Substance) ELISA Kit	Porcine TBARS	Competitive	400 nmol/mL	3.55 nmol/mL	6.25-400 nmol/ml
FY-EZe1667	Zebrafish IgM(Immunoglobulin M) ELISA Kit	Zebrafish IgM	Competitive	50 µg/mL	0.21 µg/mL	0.78-50 µg/mL
FY-ECH1668	Chicken Testo(Testosterone) ELISA Kit	Chicken Testo	Competitive	4000 pg/mL	28.7 pg/mL	62.5-4000 pg/mL
FY-EH1669	Human CDK8(Cyclin Dependent Kinase 8) ELISA Kit	CDK8	P49336 Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EG1670	Goat T3(Triiodothyronine) ELISA Kit	Goat T3	Competitive	10 ng/mL	0.047 ng/mL	0.16-10 ng/mL
FY-EG1671	Goat T4(Thyroxine) ELISA Kit	Goat T4	Competitive	300 ng/mL	1.81 ng/mL	4.69-300 ng/mL
FY-EH1672	Human AMY1(Amylase Alpha 1, Salivary) ELISA Kit	AMY1	P04745 Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER1673	Rat AMY1(Amylase Alpha 1, Salivary) ELISA Kit	AMY1	Q99N59 Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EU1674	KYNA(Kynurenic Acid)ELISA Kit	KYNA	Competitive	500 pmol/mL	3.9 pmol/mL	7.8pmol/ml-500pmol/ml
FY-EU1675	LPS(Lipopolysaccharides) ELISA Kit	LPS	Competitive	100ng/mL	0.68 ng/mL	1.57-100 ng/mL
FY-EH1676	Human IAPP(Islet Amyloid Polypeptide) ELISA Kit	IAPP	P10997 Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM1677	Mouse IAPP(Islet Amyloid Polypeptide) ELISA Kit	IAPP	P12968 Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER1678	Rat IAPP(Islet Amyloid Polypeptide) ELISA Kit	IAPP	P12969 Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1679	Human CDKN3(Cyclin Dependent Kinase Inhibitor 3) ELISA Kit	CDKN3	Q16667 Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1680	Human ANG(Angiogenin) ELISA Kit	ANG	P03950 Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM1681	Mouse ANG(Angiogenin) ELISA Kit	ANG	P21570 Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-ER1682	Rat ANG(Angiogenin) ELISA Kit	ANG	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1683	Human ANG II R1(Angiotensin II Receptor 1) ELISA Kit	ANG II R1	P30556 Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1684	Rat ANG II R1(Angiotensin II Receptor 1) ELISA Kit	ANG II R1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1685	Human ANG II R1-Ab(Angiotensin II Receptor 1 Antibody) ELISA Kit	ANG II R1-Ab	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1686	Human ANG II R2(Angiotensin II Receptor 2) ELISA Kit	ANG II R2	P50052 Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1687	Human ANGPT4(Angiopietin 4) ELISA Kit	ANGPT4	Q9Y264 Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1688	Human ANGPTL1(Angiopietin Like Protein 1) ELISA Kit	ANGPTL1	O95841 Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1689	Human ANGPTL4(Angiopietin Like Protein 4) ELISA Kit	ANGPTL4	Q9BY76 Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM1690	Mouse ANGPTL4(Angiopietin Like Protein 4) ELISA Kit	ANGPTL4	Q9Z1P8 Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1691	Human ANGPTL6(Angiopietin Like Protein 6) ELISA Kit	ANGPTL6	Q8NI99 Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1692	Human ANG-R-Tie1(Angiopietin Receptor Tie1) ELISA Kit	ANG-R-Tie1	P35590 Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1693	Mouse ANG-R-Tie1(Angiopietin Receptor Tie1) ELISA Kit	ANG-R-Tie1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1694	Human ANG-R-Tie2(Angiopietin Receptor Tie2) ELISA Kit	ANG-R-Tie2	Q02763 Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH1695	Human Ank-B(Ankyrin B) ELISA Kit	Ank-B	Q01484	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1696	Human MAG(Myelin Associated Glycoprotein) ELISA Kit	MAG	P20916	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1697	Rat MAG(Myelin Associated Glycoprotein) ELISA Kit	MAG	P07722	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1698	Human CDNF(Cerebral Dopamine Neurotrophic Factor) ELISA Kit	CDNF	Q49AH0	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1699	Human FGFR1(Fibroblast Growth Factor Receptor 1) ELISA Kit	FGFR1	P11362	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1700	Mouse FGFR1(Fibroblast Growth Factor Receptor 1) ELISA Kit	FGFR1	P16092	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1701	Rat FGFR1(Fibroblast Growth Factor Receptor 1) ELISA Kit	FGFR1	Q04589	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1702	Human CNN1(Calponin 1, Basic, Smooth Muscle) ELISA Kit	CNN1	P51911	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1703	Human CNTNAP2(Contactin Associated Protein Like Protein 2) ELISA Kit	CNTNAP2	Q9UHC6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1704	Human ACHAP(Acetylcholinesterase Associated Protein) ELISA Kit	ACHAP	O60888	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1705	Human ACO1(Aconitase 1) ELISA Kit	ACO1	P21399	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1706	Human CPE(Carboxypeptidase E) ELISA Kit	CPE	P16870	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM1707	Mouse CPE(Carboxypeptidase E) ELISA Kit	CPE	Q00493	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1708	Human GM1-Ab(Ganglioside M1 Antibody) ELISA Kit	GM1-Ab		Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH1709	Human CG α (Chorionic Gonadotropin Alpha Polypeptide) ELISA Kit	CG α	P01215	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1710	Human FGFR3(Fibroblast Growth Factor Receptor 3) ELISA Kit	FGFR3	P22607	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1711	Human DPP10(Dipeptidyl Peptidase X) ELISA Kit	DPP10	Q8N608	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1712	Human CPR(Cytochrome P450 Reductase) ELISA Kit	CPR	P16435	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1713	Human CRLF1(Cytokine Receptor Like Factor 1) ELISA Kit	CRLF1	O75462	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1714	Human CSF2R β (Colony Stimulating Factor 2 Receptor Beta) ELISA Kit	CSF2R β	P32927	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1715	Human CA3(Carbonic Anhydrase III, Muscle Specific) ELISA Kit	CA3	P07451	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH1716	Human IR-Ab(Insulin Receptor Antibody) ELISA Kit	IR-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1717	Human CYP1A2(Cytochrome P450, family 1, subfamily A, polypeptide 2) ELISA Kit	CYP1A2	P05177	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1718	Mouse CYP1A2(Cytochrome P450, family 1, subfamily A, polypeptide 2) ELISA Kit	CYP1A2	P00186	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1719	Rat CYP1A2(Cytochrome P450, family 1, subfamily A, polypeptide 2) ELISA Kit	CYP1A2	P04799	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1720	Human Keratin-Ab(Keratin Antibody) ELISA Kit	Keratin-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1721	Human DCLK1(Doublecortin Like Kinase 1) ELISA Kit	DCLK1	O15075	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1722	Mouse DCLK1(Doublecortin Like Kinase 1) ELISA Kit	DCLK1	Q9JLM8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1723	Human DPYD(Dihydropyrimidine Dehydrogenase) ELISA Kit	DPYD	Q12882	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1724	Human DSE(Dermatan Sulfate Epimerase) ELISA Kit	DSE	Q9UL01	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1725	Human EPYC(Epiphycan) ELISA Kit	EPYC	Q99645	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1726	Human FADS2(Fatty Acid Desaturase 2) ELISA Kit	FADS2	O95864	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1727	Human FAIM3(Fas Apoptotic Inhibitory Molecule 3) ELISA Kit	FAIM3	O60667	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1728	Human FASN(Fatty Acid Synthase) ELISA Kit	FASN	P49327	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1729	Mouse FASN(Fatty Acid Synthase) ELISA Kit	FASN	P19096	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1730	Rat FASN(Fatty Acid Synthase) ELISA Kit	FASN	P12785	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1731	Human FBLN4(Fibulin 4) ELISA Kit	FBLN4	O95967	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH1732	Human CYP24A1(Cytochrome P450, family 24, subfamily A, polypeptide 1) ELISA Kit	CYP24A1	Q07973	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1733	Human FCN1(Ficolin 1) ELISA Kit	FCN1	O00602	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1734	Human MAG-Ab(Myelin Associated Glycoprotein Antibody) ELISA Kit	MAG-Ab		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1735	Human Fc γ BP(Fc Fragment of IgG Binding Protein) ELISA Kit	Fc γ BP	Q9Y6R7	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EH1736	Human FcγR3A(Fc Fragment of IgG Low Affinity IIIa Receptor) ELISA Kit	FcγR3A	P08637	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH1737	Human CYP1B1(Cytochrome P450, family 1, subfamily B, polypeptide 1) ELISA Kit	CYP1B1	Q16678	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1738	Rat CYP1B1(Cytochrome P450, family 1, subfamily B, polypeptide 1) ELISA Kit	CYP1B1	Q64678	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH1739	Human FETUA(Fetuin A) ELISA Kit	FETUA	P02765	Sandwich	600ng/mL	5.63 ng/mL	9.38-600ng/mL
FY-EM1740	Mouse FETUA(Fetuin A) ELISA Kit	FETUA	P29699	Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-ER1741	Rat FETUA(Fetuin A) ELISA Kit	FETUA	P24090	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH1742	Human FGγ(Fibrinogen Gamma chain) ELISA Kit	FGγ	P02679	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH1743	Human FLOT2(Flotillin 2) ELISA Kit	FLOT2	Q14254	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH1744	Human GC(Glucagon) ELISA Kit	GC	P01275	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1745	Mouse GC(Glucagon) ELISA Kit	GC	P55095	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1746	Rat GC(Glucagon) ELISA Kit	GC	G3V6P5	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1747	Human GLRX3(Glutaredoxin 3) ELISA Kit	GLRX3	O76003	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1748	Human GIF(Gastric Intrinsic Factor) ELISA Kit	GIF	P27352	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH1749	Human C4BPβ(C4 Binding Protein Beta) ELISA Kit	C4BPβ	P20851	Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-EH1750	Human CA4(Carbonic Anhydrase IV) ELISA Kit	CA4	P22748	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH1751	Human C4BPα(C4 Binding Protein Alpha) ELISA Kit	C4BPα	P04003	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1752	Human GPLD1(Glycosylphosphatidylinositol Specific Phospholipase D1) ELISA Kit	GPLD1	P80108	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH1753	Human PRELP(Proline Arginine Rich End Leucine Rich Repeat Protein) ELISA Kit	PRELP	P51888	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1754	Human PS-Ab(Phosphatidyl Serine Antibody) ELISA Kit	PS-Ab		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1755	Human TARS(Threonyl tRNA Synthetase, cytoplasmic) ELISA Kit	TARS	P26639	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1756	Mouse TARS(Threonyl tRNA Synthetase, cytoplasmic) ELISA Kit	TARS	Q9D0R2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH1757	Human GTF2A1(General Transcription Factor II A, Polypeptide 1) ELISA Kit	GTF2A1	P52655	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1758	Human GUCA2A(Guanylate Cyclase Activator 2A) ELISA Kit	GUCA2A	Q02747	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1759	Human INSL3(Insulin Like Protein 3) ELISA Kit	INSL3	P51460	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1760	Human INSL5(Insulin Like Protein 5) ELISA Kit	INSL5	Q9Y5Q6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1761	Mouse INSL5(Insulin Like Protein 5) ELISA Kit	INSL5	Q9WUG6	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1762	Human GUCA2B(Guanylate Cyclase Activator 2B) ELISA Kit	GUCA2B	Q16661	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1763	Human GZMK(Granzyme K) ELISA Kit	GZMK	P49863	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1764	Rat GZMK(Granzyme K) ELISA Kit	GZMK		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1765	Human ALOXE3(Arachidonate Lipooxygenase 3) ELISA Kit	ALOXE3	Q9BYJ1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1766	Human CYP26A1(Cytochrome P450, family 26, subfamily A, polypeptide 1) ELISA Kit	CYP26A1	O43174	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1767	Human HCCS(Holocytochrome C Synthase) ELISA Kit	HCCS	P53701	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1768	Human CAMK2α(Calcium/Calmodulin-Dependent Protein Kinase II Alpha) ELISA Kit	CAMK2α	Q9UQM7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1769	Human BAG3(Bcl2 Associated Athanogene 3) ELISA Kit	BAG3	O95817	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1770	Human ARSB(Arylsulfatase B) ELISA Kit	ARSB	P15848	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1771	Human ANXA5(Annexin A5) ELISA Kit	ANXA5	P08758	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1772	Mouse ANXA5(Annexin A5) ELISA Kit	ANXA5	P48036	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1773	Rat ANXA5(Annexin A5) ELISA Kit	ANXA5	P14668	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1774	Human HTRA1(HtrA Serine Peptidase 1) ELISA Kit	HTRA1	Q92743	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1775	Human IGF1R(Insulin Like Growth Factor 1 Receptor) ELISA Kit	IGF1R	P08069	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1776	Human HYOU1(Hypoxia Up Regulated 1) ELISA Kit	HYOU1	Q9Y4L1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-ER1777	Rat HYOU1(Hypoxia Up Regulated 1) ELISA Kit	HYOU1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1778	Human ICAM4(Intercellular Adhesion Molecule 4) ELISA Kit	ICAM4	Q14773	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH1779	Human IDE(Insulin Degrading Enzyme) ELISA Kit	IDE	P14735	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER1780	Rat IDE(Insulin Degrading Enzyme) ELISA Kit	IDE	P35559	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1781	Human Surv-Ab(Survivin Antibody) ELISA Kit	Surv-Ab		Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH1782	Human LIPH(Lipase H) ELISA Kit	LIPH	Q8WWY8	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1783	Human LTβR(Lymphotoxin Beta Receptor) ELISA Kit	LTβR	P36941	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1784	Human TPO-Ab(Thyroid Peroxidase Antibody) ELISA Kit	TPO-Ab		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1785	Mouse TPO-Ab(Thyroid Peroxidase Antibody) ELISA Kit	TPO-Ab		Sandwich	100IU/mL	0.94 IU/mL	1.57-100IU/mL
FY-ER1786	Rat TPO-Ab(Thrombopoietin Antibody) ELISA Kit	TPO-Ab		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH1787	Human IGF2BP3(Insulin Like Growth Factor 2 mRNA Binding Protein 3) ELISA Kit	IGF2BP3	O00425	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1788	Human TTN(Titin) ELISA Kit	TTN	Q8WZ42	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1789	Human IGF2R(Insulin Like Growth Factor 2 Receptor) ELISA Kit	IGF2R	P11717	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH1790	Human IGFALS(Insulin Like Growth Factor Binding Protein, Acid Labile Subunit) ELISA Kit	IGFALS	P35858	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1791	Human IGFBP-7(Insulin Like Growth Factor Binding Protein 7) ELISA Kit	IGFBP-7	Q16270	Sandwich	60ng/mL	0.56 ng/mL	0.94-60ng/mL
FY-EH1792	Human ANXA2(Annexin A2) ELISA Kit	ANXA2	P07355	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM1793	Mouse ANXA2(Annexin A2) ELISA Kit	ANXA2	P07356	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1794	Human INSR(Insulin Receptor) ELISA Kit	INSR	P06213	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM1795	Mouse INSR(Insulin Receptor) ELISA Kit	INSR	P15208	Sandwich	60ng/mL	0.56 ng/mL	0.94-60ng/mL
FY-ER1796	Rat INSR(Insulin Receptor) ELISA Kit	INSR	P15127	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1797	Human APLN(Apelin) ELISA Kit	APLN	Q9ULZ1	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM1798	Mouse APLN(Apelin) ELISA Kit	APLN	Q9R0R4	Competitive	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-ER1799	Rat APLN(Apelin) ELISA Kit	APLN	Q9R0R3	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1800	Human ApoA1-Ab(Apolipoprotein A1 Antibody) ELISA Kit	ApoA1-Ab		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH1801	Human LCN1(Lipocalin 1) ELISA Kit	LCN1	P31025	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH1802	Human ApoA4(Apolipoprotein A4) ELISA Kit	ApoA4	P06727	Sandwich	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-ER1803	Rat ApoA4(Apolipoprotein A4) ELISA Kit	ApoA4	P02651	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1804	Human ApoB(Apolipoprotein B) ELISA Kit	ApoB	P04114(Apo	Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-EM1805	Mouse ApoB(Apolipoprotein B) ELISA Kit	ApoB	E9Q414	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1806	Rat ApoB(Apolipoprotein B) ELISA Kit	ApoB	Q7TMA5	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EH1807	Human ApoC1(Apolipoprotein C1) ELISA Kit	ApoC1	P02654	Sandwich	600ng/mL	5.63 ng/mL	9.38-600ng/mL
FY-EM1808	Mouse ApoC1(Apolipoprotein C1) ELISA Kit	ApoC1	P34928	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-ER1809	Rat ApoC1(Apolipoprotein C1) ELISA Kit	ApoC1	P19939	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH1810	Human ApoC2(Apolipoprotein C2) ELISA Kit	ApoC2	P02655	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER1811	Rat ApoC2(Apolipoprotein C2) ELISA Kit	ApoC2	G3V8D4	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH1812	Human ApoC4(Apolipoprotein C4) ELISA Kit	ApoC4	P55056	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM1813	Mouse ApoC4(Apolipoprotein C4) ELISA Kit	ApoC4	Q61268	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1814	Human ApoD(Apolipoprotein D) ELISA Kit	ApoD	P05090	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH1815	Human ApoE(Apolipoprotein E) ELISA Kit	ApoE	P02649	Sandwich	1500ng/mL	14.06 ng/mL	23.44-1500ng/mL
FY-EM1816	Mouse ApoE(Apolipoprotein E) ELISA Kit	ApoE	P08226	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER1817	Rat ApoE(Apolipoprotein E) ELISA Kit	ApoE	P02650	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL

FY-EH1818	Human ApoH(Apolipoprotein H) ELISA Kit	ApoH	P02749	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1819	Mouse ApoH(Apolipoprotein H) ELISA Kit	ApoH	Q01339	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH1820	Human ApoM(Apolipoprotein M) ELISA Kit	ApoM	O95445	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1821	Mouse ApoM(Apolipoprotein M) ELISA Kit	ApoM	Q9Z1R3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1822	Human CA6(Carbonic Anhydrase VI) ELISA Kit	CA6	P23280	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1823	Human AATF(Apoptosis Antagonizing Transcription Factor) ELISA Kit	AATF	Q9NY61	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1824	Mouse AATF(Apoptosis Antagonizing Transcription Factor) ELISA Kit	AATF	Q9JKX4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1825	Human LIPC(Lipase, Hepatic) ELISA Kit	LIPC	P11150	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1826	Mouse LIPC(Lipase, Hepatic) ELISA Kit	LIPC	P27656	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER1827	Rat LIPC(Lipase, Hepatic) ELISA Kit	LIPC		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1828	Human AIF(Apoptosis Inducing Factor) ELISA Kit	AIF	O95831	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1829	Mouse AIF(Apoptosis Inducing Factor) ELISA Kit	AIF	O70200	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH1830	Human LIPF(Lipase, Gastric) ELISA Kit	LIPF	P07098	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1831	Human APAF1(Apoptosis Protease Activating Factor 1) ELISA Kit	APAF1	O14727	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1832	Mouse APAF1(Apoptosis Protease Activating Factor 1) ELISA Kit	APAF1	O88879	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH1833	Human ASK1(Apoptosis Signal Regulating Kinase 1) ELISA Kit	ASK1	Q99683	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1834	Mouse ASK1(Apoptosis Signal Regulating Kinase 1) ELISA Kit	ASK1		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1835	Human AQP-1(Aquaporin 1) ELISA Kit	AQP-1	P29972	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EM1836	Mouse AQP-1(Aquaporin 1) ELISA Kit	AQP-1	Q02013	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1837	Human AQP-2(Aquaporin 2) ELISA Kit	AQP-2	P41181	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1838	Mouse AQP-2(Aquaporin 2) ELISA Kit	AQP-2	P56402	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER1839	Rat AQP-2(Aquaporin 2) ELISA Kit	AQP-2	P34080	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1840	Human AQP-3(Aquaporin 3) ELISA Kit	AQP-3	Q92482	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1841	Mouse AQP-3(Aquaporin 3) ELISA Kit	AQP-3	Q8R2N1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1842	Rat AQP-3(Aquaporin 3) ELISA Kit	AQP-3		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1843	Human AQP-4(Aquaporin 4) ELISA Kit	AQP-4	P55087	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1844	Mouse AQP-4(Aquaporin 4) ELISA Kit	AQP-4	P55088	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1845	Rat AQP-4(Aquaporin 4) ELISA Kit	AQP-4	P47863	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1846	Human AR(Amphiregulin) ELISA Kit	AR	P15514	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM1847	Mouse AR(Amphiregulin) ELISA Kit	AR	P31955	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1848	Rat AR(Amphiregulin) ELISA Kit	AR	P24338	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1849	Human AVPR1A(Arginine Vasopressin Receptor 1A) ELISA Kit	AVPR1A	P37288	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM1850	Mouse AVPR1A(Arginine Vasopressin Receptor 1A) ELISA Kit	AVPR1A	Q62463	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER1851	Rat AVPR1A(Arginine Vasopressin Receptor 1A) ELISA Kit	AVPR1A	P30560	Sandwich	30µg/mL	0.28 µg/mL	0.47-30µg/mL
FY-EH1852	Human KLK6(Kallikrein 6) ELISA Kit	KLK6	Q92876	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1853	Human KLK7(Kallikrein 7) ELISA Kit	KLK7	P49862	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1854	Human ARHGDIα(Rho GDP Dissociation Inhibitor Alpha) ELISA Kit	ARHGDIα	P52565	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1855	Human ARRβ1(Arrestin Beta 1) ELISA Kit	ARRβ1	P49407	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1856	Human ARRβ2(Arrestin Beta 2) ELISA Kit	ARRβ2	P32121	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1857	Human MANF(Mesencephalic Astrocyte Derived Neurotrophic Factor) ELISA Kit	MANF	P55145	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1858	Human MAP2K1(Mitogen Activated Protein Kinase Kinase 1) ELISA Kit	MAP2K1	Q02750	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH1859	Human AhR(Aryl Hydrocarbon Receptor) ELISA Kit	AhR	P35869	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1860	Mouse AhR(Aryl Hydrocarbon Receptor) ELISA Kit	AhR	P30561	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER1861	Rat AhR(Aryl Hydrocarbon Receptor) ELISA Kit	AhR	P41738	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1862	Human ARSA(Arylsulfatase A) ELISA Kit	ARSA	P15289	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1863	Human ASGR2(Asialoglycoprotein Receptor 2) ELISA Kit	ASGR2	P07307	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1864	Human ASGR1(Asialoglycoprotein Receptor 1) ELISA Kit	ASGR1	P07306	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1865	Mouse ASGR1(Asialoglycoprotein Receptor 1) ELISA Kit	ASGR1	P34927	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1866	Rat ASGR1(Asialoglycoprotein Receptor 1) ELISA Kit	ASGR1		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH1867	Human ASM(Acid Sphingomyelinase) ELISA Kit	ASM	P17405	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1868	Rat ASM(Acid Sphingomyelinase) ELISA Kit	ASM		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH1869	Human MAPK14(Mitogen Activated Protein Kinase 14) ELISA Kit	MAPK14	Q16539	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM1870	Mouse MAPK14(Mitogen Activated Protein Kinase 14) ELISA Kit	MAPK14		Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-ER1871	Rat MAPK14(Mitogen Activated Protein Kinase 14) ELISA Kit	MAPK14	P70618	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH1872	Human AST(Aspartate Aminotransferase) ELISA Kit	AST	P17174	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1873	Mouse AST(Aspartate Aminotransferase) ELISA Kit	AST	P05201	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1874	Rat AST(Aspartate Aminotransferase) ELISA Kit	AST	P13221	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1875	Human ASPN(Asporin) ELISA Kit	ASPN	Q9BXN1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1876	Human ATA(Angiotensinase A) ELISA Kit	ATA	Q07075	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1877	Human ATM(Ataxia Telangiectasia Mutated) ELISA Kit	ATM	Q13315	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1878	Mouse ATM(Ataxia Telangiectasia Mutated) ELISA Kit	ATM	Q62388	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER1879	Rat ATM(Ataxia Telangiectasia Mutated) ELISA Kit	ATM		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1880	Human ATXN1(Ataxin 1) ELISA Kit	ATXN1	P54253	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1881	Mouse ATXN1(Ataxin 1) ELISA Kit	ATXN1	P54254	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1882	Rat ATXN1(Ataxin 1) ELISA Kit	ATXN1	Q63540	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1883	Human ATF1(Activating Transcription Factor 1) ELISA Kit	ATF1	P18846	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1884	Human ATF2(Activating Transcription Factor 2) ELISA Kit	ATF2	P15336	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1885	Human ATF3(Activating Transcription Factor 3) ELISA Kit	ATF3	P18847	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1886	Human ATF4(Activating Transcription Factor 4) ELISA Kit	ATF4	P18848	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1887	Human ATF6(Activating Transcription Factor 6) ELISA Kit	ATF6	P18850	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1888	Mouse ATF6(Activating Transcription Factor 6) ELISA Kit	ATF6	F6VAN0	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1889	Rat ATF6(Activating Transcription Factor 6) ELISA Kit	ATF6	G3V909	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH1890	Human ANP(Atrial Natriuretic Peptide) ELISA Kit	ANP	P01160(124-	Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM1891	Mouse ANP(Atrial Natriuretic Peptide) ELISA Kit	ANP	P05125	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1892	Rat ANP(Atrial Natriuretic Peptide) ELISA Kit	ANP		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1893	Human MGL(Monoacylglycerol Lipase) ELISA Kit	MGL	Q99685	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1894	Human AIRE(Autoimmune Regulator) ELISA Kit	AIRE	O43918	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1895	Human AXIN2(Axis Inhibition Protein 2) ELISA Kit	AXIN2	Q9Y2T1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1896	Mouse AXIN2(Axis Inhibition Protein 2) ELISA Kit	AXIN2	O88566	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1897	Rat AXIN2(Axis Inhibition Protein 2) ELISA Kit	AXIN2		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1898	Human AZU1(Azurocidin 1) ELISA Kit	AZU1	P20160	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1899	Human MCRS1(Microspherule Protein 1) ELISA Kit	MCRS1	Q96EZ8	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EH1900	Human BAFFR(B-Cell Activation Factor Receptor) ELISA Kit	BAFFR	Q96RJ3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1901	Mouse BAFFR(B-Cell Activation Factor Receptor) ELISA Kit	BAFFR	Q9D8D0	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER1902	Rat BAFFR(B-Cell Activation Factor Receptor) ELISA Kit	BAFFR		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1903	Human BCL2L1/BCL-X(Bcl-2 Like Protein 1) ELISA Kit	BCL2L1/BCL-X	Q07817	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1904	Human BPI(Bactericidal/Permeability Increasing Protein) ELISA Kit	BPI	P17213	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1905	Human FGF6(Fibroblast Growth Factor 6) ELISA Kit	FGF6	P10767	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EM1906	Mouse FGF6(Fibroblast Growth Factor 6) ELISA Kit	FGF6	P21658	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1907	Rat FGF6(Fibroblast Growth Factor 6) ELISA Kit	FGF6	Q8R5L5	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1908	Human LDHA(Lactate Dehydrogenase A) ELISA Kit	LDHA	P00338	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER1909	Rat LDHA(Lactate Dehydrogenase A) ELISA Kit	LDHA	P04642	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1910	Human Bcl-10(B-cell Lymphoma/Leukemia 10) ELIS Kit ELISA Kit	Bcl-10	O95999	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM1911	Mouse Bcl-10(B-cell Lymphoma/Leukemia 10) ELIS Kit ELISA Kit	Bcl-10	Q9Z0H7	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER1912	Rat Bcl-10(B-cell Lymphoma/Leukemia 10) ELIS Kit ELISA Kit	Bcl-10		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1913	Human BAK1(Bcl2 Antagonist/Killer 1) ELISA Kit	BAK1	Q16611	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM1914	Mouse BAK1(BCL2 Antagonist/Killer 1) ELISA Kit	BAK1	O08734	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1915	Rat BAK1(BCL2 Antagonist/Killer 1) ELISA Kit	BAK1	F7EQF3	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH1916	Human BAD(Bcl2 Associated Death Promoter) ELISA Kit	BAD	Q92934	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EM1917	Mouse BAD(Bcl2 Associated Death Promoter) ELISA Kit	BAD	Q61337	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1918	Rat BAD(Bcl2 Associated Death Promoter) ELISA Kit	BAD		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH1919	Human BAX(Bcl-2 Associated X Protein) ELISA Kit	BAX	Q07812	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM1920	Mouse BAX(Bcl-2 Associated X Protein) ELISA Kit	BAX	Q07813	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1921	Rat BAX(Bcl-2 Associated X Protein) ELISA Kit	BAX		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1922	Human BMF(Bcl-2 Modifying Factor) ELISA Kit	BMF	Q96LC9	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EM1923	Mouse BMF(Bcl-2 Modifying Factor) ELISA Kit	BMF	Q91ZE9	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER1924	Rat BMF(BCL-2 Modifying Factor) ELISA Kit	BMF		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1925	Human BECN1(Beclin 1) ELISA Kit	BECN1	Q14457	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1926	Rat BECN1(Beclin 1) ELISA Kit	BECN1	Q91XJ1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1927	Human SEMA3B(Semaphorin 3B) ELISA Kit	SEMA3B	Q13214	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1928	Human β -EP(Beta-Endorphin) ELISA Kit	β -EP	P01189(237-	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM1929	Mouse β -EP(Beta-Endorphin) ELISA Kit	β -EP	P01193(205-	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1930	Rat β -EP(Beta-Endorphin) ELISA Kit	β -EP		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1931	Human BACE2(Beta-Site APP Cleaving Enzyme 2) ELISA Kit	BACE2	Q9Y5Z0	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1932	Rat BACE2(Beta-Site APP Cleaving Enzyme 2) ELISA Kit	BACE2	Q6IE75	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH1933	Human BID(BH3 Interacting Domain Death Agonist) ELISA Kit	BID	P55957	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1934	Mouse BID(BH3 Interacting Domain Death Agonist) ELISA Kit	BID	P70444	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1935	Rat BID(BH3 Interacting Domain Death Agonist) ELISA Kit	BID	Q9JLT6	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1936	Human MMRN1(Multimerin 1) ELISA Kit	MMRN1	Q13201	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1937	Human Bcl2L11(Bcl-2 Like Protein 11) ELISA Kit	Bcl2L11	O43521	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1938	Human BMP-10(Bone Morphogenetic Protein 10) ELISA Kit	BMP-10	O95393	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1939	Human BMPR2(Bone Morphogenetic Protein Receptor II) ELISA Kit	BMPR2	Q13873	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1940	Mouse BMPR2(Bone Morphogenetic Protein Receptor II) ELISA Kit	BMPR2	O35607	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-ER1941	Rat BMPR2(Bone Morphogenetic Protein Receptor II) ELISA Kit	BMPR2		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH1942	Human MOG(Myelin Oligodendrocyte Glycoprotein) ELISA Kit	MOG	Q16653	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1943	Mouse MOG(Myelin Oligodendrocyte Glycoprotein) ELISA Kit	MOG	Q61885	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1944	Human MT1E(Metallothionein 1E) ELISA Kit	MT1E	P04732	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1945	Human BK(Bradykinin) ELISA Kit	BK	P01042(381-	Competitive	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM1946	Mouse BK(Bradykinin) ELISA Kit	BK	O08677	Competitive	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER1947	Rat BK(Bradykinin) ELISA Kit	BK		Competitive	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH1948	Human MSMβ(Microseminoprotein Beta) ELISA Kit	MSMβ	P08118	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1949	Human BRAK/CXCL14(Breast and Kidney Expressed Chemokine) ELISA Kit	BRAK/CXCL14	Q95715	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1950	Mouse BRAK/CXCL14(Breast and Kidney Expressed Chemokine) ELISA Kit	BRAK/CXCL14	Q9WUQ5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1951	Rat BRAK/CXCL14(Breast and Kidney Expressed Chemokine) ELISA Kit	BRAK/CXCL14		Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EH1952	Human BCAN(Brevican) ELISA Kit	BCAN	Q96GW7	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1953	Human MT1M(Metallothionein 1M) ELISA Kit	MT1M	Q8N339	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH1954	Human MT3(Metallothionein 3) ELISA Kit	MT3	P25713	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1955	Human CNP(C-type Natriuretic Peptide) ELISA Kit	CNP	P23582	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM1956	Mouse CNP(C-type natriuretic peptide) ELISA Kit	CNP	Q61839	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER1957	Rat CNP(C-Type Natriuretic Peptide) ELISA Kit	CNP	P55207	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1958	Human C1q-Ab(Complement 1q Antibody) ELISA Kit	C1q-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1959	Human MUC13(Mucin 13) ELISA Kit	MUC13	Q9H3R2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1960	Human CA724(Tumor Marker CA724) ELISA Kit	CA724		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1961	Mouse CA724(Tumor Marker CA724) ELISA Kit	CA724		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1962	Human CNR-1(Cadherin Related Neuronal Receptor 1) ELISA Kit	CNR-1	P21554	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1963	Human MUC1(Mucin 1) ELISA Kit	MUC1	P15941	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1964	Mouse MUC1(Mucin 1) ELISA Kit	MUC1	Q02496 (21	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1965	Rat MUC1(Mucin 1) ELISA Kit	MUC1	Q35770	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1966	Human CABIN1(Calcineurin Binding Protein 1) ELISA Kit	CABIN1	Q9Y6J0	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1967	Mouse CABIN1(Calcineurin Binding Protein 1) ELISA Kit	CABIN1	B9EKC5	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER1968	Rat CABIN1(Calcineurin Binding Protein 1) ELISA Kit	CABIN1	O88480	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1969	Human CTR(Calcitonin Receptor) ELISA Kit	CTR	P30988	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1970	Mouse CTR(Calcitonin Receptor) ELISA Kit	CTR	Q60755	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1971	Human CACYBP(Calcyclin Binding Protein) ELISA Kit	CACYBP	Q9HB71	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM1972	Mouse CACYBP(Calcyclin Binding Protein) ELISA Kit	CACYBP	Q9CXW3	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1973	Human CALD(Caldesmon) ELISA Kit	CALD	Q05682	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1974	Rat CALD(Caldesmon) ELISA Kit	CALD		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1975	Human MSR1(Macrophage Scavenger Receptor 1) ELISA Kit	MSR1	P21757	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1976	Human CAM(Calmodulin) ELISA Kit	CAM	P0DP23(CA	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EM1977	Mouse CAM(Calmodulin) ELISA Kit	CAM	P0DP26(CA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1978	Rat CAM(Calmodulin) ELISA Kit	CAM		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1979	Human CNX(Calnexin) ELISA Kit	CNX	P27824	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM1980	Mouse CNX(Calnexin) ELISA Kit	CNX	P35564	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1981	Rat CNX(Calnexin) ELISA Kit	CNX	P35565	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH1982	Human CRT(Calreticulin) ELISA Kit	CRT	P27797	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1983	Mouse CRT(Calreticulin) ELISA Kit	CRT	P14211	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1984	Rat CRT(Calreticulin) ELISA Kit	CRT	P18418	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH1985	Human CR(Calretinin) ELISA Kit	CR	P22676	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM1986	Mouse CR(Calretinin) ELISA Kit	CR		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER1987	Rat CR(Calretinin) ELISA Kit	CR	P47728	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH1988	Human CALU(Calumenin) ELISA Kit	CALU	O43852	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1989	Mouse CALU(Calumenin) ELISA Kit	CALU	O35887	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER1990	Rat CALU(Calumenin) ELISA Kit	CALU	O35783	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH1991	Human CAM-Ab(Calmodulin Antibody) ELISA Kit	CAM-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH1992	Human MUC17(Mucin 17) ELISA Kit	MUC17	Q685J3	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH1993	Human MUC2(Mucin 2) ELISA Kit	MUC2	Q02817	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER1994	Rat MUC2(Mucin 2) ELISA Kit	MUC2	Q62635	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1995	Human MUC20(Mucin 20) ELISA Kit	MUC20	Q8N307	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH1996	Human CNR1(Cannabinoid Receptor 1, Brain) ELISA Kit	CNR1	P21554	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM1997	Mouse CNR1(Cannabinoid Receptor 1, Brain) ELISA Kit	CNR1	P47746	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH1998	Human CA1(Carbonic Anhydrase I) ELISA Kit	CA1	P00915	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM1999	Mouse CA1(Carbonic Anhydrase I) ELISA Kit	CA1	P13634	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2000	Rat CA1(Carbonic Anhydrase I) ELISA Kit	CA1		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2001	Human CA2(Carbonic Anhydrase II) ELISA Kit	CA2	P00918	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM2002	Mouse CA2(Carbonic Anhydrase II) ELISA Kit	CA2	P00920	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER2003	Rat CA2(Carbonic Anhydrase II) ELISA Kit	CA2	P27139	Sandwich	180ng/mL	1.69 ng/mL	2.82-180ng/mL
FY-EH2004	Human ANKRD1(Ankyrin Repeat Domain 1) ELISA Kit	ANKRD1	Q15327	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2005	Human cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit	cTnT/TNNT2	P45379	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2006	Mouse cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit	cTnT/TNNT2	P50752	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2007	Rat cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit	cTnT/TNNT2	P50753	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2008	Human IL-7(Interleukin 7) ELISA Kit	IL-7	P13232	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM2009	Mouse IL-7(Interleukin 7) ELISA Kit	IL-7	P10168	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2010	Rat IL-7(Interleukin 7) ELISA Kit	IL-7	P56478	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2011	Human CLS(Cardiolipin Synthase) ELISA Kit	CLS	Q9UJA2	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2012	Mouse CLS(Cardiolipin Synthase) ELISA Kit	CLS	Q80ZM8	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER2013	Rat CLS(Cardiolipin Synthase) ELISA Kit	CLS		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2014	Human CLCF1(Cardiotrophin Like Cytokine Factor 1) ELISA Kit	CLCF1	Q9UBD9	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2015	Mouse CLCF1(Cardiotrophin Like Cytokine Factor 1) ELISA Kit	CLCF1	Q9QZM3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2016	Rat CLCF1(Cardiotrophin Like Cytokine Factor 1) ELISA Kit	CLCF1	A0A8I6GGY	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2017	Human PAR1(Protease Activated Receptor 1) ELISA Kit	PAR1	P25116	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2018	Human COMP(Cartilage Oligomeric Matrix Protein) ELISA Kit	COMP	P49747	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM2019	Mouse COMP(Cartilage Oligomeric Matrix Protein) ELISA Kit	COMP	Q9R0G6	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER2020	Rat COMP(Cartilage Oligomeric Matrix Protein) ELISA Kit	COMP		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2021	Human CASP14(Caspase 14) ELISA Kit	CASP14	P31944	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2022	Mouse CASP14(Caspase 14) ELISA Kit	CASP14	O89094	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-ER2023	Rat CASP14(Caspase 14) ELISA Kit	CASP14		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2024	Human CASP8(Caspase 8) ELISA Kit	CASP8	Q14790	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2025	Mouse CASP8(Caspase 8) ELISA Kit	CASP8	O89110	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2026	Rat CASP8(Caspase 8) ELISA Kit	CASP8	Q9JHX4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2027	Human CASP4(Caspase 4) ELISA Kit	CASP4	P49662	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2028	Mouse CASP4(Caspase 4) ELISA Kit	CASP4	P70343	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER2029	Rat CASP4(Caspase 4) ELISA Kit	CASP4		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2030	Human CASP7(Caspase 7) ELISA Kit	CASP7	P55210	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2031	Mouse CASP7(Caspase 7) ELISA Kit	CASP7	P97864	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2032	Rat CASP7(Caspase 7) ELISA Kit	CASP7		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2033	Human CAD(Caspase Activated Deoxyribonuclease) ELISA Kit	CAD	O76075	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2034	Mouse CAD(Caspase Activated Deoxyribonuclease) ELISA Kit	CAD	O54788	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2035	Human CASP 9(Caspase 9) ELISA Kit	CASP 9	P55211	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2036	Mouse CASP9(Caspase 9) ELISA Kit	CASP9	Q8C3Q9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2037	Rat CASP9(Caspase 9) ELISA Kit	CASP9		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2038	Human CTNN α 1(Catenin, Alpha 1) ELISA Kit	CTNN α 1	P35221	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2039	Mouse CTNN α 1(Catenin, Alpha 1) ELISA Kit	CTNN α 1	P26231	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER2040	Rat CTNN α 1(Catenin, Alpha 1) ELISA Kit	CTNN α 1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2041	Human CTNN β 1(Catenin, Beta 1) ELISA Kit	CTNN β 1	P35222	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2042	Mouse CTNN β 1(Catenin, Beta 1) ELISA Kit	CTNN β 1	Q02248	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2043	Rat CTNN β 1(Catenin, Beta 1) ELISA Kit	CTNN β 1	Q9WU82	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2044	Human CTSD(Cathepsin D) ELISA Kit	CTSD	P07339	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-EM2045	Mouse CTSD(Cathepsin D) ELISA Kit	CTSD	P18242	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2046	Rat CTSD(Cathepsin D) ELISA Kit	CTSD	P24268	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2047	Human CTSK(Cathepsin K) ELISA Kit	CTSK	P43235	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2048	Mouse CTSK(Cathepsin K) ELISA Kit	CTSK	P55097	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2049	Rat CTSK(Cathepsin K) ELISA Kit	CTSK	O35186	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2050	Human CTSL(Cathepsin L) ELISA Kit	CTSL	P07711	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2051	Mouse CTSL(Cathepsin L) ELISA Kit	CTSL	P06797	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2052	Rat CTSL(Cathepsin L) ELISA Kit	CTSL		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2053	Human CDX2(Caudal Type Homeobox Transcription Factor 2) ELISA Kit	CDX2	Q99626	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2054	Mouse CDX2(Caudal Type Homeobox Transcription Factor 2) ELISA Kit	CDX2	P43241	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER2055	Rat CDX2(Caudal Type Homeobox Transcription Factor 2) ELISA Kit	CDX2		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2056	Human C/EBP α (CCAAT/Enhancer Binding Protein Alpha) ELISA Kit	C/EBP α	P49715	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2057	Human C/EBP β (CCAAT/Enhancer Binding Protein Beta) ELISA Kit	C/EBP β	P17676	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2058	Human C/EBP δ (CCAAT/Enhancer Binding Protein, Delta) ELISA Kit	C/EBP δ	P49716	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER2059	Rat C/EBP δ (CCAAT/Enhancer Binding Protein, Delta) ELISA Kit	C/EBP δ		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2060	Human C/EBP ϵ (CCAAT/Enhancer Binding Protein, Epsilon) ELISA Kit	C/EBP ϵ	Q15744	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2061	Human C/EBP γ (CCAAT/Enhancer Binding Protein Gamma) ELISA Kit	C/EBP γ	P53567	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2062	Human CD59(Protectin) ELISA Kit	CD59	P13987	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2063	Human CENPA(Centromere Protein A) ELISA Kit	CENPA	P49450	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EM2064	Mouse CENPA(Centromere Protein A) ELISA Kit	CENPA	O35216	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2065	Rat CENPA(Centromere Protein A) ELISA Kit	CENPA		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2066	Human CA5B(Carbonic Anhydrase V B, Mitochondrial) ELISA Kit	CA5B	Q9Y2D0	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2067	Human CFL1(Cofilin 1, Non-Muscle) ELISA Kit	CFL1	P23528	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2068	Mouse CFL1(Cofilin 1, Non-Muscle) ELISA Kit	CFL1	P18760	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2069	Rat CFL1(Cofilin 1, Non-Muscle) ELISA Kit	CFL1	P45592	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2070	Human FOS(V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog) ELISA Kit	FOS	P01100	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2071	Mouse FOS(V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog) ELISA Kit	FOS	P01101	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2072	Rat FOS(V-Fos FBJ Murine Osteosarcoma Viral Oncogene Homolog) ELISA Kit	FOS	P12841	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2073	Human MUC21(Mucin 21) ELISA Kit	MUC21	Q5SSG8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2074	Human CHEM(Chemerin) ELISA Kit	CHEM	Q99969	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2075	Mouse CHEM(Chemerin) ELISA Kit	CHEM	Q9DD06	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-ER2076	Rat CHEM(Chemerin) ELISA Kit	CHEM		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2077	Human CCL14(Chemokine C-C-Motif Ligand 14) ELISA Kit	CCL14	Q16627	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2078	Human CCL16(Chemokine C-C-Motif Ligand 16) ELISA Kit	CCL16	O15467	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2079	Human CCL3L1(Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit	CCL3L1	P16619	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM2080	Mouse CCL3L1(Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit	CCL3L1		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2081	Human CCR8(Chemokine C-C-Motif Receptor 8) ELISA Kit	CCR8	P51685	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2082	Mouse CCR8(Chemokine C-C-Motif Receptor 8) ELISA Kit	CCR8	P56484	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2083	Rat CCR8(Chemokine C-C-Motif Receptor 8) ELISA Kit	CCR8		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2084	Human TERT(Telomerase Reverse Transcriptase) ELISA Kit	TERT	O14746	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2085	Human CKLFSF4(Chemokine Like Factor Superfamily 4) ELISA Kit	CKLFSF4	Q8IZR5	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2086	Mouse CKLFSF4(Chemokine Like Factor Superfamily 4) ELISA Kit	CKLFSF4	Q8CJ61	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2087	Rat CKLFSF4(Chemokine Like Factor Superfamily 4) ELISA Kit	CKLFSF4		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2088	Human CKLFSF5(Chemokine Like Factor Superfamily 5) ELISA Kit	CKLFSF5	Q96DZ9	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2089	Mouse CKLFSF5(Chemokine Like Factor Superfamily 5) ELISA Kit	CKLFSF5	Q9D6G9	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2090	Rat CKLFSF5(Chemokine Like Factor Superfamily 5) ELISA Kit	CKLFSF5		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2091	Human CKLFSF6(Chemokine Like Factor Superfamily 6) ELISA Kit	CKLFSF6	Q9NX76	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2092	Mouse CKLFSF6(Chemokine Like Factor Superfamily 6) ELISA Kit	CKLFSF6	Q9CZ69	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2093	Rat CKLFSF6(Chemokine Like Factor Superfamily 6) ELISA Kit	CKLFSF6		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2094	Human MUC4(Mucin 4) ELISA Kit	MUC4	Q99102	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2095	Human MUC6(Mucin 6) ELISA Kit	MUC6	Q6W4X9	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2096	Human MYBPC3(Myosin Binding Protein C, Cardiac) ELISA Kit	MYBPC3	Q14896	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2097	Human MYH10(Myosin Heavy Chain 10, Non Muscle) ELISA Kit	MYH10	P35580	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2098	Human MYH9(Myosin Heavy Chain 9, Non Muscle) ELISA Kit	MYH9	P35579	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2099	Human MYL9(Myosin Light Chain 9, Regulatory) ELISA Kit	MYL9	P24844	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2100	Human CCKAR(Cholecystokinin A Receptor) ELISA Kit	CCKAR	P32238	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2101	Mouse CCKAR(Cholecystokinin A Receptor) ELISA Kit	CCKAR	O08786	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2102	Rat CCKAR(Cholecystokinin A Receptor) ELISA Kit	CCKAR	P30551	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2103	Human CETP(Cholesteryl Ester Transfer Protein) ELISA Kit	CETP	P11597	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2104	Mouse CETP(Cholesteryl Ester Transfer Protein) ELISA Kit	CETP		Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL

FY-ER2105	Rat CETP(Cholesteryl Ester Transfer Protein) ELISA Kit	CETP		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2106	Human CENPF(Centromere Protein F) ELISA Kit	CENPF	P49454	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2107	Human CHR2(Cholinergic Receptor, Muscarinic 2) ELISA Kit	CHR2	P08172	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2108	Mouse CHR2(Cholinergic Receptor, Muscarinic 2) ELISA Kit	CHR2	Q9ERZ4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2109	Rat CHR2(Cholinergic Receptor, Muscarinic 2) ELISA Kit	CHR2	P10980	Sandwich	4000pmol/mL	37.5 pmol/mL	62.5-4000pmol/mL
FY-EH2110	Human CA12(Carbonic Anhydrase XII) ELISA Kit	CA12	O43570	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2111	Human CBX3(Chromobox Homolog 3) ELISA Kit	CBX3	Q13185	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2112	Mouse CBX3(Chromobox Homolog 3) ELISA Kit	CBX3	P23198	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2113	Human CHD5(Chromodomain Helicase DNA Binding Protein 5) ELISA Kit	CHD5	Q8TD10	Sandwich	1200pg/mL	11.25 pg/mL	18.75-1200pg/mL
FY-EM2114	Mouse CHD5(Chromodomain Helicase DNA Binding Protein 5) ELISA Kit	CHD5	A2A8L1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2115	Rat CHD5(Chromodomain Helicase DNA Binding Protein 5) ELISA Kit	CHD5		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2116	Human NEF3(Neurofilament 3) ELISA Kit	NEF3	P07197	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM2117	Mouse NEF3(Neurofilament 3) ELISA Kit	NEF3	P08553	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2118	Rat NEF3(Neurofilament 3) ELISA Kit	NEF3		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2119	Human NEFL(Neurofilament, Light Polypeptide) ELISA Kit	NEFL	P07196	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM2120	Mouse NEFL(Neurofilament, Light Polypeptide) ELISA Kit	NEFL	P08551	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2121	Rat NEFL(Neurofilament, Light Polypeptide) ELISA Kit	NEFL	P19527	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2122	Human NOS1/nNOS(Nitric Oxide Synthase 1, Neuronal) ELISA Kit	NOS1/nNOS	P29475	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2123	Mouse NOS1/nNOS(Nitric Oxide Synthase 1, Neuronal) ELISA Kit	NOS1/nNOS	Q9Z0J4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2124	Rat NOS1/nNOS(Nitric Oxide Synthase 1, Neuronal) ELISA Kit	NOS1/nNOS	P29476	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2125	Human C-jun(Proto-oncogene c-Jun) ELISA Kit	C-jun	P05412	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER2126	Rat C-jun(Proto-oncogene c-Jun) ELISA Kit	C-jun		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2127	Human CC16(Clara Cell 16kD protein) ELISA Kit	CC16	P11684	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2128	Human CLDN1(Claudin 1) ELISA Kit	CLDN1	O95832	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2129	Human CD30(Cluster of Differentiation 30) ELISA Kit	CD30	P28908	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2130	Mouse CD30(Cluster of differentiation 30) ELISA Kit	CD30	Q60846	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-ER2131	Rat CD30(Cluster of differentiation 30) ELISA Kit	CD30		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2132	Human CD97(Cluster of Differentiation 97) ELISA Kit	CD97	P48960	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2133	Mouse CD97(Cluster of Differentiation 97) ELISA Kit	CD97	Q9Z0M6	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER2134	Rat CD97(Cluster of Differentiation 97) ELISA Kit	CD97		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2135	Human CD99(Cluster of Differentiation 99) ELISA Kit	CD99	P14209	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2136	Human CD4(Cluster of Differentiation 4) ELISA Kit	CD4	P01730	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2137	Mouse CD4(Cluster Of Differentiation 4) ELISA Kit	CD4	P06332	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER2138	Rat CD4(Cluster of Differentiation 4) ELISA Kit	CD4	P05540	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2139	Human NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit	NOS2/iNOS	P35228	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM2140	Mouse NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit	NOS2/iNOS	P29477	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2141	Rat NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit	NOS2/iNOS	Q06518	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2142	Human CLDN3(Claudin 3) ELISA Kit	CLDN3	O15551	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2143	Human NOS3/eNOS(Nitric Oxide Synthase 3, Endothelial) ELISA Kit	NOS3/eNOS	P29474	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2144	Mouse NOS3/eNOS(Nitric Oxide Synthase 3, Endothelial) ELISA Kit	NOS3/eNOS	P70313	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2145	Rat NOS3/eNOS(Nitric Oxide Synthase 3, Endothelial) ELISA Kit	NOS3/eNOS	Q62600	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-EH2146	Human F7(Coagulation Factor VII) ELISA Kit	F7	P08709	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2147	Mouse F7(Coagulation Factor VII) ELISA Kit	F7	P70375	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2148	Rat F7(Coagulation Factor VII) ELISA Kit	F7	Q8K3U6	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2149	Human F9(Coagulation Factor IX) ELISA Kit	F9	P00740	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EM2150	Mouse F9(Coagulation Factor IX) ELISA Kit	F9		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2151	Rat F9(Coagulation Factor IX) ELISA Kit	F9	P16296	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2152	Human F10(Coagulation Factor X) ELISA Kit	F10	P00742	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EM2153	Mouse F10(coagulation factor X) ELISA Kit	F10	O88947	Sandwich	20µg/mL	0.19 µg/mL	0.32-20µg/mL
FY-ER2154	Rat F10(Coagulation Factor X) ELISA Kit	F10	Q63207	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2155	Human F13(Coagulation Factor X III) ELISA Kit	F13	P00488(F13)	Sandwich	1600pg/mL	15 pg/mL	25-1600pg/mL
FY-ER2156	Rat F13(Coagulation Factor X III) ELISA Kit	F13	O08619	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2157	Human F12(Coagulation Factor XII) ELISA Kit	F12	P00748	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-EM2158	Mouse F12(Coagulation Factor XII) ELISA Kit	F12		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2159	Rat F12(Coagulation Factor XII) ELISA Kit	F12		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2160	Human F5(Coagulation Factor V) ELISA Kit	F5	P12259	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2161	Mouse F5(Coagulation Factor V) ELISA Kit	F5	O88783	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER2162	Rat F5(Coagulation Factor V) ELISA Kit	F5		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2163	Human F11(Coagulation Factor XI) ELISA Kit	F11	P03951	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2164	Rat F11(Coagulation Factor XI) ELISA Kit	F11		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2165	Human F13A1(Coagulation Factor X III A1 Polypeptide) ELISA Kit	F13A1	P00488	Sandwich	1200pg/mL	11.25 pg/mL	18.75-1200pg/mL
FY-EH2166	Human F13B(Coagulation Factor X III B Polypeptide) ELISA Kit	F13B	P05160	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM2167	Mouse F13B(Coagulation Factor X III B Polypeptide) ELISA Kit	F13B	Q07968	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2168	Human CFL2(Cofilin 2, Muscle) ELISA Kit	CFL2	Q9Y281	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM2169	Mouse CFL2(Cofilin 2, Muscle) ELISA Kit	CFL2	P45591	Sandwich	1600pg/mL	15 pg/mL	25-1600pg/mL
FY-ER2170	Rat CFL2(Cofilin 2, Muscle) ELISA Kit	CFL2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2171	Human COL15(Collagen Type X V) ELISA Kit	COL15	P39059(COI	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2172	Human COL6(Collagen Type VI) ELISA Kit	COL6	P12109(COI	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH2173	Human COL1α2(Collagen Type I Alpha 2) ELISA Kit	COL1α2	P08123	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2174	Mouse COL1α2(Collagen Type I Alpha 2) ELISA Kit	COL1α2	Q01149	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2175	Rat COL1α2(Collagen Type I Alpha 2) ELISA Kit	COL1α2	P02466	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2176	Human COL2α1(Collagen Type II Alpha 1) ELISA Kit	COL2α1	P02458	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER2177	Rat COL2α1(Collagen Type II Alpha 1) ELISA Kit	COL2α1		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2178	Human COL3α1(Collagen Type III Alpha 1) ELISA Kit	COL3α1	P02461	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER2179	Rat COL3α1(Collagen Type III Alpha 1) ELISA Kit	COL3α1		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2180	Human COL4α1(Collagen Type IV Alpha 1) ELISA Kit	COL4α1	P02462	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2181	Mouse COL4α1(Collagen Type IV Alpha 1) ELISA Kit	COL4α1	P02463	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2182	Rat COL4α1(Collagen Type IV Alpha 1) ELISA Kit	COL4α1		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2183	Human COL6α1(Collagen Type VI Alpha 1) ELISA Kit	COL6α1	P12109	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2184	Human COL7(Collagen Type VII) ELISA Kit	COL7	Q02388(COI	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2185	Human COL8α1(Collagen Type VIII Alpha 1) ELISA Kit	COL8α1	P27658	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2186	Human COL11α1(Collagen Type XI Alpha 1) ELISA Kit	COL11α1	P12107	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL

FY-EH2187	Human COL17(Collagen Type XVII) ELISA Kit	COL17	Q9UMD9(C1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2188	Human Collagenase I ELISA Kit	Collagenase I	P03956(MM	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2189	Human NPYR1(Neuropeptide Y Receptor Y1) ELISA Kit	NPYR1	P25929	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2190	Human CLL1(Collectin Liver 1) ELISA Kit	CLL1	Q9Y6Z7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2191	Human CCSA-4(Colon Cancer Specific Antigen-4) ELISA Kit	CCSA-4		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2192	Human CLDN9(Claudin 9) ELISA Kit	CLDN9	O95484	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2193	Human NEP(Nepriylisin) ELISA Kit	NEP	P08473	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2194	Human CLDN4(Claudin 4) ELISA Kit	CLDN4	O14493	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2195	Human C7(Complement Component 7) ELISA Kit	C7	P10643	Sandwich	800pg/mL	7.5 pg/mL	12.5-800pg/mL
FY-EH2196	Human C9(Complement Component 9) ELISA Kit	C9	P02748	Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-ER2197	Rat C9(Complement Component 9) ELISA Kit	C9		Sandwich	600ng/mL	5.63 ng/mL	9.38-600ng/mL
FY-EH2198	Human CFP(Complement Factor P) ELISA Kit	CFP	P27918	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2199	Human C3a(Complement Component 3a) ELISA Kit	C3a	P01024	(C3 Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2200	Mouse C3a(Complement Component 3a) ELISA Kit	C3a	P01027	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2201	Rat C3a(Complement Component 3a) ELISA Kit	C3a		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2202	Human C3b(Complement Fragment 3b) ELISA Kit	C3b	P01024	(C3 Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2203	Human C3c(Complement Fragment 3c) ELISA Kit	C3c	P01024	(C3 Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2204	Human NPYR2(Neuropeptide Y Receptor Y2) ELISA Kit	NPYR2	P49146	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2205	Human C5b(Complement Fragment 5b) ELISA Kit	C5b	P01031	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2206	Human NPYR5(Neuropeptide Y Receptor Y5) ELISA Kit	NPYR5	Q15761	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2207	Human CTGF(Connective Tissue Growth Factor) ELISA Kit	CTGF	P29279	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2208	Mouse CTGF(Connective Tissue Growth Factor) ELISA Kit	CTGF	P29268	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2209	Rat CTGF(Connective Tissue Growth Factor) ELISA Kit	CTGF	Q9R1E9	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2210	Human CX26(Connexin 26) ELISA Kit	CX26	P29033	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2211	Mouse CX26(Connexin 26) ELISA Kit	CX26	Q00977	Sandwich	1500pg/mL	14.06 pg/mL	23.44-1500pg/mL
FY-ER2212	Rat CX26(Connexin 26) ELISA Kit	CX26		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2213	Human CHGB(Chromogranin B) ELISA Kit	CHGB	P05060	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2214	Human DEFa1(Defensin Alpha 1, Neutrophil) ELISA Kit	DEFa1	P59665	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2215	Human CTX I (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit	CTX I		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2216	Mouse CTX I (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit	CTX I		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2217	Rat CTX I (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit	CTX I		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2218	Human NTX I (Cross Linked N-telopeptide of Type I Collagen) ELISA Kit	NTX I		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM2219	Mouse NTX I (Cross Linked N-Telopeptide of Type I Collagen) ELISA Kit	NTX I		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2220	Rat NTX I (Cross Linked N-Telopeptide of Type I Collagen) ELISA Kit	NTX I		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2221	Human CTX II (Cross Linked C-telopeptide of Type II Collagen) ELISA Kit	CTX II		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2222	Rat CTX II (Cross Linked C-telopeptide of Type II Collagen) ELISA Kit	CTX II		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2223	Human LDHB(Lactate Dehydrogenase B) ELISA Kit	LDHB	P07195	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER2224	Rat LDHB(Lactate Dehydrogenase B) ELISA Kit	LDHB	P42123	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2225	Human HMGCS1(3-Hydroxy-3-Methylglutaryl-Coenzyme A Synthase 1) ELISA Kit	HMGCS1	Q01581	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2226	Human PIII CP(Procollagen III C-Terminal ProPeptide) ELISA Kit	PIII CP	P02461	Sandwich	800pg/mL	7.5 pg/mL	12.5-800pg/mL
FY-EM2227	Mouse PIII CP(Procollagen III C-Terminal ProPeptide) ELISA Kit	PIII CP		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-EH2228	Human DKK4(Dickkopf Related Protein 4) ELISA Kit	DKK4	Q9UBT3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2229	Human CPP(Copeptin) ELISA Kit	CPP	P01185(126-	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2230	Human CLDN11(Claudin 11) ELISA Kit	CLDN11	O75508	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2231	Human CXCR1(CXC-Chemokine Receptor 1) ELISA Kit	CXCR1	P25024	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2232	Human SCG2(Secretogranin II) ELISA Kit	SCG2	P13521	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH2233	Human MYCBP(C-Myc Binding Protein) ELISA Kit	MYCBP	Q99417	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2234	Human MYC(Myc proto-oncogene protein) ELISA Kit	MYC	P01106	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2235	Mouse MYC(Myc proto-oncogene protein) ELISA Kit	MYC	P01108	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2236	Rat MYC(Myc proto-oncogene protein) ELISA Kit	MYC	P09416	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2237	Human CLEC4E(C-Type Lectin Domain Family 4 Member E) ELISA Kit	CLEC4E	Q9ULY5	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2238	Human LDHC(Lactate Dehydrogenase C) ELISA Kit	LDHC	P07864	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER2239	Rat LDHC(Lactate Dehydrogenase C) ELISA Kit	LDHC	P19629	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2240	Human EphA2(Tyrosine protein kinase receptor A2) ELISA Kit	EphA2	P29317	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2241	Human DKK2(Dickkopf Related Protein 2) ELISA Kit	DKK2	Q9UBU2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2242	Mouse DKK2(Dickkopf Related Protein 2) ELISA Kit	DKK2	Q9QYZ8	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER2243	Rat DKK2(Dickkopf Related Protein 2) ELISA Kit	DKK2		Sandwich	10IU/mL	0.1 IU/mL	0.16-10IU/mL
FY-EH2244	Human D-LDH(D-Lactate Dehydrogenase) ELISA Kit	D-LDH	Q86WU2	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2245	Mouse D-LDH(D-Lactate Dehydrogenase) ELISA Kit	D-LDH		Sandwich	300mIU/mL	2.81 mIU/mL	4.69-300mIU/mL
FY-ER2246	Rat D-LDH(D-Lactate Dehydrogenase) ELISA Kit	D-LDH		Sandwich	10pmol/mL	0.1 pmol/mL	0.16-10pmol/mL
FY-EH2247	Human PA2G4(Proliferation Associated Protein 2G4) ELISA Kit	PA2G4	Q9UQ80	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2248	Human COL1 α 1(Collagen Type I Alpha 1) ELISA Kit	COL1 α 1	P02452	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2249	Rat COL1 α 1(Collagen Type I Alpha 1) ELISA Kit	COL1 α 1	P02454	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2250	Human COL6 α 3(Collagen Type VI Alpha 3) ELISA Kit	COL6 α 3	P12111	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2251	Human COL13(Collagen Type X III) ELISA Kit	COL13	Q5TAT6(CC	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2252	Human COL14(Collagen Type X IV) ELISA Kit	COL14	Q05707(COI	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2253	Human CTX-III(Cross Linked C-telopeptide of Type III Collagen) ELISA Kit	CTX-III		Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2254	Human GlyRS-Ab(Glycyl-tRNA Synthetase Antibody) ELISA Kit	GlyRS-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2255	Human EFNA1(Ephrin A1) ELISA Kit	EFNA1	P20827	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2256	Human GATA3(GATA Binding Protein 3) ELISA Kit	GATA3	P23771	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2257	Mouse GATA3(GATA Binding Protein 3) ELISA Kit	GATA3	P23772	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2258	Rat GATA3(GATA Binding Protein 3) ELISA Kit	GATA3	Q99NH5	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2259	Human GATA4(GATA Binding Protein 4) ELISA Kit	GATA4	P43694	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2260	Mouse GATA4(GATA Binding Protein 4) ELISA Kit	GATA4	Q08369	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2261	Rat GATA4(GATA Binding Protein 4) ELISA Kit	GATA4		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2262	Human GN β 1(G Protein Beta 1) ELISA Kit	GN β 1	P62873	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2263	Human HHIP(Hedgehog Interacting Protein) ELISA Kit	HHIP	Q96QV1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2264	Mouse HHIP(Hedgehog Interacting Protein) ELISA Kit	HHIP	Q7TN16	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2265	Rat HHIP(Hedgehog Interacting Protein) ELISA Kit	HHIP	D3ZGL3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2266	Human H-Cadherin(Cadherin, Heart) ELISA Kit	H-Cadherin	P55290	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2267	Human NRG-4(Neuregulin 4) ELISA Kit	NRG-4	Q8WWG1	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2268	Mouse NRG-4(Neuregulin 4) ELISA Kit	NRG-4	Q9WTX4	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EH2269	Human PTPRK(Protein Tyrosine Phosphatase Receptor Type K) ELISA Kit	PTPRK	Q15262	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2270	Human SELL(L-Selectin) ELISA Kit	SELL	P14151	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2271	Mouse SELL(L-Selectin) ELISA Kit	SELL	P18337	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2272	Rat SELL(L-Selectin) ELISA Kit	SELL	P30836	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2273	Human MXD1(MAX Dimerization Protein 1) ELISA Kit	MXD1	Q05195	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2274	Human HBμ(Hemoglobin, Mu) ELISA Kit	HBμ	Q6B0K9	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2275	Human NAGase(N-Acetyl Beta-D-Glucosaminidase) ELISA Kit	NAGase	O60502	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2276	Mouse NAGase(N-Acetyl Beta-D-Glucosaminidase) ELISA Kit	NAGase	Q9EQQ9	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER2277	Rat NAGase(N-Acetyl Beta-D-Glucosaminidase) ELISA Kit	NAGase		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2278	Human NRG-2(Neuregulin 2) ELISA Kit	NRG-2	O14511	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2279	Rat NRG-2(Neuregulin 2) ELISA Kit	NRG-2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2280	Human NRG-3(Neuregulin 3) ELISA Kit	NRG-3	P56975	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2281	Human N-MID-OT(N-MID Osteocalcin) ELISA Kit	N-MID-OT	P02818	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2282	Human PTPRN(Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit	PTPRN	Q16849	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2283	Mouse PTPRN(Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit	PTPRN	Q60673	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER2284	Rat PTPRN(Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit	PTPRN		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2285	Human GAS(N-Acetylgalactosamine 6-Sulfatase) ELISA Kit	GAS	P34059	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2286	Mouse GAS(N-Acetylgalactosamine 6-Sulfatase) ELISA Kit	GAS	Q571E4	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER2287	Rat GAS(N-Acetylgalactosamine 6-Sulfatase) ELISA Kit	GAS		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2288	Human PTPRO(Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit	PTPRO	Q92729	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2289	Mouse PTPRO(Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit	PTPRO	E9Q612	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER2290	Rat PTPRO(Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit	PTPRO		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2291	Human P27(P27 Protein) ELISA Kit	P27	P46527	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2292	Human TP53(Tumor Protein p53) ELISA Kit	TP53	P04637	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM2293	Mouse TP53(Tumor Protein p53) ELISA Kit	TP53	P02340	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER2294	Rat TP53(Tumor Protein p53) ELISA Kit	TP53	P10361	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2295	Human PDAP1(Pdgfa Associated Protein 1) ELISA Kit	PDAP1	Q13442	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2296	Mouse PDAP1(PDGFA Associated Protein 1) ELISA Kit	PDAP1	Q3UHX2	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2297	Rat PDAP1(Pdgfa Associated Protein 1) ELISA Kit	PDAP1		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2298	Human PDCN(Podocin) ELISA Kit	PDCN	Q9NP85	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2299	Human Preptin ELISA Kit	Preptin	P01344(93~1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2300	Human P-Cadherin(Cadherin, Placental) ELISA Kit	P-Cadherin	P22223	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2301	Mouse P-Cadherin(Cadherin, Placental) ELISA Kit	P-Cadherin	P10287	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH2302	Human P-GP(Permeability Glycoprotein) ELISA Kit	P-GP	P08183	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2303	Human RAPGEF1(Rap Guanine Nucleotide Exchange Factor 1) ELISA Kit	RAPGEF1	Q13905	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2304	Human Rac-GAP1(Rac-GTPase Activating Protein 1) ELISA Kit	Rac-GAP1	Q9H0H5	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2305	Human RHOB(Ras Homolog Gene Family, Memmber B) ELISA Kit	RHOB	P62745	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2306	Human RCOR3(Rest Corepressor 3) ELISA Kit	RCOR3	Q9P2K3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2307	Human RHOA(Ras Homolog Gene Family, Memmber A) ELISA Kit	RHOA	P61586	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM2308	Mouse RHOA(Ras Homolog Gene Family, Memmber A) ELISA Kit	RHOA	Q9QU10	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2309	Rat RHOA(Ras Homolog Gene Family, Member A) ELISA Kit	RHOA		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH2310	Human Rock-2(Rho Associated Coiled Coil Containing Protein Kinase 2) ELISA Kit	Rock-2	O75116	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2311	Human RND1(Rho Family GTPase 1) ELISA Kit	RND1	Q92730	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2312	Human RUNX2(Runt-Related Transcription Factor 2) ELISA Kit	RUNX2	Q13950	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2313	Mouse RUNX2(Runt Related Transcription Factor 2) ELISA Kit	RUNX2	Q08775	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2314	Rat RUNX2(Runt Related Transcription Factor 2) ELISA Kit	RUNX2	Q9Z2J9	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2315	Human COL9α3(Collagen Type IX Alpha 3) ELISA Kit	COL9α3	Q14050	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2316	Human Salusin α(Salusin Alpha) ELISA Kit	Salusin α	Q8N2E6(21<	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2317	Human SHC/SLP-76(SH2-Containing Protein) ELISA Kit	SHC/SLP-76	Q13094	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2318	Human Slit2(Slit Homolog 2) ELISA Kit	Slit2	O94813	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2319	Mouse Slit2(Slit Homolog 2) ELISA Kit	Slit2	Q9R1B9	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2320	Rat SLIT2(Slit Homolog 2) ELISA Kit	SLIT2	Q9WVC1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2321	Human Smad1/MADH1(Mothers Against Decapentaplegic Homolog 1) ELISA Kit	Smad1/MADH1	Q15797	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2322	Mouse Smad1/MADH1(Mothers Against Decapentaplegic Homolog 1) ELISA Kit	Smad1/MADH1	P70340	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2323	Rat Smad1/MADH1(Mothers Against Decapentaplegic Homolog 1) ELISA Kit	Smad1/MADH1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2324	Human Smad4/MADH4(Mothers Against Decapentaplegic Homolog 4) ELISA Kit	Smad4/MADH4	Q13485	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2325	Mouse Smad4/MADH4(Mothers Against Decapentaplegic Homolog 4) ELISA Kit	Smad4/MADH4	P97471	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2326	Rat Smad4/MADH4(Mothers Against Decapentaplegic Homolog 4) ELISA Kit	Smad4/MADH4	O70437	Sandwich	2000pg/ml	18.75 pg/ml	31.25-2000pg/ml
FY-EH2327	Human Smad7/MADH7(Mothers Against Decapentaplegic Homolog 7) ELISA Kit	Smad7/MADH7	O15105	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2328	Mouse Smad7/MADH7(Mothers Against Decapentaplegic Homolog 7) ELISA Kit	Smad7/MADH7	O35253	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2329	Rat Smad7/MADH7(Mothers Against Decapentaplegic Homolog 7) ELISA Kit	Smad7/MADH7	O88406	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2330	Human Sp100(Sp100 Nuclear Antigen) ELISA Kit	Sp100	P23497	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2331	Human PIAS1(Protein Inhibitor of Activated STAT 1) ELISA Kit	PIAS1	O75925	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2332	Mouse PIAS1(Protein Inhibitor Of Activated STAT 1) ELISA Kit	PIAS1	O88907	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2333	Rat PIAS1(Protein Inhibitor Of Activated STAT 1) ELISA Kit	PIAS1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2334	Human PIAS3(Protein Inhibitor of Activated STAT3) ELISA Kit	PIAS3	Q9Y6X2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2335	Human PCDHβ16(Protocadherin Beta 16) ELISA Kit	PCDHβ16	Q9NRJ7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2336	Human COL5α2(Collagen Type V Alpha 2) ELISA Kit	COL5α2	P05997	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH2337	Human SAG(S Antigen) ELISA Kit	SAG	P10523	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2338	Human THP(Tamm–Horsfall Glycoprotein) ELISA Kit	THP	P07911	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2339	Mouse THP(Tamm–Horsfall Glycoprotein) ELISA Kit	THP	Q91X17	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER2340	Rat THP(Tamm–Horsfall Glycoprotein) ELISA Kit	THP	P27590	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2341	Human TBP(TATA Binding Protein) ELISA Kit	TBP	P20226	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2342	Mouse TBP(TATA Binding Protein) ELISA Kit	TBP	P29037	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2343	Rat TBP(Tata Binding Protein) ELISA Kit	TBP		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2344	Human MAPτ(Microtubule Associated Protein Tau/Tau Protein) ELISA Kit	MAPτ	P10636	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM2345	Mouse MAPτ(Microtubule Associated Protein Tau/Tau Protein) ELISA Kit	MAPτ	P10637	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2346	Rat MAPτ(Microtubule Associated Protein Tau/Tau Protein) ELISA Kit	MAPτ	P19332	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2347	Human TIEG1(TGF Beta Inducible Early Response Gene 1) ELISA Kit	TIEG1	Q13118	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2348	Mouse TIEG1(TGF Beta Inducible Early Response Gene 1) ELISA Kit	TIEG1	O89091	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER2349	Rat TIEG1(TGF Beta Inducible Early Response Gene 1) ELISA Kit	TIEG1		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2350	Human TLR9(Toll-Like Receptor 9) ELISA Kit	TLR9	Q9NR96	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EM2351	Mouse TLR9(Toll-like Receptor 9) ELISA Kit	TLR9	Q9EQU3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2352	Rat TLR9(Toll-Like Receptor 9) ELISA Kit	TLR9		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2353	Human PAR3(Protease Activated Receptor 3) ELISA Kit	PAR3	O00254	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2354	Human β -CTx(Beta Crosslaps) ELISA Kit	β -CTx		Competitive	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EM2355	Mouse β -CTx(Beta Crosslaps) ELISA Kit	β -CTx		Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER2356	Rat β -CTx(Beta Crosslaps) ELISA Kit	β -CTx		Competitive	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH2357	Human XCR1(XC-Chemokine Receptor 1) ELISA Kit	XCR1	P46094	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2358	Human PCDH15(Protocadherin 15) ELISA Kit	PCDH15	Q96QU1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2359	Mouse PCDH15(Protocadherin 15) ELISA Kit	PCDH15	Q99PJ1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2360	Human α 1BG(Alpha-1-B-Glycoprotein) ELISA Kit	α 1BG	P04217	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2361	Human α 1-AT(Alpha 1-Antitrypsin) ELISA Kit	α 1-AT	P01009	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM2362	Mouse α 1-AT(Alpha 1-Antitrypsin) ELISA Kit	α 1-AT		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER2363	Rat α 1-AT(Alpha 1-Antitrypsin) ELISA Kit	α 1-AT	P17475	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH2364	Human PLC γ 2(Phospholipase C Gamma 2, Phosphatidylinositol Specific) ELISA Kit	PLC γ 2	P16885	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2365	Human α 2-AP(α 2-Antiplasmin) ELISA Kit	α 2-AP	P08697	Sandwich	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-EM2366	Mouse α 2-AP(α 2-Antiplasmin) ELISA Kit	α 2-AP		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER2367	Rat α 2-AP(α 2-Antiplasmin) ELISA Kit	α 2-AP		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2368	Human IDU α (Alpha-L-Iduronidase) ELISA Kit	IDU α	P35475	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2369	Human NAGLU(N-Acetyl Alpha-D-Glucosaminidase) ELISA Kit	NAGLU	P54802	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2370	Human GL α (Galactosidase Alpha) ELISA Kit	GL α	P06280	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2371	Human α Actin(Alpha-Actin) ELISA Kit	α Actin		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2372	Human α -La(Alpha-Lactalbumin) ELISA Kit	α -La	P00709	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2373	Human SNC α (Synuclein Alpha) ELISA Kit	SNC α	P37840	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM2374	Mouse SNC α (Synuclein Alpha) ELISA Kit	SNC α	O55042	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2375	Rat SNC α (Synuclein Alpha) ELISA Kit	SNC α	P37377	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2376	Human PDGFR β (Platelet Derived Growth Factor Receptor Beta) ELISA Kit	PDGFR β	P09619	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2377	Human P II CP(Procollagen II C-Terminal ProPeptide) ELISA Kit	P II CP		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2378	Rat P II CP(Procollagen II C-Terminal ProPeptide) ELISA Kit	P II CP		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2379	Human PCDH γ A2(Protocadherin Gamma A2) ELISA Kit	PCDH γ A2	Q9Y5H1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2380	Human MAN1A1(Mannosidase Alpha Class 1A Member 1) ELISA Kit	MAN1A1	P33908	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2381	Human β -LPH(β -Lipotropic Hormone) ELISA Kit	β -LPH		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2382	Human DEFB1(Defensin Beta 1) ELISA Kit	DEFB1	P60022	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2383	Human DEFB2/DEFB2(Defensin Beta 2) ELISA Kit	DEFB2/DEFB2	O15263	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2384	Human β Manase(β Mannosidase) ELISA Kit	β Manase	O00462	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2385	Mouse β Manase(β mannosidase) ELISA Kit	β Manase	Q8K214	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH2386	Human CYSLTR1(Cysteinyl Leukotriene Receptor 1) ELISA Kit	CYSLTR1	Q9Y271	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2387	Human IRAK3(Interleukin 1 Receptor Associated Kinase 3) ELISA Kit	IRAK3	Q9Y616	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2388	Human GUS β (β -glucuronidase) ELISA Kit	GUS β	P08236	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2389	Mouse GUS β (β -glucuronidase) ELISA Kit	GUS β		Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER2390	Rat GUS β (β -glucuronidase) ELISA Kit	GUS β		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2391	Human SYCN(Syncollin) ELISA Kit	SYCN	Q0VAF6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH2392	Human TUBβ(β-Tubulin) ELISA Kit	TUBβ	P07437	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2393	Human γ-ECS(γ-Glutamyl Systeine Synthetase) ELISA Kit	γ-ECS	P48506	Sandwich	15ng/mL	0.14 ng/mL	0.24-15ng/mL
FY-EH2394	Human PLA2G10(Phospholipase A2, Group X) ELISA Kit	PLA2G10	O15496	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH2395	Human γGT1(Gamma Glutamyltransferase 1) ELISA Kit	γGT1	P19440	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER2396	Rat γGT1(Gamma Glutamyltransferase 1) ELISA Kit	γGT1		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2397	Human POMC(Pro-Opiomelanocortin) ELISA Kit	POMC	P01189	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2398	Mouse POMC(Pro-Opiomelanocortin) ELISA Kit	POMC	P01193	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-EH2399	Human OXA(Orexin A) ELISA Kit	OXA	O43612(34~	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2400	Mouse OXA(Orexin A) ELISA Kit	OXA	O55241	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER2401	Rat OXA(Orexin A) ELISA Kit	OXA		Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2402	Human PIK3Cβ(Phosphoinositide-3-Kinase Catalytic Beta Polypeptide) ELISA Kit	PIK3Cβ	P42338	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2403	Human OIT3(Oncoprotein Induced Transcript 3) ELISA Kit	OIT3	Q8WWZ8	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2404	Mouse OIT3(Oncoprotein Induced Transcript 3) ELISA Kit	OIT3	Q8R4V5	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2405	Rat OIT3(Oncoprotein Induced Transcript 3) ELISA Kit	OIT3		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2406	Human PIK3Cσ(Phosphoinositide-3-Kinase Catalytic Delta Polypeptide) ELISA Kit	PIK3Cσ	O00329	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2407	Human OCT2(Octamer-binding protein 2) ELISA Kit	OCT2	P09086	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM2408	Mouse OCT2(Octamer-binding protein 2) ELISA Kit	OCT2	Q00196	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2409	Human OCT4(Octamer-binding protein 4) ELISA Kit	OCT4	Q01860	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2410	Rat OCT4(Octamer-binding protein 4) ELISA Kit	OCT4		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2411	Human PLA2G12B(Phospholipase A2, Group XII B) ELISA Kit	PLA2G12B	Q9BX93	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2412	Human IL-2Rα(Interleukin 2 Receptor Alpha) ELISA Kit	IL-2Rα	P01589	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EH2413	Human IL-2Rβ(Interleukin 2 Receptor Beta) ELISA Kit	IL-2Rβ	P14784	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER2414	Rat IL-2Rβ(Interleukin 2 Receptor Beta) ELISA Kit	IL-2Rβ		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2415	Human NFIL3(Nuclear Factor, Interleukin 3 Regulated) ELISA Kit	NFIL3	Q16649	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2416	Mouse NFIL3(Nuclear Factor, Interleukin 3 Regulated) ELISA Kit	NFIL3	O08750	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2417	Rat NFIL3(Nuclear Factor, Interleukin 3 Regulated) ELISA Kit	NFIL3		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2418	Human LILRB4(Leukocyte Immunoglobulin Like Receptor Subfamily B Member 4) ELISA Kit	LILRB4	Q8NHJ6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2419	Human LIFR(Leukemia Inhibitory Factor Receptor) ELISA Kit	LIFR	P42702	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2420	Mouse LIFR(Leukemia Inhibitory Factor Receptor) ELISA Kit	LIFR	P42703	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2421	Rat LIFR(Leukemia Inhibitory Factor Receptor) ELISA Kit	LIFR	O70535	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2422	Human GARS(Glycyl tRNA Synthetase) ELISA Kit	GARS	P41250	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2423	Human EPHA1(Ephrin Type A Receptor 1) ELISA Kit	EPHA1	P21709	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2424	Human SPOP(Speckle Type POZ Protein) ELISA Kit	SPOP	O43791	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2425	Mouse SPOP(Speckle Type POZ Protein) ELISA Kit	SPOP	Q6ZWS8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2426	Human GP4/CD36(Platelet Membrane Glycoprotein IV) ELISA Kit	GP4/CD36	P16671	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2427	Mouse GP4/CD36(Platelet Membrane Glycoprotein IV) ELISA Kit	GP4/CD36	Q08857	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2428	Rat GP4/CD36(Platelet Membrane Glycoprotein IV) ELISA Kit	GP4/CD36		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2429	Human PLA2G4D(Phospholipase A2, Group IV D) ELISA Kit	PLA2G4D	Q86XP0	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2430	Human CYSLTR2(Cysteiny Leukotriene Receptor 2) ELISA Kit	CYSLTR2	Q9NS75	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2431	Mouse CYSLTR2(Cysteiny Leukotriene Receptor 2) ELISA Kit	CYSLTR2	Q920A1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2432	Rat CYSLTR2(Cysteiny Leukotriene Receptor 2) ELISA Kit	CYSLTR2	Q924T9	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EH2433	Human GAL1(Galectin 1) ELISA Kit	GAL1	P09382	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2434	Mouse GAL1(Galectin 1) ELISA Kit	GAL1	P16045	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER2435	Rat GAL1(Galectin 1) ELISA Kit	GAL1	P11762	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2436	Human GAL12(Galectin 12) ELISA Kit	GAL12	Q96DT0	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2437	Mouse GAL12(Galectin 12) ELISA Kit	GAL12	Q91VD1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2438	Human GAL14(Galectin 14) ELISA Kit	GAL14	Q8TCE9	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2439	Human GAL2(Galectin 2) ELISA Kit	GAL2	P05162	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2440	Mouse GAL2(Galectin 2) ELISA Kit	GAL2	Q9CQW5	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2441	Rat GAL2(Galectin 2) ELISA Kit	GAL2		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2442	Human EPHA3(Ephrin Type A Receptor 3) ELISA Kit	EPHA3	P29320	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2443	Human GAL4(Galectin 4) ELISA Kit	GAL4	P56470	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2444	Mouse GAL4(Galectin 4) ELISA Kit	GAL4	Q8K419	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2445	Rat GAL4(Galectin 4) ELISA Kit	GAL4	P38552	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2446	Human GAL8(Galectin 8) ELISA Kit	GAL8	O00214	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2447	Mouse GAL8(Galectin 8) ELISA Kit	GAL8	Q9JL15	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2448	Rat GAL8(Galectin 8) ELISA Kit	GAL8		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2449	Human GAL7(Galectin 7) ELISA Kit	GAL7	P47929	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EM2450	Mouse GAL7(Galectin 7) ELISA Kit	GAL7	O54974	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2451	Rat GAL7(Galectin 7) ELISA Kit	GAL7		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2452	Human GAL9(Galectin 9) ELISA Kit	GAL9	O00182	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM2453	Mouse GAL9(Galectin 9) ELISA Kit	GAL9	O08573	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER2454	Rat GAL9(Galectin 9) ELISA Kit	GAL9	P97840	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH2455	Human UBC(Urinary Bladder Cancer Antigen) ELISA Kit	UBC	P0CG48	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2456	Human PLA2R1(Phospholipase A2 Receptor 1) ELISA Kit	PLA2R1	Q13018	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2457	Human PLB(Phospholipase B) ELISA Kit	PLB	Q6P1J6	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2458	Human NP(Nandrolone Phenylpropionate) ELISA Kit	NP		Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2459	Human EFNA3(Ephrin A3) ELISA Kit	EFNA3	P52797	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2460	Human PLD2(Phospholipase D2) ELISA Kit	PLD2	O14939	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2461	Human PDXK(Pyridoxal Kinase) ELISA Kit	PDXK	O00764	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2462	Human OCLN(Occludin) ELISA Kit	OCLN	Q16625	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2463	Rat OCLN(Occludin) ELISA Kit	OCLN	Q6P6T5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2464	Human SEMA7A(Semaphorin 7A) ELISA Kit	SEMA7A	O75326	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2465	Mouse SEMA7A(Semaphorin 7A) ELISA Kit	SEMA7A	Q9QUR8	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2466	Rat SEMA7A(Semaphorin 7A) ELISA Kit	SEMA7A		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2467	Human PLVAP(Plasmalemma Vesicle Associated Protein) ELISA Kit	PLVAP	Q9BX97	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2468	Human HAGH(Hydroxyacylglyuthathione Hydrolase) ELISA Kit	HAGH	Q16775	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2469	Human SPJ(Surfactant Associated Protein J) ELISA Kit	SPJ	P56856(CLE	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2470	Human GPRC5A(G Protein Coupled Receptor, Family C, Group 5, Member A) ELIS.	GPRC5A	Q8NFFJ5	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2471	Human FABP5(Fatty Acid Binding Protein 5, Epidermal) ELISA Kit	FABP5	Q01469	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2472	Mouse FABP5(Fatty Acid Binding Protein 5, Epidermal) ELISA Kit	FABP5	Q05816	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2473	Human PDK4(Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit	PDK4	Q16654	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL

FY-EM2474	Mouse PDK4(Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit	PDK4	O70571	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2475	Human PDP(Pyruvate Dehydrogenase Phosphatase) ELISA Kit	PDP	Q9P0J1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2476	Mouse PDP(Pyruvate Dehydrogenase Phosphatase) ELISA Kit	PDP	Q3UV70	Sandwich	200mIU/mL	1.88 mIU/mL	3.13-200mIU/mL
FY-ER2477	Rat PDP(Pyruvate Dehydrogenase Phosphatase) ELISA Kit	PDP		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2478	Human PPOX(Protoporphyrinogen Oxidase) ELISA Kit	PPOX	P50336	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2479	Human PRB1(Basic Salivary Proline Rich Protein 1) ELISA Kit	PRB1	P04280	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2480	Human FOXC1(Forkhead Box Protein C1) ELISA Kit	FOXC1	Q12948	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2481	Human FOXO1(Forkhead Box Protein O1) ELISA Kit	FOXO1	Q12778	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2482	Mouse FOXO1(Forkhead Box Protein O1) ELISA Kit	FOXO1	Q9R1E0	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2483	Human FOXO3(Forkhead Box Protein O3) ELISA Kit	FOXO3	O43524	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2484	Human FOXO4(Forkhead Box Protein O4) ELISA Kit	FOXO4	P98177	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2485	Human FOXP3(Forkhead Box Protein P3) ELISA Kit	FOXP3	Q9BZS1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2486	Mouse FOXP3(Forkhead Box Protein P3) ELISA Kit	FOXP3	Q99JB6	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2487	Rat FOXP3(Forkhead Box Protein P3) ELISA Kit	FOXP3	D3ZKI1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2488	Human LATS(Long Acting Thyroid Stimulator) ELISA Kit	LATS		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2489	Human EFNA4(Ephrin A4) ELISA Kit	EFNA4	P52798	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2490	Human TFF3(Trefoil Factor 3, Intestinal) ELISA Kit	TFF3	Q07654	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM2491	Mouse TFF3(Trefoil Factor 3, Intestinal) ELISA Kit	TFF3	Q62395	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER2492	Rat TFF3(Trefoil Factor 3, Intestinal) ELISA Kit	TFF3	Q03191	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH2493	Human TAGLN(Transgelin) ELISA Kit	TAGLN	Q01995	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH2494	Human TEP1(Telomerase Associated Protein 1) ELISA Kit	TEP1	Q99973	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2495	Human PRB2(Basic Salivary Proline Rich Protein 2) ELISA Kit	PRB2	P02812	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2496	Human SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit	SOD1	P00441	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2497	Mouse SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit	SOD1	P08228	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2498	Rat SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit	SOD1	P07632	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2499	Human OGP(Osteogenic Growth Peptide) ELISA Kit	OGP		Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2500	Human FGF23(Fibroblast Growth Factor 23) ELISA Kit	FGF23	Q9GZV9	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM2501	Mouse FGF23(Fibroblast Growth Factor 23) ELISA Kit	FGF23	Q9EPC2	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2502	Rat FGF23(Fibroblast Growth Factor 23) ELISA Kit	FGF23	Q8VI82	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2503	Human FGF9(Fibroblast Growth Factor 9) ELISA Kit	FGF9	P31371	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2504	Mouse FGF9(Fibroblast Growth Factor 9) ELISA Kit	FGF9	P54130	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER2505	Rat FGF9(Fibroblast Growth Factor 9) ELISA Kit	FGF9	P36364	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2506	Human FGFBP1(Fibroblast Growth Factor Binding Protein 1) ELISA Kit	FGFBP1	Q14512	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2507	Mouse FGFBP1(Fibroblast Growth Factor Binding Protein 1) ELISA Kit	FGFBP1	O70514	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2508	Human FSHR(Follicle Stimulating Hormone Receptor) ELISA Kit	FSHR	P23945	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2509	Human TGM1(Transglutaminase 1, Keratinocyte) ELISA Kit	TGM1	P22735	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2510	Human PACAPR(Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) ELISA Kit	PACAPR	P41586	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2511	Mouse PACAPR1(Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) ELISA Kit	PACAPR1	P70205	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2512	Human EPHA10(Ephrin Type A Receptor 10) ELISA Kit	EPHA10	Q5JZY3	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2513	Human ERβ(Estrogen Receptor Beta) ELISA Kit	ERβ	Q92731	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2514	Mouse ERβ(Estrogen Receptor Beta) ELISA Kit	ERβ	O08537	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EH2515	Human ER α (Estrogen Receptor Alpha) ELISA Kit	ER α	P03372	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2516	Mouse ER α (Estrogen Receptor Alpha) ELISA Kit	ER α	P19785	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2517	Human EFNA5(Ephrin A5) ELISA Kit	EFNA5	P52803	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2518	Human MAPKAPK3(Mitogen Activated Protein Kinase Activated Protein Kinase 3)	MAPKAPK3	Q16644	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2519	Human TSH-Ab(Thyroid Stimulating Hormone Antibody) ELISA Kit	TSH-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2520	Mouse TSH-Ab(Thyroid Stimulating Hormone Antibody) ELISA Kit	TSH-Ab		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2521	Human TEF(Thyrotroph Embryonic Factor) ELISA Kit	TEF	Q10587	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2522	Mouse TEF(Thyrotroph Embryonic Factor) ELISA Kit	TEF	Q9JLC6	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2523	Rat TEF(Thyrotroph Embryonic Factor) ELISA Kit	TEF		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2524	Human TRH(Thyrotropin Releasing Hormone) ELISA Kit	TRH		Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM2525	Mouse TRH(Thyrotropin Releasing Hormone) ELISA Kit	TRH	Q62361	Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER2526	Rat TRH(Thyrotropin Releasing Hormone) ELISA Kit	TRH	P01150	Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH2527	Human GHRH(Growth Hormone Releasing Hormone) ELISA Kit	GHRH	P01286	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM2528	Mouse GHRH(Growth Hormone Releasing Hormone) ELISA Kit	GHRH	P16043	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2529	Rat GHRH(Growth Hormone Releasing Hormone) ELISA Kit	GHRH	P09916	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2530	Human δ SIP- α (Delta Sleep Inducing Peptide α) ELISA Kit	δ SIP- α		Competitive	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH2531	Human SCTR(Secretin Receptor) ELISA Kit	SCTR	P47872	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2532	Human SCT(Secretin) ELISA Kit	SCT	P09683	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2533	Mouse SCT(Secretin) ELISA Kit	SCT	Q08535	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2534	Rat SCT(Secretin) ELISA Kit	SCT	P11384	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2535	Human FABP7(Fatty Acid Binding Protein 7, Brain) ELISA Kit	FABP7	O15540	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2536	Human MCP-4/CCL13(Monocyte Chemotactic Protein 4) ELISA Kit	MCP-4/CCL13	Q99616	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2537	Human GPER1(G Protein Coupled Estrogen Receptor 1) ELISA Kit	GPER1	Q99527	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2538	Human ELN(Elastin) ELISA Kit	ELN	P15502	Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-EM2539	Mouse ELN(Elastin) ELISA Kit	ELN	P54320	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2540	Rat ELN(Elastin) ELISA Kit	ELN	Q99372	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2541	Human ELA4(Elastase 4) ELISA Kit	ELA4	Q99895	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM2542	Mouse ELA4(Elastase 4) ELISA Kit	ELA4	Q3SYP2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2543	Rat ELA4(Elastase 4) ELISA Kit	ELA4		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2544	Human PROC(Protein C) ELISA Kit	PROC	P04070	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2545	Mouse PROC(Protein C) ELISA Kit	PROC	P33587	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-ER2546	Rat PROC(Protein C) ELISA Kit	PROC	P31394	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2547	Human PDIA2(Protein Disulfide Isomerase A2) ELISA Kit	PDIA2	Q13087	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2548	Human PDIA3(Protein Disulfide Isomerase A3) ELISA Kit	PDIA3	P30101	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2549	Mouse PDIA3(Protein Disulfide Isomerase A3) ELISA Kit	PDIA3	P27773	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER2550	Rat PDIA3(Protein Disulfide Isomerase A3) ELISA Kit	PDIA3		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2551	Human PRSS2(Protease, Serine, 2) ELISA Kit	PRSS2	P07478	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2552	Human PKBG(Protein Kinase B Gamma) ELISA Kit	PKBG	Q9Y243	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2553	Mouse PKBg(Protein Kinase B Gamma) ELISA Kit	PKBg	Q9WUA6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2554	Rat PKBG(Protein Kinase B Gamma) ELISA Kit	PKBG		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2555	Human PICK1(Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit	PICK1	Q9NRD5	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL

FY-EM2556	Mouse PICK1(Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit	PICK1	Q62083	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2557	Human PKCβ1(Protein Kinase C Beta 1) ELISA Kit	PKCβ1	P05771	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2558	Human PKCγ(Protein Kinase C Gamma) ELISA Kit	PKCγ	P05129	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2559	Human PKCδ1(Protein Kinase C Delta) ELISA Kit	PKCδ1	Q05655	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2560	Human PKCε(Protein Kinase C Epsilon) ELISA Kit	PKCε	Q02156	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2561	Human PKCζ(Protein Kinase C Zeta) ELISA Kit	PKCζ	Q05513	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2562	Human PKCθ(Protein Kinase C Theta) ELISA Kit	PKCθ	Q04759	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2563	Human PRSS23(Protease, Serine 23) ELISA Kit	PRSS23	O95084	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2564	Human PKD1(Protein Kinase D1) ELISA Kit	PKD1	Q15139	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2565	Mouse PKD1(Protein Kinase D1) ELISA Kit	PKD1	Q62101	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2566	Human PKN1(Protein Kinase N1) ELISA Kit	PKN1	Q16512	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2567	Mouse PKN1(Protein Kinase N1) ELISA Kit	PKN1	P70268	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2568	Human PKIα(Protein Kinase Inhibitor Alpha) ELISA Kit	PKIα	P61925	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2569	Human PKIβ(Protein Kinase Inhibitor Beta) ELISA Kit	PKIβ	Q9C010	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2570	Human PKIγ(Protein Kinase Inhibitor Gamma) ELISA Kit	PKIγ	Q9Y2B9	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2571	Human PTPRJ(Protein Tyrosine Phosphatase Receptor Type J) ELISA Kit	PTPRJ	Q12913	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2572	Human PPP1R1A(Protein Phosphatase 1 Regulatory Subunit A) ELISA Kit	PPP1R1A	Q13522	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2573	Human PRSS1(Protease, Serine, 1) ELISA Kit	PRSS1	P07477	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2574	Mouse PRSS1(Protease, Serine, 1) ELISA Kit	PRSS1	Q9Z1R9	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER2575	Rat PRSS1(Protease, Serine, 1) ELISA Kit	PRSS1	P00762	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2576	Human PRSS12(Protease, Serine 12) ELISA Kit	PRSS12	P56730	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2577	Human PR3-Ab(Proteinase 3 Antibody) ELISA Kit	PR3-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2578	Human PAR2(Protease Activated Receptor 2) ELISA Kit	PAR2	P55085	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2579	Mouse PAR2(Protease Activated Receptor 2) ELISA Kit	PAR2	P55086	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2580	Rat PAR2(Protease Activated Receptor 2) ELISA Kit	PAR2	Q63645	Sandwich	500nIU/mL	4.69 nIU/mL	7.82-500nIU/mL
FY-EH2581	Human PCI(Protein C Inhibitor) ELISA Kit	PCI	P05154	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH2582	Human PDI(Protein Disulfide Isomerase) ELISA Kit	PDI	P07237	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2583	Mouse PDI(Protein Disulfide Isomerase) ELISA Kit	PDI	P09103	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2584	Rat PDI(Protein Disulfide Isomerase) ELISA Kit	PDI		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2585	Human PDIA5(Protein Disulfide Isomerase A5) ELISA Kit	PDIA5	Q14554	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2586	Human NGF(Nerve Growth Factor) ELISA Kit	NGF	P01138	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM2587	Mouse NGF(Nerve growth factor) ELISA Kit	NGF	P01139	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2588	Rat NGF(Nerve Growth Factor) ELISA Kit	NGF	P25427	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-EH2589	Human PRSS33(Protease, Serine 33) ELISA Kit	PRSS33	Q8NF86	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2590	Human LDL(Low Density Lipoprotein) ELISA Kit	LDL	P01130	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EM2591	Mouse LDL(Low Density Lipoprotein) ELISA Kit	LDL		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-ER2592	Rat LDL(Low Density Lipoprotein) ELISA Kit	LDL		Sandwich	10μg/mL	0.1 μg/mL	0.16-10μg/mL
FY-EH2593	Human LDLR(Low Density Lipoprotein Receptor) ELISA Kit	LDLR	P01130	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2594	Mouse LDLR(Low Density Lipoprotein Receptor) ELISA Kit	LDLR	P35951	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2595	Rat LDLR(Low Density Lipoprotein Receptor) ELISA Kit	LDLR	P35952	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2596	Human LRP6(Low Density Lipoprotein Receptor Related Protein) ELISA Kit	LRP6	O75581	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EH2597	Human RETN(Resistin) ELISA Kit	RETN	Q9HD89	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2598	Mouse RETN(Resistin) ELISA Kit	RETN	Q99P87	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2599	Rat RETN(Resistin) ELISA Kit	RETN		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2600	Human RPS6Kβ1(Ribosomal Protein S6 Kinase Beta 1) ELISA Kit	RPS6Kβ1	P23443	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2601	Human ETFβ(Electron Transfer Flavoprotein Beta Polypeptide) ELISA Kit	ETFβ	P38117	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2602	Human APP(Amyloid Precursor Protein) ELISA Kit	APP	P05067	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER2603	Rat APP(Amyloid Precursor Protein) ELISA Kit	APP	P08592	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2604	Human M-AChR M1(Muscarinic Acetylcholine Receptor M1) ELISA Kit	M-AChR M1	P11229	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2605	Human M-AChR M3(Muscarinic Acetylcholine Receptor M3) ELISA Kit	M-AChR M3	P20309	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2606	Human D1R/DRD1(Dopamine Receptor D1) ELISA Kit	D1R/DRD1	P21728	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2607	Human D2R/DRD2(Dopamine Receptor D2) ELISA Kit	D2R/DRD2	P14416	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2608	Human D3R/DRD3(Dopamine Receptor D3) ELISA Kit	D3R/DRD3	P35462	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2609	Human D4R/DRD4(Dopamine Receptor D4) ELISA Kit	D4R/DRD4	P21917	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2610	Human D5R/DRD5(Dopamine Receptor D5) ELISA Kit	D5R/DRD5	P21918	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2611	Human DDC(Dopamine Decarboxylase) ELISA Kit	DDC	P20711	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2612	Human DAT(Dopamine Transporter) ELISA Kit	DAT	Q01959	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2613	Mouse DAT(Dopamine Transporter) ELISA Kit	DAT	Q61327	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER2614	Rat DAT(Dopamine Transporter) ELISA Kit	DAT	P23977	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2615	Human DT(Dopachrome Tautomerase) ELISA Kit	DT	P40126	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2616	Human PIGR/SC(Polymeric Immunoglobulin Receptor/Membrane Secretory Comp	PIGR/SC	P01833	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2617	Human PTN(Pleiotrophin) ELISA Kit	PTN	P21246	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2618	Mouse PTN(Pleiotrophin) ELISA Kit	PTN	P63089	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2619	Rat PTN(Pleiotrophin) ELISA Kit	PTN	P63090	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2620	Human DAO(Diamine Oxidase) ELISA Kit	DAO	P19801	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM2621	Mouse DAO(Diamine Oxidase) ELISA Kit	DAO	Q8JZQ5	Sandwich	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-ER2622	Rat DAO(Diamine Oxidase) ELISA Kit	DAO	P36633	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2623	Human UGT2B4(Uridine Diphosphate Glucuronosyltransferase 2 Family Polypeptid	UGT2B4	P06133	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2624	Human DLD(Dihydropolyl Dehydrogenase) ELISA Kit	DLD	P09622	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EM2625	Mouse DLD(Dihydropolyl Dehydrogenase) ELISA Kit	DLD	O08749	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2626	Rat DLD(Dihydropolyl Dehydrogenase) ELISA Kit	DLD	Q6P6R2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2627	Human DLAT(Dihydropolyl Transacetylase) ELISA Kit	DLAT	P10515	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2628	Mouse DLAT(Dihydropolyl Transacetylase) ELISA Kit	DLAT	Q8BMF4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2629	Rat DLAT(Dihydropolyl Transacetylase) ELISA Kit	DLAT		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2630	Human DPYSL2(Dihydropyrimidinase Like Protein 2) ELISA Kit	DPYSL2	Q16555	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2631	Human DNM1(Dynamin 1) ELISA Kit	DNM1	Q05193	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2632	Mouse DNM1(Dynamin 1) ELISA Kit	DNM1	P39053	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER2633	Rat DNM1(Dynamin 1) ELISA Kit	DNM1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2634	Human DNM2(Dynamin 2) ELISA Kit	DNM2	P50570	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2635	Human Ub(Ubiquitin) ELISA Kit	Ub	G5EFF7	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2636	Mouse Ub(Ubiquitin) ELISA Kit	Ub		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2637	Human UBD(Ubiquitin D) ELISA Kit	UBD	O15205	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL

FY-EH2638	Human UBE1L(Ubiquitin Activating Enzyme E1 Like Protein) ELISA Kit	UBE1L	P41226	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH2639	Human UBE2C(Ubiquitin Conjugating Enzyme E2C) ELISA Kit	UBE2C	O00762	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2640	Human UBE3A(Ubiquitin Protein Ligase E3A) ELISA Kit	UBE3A	Q05086	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2641	Human NNE(Non-Neuronal Enolase) ELISA Kit	NNE	P06733	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2642	Rat NNE(Non-Neuronal Enolase) ELISA Kit	NNE	P04764	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2643	Human SDC3(Syndecan 3) ELISA Kit	SDC3	O75056	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2644	Human FBLN5(Fibulin 5) ELISA Kit	FBLN5	Q9UBX5	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2645	Human REG1β(Regenerating Islet Derived Protein 1 Beta) ELISA Kit	REG1β	P48304	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2646	Human SPC(Pulmonary Surfactant Associated Protein C) ELISA Kit	SPC	P11686	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2647	Rat SPC(Pulmonary Surfactant Associated Protein C) ELISA Kit	SPC	P11685	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2648	Human PARC/CCL18(Pulmonary Activation Regulated Chemokine) ELISA Kit	PARC/CCL18	P55774	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2649	Human LRP(Lung Resistance-Related Protein) ELISA Kit	LRP	Q14764	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2650	Human PRSS8(Protease, Serine 8) ELISA Kit	PRSS8	Q16651	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2651	Human SPARCL1(SPARC Like Protein 1) ELISA Kit	SPARCL1	Q14515	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2652	Mouse SPARCL1(SPARC Like Protein 1) ELISA Kit	SPARCL1	P70663	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2653	Human SEMA3A(Semaphorin 3A) ELISA Kit	SEMA3A	Q14563	Sandwich	10ng/mL	0.10 ng/mL	0.16-10ng/mL
FY-ER2654	Rat SEMA3A(Semaphorin 3A) ELISA Kit	SEMA3A	Q63548	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2655	Human slgA(Secretory Immunoglobulin A) ELISA Kit	slgA		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2656	Mouse slgA(Secretory Immunoglobulin A) ELISA Kit	slgA		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2657	Rat slgA(Secretory Immunoglobulin A) ELISA Kit	slgA		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2658	Human SLPI(Secretory Leukocyte Protease Inhibitor) ELISA Kit	SLPI	P03973	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER2659	Rat SLPI(Secretory Leukocyte Protease Inhibitor) ELISA Kit	SLPI		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2660	Human SAP(Serum Amyloid P Component) ELISA Kit	SAP	P02743	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-ER2661	Rat SAP(Serum Amyloid P Component) ELISA Kit	SAP		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2662	Human PAR4(Protease Activated Receptor 4) ELISA Kit	PAR4	Q96R10	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2663	Human GPR120(G Protein Coupled Receptor 120) ELISA Kit	GPR120	Q5NUL3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM2664	Mouse GPR120(G Protein Coupled Receptor 120) ELISA Kit	GPR120	Q7TMA4	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2665	Human REG3α(Regenerating Islet Derived Protein 3 Alpha) ELISA Kit	REG3α	Q06141	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2666	Human Cyr61(Cysteine Rich Protein, Angiogenic Inducer 61) ELISA Kit	Cyr61	O00622	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2667	Human REG3γ(Regenerating Islet Derived Protein 3 Gamma) ELISA Kit	REG3γ	Q6UW15	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2668	Human S100A8(S100 Calcium Binding Protein A8) ELISA Kit	S100A8	P05109	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM2669	Mouse S100A8(S100 Calcium Binding Protein A8) ELISA Kit	S100A8	P27005	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2670	Human S100A9(S100 Calcium Binding Protein A9) ELISA Kit	S100A9	P06702	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2671	Mouse S100A9(S100 Calcium Binding Protein A9) ELISA Kit	S100A9	P31725	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER2672	Rat S100A9(S100 Calcium Binding Protein A9) ELISA Kit	S100A9	P50116	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2673	Human GPR37(G Protein Coupled Receptor 37) ELISA Kit	GPR37	O15354	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2674	Human S100A11(S100 Calcium Binding Protein A11) ELISA Kit	S100A11	P31949	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2675	Mouse S100A11(S100 Calcium Binding Protein A11) ELISA Kit	S100A11	P50543	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER2676	Rat S100A11(S100 Calcium Binding Protein A11) ELISA Kit	S100A11	Q6B345	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2677	Human S100A12(S100 Calcium Binding Protein A12) ELISA Kit	S100A12	P80511	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2678	Human S100A13(S100 Calcium Binding Protein A13) ELISA Kit	S100A13	Q99584	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EH2679	Human S100A6(S100 Calcium Binding Protein A6) ELISA Kit	S100A6	P06703	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2680	Mouse S100A6(S100 Calcium Binding Protein A6) ELISA Kit	S100A6	P14069	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2681	Rat S100A6(S100 Calcium Binding Protein A6) ELISA Kit	S100A6	P05964	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2682	Human S100A7(S100 Calcium Binding Protein A7) ELISA Kit	S100A7	P31151	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2683	Human S100B(S100 Calcium Binding Protein B) ELISA Kit	S100B	P04271	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2684	Mouse S100B(S100 Calcium Binding Protein B) ELISA Kit	S100B	P50114	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER2685	Rat S100B(S100 Calcium Binding Protein B) ELISA Kit	S100B	P04631	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2686	Human SDC1(Syndecan 1) ELISA Kit	SDC1	P18827	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2687	Mouse SDC1(Syndecan 1) ELISA Kit	SDC1	P18828	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2688	Rat SDC1(Syndecan 1) ELISA Kit	SDC1	P26260	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2689	Human SCF(Stem Cell Factor) ELISA Kit	SCF	P21583	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2690	Mouse SCF(Stem Cell Factor) ELISA Kit	SCF	P20826	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2691	Rat SCF(Stem Cell Factor) ELISA Kit	SCF	P21581	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2692	Human SCFR(Stem Cell Growth Factor Receptor) ELISA Kit	SCFR		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2693	Mouse SCFR(Stem Cell Growth Factor Receptor) ELISA Kit	SCFR	P05532	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2694	Rat SCFR(Stem Cell Growth Factor Receptor) ELISA Kit	SCFR		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2695	Human GAL(Galanin) ELISA Kit	GAL	P22466	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2696	Rat GAL(Galanin) ELISA Kit	GAL	P10683	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2697	Human GALR1(Galanin Receptor 1) ELISA Kit	GALR1	P47211	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2698	Mouse GALR1(Galanin Receptor 1) ELISA Kit	GALR1	P56479	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2699	Rat GALR1(Galanin Receptor 1) ELISA Kit	GALR1		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2700	Human PGAM1(Phosphoglycerate Mutase 1) ELISA Kit	PGAM1	P18669	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2701	Human PK(Pyruvate Kinase, Liver And RBC) ELISA Kit	PK	P30613	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2702	Mouse PK(Pyruvate Kinase, Liver And RBC) ELISA Kit	PK	P53657	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER2703	Rat PK(Pyruvate Kinase, Liver And RBC) ELISA Kit	PK	P12928	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2704	Human PYGL(Glycogen Phosphorylase, Liver) ELISA Kit	PYGL	P06737	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2705	Mouse PYGL(Glycogen Phosphorylase, Liver) ELISA Kit	PYGL	Q9ET01	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2706	Rat PYGL(Glycogen Phosphorylase, Liver) ELISA Kit	PYGL	P09811	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2707	Human SFRP2(Secreted Frizzled Related Protein 2) ELISA Kit	SFRP2	Q96HF1	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH2708	Human GP-73(Golgi Protein 73) ELISA Kit	GP-73	Q8NBJ4	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2709	Human SIAE(Sialic Acid Acetyltransferase) ELISA Kit	SIAE	Q9HAT2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2710	Human ST2(Syntenin 2) ELISA Kit	ST2	Q9H190	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2711	Human SIGLEC1(Sialic Acid Binding Immunoglobulin Like Lectin 1) ELISA Kit	SIGLEC1	Q9BZZ2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2712	Human SIGLEC10(Sialic Acid Binding Immunoglobulin Like Lectin 10) ELISA Kit	SIGLEC10	Q96LC7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2713	Human SIGLEC2(Sialic Acid Binding Immunoglobulin Like Lectin 2) ELISA Kit	SIGLEC2	P20273	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2714	Human PAD14(Peptidyl Arginine Deiminase Type IV) ELISA Kit	PAD14	Q9UM07	Sandwich	12ng/mL	0.11 ng/mL	0.19-12ng/mL
FY-EH2715	Human SIGLEC8(Sialic Acid Binding Immunoglobulin Like Lectin 8) ELISA Kit	SIGLEC8	Q9NYZ4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2716	Human Smad3/MADH3(Mothers Against Decapentaplegic Homolog 3) ELISA Kit	Smad3/MADH3	P84022	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER2717	Rat Smad3/MADH3(Mothers Against Decapentaplegic Homolog 3) ELISA Kit	Smad3/MADH3		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2718	Human OFQ/N(Orphanin FQ/Nociceptin) ELISA Kit	OFQ/N	Q13519	Sandwich	300pg/mL	2.81 pg/mL	4.69-300pg/mL
FY-EM2719	Mouse OFQ/N(Orphanin FQ/Nociceptin) ELISA Kit	OFQ/N		Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL

FY-ER2720	Rat OFQ/N(Orphanin FQ/Nociceptin) ELISA Kit	OFQ/N		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2721	Human GCLM(Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit	GCLM	P48507	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2722	Mouse GCLM(Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit	GCLM	O09172	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2723	Rat GCLM(Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit	GCLM	P48508	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2724	Human GSTθ1/GSTt1(Glutathione S Transferase Theta 1) ELISA Kit	GSTθ1/GSTt1	P30711	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2725	Human GSTθ2/GSTt2(Glutathione S Transferase Theta 2) ELISA Kit	GSTθ2/GSTt2	P0CG29	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH2726	Human GSTω1(Glutathione S Transferase Omega 1) ELISA Kit	GSTω1	P78417	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2727	Human HDAC2(Histone Deacetylase 2) ELISA Kit	HDAC2	Q92769	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2728	Mouse HDAC2(Histone Deacetylase 2) ELISA Kit	HDAC2	P70288	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2729	Rat HDAC2(Histone Deacetylase 2) ELISA Kit	HDAC2	B1WB58	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2730	Human C1QTNF1(C1q and Tumor Necrosis Factor Related Protein 1) ELISA Kit	C1QTNF1	Q9BXJ1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM2731	Mouse GR(glucocorticoid receptor) ELISA Kit	GR	P06537	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2732	Rat GR(Glucocorticoid Receptor) ELISA Kit	GR	P70619	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2733	Human GSTpi(Glutathione S Transferases Pi) ELISA Kit	GSTpi	P09211	Sandwich	150ng/mL	1.41 ng/mL	2.35-150ng/mL
FY-EH2734	Human OPG(Osteoprotegerin) ELISA Kit	OPG	O00300	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2735	Mouse OPG(Osteoprotegerin) ELISA Kit	OPG	O08712	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER2736	Rat OPG(Osteoprotegerin) ELISA Kit	OPG	O08727	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2737	Human OC/BGP(Osteocalcin) ELISA Kit	OC/BGP	P02818	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EM2738	Mouse OC/BGP(Osteocalcin) ELISA Kit	OC/BGP	P86546	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER2739	Rat OC/BGP(Osteocalcin) ELISA Kit	OC/BGP	P04640	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2740	Human ELA3B(Elastase 3B) ELISA Kit	ELA3B	P08861	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2741	Human MUSK(Muscle Skeletal Receptor Tyrosine Kinase) ELISA Kit	MUSK	O15146	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2742	Human GSTα5(Glutathione S Transferase Alpha 5) ELISA Kit	GSTα5	Q7RTV2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2743	Human OPN(Osteopontin) ELISA Kit	OPN	P10451	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2744	Mouse OPN(Osteopontin) ELISA Kit	OPN	P10923	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2745	Rat OPN(Osteopontin) ELISA Kit	OPN	P08721	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2746	Human MPIF-1/CCL23(Myeloid Progenitor Inhibitory Factor 1) ELISA Kit	MPIF-1/CCL23	P55773	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2747	Human SOX9(Sex Determining Region Y Box Protein 9) ELISA Kit	SOX9	P48436	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2748	Human GREM1(Gremlin 1) ELISA Kit	GREM1	O60565	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2749	Mouse GREM1(GremLin 1) ELISA Kit	GREM1	O70326	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER2750	Rat GREM1(GremLin 1) ELISA Kit	GREM1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2751	Human ON(Osteonectin) ELISA Kit	ON	P09486	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2752	Mouse ON(Osteonectin) ELISA Kit	ON	P07214	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER2753	Rat ON(Osteonectin) ELISA Kit	ON	P09486	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2754	Human SF1(Steroidogenic Factor 1) ELISA Kit	SF1	Q13285	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2755	Mouse SF1(Steroidogenic Factor 1) ELISA Kit	SF1	P33242	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2756	Rat SF1(Steroidogenic Factor 1) ELISA Kit	SF1		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2757	Human PAPP2(Pappalysin 2) ELISA Kit	PAPP2	Q9BXP8	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2758	Human P-CK(Pan-Cytokeratin) ELISA Kit	P-CK		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2759	Human ALDOC(Aldolase C, Fructose Bisphosphate) ELISA Kit	ALDOC	P09972	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2760	Human PMP3(Peroxisomal Membrane Protein 3) ELISA Kit	PMP3	P28328	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL

FY-EM2761	Mouse PMP3(Peroxisomal Membrane Protein 3) ELISA Kit	PMP3	P55098	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2762	Human PPAR γ C1 α (Peroxisome Proliferator Activated Receptor Gamma Coactivator	PPAR γ C1 α	Q9UBK2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2763	Mouse PPAR γ C1 α (Peroxisome Proliferator Activated Receptor Gamma Coactivator	PPAR γ C1 α	O70343	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2764	Human PPAR- α (Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit	PPAR- α	Q07869	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2765	Mouse PPAR- α (Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit	PPAR- α	P23204	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2766	Rat PPAR- α (Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit	PPAR- α	P37230	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2767	Human PPAR- γ (Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit	PPAR- γ	P37231	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2768	Mouse PPAR- γ (Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit	PPAR- γ	P37238	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2769	Rat PPAR- γ (Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit	PPAR- γ	O88275	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2770	Human NXF1(Nuclear RNA Export Factor 1) ELISA Kit	NXF1	Q9UBU9	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2771	Mouse NXF1(Nuclear RNA Export Factor 1) ELISA Kit	NXF1	Q99JX7	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2772	Rat NXF1(Nuclear RNA Export Factor 1) ELISA Kit	NXF1		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2773	Human NMP-22(Nuclear Matrix Protein 22) ELISA Kit	NMP-22	Q14980	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2774	Human NUP107(Nucleoporin 107kDa) ELISA Kit	NUP107	P57740	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2775	Mouse NUP107(Nucleoporin 107kDa) ELISA Kit	NUP107	Q8BH74	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2776	Rat NUP107(Nucleoporin 107kDa) ELISA Kit	NUP107		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2777	Human NUP133(Nucleoporin 133kDa) ELISA Kit	NUP133	Q8WUM0	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2778	Mouse NUP133(Nucleoporin 133kDa) ELISA Kit	NUP133	Q8R0G9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2779	Human NUP155(Nucleoporin 155kDa) ELISA Kit	NUP155	O75694	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2780	Mouse NUP155(Nucleoporin 155kDa) ELISA Kit	NUP155	Q99P88	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2781	Rat NUP155(Nucleoporin 155kDa) ELISA Kit	NUP155		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2782	Human NUP188(Nucleoporin 188kDa) ELISA Kit	NUP188	Q5SRE5	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2783	Mouse NUP188(Nucleoporin 188kDa) ELISA Kit	NUP188	Q6ZQH8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2784	Rat NUP188(Nucleoporin 188kDa) ELISA Kit	NUP188		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2785	Human NUP205(Nucleoporin 205kDa) ELISA Kit	NUP205	Q92621	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2786	Mouse NUP205(Nucleoporin 205kDa) ELISA Kit	NUP205	E9Q880	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2787	Rat NUP205(Nucleoporin 205kDa) ELISA Kit	NUP205		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2788	Human NUP214(Nucleoporin 214kDa) ELISA Kit	NUP214	P35658	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2789	Mouse NUP214(Nucleoporin 214kDa) ELISA Kit	NUP214	Q80U93	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2790	Rat NUP214(Nucleoporin 214kDa) ELISA Kit	NUP214		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2791	Human NUP35(Nucleoporin 35kDa) ELISA Kit	NUP35	Q8NFB5	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2792	Mouse NUP35(Nucleoporin 35kDa) ELISA Kit	NUP35	Q8R4R6	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER2793	Rat NUP35(Nucleoporin 35kDa) ELISA Kit	NUP35		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2794	Human NUP37(Nucleoporin 37kDa) ELISA Kit	NUP37	Q8NFB4	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2795	Mouse NUP37(Nucleoporin 37kDa) ELISA Kit	NUP37	Q9CWU9	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2796	Rat NUP37(Nucleoporin 37kDa) ELISA Kit	NUP37		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2797	Human NUP50(Nucleoporin 50kDa) ELISA Kit	NUP50	Q9UKX7	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2798	Mouse NUP50(Nucleoporin 50kDa) ELISA Kit	NUP50	Q9JIH2	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2799	Human NUP85(Nucleoporin 85kDa) ELISA Kit	NUP85	Q9BW27	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2800	Mouse NUP85(Nucleoporin 85kDa) ELISA Kit	NUP85	Q8R480	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2801	Rat NUP85(Nucleoporin 85kDa) ELISA Kit	NUP85		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH2802	Human NUP98(Nucleoporin 98kDa) ELISA Kit	NUP98	P52948(1~8)	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2803	Mouse NUP98(Nucleoporin 98kDa) ELISA Kit	NUP98	Q6PFD9(1~5)	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2804	Human NRIP1(Nuclear Receptor Interacting Protein 1) ELISA Kit	NRIP1	P48552	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2805	Mouse NRIP1(Nuclear Receptor Interacting Protein 1) ELISA Kit	NRIP1	Q8CBD1	Sandwich	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-ER2806	Rat NRIP1(Nuclear Receptor Interacting Protein 1) ELISA Kit	NRIP1	G3V672	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2807	Human RNASE3/ECP(Ribonuclease A3/Eosinophil Cationic Protein) ELISA Kit	RNASE3/ECP	P12724	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-EH2808	Human RI(Ribonuclease Inhibitor) ELISA Kit	RI	P13489	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH2809	Human RNASE(Ribonuclease) ELISA Kit	RNASE	P07998	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2810	Human RPL10(Ribosomal Protein L10) ELISA Kit	RPL10	P27635	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM2811	Mouse RPL10(Ribosomal Protein L10) ELISA Kit	RPL10	Q6ZWW3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2812	Rat RPL10(Ribosomal Protein L10) ELISA Kit	RPL10		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2813	Human RPL26(Ribosomal Protein L26) ELISA Kit	RPL26	P61254	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM2814	Mouse RPL26(Ribosomal Protein L26) ELISA Kit	RPL26	P61255	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2815	Rat RPL26(Ribosomal Protein L26) ELISA Kit	RPL26		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2816	Human MART/Melan-A(Melanoma Marker A) ELISA Kit	MART/Melan-A	Q16655	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2817	Human MAGE(Melanoma-associated Antigen) ELISA Kit	MAGE	P43364	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2818	Human MLZE(Melanoma Derived Leucine Zipper Extra Nuclear Factor) ELISA Kit	MLZE	Q9BYG8	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2819	Human THRSP(Thyroid Hormone Responsive) ELISA Kit	THRSP	Q92748	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2820	Human EPOR(Erythropoietin Receptor) ELISA Kit	EPOR	P19235	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2821	Mouse EPOR(Erythropoietin receptor) ELISA Kit	EPOR	P14753	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2822	Human EXT1(Exostoses 1) ELISA Kit	EXT1	Q16394	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2823	Human GSTk1(Glutathione S Transferase Kappa 1) ELISA Kit	GSTk1	Q9Y2Q3	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH2824	Human LTA4H(Leukotriene A4 Hydrolase) ELISA Kit	LTA4H	P09960	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2825	Human GSTα4(Glutathione S Transferase Alpha 4) ELISA Kit	GSTα4	O15217	Sandwich	200mIU/mL	1.88 mIU/mL	3.13-200mIU/mL
FY-EH2826	Human SNCγ(Synuclein Gamma) ELISA Kit	SNCγ	O76070	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2827	Human TXLNa(Taxilin Alpha) ELISA Kit	TXLNa	P40222	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2828	Mouse TXLNa(Taxilin Alpha) ELISA Kit	TXLNa	Q6PAM1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2829	Rat TXLNa(Taxilin Alpha) ELISA Kit	TXLNa		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2830	Human TXLNb(Taxilin Beta) ELISA Kit	TXLNb	Q8N3L3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2831	Mouse TXLNb(Taxilin Beta) ELISA Kit	TXLNb	Q8VBT1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2832	Rat TXLNB(Taxilin Beta) ELISA Kit	TXLNB		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2833	Human Sec23b(Protein transport protein Sec23B) ELISA Kit	Sec23b	Q15437	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2834	Human RNF4(Ring Finger Protein 4) ELISA Kit	RNF4	P78317	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2835	Mouse RNF4(Ring Finger Protein 4) ELISA Kit	RNF4	Q9QZS2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2836	Rat RNF4(Ring Finger Protein 4) ELISA Kit	RNF4		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2837	Human FABP6(Fatty Acid Binding Protein 6, Ileal) ELISA Kit	FABP6	P51161	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2838	Mouse FABP6(Fatty Acid Binding Protein 6, Ileal) ELISA Kit	FABP6	P51162	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2839	Human TMPRSS2(Transmembrane Protease, Serine 2) ELISA Kit	TMPRSS2	O15393	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2840	Human TNNI1(Troponin I Type 1, Slow Skeletal) ELISA Kit	TNNI1	P19237	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-EH2841	Human FSCN(Fascin) ELISA Kit	FSCN	Q16658	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2842	Mouse FSCN(Fascin) ELISA Kit	FSCN	Q61553	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH2843	Human TNNC2(Troponin C Type 2) ELISA Kit	TNNC2	P02585	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2844	Human TNC(Tenascin C) ELISA Kit	TNC	P24821	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2845	Mouse TNC(Tenascin C) ELISA Kit	TNC	Q80YX1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2846	Rat TNC(Tenascin C) ELISA Kit	TNC	A0A0G2K1I	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2847	Human TNR(Tenascin R) ELISA Kit	TNR	Q92752	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2848	Human TNX(Tenascin X) ELISA Kit	TNX	P22105	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2849	Mouse TNX(Tenascin X) ELISA Kit	TNX	Q3V5L4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2850	Rat TNX(Tenascin X) ELISA Kit	TNX		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2851	Human IgA1(Immunoglobulin A1) ELISA Kit	IgA1	P01876	Sandwich	3000ng/mL	28.13 ng/mL	46.88-3000ng/mL
FY-EH2852	Human MYL1(Myosin Light Chain 1, Alkali, Fast Skeletal) ELISA Kit	MYL1	P05976	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH2853	Human MLCK/MYLK(Myosin Light Chain Kinase) ELISA Kit	MLCK/MYLK	Q15746	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2854	Mouse MLCK/MYLK(Myosin light chain kinase) ELISA Kit	MLCK/MYLK	Q6PDN3	Sandwich	1500pg/mL	14.06 pg/mL	23.44-1500pg/mL
FY-EH2855	Human MYH1(Myosin Heavy Chain 1, Skeletal Muscle, Adult) ELISA Kit	MYH1	P12882	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2856	Human FABP3(Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit	FABP3	P05413	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2857	Mouse FABP3(Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit	FABP3	P11404	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER2858	Rat FABP3(Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit	FABP3	P07483	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2859	Human PYGM(Glycogen Phosphorylase, Muscle) ELISA Kit	PYGM	P11217	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2860	Mouse PYGM(Glycogen Phosphorylase, Muscle) ELISA Kit	PYGM	Q9WUB3	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2861	Rat PYGM(Glycogen Phosphorylase, Muscle) ELISA Kit	PYGM	P09812	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2862	Human CKM(Creatine Kinase, Muscle) ELISA Kit	CKM	P06732	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2863	Mouse CKM(Creatine Kinase, Muscle) ELISA Kit	CKM	P07310	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER2864	Rat CKM(Creatine Kinase, Muscle) ELISA Kit	CKM	P00564	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2865	Human MYOG(Myogenin) ELISA Kit	MYOG	P15173	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2866	Human DMD(Dystrophin) ELISA Kit	DMD	P11532	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2867	Mouse DMD(Dystrophin) ELISA Kit	DMD	P11531	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2868	Rat DMD(Dystrophin) ELISA Kit	DMD	P11530	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2869	Human MSTN(Myostatin) ELISA Kit	MSTN	O14793	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2870	Human PUMA(P53 Upregulated Modulator of Apoptosis) ELISA Kit	PUMA	Q9BXH1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2871	Rat PUMA(P53 Upregulated Modulator of Apoptosis) ELISA Kit	PUMA	Q80ZG6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2872	Human MGP(Matrix Gla Protein) ELISA Kit	MGP	P08493	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2873	Human MMP-11(Matrix Metalloproteinase 11) ELISA Kit	MMP-11	P24347	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2874	Mouse MMP-11(Matrix Metalloproteinase 11) ELISA Kit	MMP-11	Q02853	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2875	Rat MMP-11(Matrix Metalloproteinase 11) ELISA Kit	MMP-11		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2876	Human MGST1(Microsomal Glutathione S Transferase 1) ELISA Kit	MGST1	P10620	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH2877	Human MMP-2(Matrix Metalloproteinase 2) ELISA Kit	MMP-2	P08253	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2878	Mouse MMP-2(Matrix Metalloproteinase 2) ELISA Kit	MMP-2	P33434	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2879	Rat MMP-2(Matrix Metalloproteinase 2) ELISA Kit	MMP-2	P33436	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2880	Human MMP-3(Matrix Metalloproteinase 3) ELISA Kit	MMP-3	P08254	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2881	Mouse MMP-3(Matrix Metalloproteinase 3) ELISA Kit	MMP-3	P28862	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2882	Rat MMP-3(Matrix Metalloproteinase 3) ELISA Kit	MMP-3	P03957	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2883	Human MMP-17(Matrix Metalloproteinase 17) ELISA Kit	MMP-17	Q9ULZ9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EM2884	Mouse MMP-17(Matrix Metalloproteinase 17) ELISA Kit	MMP-17	Q9R0S3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2885	Rat MMP-17(Matrix Metalloproteinase 17) ELISA Kit	MMP-17		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2886	Human MMP-24(Matrix Metalloproteinase 24) ELISA Kit	MMP-24	Q9Y5R2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2887	Mouse MMP-24(Matrix Metalloproteinase 24) ELISA Kit	MMP-24	Q9R0S2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2888	Rat MMP-24(Matrix Metalloproteinase 24) ELISA Kit	MMP-24		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2889	Human MMP-7(Matrix Metalloproteinase 7) ELISA Kit	MMP-7	P09237	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2890	Mouse MMP-7(Matrix Metalloproteinase 7) ELISA Kit	MMP-7	Q10738	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER2891	Rat MMP-7(Matrix Metalloproteinase 7) ELISA Kit	MMP-7	P50280	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2892	Human MMP-8(Matrix Metalloproteinase 8) ELISA Kit	MMP-8	P22894	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2893	Mouse MMP-8(Matrix Metalloproteinase 8) ELISA Kit	MMP-8	O70138	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2894	Rat MMP-8(Matrix Metalloproteinase 8) ELISA Kit	MMP-8	O88766	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2895	Human ELA2B(Elastase 2B) ELISA Kit	ELA2B	P08218	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2896	Human TIMP-2(Tissue Inhibitors of Metalloproteinase 2) ELISA Kit	TIMP-2	P16035	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2897	Mouse TIMP-2(Tissue Inhibitors of Metalloproteinase 2) ELISA Kit	TIMP-2	P25785	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2898	Rat TIMP-2(Tissue Inhibitors of Metalloproteinase 2) ELISA Kit	TIMP-2	P30121	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2899	Human TIMP-3(Tissue Inhibitors of Metalloproteinase 3) ELISA Kit	TIMP-3	P35625	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2900	Mouse TIMP-3(Tissue Inhibitors of Metalloproteinase 3) ELISA Kit	TIMP-3	P39876	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2901	Rat TIMP-3(Tissue Inhibitors of Metalloproteinase 3) ELISA Kit	TIMP-3	P48032	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2902	Human TIMP-4(Tissue Inhibitors of Metalloproteinase 4) ELISA Kit	TIMP-4	Q99727	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM2903	Mouse TIMP-4(Tissue Inhibitors of Metalloproteinase 4) ELISA Kit	TIMP-4	Q9JHB3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER2904	Rat TIMP-4(Tissue Inhibitors of Metalloproteinase 4) ELISA Kit	TIMP-4	P81556	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2905	Human LGALS3BP(Lectin Galactoside Binding, Soluble 3 Binding Protein) ELISA Kit	LGALS3BP	Q08380	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH2906	Human TNNI2(Troponin I Type 2, Fast Skeletal) ELISA Kit	TNNI2	P48788	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM2907	Mouse TNNI2(Troponin I Type 2, Fast Skeletal) ELISA Kit	TNNI2	P13412	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2908	Human TP/TXA2R(Thromboxane Receptor) ELISA Kit	TP/TXA2R	P21731	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2909	Human SDF-2(Stromal Cell Derived Factor 2) ELISA Kit	SDF-2	Q99470	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2910	Mouse SDF-2(Stromal Cell Derived Factor 2) ELISA Kit	SDF-2	Q9DCT5	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2911	Rat SDF-2(Stromal Cell Derived Factor 2) ELISA Kit	SDF-2		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2912	Human SDF2L1(Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit	SDF2L1	Q9HCN8	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2913	Mouse SDF2L1(Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit	SDF2L1	Q9ESP1	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2914	Rat SDF2L1(Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit	SDF2L1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2915	Human SDF-4(Stromal Cell Derived Factor 4) ELISA Kit	SDF-4	Q9BRK5	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2916	Mouse SDF-4(Stromal Cell Derived Factor 4) ELISA Kit	SDF-4	Q61112	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2917	Rat SDF-4(Stromal Cell Derived Factor 4) ELISA Kit	SDF-4		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2918	Human TDGF1(Teratocarcinoma Derived Growth Factor 1) ELISA Kit	TDGF1	P13385	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2919	Rat TDGF1(Teratocarcinoma Derived Growth Factor 1) ELISA Kit	TDGF1		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2920	Human KLK8(Kallikrein 8) ELISA Kit	KLK8	O60259	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2921	Human VLDLR(Very Low Density Lipoprotein Receptor) ELISA Kit	VLDLR	P98155	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH2922	Human GAL3(Galectin 3) ELISA Kit	GAL3	P17931	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2923	Mouse GAL3(Galectin 3) ELISA Kit	GAL3	P16110	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2924	Rat GAL3(Galectin 3) ELISA Kit	GAL3	P08699	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL

FY-EH2925	Human FPR2(Formyl Peptide Receptor 2) ELISA Kit	FPR2	P25090	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2926	Mouse FPR2(Formyl Peptide Receptor 2) ELISA Kit	FPR2	O88536	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2927	Human PTHR2(Parathyroid Hormone Receptor 2) ELISA Kit	PTHR2	P49190	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2928	Mouse PTHR2(Parathyroid Hormone Receptor 2) ELISA Kit	PTHR2	Q91V95	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER2929	Rat PTHR2(Parathyroid Hormone Receptor 2) ELISA Kit	PTHR2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2930	Human PTHrP(Parathyroid Hormone Related Protein) ELISA Kit	PTHrP	P12272	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2931	Rat PTHrP(Parathyroid Hormone Related Protein) ELISA Kit	PTHrP	P13085	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2932	Human TPP1(Tripeptidyl Peptidase I) ELISA Kit	TPP1	O14773	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2933	Human CD276(Cluster of Differentiation 276) ELISA Kit	CD276	Q5ZPR3	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2934	Human TBG(Thyroxine Binding Globulin) ELISA Kit	TBG	P05543	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM2935	Mouse TBG(Thyroxine Binding Globulin) ELISA Kit	TBG	P61939	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2936	Human PAH(Phenylalanine Hydroxylase) ELISA Kit	PAH	P00439	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2937	Human GJβ1(Gap Junction Protein Beta 1) ELISA Kit	GJβ1	P08034	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2938	Human PCT(Procalcitonin) ELISA Kit	PCT	P01258	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM2939	Mouse PCT(Procalcitonin) ELISA Kit	PCT	P70160	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2940	Rat PCT(Procalcitonin) ELISA Kit	PCT		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2941	Human PCPE1(Procollagen C Proteinase Enhancer 1) ELISA Kit	PCPE1	Q15113	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH2942	Human GDNF(Glial Cell Line Derived Neurotrophic Factor) ELISA Kit	GDNF	P39905	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2943	Mouse GDNF(Glial Cell Line Derived Neurotrophic Factor) ELISA Kit	GDNF	P48540	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2944	Rat GDNF(Glial Cell Line Derived Neurotrophic Factor) ELISA Kit	GDNF	Q07731	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH2945	Human CDSN(Corneodesmosin) ELISA Kit	CDSN	Q15517	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2946	Human CNTN1(Contactin 1) ELISA Kit	CNTN1	Q12860	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2947	Mouse CNTN1(Contactin 1) ELISA Kit	CNTN1	P12960	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2948	Rat CNTN1(Contactin 1) ELISA Kit	CNTN1	Q63198	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2949	Human CNTN2(Contactin 2) ELISA Kit	CNTN2	Q02246	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2950	Mouse CNTN2(Contactin 2) ELISA Kit	CNTN2	Q61330	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2951	Rat CNTN2(Contactin 2) ELISA Kit	CNTN2		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2952	Human CNTN3(Contactin 3) ELISA Kit	CNTN3	Q9P232	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM2953	Mouse CNTN3(Contactin 3) ELISA Kit	CNTN3	Q07409	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2954	Rat CNTN3(Contactin 3) ELISA Kit	CNTN3		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2955	Human CNTNAP1(Contactin Associated Protein 1) ELISA Kit	CNTNAP1	P78357	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2956	Human CNTNAP4(Contactin Associated Protein 4) ELISA Kit	CNTNAP4	Q9C0A0	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2957	Human ELMO2(Engulfment And Cell Motility 2) ELISA Kit	ELMO2	Q96JJ3	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM2958	Mouse ELMO2(Engulfment And Cell Motility 2) ELISA Kit	ELMO2	Q8BHL5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2959	Human Des(Desmin) ELISA Kit	Des	P17661	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2960	Mouse Des(Desmin) ELISA Kit	Des	P31001	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER2961	Rat Des(Desmin) ELISA Kit	Des	P48675	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2962	Human SSRP1(Structure Specific Recognition Protein 1) ELISA Kit	SSRP1	Q08945	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM2963	Mouse SSRP1(Structure Specific Recognition Protein 1) ELISA Kit	SSRP1	Q08943	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2964	Rat SSRP1(Structure Specific Recognition Protein 1) ELISA Kit	SSRP1		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH2965	Human HPR(Haptoglobin Related Protein) ELISA Kit	HPR	P00739	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL

FY-EH2966	Human TNNT1(Troponin T Type 1, Slow Skeletal) ELISA Kit	TNNT1	P13805	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2967	Rat TNNT1(Troponin T Type 1, Slow Skeletal) ELISA Kit	TNNT1	Q7TNB2	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2968	Human TNPO2(Transportin 2) ELISA Kit	TNPO2	O14787	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2969	Human UCP-2(Uncoupling Protein 2, Mitochondrial) ELISA Kit	UCP-2	P55851	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2970	Mouse UCP-2(Uncoupling Protein 2, Mitochondrial) ELISA Kit	UCP-2	P70406	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2971	Human HLTf(Helicase Like Transcription Factor) ELISA Kit	HLTF	Q14527	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM2972	Mouse HLTf(Helicase Like Transcription Factor) ELISA Kit	HLTF	Q6PCN7	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER2973	Rat HLTf(Helicase Like Transcription Factor) ELISA Kit	HLTF		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2974	Human MT2(Metallothionein 2) ELISA Kit	MT2	P02795	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER2975	Rat MT2(Metallothionein 2) ELISA Kit	MT2		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2976	Human TJP1(Tight junction protein ZO-1) ELISA Kit	TJP1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2977	Mouse TJP1(Tight junction protein ZO-1) ELISA Kit	TJP1	P39447	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER2978	Rat TJP1(Tight junction protein ZO-1) ELISA Kit	TJP1	A0A0G2K2F	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2979	Human CLC(Charcot Leyden Crystal Protein) ELISA Kit	CLC	Q05315	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2980	Human ACVB(Activin B) ELISA Kit	ACVB		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH2981	Human CRYαB(Crystallin Alpha B) ELISA Kit	CRYαB	P02511	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM2982	Mouse CRYαB(Crystallin Alpha B) ELISA Kit	CRYαB	P23927	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2983	Human Sp17(Sperm Protein 17) ELISA Kit	Sp17	Q15506	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2984	Human MSP(Macrophage Stimulating Protein) ELISA Kit	MSP	P26927	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2985	Human MCSFR(Macrophage Colony Stimulating Factor Receptor) ELISA Kit	MCSFR	P07333	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH2986	Human MIP-5(Macrophage Inflammatory Protein 5) ELISA Kit	MIP-5	Q16663	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2987	Human FZD8(Frizzled Homolog 8) ELISA Kit	FZD8	Q9H461	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2988	Human HAAO(3-Hydroxyanthranilate 3,4-Dioxygenase) ELISA Kit	HAAO		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2989	Human CRHBP(Corticotropin Releasing Hormone Binding Protein) ELISA Kit	CRHBP	P24387	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2990	Human αHSP(Alpha-Hemoglobin Stabilizing Protein) ELISA Kit	αHSP	Q9NZD4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2991	Human MMP-28(Matrix Metalloproteinase 28) ELISA Kit	MMP-28	Q9H239	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2992	Human SPTAN1(Alpha-Fodrin) ELISA Kit	SPTAN1	Q13813	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH2993	Human AZGP1(Zinc-alpha-2-glycoprotein) ELISA Kit	AZGP1	P25311	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM2994	Mouse AZGP1(Zinc-alpha-2-glycoprotein) ELISA Kit	AZGP1	Q64726	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH2995	Human MSLN(Mesothelin) ELISA Kit	MSLN	Q13421	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH2996	Human SOST(Sclerostin) ELISA Kit	SOST	Q9BQB4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM2997	Mouse SOST(Sclerostin) ELISA Kit	SOST	Q99P68	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER2998	Rat SOST(Sclerostin) ELISA Kit	SOST	Q99P67	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH2999	Human IRF5(Interferon Regulatory Factor 5) ELISA Kit	IRF5	Q13568	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3000	Mouse IRF5(Interferon Regulatory Factor 5) ELISA Kit	IRF5	P56477	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3001	Human SIRT1(Sirtuin 1) ELISA Kit	SIRT1	Q96EB6	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3002	Mouse SIRT1(Sirtuin 1) ELISA Kit	SIRT1	Q923E4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3003	Rat SIRT1(Sirtuin 1) ELISA Kit	SIRT1	A0A0G2JZ7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3004	Human VWFCP/ADAMTS13(Von Willebrand Factor Cleaving Protease) ELISA Kit	VWFCP/ADAM	Q76LX8	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3005	Human MTR(5-Methyltetrahydrofolate Homocysteine Methyltransferase) ELISA Kit	MTR	Q99707	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3006	Human PGF(Placental Growth Factor) ELISA Kit	PGF	P49763	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-EM3007	Mouse PGF(Placental Growth Factor) ELISA Kit	PGF	P49764	Sandwich	800pg/mL	7.5 pg/mL	12.5-800pg/mL
FY-ER3008	Rat PGF(Placental Growth Factor) ELISA Kit	PGF	Q63434	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3009	Human FOXC2(Forkhead Box Protein C2) ELISA Kit	FOXC2	Q99958	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3010	Human ESM1(Endothelial Cell Specific Molecule 1) ELISA Kit	ESM1	Q9NQ30	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3011	Mouse ESM1(Endothelial Cell Specific Molecule 1) ELISA Kit	ESM1	Q9QYY7	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3012	Rat ESM1(Endothelial Cell Specific Molecule 1) ELISA Kit	ESM1	P97682	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3013	Human Rac1(Ras Related C3 Botulinum Toxin Substrate 1) ELISA Kit	Rac1	P63000	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3014	Human LILRB2(Leukocyte Immunoglobulin Like Receptor Subfamily B Member 2) ELISA Kit	LILRB2	Q8N423	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3015	Human OTR(Oxytocin Receptor) ELISA Kit	OTR	P30559	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3016	Human PDIA6(Protein Disulfide Isomerase A6) ELISA Kit	PDIA6	Q15084	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3017	Human HDL(High Density Lipoprotein) ELISA Kit	HDL		Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-EM3018	Mouse HDL(High Density Lipoprotein) ELISA Kit	HDL		Sandwich	200nmol/mL	1.88 nmol/mL	3.13-200nmol/mL
FY-ER3019	Rat HDL(High Density Lipoprotein) ELISA Kit	HDL		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3020	Human Pro-MMP-10(Pro-Matrix Metalloproteinase 10) ELISA Kit	Pro-MMP-10	P09238	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3021	Human Pro-MMP-13(Pro-Matrix Metalloproteinase 13) ELISA Kit	Pro-MMP-13	P45452(MM	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3022	Human SIRPA(Signal Regulatory Protein Alpha) ELISA Kit	SIRPA	P78324	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3023	Human PDGF-AB(Platelet Derived Growth Factor AB) ELISA Kit	PDGF-AB		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3024	Mouse PDGF-AB(Platelet Derived Growth Factor AB) ELISA Kit	PDGF-AB		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3025	Human PGRN(Progranulin) ELISA Kit	PGRN	P28799(18~2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3026	Mouse PGRN(progranulin) ELISA Kit	PGRN	P28798	Sandwich	5ng/mL	0.05 ng/mL	0.08-5ng/mL
FY-ER3027	Rat PGRN(progranulin) ELISA Kit	PGRN		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3028	Human RLN2(Relaxin-2) ELISA Kit	RLN2	P04090	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3029	Human RBP1(Retinol Binding Protein 1, Cellular) ELISA Kit	RBP1	P09455	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3030	Human Surv(Survivin) ELISA Kit	Surv	O15392	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3031	Human RBP3(Retinol Binding Protein 3, Interstitial) ELISA Kit	RBP3	P10745	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3032	Rat TPO(Thyroid Peroxidase) ELISA Kit	TPO	P49745	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3033	Human TSP-1(Thrombospondin-1) ELISA Kit	TSP-1	P07996	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EM3034	Mouse TSP-1(Thrombospondin-1) ELISA Kit	TSP-1	P35441	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ER3035	Rat TSP-1(Thrombospondin-1) ELISA Kit	TSP-1	Q71SA3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3036	Human TSP-2(Thrombospondin-2) ELISA Kit	TSP-2	P35442	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3037	Human IL7R(Interleukin 7 Receptor) ELISA Kit	IL7R	P16871	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3038	Human TRAIL/TNFSF10(Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit	TRAIL/TNFSF10	P50591	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3039	Mouse TRAIL/TNFSF10(Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit	TRAIL/TNFSF10	P50592	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3040	Rat TRAIL/TNFSF10(Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit	TRAIL/TNFSF10		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3041	Human PLD1(Phospholipase D1) ELISA Kit	PLD1	Q13393	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3042	Human TREM-1(Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit	TREM-1	Q9NP99	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3043	Mouse TREM-1(Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit	TREM-1	Q9JKE2	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER3044	Rat TREM-1(Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit	TREM-1		Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3045	Human TSLP(Thymic Stromal Lymphopoietin) ELISA Kit	TSLP	Q969D9	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3046	Mouse TSLP(Thymic Stromal Lymphopoietin) ELISA Kit	TSLP	Q9JIE6	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3047	Rat TSLP(Thymic Stromal Lymphopoietin) ELISA Kit	TSLP		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-EH3048	Human PLAUR/uPAR(Plasminogen Activator, Urokinase Receptor) ELISA Kit	PLAUR/uPAR	Q03405	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3049	Mouse PLAUR/uPAR(Plasminogen Activator, Urokinase Receptor) ELISA Kit	PLAUR/uPAR	P35456	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER3050	Rat PLAUR/uPAR(Plasminogen Activator, Urokinase Receptor) ELISA Kit	PLAUR/uPAR	P49616	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3051	Human VEGF-C(Vascular Endothelial Growth Factor C) ELISA Kit	VEGF-C	P49767	Sandwich	7000pg/mL	65.63 pg/mL	109.38-7000pg/mL
FY-EM3052	Mouse VEGF-C(Vascular Endothelial Growth Factor C) ELISA Kit	VEGF-C	P97953	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3053	Rat VEGF-C(Vascular Endothelial Growth Factor C) ELISA Kit	VEGF-C	O35757	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3054	Human VEGF-D(Vascular Endothelial Growth Factor D) ELISA Kit	VEGF-D	O43915	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3055	Mouse VEGF-D(Vascular Endothelial Growth Factor D) ELISA Kit	VEGF-D	P97946	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3056	Rat VEGF-D(Vascular Endothelial Growth Factor D) ELISA Kit	VEGF-D	O35251	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3057	Human VEGFR-2/KDR(Vascular Endothelial Growth Factor Receptor 2) ELISA Kit	VEGFR-2/KDR		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM3058	Mouse VEGFR-2/KDR(Vascular Endothelial Growth Factor Receptor 2) ELISA Kit	VEGFR-2/KDR	P35918	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3059	Rat VEGFR-2/KDR(Vascular Endothelial Growth Factor Receptor 2) ELISA Kit	VEGFR-2/KDR	O08775	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-EH3060	Human CD34(Cluster Of Differentiation 34) ELISA Kit	CD34	P28906	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3061	Human CYB5A(Cytochrome b5, type A) ELISA Kit	CYB5A	P00167	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3062	Human sOX40L(soluble OX40L) ELISA Kit	sOX40L	P23510	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3063	Human BRAF(B-Raf Proto Oncogene Serine/Threonine Protein Kinase) ELISA Kit	BRAF	P15056	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3064	Human GzmA(Granzyme A) ELISA Kit	GzmA	P12544	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3065	Mouse GzmA(Granzyme A) ELISA Kit	GzmA	P11032	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3066	Rat GzmA(Granzyme A) ELISA Kit	GzmA		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3067	Human GzmB(Granzyme B) ELISA Kit	GzmB	P10144	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3068	Mouse GzmB(Granzyme B) ELISA Kit	GzmB	P04187	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3069	Rat GzmB(Granzyme B) ELISA Kit	GzmB	P18291	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3070	Human GNLY(Granulysin) ELISA Kit	GNLY	P22749	Sandwich	15ng/mL	0.14 ng/mL	0.24-15ng/mL
FY-EH3071	Human S100P(S100 Calcium Binding Protein P) ELISA Kit	S100P	P25815	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3072	Human CR2/CD21(Complement Receptor 2) ELISA Kit	CR2/CD21	P20023	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3073	Human cADPRH/CD38(Cyclic ADP Ribose Hydrolase) ELISA Kit	cADPRH/CD38	P28907	Sandwich	60ng/mL	0.56 ng/mL	0.94-60ng/mL
FY-ER3074	Rat cADPRH/CD38(Cyclic ADP Ribose Hydrolase) ELISA Kit	cADPRH/CD38	Q64244	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3075	Human CLDN5(Claudin 5) ELISA Kit	CLDN5	O00501	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3076	Rat CLDN5(Claudin 5) ELISA Kit	CLDN5	Q9JKD6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3077	Human GRIN2B(Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 2B) ELISA Kit	GRIN2B	Q13224	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3078	Human GRIN1(Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 1) ELISA Kit	GRIN1	Q05586	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3079	Human GRIN2A(Glutamate Receptor, Ionotropic, N-Methyl-D-Aspartate 2A) ELISA Kit	GRIN2A	Q12879	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3080	Human FXN(Frataxin) ELISA Kit	FXN	Q16595	Sandwich	1500pg/mL	14.06 pg/mL	23.44-1500pg/mL
FY-EH3081	Human LOX-1(Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit	LOX-1	P78380	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3082	Rat LOX-1(Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit	LOX-1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3083	Human ACTN4(Actinin Alpha 4) ELISA Kit	ACTN4	O43707	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3084	Human sAC(Soluble Adenylate Cyclase) ELISA Kit	sAC	Q96PN6	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3085	Human sFMC(Soluble Fibrin Monomer Complex) ELISA Kit	sFMC		Sandwich	100µg/mL	0.94 µg/mL	1.57-100µg/mL
FY-EH3086	Human PECAM1/CD31(Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit	PECAM1/CD31	P16284	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3087	Mouse PECAM1/CD31(Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit	PECAM1/CD31		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3088	Rat PECAM1/CD31(Platelet/Endothelial Cell Adhesion Molecule 1) ELISA Kit	PECAM1/CD31	Q3SWT0	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL

FY-EH3089	Human COL18α1(Collagen Type XVIII Alpha 1) ELISA Kit	COL18α1	P39060	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3090	Human MITF(Microphthalmia Associated Transcription Factor) ELISA Kit	MITF	O75030	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3091	Human BMP-5(Bone Morphogenetic Protein 5) ELISA Kit	BMP-5	P22003	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER3092	Rat BMP-5(Bone Morphogenetic Protein 5) ELISA Kit	BMP-5		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3093	Human IL-12 p35(Interleukin 12 p35) ELISA Kit	IL-12 p35		Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM3094	Mouse IL-12 p35(Interleukin 12 p35) ELISA Kit	IL-12 p35	P43431	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3095	Rat IL-12 p35(Interleukin 12 p35) ELISA Kit	IL-12 p35		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3096	Human IL-25(Interleukin 25) ELISA Kit	IL-25	Q9H293	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3097	Mouse IL-25(Interleukin 25) ELISA Kit	IL-25	Q8VHH8	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3098	Rat IL-25(Interleukin 25) ELISA Kit	IL-25	D3ZLB1	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3099	Human IL-34(Interleukin 34) ELISA Kit	IL-34	Q6ZMJ4	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EM3100	Mouse IL-34(Interleukin 34) ELISA Kit	IL-34	Q8R1R4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3101	Rat IL-34(Interleukin 34) ELISA Kit	IL-34	Q4KM46	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3102	Human RLN1(Relaxin-1) ELISA Kit	RLN1	P04808	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3103	Human RLN3(Relaxin-3) ELISA Kit	RLN3	Q8WXF3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3104	Human Tyk-2(Tyrosine Kinase 2) ELISA Kit	Tyk-2	P29597	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3105	Human TH(Tyrosine Hydroxylase) ELISA Kit	TH	P07101	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3106	Mouse TH(Tyrosine Hydroxylase) ELISA Kit	TH	P24529	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3107	Rat TH(Tyrosine Hydroxylase) ELISA Kit	TH		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3108	Human mTOR(Mammalian Target of Rapamycin) ELISA Kit	mTOR	P42345	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3109	Human PSPN(Persephin) ELISA Kit	PSPN	O60542	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EH3110	Human TGFβI(Transforming Growth Factor Beta Induced Protein) ELISA Kit	TGFβI	Q15582	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM3111	Mouse TGFβI(Transforming Growth Factor Beta Induced) ELISA Kit	TGFβI	P82198	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3112	Human TSC22(Transforming Growth Factor Beta Stimulated Protein Clone 22) ELIS	TSC22	Q15714	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3113	Rat TSC22(Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA K	TSC22		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3114	Human HLA-A(Leukocyte Antigen A) ELISA Kit	HLA-A	P04439	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3115	Human PKCη(Protein Kinase C Eta) ELISA Kit	PKCη	P24723	Sandwich	200mIU/mL	1.88 mIU/mL	3.13-200mIU/mL
FY-EH3116	Human UCP-1(Uncoupling Protein 1, Mitochondrial) ELISA Kit	UCP-1	P25874	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3117	Mouse UCP-1(Uncoupling Protein 1, Mitochondrial) ELISA Kit	UCP-1	P12242	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3118	Rat UCP-1(Uncoupling Protein 1, Mitochondrial) ELISA Kit	UCP-1	P04633	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3119	Human IL18BP(Interleukin 18 Binding Protein) ELISA Kit	IL18BP	O95998	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3120	Mouse IL18BP(Interleukin 18 Binding Protein) ELISA Kit	IL18BP	Q9Z0M9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3121	Rat IL18BP(Interleukin 18 Binding Protein) ELISA Kit	IL18BP	F7F5P4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3122	Human TYR(Tyrosinase) ELISA Kit	TYR	P14679	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM3123	Mouse TYR(Tyrosinase) ELISA Kit	TYR	P11344	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3124	Rat TYR(Tyrosinase) ELISA Kit	TYR		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3125	Human NFκB2(Nuclear Factor Kappa B2) ELISA Kit	NFκB2	Q00653	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3126	Mouse NFκB2(Nuclear Factor Kappa B2) ELISA Kit	NFκB2	Q9WTK5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3127	Rat NFκB2(Nuclear Factor Kappa B2) ELISA Kit	NFκB2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3128	Human FGL1(Fibrinogen Like Protein 1) ELISA Kit	FGL1	Q08830	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3129	Human FBLN2(Fibulin 2) ELISA Kit	FBLN2	P98095	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EH3130	Human FBLN3(Fibulin 3) ELISA Kit	FBLN3	Q12805	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3131	Human GMCSF-Ab(Granulocyte Macrophage Colony Stimulating Factor Antibody) ELISA Kit	GMCSF-Ab		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3132	Human PRG4(Proteoglycan 4) ELISA Kit	PRG4	Q92954	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH3133	Human FBLN1(Fibulin 1) ELISA Kit	FBLN1	P23142	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3134	Human CX40(Connexin 40) ELISA Kit	CX40	P36382	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3135	Mouse CX40(Connexin 40) ELISA Kit	CX40	Q01231	Sandwich	1500pg/mL	14.06 pg/mL	23.44-1500pg/mL
FY-ER3136	Rat CX40(Connexin 40) ELISA Kit	CX40		Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-EH3137	Human CX43(Connexin 43) ELISA Kit	CX43	P17302	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3138	Mouse CX43(Connexin 43) ELISA Kit	CX43	P23242	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3139	Rat CX43(Connexin 43) ELISA Kit	CX43	P08050	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3140	Human SDC4(Syndecan 4) ELISA Kit	SDC4	P31431	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM3141	Mouse SDC4(Syndecan 4) ELISA Kit	SDC4	O35988	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3142	Rat SDC4(Syndecan 4) ELISA Kit	SDC4	P34901	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3143	Human LAP(Leucine Aminopeptidase) ELISA Kit	LAP	P28838	Sandwich	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-EH3144	Human LTβ(Lymphotoxin Beta) ELISA Kit	LTβ	Q06643	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3145	Human LFA-1(Lymphocyte Function Associated Antigen 1) ELISA Kit	LFA-1	P20701	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3146	Human LFA-2/CD2(Lymphocyte Function Associated Antigen 2) ELISA Kit	LFA-2/CD2	P06729	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3147	Human Lptn/XCL1(Lymphotoctin) ELISA Kit	Lptn/XCL1	P47992	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3148	Mouse Lptn/XCL1(Lymphotoctin) ELISA Kit	Lptn/XCL1	P47993	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER3149	Rat Lptn/XCL1(Lymphotoctin) ELISA Kit	Lptn/XCL1	P51672	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3150	Human PGAM2(Phosphoglycerate Mutase 2) ELISA Kit	PGAM2	P15259	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3151	Human PFKM(Phosphofructokinase, Muscle) ELISA Kit	PFKM	P08237	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3152	Human pIκBα(Phospho-Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit	pIκBα	P25963	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3153	Mouse pIκBα(Phospho-Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit	pIκBα		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3154	Rat pIκBα(Phospho-Inhibitory Subunit of NF-KBα) ELISA Kit	pIκBα		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3155	Human PGM(Phosphoglucomutase) ELISA Kit	PGM	P36871	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3156	Human PCK1(Phosphoenolpyruvate Carboxykinase 1) ELISA Kit	PCK1	P35558	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3157	Mouse PCK1(Phosphoenolpyruvate Carboxykinase 1) ELISA Kit	PCK1	Q9Z2V4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3158	Rat PCK1(Phosphoenolpyruvate Carboxykinase 1) ELISA Kit	PCK1	P07379	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3159	Human Integrin beta-1(ITGB1) ELISA kit	Integrin beta-1	P05556	Sandwich	4000pg/mL	15.6 pg/mL	62.5-4000 pg/mL
FY-EH3160	Human PLSCR1(Phospholipid Scramblase 1) ELISA Kit	PLSCR1	O15162	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3161	Human PLSCR2(Phospholipid Scramblase 2) ELISA Kit	PLSCR2	Q9NRY7	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3162	Human PI3K2α(Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit	PI3K2α	O00443	Sandwich	200mIU/mL	1.88 mIU/mL	3.13-200mIU/mL
FY-EM3163	Mouse PI3K2α(Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit	PI3K2α	Q61194	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER3164	Rat PI3K2α(Phosphoinositide 3 Kinase 2 Alpha) ELISA Kit	PI3K2α		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3165	Human GPC1(Glypican 1) ELISA Kit	GPC1	P35052	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3166	Human GPC2(Glypican 2) ELISA Kit	GPC2	Q8N158	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3167	Human GPC3(Glypican 3) ELISA Kit	GPC3	P51654	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3168	Mouse GPC3(Glypican 3) ELISA Kit	GPC3	Q8CFZ4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3169	Rat GPC3(Glypican 3) ELISA Kit	GPC3	P13265	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3170	Human GPC4(Glypican 4) ELISA Kit	GPC4	O75487	Sandwich	10ng/mL	0.10 ng/mL	0.16-10ng/mL

FY-EM3171	Mouse GPC4(Glypican 4) ELISA Kit	GPC4	P51655	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3172	Human TRUB1(Probable tRNA Pseudouridine Synthase 1) ELISA Kit	TRUB1	Q8WWH5	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3173	Human PLTP(Phospholipid Transfer Protein) ELISA Kit	PLTP	P55058	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3174	Human PLCγ1(Phospholipase C Gamma 1) ELISA Kit	PLCγ1	P19174	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3175	Human FUR(Furin) ELISA Kit	FUR	P09958	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3176	Human SCCA1(Squamous Cell Carcinoma Antigen 1) ELISA Kit	SCCA1	P29508	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3177	Human SCCA2(Squamous Cell Carcinoma Antigen 2) ELISA Kit	SCCA2	P48594	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3178	Mouse SCCA2(Squamous Cell Carcinoma Antigen 2) ELISA Kit	SCCA2	Q9Z2G2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3179	Human pPLA2(Phospholipase A2, Pancreas) ELISA Kit	pPLA2	P04054	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3180	Human DS(Dermatan Sulfate) ELISA Kit	DS		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3181	Human Trx(Thioredoxin) ELISA Kit	Trx	P10599	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3182	Mouse Trx(Thioredoxin) ELISA Kit	Trx	P10639	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER3183	Rat TRX(Thioredoxin) ELISA Kit	TRX	P11232	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3184	Human TrxR1(Thioredoxin Reductase 1) ELISA Kit	TrxR1	Q16881	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3185	Mouse TrxR1(Thioredoxin Reductase 1) ELISA Kit	TrxR1	Q9JMH6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3186	Rat TrxR1(Thioredoxin Reductase 1) ELISA Kit	TrxR1	O89049	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3187	Human TBP2(Thioredoxin Binding Protein 2) ELISA Kit	TBP2	Q9H3M7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3188	Mouse TBP2(Thioredoxin Binding Protein 2) ELISA Kit	TBP2	Q8BG60	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3189	Human IAP(Integrin Associated Protein) ELISA Kit	Human IAP	Q08722	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3190	Human MCL1(Myeloid Cell Leukemia Sequence 1) ELISA Kit	Human MCL1	Q07820	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3191	Human MPA(Mycophenolic Acid) ELISA Kit	MPA		Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3192	Human FcαR/CD89(Receptor I for the Fc Fragment of IgA) ELISA Kit	FcαR/CD89	P24071	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3193	Human FcεR I (Receptor I for the Fc Fragment of IgE) ELISA Kit	FcεR I	P12319(FCE	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3194	Human FcγR II/CD32(Receptor II for the Fc Fragment of IgG) ELISA Kit	FcγR II/CD32	P31994	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-EH3195	Human FcγR3B(Fc Fragment of IgG Low Affinity IIIb Receptor) ELISA Kit	FcγR3B	O75015	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3196	Human Fcγ(Fc Fragment of IgG) ELISA Kit	Fcγ	P48200, 错	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3197	Mouse Fcγ(Fc Fragment of IgG) ELISA Kit	Fcγ		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH3198	Human MCP/CD46(Membrane Cofactor Protein) ELISA Kit	MCP/CD46	P15529	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3199	Human TSP-4(Thrombospondin-4) ELISA Kit	TSP-4	P35443	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3200	Human TDT(Terminal Deoxynucleotidyl Transferase) ELISA Kit	TDT	P04053	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3201	Human SMN1(Survival of Motor Neuron 1, Telomeric) ELISA Kit	SMN1	Q16637	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3202	Mouse SMN1(Survival Of Motor Neuron 1, Telomeric) ELISA Kit	SMN1	P97801	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER3203	Rat SMN1(Survival of Motor Neuron 1, Telomeric) ELISA Kit	SMN1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3204	Human ETRA(Endothelin Receptor A) ELISA Kit	ETRA	P25101	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3205	Mouse ETRA(Endothelin Receptor A) ELISA Kit	ETRA	Q61614	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3206	Human ECE1(Endothelin Converting Enzyme 1) ELISA Kit	ECE1	P42892	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3207	Mouse ECE1(Endothelin Converting Enzyme 1) ELISA Kit	ECE1	Q4PZA2	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3208	Human REPIN1(Replication Initiator 1) ELISA Kit	REPIN1	Q9BWE0	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3209	Human EL/LIPG(Endothelial Lipase) ELISA Kit	EL/LIPG	Q9Y5X9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3210	Human VASPIN(Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit	VASPIN	Q8IW75	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3211	Mouse VASPIN(Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit	VASPIN	Q7TMF5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-ER3212	Rat VASPIN(Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit	VASPIN	Q8R4Z1	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EH3213	Human VF(Visfatin) ELISA Kit	VF	P43490	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3214	Mouse VF(Visfatin) ELISA Kit	VF	Q99KQ4	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER3215	Rat VF(Visfatin) ELISA Kit	VF	Q80Z29	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3216	Human ERN1(Endoplasmic Reticulum to Nucleus Signalling 1) ELISA Kit	ERN1	O75460	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3217	Mouse ERN1(Endoplasmic Reticulum To Nucleus Signalling 1) ELISA Kit	ERN1	Q9EQY0	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ER3218	Rat ERN1(Endoplasmic Reticulum To Nucleus Signalling 1) ELISA Kit	ERN1		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3219	Human CFTR(Cystic Fibrosis Transmembrane Conductance Regulator) ELISA Kit	CFTR	P13569	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3220	Human NGB(Neuroglobin) ELISA Kit	NGB	Q9NPG2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3221	Mouse NGB(Neuroglobin) ELISA Kit	NGB	Q9ER97	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER3222	Rat NGB(Neuroglobin) ELISA Kit	NGB	Q99JA8	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3223	Human MAdCAM1(Mucosal Addressin Cell Adhesion Molecule 1) ELISA Kit	MAdCAM1	Q13477	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3224	Human FAK(Focal Adhesion Kinase) ELISA Kit	FAK	Q05397	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3225	Mouse FAK(Focal Adhesion Kinase) ELISA Kit	FAK	P34152	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER3226	Rat FAK(Focal Adhesion Kinase) ELISA Kit	FAK	O35346	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3227	Human OTC(Ornithine Transcarbamylase) ELISA Kit	Human OTC		Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-EM3228	Mouse OTC(Ornithine Transcarbamylase) ELISA Kit	OTC		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3229	Human ODC(Ornithine Decarboxylase) ELISA Kit	ODC	P11926	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3230	Mouse ODC(Ornithine Decarboxylase) ELISA Kit	ODC	P00860	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3231	Rat ODC(Ornithine Decarboxylase) ELISA Kit	ODC	P09057	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3232	Human GBP4(Guanylate Binding Protein 4) ELISA Kit	GBP4	Q96PP9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3233	Human UGCG(UDP-glucose ceramide glucosyltransferase) ELISA Kit	UGCG	Q16739	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3234	Human PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit	PLAU/uPA	P00749	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3235	Mouse PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit	PLAU/uPA	P06869	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3236	Rat PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit	PLAU/uPA		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3237	Human UPP1(Uridine Phosphorylase 1) ELISA Kit	UPP1	Q16831	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3238	Human ETRB(Endothelin Receptor B) ELISA Kit	ETRB	P24530	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3239	Mouse GS(Glutamine synthetase) ELISA Kit	GS	P15105	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3240	Rat GS(Glutamine synthetase) ELISA Kit	GS	P09606	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3241	Human TAT(Thrombin/Antithrombin Complex) ELISA Kit	TAT		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3242	Mouse TAT(Thrombin/Antithrombin Complex) ELISA Kit	TAT		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3243	Rat TAT(Thrombin/Antithrombin Complex) ELISA Kit	TAT		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3244	Human PT(Prothrombin) ELISA Kit	PT	P00734	Sandwich	600ng/mL	5.63 ng/mL	9.38-600ng/mL
FY-EH3245	Human F1+2(Prothrombin Fragment 1+2) ELISA Kit	F1+2		Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-ER3246	Rat F1+2(Prothrombin Fragment 1+2) ELISA Kit	F1+2		Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-EH3247	Human KERA(Keratocan) ELISA Kit	KERA	O60938	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3248	Human VCL(Vinculin) ELISA Kit	VCL	P18206	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3249	Human EED(Embryonic Ectoderm Development) ELISA Kit	EED	O75530	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3250	Mouse EED(Embryonic Ectoderm Development) ELISA Kit	EED	Q921E6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3251	Rat EED(Embryonic Ectoderm Development) ELISA Kit	EED		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3252	Human DPT(Dermatopontin) ELISA Kit	DPT	Q07507	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH3253	Human SYK(Spleen Tyrosine Kinase) ELISA Kit	SYK	P43405	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3254	Human PDPN(Podoplanin) ELISA Kit	PDPN	Q86YL7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3255	Mouse PDPN(Podoplanin) ELISA Kit	PDPN	Q62011	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3256	Rat PDPN(Podoplanin) ELISA Kit	PDPN	Q64294	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3257	Human OSF(Osteoclast Stimulating Factor 1) ELISA Kit	OSF	Q92882	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3258	Rat OSF(Osteoclast Stimulating Factor 1) ELISA Kit	OSF		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3259	Human OSCAR(Osteoclast Associated Receptor) ELISA Kit	OSCAR	Q8IYS5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3260	Mouse OSCAR(Osteoclast Associated Receptor) ELISA Kit	OSCAR	Q8VBT3	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3261	Rat OSCAR(Osteoclast Associated Receptor) ELISA Kit	OSCAR		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3262	Human G6PD(Glucose 6 Phosphate Dehydrogenase) ELISA Kit	G6PD	P11413	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3263	Mouse G6PD(Glucose 6 Phosphate Dehydrogenase) ELISA Kit	G6PD	Q8CFX1	Sandwich	100µg/mL	0.94 µg/mL	1.57-100µg/mL
FY-ER3264	Rat G6PD(Glucose 6 Phosphate Dehydrogenase) ELISA Kit	G6PD	P05370	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3265	Human GPI(Glucose 6 Phosphate Isomerase) ELISA Kit	GPI	P06744	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3266	Mouse GPI(Glucose 6 Phosphate Isomerase) ELISA Kit	GPI	P06745	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3267	Human GCK(Glucokinase) ELISA Kit	GCK	P35557	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM3268	Mouse GCK(Glucokinase) ELISA Kit	GCK	P52792	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3269	Rat GCK(Glucokinase) ELISA Kit	GCK	F7EZX6	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3270	Human GKRP(Glucokinase Regulatory Protein) ELISA Kit	GKRP	Q14397	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3271	Human GLUT1(Glucose Transporter 1) ELISA Kit	GLUT1	P11166	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3272	Mouse GLUT1(Glucose Transporter 1) ELISA Kit	GLUT1	P17809	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3273	Rat GLUT1(Glucose Transporter 1) ELISA Kit	GLUT1	P11167	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH3274	Human GLUT2(Glucose Transporter 2) ELISA Kit	GLUT2	P11168	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM3275	Mouse GLUT2(Glucose Transporter 2) ELISA Kit	GLUT2	P14246	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3276	Rat GLUT2(Glucose Transporter 2) ELISA Kit	GLUT2	P12336	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3277	Human GLUT3(Glucose Transporter 3) ELISA Kit	GLUT3	P11169	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3278	Mouse GLUT3(Glucose Transporter 3) ELISA Kit	GLUT3	P32037	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER3279	Rat GLUT3(Glucose Transporter 3) ELISA Kit	GLUT3	Q07647	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3280	Human GLUT4(Glucose Transporter 4) ELISA Kit	GLUT4	P14672	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3281	Mouse GLUT4(Glucose Transporter 4) ELISA Kit	GLUT4	P14142	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3282	Rat GLUT4(Glucose Transporter 4) ELISA Kit	GLUT4	P19357	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3283	Human GTF2B(General Transcription Factor II B) ELISA Kit	GTF2B	Q00403	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3284	Human GTF2F1(General Transcription Factor II F, Polypeptide 1) ELISA Kit	GTF2F1	P35269	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3285	Human PK2(Prokineticin 2) ELISA Kit	PK2	Q9HC23	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3286	Mouse PK2(Prokineticin 2) ELISA Kit	PK2	Q9QXU7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3287	Rat PK2(Prokineticin 2) ELISA Kit	PK2	Q8R413	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3288	Human PKR1(Prokineticin Receptor 1) ELISA Kit	PKR1	Q8TCW9	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM3289	Mouse PKR1(Prokineticin Receptor 1) ELISA Kit	PKR1	Q9JKL1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3290	Rat PKR1(Prokineticin Receptor 1) ELISA Kit	PKR1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3291	Human PKR2(Prokineticin Receptor 2) ELISA Kit	PKR2	Q8NFI6	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3292	Human EP3(Prostaglandin E Receptor 3) ELISA Kit	EP3	P43115	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3293	Human PGF2α(Prostaglandin F2α) ELISA Kit	PGF2α		Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL

FY-EM3294	Mouse PGF2 α (Prostaglandin F2 Alpha) ELISA Kit	PGF2 α		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3295	Rat PGF2 α (Prostaglandin F2 Alpha) ELISA Kit	PGF2 α		Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3296	Human FP(Prostaglandin F Receptor) ELISA Kit	FP	P43088	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER3297	Rat FP(Prostaglandin F Receptor) ELISA Kit	FP		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3298	Human PTGIS(Prostaglandin I Synthase) ELISA Kit	PTGIS	Q16647	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3299	Human PTGES2(Prostaglandin E Synthase 2) ELISA Kit	PTGES2	Q9H7Z7	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM3300	Mouse PTGES2(Prostaglandin E Synthase 2) ELISA Kit	PTGES2	Q8BWM0	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3301	Rat PTGES2(Prostaglandin E Synthase 2) ELISA Kit	PTGES2		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3302	Human PTGS1(Prostaglandin Endoperoxide Synthase 1) ELISA Kit	PTGS1	P23219	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM3303	Mouse PTGS1(Prostaglandin Endoperoxide Synthase 1) ELISA Kit	PTGS1	P22437	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER3304	Rat PTGS1(Prostaglandin Endoperoxide Synthase 1) ELISA Kit	PTGS1	Q63921	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3305	Human PTGS2/COX-2(Prostaglandin Endoperoxide Synthase 2) ELISA Kit	PTGS2/COX-2	P35354	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3306	Mouse PTGS2/COX-2(Prostaglandin Endoperoxide Synthase 2) ELISA Kit	PTGS2/COX-2	Q05769	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3307	Rat PTGS2/COX-2(Prostaglandin Endoperoxide Synthase 2) ELISA Kit	PTGS2/COX-2	P35355	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH3308	Human Dyn(Big Dynorphin) ELISA Kit	Dyn	P01213	Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3309	Human DSP(Desmoplakin) ELISA Kit	DSP	P15924	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3310	Human ERCC1(Excision Repair Cross Complementing Rodent Repair Deficiency Cc	ERCC1	P07992	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3311	Human KLK1(Kallikrein 1) ELISA Kit	KLK1	P06870	Sandwich	1200pg/mL	11.25 pg/mL	18.75-1200pg/mL
FY-EH3312	Human LXR α (Liver X Receptor Alpha) ELISA Kit	LXR α	Q13133	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3313	Rat LXR α (Liver X Receptor Alpha) ELISA Kit	LXR α		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3314	Human HSP β 6/HSP-20(Heat Shock Protein Beta 6) ELISA Kit	HSP β 6/HSP-20	O14558	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3315	Human HSP-27/HSPB1(Heat Shock Protein 27) ELISA Kit	HSP-27/HSPB1	P04792	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM3316	Mouse HSP-27/HSPB1(Heat Shock Protein 27) ELISA Kit	HSP-27/HSPB1	P14602	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER3317	Rat HSP-27/HSPB1(Heat Shock Protein 27) ELISA Kit	HSP-27/HSPB1		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3318	Human HSP-40(Heat Shock Protein 40) ELISA Kit	HSP-40	P25686	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3319	Mouse HSP-40(Heat Shock Protein 40) ELISA Kit	HSP-40	Q9QYJ3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER3320	Rat HSP-40(Heat Shock Protein 40) ELISA Kit	HSP-40		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3321	Human HSP-60/HSPD1(Heat Shock Protein 60) ELISA Kit	HSP-60/HSPD1	P10809	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3322	Human HSP-90(Heat Shock Protein 90) ELISA Kit	HSP-90	P08238	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM3323	Mouse HSP-90(Heat Shock Protein 90) ELISA Kit	HSP-90	P07901(Hsp'	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3324	Rat HSP-90(Heat Shock Protein 90) ELISA Kit	HSP-90		Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH3325	Human HSP gp96(Heat Shock Protein glycoprotein 96) ELISA Kit	HSP gp96		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3326	Mouse HSP gp96(Heat Shock Protein glycoprotein 96) ELISA Kit	HSP gp96	P08113	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER3327	Rat HSP GP96(Heat Shock Protein Glycoprotein 96) ELISA Kit	HSP GP96		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3328	Human HSF2(Heat Shock Transcription Factor 2) ELISA Kit	HSF2	Q03933	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3329	Mouse HSF2(Heat Shock Transcription Factor 2) ELISA Kit	HSF2	P38533	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER3330	Rat HSF2(Heat Shock Transcription Factor 2) ELISA Kit	HSF2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3331	Human PSG1/SP1(Pregnancy Specific Beta-1-Glycoprotein 1) ELISA Kit	PSG1/SP1	P11464	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3332	Human LZM(Lysozyme) ELISA Kit	LZM	P61626	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EM3333	Mouse LZM(Lysozyme) ELISA Kit	LZM	P17897(LYZ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH3334	Human TERF2(Telomeric Repeat Binding Factor 2) ELISA Kit	TERF2	Q15554	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH3335	Human TFF1(Trefoil Factor 1) ELISA Kit	TFF1	P04155	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH3336	Human TFF2(Trefoil Factor 2) ELISA Kit	TFF2	Q03403	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3337	Mouse TFF2(Trefoil Factor 2) ELISA Kit	TFF2	Q03404	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER3338	Rat TFF2(Trefoil Factor 2) ELISA Kit	TFF2		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3339	Human PEDF(Pigment Epithelium Derived Factor) ELISA Kit	PEDF	P36955	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3340	Mouse PEDF(Pigment Epithelium Derived Factor) ELISA Kit	PEDF	P97298	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3341	Rat PEDF(Pigment Epithelium Derived Factor) ELISA Kit	PEDF		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3342	Human SDH(Sorbitol Dehydrogenase) ELISA Kit	SDH	Q00796	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM3343	Mouse SDH(Sorbitol Dehydrogenase) ELISA Kit	SDH		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER3344	Rat SDH(Sorbitol Dehydrogenase) ELISA Kit	SDH	P27867	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3345	Human CCR2(Chemokine C-C-Motif Receptor 2) ELISA Kit	CCR2	P41597	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3346	Human Ep-CAM/CD326(Epithelial Cell Adhesion Molecule) ELISA Kit	Ep-CAM/CD326	P16422	Sandwich	100IU/mL	0.94 IU/mL	1.57-100IU/mL
FY-EH3347	Human NKA(Neurokinin A) ELISA Kit	NKA	P20366	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3348	Human NKB(Neurokinin B) ELISA Kit	NKB	Q9UHF0(81	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3349	Human GMFB(Glia Maturation Factor, Beta) ELISA Kit	GMFB	P60983	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3350	Mouse GMFB(Glia Maturation Factor, Beta) ELISA Kit	GMFB	Q9CQI3	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3351	Rat GMFB(Glia Maturation Factor, Beta) ELISA Kit	GMFB		Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH3352	Human CL(Cardiolipin) ELISA Kit	CL		Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3353	Human HK1(Hexokinase 1) ELISA Kit	HK1	P19367	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3354	Mouse HK1(Hexokinase 1) ELISA Kit	HK1	P17710	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3355	Rat HK1(Hexokinase 1) ELISA Kit	HK1	P05708	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3356	Human P2(Myelin Protein 2) ELISA Kit	P2	P02689	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3357	Human NPY(Neuropeptide Y) ELISA Kit	NPY	P01303	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3358	Mouse NPY(Neuropeptide Y) ELISA Kit	NPY	P57774	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3359	Rat NPY(Neuropeptide Y) ELISA Kit	NPY	P07808	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3360	Human NCAM/CD56(Neural Cell Adhesion Molecule) ELISA Kit	NCAM/CD56	P13591	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM3361	Mouse NCAM/CD56(Neural Cell Adhesion Molecule) ELISA Kit	NCAM/CD56	P13595	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3362	Rat NCAM/CD56(Neural Cell Adhesion Molecule) ELISA Kit	NCAM/CD56		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3363	Human L1CAM/CD171(L1-Cell Adhesion Molecule) ELISA Kit	L1CAM/CD171	P32004	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EH3364	Human NTRK2(Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit	NTRK2	Q16620	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EM3365	Mouse NTRK2(Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit	NTRK2	P15209	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER3366	Rat NTRK2(Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit	NTRK2	Q63604	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EH3367	Human NTRK3(Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit	NTRK3	Q16288	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM3368	Mouse NTRK3(Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit	NTRK3	Q6VNS1	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER3369	Rat NTRK3(Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit	NTRK3		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3370	Human NAIP(Neuronal Apoptosis Inhibitory Protein) ELISA Kit	NAIP	Q13075	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3371	Human NPHN(Nephrin) ELISA Kit	NPHN	O60500	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3372	Mouse NPHN(Nephrin) ELISA Kit	NPHN	Q9QZS7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3373	Rat NPHN(Nephrin) ELISA Kit	NPHN	Q9R044	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3374	Human GTEAGE(Glomerular Tissue Extract Advanced Glycosylation End Products) ELISA Kit	GTEAGE		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3375	Human GROγ/CXCL3(Growth Regulated Oncogene Gamma) ELISA Kit	GROγ/CXCL3	P19876	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL

FY-EM3376	Mouse GRO γ /CXCL3(Growth Regulated Oncogene Gamma) ELISA Kit	GRO γ /CXCL3	Q6W5C0	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER3377	Rat GRO γ /CXCL3(Growth Regulated Oncogene Gamma) ELISA Kit	GRO γ /CXCL3	Q10746	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3378	Human GDF1(Growth Differentiation Factor 1) ELISA Kit	GDF1	P27539	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3379	Mouse GDF1(Growth Differentiation Factor 1) ELISA Kit	GDF1	P20863	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3380	Rat GDF1(Growth Differentiation Factor 1) ELISA Kit	GDF1	Q1HI69	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3381	Human GDF10(Growth Differentiation Factor 10) ELISA Kit	GDF10	P55107	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3382	Mouse GDF10(Growth Differentiation Factor 10) ELISA Kit	GDF10	P97737	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3383	Human GDF11(Growth Differentiation Factor 11) ELISA Kit	GDF11	O95390	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3384	Mouse GDF11(Growth Differentiation Factor 11) ELISA Kit	GDF11	Q9Z1W4	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3385	Rat GDF11(Growth Differentiation Factor 11) ELISA Kit	GDF11	Q9Z217	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3386	Human GDF3(Growth Differentiation Factor 3) ELISA Kit	GDF3	Q9NR23	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM3387	Mouse GDF3(Growth Differentiation Factor 3) ELISA Kit	GDF3	Q07104	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER3388	Rat GDF3(Growth Differentiation Factor 3) ELISA Kit	GDF3	A2SY89	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3389	Human GDF6(Growth Differentiation Factor 6) ELISA Kit	GDF6	Q6KF10	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3390	Mouse GDF6(Growth Differentiation Factor 6) ELISA Kit	GDF6	P43028	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3391	Rat GDF6(Growth Differentiation Factor 6) ELISA Kit	GDF6		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3392	Human GDF7(Growth Differentiation Factor 7) ELISA Kit	GDF7	Q7Z4P5	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3393	Mouse GDF7(Growth Differentiation Factor 7) ELISA Kit	GDF7	P43029	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3394	Human GDF9(Growth Differentiation Factor 9) ELISA Kit	GDF9	O60383	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3395	Mouse GDF9(Growth Differentiation Factor 9) ELISA Kit	GDF9	Q07105	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3396	Rat GDF9(Growth Differentiation Factor 9) ELISA Kit	GDF9	Q9Z0X3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3397	Human HK2(Hexokinase 2) ELISA Kit	HK2	P52789	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH3398	Human GH2(Growth Hormone 2) ELISA Kit	GH2	P01242	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3399	Human OSBPL1A(Oxysterol Binding Protein Like Protein 1A) ELISA Kit	OSBPL1A	Q9BXW6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3400	Human GHRL(Ghrelin) ELISA Kit	GHRL	Q9UBU3	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3401	Mouse GHRL(Ghrelin) ELISA Kit	GHRL	Q9EQX0(24	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3402	Rat GHRL(Ghrelin) ELISA Kit	GHRL		Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3403	Human GHITM(Growth Hormone Inducible Transmembrane Protein) ELISA Kit	GHITM	Q9H3K2	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3404	Human GAP43(Growth Associated Protein 43) ELISA Kit	GAP43	P17677	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3405	Rat GAP43(Growth Associated Protein 43) ELISA Kit	GAP43		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3406	Human SST(Somatostatin) ELISA Kit	SST	P61278	Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM3407	Mouse SST(Somatostatin) ELISA Kit	SST	P60041	Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER3408	Rat SST(Somatostatin) ELISA Kit	SST	P60042	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3409	Human SSTR1(Somatostatin Receptor 1) ELISA Kit	SSTR1	P30872	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3410	Mouse SSTR1(Somatostatin Receptor 1) ELISA Kit	SSTR1	P30873	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3411	Human Grb10(Growth Factor Receptor Bound Protein 10) ELISA Kit	Grb10	Q13322	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3412	Mouse Grb10(Growth Factor Receptor Bound Protein 10) ELISA Kit	Grb10	Q60760	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3413	Rat GRB10(Growth Factor Receptor Bound Protein 10) ELISA Kit	GRB10		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3414	Human Grb14(Growth Factor Receptor Bound Protein 14) ELISA Kit	Grb14	Q14449	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3415	Mouse Grb14(Growth Factor Receptor Bound Protein 14) ELISA Kit	Grb14	Q9JLM9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3416	Rat GRB14(Growth Factor Receptor Bound Protein 14) ELISA Kit	GRB14		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH3417	Human Grb2(Growth Factor Receptor Bound Protein 2) ELISA Kit	Grb2	P62993	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3418	Mouse Grb2(Growth Factor Receptor Bound Protein 2) ELISA Kit	Grb2	Q60631	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3419	Rat GRB2(Growth Factor Receptor Bound Protein 2) ELISA Kit	GRB2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3420	Human Grb7(Growth Factor Receptor Bound Protein 7) ELISA Kit	Grb7	Q14451	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3421	Mouse Grb7(Growth Factor Receptor Bound Protein 7) ELISA Kit	Grb7	Q03160	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3422	Human OXB(Orexin B) ELISA Kit	OXB	O43612(70~	Competitive	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-EM3423	Mouse OXB(Orexin B) ELISA Kit	OXB	O55241	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3424	Rat OXB(Orexin B) ELISA Kit	OXB		Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3425	Human RB1(Retinoblastoma Protein 1) ELISA Kit	RB1	P06400	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3426	Human CYPB(Cyclophilin B) ELISA Kit	CYPB	P23284	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3427	Mouse CYPB(Cyclophilin B) ELISA Kit	CYPB	P24369	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-ER3428	Rat CYPB(Cyclophilin B) ELISA Kit	CYPB	P24368	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3429	Human CYPD(Cyclophilin D) ELISA Kit	CYPD	P30405	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3430	Mouse CYPD(Cyclophilin D) ELISA Kit	CYPD	Q99KR7	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3431	Rat CYPD(Cyclophilin D) ELISA Kit	CYPD	Q6DGG0	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3432	Human NRP2(Neuropilin 2) ELISA Kit	NRP2	O60462	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3433	Rat NRP2(Neuropilin 2) ELISA Kit	NRP2	O35276	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3434	Human CS(Citrate Synthase) ELISA Kit	CS	O75390	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-ER3435	Rat CS(Citrate Synthase) ELISA Kit	CS	Q8VHF5	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3436	Human OGG1(8-Oxoguanine Glycosylase 1) ELISA Kit	OGG1	O15527	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3437	Human ECFR(Eosinophil Chemotactic Factor Receptor) ELISA Kit	ECFR	P51677	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3438	Mouse ECFR(Eosinophil chemotactic factor receptor) ELISA Kit	ECFR	P51678	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3439	Human EAR1(Eosinophil Associated Ribonuclease A family Member 1) ELISA Kit	EAR1	Q53ZF0	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3440	Human EAR2(Eosinophil Associated Ribonuclease A family Member 2) ELISA Kit	EAR2	Q5WN11	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3441	Human EAR3(Eosinophil Associated Ribonuclease A family Member 3) ELISA Kit	EAR3	P10589	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3442	Human NOG(Noggin) ELISA Kit	NOG	Q13253	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH3443	Human PLN(Phospholamban) ELISA Kit	PLN	P26678	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3444	Human RIPK2(Receptor TNFRSF Interacting Serine Threonine Kinase 2) ELISA Kit	RIPK2	O43353	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3445	Human TNK2(Tyrosine Kinase, Non Receptor 2) ELISA Kit	TNK2	Q07912	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3446	Human DAF(Decay Accelerating Factor) ELISA Kit	DAF	P08174	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3447	Rat DAF(Decay Accelerating Factor) ELISA Kit	DAF		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3448	Human TNNI3K(Serine/threonine-protein kinase TNNI3K) ELISA Kit	TNNI3K	Q59H18	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3449	Human MAP2K6/MKK6(Mitogen Activated Protein Kinase Kinase 6) ELISA Kit	MAP2K6/MKK6	P52564	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3450	Human RXFP1(Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit	RXFP1	Q9HBX9	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3451	Mouse RXFP1(Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit	RXFP1	Q6R617	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3452	Human TACR1(Tachykinin Receptor 1) ELISA Kit	TACR1	P25103	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3453	Mouse TACR1(Tachykinin Receptor 1) ELISA Kit	TACR1	P30548	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER3454	Rat TACR1(Tachykinin Receptor 1) ELISA Kit	TACR1	P14600	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH3455	Human TACR2(Tachykinin Receptor 2) ELISA Kit	TACR2	P21452	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3456	Mouse TACR2(Tachykinin Receptor 2) ELISA Kit	TACR2	P30549	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3457	Rat TACR2(Tachykinin Receptor 2) ELISA Kit	TACR2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH3458	Human TACR3(Tachykinin Receptor 3) ELISA Kit	TACR3	P29371	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3459	Human PLP1(Proteolipid Protein 1, Myelin) ELISA Kit	PLP1	P60201	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3460	Human MyD88(Myeloid Differentiation Factor 88) ELISA Kit	MyD88	Q99836	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3461	Mouse MyD88(Myeloid Differentiation Factor 88) ELISA Kit	MyD88	P22366	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3462	Human POT1(Protection of Telomeres Protein 1) ELISA Kit	POT1	Q9NUX5	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3463	Human HBF(Fetal Hemoglobin) ELISA Kit	HBF		Sandwich	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-EH3464	Human PR3(Proteinase 3) ELISA Kit	PR3	P24158	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3465	Human PL(Placental Lactogen) ELISA Kit	PL	P0DML2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3466	Human PP13(Placental Protein 13) ELISA Kit	PP13	Q9UHV8	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3467	Human PLAP/ALPP(Placental Alkaline Phosphatase) ELISA Kit	PLAP/ALPP	P05187	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3468	Human CPT1A(Carnitine Palmitoyltransferase 1A, Liver) ELISA Kit	CPT1A	P50416	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3469	Human CPT1B(Carnitine Palmitoyltransferase 1B, Muscle) ELISA Kit	CPT1B	Q92523	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3470	Human OSBPL8(Oxysterol Binding Protein Like Protein 8) ELISA Kit	OSBPL8	Q9BZF1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3471	Human CFHR3(Complement Factor H Related Protein 3) ELISA Kit	CFHR3	Q02985	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3472	Human PGRP-S(Peptidoglycan Recognition Protein-S) ELISA Kit	PGRP-S	O75594	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3473	Human PPIL1(Peptidylprolyl Isomerase/Cyclophilin Like 1) ELISA Kit	PPIL1	Q9Y3C6	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3474	Human RAP1A(Ras-related protein Rap-1A) ELISA Kit	RAP1A	P62834	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3475	Human GP5(Glycoprotein V, Platelet) ELISA Kit	GP5	P40197	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EM3476	Mouse GP5(Glycoprotein V, Platelet) ELISA Kit	GP5	O08742	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER3477	Rat GP5(Glycoprotein V, Platelet) ELISA Kit	GP5		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3478	Human GP6(Glycoprotein VI, Platelet) ELISA Kit	GP6	Q9HCN6	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3479	Human GlyCAM1(Glycosylation Dependent Cell Adhesion Molecule 1) ELISA Kit	GlyCAM1	Q8IVK1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3480	Human GR α (Glucocorticoid Receptor Alpha) ELISA Kit	GR α	P04150 (1~	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3481	Mouse GR α (Glucocorticoid Receptor Alpha) ELISA Kit	GR α	P06537	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3482	Rat GR α (Glucocorticoid Receptor Alpha) ELISA Kit	GR α		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3483	Human GR β (Glucocorticoid Receptor Beta) ELISA Kit	GR β	P04150	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM3484	Mouse GR β (Glucocorticoid Receptor Beta) ELISA Kit	GR β	P06537	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3485	Rat GR β (Glucocorticoid Receptor Beta) ELISA Kit	GR β		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3486	Human GSK3 β (Glycogen Synthase Kinase 3 Beta) ELISA Kit	GSK3 β	P49841	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM3487	Mouse GSK3 β (Glycogen Synthase Kinase 3 Beta) ELISA Kit	GSK3 β	Q9WV60	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3488	Human GSK3 α (Glycogen Synthase Kinase 3 Alpha) ELISA Kit	GSK3 α	P49840	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3489	Mouse GSK3 α (Glycogen Synthase Kinase 3 Alpha) ELISA Kit	GSK3 α	Q2NL51	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER3490	Rat GSK3 α (Glycogen Synthase Kinase 3 Alpha) ELISA Kit	GSK3 α	P18265	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3491	Human A-GHRL(Acylated Ghrelin) ELISA Kit	A-GHRL	Q9UBU3	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3492	Human GLTP(Glycolipid Transfer Protein) ELISA Kit	Human GLTP	Q9N2D2	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3493	Mouse GLTP(Glycolipid Transfer Protein) ELISA Kit	GLTP	Q9JL62	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3494	Rat GLTP(Glycolipid Transfer Protein) ELISA Kit	GLTP		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3495	Human CFHR5(Complement Factor H Related Protein 5) ELISA Kit	CFHR5	Q9BXR6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3496	Human FTL(Ferritin, Light Polypeptide) ELISA Kit	FTL	P02792	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3497	Human FTH(Ferritin, Heavy Polypeptide) ELISA Kit	FTH	P02794	Sandwich	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-EM3498	Mouse FTH(Ferritin, Heavy Polypeptide) ELISA Kit	FTH	P09528	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL

FY-ER3499	Rat FTH(Ferritin, Heavy Polypeptide) ELISA Kit	FTH	P19132	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH3500	Human PIR(Pirin) ELISA Kit	PIR	O00625	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3501	Mouse PIR(Pirin) ELISA Kit	PIR	Q9D711	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3502	Human MRP1(Multidrug Resistance Associated Protein 1) ELISA Kit	MRP1	P33527	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3503	Human SYNPO(Synaptopodin) ELISA Kit	SYNPO	Q8N3V7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3504	Human SYP(Synaptophysin) ELISA Kit	SYP	P08247	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3505	Mouse SYP(Synaptophysin) ELISA Kit	SYP	Q62277	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER3506	Rat SYP(Synaptophysin) ELISA Kit	SYP	P07825	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3507	Human SAP25(Synaptosome Associated Protein 25) ELISA Kit	SAP25	P60880	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3508	Human MT(Melatonin) ELISA Kit	MT		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3509	Mouse MT(Melatonin) ELISA Kit	MT		Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER3510	Rat MT(Melatonin) ELISA Kit	MT		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3511	Human DHEA-S(Dehydroepiandrosterone Sulfate) ELISA Kit	DHEA-S		Competitive	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EM3512	Mouse DHEA-S(Dehydroepiandrosterone sulfat) ELISA Kit	DHEA-S		Competitive	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-ER3513	Rat DHEA-S(Dehydroepiandrosterone Sulfate) ELISA Kit	DHEA-S		Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3514	Human DPD(Deoxyypyridinoline) ELISA Kit	DPD		Competitive	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM3515	Mouse DPD(Deoxyypyridinoline) ELISA Kit	DPD		Competitive	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER3516	Rat DPD(Deoxyypyridinoline) ELISA Kit	DPD	O89000	Competitive	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3517	Human DNase- I (Deoxyribonuclease I) ELISA Kit	DNase- I	P24855	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3518	Human TERF1(Telomeric Repeat Binding Factor 1) ELISA Kit	TERF1	P54274	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3519	Human SA(Sialic Acid) ELISA Kit	SA		Competitive	500µg/mL	4.69 µg/mL	7.82-500µg/mL
FY-EH3520	Human SIGLEC12(Sialic Acid Binding Ig Like Lectin 12) ELISA Kit	SIGLEC12	Q96PQ1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3521	Mouse SIGLEC12(Sialic Acid Binding Ig Like Lectin 12) ELISA Kit	SIGLEC12	Q91Y57	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3522	Human SIGLEC3(Sialic Acid Binding Ig Like Lectin 3) ELISA Kit	SIGLEC3	P20138	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3523	Mouse SIGLEC3(Sialic Acid Binding Ig Like Lectin 3) ELISA Kit	SIGLEC3	Q63994	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3524	Human ITLN1(Intelectin 1/Omentin) ELISA Kit	ITLN1	Q8WWA0	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM3525	Mouse ITLN1(Intelectin 1/Omentin) ELISA Kit	ITLN1	O88310(Itn:	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3526	Rat ITLN1(Intelectin 1/Omentin) ELISA Kit	ITLN1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3527	Human TUBa3C(Tubulin Alpha 3c) ELISA Kit	TUBa3C	P0DPH7	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3528	Human CFHR4(Complement Factor H Related Protein 4) ELISA Kit	CFHR4	Q92496	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3529	Human STMN1(Stathmin 1) ELISA Kit	STMN1	P16949	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3530	Human DCX(Doublecortin) ELISA Kit	DCX	O43602	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3531	Human VDR(Vitamin D Receptor) ELISA Kit	VDR	P11473	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3532	Human CFHR2(Complement Factor H Related Protein 2) ELISA Kit	CFHR2	P36980	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3533	Human UST2(Urotensin 2) ELISA Kit	UST2	O95399	Competitive	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3534	Rat UST2(Urotensin 2) ELISA Kit	UST2		Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3535	Human CFHR1(Complement Factor H Related Protein 1) ELISA Kit	CFHR1	Q03591	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3536	Human PP(Pepsin) ELISA Kit	PP	P01298	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3537	Mouse PP(Pepsin) ELISA Kit	PP		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER3538	Rat PP(Pepsin) ELISA Kit	PP	P06303	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3539	Human PGA(Pepsinogen A) ELISA Kit	PGA	P0DJD9	Sandwich	800pg/mL	7.5 pg/mL	12.5-800pg/mL

FY-EM3540	Mouse PGA(Pepsinogen A) ELISA Kit	PGA		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3541	Human PGC(Pepsinogen C) ELISA Kit	PGC	P20142	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3542	Mouse PGC(Pepsinogen C) ELISA Kit	PGC	Q9D7R7	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER3543	Rat PGC(Pepsinogen C) ELISA Kit	PGC	P80020	Sandwich	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-EH3544	Human GKN1(Gastrokine 1) ELISA Kit	GKN1	Q9NS71	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3545	Human GKN2(Gastrokine 2) ELISA Kit	GKN2	Q86XP6	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3546	Mouse GKN2(Gastrokine 2) ELISA Kit	GKN2	Q9CQS6	Sandwich	5ng/mL	0.05 ng/mL	0.08-5ng/mL
FY-ER3547	Rat GKN2(Gastrokine 2) ELISA Kit	GKN2	Q29TV8	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3548	Human MTL(Motilin) ELISA Kit	MTL	P12872	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3549	Rat MTL(Motilin) ELISA Kit	MTL		Competitive	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EH3550	Human GT-34(Gastrin 34) ELISA Kit	GT-34	P01350(59~)	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3551	Mouse GT-34(Gastrin 34) ELISA Kit	GT-34		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3552	Rat GT-34(Gastrin 34) ELISA Kit	GT-34		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3553	Human ProGRP(Pro-Gastrin Releasing Peptide) ELISA Kit	ProGRP	P07492(54~)	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3554	Human GT(Gastrin) ELISA Kit	GT	P01350	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3555	Mouse GT(Gastrin) ELISA Kit	GT	P48757(76 ~)	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3556	Rat GT(Gastrin) ELISA Kit	GT	P04563	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3557	Human GIP(Gastric Inhibitory Polypeptide) ELISA Kit	GIP	P09681	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM3558	Mouse GIP(Gastric Inhibitory Polypeptide) ELISA Kit	GIP	P48756	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER3559	Rat GIP(Gastric Inhibitory Polypeptide) ELISA Kit	GIP	Q06145	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3560	Human MFGE8(Milk Fat Globule EGF Factor 8) ELISA Kit	MFGE8	Q08431	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3561	Human MMP23B(Matrix Metalloproteinase 23B) ELISA Kit	MMP23B	O75900	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3562	Human CK-13/KRT13(Cytokeratin 13) ELISA Kit	CK-13/KRT13	P13646	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3563	Human CK-19/KRT19(Cytokeratin 19) ELISA Kit	CK-19/KRT19	P08727	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3564	Human CK-20/KRT20(Cytokeratin 20) ELISA Kit	CK-20/KRT20		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3565	Human CYFRA21-1(Cytokeratin Fragment Antigen 21-1) ELISA Kit	CYFRA21-1	P08727	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3566	Human Bcl11A(B-Cell Lymphoma 11A) ELISA Kit	Bcl11A	Q9H165	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3567	Human CVL/EZR(Cytovillin/Ezrin) ELISA Kit	CVL/EZR	P15311	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3568	Mouse CVL/EZR(Cytovillin/Ezrin) ELISA Kit	CVL/EZR	P26040	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3569	Human CYB561(Cytochrome b-561) ELISA Kit	CYB561	P49447	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3570	Mouse CYB561(Cytochrome b-561) ELISA Kit	CYB561	Q60720	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER3571	Rat CYB561(Cytochrome b-561) ELISA Kit	CYB561		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3572	Human CYP21B(Cytochrome P450 21B) ELISA Kit	CYP21B	P08686	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3573	Human ECM1(Extracellular Matrix Protein 1) ELISA Kit	ECM1	Q16610	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3574	Mouse ECM1(Extracellular Matrix Protein 1) ELISA Kit	ECM1	Q61508	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-EH3575	Human MEPE(Matrix Extracellular Phosphoglycoprotein) ELISA Kit	MEPE	Q9NQ76	Sandwich	60ng/mL	0.56 ng/mL	0.94-60ng/mL
FY-EH3576	Human ERK1(Extracellular Signal Regulated Kinase 1) ELISA Kit	ERK1	P27361	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3577	Mouse ERK1(Extracellular Signal Regulated Kinase 1) ELISA Kit	ERK1	Q63844	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER3578	Rat ERK1(Extracellular Signal Regulated Kinase 1) ELISA Kit	ERK1	P21708	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3579	Human ERK2(Extracellular Signal Regulated Kinase 2) ELISA Kit	ERK2	P28482	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-EM3580	Mouse ERK2(Extracellular Signal Regulated Kinase 2) ELISA Kit	ERK2	P63085	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-ER3581	Rat ERK2(Extracellular Signal Regulated Kinase 2) ELISA Kit	ERK2	P63086	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3582	Human ERK3(Extracellular Signal Regulated Kinase 3) ELISA Kit	ERK3	Q16659	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3583	Human SOCS3(Suppressors of Cytokine Signaling 3) ELISA Kit	SOCS3	O14543	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3584	Mouse SOCS3(Suppressors Of Cytokine Signaling 3) ELISA Kit	SOCS3	O35718	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3585	Rat SOCS3(Suppressors Of Cytokine Signaling 3) ELISA Kit	SOCS3	O88583	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3586	Human CYTH1(Cytohesin 1) ELISA Kit	CYTH1	Q15438	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3587	Mouse CYTH1(Cytohesin 1) ELISA Kit	CYTH1	Q9QX11	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3588	Rat CYTH1(Cytohesin 1) ELISA Kit	CYTH1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3589	Human CYTIP(Cytohesin Interacting Protein) ELISA Kit	CYTIP	O60759	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3590	Mouse CYTIP(Cytohesin Interacting Protein) ELISA Kit	CYTIP	Q91VY6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3591	Rat CYTIP(Cytohesin Interacting Protein) ELISA Kit	CYTIP		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3592	Human FLNB(Filamin B, Beta) ELISA Kit	Human FLNB		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3593	Human FLNC(Filamin C, Gamma) ELISA Kit	FLNC	Q14315	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3594	Human FMOD(Fibromodulin) ELISA Kit	FMOD	Q06828	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3595	Mouse FMOD(Fibromodulin) ELISA Kit	FMOD	P50608	Sandwich	100pg/mL	0.94 pg/mL	1.57-100pg/mL
FY-ER3596	Rat FMOD(Fibromodulin) ELISA Kit	FMOD	P50609	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3597	Human FLRT2(Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit	FLRT2	O43155	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3598	Mouse FLRT2(Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit	FLRT2	Q8BLU0	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER3599	Rat FLRT2(Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit	FLRT2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3600	Human PAP(Plasmin-Antiplasmin Complex) ELISA Kit	PAP		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM3601	Mouse PAP(Plasmin-Antiplasmin Complex) ELISA Kit	PAP		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER3602	Rat PAP(Plasmin-Antiplasmin Complex) ELISA Kit	PAP		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH3603	Human Plg(Plasminogen) ELISA Kit	Plg	P00747	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER3604	Rat Plg(Plasminogen) ELISA Kit	Plg	Q01177	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH3605	Human PAI1(Plasminogen Activator Inhibitor 1) ELISA Kit	PAI1	P05121	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3606	Mouse PAI1(Plasminogen Activator Inhibitor 1) ELISA Kit	PAI1	P22777	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER3607	Rat PAI1(Plasminogen Activator Inhibitor 1) ELISA Kit	PAI1		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3608	Human PAI2(Plasminogen Activator Inhibitor 2) ELISA Kit	PAI2	P05120	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM3609	Mouse PAI2(Plasminogen Activator Inhibitor 2) ELISA Kit	PAI2	P12388	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER3610	Rat PAI2(Plasminogen Activator Inhibitor 2) ELISA Kit	PAI2	P29524	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3611	Human tPA(Tissue-type Plasminogen Activator) ELISA Kit	tPA	P00750	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM3612	Mouse tPA(Tissue-type Plasminogen Activator) ELISA Kit	tPA	P11214	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ER3613	Rat tPA(Tissue-type Plasminogen Activator) ELISA Kit	tPA	P19637	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3614	Human TAFI(Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit	TAFI	Q96IY4	Sandwich	320ng/mL	3 ng/mL	5-320ng/mL
FY-EM3615	Mouse TAFI(Thrombin activatable fibrinolysis inhibitor) ELISA Kit	TAFI	Q9JHH6	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER3616	Rat TAFI(Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit	TAFI	Q9EQV9	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3617	Human FBL(Fibrillarlin) ELISA Kit	FBL	P22087	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM3618	Mouse FBL(Fibrillarlin) ELISA Kit	FBL		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3619	Human GHR(Growth Hormone Receptor) ELISA Kit	GHR	P10912	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM3620	Mouse GHR(Growth Hormone Receptor) ELISA Kit	GHR	P16882	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3621	Human FGα(Fibrinogen Alpha) ELISA Kit	FGα	P02671	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL

FY-EM3622	Mouse FG α (Fibrinogen Alpha) ELISA Kit	FG α	E9PV24	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-ER3623	Rat FG α (Fibrinogen Alpha) ELISA Kit	FG α	P06399	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3624	Human FG β (Fibrinogen Beta) ELISA Kit	FG β	P02675	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM3625	Mouse FG β (Fibrinogen Beta) ELISA Kit	FG β	Q8K0E8	Sandwich	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-EH3626	Human FGL2(Fibrinogen Like Protein 2) ELISA Kit	FGL2	Q14314	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-EM3627	Mouse FGL2(Fibrinogen Like Protein 2) ELISA Kit	FGL2	P12804	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3628	Human SOD2(Superoxide Dismutase 2, Mitochondrial) ELISA Kit	SOD2		Sandwich	20ng/mL	0.16 ng/mL	0.31-20 ng/mL
FY-EM3629	Mouse SOD2(Superoxide Dismutase 2, Mitochondrial) ELISA Kit	SOD2	P09671	Competitive	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3630	Rat SOD2(Superoxide Dismutase 2, Mitochondrial) ELISA Kit	SOD2	P07895	Competitive	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3631	Human SNRPC(Small Nuclear Ribonucleoprotein Polypeptide C) ELISA Kit	SNRPC	P09234	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3632	Human SBEM(Small Breast Epithelial Mucin) ELISA Kit	SBEM	Q96DR8	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3633	Human TINF2(TERF1 Interacting Nuclear Factor 2) ELISA Kit	TINF2	Q9BSI4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3634	Human TSHR-Ab(Thyroid Stimulating Hormone Receptor Antibody) ELISA Kit	TSHR-Ab		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM3635	Mouse TSHR-Ab(Thyroid Stimulating Hormone Receptor Antibody) ELISA Kit	TSHR-Ab		Sandwich	100IU/mL	0.94 IU/mL	1.57-100IU/mL
FY-EH3636	Human SRP9(Signal Recognition Particle 9kDa) ELISA Kit	SRP9	P49458	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM3637	Mouse SRP9(Signal Recognition Particle 9kDa) ELISA Kit	SRP9	P49962	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3638	Human NPTX1(Neuronal Pentraxin 1) ELISA Kit	NPTX1	Q15818	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3639	Human STAT1(Signal Transducer and Activator of Transcription 1) ELISA Kit	STAT1	P42224	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3640	Mouse STAT1(Signal Transducer and Activator of Transcription 1) ELISA Kit	STAT1	P42225	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3641	Human STAT2(Signal Transducer and Activator of Transcription 2) ELISA Kit	STAT2	P52630	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3642	Mouse STAT2(Signal Transducer and Activator of Transcription 2) ELISA Kit	STAT2	Q9WVL2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3643	Human STAT4(Signal Transducer and Activator of Transcription 4) ELISA Kit	STAT4	Q14765	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3644	Mouse STAT4(Signal Transducer and Activator of Transcription 4) ELISA Kit	STAT4	P42228	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3645	Human STAT6(Signal Transducer and Activator of Transcription 6) ELISA Kit	STAT6	P42226	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3646	Mouse STAT6(Signal Transducer and Activator of Transcription 6) ELISA Kit	STAT6	P52633	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3647	Human SOX1(Sex Determining Region Y Box Protein 1) ELISA Kit	SOX1	O00570	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3648	Human SOX18(Sex Determining Region Y Box Protein 18) ELISA Kit	SOX18	P35713	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3649	Human SOX3(Sex Determining Region Y Box Protein 3) ELISA Kit	SOX3	P41225	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM3650	Mouse SOX3(Sex Determining Region Y Box Protein 3) ELISA Kit	SOX3	P53784	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3651	Human TS/TYMS(Thymidylate Synthetase) ELISA Kit	TS/TYMS	P04818	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3652	Human TLA(Thymus Leukemia Antigen) ELISA Kit	TLA		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3653	Human TP(Thymidine Phosphorylase) ELISA Kit	TP	P19971	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3654	Human TMS β 10(Thymosin Beta 10) ELISA Kit	TMS β 10	P63313	Competitive	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EH3655	Human TP5(Thymopentin) ELISA Kit	TP5		Competitive	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3656	Human NPTX2(Neuronal Pentraxin 2) ELISA Kit	NPTX2	P47972	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3657	Human VIP(Vasoactive Intestinal Peptide) ELISA Kit	VIP	P01282	Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM3658	Mouse VIP(Vasoactive Intestinal Peptide) ELISA Kit	VIP	P32648(125-	Competitive	800pg/mL	7.5 pg/mL	12.5-800pg/mL
FY-ER3659	Rat VIP(Vasoactive Intestinal Peptide) ELISA Kit	VIP	P01283	Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3660	Human VEGF165(Vascular Endothelial Growth Factor 165) ELISA Kit	VEGF165	P15692	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3661	Human GST μ 2/GSTm2(Glutathione S Transferase Mu 2, Muscle) ELISA Kit	GST μ 2/GSTm2	P28161	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3662	Human SUMO2(Small Ubiquitin Related Modifier 2) ELISA Kit	SUMO2	P61956	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EH3663	Human IDO(Indoleamine-2,3-Dioxygenase) ELISA Kit	IDO	P14902	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3664	Human OGN(Osteoglycin) ELISA Kit	OGN	P20774	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3665	Human VEGF-B(Vascular Endothelial Cell Growth Factor B) ELISA Kit	VEGF-B	P49765	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3666	Mouse VEGF-B(Vascular Endothelial cell Growth Factor B) ELISA Kit	VEGF-B	P49766	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3667	Rat VEGF-B(Vascular Endothelial Cell Growth Factor B) ELISA Kit	VEGF-B		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3668	Human VEGFR3/Flt4(Vascular Endothelial Cell Growth Factor Receptor 3) ELISA K	VEGFR3/Flt4	P35916	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3669	Mouse VEGFR3/Flt4(Vascular Endothelial Cell Growth Factor Receptor 3) ELISA K	VEGFR3/Flt4	P35917	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3670	Human ES-Ab(Endostatin Antibody) ELISA Kit	ES-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3671	Human VWF(Von Willebrand Factor) ELISA Kit	VWF	P04275	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM3672	Mouse VWF(Von Willebrand Factor) ELISA Kit	VWF	Q8CIZ8	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-ER3673	Rat VWF(Von Willebrand Factor) ELISA Kit	VWF	Q62935	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3674	Human HbC(Hemoglobin C) ELISA Kit	HbC		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3675	Human HBβ(Hemoglobin Beta) ELISA Kit	HBβ	P68871	Sandwich	100μg/mL	0.94 μg/mL	1.57-100μg/mL
FY-EH3676	Human HO1(Heme Oxygenase 1) ELISA Kit	HO1	P09601	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3677	Mouse HO1(Heme Oxygenase 1) ELISA Kit	HO1	P14901	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER3678	Rat HO1(Heme Oxygenase 1) ELISA Kit	HO1	P06762	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3679	Human HO2(Heme Oxygenase 2) ELISA Kit	HO2	P30519	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3680	Mouse HO2(Heme Oxygenase 2) ELISA Kit	HO2	O70252	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER3681	Rat HO2(Heme Oxygenase 2) ELISA Kit	HO2	P23711	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH3682	Human SGK2(Serum/Glucocorticoid Regulated Kinase 2) ELISA Kit	SGK2	Q9HBY8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3683	Human KRAS(V-Ki-Ras2 Kirsten Rat Sarcoma Viral Oncogene Homolog) ELISA Ki	KRAS	P01116	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3684	Human SEMA4C(Semaphorin 4C) ELISA Kit	SEMA4C	Q9C0C4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3685	Human SERT(Serotonin Transporter) ELISA Kit	SERT	P31645	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3686	Rat SERT(Serotonin Transporter) ELISA Kit	SERT	P31652	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3687	Human TPP(Thrombus Precursor Protein) ELISA Kit	TPP		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3688	Mouse TpP(Thrombus Precursor Protein) ELISA Kit	TpP		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER3689	Rat TPP(Thrombus Precursor Protein) ELISA Kit	TPP		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH3690	Human TXB2(Thromboxane B2) ELISA Kit	TXB2		Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM3691	Mouse TXB2(Thromboxane B2) ELISA Kit	TXB2		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3692	Rat TXB2(Thromboxane B2) ELISA Kit	TXB2		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3693	Human STC2(Stanniocalcin 2) ELISA Kit	STC2	O76061	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3694	Human FDP(Fibrinogen Degradation Product) ELISA Kit	FDP		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM3695	Mouse FDP(Fibrinogen Degradation Product) ELISA Kit	FDP		Sandwich	10μg/mL	0.1 μg/mL	0.16-10μg/mL
FY-ER3696	Rat FDP(Fibrinogen Degradation Product) ELISA Kit	FDP		Sandwich	1500pg/mL	14.06 pg/mL	23.44-1500pg/mL
FY-EH3697	Human FPA(Fibrinopeptide A) ELISA Kit	FPA	P02671	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3698	Human FPB(Fibrinopeptide B) ELISA Kit	FPB	P02675(31~	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3699	Human SSTR2(Somatostatin Receptor 2) ELISA Kit	SSTR2	P30874	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH3700	Human PHLDA2(Pleckstrin Homology Like Domain Family A, Member 2) ELISA K	PHLDA2	Q53GA4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3701	Mouse PHLDA2(Pleckstrin Homology Like Domain Family A, Member 2) ELISA Ki	PHLDA2	O08969	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER3702	Rat PHLDA2(Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit	PHLDA2		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3703	Human PAFAHIB3(Platelet Activating Factor Acetylhydrolase Ib3) ELISA Kit	PAFAHIB3	Q15102	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL

FY-EH3704	Human HAase(Hyaluronidase) ELISA Kit	HAase	Q12794	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM3705	Mouse HAase(Hyaluronidase) ELISA Kit	HAase	Q91ZJ9	Sandwich	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-ER3706	Rat HAase(Hyaluronidase) ELISA Kit	HAase		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3707	Human GP II b/IIIa(Platelet Membrane Glycoprotein II b/IIIa complex) ELISA Kit	GP II b/IIIa	P08514	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3708	Mouse GP II b/IIIa(Platelet Membrane Glycoprotein II b/IIIa complex) ELISA Kit	GP II b/IIIa	O54890	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER3709	Rat GP II b/IIIa(Platelet Membrane Glycoprotein II b/IIIa complex) ELISA Kit	GP II b/IIIa		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3710	Human ITGB3(Integrin Beta 3) ELISA Kit	ITGB3	P05106	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3711	Mouse ITGB3(Integrin Beta 3) ELISA Kit	ITGB3	O54890	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3712	Human PKP1(Plakophilin 1) ELISA Kit	PKP1	Q13835	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3713	Human PKP2(Plakophilin 2) ELISA Kit	PKP2	Q99959	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM3714	Mouse PKP2(Plakophilin 2) ELISA Kit	PKP2	Q9CQ73	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3715	Human ANGPTL8(Angiopoietin Like Protein 8) ELISA Kit	ANGPTL8	Q6UXH0	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH3716	Human PFKP(Phosphofructokinase, Platelet) ELISA Kit	PFKP	Q01813	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3717	Human PDGFR α (Platelet Derived Growth Factor Receptor Alpha) ELISA Kit	PDGFR α	P16234	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3718	Human PF4V1(Platelet Factor 4 Variant 1) ELISA Kit	PF4V1	P10720	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3719	Human GYPE(Glycophorin E) ELISA Kit	GYPE	P15421	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3720	Human SSTR4(Somatostatin Receptor 4) ELISA Kit	SSTR4	P31391	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3721	Human NOX1(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit	NOX1	Q9Y5S8	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3722	Mouse NOX1(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit	NOX1	Q8CIZ9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3723	Human EloA(Elongin A) ELISA Kit	EloA	Q14241	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3724	Mouse ELOA(Elongin A) ELISA Kit	ELOA	Q8CB77	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3725	Rat ELOA(Elongin A) ELISA Kit	ELOA		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3726	Human OxLDL-Ab(Oxidized Low Density Lipoprotein Antibody) ELISA Kit	OxLDL-Ab		Sandwich	60mIU/mL	0.56 mIU/mL	0.94-60mIU/mL
FY-EM3727	Mouse OxLDL-Ab(Oxidized Low Density Lipoprotein Antibody) ELISA Kit	OxLDL-Ab		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-ER3728	Rat OxLDL-Ab(Oxidized Low Density Lipoprotein Antibody) ELISA Kit	OxLDL-Ab		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3729	Human SSTR3(Somatostatin Receptor 3) ELISA Kit	SSTR3	P32745	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3730	Human PI(Proinsulin) ELISA Kit	PI	P01308	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3731	Mouse PI(Proinsulin) ELISA Kit	PI		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3732	Rat PI(Proinsulin) ELISA Kit	PI		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3733	Human AMY2(Amylase Alpha 2, Pancreatic) ELISA Kit	AMY2	P04746	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER3734	Rat AMY2(Amylase Alpha 2, Pancreatic) ELISA Kit	AMY2		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3735	Human GC(Glycocalicin) ELISA Kit	GC	P07359	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3736	Human GLP-2(Glucagon Like Peptide 2) ELISA Kit	GLP-2	P01275(146-	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3737	Mouse GLP-2(Glucagon Like Peptide 2) ELISA Kit	GLP-2	P55095(146-	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3738	Human JAK2(Janus Kinase 2) ELISA Kit	JAK2	O60674	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-ER3739	Rat JAK2(Janus Kinase 2) ELISA Kit	JAK2	Q62689	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3740	Human PTF1 α (Pancreas Specific Transcription Factor 1 α) ELISA Kit	PTF1 α	Q7RTS3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3741	Human REG1 α (Regenerating Islet Derived Protein 1 Alpha) ELISA Kit	REG1 α	P05451	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3742	Human REG4(Regenerating Islet Derived Protein 4) ELISA Kit	REG4	Q9BYZ8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3743	Mouse REG4(Regenerating Islet Derived Protein 4) ELISA Kit	REG4	Q9D8G5	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-ER3744	Rat REG4(Regenerating Islet Derived Protein 4) ELISA Kit	REG4		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH3745	Human SULT1A1(Sulfotransferase 1A1) ELISA Kit	SULT1A1	P50225	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3746	Human PTPN11(Protein Tyrosine Phosphatase, Non Receptor Type 11) ELISA Kit	PTPN11	Q06124	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3747	Human OSM(Oncostatin M) ELISA Kit	OSM	P13725	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3748	Mouse OSM(Oncostatin M) ELISA Kit	OSM	P53347	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3749	Rat OSM(Oncostatin M) ELISA Kit	OSM	Q65Z15	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3750	Human DCN(Decorin) ELISA Kit	DCN	P07585	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EM3751	Mouse DCN(Decorin) ELISA Kit	DCN	P28654	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER3752	Rat DCN(Decorin) ELISA Kit	DCN	Q01129	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3753	Human PHB(Prohibitin) ELISA Kit	PHB	P35232	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3754	Mouse PHB(Prohibitin) ELISA Kit	PHB	P67778	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER3755	Rat PHB(Prohibitin) ELISA Kit	PHB	P67779, Q5	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3756	Human CRY1(Cryptochrome 1) ELISA Kit	CRY1	Q16526	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3757	Mouse CRY1(Cryptochrome 1) ELISA Kit	CRY1	P97784	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3758	Rat CRY1(Cryptochrome 1) ELISA Kit	CRY1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3759	Human SEMA4D(Semaphorin 4D) ELISA Kit	SEMA4D	Q92854	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3760	Mouse SEMA4D(Semaphorin 4D) ELISA Kit	SEMA4D	O09126	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3761	Human SSTR5(Somatostatin Receptor 5) ELISA Kit	SSTR5	P35346	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3762	Human KNG1(Kininogen 1) ELISA Kit	KNG1	P01042	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3763	Human S100A2(S100 Calcium Binding Protein A2) ELISA Kit	S100A2	P29034	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3764	Human VAP-1(Vascular Adhesion Protein 1) ELISA Kit	VAP-1	Q16853	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3765	Mouse VAP-1(Vascular Adhesion Protein 1) ELISA Kit	VAP-1	O70423	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3766	Rat VAP-1(Vascular Adhesion Protein 1) ELISA Kit	VAP-1	O08590	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3767	Human PTPN6(Protein Tyrosine Phosphatase, Non Receptor Type 6) ELISA Kit	PTPN6	P29350	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3768	Human MGMT(O-6-Methylguanine DNA Methyltransferase) ELISA Kit	MGMT	P16455	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3769	Human PCDH1(Protocadherin 1) ELISA Kit	PCDH1	Q08174	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3770	Mouse PCDH1(Protocadherin 1) ELISA Kit	PCDH1	Q8CFX3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3771	Rat PCDH1(Protocadherin 1) ELISA Kit	PCDH1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3772	Human PCDHβ2(Protocadherin Beta 2) ELISA Kit	PCDHβ2	Q9Y5E7	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3773	Human FBN1(Fibrillin 1) ELISA Kit	FBN1	P35555	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3774	Mouse FBN1(Fibrillin 1) ELISA Kit	FBN1	Q61554	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3775	Human FBN2(Fibrillin 2) ELISA Kit	FBN2	P35556	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3776	Human FBN3(Fibrillin 3) ELISA Kit	FBN3	Q75N90	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3777	Human PAQR7(Progesterin and AdipoQ Receptor Family Member VII) ELISA Kit	PAQR7	Q86WK9	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3778	Human PGR(Progesterone Receptor) ELISA Kit	PGR	P06401	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3779	Mouse PGR(Progesterone receptor) ELISA Kit	PGR		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3780	Rat PGR(Progesterone Receptor) ELISA Kit	PGR	Q63449	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3781	Human sCD44(Soluble Cluster of Differentiation 44) ELISA Kit	sCD44	P16070	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3782	Human PS2(Presenilin 2) ELISA Kit	PS2	P49810	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3783	Human EGR1(Early Growth Response Protein 1) ELISA Kit	EGR1	P18146	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3784	Mouse EGR1(Early Growth Response Protein 1) ELISA Kit	EGR1	P08046	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3785	Rat EGR1(Early Growth Response Protein 1) ELISA Kit	EGR1		Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL

FY-EH3786	Human EGR2(Early Growth Response Protein 2) ELISA Kit	EGR2	P11161	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3787	Mouse EGR2(Early Growth Response Protein 2) ELISA Kit	EGR2	P08152	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3788	Rat EGR2(Early Growth Response Protein 2) ELISA Kit	EGR2		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3789	Human EGR3(Early Growth Response Protein 3) ELISA Kit	EGR3	Q06889	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3790	Mouse EGR3(Early Growth Response Protein 3) ELISA Kit	EGR3	P43300	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER3791	Rat EGR3(Early Growth Response Protein 3) ELISA Kit	EGR3		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3792	Human EGR4(Early Growth Response Protein 4) ELISA Kit	EGR4	Q05215	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3793	Mouse EGR4(Early Growth Response Protein 4) ELISA Kit	EGR4	Q9WUF2	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3794	Human NFYC(Nuclear Transcription Factor Y Gamma) ELISA Kit	NFYC	Q13952	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3795	Human MUC5AC(Mucin 5 Subtype AC) ELISA Kit	MUC5AC	P98088	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3796	Mouse MUC5AC(Mucin-5 subtype AC) ELISA Kit	MUC5AC	E9PWB6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3797	Human MUC5B(Mucin 5 Subtype B) ELISA Kit	MUC5B	Q9HC84	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM3798	Mouse MUC5B(Mucin-5 Subtype B) ELISA Kit	MUC5B	E9Q513	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3799	Human MUC7(Mucin 7) ELISA Kit	MUC7	Q8TAX7	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3800	Human eIF3a(Eukaryotic Translation Initiation Factor 3a) ELISA Kit	eIF3a	Q14152	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3801	Mouse EIF3A(Eukaryotic Translation Initiation Factor 3A) ELISA Kit	EIF3A	P23116	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3802	Human PC4(Positive Cofactor 4) ELISA Kit	PC4	P53999	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM3803	Mouse PC4(activated RNA polymerase II transcription cofactor 4) ELISA Kit	PC4	P11031	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3804	Rat PC4(Positive Cofactor 4) ELISA Kit	PC4		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3805	Human IL411(Interleukin 4 Induced Protein 1) ELISA Kit	IL411	Q96RQ9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3806	Human GYS1(Glycogen Synthase 1, Muscle) ELISA Kit	GYS1	P13807	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH3807	Human GYS2(Glycogen Synthase 2, Liver) ELISA Kit	GYS2	P54840	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER3808	Rat GYS2(Glycogen Synthase 2, Liver) ELISA Kit	GYS2		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3809	Human OSMR(Oncostatin M Receptor) ELISA Kit	OSMR	Q99650	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3810	Rat OSMR(Oncostatin M Receptor) ELISA Kit	OSMR		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3811	Human MK(Midkine) ELISA Kit	MK	P21741	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3812	Human PON1(Paraoxonase 1) ELISA Kit	PON1	P27169	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3813	Rat PON1(Paraoxonase 1) ELISA Kit	PON1	P55159	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH3814	Human HNP1-3(Neutrophil Peptide 1-3) ELISA Kit	HNP1-3		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-FN13815	Human PFN1(Profilin 1) ELISA Kit	PFN1	P07737	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3816	Human ADAM22(A Disintegrin And Metalloprotease 22) ELISA Kit	ADAM22	Q9P0K1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3817	Human TP63(Tumor Protein p63) ELISA Kit	TP63	Q9H3D4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3818	Human TNF-β(Tumor Necrosis Factor Beta) ELISA Kit	TNF-β	P01374	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3819	Mouse TNF-β(Tumor Necrosis Factor Beta) ELISA Kit	TNF-β	P09225	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3820	Rat TNF-β(Tumor Necrosis Factor Beta) ELISA Kit	TNF-β	Q06332	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3821	Human TL1A(Tumor Necrosis Factor Related Ligand 1A) ELISA Kit	TL1A	O95150	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3822	Human CCT8(T-complex protein 1 subunit theta) ELISA Kit	CCT8	P50990	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3823	Human CCT8L1(T-complex 1 protein subunit theta-like 1) ELISA Kit	CCT8L1	A6NM43	Sandwich	1500pg/mL	14.06 pg/mL	23.44-1500pg/mL
FY-EH3824	Human CCT5(T-complex 1 protein subunit epsilon) ELISA Kit	CCT5	P48643	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3825	Human CCT6A(T-complex 1 protein subunit zeta) ELISA Kit	CCT6A	P40227	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3826	Human TSTA(Tumor Specific Transplantation Antigen) ELISA Kit	TSTA		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL

FY-EH3827	Human RAG-1(Recombination Activating Gene 1) ELISA Kit	RAG-1	P15918	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3828	Human RAG-2(Recombination Activating Gene 2) ELISA Kit	RAG-2	P55895	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3829	Human CDK1(Cyclin Dependent Kinase 1) ELISA Kit	CDK1	P06493	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3830	Human CDK2(Cyclin Dependent Kinase 2) ELISA Kit	CDK2	P24941	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3831	Human CDK4(Cyclin Dependent Kinase 4) ELISA Kit	CDK4	P11802	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3832	Mouse CDK4(Cyclin Dependent Kinase 4) ELISA Kit	CDK4		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER3833	Rat CDK4(Cyclin Dependent Kinase 4) ELISA Kit	CDK4		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3834	Human CDK5(Cyclin Dependent Kinase 5) ELISA Kit	CDK5	Q00535	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3835	Mouse CDK5(Cyclin Dependent Kinase 5) ELISA Kit	CDK5	P49615	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3836	Rat CDK5(Cyclin Dependent Kinase 5) ELISA Kit	CDK5	Q03114	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3837	Human CDK7(Cyclin Dependent Kinase 7) ELISA Kit	CDK7	P50613	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH3838	Human CDK9(Cyclin Dependent Kinase 9) ELISA Kit	CDK9	P50750	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM3839	Mouse CDK9(Cyclin Dependent Kinase 9) ELISA Kit	CDK9	Q99J95	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER3840	Rat CDK9(Cyclin Dependent Kinase 9) ELISA Kit	CDK9		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3841	Human CDKN2A(Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit	CDKN2A	P42771	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM3842	Mouse CDKN2A(Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit	CDKN2A	P51480	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ER3843	Rat CDKN2A(Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit	CDKN2A	Q9R0Z3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3844	Human Ntn1(Netrin 1) ELISA Kit	Ntn1	O95631	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3845	Human Ntn4(Netrin 4) ELISA Kit	Ntn4	Q9HB63	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3846	Human MHCC/HLA-C(Major Histocompatibility Complex Class I C) ELISA Kit	MHCC/HLA-C		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3847	Human MHCE/HLA-E(Major Histocompatibility Complex Class I E) ELISA Kit	MHCE/HLA-E		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3848	Human AMACR(Alpha-methylacyl-CoA racemase) ELISA Kit	AMACR	Q9UHK6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3849	Rat AMACR(Alpha-methylacyl-CoA racemase) ELISA Kit	AMACR		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3850	Human CCNA1(Cyclin-A1) ELISA Kit	CCNA1	P78396	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3851	Mouse CCNA1(Cyclin-A1) ELISA Kit	CCNA1	Q61456	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER3852	Rat CCNA1(Cyclin-A1) ELISA Kit	CCNA1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3853	Human RELN(Reelin) ELISA Kit	RELN	P78509	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3854	Mouse RELN(Reelin) ELISA Kit	RELN	Q60841	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3855	Rat RELN(Reelin) ELISA Kit	RELN	P58751	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3856	Human TGM2(Transglutaminase 2, Tissue) ELISA Kit	TGM2	P21980	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3857	Mouse TGM2(Transglutaminase 2, Tissue) ELISA Kit	TGM2	P21981	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER3858	Rat TGM2(Transglutaminase 2, Tissue) ELISA Kit	TGM2	Q9WVJ6	Sandwich	10mIU/mL	0.1 mIU/mL	0.16-10mIU/mL
FY-EH3859	Human TGM3(Transglutaminase 3, Epidermal) ELISA Kit	TGM3	Q08188	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3860	Human IL-27(Interleukin 27) ELISA Kit	IL-27	Q8NEV9	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3861	Mouse IL-27(Interleukin 27) ELISA Kit	IL-27	Q8K316	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3862	Rat IL-27(Interleukin 27) ELISA Kit	IL-27	D3ZPW3	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3863	Human EDN(Eosinophil-Derived Neurotoxin) ELISA Kit	EDN	P10153	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3864	Human TTR(Transthyretin) ELISA Kit	TTR	P02766	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EM3865	Mouse TTR(Transthyretin) ELISA Kit	TTR	P07309	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-ER3866	Rat TTR(Transthyretin) ELISA Kit	TTR	P02767	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3867	Human TCF20(Transcription Factor 20) ELISA Kit	TCF20	Q9UGU0	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-EM3868	Mouse TCF20(Transcription Factor 20) ELISA Kit	TCF20	Q9EPQ8	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3869	Rat TCF20(Transcription factor 20) ELISA Kit	TCF20		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3870	Human GAPDH(Glyceraldehyde-3-Phosphate Dehydrogenase) ELISA Kit	GAPDH	P04406	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM3871	Mouse GAPDH(Glyceraldehyde-3-Phosphate Dehydrogenase) ELISA Kit	GAPDH	P16858	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3872	Human TFR2(Transferrin Receptor 2) ELISA Kit	TFR2	Q9UP52	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3873	Mouse TFR2(Transferrin Receptor 2) ELISA Kit	TFR2	Q9JKX3	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER3874	Rat TFR2(Transferrin Receptor 2) ELISA Kit	TFR2	B2GUY2	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3875	Human TNPO1(Transportin 1) ELISA Kit	TNPO1	Q92973	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM3876	Mouse TNPO1(Transportin 1) ELISA Kit	TNPO1	Q8BFY9	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER3877	Rat TNPO1(Transportin 1) ELISA Kit	TNPO1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3878	Human Pax(Paxillin) ELISA Kit	Pax	P49023	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3879	Mouse Pax(Paxillin) ELISA Kit	Pax	Q8V136	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3880	Human CD200R1(CD200 Receptor 1) ELISA Kit	CD200R1	Q8TD46	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH3881	Human MAOB(Amine oxidase B) ELISA Kit	MAOB	P27338	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM3882	Mouse MAOB(Amine oxidase B) ELISA Kit	MAOB	Q8BW75	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3883	Rat MAOB(Amine oxidaseB) ELISA Kit	MAOB	P19643	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3884	Human FPN(Ferroportin) ELISA Kit	FPN	Q9NP59	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3885	Mouse FPN(Ferroportin) ELISA Kit	FPN	Q9JHI9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3886	Rat FPN(Ferroportin) ELISA Kit	FPN	Q923U9	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3887	Human HEPH(Hephaestin) ELISA Kit	HEPH	Q9BQS7	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3888	Human CALP(Calprotectin) ELISA Kit	CALP	P05109(S10	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM3889	Mouse CALP(Calprotectin) ELISA Kit	CALP	P27005(S10	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER3890	Rat CALP(Calprotectin) ELISA Kit	CALP		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH3891	Human 11,12-DHETs(11,12-dihydroxyeicosatrienoic acids) ELISA Kit	11,12-DHETs		Competitive	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3892	Human PCX(Podocalyxin) ELISA Kit	PCX	O00592	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3893	Mouse PCX(Podocalyxin) ELISA Kit	PCX	Q9R0M4	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3894	Rat PCX(Podocalyxin) ELISA Kit	PCX	Q9WTQ2	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3895	Human TPS(Tissue Polypeptide Specific Antigen) ELISA Kit	TPS		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3896	Human SCARB1(Scavenger receptor class B member 1) ELISA Kit	SCARB1	Q8WTV0	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM3897	Mouse SCARB1(Scavenger receptor class B member 1) ELISA Kit	SCARB1	Q61009	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3898	Rat SCARB1(Scavenger receptor class B member 1) ELISA Kit	SCARB1		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3899	Human HPA(Heparanase) ELISA Kit	HPA	Q9Y251	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3900	Mouse HPA(Heparanase) ELISA Kit	HPA	Q6YGZ1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER3901	Rat HPA(Heparanase) ELISA Kit	HPA	Q71RP1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH3902	Human TFPI2(Tissue Factor Pathway Inhibitor 2) ELISA Kit	TFPI2	P48307	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM3903	Mouse TFPI2(Tissue Factor Pathway Inhibitor 2) ELISA Kit	TFPI2	O35536	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3904	Human proPSA(precursor of Prostate Specific Antigen) ELISA Kit	proPSA		Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-EH3905	Human UBA52(Ubiquitin A 52 Residue Ribosomal Protein Fusion Product 1) ELISA Kit	UBA52	P62987	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3906	Human PON2(Paraoxonase 2) ELISA Kit	PON2	Q15165	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3907	Human NES1(Nesfatin 1) ELISA Kit	NES1		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3908	Mouse NES1(Nesfatin 1) ELISA Kit	NES1		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL

FY-ER3909	Rat NES1(Nesfatin 1) ELISA Kit	NES1		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3910	Human PON3(Paraoxonase 3) ELISA Kit	PON3	Q15166	Sandwich	15ng/mL	0.14 ng/mL	0.24-15ng/mL
FY-EH3911	Human DDT(D-Dopachrome Tautomerase) ELISA Kit	DDT	P30046	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3912	Human sC5b-9(Soluble Terminal Complement Complex C5b-9) ELISA Kit	Human sC5b-9		Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EM3913	Mouse sC5b-9(Soluble Terminal Complement Complex C5b-9) ELISA Kit	Mouse sC5b-9		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3914	Rat sC5b-9(Soluble Terminal Complement Complex C5b-9) ELISA Kit	Rat sC5b-9		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3915	Human UCHL1(Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit	UCHL1	P09936	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM3916	Mouse UCHL1(Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit	UCHL1	Q9R0P9	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER3917	Rat UCHL1(Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit	UCHL1	Q00981	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3918	Human UBR4(Ubiquitin Protein Ligase E3 Component N-Recognin 4) ELISA Kit	UBR4	Q5T4S7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3919	Human UBE2A(Ubiquitin Conjugating Enzyme E2A) ELISA Kit	UBE2A	P49459	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3920	Human SOD3(Superoxide Dismutase 3, Extracellular) ELISA Kit	SOD3	P08294	Competitive	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EM3921	Mouse SOD3(Superoxide Dismutase 3, Extracellular) ELISA Kit	SOD3	O09164	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3922	Rat SOD3(Superoxide Dismutase 3, Extracellular) ELISA Kit	SOD3	Q08420	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3923	Human BCAR1(Breast Cancer Anti Estrogen Resistance 1) ELISA Kit	BCAR1	P56945	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3924	Human PLA2G2D(Phospholipase A2, Group II D) ELISA Kit	PLA2G2D	Q9UNK4	Sandwich	60ng/mL	0.56 ng/mL	0.94-60ng/mL
FY-EM3925	Mouse PLA2G2D(Phospholipase A2, Group II D) ELISA Kit	PLA2G2D	Q9WVF6	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3926	Human UBAP2(Ubiquitin Associated Protein 2) ELISA Kit	UBAP2	Q5T6F2	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3927	Human MASP1(Mannan Associated Serine Protease 1) ELISA Kit	MASP1	P48740	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EH3928	Human MASP2(Mannan Associated Serine Protease 2) ELISA Kit	MASP2	O00187	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3929	Human USP8(Ubiquitin Specific Peptidase 8) ELISA Kit	USP8	P40818	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3930	Human DPP6(Dipeptidyl Peptidase VI) ELISA Kit	DPP6	P42658	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH3931	Human MRC2(Mannose Receptor C Type 2) ELISA Kit	MRC2	Q9UBG0	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3932	Human MRC1(Mannose Receptor C Type 1) ELISA Kit	MRC1	P22897	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3933	Rat MRC1(Mannose Receptor C Type 1) ELISA Kit	MRC1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3934	Human GADL1(Glutamate Decarboxylase Like Protein 1) ELISA Kit	GADL1	Q6ZQY3	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3935	Human GAD1(Glutamate Decarboxylase 1) ELISA Kit	GAD1	Q99259	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3936	Human PPAR- δ (Peroxisome Proliferator Activated Receptor Delta) ELISA Kit	PPAR- δ	Q03181	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3937	Human NUP62(Nucleoporin 62kDa) ELISA Kit	NUP62	P37198	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3938	Human PCNA(Proliferating Cell Nuclear Antigen) ELISA Kit	PCNA	P12004	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3939	Mouse PCNA(Proliferating Cell Nuclear Antigen) ELISA Kit	PCNA	P17918	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3940	Rat PCNA(Proliferating Cell Nuclear Antigen) ELISA Kit	PCNA	P04961	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3941	Human RAN κ (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit	RAN κ	Q9Y6Q6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3942	Mouse RAN κ (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit	RAN κ	O35235	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3943	Human IL-33(Interleukin 33) ELISA Kit	IL-33	O95760	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3944	Mouse IL-33(Interleukin 33) ELISA Kit	IL-33	Q8BVZ5	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3945	Human MCAM(Melanoma Cell Adhesion Molecule) ELISA Kit	MCAM	P43121	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH3946	Human MIA1(Melanoma Inhibitory Activity Protein 1) ELISA Kit	MIA1	Q16674	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3947	Human MFI2(Melanotransferrin) ELISA Kit	MFI2	P08582	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH3948	Human SPT β N4(Spectrin Beta, Non Erythrocytic 4) ELISA Kit	SPT β N4	Q9H254	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3949	Human I κ B α (Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit	I κ B α	P25963	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EM3950	Mouse IκBα(Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit	IκBα	Q9Z1E3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3951	Human HNRPA2B1(Heterogeneous Nuclear Ribonucleoprotein A2/B1) ELISA Kit	HNRPA2B1	P22626	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3952	Human HNRPU(Heterogeneous Nuclear Ribonucleoprotein U) ELISA Kit	HNRPU	Q00839	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3953	Human RARα(Retinoic Acid Receptor Alpha) ELISA Kit	RARα	P10276	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3954	Human RNASEN(Ribonuclease III, Nuclear) ELISA Kit	RNASEN	Q9NRR4	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3955	Human MCSP(Melanoma Associated Chondroitin Sulfate Proteoglycan) ELISA Kit	MCSP	Q6UVK1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3956	Human αMSH(Alpha-Melanocyte Stimulating Hormone) ELISA Kit	αMSH	P01189(138-	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3957	Mouse αMSH(Alpha-Melanocyte Stimulating Hormone) ELISA Kit	αMSH	P01193(124-	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3958	Human EP(Erythrocyte Protoporphyrin) ELISA Kit	EP		Competitive	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER3959	Rat EP(Erythrocyte Protoporphyrin) ELISA Kit	EP		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3960	Human THRα(Thyroid Hormone Receptor Alpha) ELISA Kit	THRα	P10827	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3961	Human NFRκB(Nuclear Factor Related To Kappa B Binding Protein) ELISA Kit	NFRκB	Q6P4R8	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3962	Human TNKS1(Tankyrase 1) ELISA Kit	TNKS1	O95271	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3963	Human TNKS2(Tankyrase 2) ELISA Kit	TNKS2	Q9H2K2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3964	Human PARP4(Poly ADP Ribose Polymerase 4) ELISA Kit	PARP4	Q9UKK3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3965	Human ME20M(Melanoma Associated ME20) ELISA Kit	ME20M	P40967	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3966	Human PXR(Pregnane X Receptor) ELISA Kit	PXR	O75469	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3967	Human ALOX15B(Arachidonate-15-Lipoxygenase, Type B) ELISA Kit	ALOX15B	O15296	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3968	Human RELB(V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit	RELB	Q01201	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM3969	Mouse RELB(V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit	RELB	Q04863	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER3970	Rat RELB(V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit	RELB		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3971	Human IL-24(Interleukin 24) ELISA Kit	IL-24	Q13007	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3972	Mouse IL-24(Interleukin 24) ELISA Kit	IL-24	Q925S4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3973	Rat IL-24(Interleukin 24) ELISA Kit	IL-24	Q9JI24	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3974	Human AIM2(Absent In Melanoma 2) ELISA Kit	AIM2	O14862	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3975	Human TSPAN30/CD63(Tetraspanin 30) ELISA Kit	TSPAN30/CD63	P08962	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3976	Human BANP(BTG3 Associated Nuclear Protein) ELISA Kit	BANP	Q8N9N5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3977	Human PRKG1(Protein Kinase, cGMP Dependent Type I) ELISA Kit	PRKG1	Q13976	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3978	Human PDE4D(Phosphodiesterase 4D, cAMP Specific) ELISA Kit	PDE4D	Q08499	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3979	Human CKB(Creatine Kinase, Brain) ELISA Kit	CKB	P12277	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3980	Human LOXL4(Lysyl Oxidase Like Protein 4) ELISA Kit	LOXL4	Q96JB6	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH3981	Human TNFRSF1B(Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit	TNFRSF1B	P20333	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM3982	Mouse TNFRSF1B(Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit	TNFRSF1B	P25119	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER3983	Rat TNFRSF1B(Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit	TNFRSF1B	Q80WY6	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH3984	Human LL-37(Antibacterial Protein LL-37) ELISA Kit	LL-37	P49913(134-	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH3985	Human Tα1(Thymosin-α1) ELISA Kit	Tα1	P06454(2~2)	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EM3986	Mouse Tα1(Thymosin-α1) ELISA Kit	Tα1	P26350(2~2)	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER3987	Rat Tα1(Thymosin-α1) ELISA Kit	Tα1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH3988	Human PACAP-38(Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit	PACAP-38	P18509(132-	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3989	Mouse PACAP-38(Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit	PACAP-38	O70176(131-	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3990	Rat PACAP-38(Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit	PACAP-38		Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL

FY-EH3991	Human PACAP-27(Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA K	PACAP-27	P18509(132-	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3992	Mouse PACAP-27(Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Ki	PACAP-27	O70176(131-	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER3993	Rat PACAP-27(Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit	PACAP-27		Competitive	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH3994	Human IL-35(Interleukin 35) ELISA Kit	IL-35	P29459(IL12	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM3995	Mouse IL-35(Interleukin 35) ELISA Kit	IL-35		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER3996	Rat IL-35(Interleukin 35) ELISA Kit	IL-35		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3997	Human TK1(Thymidine Kinase 1, Soluble) ELISA Kit	TK1	P04183	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH3998	Human 12-LO(Arachidonate 12-Lipoxygenase) ELISA Kit	12-LO	P18054	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM3999	Mouse 12-LO(Arachidonate 12-Lipoxygenase) ELISA Kit	12-LO	P39655	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER4000	Rat 12-LO(Arachidonate 12-Lipoxygenase) ELISA Kit	12-LO		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4001	Human VIL1(Villin 1) ELISA Kit	VIL1	P09327	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4002	Human PI3/Elafin(Peptidase Inhibitor 3, Skin Derived) ELISA Kit	PI3/Elafin	P19957	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4003	Human TTPα(Alpha Tocopherol Transfer Protein) ELISA Kit	TTPα		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4004	Human IL-21(Interleukin 21) ELISA Kit	IL-21	Q9HBE4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM4005	Mouse IL-21(Interleukin 21) ELISA Kit	IL-21	Q9ES17	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4006	Rat IL-21(Interleukin 21) ELISA Kit	IL-21	A3QPB9	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4007	Human SSADH(Succinate-semialdehyde dehydrogenase, mitochondrial) ELISA Kit	SSADH	P51649	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4008	Human MMP-15(Matrix Metalloproteinase 15) ELISA Kit	MMP-15	P51511	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4009	Human MMP-14(Matrix Metalloproteinase 14) ELISA Kit	MMP-14	P50281	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM4010	Mouse MMP-14(Matrix Metalloproteinase 14) ELISA Kit	MMP-14	P53690	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4011	Rat MMP-14(Matrix Metalloproteinase 14) ELISA Kit	MMP-14	Q10739	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4012	Human VLDL(Very Low Density Lipoprotein) ELISA Kit	VLDL		Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-ER4013	Rat VLDL(Very Low Density Lipoprotein) ELISA Kit	VLDL		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4014	Human LCAT(Lecithin-Cholesterol Acyltransferase) ELISA Kit	LCAT	P04180	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4015	Mouse LCAT(Lecithin-Cholesterol Acyltransferase) ELISA Kit	LCAT	P16301	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER4016	Rat LCAT(Lecithin-Cholesterol Acyltransferase) ELISA Kit	LCAT	P18424	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4017	Human CYP2E1(Cytochrome P450, family 2, subfamily E, polypeptide 1) ELISA Kit	CYP2E1	P05181	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4018	Mouse CYP2E1(Cytochrome P450, family 2, subfamily E, polypeptide 1) ELISA Kit	CYP2E1	Q05421	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER4019	Rat CYP2E1(Cytochrome P450, family 2, subfamily E, polypeptide 1) ELISA Kit	CYP2E1	P05182	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4020	Human MICB(Major Histocompatibility Complex Class I Chain-related gene B) ELISA Kit	MICB	Q29980	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4021	Human TNNC1(Troponin C Type 1) ELISA Kit	TNNC1	P63316	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4022	Human GALT(Galactose-1-Phosphate Uridyltransferase) ELISA Kit	GALT	P07902	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4023	Human FSHB(Follitropin Subunit Beta) ELISA Kit	FSHB	P01225	Sandwich	50mIU/mL	0.47 mIU/mL	0.79-50mIU/mL
FY-ER4024	Rat FSHB(Follitropin Subunit Beta) ELISA Kit	FSHB		Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH4025	Human LCE3B(Late Cornified Envelope Protein 3B) ELISA Kit	Human LCE3B	Q5TA77	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4026	Human CUBN (Cubilin) ELISA Kit	Human CUBN	O60494	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4027	Human CBG(Corticosteroid Binding Globulin) ELISA Kit	CBG	P08185	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-ER4028	Rat CBG(Corticosteroid Binding Globulin) ELISA Kit	CBG	P31211	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH4029	Human PTK6(Protein Tyrosine Kinase 6) ELISA Kit	PTK6	Q13882	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4030	Human TNFRSF5(Tumor Necrosis Factor Receptor Superfamily, Member 5) ELISA Kit	TNFRSF5	P25942	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4031	Human HIF-2α(Hypoxia Inducible Factor 2 Alpha) ELISA Kit	HIF-2α	Q99814	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EM4032	Mouse HIF-2 α (Hypoxia Inducible Factor 2 Alpha) ELISA Kit	HIF-2 α	P97481	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4033	Human FCRL3(Fc receptor-like protein 3) ELISA Kit	FCRL3	Q96P31	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4034	Human GPS1(G Protein Pathway Suppressor 1) ELISA Kit	GPS1	Q13098	Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-EH4035	Human CYGB(Cytoglobin) ELISA Kit	CYGB	Q8WWM9	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM4036	Mouse CYGB(Cytoglobin) ELISA Kit	CYGB	Q9CX80	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4037	Rat CYGB(Cytoglobin) ELISA Kit	CYGB	Q921A4	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4038	Human CDKN1A(Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit	CDKN1A	P38936	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM4039	Mouse CDKN1A(Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit	CDKN1A		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4040	Rat CDKN1A(Cyclin Dependent Kinase Inhibitor 1A) ELISA Kit	CDKN1A	Q64315	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4041	Human CASP6(Caspase 6) ELISA Kit	CASP6	P55212	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4042	Human Itgam/CD11b(Integrin alpha-M) ELISA Kit	Itgam/CD11b	P11215	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM4043	Mouse Itgam/CD11b(Integrin alpha-M) ELISA Kit	Itgam/CD11b	P05555	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4044	Human CLSTN1(Calsyntenin 1) ELISA Kit	CLSTN1	O94985	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4045	Human SIRT2(Sirtuin 2) ELISA Kit	SIRT2	Q81XJ6	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4046	Human TNFSF13(Tumor Necrosis Factor Ligand Superfamily, Member 13) ELISA Kit	TNFSF13	O75888	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4047	Human SKIL(SKI Like Oncogene) ELISA Kit	SKIL	P12757	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4048	Human DSG3(Desmogleins 3) ELISA Kit	DSG3	P32926	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4049	Human CD200(Cluster of Differentiation 200) ELISA Kit	CD200	P41217	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4050	Human DMBT1(Deleted In Malignant Brain Tumors 1) ELISA Kit	DMBT1	Q9UGM3	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4051	Human PLXNA2(Plexin-A2) ELISA Kit	PLXNA2	O75051	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4052	Human TGM6(Transglutaminase-6) ELISA Kit	TGM6	O95932	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH4053	Human MAP1B(Microtubule Associated Protein Light Chain 1B) ELISA Kit	MAP1B	P46821(2207	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4054	Human DYNC2H1(Cytoplasmic Dynein 2 heavy Chain 1) ELISA Kit	DYNC2H1	Q8NCM8	Competitive	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4055	Human SI(Sucrase Isomaltase) ELISA Kit	SI	P14410	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4056	Human CART(Cocaine And Amphetamine Regulated Transcript) ELISA Kit	CART	Q16568	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4057	Human DUOX2(Dual Oxidase 2) ELISA Kit	DUOX2	Q9NRD8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4058	Human GGT2(Gamma-glutamyltranspeptidase 2) ELISA Kit	GGT2	P36268	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH4059	Human s β APP(Soluble β Amyloid Precursor Protein) ELISA Kit	s β APP	P05067(18~t	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4060	Human sTNFRSF12A(Soluble Tumor Necrosis Factor Receptor Superfamily, Member 12A) ELISA Kit	sTNFRSF12A	Q9NP84	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4061	Human GLP1R(Glucagon-like Peptide 1 receptor) ELISA Kit	GLP1R	P43220	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4062	Human AR(Aldose reductase) ELISA Kit	AR	P15121	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4063	Human MAGED2(Melanoma-associated antigen D2) ELISA Kit	MAGED2		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4064	Human SYN1(Synapsin 1) ELISA Kit	Human SYN1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4065	Mouse SYN1(Synapsin 1) ELISA Kit	SYN1	O88935	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4066	Human BTB(Biotinidase) ELISA Kit	Human BTB		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4067	Human HABP2(Hyaluronan-binding protein 2) ELISA Kit	HABP2	Q14520	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4068	Human SRGN(Serglycin) ELISA Kit	SRGN	P10124	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4069	Human SAT1(Spermidine/Spermine N1-Acetyltransferase 1) ELISA Kit	SAT1	P21673	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4070	Human RNLS(Renalse) ELISA Kit	RNLS	Q5VYX0	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH4071	Human Cadm4(Cell adhesion molecule 4) ELISA Kit	Cadm4		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4072	Human KEAP1(Kelch-like ECH-associated protein 1) ELISA Kit	Human KEAP1	Q14145	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EH4073	Human NEO1(Neogenin 1) ELISA Kit	NEO1	Q92859	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4074	Human CD27BP(CD27 Binding Protein) ELISA Kit	CD27BP		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4075	Human UCN3(Urocortin 3) ELISA Kit	UCN3	Q969E3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4076	Human IRF1(Interferon Regulatory Factor 1) ELISA Kit	IRF1	P10914	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4077	Human ALPI(Intestinal alkaline phosphatase 1) ELISA Kit	ALPI	P09923	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4078	Rat ALPI(Intestinal alkaline phosphatase 1) ELISA Kit	ALPI	P15693	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4079	Human ADAM33 (A Disintegrin And Metalloprotease 33) ELISA Kit	Human ADAM33	Q9BZ11	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4080	Human SRD5a1(Steroid 5 Alpha Reductase 1) ELISA Kit	SRD5a1	P18405	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4081	Human VASP(Vasodilator Stimulated Phosphoprotein) ELISA Kit	VASP	P50552	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4082	Rat VASP(Vasodilator Stimulated Phosphoprotein) ELISA Kit	VASP	F7EWC1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4083	Human DAG1(Dystrophin Associated Glycoprotein 1) ELISA Kit	DAG1	Q14118	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4084	Human HABP1/C1QBP(Hyaluronan Binding Protein 1) ELISA Kit	HABP1/C1QBP	Q07021	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4085	Mouse HABP1/C1QBP(Hyaluronan Binding Protein 1) ELISA Kit	HABP1/C1QBP	O35658	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4086	Rat HABP1/C1QBP(Hyaluronan Binding Protein 1) ELISA Kit	HABP1/C1QBP		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4087	Human SOCS1(Suppressor of Cytokine Signaling 1) ELISA Kit	SOCS1	O15524	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4088	Human SRD5a2(Steroid 5 Alpha Reductase 2) ELISA Kit	SRD5a2	P31213	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4089	Human CRP-Ab(C Reactive Protein Antibody) ELISA Kit	CRP-Ab		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH4090	Human IL-37(Interleukin 37) ELISA Kit	IL-37	Q9NZH6	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4091	Human CXCR2(CXC-Chemokine Receptor 2) ELISA Kit	CXCR2	P25025	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4092	Human MEP1α(Meprin A Alpha) ELISA Kit	MEP1α	Q16819	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4093	Human SGPP1(Sphingosine-1-Phosphate Phosphatase 1) ELISA Kit	SGPP1	Q9BX95	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4094	Human SPHK1(Sphingosine Kinase 1) ELISA Kit	SPHK1	Q9NYA1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4095	Human suPAR(Soluble Urokinase-type Plasminogen Activator Receptor) ELISA Kit	suPAR	Q03405	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4096	Human KLF4(Krueppel-Like Factor 4) ELISA Kit	KLF4	O43474	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4097	Human IL-11(Interleukin 11) ELISA Kit	IL-11	P20809	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4098	Mouse IL-11(Interleukin 11) ELISA Kit	IL-11	P47873	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4099	Rat IL-11(Interleukin 11) ELISA Kit	IL-11	Q99MF5	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4100	Human sdLDL(Small Dense Low Density Lipoprotein) ELISA Kit	sdLDL		Sandwich	20nmol/mL	0.19 nmol/mL	0.32-20nmol/mL
FY-EH4101	Human TSHR(Thyroid Stimulating Hormone Receptor) ELISA Kit	TSHR	P16473	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4102	Human TERF2IP(Telomeric Repeat Binding Factor 2 Interacting Protein) ELISA Kit	TERF2IP	Q9NYB0	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4103	Human DBT(Dihydro-lipoamide branched chain transacylase E2) ELISA Kit	DBT	P11182	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4104	Human PMEPA1(Prostate Transmembrane Protein androgen Induced 1) ELISA Kit	PMEPA1	Q969W9	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4105	Human RGC32(Response Gene To Complement 32) ELISA Kit	Human RGC32	Q9H4X1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4106	Human KAL1(Anosmin-1) ELISA Kit	KAL1		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4107	Human AD(Adropin) ELISA Kit	AD	Q6UWT2	Sandwich	800pg/mL	7.5 pg/mL	12.5-800pg/mL
FY-ER4108	Rat AD(Adropin) ELISA Kit	AD	D3ZUG5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4109	Human Salusin β(Salusin Beta) ELISA Kit	Salusin β		Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4110	Human sENG(Soluble Endoglin) ELISA Kit	sENG	P17813	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4111	Human c-MET/HGFR(hepatocyte growth factor receptor) ELISA Kit	c-MET/HGFR	P08581	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4112	Human RECK(Reversion-inducing-cysteine-rich protein with kazal motifs) ELISA Kit	RECK	O95980	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4113	Human TPH1(Tryptophan Hydroxylase 1) ELISA Kit	TPH1	P17752	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-ER4114	Rat TPH1(Tryptophan Hydroxylase 1) ELISA Kit	TPH1	P09810	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4115	Human pMAPT/pTAU(phosphorylated microtubule-associated protein tau) ELISA Kit	pMAPT/pTAU	P10636	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4116	Mouse pMAPT/pTAU(phosphorylated microtubule-associated protein tau) ELISA Kit	pMAPT/pTAU		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4117	Rat pMAPT/pTAU(phosphorylated microtubule-associated protein tau) ELISA Kit	pMAPT/pTAU		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4118	Human HDAC1(Histone Deacetylase 1) ELISA Kit	HDAC1	Q13547	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4119	Mouse HDAC1(Histone Deacetylase 1) ELISA Kit	HDAC1	O09106	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4120	Rat HDAC1(Histone Deacetylase 1) ELISA Kit	HDAC1	Q4QQW4	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4121	Human SCUBE1(Signal Peptide, CUB and EGF-like Domain-containing Protein 1) ELISA Kit	SCUBE1	Q8IWY4	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER4122	Rat SCUBE1(Signal Peptide, CUB and EGF-like Domain-containing Protein 1) ELISA Kit	SCUBE1		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4123	Human GPX1(Glutathione Peroxidase 1) ELISA Kit	GPX1	P07203	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4124	Mouse GPX1(Glutathione Peroxidase 1) ELISA Kit	GPX1	P11352	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER4125	Rat GPX1(Glutathione Peroxidase 1) ELISA Kit	GPX1	P04041	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4126	Human PSMA(Prostate specific membrane antigen) ELISA Kit	PSMA	Q04609	Sandwich	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-EH4127	Human NID2(Nidogen 2) ELISA Kit	NID2	Q14112	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4128	Human MAC(Membrane Attack Complex) ELISA Kit	MAC		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4129	Human HPF4-Ab(Heparin/Platelet Factor 4 Antibody) ELISA Kit	HPF4-Ab		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM4130	Mouse HPF4-Ab(Heparin/Platelet Factor 4 Antibody) ELISA Kit	HPF4-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4131	Rat HPF4-Ab(Heparin/Platelet Factor 4 Antibody) ELISA Kit	HPF4-Ab		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4132	Human GRIA1(Glutamate receptor 1) ELISA Kit	GRIA1	P42261	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4133	Rat GRIA1(Glutamate receptor 1) ELISA Kit	GRIA1	P19490	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4134	Human MeCP2(Methyl CpG Binding Protein 2) ELISA Kit	MeCP2	P51608	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4135	Human MBD1(Methyl CpG Binding Domain Protein 1) ELISA Kit	MBD1	Q9UIS9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4136	Human MBD2(Methyl CpG Binding Domain Protein 2) ELISA Kit	MBD2	Q9UBB5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4137	Human DNMT1(DNA Methyltransferase 1) ELISA Kit	DNMT1	P26358	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4138	Rat DNMT1(DNA Methyltransferase 1) ELISA Kit	DNMT1	Q9Z330	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4139	Human DNMT2(DNA Methyltransferase 2) ELISA Kit	DNMT2	Q14717	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4140	Human DNMT3A(DNA Methyltransferase 3A) ELISA Kit	DNMT3A	Q9Y6K1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4141	Rat DNMT3A(DNA Methyltransferase 3A) ELISA Kit	DNMT3A		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4142	Human DNMT3B(DNA Methyltransferase 3B) ELISA Kit	DNMT3B	Q9UBC3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4143	Rat DNMT3B(DNA Methyltransferase 3B) ELISA Kit	DNMT3B		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4144	Human ROMO1(Reactive oxygen species modulator 1) ELISA Kit	ROMO1	P60602	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4145	Human CTSS(Cathepsin S) ELISA Kit	CTSS	P25774	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4146	Human Ki67P(Ki-67 Protein) ELISA Kit	Ki67P	P46013	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4147	Human HE4(epididymal protein 4) ELISA Kit	HE4	Q14508	Sandwich	50pmol/mL	0.47 pmol/mL	0.79-50pmol/mL
FY-EH4148	Human FGF18(Fibroblast Growth Factor 18) ELISA Kit	FGF18	O76093	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4149	Human TACSTD2(Tumor-associated Calcium Signal Transducer 2) ELISA Kit	TACSTD2	P09758	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH4150	Human PEBP4(Phosphatidylethanolamine Binding Protein 4) ELISA Kit	PEBP4	Q96S96	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4151	Human MR-ProADM(Mid-regional pro-adrenomedullin) ELISA Kit	MR-ProADM		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4152	Human SFRP4(secreted frizzled-related protein 4) ELISA Kit	SFRP4	Q6FHJ7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4153	Human FCN3(Ficolin 3) ELISA Kit	FCN3	O75636	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4154	Human KL(Klotho) ELISA Kit	KL	Q9UEF7	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-ER4155	Rat KL(Klotho) ELISA Kit	KL	Q9Z2Y9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4156	Human ARL1(Aldose Reductase Like Protein 1) ELISA Kit	ARL1	O60218	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4157	Human PTP1B(protein tyrosine phosphatase 1B) ELISA Kit	PTP1B	P18031	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4158	Human BIRC7(Baculoviral IAP repeat-containing protein 7) ELISA Kit	BIRC7	Q96CA5	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH4159	Human C3d(Complement Fragment 3d) ELISA Kit	C3d	P01024	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH4160	Human FGF8(Fibroblast Growth Factor 8) ELISA Kit	FGF8	P55075	Sandwich	800pg/mL	7.5 pg/mL	12.5-800pg/mL
FY-EH4161	Human FGF10(Fibroblast Growth Factor 10) ELISA Kit	FGF10	O15520	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER4162	Rat FGF10(Fibroblast Growth Factor 10) ELISA Kit	FGF10	P70492	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4163	Human ABCC5(ATP Binding Cassette Transporter C5) ELISA Kit	ABCC5	O15440	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4164	Human LAG3(Lymphocyte Activation Gene 3) ELISA Kit	LAG3	P18627	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4165	Human HGS(Hepatocyte growth factor-regulated tyrosine kinase substrate) ELISA Kit	HGS	O14964	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4166	Human LEFTY1(Left/Right Determination Factor 1) ELISA Kit	LEFTY1	O75610	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4167	Human SCGB2A2(Mammaglobin-A) ELISA Kit	SCGB2A2	Q13296	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4168	Human BLVRA(Biliverdin Reductase A) ELISA Kit	BLVRA	P53004	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4169	Human BPN(biopyrrin) ELISA Kit	BPN		Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4170	Human SLC16A8(Solute Carrier Family 16, Member 8) ELISA Kit	SLC16A8	O95907	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4171	Human IER3(Immediate Early Response 3) ELISA Kit	IER3	P46695	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4172	Human ADAMTS9(A Disintegrin And Metalloproteinase With Thrombospondin 9) ELISA Kit	ADAMTS9	Q9P2N4	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4173	Human HES1(Transcription factor HES-1) ELISA Kit	HES1	Q14469	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4174	Human CSNK2A1(Casein Kinase 2, alpha 1 Polypeptide) ELISA Kit	CSNK2A1	P68400	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4175	Human CSNK2B(Casein Kinase 2, beta) ELISA Kit	CSNK2B	P67870	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4176	Rat CSNK2B(Casein Kinase 2, beta) ELISA Kit	CSNK2B		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4177	Human SMRP(Soluble Mesothelin Related Peptide) ELISA Kit	SMRP		Sandwich	20pmol/mL	0.19 pmol/mL	0.32-20pmol/mL
FY-EH4178	Human CCR4(Chemokine C-C-Motif Receptor 4) ELISA Kit	Human CCR4	P51679	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4179	Human Pro-ADM(Proadrenomedullin) ELISA Kit	Pro-ADM	P35318	Sandwich	50pmol/mL	0.47 pmol/mL	0.79-50pmol/mL
FY-ER4180	Rat Pro-ADM(Proadrenomedullin) ELISA Kit	Pro-ADM		Sandwich	100pmol/mL	0.94 pmol/mL	1.57-100pmol/mL
FY-EH4181	Human CECR1(Cat Eye Syndrome Chromosome Region, Candidate 1) ELISA Kit	CECR1	Q9NZK5	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4182	Human PROM1(Prominin-1) ELISA Kit	PROM1	O43490	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4183	Human HEPc25(Hepcidin 25) ELISA Kit	HEPC25	P81172(60~80)	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4184	Human HOXA10(Homeobox protein Hox-A10) ELISA Kit	HOXA10	P31260	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4185	Human HOXB4(Homeobox protein Hox-B4) ELISA Kit	HOXB4	P17483	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4186	Human MC5R(Melanocortin 5 Receptor) ELISA Kit	MC5R	P33032	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4187	Human QSOX1(Sulfhydryl oxidase 1) ELISA Kit	Human QSOX1	O00391	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4188	Human ATRN(Attractin) ELISA Kit	ATRN	O75882	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4189	Human SCG3(Secretogranin III) ELISA Kit	Human SCG3	Q8WXD2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4190	Human IL-28A(Interleukin 28A) ELISA Kit	IL-28A	Q8IZJ0	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4191	Human IL-28B(Interleukin 28B) ELISA Kit	IL-28B	Q8IZI9	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4192	Human IL-29(Interleukin 29) ELISA Kit	IL-29	Q8IU54	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4193	Human C4d(Complement Fragment 4d) ELISA Kit	C4d		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4194	Human VNN1(Vanin 1) ELISA Kit	VNN1	O95497	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4195	Human ANXA1(Annexin A1) ELISA Kit	ANXA1	P04083	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EH4196	Human BCL2A1(Bcl-2-related protein A1) ELISA Kit	Human BCL2A1	Q16548	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4197	Human TIMD4(T-cell immunoglobulin and mucin domain-containing protein 4) ELISA Kit	Human TIMD4	Q96H15	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4198	Human VSNL1(Visinin Like Protein 1) ELISA Kit	VSNL1	P62760	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4199	Human MC-Ab(Melanocyte Antibody) ELISA Kit	MC-Ab		Sandwich	20IU/mL	0.19 IU/mL	0.32-20IU/mL
FY-EH4200	Human ANXA3(Annexin A3) ELISA Kit	ANXA3	P12429	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4201	Human AKAP12(A Kinase Anchor Protein 12) ELISA Kit	AKAP12	Q02952	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4202	Human C1QTNF2(C1q and Tumor Necrosis Factor Related Protein 2) ELISA Kit	C1QTNF2	Q9BXJ5	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4203	Human C1QTNF3(C1q and Tumor Necrosis Factor Related Protein 3) ELISA Kit	C1QTNF3	Q9BXJ4	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4204	Human C1QTNF9(C1q and Tumor Necrosis Factor Related Protein 9) ELISA Kit	C1QTNF9	P0C862	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4205	Human FAP α (Fibroblast Activation Protein Alpha) ELISA Kit	Human FAP α	Q12884	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4206	Human GDN(Glia Derived Nexin) ELISA Kit	GDN	P07093	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4207	Mouse GDN(Glia Derived Nexin) ELISA Kit	GDN	Q07235	Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-ER4208	Rat GDN(Glia Derived Nexin) ELISA Kit	GDN		Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-EH4209	Human PLOD3(Procollagen Lysine-2-Oxoglutarate-5-Dioxygenase 3) ELISA Kit	Human PLOD3	O60568	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4210	Human CLIP(Corticotropin Like Intermediate Lobe Peptide) ELISA Kit	CLIP		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4211	Human SESN1(Sestrin 1) ELISA Kit	SESN1	Q9Y6P5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4212	Mouse SESN1(Sestrin 1) ELISA Kit	SESN1	P58006	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4213	Human LAMA4(Laminin Alpha 4) ELISA Kit	LAMA4	Q16363	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4214	Human CSNK2A2(Casein Kinase 2, alpha 2 Polypeptide) ELISA Kit	CSNK2A2	P19784	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4215	Human STUB1(E3 ubiquitin-protein ligase CHIP) ELISA Kit	Human STUB1	Q9UNE7	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4216	Mouse STUB1(E3 ubiquitin-protein ligase CHIP) ELISA Kit	STUB1	Q9WUD1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4217	Human pNF-H (Phosphorylated Neurofilament H) ELISA Kit	Human pNF-H	P12036	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4218	Human DHFR(Dihydrofolate Reductase) ELISA Kit	DHFR	P00374	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4219	Human sFAS/sAPO-1(Soluble Factor-Related Apoptosis) ELISA Kit	sFAS/sAPO-1	P25445	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4220	Rat sFAS/sAPO-1(Soluble Factor-Related Apoptosis) ELISA Kit	sFAS/sAPO-1		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4221	Human ADAMTS1(A Disintegrin and Metalloproteinase with Thrombospondin Motif Type 1) ELISA Kit	Human ADAMTS1	Q9UHI8	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4222	Human SPINK1(Serine Peptidase Inhibitor Kazal Type 1) ELISA Kit	Human SPINK1	P00995	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4223	Human WISP1(WNT1 Inducible Signaling Pathway Protein 1) ELISA Kit	WISP1	O95388	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4224	Human SEMG1(Semenogelin I) ELISA Kit	SEMG1	P04279	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4225	Human SFRP5(secreted frizzled-related protein 5) ELISA Kit	SFRP5	Q5T4F7	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4226	Human TMPRSS6(Transmembrane Protease, Serine 6) ELISA Kit	TMPRSS6	Q8IU80	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4227	Human C1QTNF12(C1q and Tumor Necrosis Factor Related Protein 12) ELISA Kit	Human C1QTNF12	Q5T7M4	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4228	Human SERPINA4(Kallistatin) ELISA Kit	SERPINA4	P29622	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH4229	Human ADAMTS7(A Disintegrin And Metalloproteinase With Thrombospondin Motif Type 7) ELISA Kit	Human ADAMTS7	Q9UKP4	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4230	Human JAK1(Janus Kinase 1) ELISA Kit	JAK1	P23458	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4231	Human CAPN1(Calpain 1) ELISA Kit	CAPN1	P07384	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4232	Rat CAPN1(calpain 1) ELISA Kit	CAPN1	P97571	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4233	Human IRS1(Insulin Receptor Substrate 1) ELISA Kit	IRS1	P35568	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4234	Rat IRS1(Insulin Receptor Substrate 1) ELISA Kit	IRS1	P35570	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4235	Human CAMK2N2(Calcium/Calmodulin Dependent Protein Kinase II Inhibitor 2) ELISA Kit	Human CAMK2N2	Q96S95	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4236	Human BPDE(Benzo Pyrene Dihydrodiol Epoxide) ELISA Kit	Human BPDE		Competitive	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL

FY-EH4237	Human ADAMTS12(A Disintegrin And Metalloproteinase With Thrombospondin 12) ELISA Kit	ADAMTS12	P58397	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4238	Human sRANKL(Soluble Receptor Activator of Nuclear factor-kB Ligand) ELISA Kit	sRANKL	O14788(140	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4239	Human JAK3(Janus Kinase 3) ELISA Kit	JAK3	P52333	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4240	Human GLDL(Glycated Low Density Lipoprotein) ELISA Kit	GLDL		Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH4241	Human ITLN2(Intelectin 2) ELISA Kit	Human ITLN2	Q8WWU7	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH4242	Human TPMT(Thiopurine S-Methyltransferase) ELISA Kit	TPMT	P51580	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4243	Human CYP2A6(Cytochrome P450 2A6) ELISA Kit	CYP2A6	P11509	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4244	Human BCHE(Butyrylcholinesterase) ELISA Kit	BCHE	P06276	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4245	Human SCUBE3(Signal Peptide, CUB and EGF-like Domain-containing Protein 3) ELISA Kit	SCUBE3	Q81X30	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4246	Human CRY2(Cryptochrome 2) ELISA Kit	CRY2	Q49AN0	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4247	Human RCN1(Reticulocalbin 1) ELISA Kit	RCN1	Q15293	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4248	Human GLDC(Glycine Dehydrogenase) ELISA Kit	Human GLDC		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4249	Human CD5L(CD5 Antigen-like) ELISA Kit	CD5L	O43866	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4250	Human COX2(Cytochrome C Oxidase Subunit II) ELISA Kit	COX2	P00403	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4251	Human PEPD(Peptidase D) ELISA Kit	PEPD	P12955	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4252	Human XBP1(X-Box Binding Protein 1) ELISA Kit	Human XBP1	P17861	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4253	Human XOD(Xanthine Oxidase) ELISA Kit	XOD	P47989	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4254	Rat XOD(Xanthine Oxidase) ELISA Kit	XOD		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4255	Human pro-IL1β(pro interleukin 1 beta) ELISA Kit	Human pro-IL1β	P01584	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4256	Human KRT1(Keratin 1) ELISA Kit	Human KRT1	P04264	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4257	Human ACC(Acetyl-CoA Carboxylase) ELISA Kit	Human ACC	Q13085(AC	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4258	Rat ACC(Acetyl-CoA Carboxylase) ELISA Kit	ACC		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4259	Human CRMP5(Collapsin Response Mediator Protein 5) ELISA Kit	CRMP5	Q9BPU6	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4260	Human GRP78(Glucose Regulated Protein 78) ELISA Kit	Human GRP78	P11021	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM4261	Mouse GRP78(Glucose Regulated Protein 78) ELISA Kit	GRP78	P20029	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4262	Human VCAM-1/CD106(Vascular Cell Adhesion Molecule 1) ELISA Kit	VCAM-1/CD106	P19320	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4263	Mouse VCAM-1/CD106(Vascular cell adhesion molecule 1) ELISA Kit	VCAM-1/CD106	P29533	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4264	Rat VCAM-1/CD106(Vascular Cell Adhesion Molecule 1) ELISA Kit	VCAM-1/CD106	P29534	Sandwich	800pg/mL	7.5 pg/mL	12.5-800pg/mL
FY-EH4265	Human HSL(Hormone Sensitive Lipase) ELISA Kit	Human HSL	Q05469	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4266	Human AITRL/TNFSF18(Activation-inducible TNF-related Ligand) ELISA Kit	Human AITRL/1	Q9UNG2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4267	Human ADAMTS5(A Disintegrin And Metalloproteinase With Thrombospondin 5) ELISA Kit	ADAMTS5	Q9UNA0	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4268	Human EAAT3(Excitatory Amino Acid Transporter 3) ELISA Kit	EAAT3	P43005	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER4269	Rat EAAT3(Excitatory Amino Acid Transporter 3) ELISA Kit	EAAT3		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4270	Human STIP1(Stress Induced Phosphoprotein 1) ELISA Kit	STIP1	P31948	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4271	Human MIOX(Myoinositol Oxygenase) ELISA Kit	MIOX	Q9UGB7	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4272	Rat MIOX(Myoinositol Oxygenase) ELISA Kit	MIOX		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4273	Human FAAH(Fatty-acid amide hydrolase) ELISA Kit	FAAH	O00519	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4274	Human ATGL(Adipose Triglyceride Lipase) ELISA Kit	ATGL	Q96AD5	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER4275	Rat ATGL(Adipose Triglyceride Lipase) ELISA Kit	ATGL		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4276	Human ADAMTS8(A Disintegrin And Metalloproteinase With Thrombospondin 8) ELISA Kit	ADAMTS8	Q9UP79	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4277	Human NLR4(NLR family CARD domain-containing protein 4) ELISA Kit	NLR4	Q9NPP4	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL

FY-EH4278	Human NFATC2(Nuclear Factor Of Activated T-Cells, Cytoplasmic 2) ELISA Kit	NFATC2	Q13469	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4279	Human ELA1(Elastase 1, Pancreatic) ELISA Kit	ELA1	Q9UNI1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4280	Human ADAMTS18(A Disintegrin And Metalloproteinase With Thrombospondin 18) ELISA Kit	ADAMTS18	Q8TE60	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4281	Human CST1(Cystatin-SN) ELISA Kit	CST1	P01037	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4282	Human TBX2(T-Box Protein 2) ELISA Kit	TBX2	Q13207	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4283	Rat TBX2(T-Box Protein 2) ELISA Kit	TBX2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4284	Human IREB2(Iron Responsive Element Binding Protein 2) ELISA Kit	IREB2	P48200	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4285	Human SPX(Spexin) ELISA Kit	SPX	Q9BT56	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4286	Human CTSA(Cathepsin A) ELISA Kit	CTSA	P10619	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4287	Human FGF3(Fibroblast Growth Factor 3) ELISA Kit	Human FGF3	P11487	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4288	Human STATH(Statherin) ELISA Kit	STATH	P02808	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH4289	Human UCP-3 (Uncoupling Protein 3, Mitochondrial) ELISA Kit	Human UCP-3	P55916	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4290	Human FABP12(Fatty Acid Binding Protein 12) ELISA Kit	FABP12	A6NFH5	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4291	Human TGFBR1(TGF-beta receptor type 1) ELISA Kit	TGFBR1	P36897	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4292	Human ADORA2a(Adenosine A2a Receptor) ELISA Kit	ADORA2a	P29274	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4293	Rat ADORA2a(Adenosine A2a Receptor) ELISA Kit	ADORA2a		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4294	Human FGF14 (Fibroblast Growth Factor 14) ELISA Kit	Human FGF14	Q92915	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4295	Human KISS-54(Kisspeptin-54) ELISA Kit	KISS-54	Q15726(68~	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4296	Human HPX(Hemopexin) ELISA Kit	HPX	P02790	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4297	Rat HPX(Hemopexin) ELISA Kit	HPX	P20059	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH4298	Human CHIT1(Chitotriosidase-1) ELISA Kit	CHIT1	Q13231	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4299	Human ADP-HMW(High Molecular Weight Adiponectin) ELISA Kit	ADP-HMW	Q15848	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH4300	Human EIF2αK3(Eukaryotic Translation Initiation Factor 2 Alpha Kinase 3) ELISA Kit	Human EIF2αK3	Q9NZJ5	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4301	Human Smad2/MADH2(Mothers Against Decapentaplegic Homolog 2) ELISA Kit	Smad2/MADH2	Q15796	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4302	Rat Smad2/MADH2(Mothers Against Decapentaplegic Homolog 2) ELISA Kit	Smad2/MADH2	Q70436	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4303	Human ADORA1(Adenosine A1 Receptor) ELISA Kit	ADORA1	P30542	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4304	Rat ADORA1(Adenosine A1 Receptor) ELISA Kit	ADORA1	P25099	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4305	Human VIM-Ab(Vimentin Antibody) ELISA Kit	VIM-Ab		Sandwich	50mIU/mL	0.47 mIU/mL	0.79-50mIU/mL
FY-EH4306	Human ART3(Ecto-ADP-ribosyltransferase 3) ELISA Kit	Human ART3	Q13508	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4307	Human TWIST1(Twist-related protein 1) ELISA Kit	TWIST1	Q15672	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM4308	Mouse TWIST1(Twist-related protein 1) ELISA Kit	TWIST1	P26687	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4309	Human Smad6/MADH6(Mothers Against Decapentaplegic Homolog 6) ELISA Kit	Smad6/MADH6	Q43541	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4310	Human SH3YL1(SH3 domain-containing YSC84-like protein 1) ELISA Kit	SH3YL1	Q96HL8	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4311	Mouse SH3YL1(SH3 domain-containing YSC84-like protein 1) ELISA Kit	SH3YL1	O08641	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4312	Human RXRα(Retinoid X Receptor Alpha) ELISA Kit	Human RXRα	P19793	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4313	Human CDH19(Cadherin 19) ELISA Kit	CDH19	Q96JQ0	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4314	Human ADORA3(Adenosine A3 Receptor) ELISA Kit	ADORA3	P0DMS8	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER4315	Rat ADORA3(Adenosine A3 Receptor) ELISA Kit	ADORA3	P28647	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4316	Human TGFBR2(TGF-beta receptor type 2) ELISA Kit	TGFBR2	P37173	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4317	Human FGF22(Fibroblast Growth Factor 22) ELISA Kit	FGF22	Q9HCT0	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4318	Human CD3e(Cluster of Differentiation 3e) ELISA Kit	CD3e	P07766	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EH4319	Human EIF2a(Eukaryotic Translation Initiation Factor 2 Alpha) ELISA Kit	Human EIF2a	Q9BY44	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4320	Human SAA4(Serum Amyloid A4) ELISA Kit	SAA4	P35542	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4321	Human SSB(Sjogren Syndrome Antigen B) ELISA Kit	Human SSB	P05455	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4322	Human PRKACA(cAMP-dependent protein kinase catalytic subunit alpha) ELISA Kit	PRKACA	P17612	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4323	Human PLEC(Plectin) ELISA Kit	PLEC	Q15149	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4324	Human SIK2(Serine/threonine protein kinase SIK2) ELISA Kit	SIK2	Q9H0K1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4325	Human ANKRD12(Ankyrin Repeat Domain-containing Protein 12) ELISA Kit	ANKRD12	Q6UB98	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4326	Human CD3g(Cluster of Differentiation 3g) ELISA Kit	CD3g	P09693	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4327	Human CMA1(Chymase 1, Mast Cell) ELISA Kit	CMA1		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM4328	Mouse CMA1(Chymase 1, Mast Cell) ELISA Kit	CMA1	P21844	Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-EH4329	Human HDGF(Hepatoma Derived Growth Factor) ELISA Kit	HDGF	P51858	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4330	Mouse HDGF(Hepatoma Derived Growth Factor) ELISA Kit	HDGF	P51859	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4331	Rat HDGF(Hepatoma Derived Growth Factor) ELISA Kit	HDGF	Q8VHK7	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4332	Human AgRP(Agouti Related Protein) ELISA Kit	AgRP	O00253	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM4333	Mouse AgRP(Agouti Related Protein) ELISA Kit	AgRP	P56473	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4334	Rat AgRP(Agouti Related Protein) ELISA Kit	AgRP		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4335	Human CA125(Carbohydrate Antigen 125) ELISA Kit	CA125		Sandwich	200IU/mL	1.04 IU/mL	3.13-200 IU/mL
FY-EH4336	Human MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	MCP-1	P13500	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM4337	Mouse MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	MCP-1	P10148	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4338	Rat MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	MCP-1	P14844	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4339	Human CFD(Complement Factor D) ELISA Kit	CFD	P00746	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4340	Mouse CFD(Complement Factor D) ELISA Kit	CFD	P03953	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-ER4341	Rat CFD(Complement Factor D) ELISA Kit	CFD	P32038	Sandwich	8000ng/mL	75 ng/mL	125-8000ng/mL
FY-EH4342	Human IL-8(Interleukin 8) ELISA Kit	IL-8	P10145	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH4343	Human CXCL16(Chemokine C-X-C-Motif Ligand 16) ELISA Kit	CXCL16	Q9H2A7	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4344	Mouse CXCL16(Chemokine C-X-C-Motif Ligand 16) ELISA Kit	CXCL16	Q8BSU2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4345	Rat CXCL16(Chemokine C-X-C-Motif Ligand 16) ELISA Kit	CXCL16	Q6AXU5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4346	Human ES(Endostatin) ELISA Kit	ES		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM4347	Mouse ES(Endostatin) ELISA Kit	ES	P39061(1591)	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4348	Rat ES(Endostatin) ELISA Kit	ES		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4349	Human Hepc(Hepcidin) ELISA Kit	Hepc	P81172	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM4350	Mouse Hepc(Hepcidin) ELISA Kit	Hepc	Q9EQ21	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-ER4351	Rat Hepc(Hepcidin) ELISA Kit	Hepc	Q99MH3	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4352	Human Gas6(Growth Arrest Specific Protein 6) ELISA Kit	Gas6	Q14393	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM4353	Mouse Gas6(Growth Arrest Specific Protein 6) ELISA Kit	Gas6	Q61592	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4354	Human sgp130(Soluble Glycoprotein 130) ELISA Kit	sgp130	P40189	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM4355	Mouse sgp130(Soluble Glucoprotein 130) ELISA Kit	sgp130	Q00560	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4356	Rat sgp130(soluble Glucoprotein 130) ELISA Kit	sgp130		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4357	Human LEP(Leptin) ELISA Kit	LEP	P41159	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4358	Mouse LEP(Leptin) ELISA Kit	LEP	P41160	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4359	Rat LEP(Leptin) ELISA Kit	LEP	P50596	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH4360	Human MAU(Microalbuminuria) ELISA Kit	MAU		Sandwich	50pg/mL	0.47 pg/mL	0.79-50pg/mL
FY-EM4361	Mouse MAU(Microalbuminuria) ELISA Kit	MAU		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4362	Rat MAU(Microalbuminuria) ELISA Kit	MAU		Sandwich	100µg/mL	0.94 µg/mL	1.57-100µg/mL
FY-EH4363	Human LH(Luteinizing Hormone) ELISA Kit	LH	P01229(LH)	Sandwich	10mIU/mL	0.09 mIU/mL	0.16-10mIU/mL
FY-EM4364	Mouse LH(Luteinizing Hormone) ELISA Kit	LH		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4365	Rat LH(Luteinizing Hormone) ELISA Kit	LH		Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH4366	Human C-P(C-Peptide) ELISA Kit	C-P	P01308	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM4367	Mouse C-P(C-Peptide) ELISA Kit	C-P	P01325	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4368	Rat C-P(C-peptide) ELISA Kit	C-P	P01322,P01325	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4369	Human OxLDL(Oxidized Low Density Lipoprotein) ELISA Kit	OxLDL	P78380	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM4370	Mouse OxLDL(Oxidized Low Density Lipoprotein) ELISA Kit	OxLDL		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4371	Rat OxLDL(Oxidized Low Density Lipoprotein) ELISA Kit	OxLDL		Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-EH4372	Human MMP-10(Matrix Metalloproteinase 10) ELISA Kit	MMP-10	P09238	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4373	Mouse MMP-10(Matrix Metalloproteinase 10) ELISA Kit	MMP-10	O55123	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER4374	Rat MMP-10(Matrix Metalloproteinase 10) ELISA Kit	MMP-10	P07152	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4375	Human MMP-13(Matrix Metalloproteinase 13) ELISA Kit	MMP-13	P45452	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4376	Mouse MMP-13(Matrix Metalloproteinase 13) ELISA Kit	MMP-13	P33435	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4377	Rat MMP-13(Matrix Metalloproteinase 13) ELISA Kit	MMP-13	P23097	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4378	Human MYO(Myoglobin) ELISA Kit	MYO	P02144	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4379	Mouse MYO(Myoglobin) ELISA Kit	MYO	P04247	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4380	Rat MYO(Myoglobin) ELISA Kit	MYO	Q9QZ76	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4381	Human GLP-1(Glucagon Like Peptide 1) ELISA Kit	GLP-1	P01275	Sandwich	100pg/mL	0.94 pg/mL	1.57-100pg/mL
FY-EM4382	Mouse GLP-1(Glucagon Like Peptide 1) ELISA Kit	GLP-1	P55095(92~107)	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4383	Rat GLP-1(Glucagon Like Peptide 1) ELISA Kit	GLP-1	P06883	Sandwich	100pg/mL	0.94 pg/mL	1.57-100pg/mL
FY-EH4384	Human C4(Complement Component 4) ELISA Kit	C4	P0C0L4(C4)	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM4385	Mouse C4(Complement Component 4) ELISA Kit	C4	P01029	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-ER4386	Rat C4(Complement Component 4) ELISA Kit	C4	P08649	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4387	Human TF(Transferrin) ELISA Kit	TF	P02787	Sandwich	200µg/mL	1.88 µg/mL	3.13-200µg/mL
FY-EM4388	Mouse TF(Transferrin) ELISA Kit	TF	Q921I1	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4389	Rat TF(Transferrin) ELISA Kit	TF	P12346	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4390	Human KIM-1(Kidney Injury Molecule 1) ELISA Kit	KIM-1	Q96D42	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM4391	Mouse KIM-1(Kidney Injury Molecule 1) ELISA Kit	KIM-1	Q5QNS5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4392	Rat KIM-1(Kidney Injury Molecule 1) ELISA Kit	KIM-1	O54947	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4393	Human P I CP(Procollagen I C-Terminal Propeptide) ELISA Kit	P I CP		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM4394	Mouse P I CP(Procollagen I C-Terminal Propeptide) ELISA Kit	P I CP	P11087	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4395	Rat P I CP(Procollagen I C-terminal Propeptide) ELISA Kit	P I CP		Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH4396	Human AChE(Acetylcholinesterase) ELISA Kit	AChE	P22303	Sandwich	1000IU/mL	9.38 IU/mL	15.63-1000IU/mL
FY-EM4397	Mouse AChE(Acetylcholinesterase) ELISA Kit	AChE	P21836	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4398	Rat AChE(Acetylcholinesterase) ELISA Kit	AChE	P37136	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4399	Human ALCAM(Activated Leukocyte Cell Adhesion Molecule) ELISA Kit	ALCAM	Q13740	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4400	Mouse ALCAM(Activated Leukocyte Cell Adhesion Molecule) ELISA Kit	ALCAM	Q61490	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-ER4401	Rat ALCAM(Activated Leukocyte Cell Adhesion Molecule) ELISA Kit	ALCAM		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4402	Human ALDOA(Aldolase A, Fructose Bisphosphate) ELISA Kit	ALDOA	P04075	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4403	Human ANGPTL2(Angiopoietin Like Protein 2) ELISA Kit	ANGPTL2	Q9UKU9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4404	Mouse ANGPTL2(Angiopoietin Like Protein 2) ELISA Kit	ANGPTL2	Q9R045	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4405	Human FGFR4(Fibroblast Growth Factor Receptor 4) ELISA Kit	FGFR4	P22455	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4406	Human CYP1A1(Cytochrome P450, family 1, subfamily A, polypeptide 1) ELISA Kit	CYP1A1	P04798	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4407	Rat CYP1A1(Cytochrome P450, family 1, subfamily A, polypeptide 1) ELISA Kit	CYP1A1	P00185	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4408	Human IGF-2(Insulin Like Growth Factor 2) ELISA Kit	IGF-2	P01344	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EM4409	Mouse IGF-2(Insulin Like Growth Factor 2) ELISA Kit	IGF-2	P09535	Sandwich	50pg/mL	0.47 pg/mL	0.79-50pg/mL
FY-ER4410	Rat IGF-2(Insulin Like Growth Factor 2) ELISA Kit	IGF-2	P01346	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-EH4411	Human IGFBP-2(Insulin-like Growth Factor Binding Protein 2) ELISA Kit	IGFBP-2	P18065	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4412	Mouse IGFBP-2(Insulin-like Growth Factor Binding Protein 2) ELISA Kit	IGFBP-2	P47877	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4413	Rat IGFBP-2(Insulin-like Growth Factor Binding Protein 2) ELISA Kit	IGFBP-2	P12843	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4414	Human ApoA2(Apolipoprotein A2) ELISA Kit	ApoA2	P02652	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EM4415	Mouse ApoA2(Apolipoprotein A2) ELISA Kit	ApoA2	P09813	24. Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4416	Rat ApoA2(Apolipoprotein A2) ELISA Kit	ApoA2	P04638	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH4417	Human ApoF(Apolipoprotein F) ELISA Kit	ApoF	Q13790	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM4418	Mouse ApoF(Apolipoprotein F) ELISA Kit	ApoF	Q91V80	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH4419	Human ApoA5(Apolipoprotein A5) ELISA Kit	ApoA5	Q6Q788	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER4420	Rat ApoA5(Apolipoprotein A5) ELISA Kit	ApoA5	Q9QUH3	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4421	Human bFGF/FGF2(Basic Fibroblast Growth Factor) ELISA Kit	bFGF/FGF2	P09038	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4422	Mouse bFGF/FGF2(Basic Fibroblast Growth Factor) ELISA Kit	bFGF/FGF2	P15655	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4423	Rat bFGF/FGF2(Basic Fibroblast Growth Factor) ELISA Kit	bFGF/FGF2	P13109	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4424	Human FGF4(Fibroblast Growth Factor 4) ELISA Kit	FGF4	P08620	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4425	Mouse FGF4(Fibroblast Growth Factor 4) ELISA Kit	FGF4	P11403	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4426	Human BMP-1(Bone Morphogenetic Protein 1) ELISA Kit	BMP-1	P13497	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4427	Mouse BMP-1(Bone Morphogenetic Protein 1) ELISA Kit	BMP-1	P98063	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4428	Rat BMP-1(Bone Morphogenetic Protein 1) ELISA Kit	BMP-1		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4429	Human BMP-15(Bone Morphogenetic Protein 15) ELISA Kit	BMP-15	O95972	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4430	Human BMP-6(Bone Morphogenetic Protein 6) ELISA Kit	BMP-6	P22004	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4431	Mouse BMP-6(Bone Morphogenetic Protein 6) ELISA Kit	BMP-6	P20722	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4432	Rat BMP-6(Bone Morphogenetic Protein 6) ELISA Kit	BMP-6	Q04906	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4433	Human CEA(Carcinoembryonic Antigen) ELISA Kit	CEA	P06731	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4434	Mouse CEA(Carcinoembryonic Antigen) ELISA Kit	CEA	Q3UKK2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER4435	Rat CEA(Carcinoembryonic Antigen) ELISA Kit	CEA		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4436	Human CKLF(Chemokine Like Factor) ELISA Kit	CKLF	Q9UBR5	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-EM4437	Mouse CKLF(Chemokine Like Factor) ELISA Kit	CKLF	Q9DAS1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4438	Rat CKLF(Chemokine Like Factor) ELISA Kit	CKLF		Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH4439	Human COL3(Collagen Type III) ELISA Kit	COL3	P02461(COI	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4440	Mouse COL3(Collagen Type III) ELISA Kit	COL3	P08121	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-ER4441	Rat COL3(Collagen Type III) ELISA Kit	COL3	P13941	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-EH4442	Human COL10(Collagen Type X) ELISA Kit	COL10	Q03692(COI	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4443	Mouse COL10(Collagen Type X) ELISA Kit	COL10	Q05306	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4444	Rat COL10(Collagen Type X) ELISA Kit	COL10		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4445	Human COLEC11(Collectin Sub-Family Member 11) ELISA Kit	COLEC11	Q9BWP8	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4446	Human C1INH(Complement 1 Inhibitor) ELISA Kit	C1INH	P05155	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EM4447	Mouse C1INH(Complement 1 Inhibitor) ELISA Kit	C1INH	P97290	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER4448	Rat C1INH(Complement 1 Inhibitor) ELISA Kit	C1INH		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH4449	Human C1q(Complement 1q) ELISA Kit	C1q	P02745(C1Q	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4450	Rat C1q(Complement 1q) ELISA Kit	C1q	P31721	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4451	Human C3(Complement Component 3) ELISA Kit	C3	P01024	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4452	Mouse C3(Complement Component 3) ELISA Kit	C3	P01027	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER4453	Rat C3(Complement Component 3) ELISA Kit	C3	P01026	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH4454	Human C2(Complement Component 2) ELISA Kit	C2	P06681	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4455	Rat C2(Complement Component 2) ELISA Kit	C2		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4456	Human CFB(Complement Factor B) ELISA Kit	CFB	P00751	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4457	Mouse CFB(Complement Factor B) ELISA Kit	CFB	P04186	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4458	Rat CFB(Complement Factor B) ELISA Kit	CFB		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4459	Human CFH(Complement Factor H) ELISA Kit	CFH	P08603	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM4460	Mouse CFH(Complement Factor H) ELISA Kit	CFH	P06909	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER4461	Rat CFH(Complement Factor H) ELISA Kit	CFH	G3V9R2	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-EH4462	Human CFI(Complement Factor I) ELISA Kit	CFI	P05156	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4463	Human C4b(Complement Component 4b) ELISA Kit	C4b	P0C0L5	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4464	Human α 2-M(Alpha-2 Macroglobulin) ELISA Kit	α 2-M	P01023	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-ER4465	Rat α 2-M(Alpha-2 Macroglobulin) ELISA Kit	α 2-M		Sandwich	600ng/mL	5.63 ng/mL	9.38-600ng/mL
FY-EH4466	Human GL β (Galactosidase Beta) ELISA Kit	GL β	P16278	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-ER4467	Rat GL β (Galactosidase Beta) ELISA Kit	GL β		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4468	Human MI γ /CXCL9(Monocyte Interferon Gamma Inducing Factor) ELISA Kit	MI γ /CXCL9	Q07325	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4469	Mouse MI γ /CXCL9(Monocyte Interferon Gamma Inducing Factor) ELISA Kit	MI γ /CXCL9	P18340	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4470	Rat MI γ /CXCL9(Monocyte Interferon Gamma Inducing Factor) ELISA Kit	MI γ /CXCL9	Q8K4B1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4471	Human PKM2(Pyruvate Kinase, Muscle) ELISA Kit	PKM2	P14618	Sandwich	1000IU/mL	9.38 IU/mL	15.63-1000IU/mL
FY-EM4472	Mouse PKM2(Pyruvate Kinase, Muscle) ELISA Kit	PKM2	P52480	Sandwich	50IU/mL	0.47 IU/mL	0.79-50IU/mL
FY-ER4473	Rat PKM2(Pyruvate Kinase, Muscle) ELISA Kit	PKM2	P11980	Sandwich	50IU/ml	0.47 IU/ml	0.79-50IU/ml
FY-EH4474	Human VTN(Vitronectin) ELISA Kit	VTN	P04004	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH4475	Human D β H(Dopamine- β -Hydroxylase) ELISA Kit	D β H	P09172	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH4476	Human HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit	HIF-1 α	Q16665	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM4477	Mouse HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit	HIF-1 α	Q61221	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-ER4478	Rat HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit	HIF-1 α	O35800	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH4479	Human LRG1(Leucine Rich Alpha-2-Glycoprotein 1) ELISA Kit	LRG1	P02750	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4480	Human GDH/GLDH(Glutamate dehydrogenase) ELISA Kit	GDH/GLDH	P00367	Sandwich	20mIU/mL	0.19 mIU/mL	0.32-20mIU/mL
FY-EM4481	Mouse GDH/GLDH(Glutamate dehydrogenase) ELISA Kit	GDH/GLDH		Sandwich	10mIU/mL	0.1 mIU/mL	0.16-10mIU/mL
FY-ER4482	Rat GDH/GLDH(Glutamate Dehydrogenase) ELISA Kit	GDH/GLDH		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL

FY-EH4483	Human GAD2(Glutamate Decarboxylase 2) ELISA Kit	GAD2	Q05329	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM4484	Mouse GAD2(Glutamate Decarboxylase 2, Acid) ELISA Kit	GAD2	P48320	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER4485	Rat GAD2(Glutamate Decarboxylase 2, Acid) ELISA Kit	GAD2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4486	Human GSTα1(Glutathione S Transferase Alpha 1) ELISA Kit	GSTα1	P08263	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4487	Human IgA(Immunoglobulin A) ELISA Kit	IgA	P01876(IGH)	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4488	Mouse IgA(Immunoglobulin A) ELISA Kit	IgA		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-ER4489	Rat IgA(Immunoglobulin A) ELISA Kit	IgA		Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH4490	Human TLN(Talin) ELISA Kit	TLN	Q9Y490	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-EM4491	Mouse TLN(Talin) ELISA Kit	TLN	P26039	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4492	Human MMP-1(Matrix Metalloproteinase 1) ELISA Kit	MMP-1	P03956	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4493	Mouse MMP-1(Matrix Metalloproteinase 1) ELISA Kit	MMP-1		Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-ER4494	Rat MMP-1(Matrix Metalloproteinase 1) ELISA Kit	MMP-1	P81563	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4495	Human IgA2(Immunoglobulin A2) ELISA Kit	IgA2	P01877(IGH)	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4496	Human MMP-9(Matrix Metalloproteinase 9) ELISA Kit	MMP-9	P14780	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4497	Mouse MMP-9(MatrixMetalloproteinase 9) ELISA Kit	MMP-9	P41245	Sandwich	5ng/mL	0.0468ng/mL	0.1-5ng/mL
FY-ER4498	Rat MMP-9(Matrix Metalloproteinase 9) ELISA Kit	MMP-9	P50282	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH4499	Human I-PTH(intact Parathormone) ELISA Kit	I-PTH	P01270	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-EM4500	Mouse I-PTH(intact Parathormone) ELISA Kit	I-PTH		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4501	Rat I-PTH(Intact Parathormone) ELISA Kit	I-PTH		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4502	Rat TG(Transglutaminase) ELISA Kit	TG	P06882	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4503	Human HP(Haptoglobin) ELISA Kit	HP	P00738	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EM4504	Mouse HP(Haptoglobin) ELISA Kit	HP	Q61646	Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-ER4505	Rat HP(Haptoglobin) ELISA Kit	HP	P06866	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH4506	Human MIP-3β/ELC/CCL19(Macrophage Inflammatory Protein 3β) ELISA Kit	MIP-3β/ELC/CC	Q99731	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM4507	Mouse MIP-3β/ELC/CCL19(Macrophage Inflammatory Protein 3β) ELISA Kit	MIP-3β/ELC/CC	O70460	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4508	Human HBsAg(Hepatitis B surface antigen) ELISA Kit	HBsAg	Q8JXB9	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4509	Human sST2(Soluble ST2) ELISA Kit	sST2	P78524	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4510	Human EGFR2(Epidermal Growth Factor Receptor 2) ELISA Kit	EGFR2	P04626	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4511	Mouse EGFR2(Epidermal Growth Factor Receptor 2) ELISA Kit	EGFR2	P70424	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4512	Human ApoB48(Apolipoprotein B48) ELISA Kit	ApoB48	Q0VD83	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4513	Human sTfR1(Soluble Transferrin Receptor 1) ELISA Kit	sTfR1	P02786	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH4514	Human ORM2(Orosomucoid 2) ELISA Kit	ORM2	P19652	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4515	Human IgG1(Immunoglobulin G1) ELISA Kit	IgG1	P01857	Sandwich	10μg/mL	0.1 μg/mL	0.16-10μg/mL
FY-EM4516	Mouse IgG1(Immunoglobulin G1) ELISA Kit	IgG1		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH4517	Human IgG3(Immunoglobulin G3) ELISA Kit	IgG3	P01860	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EM4518	Mouse IgG3(Immunoglobulin G3) ELISA Kit	IgG3	P03987	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH4519	Human IgG4(Immunoglobulin G4) ELISA Kit	IgG4	P01861	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EH4520	Human NT-proANP(N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit	NT-proANP	P01160(26~1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER4521	Rat NT-proANP(N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit	NT-proANP	P01161	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4522	Human BGN(Biglycan) ELISA Kit	BGN	P21810	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4523	Human NRG-1(Neuregulin 1) ELISA Kit	NRG-1	Q02297	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-EM4524	Mouse NRG-1(Neuregulin 1) ELISA Kit	NRG-1	Q6DR98	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4525	Rat NRG-1(Neuregulin 1) ELISA Kit	NRG-1	P43322	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4526	Human GFAP(Glial Fibrillary Acidic Protein) ELISA Kit	GFAP	P14136	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4527	Mouse GFAP(glial fibrillary acidic protein) ELISA Kit	GFAP	P03995	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4528	Rat GFAP(Glial Fibrillary Acidic Protein) ELISA Kit	GFAP	P47819	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4529	Human NT-4(Neurotrophin 4) ELISA Kit	NT-4	P34130	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4530	Mouse NT-4(Neurotrophin 4) ELISA Kit	NT-4	Q80VU4	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4531	Rat NT-4(Neurotrophin 4) ELISA Kit	NT-4	P34131	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4532	Human ARG2(Arginase II) ELISA Kit	ARG2	P78540	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4533	Human CYPA(Cyclophilin A) ELISA Kit	CYPA	P62937	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM4534	Mouse CYPA(Cyclophilin A) ELISA Kit	CYPA	P17742	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4535	Rat CYPA(Cyclophilin A) ELISA Kit	CYPA		Sandwich	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-EH4536	Human FETUB(Fetuin B) ELISA Kit	FETUB	Q9UGM5	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4537	Human KRT10(Keratin, Type I Cytoskeletal 10) ELISA Kit	KRT10	P13645	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH4538	Human KISS1(Kisspeptin 1) ELISA Kit	KISS1	Q15726	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4539	Mouse KISS1(Kisspeptin 1) ELISA Kit	KISS1	Q6Y4S4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4540	Rat KISS1(Kisspeptin 1) ELISA Kit	KISS1		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4541	Human SHBG(Sex Hormone Binding Globulin) ELISA Kit	SHBG	P04278	Sandwich	100pmol/mL	0.94 pmol/mL	1.57-100pmol/mL
FY-EM4542	Mouse SHBG(Sex Hormone Binding Globulin) ELISA Kit	SHBG	P97497	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4543	Rat SHBG(Sex Hormone Binding Globulin) ELISA Kit	SHBG	P08689	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4544	Human TECK/CCL25(Thymus Expressed Chemokine) ELISA Kit	TECK/CCL25	O15444	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM4545	Mouse TECK/CCL25(Thymus Expressed Chemokine) ELISA Kit	TECK/CCL25	O35903	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4546	Rat TECK/CCL25(Thymus Expressed Chemokine) ELISA Kit	TECK/CCL25		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4547	Human OLFM4(Olfactomedin 4) ELISA Kit	OLFM4	Q6UX06	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM4548	Mouse OLFM4(Olfactomedin 4) ELISA Kit	OLFM4	Q3UZZ4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4549	Human VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit	VE-Cadherin	P33151	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4550	Mouse VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit	VE-Cadherin	P55284	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4551	Rat VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit	VE-Cadherin		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4552	Human IgE(Immunoglobulin E) ELISA Kit	IgE	P01854	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM4553	Mouse IgE(Immunoglobulin E) ELISA Kit	IgE		Sandwich	8000ng/mL	75 ng/mL	125-8000ng/mL
FY-ER4554	Rat IgE(Immunoglobulin E) ELISA Kit	IgE		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4555	Human ALB(Albumin) ELISA Kit	ALB	P02768	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM4556	Mouse ALB(Albumin) ELISA Kit	ALB	P07724	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-ER4557	Rat ALB(Albumin) ELISA Kit	ALB	P02770	Sandwich	200µg/mL	1.88 µg/mL	3.13-200µg/mL
FY-EH4558	Human FG(Fibrinogen) ELISA Kit	FG	P02671(FGα	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EM4559	Mouse FG(Fibrinogen) ELISA Kit	FG		Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-ER4560	Rat FG(Fibrinogen) ELISA Kit	FG		Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-EH4561	Human LBP(Lipopolysaccharide Binding Protein) ELISA Kit	LBP	P18428	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM4562	Mouse LBP(Lipopolysaccharide Binding Protein) ELISA Kit	LBP	Q61805	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER4563	Rat LBP(Lipopolysaccharide Binding Protein) ELISA Kit	LBP	Q63313	Sandwich	120ng/mL	1.13 ng/mL	1.88-120ng/mL
FY-EH4564	Human DKK3(Dickkopf Related Protein 3) ELISA Kit	DKK3	Q9UBP4	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL

FY-EM4565	Mouse DKK3(Dickkopf Related Protein 3) ELISA Kit	DKK3	Q9QUN9	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4566	Rat DKK3(Dickkopf Related Protein 3) ELISA Kit	DKK3		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4567	Human iC3b(Complement iC3b) ELISA Kit	iC3b	P01024	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4568	Human SAA2(Serum Amyloid A2) ELISA Kit	SAA2	P0DJ19	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4569	Mouse SAA2(Serum Amyloid A2) ELISA Kit	SAA2	P05367(20~1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4570	Human ICAM-1/CD54(intercellular adhesion molecule 1) ELISA Kit	ICAM-1/CD54		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4571	Mouse ICAM-1/CD54(Intercellular Adhesion Molecule 1) ELISA Kit	ICAM-1/CD54	P13597	Sandwich	1000ng/mL	0.82 ng/mL	1.37-1000 ng/mL
FY-ER4572	Rat ICAM-1/CD54(Intercellular Adhesion Molecule 1) ELISA Kit	ICAM-1/CD54	Q00238	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4573	Human F8(Coagulation Factor VIII) ELISA Kit	F8	P00451	Sandwich	10IU/mL	0.1 IU/mL	0.16-10IU/mL
FY-EM4574	Mouse F8(Coagulation Factor VIII) ELISA Kit	F8	Q06194	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4575	Rat F8(Coagulation Factor VIII) ELISA Kit	F8		Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH4576	Human sVEGFR-1(soluble Vascular endothelial growth factor receptor 1) ELISA Kit	sVEGFR-1	P17948	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH4577	Human ucOC(Undercarboxylated Osteocalcin) ELISA Kit	ucOC		Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-ER4578	Rat ucOC(Undercarboxylated Osteocalcin)ELISA Kit	ucOC		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4579	Human ANG1(Angiopoietin 1) ELISA Kit	ANG1	Q15389	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM4580	Mouse ANG1(Angiopoietin 1) ELISA Kit	ANG1	O08538	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER4581	Rat ANG1(Angiopoietin 1) ELISA Kit	ANG1	O35460	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4582	Human Irisin ELISA Kit	Irisin		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4583	Mouse Irisin ELISA Kit	Irisin	Q8K4Z2	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4584	Rat Irisin ELISA Kit	Irisin		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4585	Human ADP/Acrp30(Adiponectin) ELISA Kit	ADP/Acrp30		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4586	Mouse ADP/Acrp30(Adiponectin) ELISA Kit	ADP/Acrp30	Q60994	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4587	Rat ADP/Acrp30(Adiponectin) ELISA Kit	ADP/Acrp30		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4588	Human TLR4(Toll-Like Receptor 4) ELISA Kit	TLR4	O00206	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4589	Mouse TLR4(Toll-Like Receptor 4) ELISA Kit	TLR4	Q9QUK6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4590	Rat TLR4(Toll-Like Receptor 4) ELISA Kit	TLR4	Q9QX05	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4591	Human IFN α (Interferon Alpha) ELISA Kit	IFN α		Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM4592	Mouse IFN α (Interferon Alpha) ELISA Kit	IFN α	P01572,P015	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4593	Rat IFN α (Interferon Alpha) ELISA Kit	IFN α	P05011	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4594	Human NGAL(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	NGAL	P80188	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-EM4595	Mouse NGAL(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	NGAL	P11672	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ER4596	Rat NGAL(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	NGAL	P30152	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4597	Human D2D(D-Dimer) ELISA Kit	D2D		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EM4598	Mouse D2D(D-Dimer) ELISA Kit	D2D		Competitive	5000ng/mL	46.88 ng/mL	78.13-5000 ng/mL
FY-ER4599	Rat D2D(D-Dimer) ELISA Kit	D2D		Competitive	5000ng/mL	46.88 ng/mL	78.13-5000 ng/mL
FY-EH4600	Human DPP4(Dipeptidyl Peptidase IV) ELISA Kit	DPP4	P27487	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM4601	Mouse DPP4(Dipeptidyl Peptidase IV) ELISA Kit	DPP4	P28843	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER4602	Rat DPP4(Dipeptidyl Peptidase IV) ELISA Kit	DPP4	P14740	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4603	Human sCD14(Soluble Cluster of Differentiation 14) ELISA Kit	sCD14	P08571	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EM4604	Mouse sCD14(Soluble Cluster of Differentiation 14) ELISA Kit	sCD14	P10810	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4605	Rat sCD14(Soluble Cluster of Differentiation 14) ELISA Kit	sCD14		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL

FY-EH4606	Human EMMPRIN/CD147(Extracellular Matrix Metalloproteinase Inducer) ELISA Kit	EMMPRIN/CD147	P35613	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4607	Human CTSB(Cathepsin B) ELISA Kit	CTSB	P07858	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM4608	Mouse CTSB(Cathepsin B) ELISA Kit	CTSB	P10605(80~)	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4609	Rat CTSB(Cathepsin B) ELISA Kit	CTSB	P00787	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4610	Human FABP1(Fatty Acid Binding Protein 1, Liver) ELISA Kit	FABP1	P07148	Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-EM4611	Mouse FABP1(Fatty Acid Binding Protein 1, Liver) ELISA Kit	FABP1	P12710	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4612	Rat FABP1(Fatty Acid Binding Protein 1, Liver) ELISA Kit	FABP1	P02692	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4613	Human AGP1(Alpha-1-Acid Glycoprotein 1) ELISA Kit	AGP1	P02763	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM4614	Mouse AGP1(Alpha-1-Acid Glycoprotein 1) ELISA Kit	AGP1	Q60590	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH4615	Human BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	BDNF	P23560	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4616	Mouse BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	BDNF	P21237	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4617	Rat BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	BDNF	P23363	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4618	Human CASP1(Caspase 1) ELISA Kit	CASP1	P29466	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM4619	Mouse CASP1(Caspase 1) ELISA Kit	CASP1	P29452	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4620	Rat CASP1(Caspase 1) ELISA Kit	CASP1	P43527	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4621	Human CASP3(Caspase 3) ELISA Kit	CASP3	P42574	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4622	Mouse CASP3(Caspase 3) ELISA Kit	CASP3	P70677	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4623	Rat CASP3(Caspase 3) ELISA Kit	CASP3	P55213	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4624	Human MIP-3α(Macrophage Inflammatory Protein 3 Alpha) ELISA Kit	MIP-3α	P78556	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4625	Mouse MIP-3α(Macrophage Inflammatory Protein 3 Alpha) ELISA Kit	MIP-3α	O89093	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER4626	Rat MIP-3α(Macrophage Inflammatory Protein 3 Alpha) ELISA Kit	MIP-3α	P97884	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4627	Human MDC/CCL22(Macrophage-Derived Chemokine) ELISA Kit	MDC/CCL22	O00626	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM4628	Mouse MDC/CCL22(Macrophage-Derived Chemokine) ELISA Kit	MDC/CCL22	O88430	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4629	Rat MDC/CCL22(Macrophage-Derived Chemokine) ELISA Kit	MDC/CCL22		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4630	Human CLU(Clusterin) ELISA Kit	CLU	P10909	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4631	Mouse CLU(Clusterin) ELISA Kit	CLU	Q06890	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-ER4632	Rat CLU(Clusterin) ELISA Kit	CLU	P05371	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4633	Human CRP(C-Reactive Protein) ELISA Kit	CRP	P02741	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-EM4634	Mouse CRP(C-Reactive-Protein) ELISA Kit	CRP		Sandwich	25ng/mL	0.02 ng/mL	0.03-25ng/mL
FY-ER4635	Rat CRP(C-Reactive Protein) ELISA Kit	CRP	P48199	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4636	Human ENA-78/CXCL5(Epithelial Neutrophil Activating Peptide 78) ELISA Kit	ENA-78/CXCL5	P42830	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4637	Human AH-Ab(Apolipoprotein Antibody) ELISA Kit	AH-Ab		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4638	Human SDF-1/CXCL12(Stromal Cell Derived Factor 1) ELISA Kit	SDF-1/CXCL12	P48061	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4639	Mouse SDF-1/CXCL12(Stromal Cell Derived Factor 1) ELISA Kit	SDF-1/CXCL12	P40224	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4640	Rat SDF-1/CXCL12(Stromal Cell Derived Factor 1) ELISA Kit	SDF-1/CXCL12	Q9QZD1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4641	Human αFP(Alpha-Fetoprotein) ELISA Kit	αFP	P02771	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4642	Mouse αFP(Alpha-Fetoprotein) ELISA Kit	αFP	P02772	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4643	Rat αFP(Alpha-Fetoprotein) ELISA Kit	αFP	P02773	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4644	Human APN(Aminopeptidase N) ELISA Kit	APN		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4645	Human Flt3L(FMS-like Tyrosine Kinase 3 Ligand) ELISA Kit	Flt3L	P49771	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4646	Mouse Flt3L(FMS-like Tyrosine Kinase 3 Ligand) ELISA Kit	Flt3L	P49772	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-EH4647	Human G-CSF(Granulocyte Colony-stimulating Factor) ELISA Kit	G-CSF	P09919	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-EM4648	Mouse G-CSF(Granulocyte Colony-stimulating Factor) ELISA Kit	G-CSF	P09920	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4649	Rat G-CSF(Granulocyte Colony Stimulating Factor) ELISA Kit	G-CSF	P48750	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4650	Human GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	GM-CSF	P04141	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM4651	Mouse GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	GM-CSF	P01587	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4652	Rat GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	GM-CSF	P48750	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4653	Human IGF-1(Insulin-like Growth Factor 1) ELISA Kit	IGF-1	P05019	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4654	Mouse IGF-1(Insulin-like Growth Factor 1) ELISA Kit	IGF-1	P05017	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-ER4655	Rat IGF-1(Insulin-like Growth Factor 1) ELISA Kit	IGF-1	P08025	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4656	Human IGFBP-3(Insulin-like Growth Factor Binding Protein 3) ELISA Kit	IGFBP-3	P17936	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM4657	Mouse IGFBP-3(Insulin-like Growth Factor Binding Protein 3) ELISA Kit	IGFBP-3	P47878	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-ER4658	Rat IGFBP-3(Insulin-Like Growth Factor Binding Protein 3) ELISA Kit	IGFBP-3	P15473	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4659	Human FGF7/KGF(Fibroblast Growth Factor 7) ELISA Kit	FGF7/KGF	P21781	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4660	Mouse FGF7/KGF(Fibroblast Growth Factor 7) ELISA Kit	FGF7/KGF	P36363	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4661	Rat FGF7/KGF(Fibroblast Growth Factor 7) ELISA Kit	FGF7/KGF	Q02195	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4662	Human IL-2(Interleukin 2) ELISA Kit	IL-2	P60568	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM4663	Mouse IL-2(Interleukin 2) ELISA Kit	IL-2	P04351	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4664	Rat IL-2(Interleukin 2) ELISA Kit	IL-2	P17108	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4665	Human IL-17A(Interleukin 17A) ELISA Kit	IL-17A	Q16552	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-EM4666	Mouse IL-17A(Interleukin 17A) ELISA Kit	IL-17A		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4667	Rat IL-17A(Interleukin 17A) ELISA Kit	IL-17A	Q61453	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4668	Human ACAT2(Acetyl Coenzyme A Acetyltransferase 2) ELISA Kit	ACAT2	Q9BWD1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4669	Human TnTf/TNNT3(Troponin T, fast skeletal muscle) ELISA Kit	TnTf/TNNT3	P45378	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4670	Mouse TnTf/TNNT3(Troponin T, fast skeletal muscle) ELISA Kit	TnTf/TNNT3	Q9QZ47	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4671	Rat TnTf/TNNT3(Troponin T Type 3, Fast Skeletal) ELISA Kit	TnTf/TNNT3		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4672	Human IFABP/FABP2(Intestinal Fatty Acid Binding Protein) ELISA Kit	IFABP/FABP2	P12104	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4673	Mouse IFABP/FABP2(Intestinal Fatty Acid Binding Protein) ELISA Kit	IFABP/FABP2	P55050	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4674	Rat IFABP/FABP2(Intestinal Fatty Acid Binding Protein) ELISA Kit	IFABP/FABP2	P02693	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4675	Human LP-a(Apolipoprotein a) ELISA Kit	LP-a	P08519	Sandwich	15ng/mL	0.14 ng/mL	0.24-15ng/mL
FY-EH4676	Human FE(Ferritin) ELISA Kit	FE	P02792(FTL)	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-EM4677	Mouse FE(Ferritin) ELISA Kit	FE		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4678	Rat FE(Ferritin) ELISA Kit	FE		Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH4679	Human IgG(Immunoglobulin G) ELISA Kit	IgG	P01860(IGH)	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4680	Mouse IgG(Immunoglobulin G) ELISA Kit	IgG	P01863	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4681	Rat IgG(Immunoglobulin G) ELISA Kit	IgG	P20762	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4682	Human HCG(Chorionic Gonadotropin) ELISA Kit	HCG	P01215(CG)	Sandwich	500mIU/mL	4.69 mIU/mL	7.82-500mIU/mL
FY-EH4683	Human P I NP(Procollagen I N-Terminal Propeptide) ELISA Kit	P I NP		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4684	Mouse P I NP(Procollagen I N-Terminal Propeptide) ELISA Kit	P I NP		Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-ER4685	Rat P I NP(Procollagen I N-Terminal Propeptide) ELISA Kit	P I NP		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4686	Human Actin-Ab(Actin Antibody) ELISA Kit	Actin-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4687	Human ALB-Ab(Albumin Antibody) ELISA Kit	ALB-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL

FY-EH4688	Human AP17(Apelin 17) ELISA Kit	AP17	Q9ULZ1(61-	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4689	Human INHA(Inhibin Alpha Chain) ELISA Kit	INHA	P05111	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4690	Mouse INHA(Inhibin Alpha Chain) ELISA Kit	INHA	Q04997	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4691	Rat INHA(Inhibin Alpha Chain) ELISA Kit	INHA	P17490	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4692	Human ACD(Adrenocortical Dysplasia Homolog) ELISA Kit	ACD	Q96AP0	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4693	Human AChR-Ab(Acetylcholine Receptor Antibody) ELISA Kit	AChR-Ab		Sandwich	5000fmol/mL	46.88 fmol/mL	78.13-5000fmol/mL
FY-ER4694	Rat AChR-Ab(Acetylcholine Receptor Antibody) ELISA Kit	AChR-Ab		Sandwich	4000fmol/mL	37.5 fmol/mL	62.5-4000fmol/mL
FY-EH4695	Human ACO-2(Aconitase 2) ELISA Kit	ACO-2	Q99798	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4696	Human ACP1(Acid Phosphatase 1) ELISA Kit	ACP1	P24666	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4697	Human ADCY3(Adenylate Cyclase 3) ELISA Kit	ADCY3	O60266	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4698	Human ADCY4(Adenylate Cyclase 4) ELISA Kit	ADCY4	Q8NFM4	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4699	Human ADCYAP1/PACAP(Adenylate Cyclase Activating Polypeptide 1, Pituitary) ELISA Kit	ADCYAP1/PACAP		Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH4700	Human AFAP1(Actin Filament Associated Protein 1) ELISA Kit	AFAP1	Q8N556	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4701	Human AIG1(Androgen Induced Protein 1) ELISA Kit	AIG1	Q9N5V5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4702	Human AKAP1(A Kinase Anchor Protein 1) ELISA Kit	AKAP1	Q92667	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4703	Human ALK(Anaplastic Lymphoma Kinase) ELISA Kit	ALK		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4704	Human Ang- II (Angiotensin II) ELISA Kit	Ang- II	P01019(34~	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4705	Mouse Ang- II (Angiotensin II) ELISA Kit	Ang- II	P11859(25~	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4706	Rat Ang- II (Angiotensin II) ELISA Kit	Ang- II		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4707	Human Ang- I (Angiotensin I) ELISA Kit	Ang- I	P01019(34~	Competitive	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EM4708	Mouse Ang- I (Angiotensin I) ELISA Kit	Ang- I	P11859(25~	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4709	Rat Ang- I (Angiotensin I) ELISA Kit	Ang- I		Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4710	Human ANGPTL3(Angiopoietin Like Protein 3) ELISA Kit	ANGPTL3	Q9Y5C1	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM4711	Mouse ANGPTL3(Angiopoietin Like Protein 3) ELISA Kit	ANGPTL3	Q9R182	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4712	Rat ANGPTL3(Angiopoietin Like Protein 3) ELISA Kit	ANGPTL3		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4713	Human Ang-III(Angiotensin III) ELISA Kit	Ang-III	P01019(35~	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4714	Human RARS(Arginyl tRNA Synthetase) ELISA Kit	RARS	P54136	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM4715	Mouse RARS(Arginyl tRNA Synthetase) ELISA Kit	RARS	Q9D0I9	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4716	Human NE-Ab(Neutrophil Elastase Antibody) ELISA Kit	NE-Ab		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH4717	Human IleRS-Ab(Isoleucyl-tRNA Synthetase Antibody) ELISA Kit	IleRS-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4718	Human P62-Ab(P62 Antibody) ELISA Kit	P62-Ab		Sandwich	8000ng/mL	75 ng/mL	125-8000ng/mL
FY-EH4719	Human PL-Ab(Phospholipid Antibody) ELISA Kit	PL-Ab		Sandwich	100µg/mL	0.94 µg/mL	1.57-100µg/mL
FY-EH4720	Human iPLA2(Phospholipase A2, Calcium Independent) ELISA Kit	iPLA2	O60733	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4721	Mouse iPLA2(Phospholipase A2, Calcium Independent) ELISA Kit	iPLA2	P97819	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4722	Human GUK1(Guanylate Kinase 1) ELISA Kit	GUK1	Q16774	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4723	Human HB(Hemoglobin) ELISA Kit	HB	P69905,P688	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4724	Mouse HB(Hemoglobin) ELISA Kit	HB		Sandwich	10000ng/mL	93.75 ng/mL	156.25-10000ng/mL
FY-EH4725	Human HRG(Histidine Rich Glycoprotein) ELISA Kit	HRG	P04196	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4726	Human HSPG1(Heparan Sulfate Proteoglycan 1) ELISA Kit	HSPG1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4727	Human IL-10Rβ(Interleukin 10 Receptor Beta) ELISA Kit	IL-10Rβ	Q08334	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4728	Rat IL-10Rβ(Interleukin 10 Receptor Beta) ELISA Kit	IL-10Rβ		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL

FY-EH4729	Human AT-III(Antithrombin III) ELISA Kit	AT-III		Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EM4730	Mouse AT-III(Antithrombin III) ELISA Kit	AT-III		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4731	Human TG-Ab(Thyroglobulin Antibody) ELISA Kit	TG-Ab		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM4732	Mouse TG-Ab(Thyroglobulin Antibody) ELISA Kit	TG-Ab		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH4733	Human IGFBP-1(Insulin-like Growth Factor Binding Protein 1) ELISA Kit	IGFBP-1	P08833	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4734	Mouse IGFBP-1(Insulin-like Growth Factor Binding Protein 1) ELISA Kit	IGFBP-1	P47876	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4735	Rat IGFBP-1(Insulin-like Growth Factor Binding Protein 1) ELISA Kit	IGFBP-1	P21743	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4736	Human IGFBP-5(Insulin-like Growth Factor Binding Protein 5) ELISA Kit	IGFBP-5	P24593	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4737	Mouse IGFBP-5(Insulin-like Growth Factor Binding Protein 5) ELISA Kit	IGFBP-5	Q07079	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-ER4738	Rat IGFBP-5(Insulin-like Growth Factor Binding Protein 5) ELISA Kit	IGFBP-5	P24594	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4739	Human IGFBP-4(Insulin-like Growth Factor Binding Protein 4) ELISA Kit	IGFBP-4	P22692	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM4740	Mouse IGFBP-4(Insulin-like Growth Factor Binding Protein 4) ELISA Kit	IGFBP-4	P47879	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER4741	Rat IGFBP-4(Insulin-like Growth Factor Binding Protein 4) ELISA Kit	IGFBP-4	P21744	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4742	Human IGFBP-6(Insulin-like Growth Factor Binding Protein 6) ELISA Kit	IGFBP-6	P24592	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4743	Mouse IGFBP-6(Insulin-like Growth Factor Binding Protein 6) ELISA Kit	IGFBP-6	P47880	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER4744	Rat IGFBP-6(Insulin-like Growth Factor Binding Protein 6) ELISA Kit	IGFBP-6	P35572	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4745	Human AP13(Apelin 13) ELISA Kit	AP13	Q9ULZ1(65-	Competitive	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-ER4746	Rat AP13(Apelin 13) ELISA Kit	AP13	Q9R0R3	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4747	Human AP36(Apelin 36) ELISA Kit	AP36	Q9ULZ1(42-	Competitive	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-ER4748	Rat AP36(Apelin 36) ELISA Kit	AP36		Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4749	Human ApoC3(Apolipoprotein C3) ELISA Kit	ApoC3	P02656	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EM4750	Mouse ApoC3(Apolipoprotein C3) ELISA Kit	ApoC3	P33622	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4751	Rat ApoC3(Apolipoprotein C3) ELISA Kit	ApoC3	P06759	Sandwich	400ng/ml	3.75 ng/ml	6.25-400ng/ml
FY-EH4752	Human ApoB100(Apolipoprotein B100) ELISA Kit	ApoB100	P04114	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-EM4753	Mouse ApoB100(Apolipoprotein B100) ELISA Kit	ApoB100		Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-ER4754	Rat ApoB100(Apolipoprotein B100) ELISA Kit	ApoB100	Q7TMA5	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-EH4755	Human ITIH4(Inter-alpha-trypsin Inhibitor Heavy Chain H4) ELISA Kit	ITIH4	Q14624	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4756	Human ARG1(Arginase 1) ELISA Kit	ARG1	P05089	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM4757	Mouse ARG1(Arginase 1) ELISA Kit	ARG1	Q61176	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER4758	Rat ARG1(Arginase 1) ELISA Kit	ARG1	P07824	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4759	Human ATXN2(Ataxin 2) ELISA Kit	ATXN2	Q99700	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM4760	Mouse ATXN2(Ataxin 2) ELISA Kit	ATXN2	O70305	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER4761	Rat ATXN2(Ataxin 2) ELISA Kit	ATXN2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4762	Human ATF5(Activating Transcription Factor 5) ELISA Kit	ATF5	Q9Y2D1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4763	Human ATF7(Activating Transcription Factor 7) ELISA Kit	ATF7	P17544	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4764	Human ABCG1(ATP Binding Cassette Transporter G1) ELISA Kit	ABCG1	P45844	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4765	Mouse ABCG1(ATP Binding Cassette Transporter G1) ELISA Kit	ABCG1	Q64343	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4766	Rat ABCG1(ATP Binding Cassette Transporter G1) ELISA Kit	ABCG1		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4767	Human ABCA1(ATP Binding Cassette Transporter A1) ELISA Kit	ABCA1	O95477	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4768	Mouse ABCA1(ATP Binding Cassette Transporter A1) ELISA Kit	ABCA1	P41233	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4769	Rat ABCA1(ATP Binding Cassette Transporter A1) ELISA Kit	ABCA1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH4770	Human ABCG2(ATP Binding Cassette Transporter G2) ELISA Kit	ABCG2	Q9UNQ0	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4771	Human Aβ1-40(Amyloid Beta 1-40) ELISA Kit	Aβ1-40	P05067	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4772	Mouse Aβ1-40(Amyloid Beta 1-40) ELISA Kit	Aβ1-40	P12023(672-	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER4773	Rat Aβ1-40(Amyloid Beta 1-40) ELISA Kit	Aβ1-40		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4774	Human Aβ1-42(Amyloid Beta 1-42) ELISA Kit	Aβ1-42	P05067	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4775	Mouse Aβ1-42(Amyloid Beta 1-42) ELISA Kit	Aβ1-42	P12023(672-	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-ER4776	Rat Aβ1-42(Amyloid Beta 1-42) ELISA Kit	Aβ1-42		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4777	Human Bcl-3(B-cell Lymphoma/Leukemia 3) ELIS Kit ELISA Kit	Bcl-3	P20749	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM4778	Mouse Bcl-3(B-cell Lymphoma/Leukemia 3) ELIS Kit ELISA Kit	Bcl-3	Q9Z2F6	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4779	Rat Bcl-3(B-cell Lymphoma/Leukemia 3) ELISA Kit	Bcl-3		Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-EH4780	Human BLNK(B-Cell Linker Protein) ELISA Kit	BLNK	Q8WV28	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4781	Human BTC(Betacellulin) ELISA Kit	BTC	P35070	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4782	Mouse BTC(Betacellulin) ELISA Kit	BTC	Q543J8	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER4783	Rat BTC(Betacellulin) ELISA Kit	BTC	Q9JJM4	Sandwich	600pg/mL	5.63 pg/mL	9.38-600pg/mL
FY-EH4784	Human BMPR1A(Bone Morphogenetic Protein Receptor 1A) ELISA Kit	BMPR1A	P36894	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4785	Mouse BMPR1A(Bone Morphogenetic Protein Receptor 1A) ELISA Kit	BMPR1A	P36895	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4786	Rat BMPR1A(Bone Morphogenetic Protein Receptor 1A) ELISA Kit	BMPR1A	Q78EA7	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4787	Human IBSP(Bone Sialoprotein 2) ELISA Kit	IBSP		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM4788	Mouse IBSP(Bone Sialoprotein 2) ELISA Kit	IBSP		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4789	Rat IBSP(Bone Sialoprotein 2) ELISA Kit	IBSP		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4790	Human BFP(Brain Finger Protein) ELISA Kit	BFP	Q9ULX5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4791	Mouse BFP(Brain Finger Protein) ELISA Kit	BFP	Q96DY5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4792	Human C4BP(C4 Binding Protein) ELISA Kit	C4BP	P04003(C4B	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4793	Human CGRP1(Calcitonin Gene Related Peptide 1) ELISA Kit	CGRP1	P06881	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4794	Mouse CGRP1(Calcitonin Gene Related Peptide 1) ELISA Kit	CGRP1	Q99JA0	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4795	Rat CGRP1(Calcitonin Gene Related Peptide 1) ELISA Kit	CGRP1		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4796	Human CaSR(Calcium Sensing Receptor) ELISA Kit	CaSR	P41180	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM4797	Mouse CaSR(Calcium Sensing Receptor) ELISA Kit	CaSR	Q9QY96	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4798	Rat CASR(Calcium Sensing Receptor) ELISA Kit	CASR	P48442	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4799	Human CAP1(Adenylate Cyclase Associated Protein 1) ELISA Kit	CAP1	Q01518	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4800	Human CA19-9(Carbohydrate antigen19-9) ELISA Kit	CA19-9		Sandwich	400IU/mL	3.75 IU/mL	6.25-400IU/mL
FY-EH4801	Human CMLC-1(Cardiac Myosin Light Chain 1) ELISA Kit	CMLC-1	P08590	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER4802	Rat CMLC-1(Cardiac Myosin Light Chain 1) ELISA Kit	CMLC-1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4803	Human TNNI3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit	TNNI3/cTn-I	P19429	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-EM4804	Mouse TNNI3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit	TNNI3/cTn-I	P48787	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER4805	Rat TNNI3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit	TNNI3/cTn-I	P23693	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4806	Human CT1(Cardiotrophin 1) ELISA Kit	CT1	Q16619	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4807	Mouse CT1(Cardiotrophln 1) ELISA Kit	CT1	Q60753	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4808	Rat CT1(Cardiotrophin 1) ELISA Kit	CT1	Q63086	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4809	Human CASP10(Caspase 10) ELISA Kit	CASP10	Q92851	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4810	Mouse CASP10(Caspase 10) ELISA Kit	CASP10		Sandwich	1500pg/mL	14.06 pg/mL	23.44-1500pg/mL

FY-EH4811	Human CASP12(Caspase 12) ELISA Kit	CASP12		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4812	Rat CASP12(Caspase 12) ELISA Kit	CASP12	Q920D5	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4813	Human Cav-1(Caveolin-1) ELISA Kit	Cav-1	Q03135	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4814	Mouse Cav-1(Caveolin-1) ELISA Kit	Cav-1	P49817	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ER4815	Rat CAV-1(Caveolin-1) ELISA Kit	CAV-1	P41350	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4816	Human CCR1(CC-Chemokine Receptor 1) ELISA Kit	CCR1	P32246	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4817	Rat CCR1(CC-Chemokine Receptor 1) ELISA Kit	CCR1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4818	Human CDC23(Cell Division Cycle Protein 23) ELISA Kit	CDC23	Q9UJX2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4819	Human CENPE(Centromere Protein E) ELISA Kit	CENPE	Q02224	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4820	Mouse CENPE(Centromere Protein E) ELISA Kit	CENPE	Q6RT24	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4821	Rat CENPE(Centromere Protein E) ELISA Kit	CENPE		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4822	Human CENPH(Centromere Protein H) ELISA Kit	CENPH	Q9H3R5	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4823	Mouse CENPH(Centromere Protein H) ELISA Kit	CENPH	Q9QYM8	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4824	Human CENPI(Centromere Protein I) ELISA Kit	CENPI	Q92674	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4825	Mouse CENPI(Centromere Protein I) ELISA Kit	CENPI	Q8K1K4	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4826	Rat CENPI(Centromere Protein I) ELISA Kit	CENPI		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4827	Human CENPB(Centromere Protein B) ELISA Kit	CENPB	P07199	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM4828	Mouse CENPB(Centromere Protein B) ELISA Kit	CENPB	P27790	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4829	Rat CENPB(Centromere Protein B) ELISA Kit	CENPB		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4830	Human CEP110(Centrosomal Protein 110kDa) ELISA Kit	CEP110	Q7Z7A1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM4831	Mouse CEP110(Centrosomal Protein 110kDa) ELISA Kit	CEP110	Q7TSH4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4832	Rat CEP110(Centrosomal Protein 110kDa) ELISA Kit	CEP110		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4833	Human CPN10(Chaperonin 10) ELISA Kit	CPN10	P61604	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH4834	Human CCL1(Chemokine C-C-Motif Ligand 1) ELISA Kit	CCL1	P22362	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4835	Mouse CCL1(Chemokine C-C-Motif Ligand 1) ELISA Kit	CCL1	P10146	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4836	Human CCL4L1(Chemokine C-C-Motif Ligand 4 Like Protein 1) ELISA Kit	CCL4L1	Q8NHW4	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4837	Human CCRL1(Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit	CCRL1	Q9NBP9	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4838	Mouse CCRL1(Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit	CCRL1	Q924I3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4839	Rat CCRL1(Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit	CCRL1		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4840	Human CX3CR1(Chemokine C-X3-C-Motif Receptor 1) ELISA Kit	CX3CR1	P49238	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4841	Mouse CX3CR1(Chemokine C-X3-C-Motif Receptor 1) ELISA Kit	CX3CR1	Q9Z0D9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4842	Rat CX3CR1(Chemokine C-X3-C-Motif Receptor 1) ELISA Kit	CX3CR1	P35411	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4843	Human CKLFSF7(Chemokine Like Factor Superfamily 7) ELISA Kit	CKLFSF7	Q96FZ5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4844	Mouse CKLFSF7(Chemokine Like Factor Superfamily 7) ELISA Kit	CKLFSF7	Q9ESD6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4845	Rat CKLFSF7(Chemokine Like Factor Superfamily 7) ELISA Kit	CKLFSF7		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4846	Human CKLFSF8(Chemokine Like Factor Superfamily 8) ELISA Kit	CKLFSF8	Q8IZV2	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4847	Mouse CKLFSF8(Chemokine Like Factor Superfamily 8) ELISA Kit	CKLFSF8	Q9CZR4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4848	Human CCKBR(Cholecystokinin B Receptor) ELISA Kit	CCKBR	P32239	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4849	Mouse CCKBR(Cholecystokinin B Receptor) ELISA Kit	CCKBR	P56481	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4850	Human CHAc(Choline acetylase) ELISA Kit	CHAc	P28329	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4851	Human PC/CPG(Choline Phosphoglyceride) ELISA Kit	PC/CPG		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL

FY-EM4852	Mouse PC/CPG(choline phosphoglyceride) ELISA Kit	PC/CPG		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4853	Human CHD3(Chromodomain Helicase DNA Binding Protein 3) ELISA Kit	CHD3	Q12873	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4854	Mouse CHD3(Chromodomain Helicase DNA Binding Protein 3) ELISA Kit	CHD3	B1AR17	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4855	Rat CHD3(Chromodomain Helicase DNA Binding Protein 3) ELISA Kit	CHD3		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4856	Human CHGA(Chromogranin A) ELISA Kit	CHGA	P10645	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM4857	Mouse CHGA(Chromogranin A) ELISA Kit	CHGA	P26339	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4858	Rat CHGA(Chromogranin A) ELISA Kit	CHGA	P10354	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4859	Human CLDN14(Claudin 14) ELISA Kit	CLDN14	O95500	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4860	Human CD109(Cluster of Differentiation 109) ELISA Kit	CD109	Q6YHK3	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4861	Mouse CD109(Cluster of Differentiation 109) ELISA Kit	CD109	Q8R422	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4862	Rat CD109(Cluster Of Differentiation 109) ELISA Kit	CD109		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4863	Mouse GRN(Granulin) ELISA Kit	GRN		Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-EH4864	Human CCSA-2(Colon Cancer Specific Antigen-2) ELISA Kit	CCSA-2		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4865	Human CCSA-3(Colon Cancer Specific Antigen-3) ELISA Kit	CCSA-3		Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH4866	Human CSF2R α (Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit	CSF2R α	P15509	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4867	Human P II NP(Procollagen II N-Terminal ProPeptide) ELISA Kit	P II NP	P02458(COI	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4868	Rat P II NP(Procollagen II N-Terminal Propeptide) ELISA Kit	P II NP		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4869	Human GCSFR(Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit	GCSFR	Q99062	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4870	Mouse GCSFR(Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit	GCSFR	P40223	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4871	Rat GCSFR(Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit	GCSFR		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4872	Human C1qR1(Complement Component 1q Receptor 1) ELISA Kit	C1qR1	Q9NPY3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4873	Mouse C1qR1(Complement Component 1q Receptor 1) ELISA Kit	C1qR1	O89103	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4874	Human C1r(Complement Component 1R) ELISA Kit	C1r	P00736	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4875	Rat C1R(Complement Component 1Q Receptor) ELISA Kit	C1R		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4876	Human C5(Complement Component 5) ELISA Kit	C5	P01031	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM4877	Mouse C5(Complement Component 5) ELISA Kit	C5	P06684	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4878	Rat C5(Complement Component 5) ELISA Kit	C5	P08650	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4879	Human C8b(Complement Component 8b) ELISA Kit	C8b	P07358	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4880	Human C8g(Complement Component 8g) ELISA Kit	C8g	P07360	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4881	Human CTAPIII(Connective Tissue Activating Peptide III) ELISA Kit	CTAPIII	P02775	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EH4882	Human CX31(Connexin 31) ELISA Kit	CX31	O75712	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4883	Mouse CX31(Connexin 31) ELISA Kit	CX31	P28231	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4884	Rat CX31(Connexin 31) ELISA Kit	CX31		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH4885	Human HELIX- II (Type II Collagen Helical Peptide) ELISA Kit	HELIX- II		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4886	Human MTHFR(Methylenetetrahydrofolate Reductase) ELISA Kit	MTHFR	P42898	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4887	Human PRKA α 1(Protein Kinase, AMP Activated Alpha 1) ELISA Kit	PRKA α 1	Q13131	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4888	Human IKBKB(Inhibitor Of Kappa-Light Polypeptide Gene Enhancer In B-Cells Kinase) ELISA Kit	IKBKB	O14920	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4889	Human PTPRB(Protein Tyrosine Phosphatase Receptor Type B) ELISA Kit	PTPRB	P23467	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4890	Human CREB1(Cyclic AMP Response Element Binding Protein 1) ELISA Kit	CREB1	P16220	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4891	Mouse CREB1(Cyclic AMP Response Element Binding Protein 1) ELISA Kit	CREB1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4892	Rat CREB1(Cyclic AMP Response Element Binding Protein 1) ELISA Kit	CREB1		Sandwich	20pmol/mL	0.19 pmol/mL	0.32-20pmol/mL

FY-EH4893	Human CXCR3(CXC-Chemokine Receptor 3) ELISA Kit	CXCR3	P49682	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4894	Mouse CXCR3(CXC-Chemokine Receptor 3) ELISA Kit	CXCR3	O88410	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4895	Rat CXCR3(CXC-Chemokine Receptor 3) ELISA Kit	CXCR3	Q9JIJ9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4896	Human CLEC2C(C-Type Lectin Domain Family 2, Member C) ELISA Kit	CLEC2C	Q07108	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM4897	Mouse CLEC2C(C-Type Lectin Domain Family 2, Member C) ELISA Kit	CLEC2C		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER4898	Rat CLEC2C(C-Type Lectin Domain Family 2, Member C) ELISA Kit	CLEC2C		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4899	Human SELE(E-Selectin) ELISA Kit	SELE	P16581	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM4900	Mouse SELE(E-selectin) ELISA Kit	SELE	Q00690	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4901	Rat SELE(E-selectin) ELISA Kit	SELE	P98105	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH4902	Human FBXO2(F-Box Protein 2) ELISA Kit	FBXO2	Q9UK22	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4903	Mouse FBXO2(F-Box Protein 2) ELISA Kit	FBXO2	Q80UW2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4904	Human Flt3(FMS-like Tyrosine Kinase 3) ELISA Kit	Flt3	P36888	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4905	Mouse Flt3(FMS-like Tyrosine Kinase 3) ELISA Kit	Flt3	Q00342	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4906	Human GATA1(GATA Binding Protein 1) ELISA Kit	GATA1	P15976	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM4907	Mouse GATA1(GATA Binding Protein 1) ELISA Kit	GATA1	P17679	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4908	Rat GATA1(GATA Binding Protein 1) ELISA Kit	GATA1	P43429	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4909	Human GATA2(GATA Binding Protein 2) ELISA Kit	GATA2	P23769	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4910	Mouse GATA2(GATA Binding Protein 2) ELISA Kit	GATA2	O09100	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4911	Rat GATA2(GATA Binding Protein 2) ELISA Kit	GATA2	Q924Y4	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4912	Human GATA5(GATA Binding Protein 5) ELISA Kit	GATA5	Q9BWX5	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4913	Mouse GATA5(GATA Binding Protein 5) ELISA Kit	GATA5	P97489	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4914	Rat GATA5(GATA Binding Protein 5) ELISA Kit	GATA5	Q5U2V0	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4915	Human GDI1(GDP Dissociation Inhibitor 1) ELISA Kit	GDI1	P31150	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4916	Human GAP(Ras Gtpase Activating Protein) ELISA Kit	GAP	P20936	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4917	Rat GAP(Ras Gtpase Activating Protein) ELISA Kit	GAP		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4918	Human RGS4(Regulator of G Protein Signaling 4) ELISA Kit	RGS4	P49798	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4919	Mouse RGS4(Regulator Of G Protein Signaling 4) ELISA Kit	RGS4	O08899	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4920	Rat RGS4(Regulator Of G Protein Signaling 4) ELISA Kit	RGS4		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4921	Human NOD9(Nucleotide-binding oligomerization domain protein 9) ELISA Kit	NOD9		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4922	Human S100A1(S100 Calcium Binding Protein A1) ELISA Kit	S100A1	P23297	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4923	Mouse S100A1(S100 Calcium Binding Protein A1) ELISA Kit	S100A1	P56565	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4924	Rat S100A1(S100 Calcium Binding Protein A1) ELISA Kit	S100A1	P04631	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4925	Human SHH(Sonic Hedgehog Protein) ELISA Kit	SHH	Q15465	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM4926	Mouse SHH(Sonic Hedgehog Protein) ELISA Kit	SHH	Q62226	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4927	Rat SHH(Sonic Hedgehog Protein) ELISA Kit	SHH	Q63673	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4928	Human TDP43(TAR DNA Binding Protein 43) ELISA Kit	TDP43	Q13148	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4929	Mouse TDP43(TAR DNA Binding Protein 43) ELISA Kit	TDP43	Q921F2	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4930	Rat TDP43(TAR DNA Binding Protein 43) ELISA Kit	TDP43		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4931	Human TAF1(Transcription initiation factor TFIID subunit 1) ELISA Kit	TAF1		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4932	Mouse TAF1(Transcription initiation factor TFIID subunit 1) ELISA Kit	TAF1		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4933	Human TAF2(Transcription initiation factor TFIID subunit 2) ELISA Kit	TAF2		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-EM4934	Mouse TAF2(Transcription initiation factor TFIID subunit 2) ELISA Kit	TAF2		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4935	Human TAF13(Transcription initiation factor TFIID subunit 13) ELISA Kit	TAF13	Q15543	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4936	Human TLR-2(Toll-like Receptor 2) ELISA Kit	TLR-2	O60603	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4937	Mouse TLR-2(Toll-like Receptor 2) ELISA Kit	TLR-2	Q9QUN7	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4938	Rat TLR-2(Toll-like Receptor 2) ELISA Kit	TLR-2	E9PTD9	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4939	Human TU3A(TU3A Protein) ELISA Kit	TU3A	O95990	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4940	Human LAT(Linker for activation of T-cell) ELISA Kit	LAT	O43561	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4941	Mouse LAT(Linker for activation of T-cell) ELISA Kit	LAT	O54957	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4942	Human TALLA-1/CD231(T-cell Acute Lymphoblastic Leukemia Antigen) ELISA Kit	TALLA-1/CD23	P41732	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4943	Human OCT1(Octamer-binding protein 1) ELISA Kit	OCT1	P14859	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4944	Mouse OCT1(Octamer-binding protein 1) ELISA Kit	OCT1	P25425	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4945	Human IL-1R1(Interleukin 1 Receptor Type I) ELISA Kit	IL-1R1	P14778	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4946	Mouse IL-1R1(Interleukin 1 Receptor Type I) ELISA Kit	IL-1R1	P13504	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-ER4947	Rat IL-1R1(Interleukin 1 Receptor Type I) ELISA Kit	IL-1R1		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH4948	Human LCA/CD45(Leukocyte Common Antigen) ELISA Kit	LCA/CD45	P08575	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4949	Human NSE(Neuron Specific Enolase) ELISA Kit	NSE	P09104	Sandwich	5000pg/mL	35.16 pg/mL	78.13-5000pg/mL
FY-EM4950	Mouse NSE(Neuron Specific Enolase) ELISA Kit	NSE	P17183	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4951	Rat NSE(Neuron Specific Enolase) ELISA Kit	NSE	P07323	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4952	Human cPLA2(Cytosolic Phospholipase A2) ELISA Kit	cPLA2	P47712	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4953	Mouse cPLA2(Cytosolic Phospholipase A2) ELISA Kit	cPLA2	P47713	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4954	Rat cPLA2(Cytosolic Phospholipase A2) ELISA Kit	cPLA2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4955	Human PDHB(Pyruvate dehydrogenase E1 component subunit Beta) ELISA Kit	PDHB		Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH4956	Human MIgA(Membrane Immunoglobulin A) ELISA Kit	MIgA		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EM4957	Mouse MIgA(Membrane Immunoglobulin A) ELISA Kit	MIgA		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4958	Human MIgD(Membrane Immunoglobulin D) ELISA Kit	MIgD		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH4959	Human MIgM(Membrane Immunoglobulin M) ELISA Kit	MIgM		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4960	Human PDHA1(Pyruvate dehydrogenase E1 component subunit alpha) ELISA Kit	PDHA1		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4961	Mouse PDHA1(Pyruvate dehydrogenase E1 component subunit alpha) ELISA Kit	PDHA1		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER4962	Rat PDHA1(Pyruvate dehydrogenase E1 component subunit alpha) ELISA Kit	PDHA1		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4963	Human VIM(Vimentin) ELISA Kit	VIM	P08670	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH4964	Human FOXP1(Forkhead Box Protein P1) ELISA Kit	FOXP1	Q9H334	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4965	Human FGFBP2(Fibroblast Growth Factor Binding Protein 2) ELISA Kit	FGFBP2	Q9BYJ0	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH4966	Human FRS2(Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit	FRS2	Q8WU20	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM4967	Mouse FRS2(Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit	FRS2	Q8C180	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH4968	Human VLA4(Very Late Appearing Antigen 4) ELISA Kit	VLA4	P13612(ITG.	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4969	Rat VLA4(Very Late Appearing Antigen 4) ELISA Kit	VLA4		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH4970	Human PRF1(Perforin 1) ELISA Kit	PRF1	P14222	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM4971	Mouse PRF1(Perforin 1) ELISA Kit	PRF1	P10820	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER4972	Rat PRF1(Perforin 1) ELISA Kit	PRF1	P35763	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4973	Human FSH(Follicle Stimulating Hormone) ELISA Kit	FSH	P01225	Sandwich	10mIU/mL	0.09 mIU/mL	0.16-10mIU/mL
FY-EM4974	Mouse FSH(Follicle Stimulating Hormone) ELISA Kit	FSH		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL

FY-ER4975	Rat FSH(Follicle Stimulating Hormone) ELISA Kit	FSH		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH4976	Human CRH(Corticotropin Releasing Hormone) ELISA Kit	CRH	P06850	Competitive	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-EM4977	Mouse CRH(Corticotropin Releasing Hormone) ELISA Kit	CRH	Q8CIT0	Competitive	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER4978	Rat CRH(Corticotropin Releasing Hormone) ELISA Kit	CRH	P01143	Competitive	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4979	Human PYGB(Glycogen Phosphorylase, Brain) ELISA Kit	PYGB	P11216	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM4980	Mouse PYGB(Glycogen Phosphorylase, Brain) ELISA Kit	PYGB	Q8CI94	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER4981	Rat PYGB(Glycogen Phosphorylase, Brain) ELISA Kit	PYGB		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH4982	Human MAOA(Type A Monoamine Oxidase) ELISA Kit	MAOA	P21397	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4983	Mouse MAOA(Monoamine Oxidase Type A) ELISA Kit	MAOA	Q64133	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER4984	Rat MAOA(Monoamine Oxidase Type A) ELISA Kit	MAOA	P21396	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4985	Human MCP-2/CCL8(Monocyte Chemotactic Protein 2) ELISA Kit	MCP-2/CCL8	P80075	Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-EM4986	Mouse MCP-2/CCL8(Monocyte Chemotactic Protein 2) ELISA Kit	MCP-2/CCL8	Q9Z121	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH4987	Human SREBP(Sterol Regulatory Element Binding Protein) ELISA Kit	SREBP	P36956(SRE	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM4988	Mouse SREBP(Sterol Regulatory Element Binding Protein) ELISA Kit	SREBP	Q9WTN3	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER4989	Rat SREBP(Sterol Regulatory Element Binding Protein) ELISA Kit	SREBP		Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-EH4990	Human EMILIN1(Elastin Microfibril Interface Located Protein 1) ELISA Kit	EMILIN1	Q9Y6C2	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM4991	Mouse EMILIN1(Elastin Microfibril Interface Located Protein 1) ELISA Kit	EMILIN1	Q99K41	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER4992	Rat EMILIN1(Elastin Microfibril Interface Located Protein 1) ELISA Kit	EMILIN1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH4993	Human PKB(Protein Kinase B) ELISA Kit	PKB	P31749	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM4994	Mouse PKB(Protein Kinase B) ELISA Kit	PKB	P31750	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER4995	Rat PKB(Protein Kinase B) ELISA Kit	PKB		Sandwich	200IU/mL	1.88 IU/mL	3.13-200IU/mL
FY-EH4996	Human PKD2(Protein Kinase D2) ELISA Kit	PKD2	Q9BZL6	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM4997	Mouse PKD2(Protein Kinase D2) ELISA Kit	PKD2	Q8BZ03	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH4998	Human LMP7/PSMB8(Low molecular mass protein 7/Proteasome subunit beta type- ξ LMP7/PSMB8			Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH4999	Human VCAN(Versican/PG-M) ELISA Kit	VCAN	P13611	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-EH5000	Human PYY(Peptide YY) ELISA Kit	PYY	P10082	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5001	Mouse PYY(Peptide YY) ELISA Kit	PYY	Q9EPS2	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5002	Rat PYY(Peptide YY) ELISA Kit	PYY	P10631	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5003	Human DPYSL3(Dihydropyrimidinase Like Protein 3) ELISA Kit	DPYSL3		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5004	Human FDFT1(Farnesyl Diphosphate Farnesyltransferase 1) ELISA Kit	FDFT1	P37268	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH5005	Human SP-A1(Pulmonary surfactant-associated protein A1) ELISA Kit	SP-A1	Q8IWL2	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5006	Mouse SP-A1(Pulmonary surfactant-associated protein A1) ELISA Kit	SP-A1	P35242	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5007	Rat SP-A1(Pulmonary surfactant-associated protein A1) ELISA Kit	SP-A1		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5008	Human SP-B(Pulmonary surfactant-associated protein B) ELISA Kit	SP-B	P07988	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5009	Rat SP-B(Pulmonary surfactant-associated protein B) ELISA Kit	SP-B	P22355	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH5010	Human SP-D(Pulmonary surfactant-associated protein D) ELISA Kit	SP-D	P35247	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5011	Mouse SP-D(Pulmonary surfactant-associated protein D) ELISA Kit	SP-D	P50404	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5012	Rat SP-D(Pulmonary surfactant-associated protein D) ELISA Kit	SP-D	P35248	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH5013	Human PRCP(Prolylcarboxypeptidase) ELISA Kit	PRCP	P42785	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5014	Rat PRCP(Prolylcarboxypeptidase) ELISA Kit	PRCP		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5015	Human RPA1(Replication protein A1) ELISA Kit	RPA1	P27694	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-EH5016	Human MBP-C/MBL2(Mannose Binding Protein C/Mannose Binding Lectin) ELISA	MBP-C/MBL2	P11226	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5017	Rat MBP-C/MBL2(Mannose Binding Protein C/Mannose Binding Lectin) ELISA Kit	MBP-C/MBL2		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5018	Human RNASEH1(Ribonuclease H1) ELISA Kit	RNASEH1	O60930	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER5019	Rat RNASEH1(Ribonuclease H1) ELISA Kit	RNASEH1		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH5020	Human RNASET2(Ribonuclease T2) ELISA Kit	RNASET2	O00584	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5021	Mouse RNASET2(Ribonuclease T2) ELISA Kit	RNASET2	C0HKG5(Rn	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5022	Rat RNASET2(Ribonuclease T2) ELISA Kit	RNASET2		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5023	Human NFKB-p105(Nuclear factor NF-kappa-B p105 subunit) ELISA Kit	NFKB-p105		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5024	Mouse NFKB-p105(Nuclear factor NF-kappa-B p105 subunit) ELISA Kit	NFKB-p105	P25799	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5025	Rat NFKB-p105(Nuclear factor NF-kappa-B p105 subunit) ELISA Kit	NFKB-p105	Q63369	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5026	Human NFKB-p65(Nuclear factor NF-kappa-B p65 subunit) ELISA Kit	NFKB-p65	Q04206	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5027	Mouse NFKB-p65(Nuclear Factor Kappa B p65) ELISA Kit	NFKB-p65	Q04207	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5028	Rat NFKB-p65(Nuclear Factor Kappa B p65) ELISA Kit	NFKB-p65		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH5029	Human β MSH(beta-Melanocyte Stimulating Hormone) ELISA Kit	β MSH		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5030	Mouse β MSH(beta-Melanocyte Stimulating Hormone) ELISA Kit	β MSH	P01193(185-	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5031	Rat β MSH(beta-Melanocyte Stimulating Hormone) ELISA Kit	β MSH		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH5032	Human GRIK2(Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit	GRIK2	Q13002	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5033	Rat GRIK2(Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit	GRIK2		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH5034	Human MUSK-Ab(Muscle Skeletal Receptor Tyrosine Kinase Antibody) ELISA Kit	MUSK-Ab		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH5035	Human CKMB(Creatine Kinase MB Isoenzyme) ELISA Kit	CKMB	P06732	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EM5036	Mouse CKMB(Creatine Kinase MB Isoenzyme) ELISA Kit	CKMB		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5037	Rat CKMB(Creatine Kinase MB Isoenzyme) ELISA Kit	CKMB		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5038	Human T4-Ab(Thyroxine Antibody) ELISA Kit	T4-Ab		Sandwich	3000ng/mL	28.13 ng/mL	46.88-3000ng/mL
FY-EM5039	Mouse T4-Ab(Thyroxine Antibody) ELISA Kit	T4-Ab		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER5040	Rat T4-Ab(Thyroxine Antibody) ELISA Kit	T4-Ab		Sandwich	100pmol/mL	0.94 pmol/mL	1.57-100pmol/mL
FY-EH5041	Human CNTN4(Contactin 4) ELISA Kit	CNTN4	Q81WV2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5042	Mouse CNTN4(Contactin 4) ELISA Kit	CNTN4	Q69Z26	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5043	Rat CNTN4(Contactin 4) ELISA Kit	CNTN4		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5044	Human UPP2(Uridine Phosphorylase 2) ELISA Kit	UPP2	O95045	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5045	Human MT1(Metallothionein 1) ELISA Kit	MT1	P04731	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5046	Mouse MT1(Metallothionein 1) ELISA Kit	MT1	P02802	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5047	Rat MT1(Metallothionein 1) ELISA Kit	MT1	P02803	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5048	Human SSFA2(Sperm Specific Antigen 2) ELISA Kit	SSFA2	P28290	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5049	Mouse SSFA2(Sperm Specific Antigen 2) ELISA Kit	SSFA2	Q922B9	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5050	Rat SSFA2(Sperm Specific Antigen 2) ELISA Kit	SSFA2		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH5051	Human PD-1(Programmed Cell Death Protein 1) ELISA Kit	PD-1	Q15116	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5052	Human PD-L1(Programmed Cell Death Protein 1 Ligand 1) ELISA Kit	PD-L1	Q9NZQ7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5053	Human GP210-Ab(Glucoprotein 210 Antibody) ELISA Kit	GP210-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5054	Mouse GP210-Ab(Glucoprotein 210 Antibody) ELISA Kit	GP210-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH5055	Human τ PK1(Tau-Protein Kinase 1) ELISA Kit	τ PK1	Q5TCY1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5056	Human TRACP-5b(Tartrate Resistant Acid Phosphatase 5b) ELISA Kit	TRACP-5b		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL

FY-EM5057	Mouse TRACP-5b(Tartrate Resistant Acid Phosphatase 5b) ELISA Kit	TRACP-5b	Q05117	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5058	Rat TRACP-5b(Tartrate Resistant Acid Phosphatase 5b) ELISA Kit	TRACP-5b		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH5059	Human HMGB-1(High Mobility Group Protein B1) ELISA Kit	HMGB-1	P09429	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5060	Mouse HMGB-1(High Mobility Group Protein B1) ELISA Kit	HMGB-1	P63158	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5061	Rat HMGB-1(High Mobility Group Protein B1) ELISA Kit	HMGB-1	P63159	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5062	Human GH-Ab(Growth Hormone Antibody) ELISA Kit	GH-Ab		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5063	Rat GH-Ab(Growth Hormone Antibody) ELISA Kit	GH-Ab		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH5064	Human GP-Ab(Glycoprotein Antibody) ELISA Kit	GP-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH5065	Human NFE2L2(Nuclear Factor, Erythroid Derived 2, Like 2) ELISA Kit	NFE2L2	Q16236	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5066	Mouse NFE2L2(Nuclear Factor, Erythroid Derived 2, Like 2) ELISA Kit	NFE2L2	Q60795	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-ER5067	Rat NFE2L2(Nuclear Factor, Erythroid Derived 2, Like 2) ELISA Kit	NFE2L2	P57790	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5068	Human DNP-Ab(Deoxyribonucleic Protein Antibody) ELISA Kit	DNP-Ab		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5069	Human PDGF-BB(Platelet Derived Growth Factor BB) ELISA Kit	PDGF-BB	P01127(PDC	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5070	Mouse PDGF-BB(Platelet Derived Growth Factor BB) ELISA Kit	PDGF-BB		Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ER5071	Rat PDGF-BB(Platelet Derived Growth Factor BB) ELISA Kit	PDGF-BB		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5072	Human RBP4(Retinol Binding Protein 4) ELISA Kit	RBP4	P02753	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5073	Mouse RBP4(Retinol Binding Protein 4) ELISA Kit	RBP4	Q00724	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER5074	Rat RBP4(Retinol Binding Protein 4) ELISA Kit	RBP4	P04916	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH5075	Human TGF-β2(Transforming Growth Factor Beta 2) ELISA Kit	TGF-β2	P61812	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5076	Mouse TGF-β2(Transforming Growth Factor Beta 2) ELISA Kit	TGF-β2	P27090(303-	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5077	Rat TGF-β2(Transforming Growth Factor Beta 2) ELISA Kit	TGF-β2		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5078	Human TAP1(Transporter associated with Antigen Processing 1) ELISA Kit	TAP1	Q03518	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5079	Human CD28(Cluster of Differentiation 28) ELISA Kit	CD28	P10747	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER5080	Rat CD28(Cluster of Differentiation 28) ELISA Kit	CD28		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5081	Human CD30L(Cluster of Differentiation 30 Ligand) ELISA Kit	CD30L	P32971	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM5082	Mouse CD30L(Cluster of Differentiation 30 Ligand) ELISA Kit	CD30L	P32972	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5083	Human CD86(T-Lymphocyte Activation Antigen CD86) ELISA Kit	CD86	P42081	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5084	Human sTNF(Soluble Tumor Necrosis Factor) ELISA Kit	sTNF		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5085	Human MHCG/HLA-G(Major Histocompatibility Complex Class I G) ELISA Kit	MHCG/HLA-G	P17693	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH5086	Human PDGF-A(Platelet Derived Growth Factor Subunit A) ELISA Kit	PDGF-A	P04085	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5087	Mouse PDGF-A(Platelet Derived Growth Factor Subunit A) ELISA Kit	PDGF-A	P20033	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5088	Rat PDGF-A(Platelet Derived Growth Factor Subunit A) ELISA Kit	PDGF-A		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5089	Human PDGF-B(Platelet Derived Growth Factor Subunit B) ELISA Kit	PDGF-B	P01127	Sandwich	1500pg/mL	14.06 pg/mL	23.44-1500pg/mL
FY-EM5090	Mouse PDGF-B(Platelet Derived Growth Factor Subunit B) ELISA Kit	PDGF-B	P31240	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5091	Rat PDGF-B(Platelet Derived Growth Factor Subunit B) ELISA Kit	PDGF-B	Q05028	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5092	Human CX37(Connexin 37) ELISA Kit	CX37	P35212	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5093	Mouse CX37(Connexin 37) ELISA Kit	CX37	P28235	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5094	Rat CX37(Connexin 37) ELISA Kit	CX37		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5095	Human LFA-3/CD58(Lymphocyte Function Associated Antigen 3) ELISA Kit	LFA-3/CD58	P19256	Sandwich	800pg/mL	7.5 pg/mL	12.5-800pg/mL
FY-EH5096	Human pERK1/2(Phospho Extracellular Signal Regulated Kinase 1/2) ELISA Kit	pERK1/2		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5097	Mouse pERK1/2(Phospho Extracellular Signal Regulated Kinase 1/2) ELISA Kit	pERK1/2		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL

FY-ER5098	Rat pERK1/2(Phospho Extracellular Signal Regulated Kinase 1/2) ELISA Kit	pERK1/2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5099	Human PLA2P(Phospholipase A2 Activating Protein) ELISA Kit	PLA2P	Q9Y263	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5100	Mouse PLA2P(Phospholipase A2 Activating Protein) ELISA Kit	PLA2P	P27612	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5101	Rat PLA2P(Phospholipase A2 Activating Protein) ELISA Kit	PLA2P		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5102	Human PLA1(Phospholipase A1) ELISA Kit	PLA1	Q53H76	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5103	Mouse PLA1(Phospholipase A1) ELISA Kit	PLA1		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER5104	Rat PLA1(Phospholipase A1) ELISA Kit	PLA1		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH5105	Human PLCβ1/PI-PLC(Phospholipase C Beta 1, Phosphoinositide Specific) ELISA K	PLCβ1/PI-PLC	Q9NQ66	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH5106	Human MS(Melatonin Sulfate) ELISA Kit	MS		Competitive	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-EH5107	Human EM2(Endomorphin 2) ELISA Kit	EM2	P35372	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5108	Human PGDS2(Prostaglandin D2 Synthase) ELISA Kit	PGDS2	P41222	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH5109	Human PIBF(Progesterone Induced Blocking Factor) ELISA Kit	PIBF	Q8WXW3	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5110	Mouse PIBF(Progesterone Induced Blocking Factor) ELISA Kit	PIBF	E9Q6K3	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH5111	Human DEFα5(Defensin Alpha 5, Paneth Cell Specific) ELISA Kit	DEFα5	Q01523	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5112	Human IgM(Immunoglobulin M) ELISA Kit	IgM	P01871	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5113	Mouse IgM(Immunoglobulin M) ELISA Kit	IgM	P01872	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER5114	Rat IgM(Immunoglobulin M) ELISA Kit	IgM		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH5115	Human hs-Cys-C(high-sensitivity Cystatin C) ELISA Kit	hs-Cys-C	P01034	Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-EH5116	Human EP2(Prostaglandin E Receptor 2) ELISA Kit	EP2	P43116	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5117	Mouse EP2(Prostaglandin E Receptor 2) ELISA Kit	EP2	Q62053	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5118	Rat EP2(Prostaglandin E Receptor 2) ELISA Kit	EP2		Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-EH5119	Human ACP3/PAP(Prostatic Acid Phosphatase) ELISA Kit	ACP3/PAP	P15309	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM5120	Mouse ACP3/PAP(Prostatic Acid Phosphatase) ELISA Kit	ACP3/PAP	Q8CE08	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5121	Rat ACP3/PAP(Prostatic Acid Phosphatase) ELISA Kit	ACP3/PAP		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EH5122	Human DSG1(Desmogleins 1) ELISA Kit	DSG1	Q02413	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH5123	Human HJV(Hemojuvelin) ELISA Kit	HJV	Q6ZVN8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5124	Mouse HJV(Hemojuvelin) ELISA Kit	HJV	Q7TQ32	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5125	Rat HJV(Hemojuvelin) ELISA Kit	HJV	Q8N7M5	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH5126	Human HSP-70/HSPA9(Heat Shock 70 kDa Protein 9) ELISA Kit	HSP-70/HSPA9	P38646	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5127	Mouse HSP-70/HSPA9(Heat Shock Protein 70) ELISA Kit	HSP-70/HSPA9	Q61696&P1	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ER5128	Rat HSP-70/HSPA9(Heat Shock Protein 70) ELISA Kit	HSP-70/HSPA9	P48721	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5129	Human HSF1(Heat Shock Factor 1) ELISA Kit	HSF1	Q00613	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5130	Mouse HSF1(Heat Shock Factor 1) ELISA Kit	HSF1	P38532	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5131	Rat HSF1(Heat Shock Factor 1) ELISA Kit	HSF1	P0DMW0,	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5132	Human WARS(Tryptophanyl tRNA Synthetase) ELISA Kit	WARS	P23381	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH5133	Human NT(Neurotensin) ELISA Kit	NT	P30990	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH5134	Human NT-3(Neurotrophin 3) ELISA Kit	NT-3	P20783	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5135	Mouse NT-3(Neurotrophin 3) ELISA Kit	NT-3	P20181	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5136	Rat NT-3(Neurotrophin 3) ELISA Kit	NT-3	P18280	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5137	Human GROβ/CXCL2(Growth Regulated Oncogene Beta) ELISA Kit	GROβ/CXCL2	P19875	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5138	Mouse GROβ/CXCL2(Growth Regulated Oncogene Beta) ELISA Kit	GROβ/CXCL2	P10889	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-ER5139	Rat GROβ/CXCL2(Growth Regulated Oncogene Beta) ELISA Kit	GROβ/CXCL2	P30348	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5140	Human NE/ELA2(Neutrophil Elastase/Elastase-2) ELISA Kit	NE/ELA2	P08246	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5141	Mouse NE/ELA2(Neutrophil Elastase/Elastase-2) ELISA Kit	NE/ELA2	Q3UP87	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5142	Rat NE/ELA2(Neutrophil Elastase/Elastase-2) ELISA Kit	NE/ELA2		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH5143	Human STK39(Serine/Threonine Protein Kinase 39) ELISA Kit	STK39	Q9UEW8	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5144	Mouse STK39(Serine/Threonine Protein Kinase 39) ELISA Kit	STK39	Q9Z1W9	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5145	Rat STK39(Serine/Threonine Protein Kinase 39) ELISA Kit	STK39		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH5146	Human DAP6(Death Associated Protein 6) ELISA Kit	DAP6	Q9UER7	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5147	Mouse DAP6(Death Associated Protein 6) ELISA Kit	DAP6	O35613	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5148	Rat DAP6(Death Associated Protein 6) ELISA Kit	DAP6		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH5149	Human MPO(Myeloperoxidase) ELISA Kit	MPO		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5150	Mouse MPO(Myeloperoxidase) ELISA Kit	MPO	62.5-4000 pg	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000 pg/mL
FY-EH5151	Human OB(Obestatin) ELISA Kit	OB	Q9UBU3(76	Competitive	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5152	Rat OB(Obestatin) ELISA Kit	OB		Competitive	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5153	Human TPPP(Tubulin Polymerization Promoting Protein) ELISA Kit	TPPP	O94811	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM5154	Mouse TPPP(Tubulin Polymerization Promoting Protein) ELISA Kit	TPPP	Q7TQD2	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5155	Rat TPPP(Tubulin Polymerization Promoting Protein) ELISA Kit	TPPP	A0A0G2K2I	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5156	Human VIPR1(Vasoactive Intestinal Peptide Receptor 1) ELISA Kit	VIPR1	P32241	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5157	Human CTLA4(Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit	CTLA4	P16410	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5158	Mouse CTLA4(Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit	CTLA4	P09793	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-ER5159	Rat CTLA4(Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit	CTLA4		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5160	Human CK-17/KRT17(Cytokeratin 17) ELISA Kit	CK-17/KRT17	Q04695	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5161	Human CK-18/KRT18(Cytokeratin 18) ELISA Kit	CK-18/KRT18	P05783	Sandwich	400mIU/mL	3.75 mIU/mL	6.25-400mIU/mL
FY-ER5162	Rat CK-18/KRT18(Cytokeratin 18) ELISA Kit	CK-18/KRT18	Q5BJY9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5163	Human CYP7A1(Cholesterol 7-Alpha-Hydroxylase) ELISA Kit	CYP7A1	P22680	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5164	Mouse CYP7A1(Cholesterol 7-Alpha-Hydroxylase) ELISA Kit	CYP7A1	Q64505	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5165	Rat CYP7A1(Cholesterol 7-Alpha-Hydroxylase) ELISA Kit	CYP7A1	P18125	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH5166	Human CCND1(Cyclin D1) ELISA Kit	CCND1	P24385	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5167	Mouse CCND1(Cyclin D1) ELISA Kit	CCND1	P25322	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5168	Rat CCND1(Cyclin D1) ELISA Kit	CCND1	P39948	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5169	Human CCND2(Cyclin D2) ELISA Kit	CCND2	P30279	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5170	Mouse CCND2(Cyclin D2) ELISA Kit	CCND2	P30280	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ER5171	Rat CCND2(Cyclin D2) ELISA Kit	CCND2		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5172	Human FLNA(Filamin A, Alpha) ELISA Kit	Human FLNA		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5173	Mouse FLNa(Filamin A) ELISA Kit	FLNa	Q8BTM8	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5174	Human Plasmin ELISA Kit	Plasmin		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5175	Human CKMT1A(Creatine kinase U-type, mitochondrial) ELISA Kit	CKMT1A	P12532	Sandwich	200mIU/mL	1.88 mIU/mL	3.13-200mIU/mL
FY-EH5176	Human SNRPD1(Small Nuclear Ribonucleoprotein D1) ELISA Kit	SNRPD1	P62314	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5177	Mouse SNRPD1(Small Nuclear Ribonucleoprotein D1) ELISA Kit	SNRPD1	P62315	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5178	Rat SNRPD1(Small Nuclear Ribonucleoprotein D1) ELISA Kit	SNRPD1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5179	Human STAT5A(Signal Transducer and Activator of Transcription 5A) ELISA Kit	STAT5A	P42229	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL

FY-EM5180	Mouse STAT5A(Signal Transducer and Activator of Transcription 5A) ELISA Kit	STAT5A	P42230	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5181	Human SOX2(Sex Determining Region Y Box Protein 2) ELISA Kit	SOX2	P48431	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5182	Rat SOX2(Sex Determining Region Y Box Protein 2) ELISA Kit	SOX2		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5183	Human TMSβ4(Thymosin Beta 4) ELISA Kit	TMSβ4	P62328	Competitive	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-ER5184	Rat TMSβ4(Thymosin Beta 4) ELISA Kit	TMSβ4		Competitive	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-EH5185	Human HbH(Hemoglobin H) ELISA Kit	HbH		Sandwich	100μg/mL	0.94 μg/mL	1.57-100μg/mL
FY-EH5186	Human SEPP1(Selenoprotein P) ELISA Kit	SEPP1	P49908	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5187	Mouse SEPP1(Selenoprotein P) ELISA Kit	SEPP1	P70274	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5188	Rat SEPP1(Selenoprotein P) ELISA Kit	SEPP1	P25236	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5189	Human SAA(Serum Amyloid A) ELISA Kit	SAA	P0DJJ8	Sandwich	25ng/mL	0.23 ng/mL	0.39-25 ng/mL
FY-EM5190	Mouse SAA(Serum Amyloid A) ELISA Kit	SAA	P05366	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-ER5191	Rat SAA(Serum Amyloid A) ELISA Kit	SAA		Sandwich	100pg/mL	0.94 pg/mL	1.57-100pg/mL
FY-EH5192	Human CR1(Complement Receptor type 1) ELISA Kit	CR1	P17927	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5193	Human BMG/β2-MG(Beta-2-Microglobulin) ELISA Kit	BMG/β2-MG	P61769	Competitive	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EM5194	Mouse BMG/β2-MG(Beta-2-Microglobulin) ELISA Kit	BMG/β2-MG		Sandwich	20ng/mL	0.02 ng/mL	0.03-20ng/mL
FY-ER5195	Rat BMG/β2-MG(Beta-2-Microglobulin) ELISA Kit	BMG/β2-MG	P07151	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5196	Human PAF(Platelet Activating Factor) ELISA Kit	PAF	Q15004	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH5197	Human PDGF-C(Platelet Derived Growth Factor C) ELISA Kit	PDGF-C	Q9NRA1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5198	Mouse PDGF-C(Platelet Derived Growth Factor C) ELISA Kit	PDGF-C	Q8C119	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5199	Rat PDGF-C(Platelet Derived Growth Factor C) ELISA Kit	PDGF-C	Q9EQX6	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5200	Human PDGF-D(Platelet Derived Growth Factor D) ELISA Kit	PDGF-D	Q9GZP0	Sandwich	800pg/mL	7.5 pg/mL	12.5-800pg/mL
FY-EM5201	Mouse PDGF-D(Platelet Derived Growth Factor D) ELISA Kit	PDGF-D	Q92517	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5202	Rat PDGF-D(Platelet Derived Growth Factor D) ELISA Kit	PDGF-D	Q9EQT1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5203	Human NADPH(Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit	NADPH		Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5204	Mouse NADPH(Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit	NADPH		Competitive	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-ER5205	Rat NADPH(Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit	NADPH		Competitive	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EH5206	Human NOX4(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit	NOX4	Q9NPH5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5207	Mouse NOX4(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit	NOX4	Q9JH18	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5208	Human NOX5(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 5) ELISA Kit	NOX5	Q96PH1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5209	Human TAP(Trypsinogen Activation Peptide) ELISA Kit	TAP		Competitive	16pmol/mL	0.15 pmol/mL	0.25-16pmol/mL
FY-EM5210	Mouse TAP(Trypsinogen Activation Peptide) ELISA Kit	TAP	P97435	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5211	Rat TAP(Trypsinogen Activation Peptide) ELISA Kit	TAP		Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5212	Human CA242(Pancreatic Carcinoma Markers-CA242) ELISA Kit	CA242		Sandwich	50IU/mL	0.47 IU/mL	0.79-50IU/mL
FY-EH5213	Human FNDC5(Fibronectin type III domain-containing protein 5) ELISA Kit	FNDC5	Q8NAU1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5214	Mouse FNDC5(Fibronectin type III domain-containing protein 5) ELISA Kit	FNDC5	Q8K4Z2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER5215	Rat FNDC5(Fibronectin type III domain-containing protein 5) ELISA Kit	FNDC5	Q8K3V5	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH5216	Human MEC/CCL28(Mucosae Associated Epithelia Chemokine) ELISA Kit	MEC/CCL28	Q9NRJ3	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5217	Human LpPLA2(Lipoprotein-associated Phospholipase A2) ELISA Kit	LpPLA2	Q13093	Sandwich	10ng/mL	0.024 ng/mL	0.156-10ng/mL
FY-EM5218	Mouse LpPLA2(Lipoprotein-associated Phospholipase A2) ELISA Kit	LpPLA2	Q60963	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5219	Rat LpPLA2(Lipoprotein-associated Phospholipase A2) ELISA Kit	LpPLA2		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EH5220	Human TACE/ADAM17(TNF α Converting Enzyme) ELISA Kit	TACE/ADAM17	P78536	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-ER5221	Rat TACE/ADAM17(TNF α Converting Enzyme) ELISA Kit	TACE/ADAM17		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5222	Human TSGF(Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit	TSGF		Sandwich	1000IU/mL	9.38 IU/mL	15.63-1000IU/mL
FY-EM5223	Mouse TSGF(Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit	TSGF		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5224	Rat TSGF(Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit	TSGF		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5225	Human TFR1(Transferrin Receptor 1) ELISA Kit	TFR1	P02786	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5226	Mouse TFR1(Transferrin Receptor 1) ELISA Kit	TFR1	Q62351	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER5227	Rat TFR1(Transferrin Receptor 1) ELISA Kit	TFR1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5228	Human CD8(Cluster of Differentiation 8) ELISA Kit	CD8	P01732(CD8	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EM5229	Mouse CD8(Cluster of Differentiation 8) ELISA Kit	CD8	P01731(Cd8	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-ER5230	Rat CD8(Cluster of Differentiation 8) ELISA Kit	CD8		Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH5231	Human HS(Heparan Sulfate) ELISA Kit	HS	N/A	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5232	Human SEPCR/PROCR(Soluble Endothelial Protein C Receptor) ELISA Kit	SEPCR/PROCR	Q9UNN8	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM5233	Mouse SEPCR/PROCR(Soluble Endothelial Protein C Receptor) ELISA Kit	SEPCR/PROCR	Q64695	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5234	Human PLA2G2A(Phospholipase A2, Group II A) ELISA Kit	PLA2G2A	P14555	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5235	Rat PLA2G2A(Phospholipase A2, Group II A) ELISA Kit	PLA2G2A		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5236	Human PS1(Presenilin 1) ELISA Kit	PS1	P49768	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5237	Mouse PS1(Presenilin 1) ELISA Kit	PS1	P49769	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5238	Rat PS1(Presenilin 1) ELISA Kit	PS1	P97887	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH5239	Human HMG-CoAR(3-Hydroxy-3-Methylglutaryl Coenzyme A Reductase) ELISA Kit	HMG-CoAR		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER5240	Rat HMG-CoAR(3-Hydroxy-3-Methylglutaryl Coenzyme A Reductase) ELISA Kit	HMG-CoAR	P51639	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5241	Human IRF3(Interferon Regulatory Factor 3) ELISA Kit	IRF3	Q14653	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5242	Human ZPI(Protein Z Dependent Protease Inhibitor) ELISA Kit	ZPI	Q9UK55	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH5243	Human CD74(HLA class II histocompatibility antigen gamma chain) ELISA Kit	CD74	P04233	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH5244	Human RTN4(Reticulon 4) ELISA Kit	RTN4	Q9NQC3	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5245	Mouse RTN4(Reticulon 4) ELISA Kit	RTN4	Q99P72	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5246	Rat RTN4(Reticulon 4) ELISA Kit	RTN4	Q9JK11	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5247	Human MAP2(Microtubule Associated Protein 2) ELISA Kit	MAP2	P11137	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5248	Human OPRM1(Mu-type opioid receptor) ELISA Kit	OPRM1	P35372	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5249	Mouse OPRM1(Mu-type opioid receptor) ELISA Kit	Mouse OPRM1		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5250	Rat OPRM1(Mu-type opioid receptor) ELISA Kit	Rat OPRM1		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5251	Human sAGER(Soluble Advanced Glycosylation End Product Specific Receptor) ELISA Kit	sAGER	Q15109	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH5252	Human IL-26(Interleukin 26) ELISA Kit	IL-26	Q9NPH9	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5253	Human ITGa2B(Integrin Alpha 2B) ELISA Kit	ITGa2B	P08514	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5254	Human ALPL(Alkaline Phosphatase, Liver/Bone/Kidney) ELISA Kit	ALPL	P05186	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5255	Mouse ALPL(Alkaline Phosphatase, Liver/Bone/Kidney) ELISA Kit	ALPL	P09242	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5256	Rat ALPL(Alkaline Phosphatase, Liver/Bone/Kidney) ELISA Kit	ALPL	P08289	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EH5257	Human INS(Insulin) ELISA Kit	INS	P01308	Sandwich	50 μ IU/mL	0.47 μ IU/mL	0.79-50 μ IU/mL
FY-EM5258	Mouse INS(Insulin) ELISA Kit	INS	P01325,P01322	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-ER5259	Rat INS(Insulin) ELISA Kit	INS	P01322,P01325	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5260	Human TSH(Thyroid Stimulating Hormone) ELISA Kit	TSH	P01222(TSH	Sandwich	20 μ IU/mL	0.19 μ IU/mL	0.32-20 μ IU/mL
FY-EM5261	Mouse TSH(Thyroid Stimulating Hormone) ELISA Kit	TSH		Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL

FY-ER5262	Rat TSH(Thyroid Stimulating Hormone) ELISA Kit	TSH		Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-EH5263	Human EPO(Erythropoietin) ELISA Kit	EPO	P01588	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EM5264	Mouse EPO(Erythropoietin) ELISA Kit	EPO	P07321	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER5265	Rat EPO(Erythropoietin) ELISA Kit	EPO	P29676	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH5266	Human F-hCGβ(Free Chorionic Gonadotropin Beta subunit) ELISA Kit	F-hCGβ		Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EH5267	Human CGβ(Chorionic Gonadotrophin Beta) ELISA Kit	CGβ	P0DN86	Sandwich	1000mIU/mL	9.38 mIU/mL	15.63-1000mIU/mL
FY-EH5268	Human TWEAK(Tumour Necrosis Factor Related Weak Inducer of Apoptosis) ELISA Kit	TWEAK	O43508	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH5269	Human hs-CRP(high-sensitivity C-Reactive Protein) ELISA Kit	hs-CRP	P02741	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5270	Mouse hs-CRP(high-sensitivity C-Reactive Protein) ELISA Kit	hs-CRP		Sandwich	5000pg/mL	4.12 pg/mL	6.86-5000pg/mL
FY-ER5271	Rat hs-CRP(high-sensitivity C-Reactive Protein) ELISA Kit	hs-CRP	P48199	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH5272	Human LTF/LF(Lactoferrin) ELISA Kit	LTF/LF	P02788	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5273	Mouse LTF/LF(Lactoferrin) ELISA Kit	LTF/LF	P08071	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER5274	Rat LTF/LF(Lactoferrin) ELISA Kit	LTF/LF		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5275	Human HSD17β10(17-Beta-Hydroxysteroid Dehydrogenase Type 10) ELISA Kit	HSD17β10	Q99714	Sandwich	1600pg/mL	15 pg/mL	25-1600pg/mL
FY-EH5276	Human SORT1(Sortilin) ELISA Kit	SORT1	Q99523	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5277	Human PD-L2(Programmed Cell Death Protein 1 Ligand 2) ELISA Kit	PD-L2		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5278	Rat PD-L2(Programmed Cell Death Protein 1 Ligand 2) ELISA Kit	PD-L2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5279	Human PIVKA- II (Protein Induced by Vitamin K Absence or Antagonist- II) ELISA Kit	PIVKA- II		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5280	Human N-ERC/Mesothelin ELISA Kit	N-ERC/Mesothelin	Q13421(37~)	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5281	Human IL-31(Interleukin 31) ELISA Kit	IL-31	Q6EBC2	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH5282	Human B3GALTL(Beta-1,3-glucosyltransferase) ELISA Kit	B3GALTL	Q6Y288	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH5283	Human CXCR4(Chemokine C-X-C-Motif Receptor 4) ELISA Kit	Human CXCR4	P61073	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH5284	Human C-ERC/Mesothelin ELISA Kit	C-ERC/Mesothelin	Q13421(296)	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH5285	Human fFN(Fetal Fibronectin) ELISA Kit	fFN	P02751	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH5286	Human CF6(Coupling Factor 6) ELISA Kit	CF6	P18859	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5287	Rat CF6(Coupling Factor 6) ELISA Kit	CF6		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH5288	Human Ang1-7(Angiotensin 1-7) ELISA Kit	Human Ang1-7	P01019(34~)	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5289	Mouse Ang1-7(Angiotensin 1-7)ELISA Kit	Ang1-7	P11859(25~)	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5290	Rat Ang1-7(Angiotensin 1-7) ELISA Kit	Ang1-7		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5291	Human γMSH(gamma-Melanocyte Stimulating Hormone) ELISA Kit	γMSH	P01189(77~)	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5292	Mouse γMSH(gamma-Melanocyte Stimulating Hormone) ELISA Kit	γMSH	P01193(77~)	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5293	Human Zonulin ELISA Kit	Zonulin	P00738	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH5294	Human OGDH(2-oxoglutarate dehydrogenase, mitochondrial) ELISA Kit	Human OGDH	Q02218	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER5295	Rat OGDH(2-oxoglutarate dehydrogenase, mitochondrial) ELISA Kit	OGDH	Q5X178	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5296	Human AMPH(Amphiphysin) ELISA Kit	AMPH	P49418	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5297	Human ICAM-3/CD50(Intercellular Adhesion Molecule 3) ELISA Kit	ICAM-3/CD50	P32942	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH5298	Human CD19(Cluster of Differentiation 19) ELISA Kit	CD19	P15391	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5299	Mouse CD19(Cluster of Differentiation 19) ELISA Kit	CD19	P25918	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5300	Human HNF1α(Hepatocyte Nuclear Factor 1 Alpha) ELISA Kit	Human HNF1α	P20823	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5301	Mouse HNF1α(Hepatocyte Nuclear Factor 1 Alpha) ELISA Kit	HNF1α	P22361	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH5302	Human CoQ10B(Coenzyme Q10 Homolog B) ELISA Kit	CoQ10B	Q9H8M1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL

FY-EH5303	Human RANκL(Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit	RANκL	O14788	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5304	Mouse RANκL(Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit	RANκL	O35235	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5305	Rat RANκL(Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit	RANκL	Q9ESE2	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5306	Human M I CA(major histocompatibility complex class I chain-related gene A) ELISA Kit	M I CA	Q29983	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5307	Human CP(Ceruloplasmin) ELISA Kit	CP	P00450	Sandwich	50μg/mL	0.47 μg/mL	0.79-50μg/mL
FY-EM5308	Mouse CP(Ceruloplasmin) ELISA Kit	CP	Q61147	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER5309	Rat CP(Ceruloplasmin) ELISA Kit	CP	P13635	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EH5310	Human NT-proBNP(N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit	NT-proBNP	P16860(NPP)	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5311	Mouse NT-proBNP(N-terminal pro-Brain Natriuretic Peptide) ELISA Kit	NT-proBNP		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5312	Rat NT-proBNP(N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit	NT-proBNP		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH5313	Human PTX 3/TSG-14(Pentraxin 3) ELISA Kit	PTX 3/TSG-14	P26022	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH5314	Human IL-6(Interleukin 6) ELISA Kit	Human IL-6	P05231	Sandwich	100pg/mL	0.94 pg/mL	1.57-100pg/mL
FY-EM5315	Mouse IL-6(Interleukin 6) ELISA Kit	IL-6	P08505	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5316	Rat IL-6(Interleukin 6) ELISA Kit	IL-6	P20607	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5317	Human IL-10(Interleukin 10) ELISA Kit	Human IL-10	P22301	Sandwich	100pg/mL	0.94 pg/mL	1.57-100pg/mL
FY-EM5318	Mouse IL-10(Interleukin 10) ELISA Kit	IL-10	P18893	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5319	Rat IL-10(Interleukin 10) ELISA Kit	IL-10	P29456	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5320	Human POSTN/OSF-2(Periostin) ELISA Kit	POSTN/OSF-2	Q15063	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5321	Mouse POSTN/OSF-2(Periostin) ELISA Kit	POSTN/OSF-2	Q62009	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5322	Rat POSTN/OSF-2(Periostin) ELISA Kit	POSTN/OSF-2	A0A816AHK	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EH5323	Human IgG2(Immunoglobulin G2) ELISA Kit	IgG2	P01859	Competitive	100μg/mL	0.94 μg/mL	1.57-100μg/mL
FY-EH5324	Human GDF2(Growth Differentiation Factor 2) ELISA Kit	GDF2	Q9UK05	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5325	Mouse GDF2(Growth Differentiation Factor 2) ELISA Kit	GDF2	Q9WV56	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-ER5326	Rat GDF2(Growth Differentiation Factor 2) ELISA Kit	GDF2	M0RB87	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5327	Human NRP1(Neuropilin 1) ELISA Kit	NRP1	O14786	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5328	Rat NRP1(Neuropilin 1) ELISA Kit	NRP1	Q9QWJ9	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5329	Human BNP(Brain Natriuretic Peptide) ELISA Kit	BNP	P16860	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-EM5330	Mouse BNP(Brain Natriuretic Peptide) ELISA Kit	BNP	P40753	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5331	Rat BNP(Brain Natriuretic Peptide) ELISA Kit	BNP		Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5332	Human DBP(Vitamin D-binding protein) ELISA Kit	DBP	P02774	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5333	Rat DBP(Vitamin D Binding Protein) ELISA Kit	DBP	P16443	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH5334	Human α-Glu(Alpha Glucosidase) ELISA Kit	α-Glu	P10253	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-EM5335	Mouse α-Glu(Alpha Glucosidase) ELISA Kit	α-Glu		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5336	Rat α-Glu(Alpha Glucosidase) ELISA Kit	α-Glu		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH5337	Human MIF(Macrophage Migration Inhibitory Factor) ELISA Kit	MIF	P14174	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5338	Mouse MIF(Macrophage Migration Inhibitory Factor) ELISA Kit	MIF	P34884	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5339	Rat MIF(Macrophage Migration Inhibitory Factor) ELISA Kit	MIF	P30904	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EH5340	Human ENG(Endoglin) ELISA Kit	ENG	P17813	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5341	Mouse ENG(Endoglin) ELISA Kit	ENG	Q63961	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER5342	Rat ENG(Endoglin) ELISA Kit	ENG		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5343	Human LEPR(Leptin Receptor) ELISA Kit	LEPR	P48357	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL

FY-EM5344	Mouse LEPR(Leptin Receptor) ELISA Kit	LEPR	P48356	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER5345	Rat LEPR(Leptin Receptor) ELISA Kit	LEPR	Q62959	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5346	Human VEGFR1/FLT1(Vascular Endothelial Growth Factor Receptor 1) ELISA Kit	VEGFR1/FLT1	P17948	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5347	Mouse VEGFR1/FLT1(Vascular Endothelial Growth Factor Receptor 1) ELISA Kit	VEGFR1/FLT1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5348	Rat VEGFR1/FLT1(Vascular Endothelial Growth Factor Receptor 1) ELISA Kit	VEGFR1/FLT1	P53767	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5349	Human sTREM-1(soluble Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit	sTREM-1	Q9NP99	Sandwich	100pg/mL	0.94 pg/mL	1.57-100pg/mL
FY-EH5350	Human ACE(Angiotensin I Converting Enzyme) ELISA Kit	ACE	P12821	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5351	Mouse ACE(Angiotensin I Converting Enzyme) ELISA Kit	ACE	P09470	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5352	Rat ACE(Angiotensin I Converting Enzyme) ELISA Kit	ACE	P47820	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5353	Human RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	RANTES	P13501	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5354	Mouse RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	RANTES	P30882	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5355	Rat RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	RANTES	P50231	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5356	Human SELP(P-Selectin) ELISA Kit	SELP	P16109	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5357	Mouse SELP(P-Selectin) ELISA Kit	SELP	Q01102	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5358	Rat SELP(P-Selectin) ELISA Kit	SELP	P98106	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EH5359	Human IL17RA(Interleukin-17 receptor A) ELISA Kit	IL17RA	Q96F46	Sandwich	20pg/mL	0.19 pg/mL	0.32-20pg/mL
FY-EH5360	Human PCSK9(Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit	PCSK9	Q8NBP7	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5361	Mouse PCSK9(Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit	PCSK9	Q80W65	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH5362	Rat PCSK9(Proprotein convertase subtilisin/kexin type 9) ELISA Kit	PCSK9	P59996	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH5363	Human PDGF-AA(Platelet Derived Growth Factor AA) ELISA Kit	PDGF-AA	P04085	Sandwich	50pg/mL	0.47 pg/mL	0.79-50pg/mL
FY-EM5364	Mouse PDGF-AA(Platelet Derived Growth Factor AA) ELISA Kit	PDGF-AA	P20033	Sandwich	50pg/mL	0.47 pg/mL	0.79-50pg/mL
FY-EH5365	Human PF4(Platelet Factor 4) ELISA Kit	Human PF4	P02776	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM5366	Mouse PF4(Platelet Factor 4) ELISA Kit	PF4	Q9Z126	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5367	Rat PF4(Platelet Factor 4) ELISA Kit	PF4	P06765	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5368	Human ACV-A(Activin A) ELISA Kit	ACV-A		Sandwich	5000pg/mL	4.63 pg/mL	78.13-5000 pg/mL
FY-EM5369	Mouse ACV-A(Activin A) ELISA Kit	ACV-A	Q04998	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5370	Rat ACV-A(Activin A) ELISA Kit	ACV-A		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5371	Mouse ANG3(Angiopoietin 3) ELISA Kit	ANG3	P97802	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5372	Mouse CHI3L1(Chitinase 3-like 1) ELISA Kit	CHI3L1	Q61362	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5373	Mouse CS-Ab(Centrosome Antibody) ELISA Kit	CS-Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5374	Mouse IgER-Ab(IgE Receptor Antibody) ELISA Kit	IgER-Ab		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM5375	Mouse OMGP-Ab(Oligodendrocyte Myelin Glycoprotein Antibody) ELISA Kit	OMGP-Ab		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5376	Mouse TR-Ab(Thrombin Receptor Antibody) ELISA Kit	TR-Ab		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5377	Mouse ACF(Apobec 1 Complementation Factor) ELISA Kit	ACF	Q5YD48	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5378	Mouse AATK(Apoptosis Associated Tyrosine Kinase) ELISA Kit	AATK	Q80YE4	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5379	Mouse AIP(Aryl Hydrocarbon Receptor Interacting Protein) ELISA Kit	AIP	O00170	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5380	Rat AIP(Aryl Hydrocarbon Receptor Interacting Protein) ELISA Kit	AIP		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5381	Mouse ARNT2(Aryl Hydrocarbon Receptor Nuclear Translocator 2) ELISA Kit	ARNT2	Q61324	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5382	Rat ARNT2(Aryl Hydrocarbon Receptor Nuclear Translocator 2) ELISA Kit	ARNT2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5383	Mouse BMP8B(Bone Morphogenetic Protein 8B) ELISA Kit	BMP8B	P55105	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5384	Mouse Ca-ATPase(Calcium ATPase) ELISA Kit	Ca-ATPase	G5E829	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL

FY-EM5385	Mouse CAST(Calpastatin) ELISA Kit	CAST	P51125	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5386	Mouse CEBPD(CCAAT/Enhancer Binding Protein Delta) ELISA Kit	CEBPD	Q00322	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5387	Mouse CROP(Cisplatin Resistance Associated Overexpressed Protein) ELISA Kit	CROP	Q5SUF2	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5388	Rat CROP(Cisplatin Resistance Associated Overexpressed Protein) ELISA Kit	CROP		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5389	Mouse CC17(Clara Cell Protein) ELISA Kit	CC17	Q06318	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5390	Rat CC17(Clara Cell Protein) ELISA Kit	CC17		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5391	Mouse CD3(Cluster of differentiation 3) ELISA Kit	CD3		Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER5392	Rat CD3(Cluster of Differentiation 3) ELISA Kit	CD3		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5393	Mouse CD40(Cluster of Differentiation 40) ELISA Kit	CD40	P27512	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5394	Mouse CD40L(Cluster of Differentiation 40 Ligand) ELISA Kit	CD40L	P27548	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-ER5395	Rat CD40L(Cluster Of Differentiation 40 Ligand) ELISA Kit	CD40L	Q9Z2V2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5396	Mouse CSF2Ra(Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit	CSF2Ra	Q00941	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5397	Mouse COL1(Collagen Type I) ELISA Kit	COL1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5398	Rat COL1(Collagen Type I) ELISA Kit	COL1		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EM5399	Mouse COL2(Collagen Type II) ELISA Kit	COL2	P28481	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER5400	Rat COL2(Collagen Type II) ELISA Kit	COL2		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5401	Mouse CDKN1B(Cyclin Dependent Kinase Inhibitor 1B) ELISA Kit	CDKN1B	P46414	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM5402	Mouse DEFA5(Defensin Alpha 5, Paneth Cell Specific) ELISA Kit	DEFA5	P28312	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5403	Mouse H3(Histone H3) ELISA Kit	H3	P68433(Hist	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5404	Mouse D1R(Dopamine D1 receptor) ELISA Kit	D1R	Q61616	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-EM5405	Mouse DRD3(Dopamine Receptor D3) ELISA Kit	DRD3	P30728	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5406	Rat DRD3(Dopamine Receptor D3) ELISA Kit	DRD3		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5407	Mouse DRD4(Dopamine Receptor D4) ELISA Kit	DRD4		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5408	Rat DRD4(Dopamine Receptor D4) ELISA Kit	DRD4	P30729	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5409	Mouse EL(Endothelial lipase) ELISA Kit	EL	Q9WVG5	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5410	Rat EL(Endothelial lipase) ELISA Kit	EL		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5411	Mouse Ear 14(Eosinophil-associated, ribonuclease A family, member 14) ELISA Kit	Ear 14	F6YUX7	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-EM5412	Mouse FRS3(Fibroblast Growth Factor Receptor SubstRate 3) ELISA Kit	FRS3	Q91WJ0	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5413	Rat FRS3(Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit	FRS3		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5414	Mouse FLI1(Friend Leukemia Virus Integration 1) ELISA Kit	FLI1	P26323	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5415	Mouse GNB1(G Protein Beta 1) ELISA Kit	GNb1	P62874	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5416	Rat GNB1(G Protein Beta 1) ELISA Kit	GNb1		Sandwich	1000μIU/mL	9.38 μIU/mL	15.63-1000μIU/mL
FY-EM5417	Mouse GLb(Galactosidase, Beta) ELISA Kit	GLb	P23780	Sandwich	10mIU/mL	0.1 mIU/mL	0.16-10mIU/mL
FY-EM5418	Mouse GAL6(Galectin 6) ELISA Kit	GAL6		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5419	Mouse gGT1(Gamma Glutamyltransferase 1) ELISA Kit	gGT1	Q60928	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5420	Mouse GJb1(Gap Junction Protein Beta 1) ELISA Kit	GJb1	P28230	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5421	Rat GJb1(Gap Junction Protein Beta 1) ELISA Kit	GJb1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5422	Mouse GTF II B(General Transcription Factor II B) ELISA Kit	GTF II B	P62915	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5423	Mouse GTF2H1(General Transcription Factor II H, Polypeptide 1) ELISA Kit	GTF2H1	Q9DBA9	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5424	Mouse GSTa1(Glutathione S Transferase Alpha 1) ELISA Kit	GSTa1	P13745	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5425	Rat GSTA1(Glutathione S Transferase Alpha 1) ELISA Kit	GSTA1	P00502	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EM5426	Mouse GSTo1(Glutathione S Transferase Omega 1) ELISA Kit	GSTo1	O09131	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5427	Rat GSTO1(Glutathione S Transferase Omega 1) ELISA Kit	GSTO1		Sandwich	2000IU/mL	18.75 IU/mL	31.25-2000IU/mL
FY-EM5428	Mouse GSTt1(Glutathione S Transferase Theta 1) ELISA Kit	GSTt1	Q64471	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5429	Mouse GSTt2(Glutathione S Transferase Theta 2) ELISA Kit	GSTt2	Q61133	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-EM5430	Mouse GDF5(Growth Differentiation Factor 5) ELISA Kit	GDF5	P43027	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5431	Rat GDF5(Growth Differentiation Factor 5) ELISA Kit	GDF5	M0RAY4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5432	Mouse HSP-20(Heat Shock Protein 20) ELISA Kit	HSP-20		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5433	Mouse HSP-60(Heat Shock Protein 60) ELISA Kit	HSP-60	P63038	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5434	Rat HSP-60(Heat Shock Protein 60) ELISA Kit	HSP-60		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5435	Mouse HC II(Heparin Cofactor II) ELISA Kit	HC II	P49182	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5436	Rat HC II(Heparin Cofactor II) ELISA Kit	HC II		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EM5437	Mouse HB-EGF(Heparin-binding Epidermal Growth Factor-like Growth Factor) ELISA Kit	HB-EGF	Q06186	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5438	Rat HB-EGF(Heparin-binding Epidermal Growth Factor-like Growth Factor) ELISA Kit	HB-EGF	Q06175	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5439	Mouse HNF1b(Hepatocyte Nuclear Factor 1 Beta) ELISA Kit	HNF1b	P27889	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM5440	Mouse HNF4a(Hepatocyte Nuclear Factor 4 Alpha) ELISA Kit	HNF4a	P49698	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ER5441	Rat HNF4a(Hepatocyte Nuclear Factor 4 Alpha) ELISA Kit	HNF4a		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5442	Mouse HNF6a(Hepatocyte Nuclear Factor 6 Alpha) ELISA Kit	HNF6a	O08755	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-EM5443	Mouse HNF6b(Hepatocyte Nuclear Factor 6 Beta) ELISA Kit	HNF6b	Q6XBJ3	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5444	Mouse HNRPA1(Heterogeneous Nuclear Ribonucleoprotein A1) ELISA Kit	HNRPA1	P49312	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-ER5445	Rat HNRPA1(Heterogeneous Nuclear Ribonucleoprotein A1) ELISA Kit	HNRPA1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5446	Mouse ID(deiodinase) ELISA Kit	ID		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5447	Mouse IgG2a(Immunoglobulin G2a) ELISA Kit	IgG2a	P01863	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-ER5448	Rat IgG2a(Immunoglobulin G2a) ELISA Kit	IgG2a	P20760	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EM5449	Mouse INHbA(Inhibin Beta A) ELISA Kit	INHbA	Q04998	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5450	Mouse INHBP(Inhibin Binding Protein) ELISA Kit	INHBP	Q7TQA1	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5451	Rat INHBP(Inhibin Binding Protein) ELISA Kit	INHBP		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5452	Mouse IFN α -Ab(Interferon- α Antibody) ELISA Kit	IFN α -Ab		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5453	Rat IFN α -Ab(Interferon α Antibody) ELISA Kit	Rat IFN α -Ab		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5454	Mouse FIL1h(Interleukin 1 Eta) ELISA Kit	FIL1h		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5455	Rat FIL1h(Interleukin 1 Eta) ELISA Kit	FIL1h		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5456	Mouse PAL(L-Phenylalanine ammonia-lyase) ELISA Kit	PAL		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5457	Mouse LHRH(Luteinizing Hormone-Releasing Hormone) ELISA Kit	LHRH	P01148(24~26)	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5458	Rat LHRH(Luteinizing Hormone-Releasing Hormone) ELISA Kit	LHRH		Competitive	200mIU/mL	1.88 mIU/mL	3.13-200mIU/mL
FY-EM5459	Mouse MHC-II/H-2 II(Major Histocompatibility Complex-II) ELISA Kit	MHC-II/H-2 II		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5460	Mouse MHC-III/H-2 III(Major Histocompatibility Complex-III) ELISA Kit	MHC-III/H-2 III		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5461	Mouse TRPV1(Transient Receptor Potential Cation Channel Subfamily V, Member 1) ELISA Kit	TRPV1	Q704Y3	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5462	Mouse CHRM1(Cholinergic Receptor, Muscarinic 1) ELISA Kit	CHRM1	P12657	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5463	Mouse NUP153(Nucleoporin 153kDa) ELISA Kit	NUP153	E9Q3G8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5464	Rat NUP153(Nucleoporin 153kDa) ELISA Kit	NUP153		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5465	Mouse NUP160(Nucleoporin 160kDa) ELISA Kit	NUP160	Q9Z0W3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5466	Rat NUP160(Nucleoporin 160kDa) ELISA Kit	NUP160		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL

FY-EM5467	Mouse PIN4(Peptidyl Prolyl Cis/Trans Isomerase NIMA Interacting 4 Protein) ELISA Kit	PIN4	Q9CWW6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5468	Mouse PLCg1(Phospholipase C Gamma 1) ELISA Kit	PLCg1	Q62077	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5469	Rat PLCg1(Phospholipase C Gamma 1) ELISA Kit	PLCg1		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5470	Mouse PACAP(Pituitary Adenylate Cyclase Activating Polypeptide) ELISA Kit	PACAP	O70176	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5471	Rat PACAP(Pituitary Adenylate Cyclase Activating Polypeptide) ELISA Kit	PACAP		Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM5472	Mouse PDGFsR- α (Platelet-Derived Growth Factor Soluble Receptor α) ELISA Kit	PDGFsR- α	P26618	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5473	Rat PDGFsR- α (Platelet-Derived Growth Factor Soluble Receptor α) ELISA Kit	PDGFsR- α		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5474	Mouse PVRL1(Poliiovirus Receptor Related Protein 1) ELISA Kit	PVRL1	Q9JKF6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5475	Rat PVRL1(Poliiovirus Receptor Related Protein 1) ELISA Kit	PVRL1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5476	Mouse PVRL3(Poliiovirus Receptor Related Protein 3) ELISA Kit	PVRL3	Q9JLB9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5477	Mouse Pro-ANP(Pro Atrial Natriuretic Peptide) ELISA Kit	Pro-ANP		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5478	Rat Pro-ANP(Pro Atrial Natriuretic Peptide) ELISA Kit	Pro-ANP		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5479	Mouse PGD2S(Prostaglandin D2 Synthase) ELISA Kit	PGD2S	O09114	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5480	Rat PGD2S(Prostaglandin D2 Synthase) ELISA Kit	PGD2S		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5481	Mouse PSA(Prostate Specific Antigen) ELISA Kit	PSA	P00756	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5482	Rat PSA(Prostate Specific Antigen) ELISA Kit	PSA		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5483	Mouse PKCb1(Protein Kinase C Beta 1) ELISA Kit	PKCb1	P68404	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5484	Rat PKCB1(Protein Kinase C Beta 1) ELISA Kit	PKCB1	P68403	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5485	Mouse PKCd(Protein Kinase C Delta) ELISA Kit	PKCd	P28867	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5486	Rat PKCD(Protein Kinase C Delta) ELISA Kit	PKCD	P09215	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5487	Mouse PKCe(Protein Kinase C Epsilon) ELISA Kit	PKCe	P16054	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5488	Rat PKCE(Protein Kinase C Epsilon) ELISA Kit	PKCE	P09216	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5489	Mouse PKCg(Protein Kinase C Gamma) ELISA Kit	PKCg	P63318	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER5490	Rat PKCG(Protein Kinase C Gamma) ELISA Kit	PKCG		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5491	Mouse PKCi(Protein Kinase C Iota) ELISA Kit	PKCi	Q62074	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5492	Rat PKCI(Protein Kinase C Iota) ELISA Kit	PKCI		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5493	Mouse PKCq(Protein Kinase C Theta) ELISA Kit	PKCq	Q02111	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5494	Rat PKCq(Protein Kinase C Theta) ELISA Kit	PKCq		Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EM5495	Mouse PKCz(Protein Kinase C Zeta) ELISA Kit	PKCz	Q02956	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5496	Rat PKCZ(Protein Kinase C Zeta) ELISA Kit	PKCZ	P09217	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5497	Mouse PKIa(Protein Kinase Inhibitor Alpha) ELISA Kit	PKIa	P63248	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5498	Rat PKIA(Protein Kinase Inhibitor Alpha) ELISA Kit	PKIA		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5499	Mouse PKIb(Protein Kinase Inhibitor Beta) ELISA Kit	PKIb	Q04758	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5500	Rat PKIB(Protein Kinase Inhibitor Beta) ELISA Kit	PKIB		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5501	Mouse PKIg(Protein Kinase Inhibitor Gamma) ELISA Kit	PKIg	O70139	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5502	Rat PKIG(Protein Kinase Inhibitor Gamma) ELISA Kit	PKIG		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5503	Mouse PRKAA1(Protein Kinase, AMP Activated Alpha 1) ELISA Kit	PRKAA1	Q5EG47	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5504	Rat PRKAA1(Protein Kinase, AMP Activated Alpha 1) ELISA Kit	PRKAA1	P54645	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5505	Mouse RAPGEF2(Rap Guanine Nucleotide Exchange Factor 2) ELISA Kit	RAPGEF2	Q8CHG7	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5506	Mouse RasGAP(Ras GTPase Activating Protein) ELISA Kit	RasGAP	F6SEU4	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5507	Mouse RNASEA(Ribonuclease A) ELISA Kit	RNASEA	P00683	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL

FY-ER5508	Rat RNASEA(Ribonuclease A) ELISA Kit	RNASEA		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5509	Mouse RNASEH(Ribonuclease H) ELISA Kit	RNASEH	O70338	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5510	Mouse SESN2(Sestrin 2) ELISA Kit	SESN2	P58043	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5511	Mouse SESN3(Sestrin 3) ELISA Kit	SESN3	Q9CYP7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5512	Mouse MYH7(Myosin Heavy Chain 7, Cardiac Muscle, Beta) ELISA Kit	MYH7	Q91Z83	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5513	Mouse IL-1sR I (Soluble Interleukin-1 Receptor I) ELISA Kit	IL-1sR I	P13504	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5514	Mouse sMHC-1(soluble myosin heavy chain 1) ELISA Kit	sMHC-1	Q5SX40	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5515	Mouse sMHC-2(soluble myosin heavy chain 2) ELISA Kit	sMHC-2	G3UW82	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5516	Mouse SDF-1 β (Stromal Cell Derived Factor 1 Beta) ELISA Kit	SDF-1 β	P40224	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5517	Mouse SYT1(Synaptotagmin 1) ELISA Kit	SYT1	P46096	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5518	Mouse TEK(Tyrosine Kinase, Endothelial) ELISA Kit	TEK	Q02858	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5519	Rat TEK(Tyrosine Kinase, Endothelial) ELISA Kit	TEK		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5520	Mouse TDG1(Teratocarcinoma Derived Growth Factor 1) ELISA Kit	TDG1	P51865	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM5521	Mouse TR(Thrombin Receptor) ELISA Kit	TR	P30558	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5522	Rat TR(Thrombin Receptor) ELISA Kit	TR		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5523	Mouse THBS2(Thrombospondin 2) ELISA Kit	THBS2	Q03350	Sandwich	20mg/mL	0.19 mg/mL	0.32-20mg/mL
FY-EM5524	Mouse THPOR-Ab(Thrombopoietin Receptor Antibody) ELISA Kit	THPOR-Ab		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5525	Mouse FBXO32(F-Box Protein 32) ELISA Kit	FBXO32	Q9CPU7	Sandwich	2500pg/mL	23.44 pg/mL	39.07-2500pg/mL
FY-EM5526	Mouse MuRF1 (muscle-specific RING-finger protein 1) ELISA Kit	MuRF1	Q38HM4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5527	Mouse TMED1(Transmembrane Emp24 Domain Containing Protein 1) ELISA Kit	TMED1	Q3V009	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5528	Rat TMED1(Transmembrane Emp24 Domain Containing Protein 1) ELISA Kit	TMED1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5529	Mouse TUBb(Tubulin Beta) ELISA Kit	TUBb	P99024	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5530	Rat TUBB(Tubulin Beta) ELISA Kit	TUBB		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5531	Mouse α 1-MG(α 1-Microglobulin) ELISA Kit	α 1-MG		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5532	Rat α 1-MG(α 1-Microglobulin) ELISA Kit	α 1-MG	Q63041	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EM5533	Mouse α 2-PI(α 2-plasmin inhibitor) ELISA Kit	α 2-PI	Q61247	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EM5534	Mouse IDUA(Alpha-L-iduronidase) ELISA Kit	IDUA	P48441	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5535	Mouse α NAG(α N-acetylglucosaminidase) ELISA Kit	α NAG		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5536	Mouse β -Glu(Beta Glucosidase) ELISA Kit	β -Glu		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5537	Mouse SEMA3C(Semaphorin-3C) ELISA Kit	SEMA3C	Q62181	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5538	Mouse ANTXR2(Anthrax Toxin Receptor 2) ELISA Kit	ANTXR2	Q6DFX2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5539	Mouse CCNE1(Cyclin-E1) ELISA Kit	CCNE1	Q61457	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5540	Rat CCNE1(Cyclin-E1) ELISA Kit	CCNE1		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5541	Mouse RAR β (Retinoic Acid Receptor Beta) ELISA Kit	RAR β	P22605	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5542	Mouse CCNE2(Cyclin-E2) ELISA Kit	CCNE2	Q9Z238	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5543	Rat CCNE2(Cyclin-E2) ELISA Kit	CCNE2		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5544	Mouse PII6(Peptidase Inhibitor 16) ELISA Kit	PII6	Q9ET66	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5545	Mouse TOP1(Topoisomerase I) ELISA Kit	TOP1	Q04750	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5546	Mouse S100A10(S100 Calcium-binding Protein A10) ELISA Kit	S100A10	P08207	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5547	Mouse G6PC(Glucose-6-Phosphatase, Catalytic) ELISA Kit	G6PC	P35576	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5548	Mouse IgG2b(Immunoglobulin G2b) ELISA Kit	IgG2b		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL

FY-ER5549	Rat IgG2b(Immunoglobulin G2b) ELISA Kit	IgG2b		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5550	Mouse IgG2c(Immunoglobulin G2c) ELISA Kit	IgG2c		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER5551	Rat IgG2c(Immunoglobulin G2c) ELISA Kit	IgG2c		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5552	Mouse ACTα2(Actin Alpha 2, Smooth Muscle) ELISA Kit	ACTα2	P62737	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER5553	Rat ACTα2(Actin Alpha 2, Smooth Muscle) ELISA Kit	ACTα2	P62738	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5554	Mouse CD163(Cluster of Differentiation 163) ELISA Kit	CD163	Q2VLH6	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5555	Rat CD163(Cluster of Differentiation 163) ELISA Kit	CD163		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5556	Mouse IL8Rα(Interleukin 8 Receptor Alpha) ELISA Kit	IL8Rα	Q810W6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5557	Mouse IL8Rβ(Interleukin 8 Receptor Beta) ELISA Kit	IL8Rβ	P35343	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5558	Mouse Slc18a2(Solute Carrier Family 18,Member 2) ELISA Kit	Slc18a2	Q8BRU6	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5559	Mouse Slc18a1(Solute Carrier Family 18,Member 1) ELISA Kit	Slc18a1	Q8R090	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5560	Mouse S1PR1(Sphingosine 1 Phosphate Receptor 1) ELISA Kit	S1PR1	O08530	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5561	Mouse S1PR2(Sphingosine 1 Phosphate Receptor 2) ELISA Kit	S1PR2	P52592	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5562	Mouse S1PR3(Sphingosine 1 Phosphate Receptor 3) ELISA Kit	S1PR3	Q9Z0U9	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5563	Mouse S1PR4(Sphingosine 1 Phosphate Receptor 4) ELISA Kit	S1PR4	Q9Z0L1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5564	Mouse S1PR5(Sphingosine 1 Phosphate Receptor 5) ELISA Kit	S1PR5	Q91X56	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5565	Mouse DEFb14(beta-defensin 14) ELISA Kit	DEFb14	Q7TNV9	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5566	Mouse AK1(Adenylate Kinase 1) ELISA Kit	AK1	Q9R0Y5	Sandwich	50mIU/mL	0.47 mIU/mL	0.79-50mIU/mL
FY-EM5567	Mouse sCD146(Soluble Cluster of Differentiation 146)ELISA Kit	sCD146	Q8R2Y2	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5568	Mouse TITF1(Thyroid transcription factor 1) ELISA Kit	TITF1	P50220	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5569	Mouse CNR2(Cannabinoid Receptor 2) ELISA Kit	CNR2	P47936	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5570	Mouse IARS(Isoleucyl tRNA Synthetase) ELISA Kit	IARS	Q8BU30	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM5571	Mouse ACV-AB(Activin AB) ELISA Kit	ACV-AB		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5572	Mouse DNMT3L(DNA[cytosine-5]-methyltransferase 3-like) ELISA Kit	DNMT3L	Q9CWR8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5573	Mouse LAMP1(Lysosomal Associated Membrane Protein 1) ELISA Kit	LAMP1	P11438	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5574	Mouse MMP-19(Matrix Metalloproteinase 19) ELISA Kit	MMP-19	Q9JHI0	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5575	Mouse AP-1(Transcription Factor AP-1) ELISA Kit	AP-1	P05627	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5576	Mouse APC(Activated Protein C) ELISA Kit	APC	P33587(199)	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5577	Mouse TBK1(TANK Binding Kinase 1) ELISA Kit	TBK1	Q9WUN2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM5578	Mouse sAPPα(soluble amyloid precursor protein α) ELISA Kit	sAPPα		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5579	Mouse CCL6(Chemokine C-C-Motif Ligand 6) ELISA Kit	CCL6	P27784	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5580	Mouse ABHD6(Monoacylglycerol Lipase ABHD6) ELISA Kit	ABHD6	Q8R2Y0	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5581	Mouse RAB35(Ras-related protein Rab-35)ELISA Kit	RAB35	Q6PHN9	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5582	Mouse LIPA(Lipase A, Lysosomal Acid) ELISA Kit	LIPA	Q9Z0M5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5583	Mouse Nrf1(Nuclear Respiratory Factor 1) ELISA Kit	Nrf1	Q9WU00	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5584	Mouse CTSG(Cathepsin G) ELISA Kit	CTSG	P28293	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5585	Mouse TLR7(Toll Like Receptor 7) ELISA Kit	TLR7	P58681	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5586	Mouse RARγ(Retinoic Acid Receptor Gamma) ELISA Kit	RARγ	P18911	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5587	Mouse ACTN1(Actinin Alpha 1) ELISA Kit	ACTN1	Q7TPR4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5588	Mouse TFAM(Transcription Factor A, Mitochondrial) ELISA Kit	TFAM	P40630	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5589	Mouse NEO(Neogenin) ELISA Kit	NEO	P97798	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL

FY-EM5590	Mouse PPAR γ C1 β (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 PPAR γ C1 β) ELISA Kit	PPAR γ C1 β	Q8VHJ7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5591	Mouse CCL12/MCP-5(Monocyte Chemotactic Protein 5) ELISA Kit	CCL12/MCP-5	Q62401	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH5592	Human AMH(Anti-Mullerian Hormone) ELISA Kit	AMH		Sandwich	5000pg/mL	4.63 pg/mL	78.13-5000 pg/mL
FY-EM5593	Mouse AMH(Anti-Mullerian Hormone) ELISA Kit	AMH	P27106	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5594	Rat AMH(Anti-Mullerian Hormone) ELISA Kit	AMH		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH5595	Human Cys-C(Cystatin C) ELISA Kit	Cys-C	P01034	Sandwich	100 ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-EM5596	Mouse Cys-C(Cystatin C) ELISA Kit	Cys-C	P21460	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5597	Rat Cys-C(Cystatin C) ELISA Kit	Cys-C	P14841	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5598	Mouse Klotho(Klotho) ELISA Kit	Klotho	O35082	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5599	Mouse AQP-0(Aquaporin 0) ELISA Kit	AQP-0	P51180	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5600	Mouse CXCL15(Chemokine C-X-C-Motif Ligand 15) ELISA Kit	CXCL15	Q9WVL7	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5601	Mouse F2(Coagulation Factor II) ELISA Kit	F2	P19221	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5602	Rat F2(Coagulation Factor II) ELISA Kit	F2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5603	Mouse CCND3(Cyclin D3) ELISA Kit	CCND3	P30282	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5604	Rat CCND3(Cyclin D3) ELISA Kit	CCND3		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5605	Mouse DNASE1(Deoxyribonuclease 1) ELISA Kit	DNASE1	P49183	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-ER5606	Rat DNASE1(Deoxyribonuclease I) ELISA Kit	DNASE1		Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-EM5607	Mouse EPCAM(Epithelial Cell Adhesion Molecule) ELISA Kit	EPCAM	Q99JW5	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5608	Mouse CXCL5(Chemokine C-X-C-Motif Ligand 5) ELISA Kit	CXCL5	P50228	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5609	Rat CXCL5(Chemokine C-X-C-Motif Ligand 5) ELISA Kit	CXCL5	P97885	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5610	Mouse HABP(Hyaluronate binding protein) ELISA Kit	HABP		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5611	Mouse ICAM-2/CD102(Intercellular Adhesion Molecule 2) ELISA Kit	ICAM-2/CD102	P35330	Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER5612	Rat ICAM-2/CD102(Intercellular Adhesion Molecule 2) ELISA Kit	ICAM-2/CD102		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5613	Mouse FIL1d(Interleukin 1 Delta) ELISA Kit	FIL1d	Q9QYY1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5614	Rat FIL1d(Interleukin 1 Delta) ELISA Kit	FIL1d		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5615	Mouse HAT1(Histone Acetyltransferase 1) ELISA Kit	HAT1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5616	Mouse ATPase(ATPase, Na ⁺ /K ⁺ Transporting) ELISA Kit	ATPase	Q8VDN2	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EM5617	Mouse PVRL2(Poliiovirus Receptor Related Protein 2) ELISA Kit	PVRL2	P32507	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5618	Rat PVRL2(Poliiovirus Receptor Related Protein 2) ELISA Kit	PVRL2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5619	Mouse SPLA2(Secreted Phospholipase A2) ELISA Kit	SPLA2	Q9Z0Y2	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5620	Rat SPLA2(Secreted Phospholipase A2) ELISA Kit	SPLA2		Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-EM5621	Mouse sCD86(Soluble Cluster of Differentiation 86) ELISA Kit	sCD86	P42082	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5622	Rat sCD86(Soluble Cluster of Differentiation 86) ELISA Kit	sCD86		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM5623	Mouse IL-2R α /CD25(Interleukin-2 Receptor alpha chain) ELISA Kit	IL-2R α /CD25		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5624	Rat IL-2R α /CD25(Interleukin-2 Receptor alpha chain) ELISA Kit	IL-2R α /CD25		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5625	Mouse IL-2R β /CD122(Interleukin-2 Receptor beta chain) ELISA Kit	IL-2R β /CD122		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5626	Mouse SNAP-25(Synapsome Associated Protein 25) ELISA Kit	SNAP-25	P60879	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER5627	Rat SNAP-25(Synapsome Associated Protein 25) ELISA Kit	SNAP-25	P60881	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5628	Mouse TAF4(Transcription initiation factor TFIID subunit 4) ELISA Kit	TAF4		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5629	Mouse THPO-Ab(Thrombopoietin Antibody) ELISA Kit	THPO-Ab		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-EM5630	Mouse TGF- α (Transforming Growth Factor Alpha) ELISA Kit	TGF- α	P48030(39~4	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-ER5631	Rat TGF- α (Transforming Growth Factor Alpha) ELISA Kit	TGF- α	P01134	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EH5632	Human TGF- β 3(Transforming Growth Factor Beta 3) ELISA Kit	TGF- β 3		Sandwich	1000pg/mL	5.47 pg/mL	15.63-1000 pg/mL
FY-EM5633	Mouse TGF- β 3(Transforming Growth Factor Beta 3) ELISA Kit	TGF- β 3	P17125(299-	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5634	Mouse TUBa3C(Tubulin Alpha 3c) ELISA Kit	TUBa3C		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5635	Mouse U-T4(Ultrasensitivity Thyroxine) ELISA Kit	U-T4		Competitive	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5636	Rat U-T4(Ultrasensitivity Thyroxine) ELISA Kit	U-T4		Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5637	Mouse β -Hex A(β -hexosaminidase A) ELISA Kit	β -Hex A	P29416	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5638	Rat β -Hex A(β -hexosaminidase A) ELISA Kit	β -Hex A		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-EH5639	Human β TG(Thromboglobulin, Beta) ELISA Kit	β TG		Sandwich	2000pg/mL	13.7 pg/mL	31.25-2000 pg/mL
FY-EM5640	Mouse β TG/PBP/CXCL7/NAP2(Thromboglobulin, Beta) ELISA Kit	β TG/PBP/CXCL7	Q9EQI5	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER5641	Rat β TG/PBP/CXCL7/NAP2(Thromboglobulin, Beta) ELISA Kit	β TG/PBP/CXCL7/NAP2		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5642	Mouse IFN- γ R1(Interferon Gamma Receptor 1) ELISA Kit	IFN- γ R1	P15261	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5643	Mouse MAP1LC3A(Microtubule-associated proteins 1A/1B light chain 3A) ELISA Kit	MAP1LC3A		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EM5644	Mouse MAP1LC3B(Microtubule-associated proteins 1A/1B light chain 3B) ELISA Kit	MAP1LC3B	Q9CQV6	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5645	Mouse LRP-1(Low-Density Lipoprotein-Receptor-Related Protein 1) ELISA Kit	LRP-1	Q91ZX7	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER5646	Rat LRP-1(Low-Density Lipoprotein-Receptor-Related Protein 1) ELISA Kit	LRP-1	G3V928	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-EM5647	Mouse Ang1-9(Angiotensin 1-9) ELISA Kit	Ang1-9	P11859(25~	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5648	Mouse IFN α 4(Interferon Alpha 4) ELISA Kit	IFN α 4	P07351	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5649	Mouse Pro-MMP-9(Pro-Matrix Metalloproteinase-9) ELISA Kit	Pro-MMP-9	P41245	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5650	Mouse gp130(Glucoprotein 130) ELISA Kit	gp130	Q00560	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER5651	Rat gp130(Glucoprotein 130) ELISA Kit	gp130	P40190	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5652	Rat hs-INS(high-sensitivity Insulin) ELISA Kit	hs-INS		Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-ER5653	Rat CAMK2(Calcium/Calmodulin-Dependent Protein Kinase- II) ELISA Kit	CAMK2		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5654	Rat CACT(Carnitine-Acylcarnitine Translocase) ELISA Kit	CACT	P97521	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5655	Rat SCD(Stearoyl Coenzyme A Desaturase) ELISA Kit	SCD		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5656	Rat CD6(Cluster of Differentiation 6) ELISA Kit	CD6		Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-ER5657	Rat COL9 α 1(Collagen Type IX Alpha 1) ELISA Kit	COL9 α 1		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5658	Rat Collagenase I ELISA Kit	Collagenase I		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5659	Rat CVL(Cytovillin) ELISA Kit	CVL		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5660	Rat D-5(Defensin-5) ELISA Kit	D-5		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER5661	Rat DRD1(Dopamine Receptor D1) ELISA Kit	DRD1	P18901	Sandwich	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-ER5662	Rat DRD2(Dopamine Receptor D2) ELISA Kit	DRD2	P61169	Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-ER5663	Rat ACS(Fatty Acyl-CoA synthetase) ELISA Kit	ACS		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5664	Rat GTF II B(General Transcription Factor II B) ELISA Kit	GTF II B		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER5665	Rat G-6-Pase/G6PT/GSD(Glucose 6 Phosphate) ELISA Kit	G-6-Pase/G6PT/GSD		Sandwich	5000mIU/mL	46.88 mIU/mL	78.13-5000mIU/mL
FY-ER5666	Rat GSTT1(Glutathione S Transferase Theta 1) ELISA Kit	GSTT1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5667	Rat GSTT2(Glutathione S Transferase Theta 2) ELISA Kit	GSTT2		Sandwich	100mIU/mL	0.94 mIU/mL	1.57-100mIU/mL
FY-ER5668	Rat GnRHR(Gonadotropin-Releasing Hormone Receptor) ELISA Kit	GnRHR		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5669	Rat HSPA4(Heat Shock 70kDa Protein 4) ELISA Kit	HSPA4		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5670	Rat HSP47(Heat Shock Protein 47) ELISA Kit	HSP47		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5671	Rat HNF1a(Hepatocyte Nuclear Factor 1 Alpha) ELISA Kit	HNF1a		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-ER5672	Rat HAT1(Histone Acetyltransferase 1) ELISA Kit	HAT1	Q5M939	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5673	Mouse histon-H2b ELISA Kit	histon-H2b		Sandwich	25ng/mL	0.24 ng/mL	0.4-25ng/mL
FY-ER5674	Rat Histon-H2b ELISA Kit	Histon-H2b		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5675	Rat IDS(Iduronate Sulfatase) ELISA Kit	IDS		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5676	Rat IKB α (Inhibitory Subunit of NF- κ B α) ELISA Kit	IKB α		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5677	Rat IFI16/p16(Interferon-Inducible Protein 16) ELISA Kit	IFI16/p16		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5678	Rat PEPCK(Phosphoenolpyruvate Carboxykinase) ELISA Kit	PEPCK		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5679	Rat MHC II /AgB II /H-1 II (Major Histocompatibility Complex II) ELISA Kit	MHC II /AgB II /H-1 II		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-ER5680	Rat MHC III/AgB III/H-1 III (Major Histocompatibility Complex III) ELISA Kit	MHC III/AgB III/H-1 III		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-EM5681	Mouse AMT(Aminomethyltransferase) ELISA Kit	AMT	Q8CFA2	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5682	Rat MAT(Matrilysin) ELISA Kit	MAT		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER5683	Rat PTEN(Phosphatase and Tensin Homolog) ELISA Kit	PTEN	O54857	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5684	Rat NOX1(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit	NOX1	Q9WV87	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5685	Rat NOX4(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit	NOX4	Q924V1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5686	Rat NUP198(Nucleoporin 198kDa) ELISA Kit	NUP198		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5687	Rat PMP3(Peroxisomal Membrane Protein 3) ELISA Kit	PMP3		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5688	Rat pACCase(Phosphorylated Acetyl-CoA Carboxylase) ELISA Kit	pACCase		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5689	Rat PICK1(Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit	PICK1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5690	Rat PKD1(Protein Kinase D1) ELISA Kit	PKD1		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5691	Rat PKD2(Protein Kinase D2) ELISA Kit	PKD2		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5692	Rat PKN1(Protein Kinase N1) ELISA Kit	PKN1		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5693	Rat GRAP2(GRB2-Related Adapter protein 2) ELISA Kit	GRAP2		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5694	Rat PDK4(Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit	PDK4		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5695	Rat STAT5B(Signal Transducer And Activator Of Transcription 5B) ELISA Kit	STAT5B		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5696	Rat SSTR1(Somatostatin Receptor 1) ELISA Kit	SSTR1	P28646	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5697	Rat ChAT(Choline Acetyltransferase) ELISA Kit	ChAT	P32738	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5698	Rat MSRA(Methionine Sulfoxide Reductase A) ELISA Kit	MSRA		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5699	Rat MSRB1(Methionine Sulfoxide Reductase B1) ELISA Kit	MSRB1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5700	Rat SERP1(Stress-associated Endoplasmic Reticulum Protein 1) ELISA Kit	SERP1		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5701	Rat ERR α (Estrogen Related Receptor Alpha) ELISA Kit	ERR α	Q5QJV7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5702	Rat HRAS(V-Ha-Ras Harvey Rat Sarcoma Viral Oncogene Homolog) ELISA Kit	HRAS		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5703	Rat XIAP(X-Linked Inhibitor Of Apoptosis Protein) ELISA Kit	XIAP	Q9R0I6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5704	Rat MDH1(Malate Dehydrogenase 1) ELISA Kit	MDH1	O88989	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5705	Rat FLAP(5-Lipoxygenase Activating Protein) ELISA Kit	FLAP		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5706	Rat GNMT(Glycine N-methyltransferase) ELISA Kit	GNMT		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5707	Rat NTSR2(Neurotensin Receptor Type 2) ELISA Kit	NTSR2		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5708	Rat HSD11 β 1(11-Beta-Hydroxysteroid Dehydrogenase Type 1) ELISA Kit	HSD11 β 1	P16232	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5709	Rat RGN(Regucalcin) ELISA Kit	RGN		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5710	Rat OXM(Oxyntomodulin) ELISA Kit	OXM		Competitive	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5711	Rat IRF4(Interferon Regulatory Factor 4) ELISA Kit	IRF4	D3ZEX6	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5712	Rat UGT1A1(UDP Glucuronosyltransferase 1 Family, Polypeptide A1) ELISA Kit	UGT1A1	Q64550	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL

FY-ER5713	Rat PEMT(Phosphatidylethanolamine N-methyltransferase) ELISA Kit	PEMT		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5714	Rat FECH(Ferrochelatase) ELISA Kit	FECH		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5715	Rat BCL2L1(Bcl-2-Like Protein 1) ELISA Kit	BCL2L1		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5716	Rat MC4R(Melanocortin Receptor 4) ELISA Kit	MC4R	P70596	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-ER5717	Rat ICD(Isocitrate Dehydrogenase) ELISA Kit	ICD		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5718	Rat GPIb(Platelet Membrane Glycoprotein Ib) ELISA Kit	GPIb		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5719	Rat GALR2(Galanin Receptor2) ELISA Kit	GALR2		Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-ER5720	Rat GALR3(Galanin Receptor3) ELISA Kit	GALR3		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5721	Rat SOD4(Superoxide Dismutase 4, Copper Chaperone) ELISA Kit	SOD4		Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5722	Rat GPIb-IX-V(Platelet Membrane Glycoprotein GPIb-IX-V Complex) ELISA Kit	GPIb-IX-V		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5723	Rat CYP2C9(Cytochrome P450, family 2, subfamily C, polypeptide 9) ELISA Kit	CYP2C9		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5724	Rat CYP2C19(Cytochrome P450, family 2, subfamily C, polypeptide 19) ELISA Kit	CYP2C19		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5725	Rat ALAD(Aminolevulinatase Delta Dehydratase) ELISA Kit	ALAD	P06214	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5726	Rat LIPE(Lipase, Hormone Sensitive) ELISA Kit	LIPE	P15304	Sandwich	5ng/mL	0.05 ng/mL	0.08-5ng/mL
FY-ER5727	Rat sCD80(Soluble Cluster of Differentiation 80) ELISA Kit	sCD80		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5728	Rat PKCα(Protein Kinase Cα) ELISA Kit	PKCα		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5729	Rat SDF-1α(Stromal Cell Derived Factor 1α) ELISA Kit	SDF-1α		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5730	Rat dGHRl(Desacyl Ghrelin) ELISA Kit	dGHRl		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5731	Rat ADRβ1(Adrenergic Receptor Beta 1) ELISA Kit	ADRβ1	P18090	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5732	Rat ADRβK1(Adrenergic Receptor Beta Kinase 1) ELISA Kit	ADRβK1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5733	Rat F-actin ELISA Kit	F-actin		Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-ER5734	Rat ADORA2b(Adenosine A2b Receptor) ELISA Kit	ADORA2b		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5735	Rat PC(Pyruvate Carboxylase) ELISA Kit	PC		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5736	Rat CAPN2(Calpain-2 catalytic subunit) ELISA Kit	CAPN2		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5737	Rat β-Hex B(β-hexosaminidase B) ELISA Kit	β-Hex B		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5738	Rat ICOS(Inducible T-cell costimulator) ELISA Kit	ICOS		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5739	Rat ICOSLG(Inducible T-Cell Co Stimulator Ligand) ELISA Kit	ICOSLG		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5740	Rat ERK1/2(Extracellular Signal Regulated Kinase 1/2) ELISA Kit	ERK1/2		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5741	Rat AZU/HBP(Azurocidin/Heparin Binding Protein) ELISA Kit	AZU/HBP		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5742	Rat CD68(Cluster of Differentiation 68) ELISA Kit	CD68		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5743	Rat sE-selectin(soluble E-selectin) ELISA Kit	sE-selectin		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER5744	Rat δSIP(Delta Sleep Inducing Peptide) ELISA Kit	δSIP		Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5745	Rat GSTP(Glutathione S-transferase P) ELISA Kit	GSTP	P04906	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5746	Rat HRH4(Histamine Receptor H4) ELISA Kit	HRH4	Q91ZY1	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5747	Rat RFTN1(Raft Linking Protein) ELISA Kit	RFTN1	Q4G023	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5748	Rat LXRβ(liver X Receptor Beta) ELISA Kit	LXRβ		Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER5749	Rat APLNR(Apelin Receptor) ELISA Kit	APLNR	Q9JHG3	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5750	Rat NEFH(Neurofilament, Heavy Polypeptide) ELISA Kit	NEFH		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5751	Rat WNT10B(Wingless-Type MMTV Integration Site Family, Member 10B) ELISA Kit	WNT10B		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5752	Rat ACTα1(Actin Alpha 1, Skeletal Muscle) ELISA Kit	ACTα1		Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-ER5753	Rat MC3R(Melanocortin Receptor 3) ELISA Kit	MC3R		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL

FY-ER5754	Rat JNK1(C-Jun N-Terminal Kinases 1) ELISA Kit	JNK1		Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-ER5755	Rat CD44(Cluster of Differentiation 44) ELISA Kit	CD44		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5756	Rat F3(Coagulation Factor III) ELISA Kit	F3		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5757	Rat CTX-III(Cross Linked C-telopeptide of Type III Collagen) ELISA Kit	CTX-III		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5758	Rat HSPG(Heparan Sulfate Proteoglycan) ELISA Kit	HSPG	P34900	Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5759	Rat INHbC(Inhibin Beta C) ELISA Kit	INHbC	Q9WUK5	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5760	Rat MHC I /AgB I /H-1 I (Major Histocompatibility Complex I) ELISA Kit	MHC I /AgB I /H-1 I		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5761	Rat REG(Regenerating Islet Derived Protein) ELISA Kit	REG		Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-ER5762	Rat TLR6(Toll-Like Receptor 6) ELISA Kit	TLR6		Sandwich	20ng/mL	0.19 ng/mL	0.32-20ng/mL
FY-ER5763	Rat TAF1A(Tata Box Binding Protein Associated Factor 1A) ELISA Kit	TAF1A		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-ER5764	Rat TXB3(Thromboxanes B3) ELISA Kit	TXB3		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5765	Rat PTX-3/TSG-14(Pentraxin 3) ELISA Kit	PTX-3/TSG-14		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5766	Rat TGF-βRII(Transforming Growth Factor Receptor II) ELISA Kit	TGF-βRII		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-ER5767	Rat ADAMTS13/VWFPC(Von Willebrand Factor Cleaving Protease) ELISA Kit	ADAMTS13/VWFPC		Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5768	Rat PFK-A(6-phosphofructokinase Type A) ELISA Kit	PFK-A		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5769	Rat Notch 2(Notch Homolog 2) ELISA Kit	Notch 2	Q9QW30	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-ER5770	Rat GFRA1(Glial Cell Line Derived Neurotrophic Factor Receptor Alpha 1) ELISA Kit	GFRA1	Q62997	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5771	Mouse MHC I/H-2 I (Major Histocompatibility Complex I) ELISA Kit	MHC I/H-2 I		Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM5772	Mouse Gal(α-Galactosyl Epitopes) ELISA Kit	Gal		Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-EM5773	Mouse KARSL(Lysyl tRNA Synthetase) ELISA Kit	KARSL	Q99MN1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EH5774	Human RET(Ret Proto-Oncogene) ELISA Kit	RET	P07949	Sandwich	100pg/mL	0.94 pg/mL	1.56-100 pg/mL
FY-EH5775	Human TM(Thrombomodulin) ELISA Kit	TM	P07204	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-EM5776	Mouse TM(Thrombomodulin) ELISA Kit	TM	P15306	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-ER5777	Rat TM(Thrombomodulin) ELISA Kit	TM		Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-EH5778	Human ANG(Angiostatin) ELISA Kit	ANG	P00747(79-4	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-EH5779	Human MYC(V-Myc Myelocytomatosis Viral Oncogene Homolog) ELISA Kit	MYC	B3CJ74	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-EH5780	Human GLO1(Glyoxalase I) ELISA Kit	GLO1	Q04760	Sandwich	100ng/mL	0.94 ng/mL	1.57-100ng/mL
FY-ER5781	Rat GLO1(Glyoxalase I) ELISA Kit	GLO1		Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EH5782	Human IAP(immunosuppressive acidic protein) ELISA Kit	IAP	Q08722	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH5783	Human PL(Pancreatic Lipase) ELISA Kit	PL	P16233	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EM5784	Mouse PL(Pancreatic Lipase) ELISA Kit	PL	Q6P8U6	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER5785	Rat PL(Pancreatic Lipase) ELISA Kit	PL	P27657	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-EH5786	Human PP(Pancreatic Polypeptide) ELISA Kit	PP	P01298	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-EM5787	Mouse PP(Pancreatic Polypeptide) ELISA Kit	PP	P10601	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5788	Human GS(Gelsolin) ELISA Kit	GS	P06396	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-EM5789	Mouse GS(Gelsolin) ELISA Kit	GS	P13020	Sandwich	500ng/mL	4.69 ng/mL	7.82-500ng/mL
FY-ER5790	Rat GS(Gelsolin) ELISA Kit	GS	P09606	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-EH5791	Human GR(Glutathione Reductase) ELISA Kit	GR	P00390	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-EM5792	Mouse GR(Glutathione Reductase) ELISA Kit	GR	P47791	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER5793	Rat GR(Glutathione Reductase) ELISA Kit	GR	P70619	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-EM5794	Mouse TF(Tissue factor) ELISA Kit	TF	P20352	Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL

FY-ER5795	Rat TF(Tissue Factor) ELISA Kit	TF	P42533	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5796	Human TPO(Thrombopoietin) ELISA Kit	TPO	P40225	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EM5797	Mouse TPO(Thrombopoietin) ELISA Kit	TPO	P40226	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER5798	Rat TPO(Thrombopoietin) ELISA Kit	TPO	P49745	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-EH5799	Human TG(Thyroglobulin) ELISA Kit	TG	P01266	Sandwich	50ng/mL	0.47 ng/mL	0.79-50ng/mL
FY-EM5800	Mouse TG(Thyroglobulin) ELISA Kit	TG	O08710	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-ER5801	Rat TG(Thyroglobulin) ELISA Kit	TG	P06882	Sandwich	500pg/mL	4.69 pg/mL	7.82-500pg/mL
FY-ER5802	Rat NLRP3 (NACHT, LRR and PYD domains-containing protein 3) ELISA kit	NLRP3	D4A523	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-ER5803	Rat ERFE(Erythroferrone) ELISA Kit	ERFE	D4AB34	Sandwich	10ng/mL	0.094 ng/mL	0.156-10ng/mL
FY-EM5804	Mouse ALT(Alanine Transaminase) ELISA Kit	ALT	Q8QZR5	Sandwich	10mU/mL	0.09 mU/mL	0.16-10mU/mL
FY-EH5805	Human AEA(Anandamide) ELISA Kit	AEA		Competitive	200ng/mL	2.47 ng/mL	2.47-200ng/mL
FY-EH5806	Human HbA1c(Glycated Hemoglobin) ELISA Kit	HbA1c		Sandwich	200ng/mL	1.48 ng/mL	2.47-200ng/mL
FY-EH5807	Human 2-AG(2-Arachidonoylglycerol) ELISA Kit2-AG	2-AG		Competitive	300ng/mL	3.7 ng/mL	3.70-300ng/mL
New							
FY-B0001	Bovine BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	BMP-2		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0002	Bovine BMP-4(Bone Morphogenetic Protein 4)ELISA kit	BMP-4		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0003	Bovine ECF/CCL11(Eosinophil Chemotactic Factor)ELISA kit	ECF/CCL11		Sandwich	4000pg/mL	37.5pg/mL	62.5-4000pg/mL
FY-B0004	Bovine IP-10/CXCL10(Interferon Gamma Induced Protein 10kDa) ELISA Kit	IP-10/CXCL10		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0005	Bovine EGF(Epidermal Growth Factor) ELISA Kit	EGF		Sandwich	10ng/mL	0.067 ng/mL	0.16-10 ng/mL
FY-B0006	Bovine EG-VEGF(Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit	EG-VEGF		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0007	Bovine ET-1(Endothelin1) ELISA Kit	ET-1		Sandwich	1000pg/mL	4.74 pg/mL	15.63-1000 pg/mL
FY-B0008	Bovine AFGF/FGF1(Acidic Fibroblast Growth Factor 1)ELISA Kit	AFGF/FGF1		Sandwich	4000pg/mL	37.5pg/mL	62.5-4000pg/mL
FY-B0009	Bovine HGF(Hepatocyte Growth Factor) ELISA Kit	HGF		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0010	Bovine IFN-β(Interferon Beta)ELISA Kit	IFN-β		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0011	Bovine IL-1α(Interleukin 1 Alpha) ELISA Kit	IL-1α		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0012	Bovine IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist)ELISA kit	IL-1ra/IL-1F3		Sandwich	2000pg/mL	18.8pg/mL	31.25-2000pg/mL
FY-B0013	Bovine LIF(Leukemia Inhibitory Factor) ELISA Kit	LIF		Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-B0014	Bovine IL-4(Interleukin 4) ELISA Kit	IL-4		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0015	Bovine IL-13(Interleukin 13) ELISA Kit	IL-13		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0016	Bovine IL-22(Interleukin 22) ELISA Kit	IL-22		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0017	Bovine IFN-γ(Interferon Gamma) ELISA Kit	IFN-γ		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0018	Bovine TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	TNF-α		Sandwich	500pg/mL	3.1pg/mL	7.81-500pg/mL
FY-B0019	Bovine VEGF-A(Vascular Endothelial Growth Factor A) ELISA Kit	VEGF-A		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0020	Bovine PRL(Prolactin) ELISA Kit	PRL		Sandwich	200ng/mL	1.39ng/mL	3.12-200ng/mL
FY-B0021	Bovine IL-1β(Interleukin 1 Beta) ELISA Kit	IL-1β		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0022	Bovine IL-12(Interleukin 12) ELISA Kit	IL-12		Sandwich	10ng/mL	0.055 ng/mL	0.16-10 ng/mL
FY-B0023	Bovine IL-12 p40(Interleukin 12 p40)ELISA Kit	IL-12 p40		Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0024	Bovine COL4(Collagen Type IV) ELISA Kit	COL4		Sandwich	5000pg/mL	46.89pg/mL	78.13-5000pg/mL
FY-B0025	Bovine FN(Fibronectin) ELISA Kit	FN		Sandwich	10ng/mL	0.1ng/mL	0.16-10ng/mL
FY-B0026	Bovine PIIIINP(N-Terminal Procollagen III Propeptide) ELISA Kit	PIIINP		Sandwich	400ng/mL	2.33 ng/mL	6.25-400 ng/mL
FY-B0027	Bovine TIMP-1(Tissue Inhibitors Of Metalloproteinase 1)ELISA kit	TIMP-1		Sandwich	150ng/mL	1.41ng/mL	2.34-150ng/mL

FY-B0028	Bovine IL-5(Interleukin 5) ELISA Kit	IL-5	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0029	Bovine IL-15(Interleukin 15) ELISA Kit	IL-15	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0030	Bovine IL-16(Interleukin 16) ELISA Kit	IL-16	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0031	Bovine IL-18(Interleukin 18) ELISA Kit	IL-18	Sandwich	1000pg/mL	6.1pg/mL	15.63-1000pg/mL
FY-B0032	Bovine ADM(Adrenomedullin) ELISA Kit	ADM	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0033	Bovine INHB(Inhibin B) ELISA Kit	INHB	Sandwich	2000pg/mL	11.5pg/mL	31.25-2000pg/mL
FY-B0034	Bovine APOA4(Apolipoprotein A4) ELISA Kit	APOA4	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-B0035	Bovine AR(Amphiregulin) ELISA kit	AR	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-B0036	Bovine AST(Aspartate Aminotransferase) ELISA Kit	AST	Sandwich	10ng/mL	0.1ng/mL	0.16-10ng/mL
FY-B0037	Bovine Bax(Bcl-2 Associated X Protein) ELISA Kit	Bax	Sandwich	20ng/mL	0.12ng/mL	0.31-20ng/mL
FY-B0038	Bovine CASP8(Caspase 8) ELISA Kit	CASP8	Sandwich	20ng/mL	0.16 ng/mL	0.32-20 ng/mL
FY-B0039	Bovine CTSD(Cathepsin D) ELISA Kit	CTSD	Sandwich	20ng/mL	0.13ng/mL	0.31-20ng/mL
FY-B0040	Bovine CCKAR(Cholecystokinin A Receptor) ELISA Kit	CCKAR	Sandwich	20ng/mL	0.12ng/mL	0.31-20ng/mL
FY-B0041	Bovine NOS2/iNOS(NitricOxide Synthase 2, Inducible) ELISA Kit	NOS2/iNOS	Sandwich	20ng/mL	0.126 ng/mL	0.32-20 ng/mL
FY-B0042	Bovine NOS3/eNOS(Nitric Oxide Synthase 3, Endothelial)ELISA Kit	NOS3/eNOS	Sandwich	1000pg/mL	6.1 pg/mL	15.63-1000 pg/mL
FY-B0043	Bovine COL1 α 2(Collagen Type I Alpha 2)ELISA kit	COL1 α 2	Sandwich	400ng/mL	3.75ng/mL	6.25-400ng/mL
FY-B0044	Bovine α -LA(Alpha-Lactalbumin) ELISA Kit	α -LA	Sandwich	100ng/mL	0.67ng/mL	1.56-100ng/mL
FY-B0045	Bovine DEFB1(Defensin Beta 1) ELISA Kit	DEFB1	Sandwich	10ng/mL	0.1ng/mL	0.16-10ng/mL
FY-B0046	Bovine DEFB2/DEFB2(Defensin Beta 2)ELISA kit	DEFB2/DEFB2	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0047	Bovine IL-2R α (Interleukin 2 Receptor Alpha)ELISA kit	IL-2R α	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0048	Bovine GAL1(Galectin 1) ELISA Kit	GAL1	Sandwich	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-B0049	Bovine GAL8(Galectin 8) ELISA Kit	GAL8	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-B0050	Bovine FGF9(Fibroblast Growth Factor 9) ELISA Kit	FGF9	Sandwich	4000pg/mL	23.2pg/mL	62.5-4000pg/mL
FY-B0051	Bovine ELN(Elastin) ELISA Kit	ELN	Sandwich	200ng/mL	1.34ng/mL	3.12-200ng/mL
FY-B0052	Bovine NGF(Nerve Growth Factor) ELISA Kit	NGF	Sandwich	10ng/mL	0.1ng/mL	0.16-10ng/mL
FY-B0053	Bovine RETN(Resistin) ELISA Kit	RETN	Sandwich	5000pg/mL	31pg/mL	78.13-5000pg/mL
FY-B0054	Bovine PIGR/SC(Polymeric Immunoglobulin Receptor/Membrane Secretory Compon	PIGR/SC	Sandwich	40ng/mL	0.38ng/mL	0.63-40ng/mL
FY-B0055	Bovine PTN(Pleiotrophin) ELISA Kit	PTN	Sandwich	2000pg/mL	11.4pg/mL	31.25-2000pg/mL
FY-B0056	Bovine sIgA(Secretory Immunoglobulin A) ELISA Kit	sIgA	Sandwich	50ng/mL	0.35ng/mL	0.78-50ng/mL
FY-B0057	Bovine S100A8(S100 Calcium Binding Protein A8) ELISA Kit	S100A8	Sandwich	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-B0058	Bovine S100A9(S100 Calcium Binding Protein A9) ELISA Kit	S100A9	Sandwich	20ng/mL	0.1ng/mL	0.31-20ng/mL
FY-B0059	Bovine SCF(Stem Cell Factor) ELISA Kit	SCF	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0060	Bovine OC/BGP(Osteocalcin) ELISA Kit	OC/BGP	Sandwich	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-B0061	Bovine OPN(Osteopontin) ELISA Kit	OPN	Sandwich	100ng/mL	0.55ng/mL	1.56-100ng/mL
FY-B0062	Bovine MMP-3(Matrix Metalloproteinase 3) ELISA Kit	MMP-3	Sandwich	10ng/mL	0.064ng/mL	0.16-10ng/mL
FY-B0063	Bovine MMP-7(Matrix Metalloproteinase 7) ELISA Kit	MMP-7	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0064	Bovine TIMP-2(Tissue Inhibitors Of Metalloproteinase 2) ELISA Kit	TIMP-2	Sandwich	4000pg/mL	26.5pg/mL	62.5-4000pg/mL
FY-B0065	Bovine TIMP-3(Tissue Inhibitors Of Metalloproteinase 3) ELISA Kit	TIMP-3	Sandwich	20ng/mL	0.12ng/mL	0.31-20ng/mL
FY-B0066	Bovine TIMP-4(Tissue Inhibitors Of Metalloproteinase 4) ELISA Kit	TIMP-4	Sandwich	20ng/mL	0.13 ng/mL	0.31-20ng/mL
FY-B0067	Bovine PTHrP(Parathyroid Hormone Related Protein) ELISA Kit	PTHrP	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0068	Bovine PCT (procalcitonin) ELISA kit	PCT	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000 pg/mL

FY-B0069	Bovine PGF(Placenta Growth Factor)ELISA kit	PGF	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0070	Bovine LOX-1(Lectin Like Oxidized Low Density Lipoprotein Receptor 1)ELISA kit	LOX-1	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-B0071	Bovine IL-12 p35(Interleukin 12 p35)ELISA Kit	IL-12 p35	Sandwich	500pg/mL	4.69pg/mL	7.81-500pg/mL
FY-B0072	Bovine TYR(Tyrosinase) ELISA Kit	TYR	Sandwich	500ng/mL	3.1ng/mL	7.81-500ng/mL
FY-B0073	Bovine VF(Visfatin) ELISA Kit	VF	Sandwich	200ng/mL	1.27ng/mL	3.13-200ng/mL
FY-B0074	Bovine OTC(Ornithine Transcarbamylase)ELISA kit	OTC	Sandwich	400ng/mL	3.75ng/mL	6.25-400ng/mL
FY-B0075	Bovine PLAU/uPA(Plasminogen Activator, Urokinase)ELISA kit	PLAU/uPA	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-B0076	Bovine DPT(Dermatopontin) ELISA Kit	DPT	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-B0077	Bovine HSP-90(Heat Shock Protein 90)ELISA kit	HSP-90	Sandwich	200ng/mL	1.88ng/mL	3.13-200ng/mL
FY-B0078	Bovine LZM(Lysozyme) ELISA Kit	LZM	Sandwich	20ng/mL	0.13ng/mL	0.31-20ng/mL
FY-B0079	Bovine PEDF(Pigment Epithelium Derived Factor) ELISA Kit	PEDF	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-B0080	Bovine SDH(Sorbitol Dehydrogenase) ELISA Kit	SDH	Sandwich	200ng/mL	1.13ng/mL	3.13-200ng/mL
FY-B0081	Bovine FTL(Ferritin, Light Polypeptide) ELISA Kit	FTL	Sandwich	50ng/mL	0.28ng/mL	0.78-50ng/mL
FY-B0082	Bovine DNase-I(Deoxyribonuclease I)ELISA kit	DNase-I	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-B0083	Bovine PGA(Pepsinogen A) ELISA Kit	PGA	Sandwich	20ng/mL	0.11ng/mL	0.31-20ng/mL
FY-B0084	Bovine Plg(Plasminogen) ELISA Kit	Plg	Sandwich	400ng/mL	3.75ng/mL	6.25-400ng/mL
FY-B0085	Bovine OSM(Oncostatin M) ELISA Kit	OSM	Sandwich	500pg/mL	4.69pg/mL	7.81-500pg/mL
FY-B0086	Bovine TNF-β(Tumor Necrosis Factor Beta)ELISA kit	TNF-β	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-B0087	Bovine LL-37(Antibacterial Protein LL-37) ELISA Kit	LL-37	Sandwich	50ng/mL	0.27 ng/mL	0.79-50 ng/mL
FY-B0088	Bovine GPX1(Glutathione Peroxidase 1) ELISA Kit	GPX1	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0089	Bovine ROMO1(Reactive Oxygen Species Modulator 1) ELISA Kit	ROMO1	Sandwich	20ng/mL	0.12ng/mL	0.31-20ng/mL
FY-B0090	Bovine QSOX1(Quiescin Q6 Sulfhydryl Oxidase 1) ELISA Kit	QSOX1	Sandwich	20ng/mL	0.11ng/mL	0.31-20ng/mL
FY-B0091	Bovine MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	MCP-1	Sandwich	2000pg/mL	13.9 pg/mL	31.25-2000 pg/mL
FY-B0092	Bovine IL-8(Interleukin 8) ELISA Kit	IL-8	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0093	Bovine Hpc(Hepcidin) ELISA Kit	Hpc	Sandwich	2000pg/mL	11.7pg/mL	31.25-2000pg/mL
FY-B0094	Bovine LEP(Leptin) ELISA Kit	LEP	Sandwich	10ng/mL	0.1ng/mL	0.16-10ng/mL
FY-B0095	Bovine MMP-13(Matrix Metalloproteinase 13) ELISA Kit	MMP-13	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-B0096	Bovine MYO(Myoglobin) ELISA Kit	MYO	Sandwich	50ng/mL	0.31 ng/mL	0.79-50 ng/mL
FY-B0097	Bovine IGFBP-2(Insulin Like Growth Factor Binding Protein 2) ELISA Kit	IGFBP-2	Sandwich	60ng/mL	0.37ng/mL	0.94-60ng/mL
FY-B0098	Bovine FGF4(Fibroblast Growth Factor 4)ELISA kit	FGF4	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0099	Bovine COL3(Collagen Type III) ELISA Kit	COL3	Sandwich	30ng/mL	0.19ng/mL	0.47-30ng/mL
FY-B0100	Bovine C3(Complement Component 3) ELISA Kit	C3	Sandwich	200ng/mL	1.87ng/mL	3.12-200ng/mL
FY-B0101	Bovine MIγ/CXCL9(Monocyte Interferon Gamma Inducing Factor) ELISA Kit	MIγ/CXCL9	Sandwich	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-B0102	Bovine GDH/GLDH(Glutamate dehydrogenase)ELISA kit	GDH/GLDH	Sandwich	500ng/mL	4.69ng/mL	7.81-500ng/mL
FY-B0103	Bovine IgA(Immunoglobulin A) ELISA Kit	IgA	Sandwich	50ng/mL	0.37ng/mL	0.78-50ng/mL
FY-B0104	Bovine MMP-1(Matrix Metalloproteinase 1) ELISA Kit	MMP-1	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-B0105	Bovine MMP-9(Matrix Metalloproteinase 9) ELISA Kit	MMP-9	Sandwich	20ng/mL	0.11ng/mL	0.31-20ng/mL
FY-B0106	Bovine Hp(Haptoglobin)ELISA kit	Hp	Sandwich	1000ng/mL	9.38ng/mL	15.63-1000ng/mL
FY-B0107	Bovine GFAP(Glial Fibrillary Acidic Protein) ELISA Kit	GFAP	Sandwich	2000pg/mL	12.8pg/mL	31.25-2000pg/mL
FY-B0108	Bovine FETUB(Fetuin B) ELISA Kit	FETUB	Sandwich	2000ng/mL	11.7ng/mL	31.25-2000ng/mL
FY-B0109	Bovine IgE(Immunoglobulin E) ELISA Kit	IgE	Sandwich	4000ng/mL	28.5 ng/mL	62.5-4000 ng/mL

FY-B0110	Bovine FG(Fibrinogen) ELISA Kit	FG	Sandwich	2500ng/mL	17ng/mL	39.06-2500ng/mL
FY-B0111	Bovine LBP(Lipopolysaccharide Binding Protein) ELISA Kit	LBP	Sandwich	40ng/mL	0.29ng/mL	0.63-40ng/mL
FY-B0112	Bovine ICAM-1/CD54(intercellular adhesion molecule 1)ELISA Kit	ICAM-1/CD54	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-B0113	Bovine ANG1(Angiopoinetin 1)ELISA kit	ANG1	Sandwich	40ng/mL	0.38ng/mL	0.63-40ng/mL
FY-B0114	Bovine ADP/Acrp30(Adiponectin)ELISA Kit	ADP/Acrp30	Sandwich	50ng/mL	0.47ng/mL	0.78-50ng/mL
FY-B0115	Bovine TLR4(Toll-Like Receptor 4) ELISA Kit	TLR4	Sandwich	40ng/mL	0.28ng/mL	0.63-40ng/mL
FY-B0116	Bovine IFN- α (Interferon Alpha) ELISA Kit	IFN- α	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0117	Bovine NGAL(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	NGAL	Sandwich	100ng/mL	0.66 ng/mL	1.57-100 ng/mL
FY-B0118	Bovine AGP(Alpha-1-Acid Glycoprotein) ELISA Kit	AGP	Sandwich	1000ng/mL	9.38ng/mL	15.63-1000ng/mL
FY-B0119	Bovine CASP3(Caspase 3) ELISA Kit	CASP3	Sandwich	10ng/mL	0.1ng/mL	0.16-10ng/mL
FY-B0120	Bovine CRP(C-Reactive Protein) ELISA Kit	CRP	Sandwich	400ng/mL	2.43ng/mL	6.25-400ng/mL
FY-B0121	Bovine α FP(Alpha-Fetoprotein)ELISA Kit	α FP	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-B0122	Bovine IGF-1(Insulin-Like Growth Factor 1) ELISA Kit	IGF-1	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-B0123	Bovine IGFBP-3(Insulin-Like Growth Factor Binding Protein 3) ELISA Kit	IGFBP-3	Sandwich	40ng/mL	0.27ng/mL	0.63-40ng/mL
FY-B0124	Bovine IL-2(Interleukin 2) ELISA Kit	IL-2	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0125	Bovine IL-17A(Interleukin 17A)ELISA Kit	IL-17A	Sandwich	500pg/mL	4.69pg/mL	7.81-500pg/mL
FY-B0126	Bovine P I NP(Procollagen I N-Terminal Propeptide) ELISA Kit	P I NP	Sandwich	200ng/mL	1.19 ng/mL	3.13-200 ng/mL
FY-B0127	Bovine IGFBP-1(Insulin-Like Growth Factor Binding Protein 1)ELISA kit	IGFBP-1	Sandwich	400ng/mL	3.75ng/mL	6.25-400ng/mL
FY-B0128	Bovine IGFBP-5(Insulin-Like Growth Factor Binding Protein 5) ELISA Kit	IGFBP-5	Sandwich	4000pg/mL	24.4pg/mL	62.5-4000pg/mL
FY-B0129	Bovine IGFBP-6(Insulin-Like Growth Factor Binding Protein 6) ELISA Kit	IGFBP-6	Sandwich	160ng/mL	0.94ng/mL	2.5-160ng/mL
FY-B0130	Bovine APOB100(Apolipoprotein B100) ELISA Kit	APOB100	Sandwich	1200ng/mL	11.25ng/mL	18.75-1200ng/mL
FY-B0131	Bovine BMPR1A(Bone Morphogenetic Protein Receptor 1A) ELISA Kit	BMPR1A	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-B0132	Bovine PIINP(Procollagen II N-Terminal Propeptide)ELISA kit	PIINP	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-B0133	Bovine S100A1(S100 Calcium Binding Protein A1)ELISA kit	S100A1	Sandwich	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-B0134	Bovine TLR-2(Toll-Like Receptor 2)ELISA kit	TLR-2	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-B0135	Bovine MCP-2/CCL8(Monocyte Chemotactic Protein 2)ELISA kit	MCP-2/CCL8	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0136	Bovine RBP4(Retinol Binding Protein 4) ELISA Kit	RBP4	Sandwich	20ng/mL	0.14ng/mL	0.31-20ng/mL
FY-B0137	Bovine TGF- β 2(Transforming Growth Factor Beta 2)ELISA kit	TGF- β 2	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0138	Bovine PDGF-B(Platelet Derived Growth Factor Subunit B) ELISA Kit	PDGF-B	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0139	Bovine IgM(Immunoglobulin M) ELISA Kit	IgM	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-B0140	Bovine MPO(Myeloperoxidase) ELISA Kit	MPO	Sandwich	100ng/mL	0.62ng/mL	1.56-100ng/mL
FY-B0141	Bovine SAA(Serum Amyloid A) ELISA Kit	SAA	Sandwich	400ng/mL	2.72ng/mL	6.25-400ng/mL
FY-B0142	Bovine LTF(Lactoferrin) ELISA Kit	LTF	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-B0143	Bovine PTX 3/TSG-14(Pentraxin 3) ELISA Kit	PTX 3/TSG-14	Sandwich	20ng/mL	0.16 ng/mL	0.32-20 ng/mL
FY-B0144	Bovine IL-6(Interleukin 6) ELISA Kit	IL-6	Sandwich	500pg/mL	3.0pg/mL	7.81-500pg/mL
FY-B0145	Bovine IL-10(Interleukin 10) ELISA Kit	IL-10	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0146	Bovine ACE(Angiotensin I Converting Enzyme) ELISA Kit	ACE	Sandwich	400ng/mL	2.15ng/mL	6.25-400ng/mL
FY-B0147	Bovine PF4(Platelet Factor 4) ELISA Kit	PF4	Sandwich	400pg/mL	2.85pg/mL	6.25-400pg/mL
FY-B0148	Bovine Cys-C(Cystatin C) ELISA Kit	Cys-C	Sandwich	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-B0149	Bovine ANG(Angiogenin) ELISA Kit	ANG	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0150	Bovine GS(Gelsolin)ELISA kit	GS	Sandwich	500ng/mL	4.69ng/mL	7.81-500ng/mL

FY-B0151	Bovine TPO(Thrombopoietin) ELISA Kit	TPO	Sandwich	500pg/mL	3.1pg/mL	7.81-500pg/mL
FY-B0152	Bovine ALT(Alanine Aminotransferase)ELISA kit	ALT	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-B0153	Bovine HSP90aB1(Heat Shock Protein 90kDa Alpha B1)ELISA kit	HSP90aB1	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-B0154	Bovine NPY(Neuropeptide Y) ELISA Kit	NPY	Competitive	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-B0155	Bovine IgG(Immunoglobulin G) ELISA Kit	IgG	Competitive	100µg/mL	0.53µg/mL	1.56-100µg/mL
FY-B0156	Bovine AMH(Anti-Mullerian Hormone) ELISA Kit	AMH	Competitive	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-B0157	Bovine INS(Insulin) ELISA Kit	INS	Competitive	10000pg/mL	93.75pg/mL	156.25-10000pg/mL
FY-B0158	Bovine LH(Luteinizing Hormone) ELISA Kit	LH	Competitive	30000pg/mL	281.25pg/mL	468.75-30000pg/mL
FY-B0159	Bovine GHRL(Ghrelin) ELISA Kit	GHRL	Competitive	10000pg/mL	93.75pg/mL	156.25-10000pg/mL
FY-B0160	Bovine FSH(Follicle Stimulating Hormone) ELISA Kit	FSH	Competitive	240ng/mL	2.25ng/mL	3.75-240ng/mL
FY-B0161	Bovine ACTH(Adrenocorticotrophic Hormone) ELISA Kit	ACTH	Competitive	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0162	Bovine VIP(Vasoactive Intestinal Peptide) ELISA Kit	VIP	Competitive	500pg/mL	4.69pg/mL	7.81-500pg/mL
FY-B0163	Bovine bFGF/FGF2(Fibroblast Growth Factor 2, Basic) ELISA Kit	bFGF/FGF2	Competitive	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0164	Bovine HB(Hemoglobin) ELISA Kit	HB	Competitive	1000µg/mL	9.38µg/mL	15.63-1000µg/mL
FY-B0165	Bovine GH(Growth Hormone) ELISA Kit	GH	Competitive	10000pg/mL	93.8pg/mL	156.25-10000pg/mL
FY-B0166	Bovine KISS1(Kisspeptin 1) ELISA Kit	KISS1	Competitive	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0167	Bovine IGFBP-4(Insulin-Like Growth Factor Binding Protein 4)ELISA kit	IGFBP-4	Competitive	6000pg/mL	56.25pg/mL	93.75-6000pg/mL
FY-B0168	Bovine Ang-II(Angiotensin II)ELISA kit	Ang-II	Competitive	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-B0169	Bovine β-EP(Beta-Endorphin) ELISA Kit	β-EP	Competitive	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0170	Bovine I-PTH(intact Parathormone)ELISA Kit	I-PTH	Competitive	800pg/mL	7.5pg/mL	12.5-800pg/mL
FY-B0171	Bovine CTXI(Cross Linked C-Telopeptide Of Type I Collagen) ELISA Kit	CTXI	Competitive	10000pg/mL	52.9 pg/mL	156.25-10000 pg/mL
FY-B0172	Bovine CPP(Copeptin) ELISA Kit	CPP	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0173	Bovine GAL3(Galectin 3) ELISA Kit	GAL3	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-B0174	Bovine FABP3(Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit	FABP3	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-B0175	Bovine MR-ProADM(Mid-regional pro-adrenomedullin) ELISA Kit	MR-ProADM	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0176	Bovine NT-proBNP(N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit	NT-proBNP	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-B0177	Bovine TnTf/TNNT3(Troponin T, fast skeletal muscle) ELISA Kit	TnTf/TNNT3	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-B0178	Bovine TNNI3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit	TNNI3/cTn-I	Sandwich	20ng/mL	0.14ng/mL	0.31-20ng/mL
FY-SH0001	Sheep ANG2(Angiopoietin 2) ELISA Kit	ANG2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0002	Sheep BAFF/CD257(B-Cell Activating Factor) ELISA Kit	BAFF/CD257	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0003	Sheep BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	BMP-2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0004	Sheep BMP-7(Bone Morphogenetic Protein 7) ELISA Kit	BMP-7	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0005	Sheep CRN(Corin) ELISA Kit	CRN	Sandwich	2500pg/mL	23.44 pg/mL	39.06-2500pg/mL
FY-SH0006	Sheep CX3CL1(Chemokine C-X3-C-Motif Ligand 1) ELISA Kit	CX3CL1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0007	Sheep GROα/CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit	GROα/CXCL1	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0008	Sheep I-TAC/CXCL11(Interferon Inducible T-cell Alpha Chemoattractant) ELISA Ki	I-TAC/CXCL11	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0009	Sheep FAS/CD95(Factor Related Apoptosis) ELISA Kit	FAS/CD95	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0010	Sheep IL-23(Interleukin 23) ELISA Kit	IL-23	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0011	Sheep IFN-γ(Interferon Gamma) ELISA Kit	IFN-γ	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0012	Sheep TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	TNF-α	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-SH0013	Sheep Cys-A(Cystatin A) ELISA Kit	Cys-A	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0014	Sheep IL-1 β (Interleukin 1 Beta) ELISA Kit	IL-1 β	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0015	Sheep IL-12(Interleukin 12) ELISA Kit	IL-12	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0016	Sheep IL-12 p40(Interleukin 12 p40) ELISA Kit	IL-12 p40	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0017	Sheep MBP(Myelin Basic Protein) ELISA Kit	MBP	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0018	Sheep GH(Growth Hormone) ELISA Kit	GH	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0019	Sheep TNFRSF1A(Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit	TNFRSF1A	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-SH0020	Sheep IL-18(Interleukin 18) ELISA Kit	IL-18	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0021	Sheep ADAM9(A Disintegrin And Metalloprotease 9) ELISA Kit	ADAM9	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0022	Sheep FABP4(Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit	FABP4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-SH0023	Sheep ALDM(Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit	ALDM	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0024	Sheep AMBP(Alpha-1-Microglobulin/Bikunin Precursor) ELISA Kit	AMBP	Sandwich	150ng/mL	1.4 ng/mL	2.34-150ng/mL
FY-SH0025	Sheep ANG II R1(Angiotensin II Receptor 1) ELISA Kit	ANG II R1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0026	Sheep ANG II R2(Angiotensin II Receptor 2) ELISA Kit	ANG II R2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0027	Sheep CSF2R β (Colony Stimulating Factor 2 Receptor Beta) ELISA Kit	CSF2R β	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0028	Sheep CYP1A2(Cytochrome P450, family 1, subfamily A, polypeptide 2) ELISA Kit	CYP1A2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0029	Sheep CYP24A1(Cytochrome P450, family 24, subfamily A, polypeptide 1) ELISA Kit	CYP24A1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0030	Sheep CYP1B1(Cytochrome P450, family 1, subfamily B, polypeptide 1) ELISA Kit	CYP1B1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0031	Sheep PRELP(Proline Arginine Rich End Leucine Rich Repeat Protein) ELISA Kit	PRELP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0032	Sheep ARSB(Arylsulfatase B) ELISA Kit	ARSB	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0033	Sheep ApoE(Apolipoprotein E) ELISA Kit	ApoE	Sandwich	800ng/mL	7.03 ng/mL	12.5-800ng/mL
FY-SH0034	Sheep MANF(Mesencephalic Astrocyte Derived Neurotrophic Factor) ELISA Kit	MANF	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0035	Sheep ATA(Angiotensinase A) ELISA Kit	ATA	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0036	Sheep BAFRR(B-Cell Activation Factor Receptor) ELISA Kit	BAFRR	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0037	Sheep BCL2L1/BCL-X(Bcl-2 Like Protein 1) ELISA Kit	BCL2L1/BCL-X	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0038	Sheep Bcl-10(B-cell Lymphoma/Leukemia 10) ELISA Kit	Bcl-10	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0039	Sheep BACE2(Beta-Site APP Cleaving Enzyme 2) ELISA Kit	BACE2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0040	Sheep BRAK/CXCL14(Breast and Kidney Expressed Chemokine) ELISA Kit	BRAK/CXCL14	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0041	Sheep MUC1(Mucin 1) ELISA Kit	MUC1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0042	Sheep cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit	cTnT/TNNT2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0043	Sheep c-fos ELISA Kit	c-fos	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-SH0044	Sheep CKLFSF5(Chemokine Like Factor Superfamily 5) ELISA Kit	CKLFSF5	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0045	Sheep CD4(Cluster of Differentiation 4) ELISA Kit	CD4	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-SH0046	Sheep F13A1(Coagulation Factor X III A1 Polypeptide) ELISA Kit	F13A1	Sandwich	600pg/mL	5.63 pg/mL	9.38-600pg/mL
FY-SH0047	Sheep COL15(Collagen Type X V) ELISA Kit	COL15	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0048	Sheep CLDN4(Claudin 4) ELISA Kit	CLDN4	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0049	Sheep GATA4(GATA Binding Protein 4) ELISA Kit	GATA4	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0050	Sheep Preptin ELISA Kit	Preptin	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0051	Sheep Sp100(Sp100 Nuclear Antigen) ELISA Kit	Sp100	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0052	Sheep β Manase(β Mannosidase) ELISA Kit	β Manase	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0053	Sheep GP4/CD36(Platelet Membrane Glycoprotein IV) ELISA Kit	GP4/CD36	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL

FY-SH0054	Sheep GAL7(Galectin 7) ELISA Kit	GAL7	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-SH0055	Sheep ELN(Elastin) ELISA Kit	ELN	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0056	Sheep PKI α (Protein Kinase Inhibitor Alpha) ELISA Kit	PKI α	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0057	Sheep PRSS12(Protease, Serine 12) ELISA Kit	PRSS12	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0058	Sheep RETN(Resistin) ELISA Kit	RETN	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0059	Sheep UBE2C(Ubiquitin Conjugating Enzyme E2C) ELISA Kit	UBE2C	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0060	Sheep PAR4(Protease Activated Receptor 4) ELISA Kit	PAR4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0061	Sheep S100B(S100 Calcium Binding Protein B) ELISA Kit	S100B	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-SH0062	Sheep SDC1(Syndecan 1) ELISA Kit	SDC1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0063	Sheep SIAE(Sialic Acid Acetyltransferase) ELISA Kit	SIAE	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0064	Sheep Smad3/MADH3(Mothers Against Decapentaplegic Homolog 3) ELISA Kit	Smad3/MADH3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-SH0065	Sheep OPG(Osteoprotegerin) ELISA Kit	OPG	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0066	Sheep OC/BGP(Osteocalcin) ELISA Kit	OC/BGP	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-SH0067	Sheep ON(Osteonectin) ELISA Kit	ON	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0068	Sheep TXLNa(Taxilin Alpha) ELISA Kit	TXLNa	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0069	Sheep TXLNb(Taxilin Beta) ELISA Kit	TXLNb	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0070	Sheep CLC(Charcot Leyden Crystal Protein) ELISA Kit	CLC	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0071	Sheep OTR(Oxytocin Receptor) ELISA Kit	OTR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0072	Sheep Pro-MMP-13(Pro-Matrix Metalloproteinase 13) ELISA Kit	Pro-MMP-13	Sandwich	2500pg/mL	23.44 pg/mL	39.06-2500pg/mL
FY-SH0073	Sheep TSP-1(Thrombospondin-1) ELISA Kit	TSP-1	Sandwich	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-SH0074	Sheep PGAM2(Phosphoglycerate Mutase 2) ELISA Kit	PGAM2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0075	Sheep TrxR1(Thioredoxin Reductase 1) ELISA Kit	TrxR1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0076	Sheep Fc ϵ R I(Receptor I for the Fc Fragment of IgE) ELISA Kit	Fc ϵ R I	Sandwich	25ng/mL	0.23 ng/mL	0.39-25ng/mL
FY-SH0077	Sheep FAK(Focal Adhesion Kinase) ELISA Kit	FAK	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0078	Sheep GCK(Glucokinase) ELISA Kit	GCK	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-SH0079	Sheep GAP43(Growth Associated Protein 43) ELISA Kit	GAP43	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0080	Sheep PL(Placental Lactogen) ELISA Kit	PL	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-SH0081	Sheep PLAP/ALPP(Placental Alkaline Phosphatase) ELISA Kit	PLAP/ALPP	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0082	Sheep SYP(Synaptophysin) ELISA Kit	SYP	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-SH0083	Sheep TUB α 3C(Tubulin Alpha 3c) ELISA Kit	TUB α 3C	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0084	Sheep DCX(Doublecortin) ELISA Kit	DCX	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0085	Sheep PGC(Pepsinogen C) ELISA Kit	PGC	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-SH0086	Sheep CK-19/KRT19(Cytokeratin 19) ELISA Kit	CK-19/KRT19	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0087	Sheep FLN β (Filamin B, Beta) ELISA Kit	FLN β	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0088	Sheep STAT2(Signal Transducer and Activator of Transcription 2) ELISA Kit	STAT2	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0089	Sheep GP II b/IIIa(Platelet Membrane Glycoprotein II b/IIIa complex) ELISA Kit	GP II b/IIIa	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-SH0090	Sheep PHB(Prohibitin) ELISA Kit	PHB	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0091	Sheep IL-27(Interleukin 27) ELISA Kit	IL-27	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0092	Sheep GADL1(Glutamate Decarboxylase Like Protein 1) ELISA Kit	GADL1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0093	Sheep PARP4(Poly ADP Ribose Polymerase 4) ELISA Kit	PARP4	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0094	Sheep TNFRSF1B(Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit	TNFRSF1B	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL

FY-SH0095	Sheep $\text{T}\alpha\text{1}$ (Thymosin- α1) ELISA Kit	Ta1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-SH0096	Sheep IL-35(Interleukin 35) ELISA Kit	IL-35	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0097	Sheep LH(Luteinizing Hormone) ELISA Kit	LH	Sandwich	100mIU/mL	0.94 mIU/mL	1.56-100mIU/mL
FY-SH0098	Sheep C4(Complement Component 4) ELISA Kit	C4	Sandwich	150ng/mL	1.4 ng/mL	2.34-150ng/mL
FY-SH0099	Sheep KIM-1(Kidney Injury Molecule 1) ELISA Kit	KIM-1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0100	Sheep C3(Complement Component 3) ELISA Kit	C3	Sandwich	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-SH0101	Sheep IgA(Immunoglobulin A) ELISA Kit	IgA	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-SH0102	Sheep KISS1(Kisspeptin 1) ELISA Kit	KISS1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-SH0103	Sheep SHBG(Sex Hormone Binding Globulin) ELISA Kit	SHBG	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-SH0104	Sheep VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit	VE-Cadherin	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-SH0105	Sheep FG(Fibrinogen) ELISA Kit	FG	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-SH0106	Sheep ICAM-1/CD54(Intercellular Adhesion Molecule 1) ELISA Kit	ICAM-1/CD54	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0107	Sheep ADP/Acrp30(Adiponectin) ELISA Kit	ADP/Acrp30	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0108	Sheep IFN- α (Interferon Alpha) ELISA Kit	IFN- α	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0109	Sheep LCN2(Lipopocalin 2) ELISA Kit	LCN2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0110	Sheep FABP1(Fatty Acid Binding Protein 1, Liver) ELISA Kit	FABP1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-SH0111	Sheep CASP1(Caspase 1) ELISA Kit	CASP1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0112	Sheep APN(Aminopeptidase N) ELISA Kit	APN	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0113	Sheep IGF-1(Insulin-like Growth Factor 1) ELISA Kit	IGF-1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-SH0114	Sheep IL-2(Interleukin 2) ELISA Kit	IL-2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0115	Sheep IL-17A(Interleukin 17A) ELISA Kit	IL-17A	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0116	Sheep TnTf/TNNT3(Troponin T, fast skeletal muscle) ELISA Kit	TnTf/TNNT3	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0117	Sheep IGFBP-6(Insulin Like Growth Factor Binding Protein 6) ELISA Kit	IGFBP-6	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-SH0118	Sheep ITIH4(Inter-alpha-trypsin Inhibitor Heavy Chain H4) ELISA Kit	ITIH4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-SH0119	Sheep A β 1-40(Amyloid Beta 1-40) ELISA Kit	A β 1-40	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0120	Sheep A β 1-42(Amyloid Beta 1-42) ELISA Kit	A β 1-42	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0121	Sheep Bcl-3(B-cell Lymphoma/Leukemia 3) ELIS Kit ELISA Kit	Bcl-3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0122	Sheep BTC(Betacellulin) ELISA Kit	BTC	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0123	Sheep CGRP1(Calcitonin Gene Related Peptide 1) ELISA Kit	CGRP1	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0124	Sheep TNNI3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit	TNNI3/cTn-I	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-SH0125	Sheep CT1(Cardiotrophin 1) ELISA Kit	CT1	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0126	Sheep CLDN14(Claudin 14) ELISA Kit	CLDN14	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0127	Sheep C1qR1(Complement Component 1q Receptor 1) ELISA Kit	C1qR1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0128	Sheep CXCR3(CXC-Chemokine Receptor 3) ELISA Kit	CXCR3	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0129	Sheep MIgD(Membrane Immunoglobulin D) ELISA Kit	MIgD	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-SH0130	Sheep PRF1(Perforin 1) ELISA Kit	PRF1	Sandwich	25ng/mL	0.23 ng/mL	0.39-25ng/mL
FY-SH0131	Sheep FSH(Follicle Stimulating Hormone) ELISA Kit	FSH	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-SH0132	Sheep FDFT1(Farnesyl Diphosphate Farnesyltransferase 1) ELISA Kit	FDFT1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0133	Sheep RNASEH1(Ribonuclease H1) ELISA Kit	RNASEH1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-SH0134	Sheep τ PK1(Tau-Protein Kinase 1) ELISA Kit	τ PK1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0135	Sheep TAP1(Transporter associated with Antigen Processing 1) ELISA Kit	TAP1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL

FY-SH0136	Sheep HSP-70/HSPA9(Heat Shock 70 kDa Protein 9) ELISA Kit	HSP-70/HSPA9	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-SH0137	Sheep GROβ/CXCL2(Growth Regulated Oncogene Beta) ELISA Kit	GROβ/CXCL2	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0138	Sheep NE/ELA2(Neutrophil Elastase/Elastase-2) ELISA Kit	NE/ELA2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0139	Sheep CTLA4(Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit	CTLA4	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0140	Sheep CCND1(Cyclin D1) ELISA Kit	CCND1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0141	Sheep SAA(Serum Amyloid A) ELISA Kit	SAA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-SH0142	Sheep EPO(Erythropoietin) ELISA Kit	EPO	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0143	Sheep CP(Ceruloplasmin) ELISA Kit	CP	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-SH0144	Sheep IL6(Interleukin 6) ELISA Kit	IL6	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-SH0145	Sheep IL-10(Interleukin 10) ELISA Kit	IL-10	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0146	Sheep VEGFR1/FLT1(Vasculoar Endothelial Growth Factor Receptor 1) ELISA Kit	VEGFR1/FLT1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-SH0147	Sheep SELP(P-Selectin) ELISA Kit	SELP	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-SH0148	Sheep FRS3(Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit	FRS3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0149	Sheep AMH(Anti-Mullerian Hormone) ELISA Kit	AMH	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0150	Sheep βTG/PBP/CXCL7/NAP2(Thromboglobulin, Beta) ELISA Kit	βTG/PBP/CXCL7/NAP2	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0151	Sheep Gal(α-Galactosyl Epitopes) ELISA Kit	Gal	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0152	Sheep GC(Glucagon) ELISA Kit	GC	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0153	Sheep SOD2(Superoxide Dismutase 2, Mitochondrial) ELISA Kit	SOD2	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0154	Sheep GHRL(Ghrelin) ELISA Kit	GHRL	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-SH0155	Sheep β-EP(Beta-Endorphin) ELISA Kit	β-EP	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0156	Sheep DHEA-S(Dehydroepiandrosterone Sulfate) ELISA Kit	DHEA-S	Competitive	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-SH0157	Sheep SST(Somatostatin) ELISA Kit	SST	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0158	Sheep PGF2α(Prostaglandin F2α) ELISA Kit	PGF2α	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-SH0159	Sheep TMSβ4(Thymosin Beta 4) ELISA Kit	TMSβ4	Competitive	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-SH0160	Sheep ADH(Antidiuretic Hormone) ELISA Kit	ADH	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0161	Sheep VIP(Vasoactive Intestinal Peptide) ELISA Kit	VIP	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0162	Sheep MT(Melatonin) ELISA Kit	MT	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-SH0163	Sheep Ang- II (Angiotensin II) ELISA Kit	Ang- II	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-SH0164	Sheep AP13(Apelin 13) ELISA Kit	AP13	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-SH0165	Sheep Ang- I (Angiotensin I) ELISA Kit	Ang- I	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ca0001	Canine BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	BMP-2	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0002	Canine BMP-4(Bone Morphogenetic Protein 4) ELISA Kit	BMP-4	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0003	Canine MIP-1α(Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	MIP-1α	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0004	Canine MIP-1β(Macrophage Inflammatory Protein 1 Beta) ELISA Kit	MIP-1β	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0005	Canine IP-10/CXCL10(Interferon Gamma Induced Protein 10kDa) ELISA Kit	IP-10/CXCL10	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0006	Canine ET-1(Endothelin 1) ELISA Kit	ET-1	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0007	Canine HGF(Hepatocyte Growth Factor) ELISA Kit	HGF	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0008	Canine IL-1α(Interleukin 1 Alpha) ELISA Kit	IL-1α	Sandwich	500pg/mL	4.69pg/mL	7.81-500pg/mL
FY-Ca0009	Canine IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit	IL-1ra/IL-1F3	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0010	Canine IL-3(Interleukin 3) ELISA Kit	IL-3	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0011	Canine IL-4(Interleukin 4) ELISA Kit	IL-4	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL

FY-Ca0012	Canine IL-13(Interleukin 13) ELISA Kit	IL-13	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0013	Canine IFN- γ (Interferon Gamma) ELISA Kit	IFN- γ	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0014	Canine TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	TNF- α	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0015	Canine REN(Renin) ELISA Kit	REN	Sandwich	5000ng/mL	46.88ng/mL	78.13-5000ng/mL
FY-Ca0016	Canine APOA1(Apolipoprotein A1) ELISA Kit	APOA1	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-Ca0017	Canine CT(Calcitonin) ELISA Kit	CT	Competitive	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0018	Canine ACTH(Adrenocorticotrophic Hormone) ELISA Kit	ACTH	Competitive	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0019	Canine IL-1 β (Interleukin 1 Beta) ELISA Kit	IL-1 β	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0020	Canine IL-12(Interleukin 12) ELISA Kit	IL-12	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0021	Canine IL-12 p40(Interleukin 12 p40) ELISA Kit	IL-12 p40	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0022	Canine GH(Growth Hormone) ELISA Kit	GH	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0023	Canine TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	TIMP-1	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0024	Canine C5a(Complement Component 5a) ELISA Kit	C5a	Sandwich	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-Ca0025	Canine IL-5(Interleukin 5) ELISA Kit	IL-5	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0026	Canine IL-18(Interleukin 18) ELISA Kit	IL-18	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0027	Canine ADM(Adrenomedullin) ELISA Kit	ADM	Competitive	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0028	Canine APOE(Apolipoprotein E) ELISA Kit	APOE	Sandwich	300ng/mL	2.81ng/mL	4.69-300ng/mL
FY-Ca0029	Canine AST(Aspartate Aminotransferase) ELISA Kit	AST	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0030	Canine ANP(Atrial Natriuretic Peptide) ELISA Kit	ANP	Competitive	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0031	Canine cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit	cTnT/TNNT2	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0032	Canine COMP(Cartilage Oligomeric Matrix Protein) ELISA Kit	COMP	Sandwich	200ng/mL	1.88ng/mL	3.13-200ng/mL
FY-Ca0033	Canine CCKAR(Cholecystokinin A Receptor) ELISA Kit	CCKAR	Sandwich	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-Ca0034	Canine NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit	NOS2/iNOS	Sandwich	50ng/mL	0.47ng/mL	0.78-50ng/mL
FY-Ca0035	Canine COL1a2(Collagen Type I Alpha 2) ELISA Kit	COL1a2	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0036	Canine C3a(Complement Component 3a) ELISA Kit	C3a	Sandwich	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-Ca0037	Canine α 1-AT(Alpha 1-Antitrypsin) ELISA Kit	α 1-AT	Sandwich	800ng/mL	7.5ng/mL	12.5-800ng/mL
FY-Ca0038	Canine SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit	SOD1	Competitive	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-Ca0039	Canine ELN(Elastin) ELISA Kit	ELN	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0040	Canine sIgA(Secretory Immunoglobulin A) ELISA Kit	sIgA	Sandwich	40ng/mL	0.38ng/mL	0.63-40ng/mL
FY-Ca0041	Canine SCF(Stem Cell Factor) ELISA Kit	SCF	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0042	Canine PYGL(Glycogen Phosphorylase, Liver) ELISA Kit	PYGL	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0043	Canine OPG(Osteoprotegerin) ELISA Kit	OPG	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0044	Canine OC/BGP(Osteocalcin) ELISA Kit	OC/BGP	Sandwich	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-Ca0045	Canine OPN(Osteopontin) ELISA Kit	OPN	Sandwich	40ng/mL	0.38ng/mL	0.63-40ng/mL
FY-Ca0046	Canine PYGM(Glycogen Phosphorylase, Muscle) ELISA Kit	PYGM	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0047	Canine CKM(Creatine Kinase, Muscle) ELISA Kit	CKM	Sandwich	50ng/mL	0.47ng/mL	0.78.13-50ng/mL
FY-Ca0048	Canine DMD(Dystrophin) ELISA Kit	DMD	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0049	Canine MMP-3(Matrix Metalloproteinase 3) ELISA Kit	MMP-3	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0050	Canine PCT(Procalcitonin) ELISA Kit	PCT	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0051	Canine TRAIL/TNFSF10(Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit	TRAIL/TNFSF10	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0052	Canine GCK(Glucokinase) ELISA Kit	GCK	Sandwich	40ng/mL	0.38ng/mL	0.63-40ng/mL

FY-Ca0053	Canine GLUT4(Glucose Transporter 4) ELISA Kit	GLUT4	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0054	Canine GAP43(Growth Associated Protein 43) ELISA Kit	GAP43	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0055	Canine GSK3 α (Glycogen Synthase Kinase 3 Alpha) ELISA Kit	GSK3 α	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0056	Canine GT(Gastrin) ELISA Kit	GT	Competitive	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0057	Canine Plg(Plasminogen) ELISA Kit	Plg	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0058	Canine SOD2(Superoxide Dismutase 2, Mitochondrial) ELISA Kit	SOD2	Competitive	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0059	Canine VWF(Von Willebrand Factor) ELISA Kit	VWF	Competitive	1000ng/mL	9.38ng/mL	15.63-1000ng/mL
FY-LpCa0060	Canine GLP-2(Glucagon Like Peptide 2) ELISA Kit	GLP-2	Competitive	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0061	Canine TNF- β (Tumor Necrosis Factor Beta) ELISA Kit	TNF- β	Sandwich	500pg/mL	4.69pg/mL	7.81-500pg/mL
FY-Ca0062	Canine TTR(Transthyretin) ELISA Kit	TTR	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0063	Canine CALP(Calprotectin) ELISA Kit	CALP	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0064	Canine NES1(Nesfatin 1) ELISA Kit	NES1	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0065	Canine SOD3(Superoxide Dismutase 3, Extracellular) ELISA Kit	SOD3	Competitive	200ng/mL	1.88ng/mL	3.13-200ng/mL
FY-Ca0066	Canine CKB(Creatine Kinase, Brain) ELISA Kit	CKB	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0067	Canine T α 1(Thymosin- α 1) ELISA Kit	T α 1	Sandwich	200ng/mL	1.88ng/mL	3.13-200ng/mL
FY-Ca0068	Canine PACAP-38(Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit	PACAP-38	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0069	Canine PACAP-27(Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit	PACAP-27	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0070	Canine ADORA1(Adenosine A1 Receptor) ELISA Kit	ADORA1	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0071	Canine MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	MCP-1	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0072	Canine IL-8(Interleukin 8) ELISA Kit	IL-8	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0073	Canine Hpcp(Hpcidin) ELISA Kit	Hpcp	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0074	Canine LEP(Leptin) ELISA Kit	LEP	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0075	Canine MAU(Microalbuminuria) ELISA Kit	MAU	Sandwich	4000pg/mL	37.5pg/mL	62.5-4000pg/mL
FY-Ca0076	Canine MYO(Myoglobin) ELISA Kit	MYO	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-Ca0077	Canine TF(Transferrin) ELISA Kit	TF	Sandwich	400pg/mL	3.75ng/mL	6.25-400pg/mL
FY-Ca0078	Canine KIM-1(Kidney Injury Molecule 1) ELISA Kit	KIM-1	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0079	Canine C3(Complement Component 3) ELISA Kit	C3	Sandwich	200pg/mL	1.88ng/mL	3.13-200pg/mL
FY-Ca0080	Canine IgA(Immunoglobulin A) ELISA Kit	IgA	Sandwich	1000ng/mL	9.38ng/mL	15.63-1000ng/mL
FY-Ca0081	Canine MMP-9(Matrix Metalloproteinase 9) ELISA Kit	MMP-9	Sandwich	200ng/mL	1.88ng/mL	3.13-200ng/mL
FY-Ca0082	Canine I-PTH(intact Parathormone) ELISA Kit	I-PTH	Sandwich	800pg/mL	7.5pg/mL	12.5-800pg/mL
FY-Ca0083	Canine HP(Haptoglobin) ELISA Kit	HP	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-Ca0084	Canine MIP-3 β /ELC/CCL19(Macrophage Inflammatory Protein 3 β) ELISA Kit	MIP-3 β /ELC/CCL19	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0085	Canine IgE(Immunoglobulin E) ELISA Kit	IgE	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-Ca0086	Canine ALB(Albumin) ELISA Kit	ALB	Competitive	20000ng/mL	187.5ng/mL	312.5-20000ng/mL
FY-Ca0087	Canine FG(Fibrinogen) ELISA Kit	FG	Sandwich	4000pg/mL	37.5pg/mL	62.5-4000pg/mL
FY-Ca0088	Canine SAA2(Serum Amyloid A2) ELISA Kit	SAA2	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0089	Canine ICAM-1/CD54(intercellular adhesion molecule 1)ELISA Kit	ICAM-1/CD54	Sandwich	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-Ca0090	Canine F8(Coagulation Factor VIII (F8) ELISA Kit	F8	Sandwich	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-Ca0091	Canine ADP/Acrp30(Adiponectin) ELISA Kit	ADP/Acrp30	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-Ca0092	Canine IFN- α (Interferon Alpha) ELISA Kit	IFN- α	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0093	Canine BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	BDNF	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL

FY-Ca0094	Canine MIP-3 α (Macrophage Inflammatory Protein 3 Alpha) ELISA Kit	MIP-3 α	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0095	Canine CLU(Clusterin) ELISA Kit	CLU	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-Ca0096	Canine CRP(C-Reactive Protein) ELISA Kit	CRP	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0097	Canine α FP(Alpha-Fetoprotein) ELISA Kit	α FP	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0098	Canine GM-CSF(Granulocyte-Macrophage Colony Stimulating Factor) ELISA Kit	GM-CSF	Sandwich	500pg/mL	4.69pg/mL	7.81-500pg/mL
FY-Ca0099	Canine IGF-1(Insulin-Like Growth Factor 1) ELISA Kit	IGF-1	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-Ca0100	Canine IL-2(Interleukin 2) ELISA Kit	IL-2	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0101	Canine IL-17A(Interleukin 17A) ELISA Kit	IL-17A	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0102	Canine FE(Ferritin) ELISA Kit	FE	Sandwich	50pg/mL	0.47ng/mL	0.78-50pg/mL
FY-Ca0103	Canine IgG(Immunoglobulin G) ELISA Kit	IgG	Sandwich	200ng/mL	1.88ng/mL	3.13-200ng/mL
FY-Ca0104	Canine Ang- II (Angiotensin II) ELISA Kit	Ang- II	Competitive	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0105	Canine Ang- I (Angiotensin I) ELISA Kit	Ang- I	Competitive	5000pg/mL	46.88pg/mL	78.13-5000pg/mL
FY-Ca0106	Canine Ang-III(Angiotensin III) ELISA Kit	Ang-III	Sandwich	4000pg/mL	37.5pg/mL	62.5-4000pg/mL
FY-Ca0107	Canine TNNI3/cTn-I(Troponin I Type 3, Cardiac) ELISA Kit	TNNI3/cTn-I	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0108	Canine Cav-1(Caveolin-1) ELISA Kit	Cav-1	Sandwich	15ng/mL	0.14ng/mL	0.24-15ng/mL
FY-Ca0109	Canine SELE(E-selectin) ELISA Kit	SELE	Sandwich	80pg/mL	0.75ng/mL	1.25-80pg/mL
FY-Ca0110	Canine MCP-2/CCL8(Monocyte Chemotactic Protein 2)ELISA Kit	MCP-2/CCL8	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0111	Canine β MSH(beta-Melanocyte Stimulating Hormone) ELISA Kit	β MSH	Sandwich	2500pg/mL	23.44pg/mL	39.06-2500pg/mL
FY-Ca0112	Canine CKMB(Creatine Kinase MB Isoenzyme) ELISA Kit	CKMB	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0113	Canine HMGB-1(High Mobility Group Protein B1) ELISA Kit	HMGB-1	Sandwich	500pg/mL	4.69pg/mL	7.81-500pg/mL
FY-Ca0114	Canine TGF- β 2(Transforming Growth Factor Beta 2) ELISA Kit	TGF- β 2	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0115	Canine NE/ELA2(Neutrophil Elastase/Elastase-2) ELISA Kit	NE/ELA2	Sandwich	100ng/mL	0.94ng/mL	1.56-100ng/mL
FY-Ca0116	Canine PS1(Presenilin 1) ELISA Kit	PS1	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0117	Canine INS(Insulin) ELISA Kit	INS	Sandwich	200ng/mL	1.88ng/mL	3.13-200ng/mL
FY-Ca0118	Canine TSH(Thyroid Stimulating Hormone) ELISA Kit	TSH	Sandwich	50pg/mL	0.47ng/mL	0.78-50pg/mL
FY-Ca0119	Canine EPO(Erythropoietin) ELISA Kit	EPO	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0120	Canine LTF/LF(Lactoferrin) ELISA Kit	LTF/LF	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0121	Canine IL-31(Interleukin 31) ELISA Kit	IL-31	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0122	Canine NT-ProBNP(N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit	NT-ProBNP	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0123	Canine IL-6(Interleukin 6) ELISA Kit	IL-6	Sandwich	4000pg/mL	37.5pg/mL	62.5-4000pg/mL
FY-Ca0124	Canine IL-10(Interleukin 10) ELISA Kit	IL-10	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0125	Canine BNP(Brain Natriuretic Peptide) ELISA Kit	BNP	Competitive	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0126	Canine ENG(Endoglin) ELISA Kit	ENG	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0127	Canine RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	RANTES	Sandwich	2000pg/mL	18.75pg/mL	31.25-2000pg/mL
FY-Ca0128	Canine CD40L(Cluster of Differentiation 40 Ligand) ELISA Kit	CD40L	Sandwich	10ng/mL	0.094ng/mL	0.16-10ng/mL
FY-Ca0129	Canine Cys-C(Cystatin C) ELISA Kit	Cys-C	Sandwich	20ng/mL	0.19ng/mL	0.31-20ng/mL
FY-Ca0130	Canine TM(Thrombomodulin) ELISA Kit	TM	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-Ca0131	Canine TPO(Thrombopoietin) ELISA Kit	TPO	Sandwich	1000pg/mL	9.38pg/mL	15.63-1000pg/mL
FY-RB0001	Rabbit ANG2(Angiopoietin 2) ELISA Kit	ANG2	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB0002	Rabbit BAFF/CD257(B-Cell Activating Factor) ELISA Kit	BAFF/CD257	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB0003	Rabbit BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	BMP-2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL

FY-RB0004	Rabbit BMP-4(Bone Morphogenetic Protein 4) ELISA Kit	BMP-4	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0005	Rabbit BMP-7(Bone Morphogenetic Protein 7) ELISA Kit	BMP-7	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0006	Rabbit E-Cad(E-Cadherin) ELISA Kit	E-Cad	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0007	Rabbit CA9(Carbonic Anhydrase IX) ELISA Kit	CA9	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0008	Rabbit MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	MIP-1 α	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0009	Rabbit MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit	MIP-1 β	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0010	Rabbit MCP-3(Monocyte Chemotactic Protein 3) ELISA Kit	MCP-3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0011	Rabbit ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit	ECF/CCL11	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0012	Rabbit TARC(Thymus Activation Regulated Chemokine) ELISA Kit	TARC	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0013	Rabbit MPIF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	MPIF2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0014	Rabbit CTACK(Cutaneous T-cell Attracting Chemokine) ELISA Kit	CTACK	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0015	Rabbit Fc ϵ R II /CD23(Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit	Fc ϵ R II /CD23	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0016	Rabbit CD40L(Cluster of Differentiation 40 Ligand) ELISA Kit	CD40L	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0017	Rabbit CD163(Cluster of Differentiation 163) ELISA Kit	CD163	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0018	Rabbit GP39(Glycoprotein 39, Cartilage) ELISA Kit	GP39	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0019	Rabbit CNTF(Ciliary Neurotrophic Factor) ELISA Kit	CNTF	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0020	Rabbit TF(Tissue Factor) ELISA Kit	TF	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0021	Rabbit CRN(Corin) ELISA Kit	CRN	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0022	Rabbit CX3CL1(Chemokine C-X3-C-Motif Ligand 1) ELISA Kit	CX3CL1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0023	Rabbit GRO α /CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit	GRO α /CXCL1	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0024	Rabbit GCP-2(Granulocyte Chemotactic Protein 2) ELISA Kit	GCP-2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0025	Rabbit IP-10/CXCL10(Interferon Gamma Induced Protein 10kDa) ELISA Kit	IP-10/CXCL10	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0026	Rabbit I-TAC/CXCL11(Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit	I-TAC/CXCL11	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0027	Rabbit BLC(B-Lymphocyte Chemoattractant) ELISA Kit	BLC	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0028	Rabbit Cyt-C(Cytochrome C) ELISA Kit	Cyt-C	Sandwich	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-RB0029	Rabbit DKK1(Dickkopf Related Protein 1) ELISA Kit	DKK1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0030	Rabbit EGF(Epidermal Growth Factor) ELISA Kit	EGF	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0031	Rabbit EGFR(Epidermal Growth Factor Receptor) ELISA Kit	EGFR	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0032	Rabbit EG-VEGF(Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit	EG-VEGF	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-RB0033	Rabbit ET-1(Endothelin 1) ELISA Kit	ET-1	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0034	Rabbit EPCR(Endothelial Protein C Receptor) ELISA Kit	EPCR	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0035	Rabbit FAS/CD95(Factor Related Apoptosis) ELISA Kit	FAS/CD95	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0036	Rabbit FASL/TNFSF6(Factor Related Apoptosis Ligand) ELISA Kit	FASL/TNFSF6	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0037	Rabbit AFGF/FGF1(Acidic Fibroblast Growth Factor 1) ELISA Kit	AFGF/FGF1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0038	Rabbit FGF19(Fibroblast Growth Factor 19) ELISA Kit	FGF19	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0039	Rabbit FGF21(Fibroblast Growth Factor 21) ELISA Kit	FGF21	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0040	Rabbit FST(Follistatin) ELISA Kit	FST	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0041	Rabbit HGF(Hepatocyte Growth Factor) ELISA Kit	HGF	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0042	Rabbit IFN- β (Interferon Beta) ELISA Kit	IFN- β	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0043	Rabbit IL-1 α (Interleukin 1 Alpha) ELISA Kit	IL-1 α	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0044	Rabbit IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit	IL-1ra/IL-1F3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL

FY-RB0045	Rabbit IL-1R2(Interleukin 1 Receptor Type II) ELISA Kit	IL-1R2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0046	Rabbit LIF(Leukemia Inhibitory Factor) ELISA Kit	LIF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0047	Rabbit TNFSF14(Tumor Necrosis Factor Ligand Superfamily, Member 14) ELISA Kit	TNFSF14	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0048	Rabbit MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit	MCSF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0049	Rabbit Pro-MMP-1(pro-Matrix Metalloproteinase 1) ELISA Kit	Pro-MMP-1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0050	Rabbit IL-3(Interleukin 3) ELISA Kit	IL-3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0051	Rabbit IL-4(Interleukin 4) ELISA Kit	IL-4	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0052	Rabbit IL-13(Interleukin 13) ELISA Kit	IL-13	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0053	Rabbit IL-22(Interleukin 22) ELISA Kit	IL-22	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0054	Rabbit IL-23(Interleukin 23) ELISA Kit	IL-23	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0055	Rabbit IFN- γ (Interferon Gamma) ELISA Kit	IFN- γ	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0056	Rabbit TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	TNF- α	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0057	Rabbit EREG(Epiregulin) ELISA Kit	EREG	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0058	Rabbit Bcl-2(B-cell Leukemia/Lymphoma 2) ELISA Kit	Bcl-2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0059	Rabbit REN(Renin) ELISA Kit	REN	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0060	Rabbit CTSC(Cathepsin C) ELISA Kit	CTSC	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0061	Rabbit ApoA1(Apolipoprotein A1) ELISA Kit	ApoA1	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB0062	Rabbit MMP-25(Matrix Metalloproteinase 25) ELISA Kit	MMP-25	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0063	Rabbit MMP-12(Matrix Metalloproteinase 12) ELISA Kit	MMP-12	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0064	Rabbit CT(Calcitonin) ELISA Kit	CT	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0066	Rabbit PAPP(A) (Pregnancy Associated Plasma Protein A) ELISA Kit	PAPP(A)	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0067	Rabbit CSTA(Cystatin A) ELISA Kit	CSTA	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0068	Rabbit IL-1 β (Interleukin 1 Beta) ELISA Kit	IL-1 β	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0069	Rabbit IL-12(Interleukin 12) ELISA Kit	IL-12	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0070	Rabbit IL-12 p40(Interleukin 12 p40) ELISA Kit	IL-12 p40	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0072	Rabbit CD42(Cluster of Differentiation 42) ELISA Kit	CD42	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0073	Rabbit FGFR2(Fibroblast Growth Factor Receptor 2) ELISA Kit	FGFR2	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0074	Rabbit MBP(Myelin Basic Protein) ELISA Kit	MBP	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0075	Rabbit TE(Telomerase) ELISA Kit	TE	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0076	Rabbit TFPI(Tissue Factor Pathway Inhibitor) ELISA Kit	TFPI	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0077	Rabbit TM(Thrombin) ELISA Kit	TM	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0078	Rabbit FSTL1(Follistatin Like Protein 1) ELISA Kit	FSTL1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0079	Rabbit GH(Growth Hormone) ELISA Kit	GH	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0080	Rabbit COL4(Collagen Type IV) ELISA Kit	COL4	Sandwich	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-RB0081	Rabbit FN(Fibronectin) ELISA Kit	FN	Sandwich	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-RB0082	Rabbit IL-9(Interleukin 9) ELISA Kit	IL-9	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0083	Rabbit PC I (Procollagen I) ELISA Kit	PC I	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0084	Rabbit PCIII(Procollagen III) ELISA Kit	PCIII	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0085	Rabbit PIII(NP(N-Terminal Procollagen III Propeptide) ELISA Kit	PIIINP	Sandwich	60pg/mL	0.56 pg/mL	0.94-60pg/mL
FY-RB0086	Rabbit TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	TIMP-1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0087	Rabbit ADCY2(Adenylate Cyclase 2) ELISA Kit	ADCY2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL

FY-RB0088	Rabbit AGGF1(Angiogenic Factor with G Patch and FHA Domains 1) ELISA Kit	AGGF1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0089	Rabbit C4a(Complement Component 4a) ELISA Kit	C4a	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0090	Rabbit C5a(Complement Component 5a) ELISA Kit	C5a	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0091	Rabbit IL-5(Interleukin 5) ELISA Kit	IL-5	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0092	Rabbit IL-6R(Interleukin 6 Receptor) ELISA Kit	IL-6R	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0093	Rabbit NCAD(Neural Cadherin) ELISA Kit	NCAD	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0094	Rabbit PDGFRL(Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit	PDGFRL	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0095	Rabbit LUM(Lumican) ELISA Kit	LUM	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0096	Rabbit PSAP(Prosaposin) ELISA Kit	PSAP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0097	Rabbit PTGES(Prostaglandin E Synthase, Microsomal) ELISA Kit	PTGES	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0098	Rabbit UCRP(Ubiquitin Cross Reactive Protein) ELISA Kit	UCRP	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0100	Rabbit 15-LO(Arachidonate 15-Lipoxygenase) ELISA Kit	15-LO	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0101	Rabbit BACE1(Beta-site APP Cleaving Enzyme 1) ELISA Kit	BACE1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0102	Rabbit TNFRSF1A(Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit	TNFRSF1A	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0103	Rabbit IL-15(Interleukin 15) ELISA Kit	IL-15	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0104	Rabbit 5-LO(Arachidonate 5-Lipoxygenase) ELISA Kit	5-LO	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0105	Rabbit 5-NT(5-Nucleotidase) ELISA Kit	5-NT	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0106	Rabbit α 1ACT(Alpha1 Antichymotrypsin) ELISA Kit	α 1ACT	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0107	Rabbit IL-16(Interleukin 16) ELISA Kit	IL-16	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0108	Rabbit AR(Androgen Receptor) ELISA Kit	AR	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0109	Rabbit ARO(Aromatase) ELISA Kit	ARO	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0110	Rabbit BLVRB(Biliverdin Reductase B) ELISA Kit	BLVRB	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0111	Rabbit IL-18(Interleukin 18) ELISA Kit	IL-18	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0112	Rabbit IL-19(Interleukin 19) ELISA Kit	IL-19	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0113	Rabbit ACTg1(Actin Gamma 1) ELISA Kit	ACTg1	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0114	Rabbit ACTN2(Actinin Alpha 2) ELISA Kit	ACTN2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0115	Rabbit ACTN3(Actinin Alpha 3) ELISA Kit	ACTN3	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0116	Rabbit ADA(Adenosine Deaminase) ELISA Kit	ADA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0117	Rabbit ADAM10(A Disintegrin And Metalloprotease 10) ELISA Kit	ADAM10	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0118	Rabbit ADAM8(A Disintegrin And Metalloprotease 8) ELISA Kit	ADAM8	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0119	Rabbit ADAM9(A Disintegrin And Metalloprotease 9) ELISA Kit	ADAM9	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0120	Rabbit ADAMTS4(ADAM with Thrombospondin Type 1 Motif 4) ELISA Kit	ADAMTS4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0121	Rabbit ADCY1(Adenylate Cyclase 1, Brain) ELISA Kit	ADCY1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0122	Rabbit ADH1(Alcohol Dehydrogenase 1) ELISA Kit	ADH1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0124	Rabbit ADIPOR1(Adiponectin Receptor 1) ELISA Kit	ADIPOR1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0125	Rabbit ADIPOR2(Adiponectin Receptor 2) ELISA Kit	ADIPOR2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0126	Rabbit ADM(Adrenomedullin) ELISA Kit	ADM	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0127	Rabbit ADRA1A(Adrenergic Receptor alpha-1A) ELISA Kit	ADRA1A	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0128	Rabbit ADRP(Adipose Differentiation Related Protein) ELISA Kit	ADRP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0129	Rabbit IL-20(Interleukin 20) ELISA Kit	IL-20	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0130	Rabbit DLK-1(Delta Like 1 Homolog) ELISA Kit	DLK-1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-RB0131	Rabbit ACE2(Angiotensin I Converting Enzyme 2) ELISA Kit	ACE2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0132	Rabbit BMP-3(Bone Morphogenetic Protein 3) ELISA Kit	BMP-3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0133	Rabbit FABP4(Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit	FABP4	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0134	Rabbit AFP-L3(Alpha-Fetoprotein Lens Culinaris Agglutinin 3) ELISA Kit	AFP-L3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0135	Rabbit FUCA(Alpha-L-Fucosidase, Tissue) ELISA Kit	FUCA	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0136	Rabbit CCNB(Cyclin-B) ELISA Kit	CCNB	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0137	Rabbit AGC(Aggrecan) ELISA Kit	AGC	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0138	Rabbit AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA	AGER	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0139	Rabbit CD3d(Cluster of Differentiation 3d) ELISA Kit	CD3d	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0140	Rabbit AGR2(Anterior Gradient Protein 2) ELISA Kit	AGR2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0141	Rabbit AGRN(Agrin) ELISA Kit	AGRN	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0142	Rabbit AGT(Angiotensinogen) ELISA Kit	AGT	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0143	Rabbit AIF1(Allograft Inflammatory Factor 1) ELISA Kit	AIF1	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0144	Rabbit ALDM(Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit	ALDM	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0145	Rabbit Alk-Smase(Alkaline Sphingomyelinase) ELISA Kit	Alk-Smase	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0146	Rabbit CDK6(Cyclin Dependent Kinase 6) ELISA Kit	CDK6	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0147	Rabbit AMBP(Alpha-1-Microglobulin/Bikunin Precursor) ELISA Kit	AMBP	Sandwich	12ng/mL	0.11 ng/mL	0.19-12ng/mL
FY-RB0148	Rabbit CDK8(Cyclin Dependent Kinase 8) ELISA Kit	CDK8	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0149	Rabbit AMY1(Amylase Alpha 1, Salivary) ELISA Kit	AMY1	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0151	Rabbit CDKN3(Cyclin Dependent Kinase Inhibitor 3) ELISA Kit	CDKN3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0152	Rabbit ANG(Angiogenin) ELISA Kit	ANG	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0153	Rabbit ANG II R1(Angiotensin II Receptor 1) ELISA Kit	ANG II R1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0154	Rabbit ANG II R2/AGTR2(Angiotensin II Receptor 2) ELISA Kit	ANG II R2/AGTR2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0155	Rabbit ANGPT4(Angiopoietin 4) ELISA Kit	ANGPT4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0156	Rabbit ANGPTL1(Angiopoietin Like Protein 1) ELISA Kit	ANGPTL1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0157	Rabbit ANGPTL4(Angiopoietin Like Protein 4) ELISA Kit	ANGPTL4	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0158	Rabbit ANGPTL6(Angiopoietin Like Protein 6) ELISA Kit	ANGPTL6	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0159	Rabbit ANG-R-Tie1(Angiopoietin Receptor Tie1) ELISA Kit	ANG-R-Tie1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0160	Rabbit ANG-R-Tie2(Angiopoietin Receptor Tie2) ELISA Kit	ANG-R-Tie2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0161	Rabbit Ank-B(Ankyrin B) ELISA Kit	Ank-B	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0162	Rabbit MAG(Myelin Associated Glycoprotein) ELISA Kit	MAG	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0163	Rabbit CDNF(Cerebral Dopamine Neurotrophic Factor) ELISA Kit	CDNF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0164	Rabbit FGFR1(Fibroblast Growth Factor Receptor 1) ELISA Kit	FGFR1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0165	Rabbit CNN1(Calponin 1, Basic, Smooth Muscle) ELISA Kit	CNN1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0166	Rabbit CNTNAP2(Contactin Associated Protein Like Protein 2) ELISA Kit	CNTNAP2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0167	Rabbit ACHAP(Acetylcholinesterase Associated Protein) ELISA Kit	ACHAP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0168	Rabbit ACO1(Aconitase 1) ELISA Kit	ACO1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0169	Rabbit CPE(Carboxypeptidase E) ELISA Kit	CPE	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0170	Rabbit CGα(Chorionic Gonadotropin Alpha Polypeptide) ELISA Kit	CGα	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0171	Rabbit FGFR3(Fibroblast Growth Factor Receptor 3) ELISA Kit	FGFR3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0172	Rabbit DPP10(Dipeptidyl Peptidase X) ELISA Kit	DPP10	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL

FY-RB0173	Rabbit CPR(Cytochrome P450 Reductase) ELISA Kit	CPR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0174	Rabbit CRLF1(Cytokine Receptor Like Factor 1) ELISA Kit	CRLF1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0175	Rabbit CSF2Rβ(Colony Stimulating Factor 2 Receptor Beta) ELISA Kit	CSF2Rβ	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0176	Rabbit CA3(Carbonic Anhydrase III, Muscle Specific) ELISA Kit	CA3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0177	Rabbit CYP1A2(Cytochrome P450, family 1, subfamily A, polypeptide 2) ELISA Kit	CYP1A2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0178	Rabbit DCLK1(Doublecortin Like Kinase 1) ELISA Kit	DCLK1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0179	Rabbit DPYD(Dihydropyrimidine Dehydrogenase) ELISA Kit	DPYD	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0180	Rabbit DSE(Dermatan Sulfate Epimerase) ELISA Kit	DSE	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0181	Rabbit EPYC(Epiphycan) ELISA Kit	EPYC	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0182	Rabbit FADS2(Fatty Acid Desaturase 2) ELISA Kit	FADS2	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0183	Rabbit FAIM3(Fas Apoptotic Inhibitory Molecule 3) ELISA Kit	FAIM3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0184	Rabbit FASN(Fatty Acid Synthase) ELISA Kit	FASN	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0185	Rabbit FBLN4(Fibulin 4) ELISA Kit	FBLN4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0186	Rabbit CYP24A1(Cytochrome P450, family 24, subfamily A, polypeptide 1) ELISA Kit	CYP24A1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0187	Rabbit FCN1(Ficolin 1) ELISA Kit	FCN1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0188	Rabbit FcγBP(Fc Fragment of IgG Binding Protein) ELISA Kit	FcγBP	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0189	Rabbit FcγR3A(Fc Fragment of IgG Low Affinity IIIa Receptor) ELISA Kit	FcγR3A	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0190	Rabbit CYP1B1(Cytochrome P450, family 1, subfamily B, polypeptide 1) ELISA Kit	CYP1B1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0191	Rabbit FETUA(Fetuin A) ELISA Kit	FETUA	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0192	Rabbit FGy(Fibrinogen Gamma chain) ELISA Kit	FGy	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0193	Rabbit FLOT2(Flotillin 2) ELISA Kit	FLOT2	Sandwich	1600pg/mL	15 pg/mL	25-1600pg/mL
FY-RB0195	Rabbit GLRX3(Glutaredoxin 3) ELISA Kit	GLRX3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0196	Rabbit GIF(Gastric Intrinsic Factor) ELISA Kit	GIF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0197	Rabbit C4BPβ(C4 Binding Protein Beta) ELISA Kit	C4BPβ	Sandwich	1200pg/mL	11.25 pg/mL	18.75-1200pg/mL
FY-RB0198	Rabbit CA4(Carbonic Anhydrase IV) ELISA Kit	CA4	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0199	Rabbit C4BPα(C4 Binding Protein Alpha) ELISA Kit	C4BPα	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0200	Rabbit GPLD1(Glycosylphosphatidylinositol Specific Phospholipase D1) ELISA Kit	GPLD1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0201	Rabbit PRELP(Proline Arginine Rich End Leucine Rich Repeat Protein) ELISA Kit	PRELP	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0202	Rabbit TARS(Threonyl tRNA Synthetase, cytoplasmic) ELISA Kit	TARS	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0203	Rabbit GTF2A1(General Transcription Factor II A, Polypeptide 1) ELISA Kit	GTF2A1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0204	Rabbit GUCA2A(Guanylate Cyclase Activator 2A) ELISA Kit	GUCA2A	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0205	Rabbit INSL3(Insulin Like Protein 3) ELISA Kit	INSL3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0206	Rabbit INSL5(Insulin Like Protein 5) ELISA Kit	INSL5	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0207	Rabbit GUCA2B(Guanylate Cyclase Activator 2B) ELISA Kit	GUCA2B	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0208	Rabbit GZMK(Granzyme K) ELISA Kit	GZMK	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0209	Rabbit CYP26A1(Cytochrome P450, family 26, subfamily A, polypeptide 1) ELISA Kit	CYP26A1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0210	Rabbit HCCS(Holocytochrome C Synthase) ELISA Kit	HCCS	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0211	Rabbit CAMK2α(Calcium/Calmodulin-Dependent Protein Kinase II Alpha) ELISA Kit	CAMK2α	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0212	Rabbit BAG3(Bcl2 Associated Athanogene 3) ELISA Kit	BAG3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0213	Rabbit ARSB(Arylsulfatase B) ELISA Kit	ARSB	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0214	Rabbit ANXA5(Annexin A5) ELISA Kit	ANXA5	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL

FY-RB0215	Rabbit HTRA1(HtrA Serine, Peptidase 1) ELISA Kit	HTRA1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0216	Rabbit IGF1R(Insulin Like Growth Factor 1 Receptor) ELISA Kit	IGF1R	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0217	Rabbit HYOU1(Hypoxia Up Regulated 1) ELISA Kit	HYOU1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0218	Rabbit ICAM4(Intercellular Adhesion Molecule 4) ELISA Kit	ICAM4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0219	Rabbit IDE(Insulin Degrading Enzyme) ELISA Kit	IDE	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0220	Rabbit LIPH(Lipase H) ELISA Kit	LIPH	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0221	Rabbit LTβR(Lymphotoxin Beta Receptor) ELISA Kit	LTβR	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0222	Rabbit IGF2BP3(Insulin Like Growth Factor 2 mRNA Binding Protein 3) ELISA Kit	IGF2BP3	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0223	Rabbit TTN(Titin) ELISA Kit	TTN	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0224	Rabbit IGF2R(Insulin Like Growth Factor 2 Receptor) ELISA Kit	IGF2R	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0225	Rabbit IGFALS(Insulin Like Growth Factor Binding Protein, Acid Labile Subunit) ELISA Kit	IGFALS	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0226	Rabbit IGFBP-7(Insulin Like Growth Factor Binding Protein 7) ELISA Kit	IGFBP-7	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0227	Rabbit ANXA2(Annexin A2) ELISA Kit	ANXA2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0228	Rabbit INSR(Insulin Receptor) ELISA Kit	INSR	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0230	Rabbit LCN1(Lipocalin 1) ELISA Kit	LCN1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0231	Rabbit ApoA4(Apolipoprotein A4) ELISA Kit	ApoA4	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0232	Rabbit ApoB(Apolipoprotein B) ELISA Kit	ApoB	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0233	Rabbit ApoC1(Apolipoprotein C1) ELISA Kit	ApoC1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0234	Rabbit ApoC2(Apolipoprotein C2) ELISA Kit	ApoC2	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB0235	Rabbit ApoC4(Apolipoprotein C4) ELISA Kit	ApoC4	Sandwich	20pg/mL	0.19 pg/mL	0.31-20pg/mL
FY-RB0236	Rabbit ApoD(Apolipoprotein D) ELISA Kit	ApoD	Sandwich	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-RB0237	Rabbit ApoE(Apolipoprotein E) ELISA Kit	ApoE	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0238	Rabbit ApoH(Apolipoprotein H) ELISA Kit	ApoH	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0239	Rabbit ApoM(Apolipoprotein M) ELISA Kit	ApoM	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0240	Rabbit CA6(Carbonic Anhydrase VI) ELISA Kit	CA6	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0241	Rabbit AATF(Apoptosis Antagonizing Transcription Factor) ELISA Kit	AATF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0242	Rabbit LIPC(Lipase, Hepatic) ELISA Kit	LIPC	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0243	Rabbit AIF(Apoptosis Inducing Factor) ELISA Kit	AIF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0244	Rabbit LIPF(Lipase, Gastric) ELISA Kit	LIPF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0245	Rabbit APAF1(Apoptosis Protease Activating Factor 1) ELISA Kit	APAF1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0246	Rabbit ASK1(Apoptosis Signal Regulating Kinase 1) ELISA Kit	ASK1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0247	Rabbit AQP-1(Aquaporin 1) ELISA Kit	AQP-1	Sandwich	300pg/mL	2.81 pg/mL	4.69-300pg/mL
FY-RB0248	Rabbit AQP-2(Aquaporin 2) ELISA Kit	AQP-2	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0249	Rabbit AQP-3(Aquaporin 3) ELISA Kit	AQP-3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0250	Rabbit AQP-4(Aquaporin 4) ELISA Kit	AQP-4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0251	Rabbit AR(Amphiregulin) ELISA Kit	AR	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0252	Rabbit AVPR1A(Arginine Vasopressin Receptor 1A) ELISA Kit	AVPR1A	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0253	Rabbit KLK6(Kallikrein 6) ELISA Kit	KLK6	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0254	Rabbit KLK7(Kallikrein 7) ELISA Kit	KLK7	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0255	Rabbit ARHGD1α(Rho GDP Dissociation Inhibitor Alpha) ELISA Kit	ARHGD1α	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0256	Rabbit ARRβ1(Arrestin Beta 1) ELISA Kit	ARRβ1	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL

FY-RB0257	Rabbit ARRβ2(Arrestin Beta 2) ELISA Kit	ARRβ2	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0258	Rabbit MANF(Mesencephalic Astrocyte Derived Neurotrophic Factor) ELISA Kit	MANF	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0259	Rabbit MAP2K1(Mitogen Activated Protein Kinase Kinase 1) ELISA Kit	MAP2K1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0260	Rabbit AhR(Aryl Hydrocarbon Receptor) ELISA Kit	AhR	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0261	Rabbit ARSA(Arylsulfatase A) ELISA Kit	ARSA	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0262	Rabbit ASGR2(Asialoglycoprotein Receptor 2) ELISA Kit	ASGR2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0263	Rabbit ASGR1(Asialoglycoprotein Receptor 1) ELISA Kit	ASGR1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0264	Rabbit ASM(Acid Sphingomyelinase) ELISA Kit	ASM	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0265	Rabbit MAPK14(Mitogen Activated Protein Kinase 14) ELISA Kit	MAPK14	Sandwich	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-RB0266	Rabbit AST(Aspartate Aminotransferase) ELISA Kit	AST	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0267	Rabbit ASPN(Asporin) ELISA Kit	ASPN	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0268	Rabbit ATA(Angiotensinase A) ELISA Kit	ATA	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0269	Rabbit ATM(Ataxia Telangiectasia Mutated) ELISA Kit	ATM	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0270	Rabbit ATXN1(Ataxin 1) ELISA Kit	ATXN1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0271	Rabbit ATF6(Activating Transcription Factor 6) ELISA Kit	ATF6	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB0273	Rabbit AIRE(Autoimmune Regulator) ELISA Kit	AIRE	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0274	Rabbit AXIN2(Axis Inhibition Protein 2) ELISA Kit	AXIN2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0275	Rabbit AZU1(Azurocidin 1) ELISA Kit	AZU1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0276	Rabbit MCRS1(Microspherule Protein 1) ELISA Kit	MCRS1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0277	Rabbit BAFR(B-Cell Activation Factor Receptor) ELISA Kit	BAFR	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0278	Rabbit BCL2L1/BCL-X(Bcl-2 Like Protein 1) ELISA Kit	BCL2L1/BCL-X	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0279	Rabbit BPI(Bactericidal/Permeability Increasing Protein) ELISA Kit	BPI	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0280	Rabbit FGF6(Fibroblast Growth Factor 6) ELISA Kit	FGF6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0281	Rabbit LDHA(Lactate Dehydrogenase A) ELISA Kit	LDHA	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0282	Rabbit Bcl-10(B-cell Lymphoma/Leukemia 10) ELISA Kit	Bcl-10	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0283	Rabbit BAK1(Bcl2 Antagonist/Killer 1) ELISA Kit	BAK1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0284	Rabbit BAD(Bcl2 Associated Death Promoter) ELISA Kit	BAD	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0285	Rabbit Bax(Bcl-2 Associated X Protein) ELISA Kit	Bax	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB0286	Rabbit BMF(Bcl-2 Modifying Factor) ELISA Kit	BMF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0287	Rabbit BECN1(Beclin 1) ELISA Kit	BECN1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0288	Rabbit SEMA3B(Semaphorin 3B) ELISA Kit	SEMA3B	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0290	Rabbit BACE2(Beta-Site APP Cleaving Enzyme 2) ELISA Kit	BACE2	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB0291	Rabbit Bid(BH3 Interacting Domain Death Agonist) ELISA Kit	Bid	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0292	Rabbit MMRN1(Multimerin 1) ELISA Kit	MMRN1	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0293	Rabbit Bcl2L11(Bcl-2 Like Protein 11) ELISA Kit	Bcl2L11	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0294	Rabbit BMP-10(Bone Morphogenetic Protein 10) ELISA Kit	BMP-10	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0295	Rabbit BMPR2(Bone Morphogenetic Protein Receptor II) ELISA Kit	BMPR2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0296	Rabbit MOG(Myelin Oligodendrocyte Glycoprotein) ELISA Kit	MOG	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0297	Rabbit MT1E(Metallothionein 1E) ELISA Kit	MT1E	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0299	Rabbit MSMβ(Microseminoprotein Beta) ELISA Kit	MSMβ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0300	Rabbit BRAK/CXCL14(Breast and Kidney Expressed Chemokine) ELISA Kit	BRAK/CXCL14	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-RB0301	Rabbit BCAN(Brevican) ELISA Kit	BCAN	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0302	Rabbit MT1M(Metallothionein 1M) ELISA Kit	MT1M	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0303	Rabbit MT3(Metallothionein 3) ELISA Kit	MT3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0305	Rabbit MUC13(Mucin 13) ELISA Kit	MUC13	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0306	Rabbit CA724(Tumor Marker CA724) ELISA Kit	CA724	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0307	Rabbit CNR-1(Cadherin Related Neuronal Receptor1) ELISA Kit	CNR-1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0308	Rabbit MUC1(Mucin 1) ELISA Kit	MUC1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0309	Rabbit CABIN1(Calcineurin Binding Protein 1) ELISA Kit	CABIN1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0310	Rabbit CTR(Calcitonin Receptor) ELISA Kit	CTR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0311	Rabbit CACYBP(Calcyclin Binding Protein) ELISA Kit	CACYBP	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0312	Rabbit CALD(Caldesmon) ELISA Kit	CALD	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0313	Rabbit MSR1(Macrophage Scavenger Receptor 1) ELISA Kit	MSR1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0314	Rabbit CAM(Calmodulin) ELISA Kit	CAM	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB0315	Rabbit CNX(Calnexin) ELISA Kit	CNX	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0316	Rabbit CRT(Calreticulin) ELISA Kit	CRT	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0317	Rabbit CR(Calretinin) ELISA Kit	CR	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0318	Rabbit CALU(Calumenin) ELISA Kit	CALU	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0319	Rabbit MUC17(Mucin 17) ELISA Kit	MUC17	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0320	Rabbit MUC2(Mucin 2) ELISA Kit	MUC2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0321	Rabbit MUC20(Mucin 20) ELISA Kit	MUC20	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0322	Rabbit CNR1(Cannabinoid Receptor 1, Brain) ELISA Kit	CNR1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0323	Rabbit CA1(Carbonic Anhydrase I) ELISA Kit	CA1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0324	Rabbit CA2(Carbonic Anhydrase II) ELISA Kit	CA2	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB0325	Rabbit ANKRD1(Ankyrin Repeat Domain 1) ELISA Kit	ANKRD1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0326	Rabbit cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit	cTnT/TNNT2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0327	Rabbit CLS(Cardiolipin Synthase) ELISA Kit	CLS	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB0328	Rabbit CLCF1(Cardiotrophin Like Cytokine Factor 1) ELISA Kit	CLCF1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0329	Rabbit PAR1(Protease Activated Receptor 1) ELISA Kit	PAR1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0330	Rabbit CASP14(Caspase 14) ELISA Kit	CASP14	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0331	Rabbit CASP4(Caspase 4) ELISA Kit	CASP4	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB0332	Rabbit CASP7(Caspase 7) ELISA Kit	CASP7	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0333	Rabbit CAD(Caspase Activated Deoxyribonuclease) ELISA Kit	CAD	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0334	Rabbit CASP9(Caspase 9) ELISA Kit	CASP9	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0335	Rabbit CTNN α 1(Catenin, Alpha 1) ELISA Kit	CTNN α 1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0336	Rabbit CTNN β 1(Catenin, Beta 1) ELISA Kit	CTNN β 1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0337	Rabbit CTS(D) (Cathepsin D) ELISA Kit	CTSD	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0338	Rabbit CTSK(Cathepsin K) ELISA Kit	CTSK	Sandwich	100fmol/mL	0.94 fmol/mL	1.56-100fmol/mL
FY-RB0339	Rabbit CTSL(Cathepsin L) ELISA Kit	CTSL	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB0340	Rabbit CDX2(Caudal Type Homeobox Transcription Factor 2) ELISA Kit	CDX2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0341	Rabbit C/EBP α (CCAAT/Enhancer Binding Protein Alpha) ELISA Kit	C/EBP α	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0342	Rabbit C/EBP β (CCAAT/Enhancer Binding Protein Beta) ELISA Kit	C/EBP β	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-RB0343	Rabbit C/EBP δ (CCAAT/Enhancer Binding Protein, Delta) ELISA Kit	C/EBP δ	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0344	Rabbit C/EBP ϵ (CCAAT/Enhancer Binding Protein, Epsilon) ELISA Kit	C/EBP ϵ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0345	Rabbit CEBP γ (CCAAT/Enhancer Binding Protein Gamma) ELISA Kit	CEBP γ	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0346	Rabbit CD59(Protectin) ELISA Kit	CD59	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0347	Rabbit CENPA(Centromere Protein A) ELISA Kit	CENPA	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0348	Rabbit CA5B(Carbonic Anhydrase VB, Mitochondrial) ELISA Kit	CA5B	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0349	Rabbit MUC21(Mucin 21) ELISA Kit	MUC21	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0350	Rabbit CCL14(Chemokine C-C-Motif Ligand 14) ELISA Kit	CCL14	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0351	Rabbit CCL16(Chemokine C-C-Motif Ligand 16) ELISA Kit	CCL16	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0352	Rabbit CCL3L1(Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit	CCL3L1	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0353	Rabbit CCR8(Chemokine C-C-Motif Receptor 8) ELISA Kit	CCR8	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0354	Rabbit TERT(Telomerase Reverse Transcriptase) ELISA Kit	TERT	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0355	Rabbit CKLFSF4(Chemokine Like Factor Superfamily 4) ELISA Kit	CKLFSF4	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0356	Rabbit CKLFSF5(Chemokine Like Factor Superfamily 5) ELISA Kit	CKLFSF5	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0357	Rabbit CKLFSF6(Chemokine Like Factor Superfamily 6) ELISA Kit	CKLFSF6	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0358	Rabbit MUC6(Mucin 6) ELISA Kit	MUC6	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0359	Rabbit MYBPC3(Myosin Binding Protein C, Cardiac) ELISA Kit	MYBPC3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0360	Rabbit MYH10(Myosin Heavy Chain 10, Non Muscle) ELISA Kit	MYH10	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0361	Rabbit MYH9(Myosin Heavy Chain 9, Non Muscle) ELISA Kit	MYH9	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0362	Rabbit MYL9(Myosin Light Chain 9, Regulatory) ELISA Kit	MYL9	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0363	Rabbit CCKAR(Cholecystokinin A Receptor) ELISA Kit	CCKAR	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0364	Rabbit CETP(Cholesteryl Ester Transfer Protein) ELISA Kit	CETP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0365	Rabbit CENPF(Centromere Protein F) ELISA Kit	CENPF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0366	Rabbit CHRM2(Cholinergic Receptor, Muscarinic 2) ELISA Kit	CHRM2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0367	Rabbit CA12(Carbonic Anhydrase XII) ELISA Kit	CA12	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0368	Rabbit CBX3(Chromobox Homolog 3) ELISA Kit	CBX3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0369	Rabbit CHD5(Chromodomain Helicase DNA Binding Protein 5) ELISA Kit	CHD5	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0370	Rabbit NEF3(Neurofilament 3) ELISA Kit	NEF3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0371	Rabbit NEFL(Neurofilament, Light Polypeptide) ELISA Kit	NEFL	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0372	Rabbit NOS1/nNOS(Nitric Oxide Synthase 1, Neuronal) ELISA Kit	NOS1/nNOS	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0373	Rabbit C-jun(Proto-oncogene c-Jun) ELISA Kit	C-jun	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0374	Rabbit CC16(Clara Cell 16KD Protein) ELISA Kit	CC16	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0375	Rabbit CLDN1(Claudin 1) ELISA Kit	CLDN1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0376	Rabbit CD30(Cluster of Differentiation 30) ELISA Kit	CD30	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0377	Rabbit CD97(Cluster of Differentiation 97) ELISA Kit	CD97	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0378	Rabbit CD99(Cluster of Differentiation 99) ELISA Kit	CD99	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0379	Rabbit CD4(Cluster of Differentiation 4) ELISA Kit	CD4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0380	Rabbit NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit	NOS2/iNOS	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0381	Rabbit CLDN3(Claudin 3) ELISA Kit	CLDN3	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0382	Rabbit NOS3/eNOS(Nitric Oxide Synthase 3, Endothelial) ELISA Kit	NOS3/eNOS	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0383	Rabbit F7(Coagulation Factor VII) ELISA Kit	F7	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL

FY-RB0384	Rabbit F9(Coagulation Factor IX) ELISA Kit	F9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0385	Rabbit F10(Coagulation Factor X) ELISA Kit	F10	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-RB0386	Rabbit F13(Coagulation Factor X III) ELISA Kit	F13	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-RB0387	Rabbit F12(Coagulation Factor XII) ELISA Kit	F12	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0388	Rabbit F11(Coagulation Factor XI) ELISA Kit	F11	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0389	Rabbit F13A1(Coagulation Factor X III A1 Polypeptide) ELISA Kit	F13A1	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-RB0390	Rabbit F13B(Coagulation Factor X III B Polypeptide) ELISA Kit	F13B	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0391	Rabbit COL15(Collagen Type X V) ELISA Kit	COL15	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0392	Rabbit COL6(Collagen Type VI) ELISA Kit	COL6	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0393	Rabbit COL1 α 2(Collagen Type I Alpha 2) ELISA Kit	COL1 α 2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0394	Rabbit COL3 α 1(Collagen Type III Alpha 1) ELISA Kit	COL3 α 1	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB0395	Rabbit COL4 α 1(Collagen Type IV Alpha 1) ELISA Kit	COL4 α 1	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0396	Rabbit COL6 α 1(Collagen Type VI Alpha 1) ELISA Kit	COL6 α 1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0397	Rabbit COL7(Collagen Type VII) ELISA Kit	COL7	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0398	Rabbit COL8 α 1(Collagen Type VIII Alpha 1) ELISA Kit	COL8 α 1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0399	Rabbit COL11 α 1(Collagen Type XI Alpha 1) ELISA Kit	COL11 α 1	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-RB0400	Rabbit COL17(Collagen Type X VII) ELISA Kit	COL17	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0401	Rabbit NPYR1(Neuropeptide Y Receptor Y1) ELISA Kit	NPYR1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0402	Rabbit CLL1(Collectin Liver 1) ELISA Kit	CLL1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0403	Rabbit CLDN9(Claudin 9) ELISA Kit	CLDN9	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0404	Rabbit CLDN4(Claudin 4) ELISA Kit	CLDN4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0405	Rabbit C7(Complement Component 7) ELISA Kit	C7	Sandwich	30pg/mL	0.28 pg/mL	0.47-30pg/mL
FY-RB0406	Rabbit CFP(Complement Factor P) ELISA Kit	CFP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0407	Rabbit C3a(Complement Component 3a) ELISA Kit	C3a	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0408	Rabbit C3c(Complement Fragment 3c) ELISA Kit	C3c	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0409	Rabbit NPYR2(Neuropeptide Y Receptor Y2) ELISA Kit	NPYR2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0410	Rabbit C5b(Complement Fragment 5b) ELISA Kit	C5b	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0411	Rabbit NPYR5(Neuropeptide Y Receptor Y5) ELISA Kit	NPYR5	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0412	Rabbit CTGF(Connective Tissue Growth Factor) ELISA Kit	CTGF	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-RB0413	Rabbit CX26(Connexin 26) ELISA Kit	CX26	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0414	Rabbit CHGB(Chromogranin B) ELISA Kit	CHGB	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0415	Rabbit DEF α 1(Defensin Alpha 1, Neutrophil) ELISA Kit	DEF α 1	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0416	Rabbit CTX I (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit	CTX I	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0417	Rabbit NTXI(Cross Linked N-telopeptide of Type I Collagen) ELISA Kit	NTXI	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0418	Rabbit CTX- II (Cross Linked C-telopeptide of Type II Collagen) ELISA Kit	CTX- II	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0419	Rabbit LDHB(Lactate Dehydrogenase B) ELISA Kit	LDHB	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0420	Rabbit HMGCS1(3-Hydroxy-3-Methylglutaryl-Coenzyme A Synthase 1) ELISA Kit	HMGCS1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0421	Rabbit PIII CP(Procollagen III C-Terminal ProPeptide) ELISA Kit	PIII CP	Sandwich	30pg/mL	0.28 pg/mL	0.47-30pg/mL
FY-RB0422	Rabbit DKK4(Dickkopf Related Protein 4) ELISA Kit	DKK4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0423	Rabbit CPP(Copeptin) ELISA Kit	CPP	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0424	Rabbit CLDN11(Claudin 11) ELISA Kit	CLDN11	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-RB0425	Rabbit CXCR1(CXC-Chemokine Receptor 1) ELISA Kit	CXCR1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0426	Rabbit SCG2(Secretogranin II) ELISA Kit	SCG2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0427	Rabbit MYCBP(C-Myc Binding Protein) ELISA Kit	MYCBP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0428	Rabbit c-myc(c-myc Oncogene Product) ELISA Kit	c-myc	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0429	Rabbit CLEC4E(C-Type Lectin Domain Family 4 Member E) ELISA Kit	CLEC4E	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0430	Rabbit LDHC(Lactate Dehydrogenase C) ELISA Kit	LDHC	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB0431	Rabbit EphA2(Tyrosine protein kinase receptor A2) ELISA Kit	EphA2	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0432	Rabbit DKK2(Dickkopf Related Protein 2) ELISA Kit	DKK2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0433	Rabbit D-LDH(D-Lactate Dehydrogenase) ELISA Kit	D-LDH	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0434	Rabbit PA2G4(Proliferation Associated Protein 2G4) ELISA Kit	PA2G4	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0435	Rabbit COL1 α 1(Collagen Type I Alpha 1) ELISA Kit	COL1 α 1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0436	Rabbit COL6 α 3(Collagen Type VI Alpha 3) ELISA Kit	COL6 α 3	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0437	Rabbit COL13(Collagen Type X III) ELISA Kit	COL13	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0438	Rabbit COL14(Collagen Type XIV) ELISA Kit	COL14	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0439	Rabbit CTX-III(Cross Linked C-telopeptide of Type III Collagen) ELISA Kit	CTX-III	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0440	Rabbit EFNA1(Ephrin A1) ELISA Kit	EFNA1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0441	Rabbit GATA3(GATA Binding Protein 3) ELISA Kit	GATA3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0442	Rabbit GATA4(GATA Binding Protein 4) ELISA Kit	GATA4	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0443	Rabbit GNB1(G Protein Beta 1) ELISA Kit	GNB1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0444	Rabbit HHIP(Hedgehog Interacting Protein) ELISA Kit	HHIP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0445	Rabbit H-Cadherin(Cadherin, Heart) ELISA Kit	H-Cadherin	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0446	Rabbit NRG-4(Neuregulin 4) ELISA Kit	NRG-4	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0447	Rabbit PTPRK(Protein Tyrosine Phosphatase Receptor Type K) ELISA Kit	PTPRK	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0448	Rabbit SELL(L-Selectin) ELISA Kit	SELL	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0449	Rabbit MXD1(MAX Dimerization Protein 1) ELISA Kit	MXD1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0450	Rabbit HB μ (Hemoglobin, Mu) ELISA Kit	HB μ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0451	Rabbit NAGase(N-Acetyl Beta-D-Glucosaminidase) ELISA Kit	NAGase	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0452	Rabbit NRG-2(Neuregulin 2) ELISA Kit	NRG-2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0453	Rabbit NRG-3(Neuregulin 3) ELISA Kit	NRG-3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0454	Rabbit N-MID-OT(N-MID Osteocalcin) ELISA Kit	N-MID-OT	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0455	Rabbit PTPRN(Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit	PTPRN	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0456	Rabbit GAS(N-Acetylgalactosamine 6-Sulfatase) ELISA Kit	GAS	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0457	Rabbit PTPRO(Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit	PTPRO	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0458	Rabbit P27(P27 Protein) ELISA Kit	P27	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0459	Rabbit TP53(Tumor Protein p53) ELISA Kit	TP53	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0460	Rabbit PDAP1(Pdgfa Associated Protein 1) ELISA Kit	PDAP1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0461	Rabbit PDCN(Podocin) ELISA Kit	PDCN	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0462	Rabbit Preptin ELISA Kit	Preptin	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0463	Rabbit P-Cadherin(Cadherin, Placental) ELISA Kit	P-Cadherin	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0464	Rabbit RAPGEF1(Rap Guanine Nucleotide Exchange Factor 1) ELISA Kit	RAPGEF1	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0465	Rabbit Rac-GAP1(Rac-GTPase Activating Protein 1) ELISA Kit	Rac-GAP1	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL

FY-RB0466	Rabbit RHOB(Ras Homolog Gene Family, Memmber B) ELISA Kit	RHOB	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0467	Rabbit RCOR3(Rest Corepressor 3) ELISA Kit	RCOR3	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0468	Rabbit RHOA(Ras Homolog Gene Family, Memmber A) ELISA Kit	RHOA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0469	Rabbit Rock-2(Rho Associated Coiled Coil Containing Protein Kinase 2) ELISA Kit	Rock-2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0470	Rabbit RND1(Rho Family GTPase 1) ELISA Kit	RND1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0471	Rabbit RUNX2(Runt-Related Transcription Factor 2) ELISA Kit	RUNX2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0472	Rabbit COL9a3(Collagen Type IX Alpha 3) ELISA Kit	COL9a3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0474	Rabbit SHC/SLP-76(SH2-Containing Protein) ELISA Kit	SHC/SLP-76	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0475	Rabbit Slit2(Slit Homolog 2) ELISA Kit	Slit2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0476	Rabbit Smad1/MADH1(Mothers Against Decapentaplegic Homolog 1) ELISA Kit	Smad1/MADH1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0477	Rabbit Smad4/MADH4(Mothers Against Decapentaplegic Homolog 4) ELISA Kit	Smad4/MADH4	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0478	Rabbit Smad7/MADH7(Mothers Against Decapentaplegic Homolog 7) ELISA Kit	Smad7/MADH7	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0479	Rabbit Sp100(Sp100 Nuclear Antigen) ELISA Kit	Sp100	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0480	Rabbit PIAS1(Protein Inhibitor of Activated STAT 1) ELISA Kit	PIAS1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0481	Rabbit PIAS3(Protein Inhibitor of Activated STAT3) ELISA Kit	PIAS3	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0482	Rabbit PCDHβ16(Protocadherin Beta 16) ELISA Kit	PCDHβ16	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0483	Rabbit COL5a2(Collagen Type V Alpha 2) ELISA Kit	COL5a2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0484	Rabbit SAG(S Antigen) ELISA Kit	SAG	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0485	Rabbit TBP(TATA Binding Protein) ELISA Kit	TBP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0486	Rabbit MAPτ(Microtubule Associated Protein Tau/Tau Protein) ELISA Kit	MAPτ	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0487	Rabbit TIEG1(TGF Beta Inducible Early Response Gene 1) ELISA Kit	TIEG1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0488	Rabbit TLR9(Toll-Like Receptor 9) ELISA Kit	TLR9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0489	Rabbit PAR3(Protease Activated Receptor 3) ELISA Kit	PAR3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0491	Rabbit XCR1(XC-Chemokine Receptor 1) ELISA Kit	XCR1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0492	Rabbit PCDH15(Protocadherin 15) ELISA Kit	PCDH15	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0493	Rabbit α1BG(Alpha-1-B-Glycoprotein) ELISA Kit	α1BG	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0494	Rabbit α1-AT(Alpha 1-Antitrypsin) ELISA Kit	α1-AT	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-RB0495	Rabbit PLCγ2(Phospholipase C Gamma 2, Phosphatidylinositol Specific) ELISA Kit	PLCγ2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0496	Rabbit α2-AP(α2-Antiplasmin) ELISA Kit	α2-AP	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB0497	Rabbit IDUα(Alpha-L-Iduronidase) ELISA Kit	IDUα	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0498	Rabbit NAGLU(N-Acetyl Alpha-D-Glucosaminidase) ELISA Kit	NAGLU	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0499	Rabbit GLα(Galactosidase Alpha) ELISA Kit	GLα	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0500	Rabbit α Actin(Alpha-Actin) ELISA Kit	α Actin	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0501	Rabbit α-La(Alpha-Lactalbumin) ELISA Kit	α-La	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0502	Rabbit SNCα(Synuclein Alpha) ELISA Kit	SNCα	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0503	Rabbit PDGFRβ(Platelet Derived Growth Factor Receptor Beta) ELISA Kit	PDGFRβ	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0504	Rabbit PCDHγA2(Protocadherin Gamma A2) ELISA Kit	PCDHγA2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0505	Rabbit MAN1A1(Mannosidase Alpha Class 1A Member 1) ELISA Kit	MAN1A1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0506	Rabbit β-LPH(β-Lipotropic Hormone) ELISA Kit	β-LPH	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0507	Rabbit β Manase(β Mannosidase) ELISA Kit	β Manase	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0508	Rabbit CYSLTR1(CysteinyI Leukotriene Receptor 1) ELISA Kit	CYSLTR1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL

FY-RB0509	Rabbit IRAK3(Interleukin 1 Receptor Associated Kinase 3) ELISA Kit	IRAK3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0510	Rabbit GUSβ(β-glucuronidase) ELISA Kit	GUSβ	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0511	Rabbit SYCN(Syncollin) ELISA Kit	SYCN	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0512	Rabbit TUBβ(β-Tubulin) ELISA Kit	TUBβ	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0513	Rabbit γ-ECS(γ-Glutamyl Systeine Synthetase) ELISA Kit	γ-ECS	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-RB0514	Rabbit PLA2G10(Phospholipase A2, Group X) ELISA Kit	PLA2G10	Sandwich	20pg/mL	0.19 pg/mL	0.31-20pg/mL
FY-RB0515	Rabbit γGT1(Gamma Glutamyltransferase 1) ELISA Kit	γGT1	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0516	Rabbit POMC(Pro-Opiomelanocortin) ELISA Kit	POMC	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0518	Rabbit PIK3Cβ(Phosphoinositide-3-Kinase Catalytic Beta Polypeptide) ELISA Kit	PIK3Cβ	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0519	Rabbit OIT3(Oncoprotein Induced Transcript 3) ELISA Kit	OIT3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0520	Rabbit PIK3Cδ(Phosphoinositide-3-Kinase Catalytic Delta Polypeptide) ELISA Kit	PIK3Cδ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0521	Rabbit OCT2(Octamer-binding protein 2) ELISA Kit	OCT2	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0522	Rabbit OCT4(Octamer-binding protein 4) ELISA Kit	OCT4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0523	Rabbit PLA2G12B(Phospholipase A2, Group XII B) ELISA Kit	PLA2G12B	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0524	Rabbit NFIL3(Nuclear Factor, Interleukin 3 Regulated) ELISA Kit	NFIL3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0525	Rabbit LILRB4(Leukocyte Immunoglobulin Like Receptor Subfamily B Member 4) ELISA Kit	LILRB4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0526	Rabbit LIFR(Leukemia Inhibitory Factor Receptor) ELISA Kit	LIFR	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0527	Rabbit GARS(Glycyl tRNA Synthetase) ELISA Kit	GARS	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0528	Rabbit EPHA1(Ephrin Type A Receptor 1) ELISA Kit	EPHA1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0529	Rabbit SPOP(Speckle Type POZ Protein) ELISA Kit	SPOP	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0530	Rabbit GP4/CD36(Platelet Membrane Glycoprotein IV) ELISA Kit	GP4/CD36	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0531	Rabbit PLA2G4D(Phospholipase A2, Group IV D) ELISA Kit	PLA2G4D	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0532	Rabbit CYSLTR2(Cysteiny Leukotriene Receptor 2) ELISA Kit	CYSLTR2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0533	Rabbit GAL1(Galectin 1) ELISA Kit	GAL1	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB0534	Rabbit GAL12(Galectin 12) ELISA Kit	GAL12	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0535	Rabbit GAL14(Galectin 14) ELISA Kit	GAL14	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0536	Rabbit GAL2(Galectin 2) ELISA Kit	GAL2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0537	Rabbit EPHA3(Ephrin Type A Receptor 3) ELISA Kit	EPHA3	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0538	Rabbit GAL4(Galectin 4) ELISA Kit	GAL4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0539	Rabbit GAL8(Galectin 8) ELISA Kit	GAL8	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0540	Rabbit GAL7(Galectin 7) ELISA Kit	GAL7	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0541	Rabbit GAL9(Galectin 9) ELISA Kit	GAL9	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB0542	Rabbit UBC(Urinary Bladder Cancer Antigen) ELISA Kit	UBC	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0543	Rabbit PLA2R1(Phospholipase A2 Receptor 1) ELISA Kit	PLA2R1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0544	Rabbit PLB(Phospholipase B) ELISA Kit	PLB	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0546	Rabbit EFNA3(Ephrin A3) ELISA Kit	EFNA3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0547	Rabbit PLD2(Phospholipase D2) ELISA Kit	PLD2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0548	Rabbit PDXK(Pyridoxal Kinase) ELISA Kit	PDXK	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0549	Rabbit OCLN(Occludin) ELISA Kit	OCLN	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0550	Rabbit SEMA7A(Semaphorin 7A) ELISA Kit	SEMA7A	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0551	Rabbit PLVAP(Plasmalemma Vesicle Associated Protein) ELISA Kit	PLVAP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-RB0552	Rabbit SPJ(Surfactant Associated Protein J) ELISA Kit	SPJ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0553	Rabbit GPRC5A(G Protein Coupled Receptor, Family C, Group 5, Member A) ELISA Kit	GPRC5A	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0554	Rabbit FABP5(Fatty Acid Binding Protein 5, Epidermal) ELISA Kit	FABP5	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0555	Rabbit PDK4(Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit	PDK4	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0556	Rabbit PDP(Pyruvate Dehydrogenase Phosphatase) ELISA Kit	PDP	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0557	Rabbit PPOX(Protoporphyrinogen Oxidase) ELISA Kit	PPOX	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0558	Rabbit PRB1(Basic Salivary Proline Rich Protein 1) ELISA Kit	PRB1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0559	Rabbit FOXO1(Forkhead Box Protein O1) ELISA Kit	FOXO1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0560	Rabbit FOXO3(Forkhead Box Protein O3) ELISA Kit	FOXO3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0561	Rabbit FOXO4(Forkhead Box Protein O4) ELISA Kit	FOXO4	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0562	Rabbit FOXP3(Forkhead Box Protein P3) ELISA Kit	FOXP3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0563	Rabbit LATS(Long Acting Thyroid Stimulator) ELISA Kit	LATS	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0564	Rabbit EFNA4(Ephrin A4) ELISA Kit	EFNA4	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0565	Rabbit TFF3(Trefoil Factor 3, Intestinal) ELISA Kit	TFF3	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB0566	Rabbit TAGLN(Transgelin) ELISA Kit	TAGLN	Sandwich	16ng/mL	0.15 ng/mL	0.25-16ng/mL
FY-RB0567	Rabbit TEP1(Telomerase Associated Protein 1) ELISA Kit	TEP1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0568	Rabbit PRB2(Basic Salivary Proline Rich Protein 2) ELISA Kit	PRB2	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0571	Rabbit FGF23(Fibroblast Growth Factor 23) ELISA Kit	FGF23	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0572	Rabbit FGF9(Fibroblast Growth Factor 9) ELISA Kit	FGF9	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0573	Rabbit FGFBP1(Fibroblast Growth Factor Binding Protein 1) ELISA Kit	FGFBP1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0574	Rabbit FSHR(Follicle Stimulating Hormone Receptor) ELISA Kit	FSHR	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0575	Rabbit TGM1(Transglutaminase 1, Keratinocyte) ELISA Kit	TGM1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0576	Rabbit PACAPR(Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) ELISA Kit	PACAPR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0577	Rabbit EPHA10(Ephrin Type A Receptor 10) ELISA Kit	EPHA10	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0578	Rabbit ERβ(Estrogen Receptor Beta) ELISA Kit	ERβ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0579	Rabbit ERα(Estrogen Receptor Alpha) ELISA Kit	ERα	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0580	Rabbit EFNA5(Ephrin A5) ELISA Kit	EFNA5	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0581	Rabbit MAPKAPK3(Mitogen Activated Protein Kinase Activated Protein Kinase 3) ELISA Kit	MAPKAPK3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0582	Rabbit TEF(Thyrotroph Embryonic Factor) ELISA Kit	TEF	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0584	Rabbit GHRH(Growth Hormone Releasing Hormone) ELISA Kit	GHRH	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0586	Rabbit SCTR(Secretin Receptor) ELISA Kit	SCTR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0588	Rabbit FABP7(Fatty Acid Binding Protein 7, Brain) ELISA Kit	FABP7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0589	Rabbit MCP-4/CCL13(Monocyte Chemotactic Protein 4) ELISA Kit	MCP-4/CCL13	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0590	Rabbit GPER1(G Protein Coupled Estrogen Receptor 1) ELISA Kit	GPER1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0591	Rabbit ELN(Elastin) ELISA Kit	ELN	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB0592	Rabbit ELA4(Elastase 4) ELISA Kit	ELA4	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0593	Rabbit PROC(Protein C) ELISA Kit	PROC	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0594	Rabbit PDIA2(Protein Disulfide Isomerase A2) ELISA Kit	PDIA2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0595	Rabbit PDIA3(Protein Disulfide Isomerase A3) ELISA Kit	PDIA3	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB0596	Rabbit PRSS2(Protease, Serine, 2) ELISA Kit	PRSS2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0597	Rabbit PKBG(Protein Kinase B Gamma) ELISA Kit	PKBG	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL

FY-RB0598	Rabbit PICK1(Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit	PICK1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0599	Rabbit PKCβ1(Protein Kinase C Beta 1) ELISA Kit	PKCβ1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0600	Rabbit PKCγ(Protein Kinase C Gamma) ELISA Kit	PKCγ	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0601	Rabbit PKCδ1(Protein Kinase C Delta) ELISA Kit	PKCδ1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0602	Rabbit PKCε(Protein Kinase C Epsilon) ELISA Kit	PKCε	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0603	Rabbit PKCζ(Protein Kinase C Zeta) ELISA Kit	PKCζ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0604	Rabbit PKCθ(Protein Kinase C Theta) ELISA Kit	PKCθ	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0605	Rabbit PRSS23(Protease, Serine 23) ELISA Kit	PRSS23	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0606	Rabbit PKD1(Protein Kinase D1) ELISA Kit	PKD1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0607	Rabbit PKN1(Protein Kinase N1) ELISA Kit	PKN1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0608	Rabbit PKIα(Protein Kinase Inhibitor Alpha) ELISA Kit	PKIα	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0609	Rabbit PKIβ(Protein Kinase Inhibitor Beta) ELISA Kit	PKIβ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0610	Rabbit PKIγ(Protein Kinase Inhibitor Gamma) ELISA Kit	PKIγ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0611	Rabbit PTPRJ(Protein Tyrosine Phosphatase Receptor Type J) ELISA Kit	PTPRJ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0612	Rabbit PPP1R1A(Protein Phosphatase 1 Regulatory Subunit A) ELISA Kit	PPP1R1A	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0613	Rabbit PRSS1(Protease, Serine, 1) ELISA Kit	PRSS1	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB0614	Rabbit PRSS12(Protease, Serine 12) ELISA Kit	PRSS12	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0615	Rabbit PAR2(Protease Activated Receptor 2) ELISA Kit	PAR2	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB0616	Rabbit PCI(Protein C Inhibitor) ELISA Kit	PCI	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0617	Rabbit PDI(Protein Disulfide Isomerase) ELISA Kit	PDI	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0618	Rabbit PDIA5(Protein Disulfide Isomerase A5) ELISA Kit	PDIA5	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0619	Rabbit PRSS33(Protease, Serine 33) ELISA Kit	PRSS33	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0620	Rabbit LDL(Low Density Lipoprotein) ELISA Kit	LDL	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0621	Rabbit LDLR(Low Density Lipoprotein Receptor) ELISA Kit	LDLR	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0622	Rabbit LRP6(Low Density Lipoprotein Receptor Related Protein) ELISA Kit	LRP6	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0623	Rabbit RETN(Resistin) ELISA Kit	RETN	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0624	Rabbit RPS6Kβ1(Ribosomal Protein S6 Kinase Beta 1) ELISA Kit	RPS6Kβ1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0625	Rabbit ETFβ(Electron Transfer Flavoprotein Beta Polypeptide) ELISA Kit	ETFβ	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0626	Rabbit APP(Amyloid Precursor Protein) ELISA Kit	APP	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0627	Rabbit M-AChR M1(Muscarinic Acetylcholine Receptor M1) ELISA Kit	M-AChR M1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0628	Rabbit M-AChR M3(Muscarinic Acetylcholine Receptor M3) ELISA Kit	M-AChR M3	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0629	Rabbit D1R/DRD1(Dopamine Receptor D1) ELISA Kit	D1R/DRD1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0630	Rabbit D2R/DRD2(Dopamine Receptor D2) ELISA Kit	D2R/DRD2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0631	Rabbit D3R/DRD3(Dopamine Receptor D3) ELISA Kit	D3R/DRD3	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0632	Rabbit D4R/DRD4(Dopamine Receptor D4) ELISA Kit	D4R/DRD4	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0633	Rabbit D5R/DRD5(Dopamine Receptor D5) ELISA Kit	D5R/DRD5	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0634	Rabbit DDC(Dopamine Decarboxylase) ELISA Kit	DDC	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0635	Rabbit DAT(Dopamine Transporter) ELISA Kit	DAT	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0636	Rabbit DT(Dopachrome Tautomerase) ELISA Kit	DT	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0637	Rabbit PIGR/SC(Polymeric Immunoglobulin Receptor/Membrane Secretory Compon	PIGR/SC	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0638	Rabbit PTN(Pleiotrophin) ELISA Kit	PTN	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL

FY-RB0639	Rabbit DAO(Diamine Oxidase) ELISA Kit	DAO	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0640	Rabbit UGT2B4(Uridine Diphosphate Glucuronosyltransferase 2 Family Polypeptide) ELISA Kit	UGT2B4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0641	Rabbit DLD(Dihydrolipoyl Dehydrogenase) ELISA Kit	DLD	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0642	Rabbit DLAT(Dihydrolipoyl Transacetylase) ELISA Kit	DLAT	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0643	Rabbit DPYSL2(Dihydropyrimidinase Like Protein 2) ELISA Kit	DPYSL2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0644	Rabbit DNM1(Dynamin 1) ELISA Kit	DNM1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0645	Rabbit DNM2(Dynamin 2) ELISA Kit	DNM2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0646	Rabbit Ub(Ubiquitin) ELISA Kit	Ub	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0647	Rabbit UBD(Ubiquitin D) ELISA Kit	UBD	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0648	Rabbit UBE1L(Ubiquitin Activating Enzyme E1 Like Protein) ELISA Kit	UBE1L	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0649	Rabbit UBE2C(Ubiquitin Conjugating Enzyme E2C) ELISA Kit	UBE2C	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0650	Rabbit UBE3A(Ubiquitin Protein Ligase E3A) ELISA Kit	UBE3A	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0651	Rabbit NNE(Non-Neuronal Enolase) ELISA Kit	NNE	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0652	Rabbit SDC3(Syndecan 3) ELISA Kit	SDC3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0653	Rabbit FBLN5(Fibulin 5) ELISA Kit	FBLN5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0654	Rabbit REG1β(Regenerating Islet Derived Protein 1 Beta) ELISA Kit	REG1β	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0655	Rabbit SPC(Pulmonary Surfactant Associated Protein C) ELISA Kit	SPC	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB0656	Rabbit LRP(Lung Resistance-Related Protein) ELISA Kit	LRP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0657	Rabbit PRSS8(Protease, Serine 8) ELISA Kit	PRSS8	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0658	Rabbit SPARCL1(SPARC Like Protein 1) ELISA Kit	SPARCL1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0659	Rabbit SEMA3A(Semaphorin 3A) ELISA Kit	SEMA3A	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0660	Rabbit sIgA(Secretory Immunoglobulin A) ELISA Kit	sIgA	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB0661	Rabbit SLPI(Secretory Leukocyte Protease Inhibitor) ELISA Kit	SLPI	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0662	Rabbit SAP(Serum Amyloid P Component) ELISA Kit	SAP	Sandwich	16ng/mL	0.15 ng/mL	0.25-16ng/mL
FY-RB0663	Rabbit PAR4(Protease Activated Receptor 4) ELISA Kit	PAR4	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0664	Rabbit GPR120(G Protein Coupled Receptor 120) ELISA Kit	GPR120	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0665	Rabbit Cyr61(Cysteine Rich Protein, Angiogenic Inducer 61) ELISA Kit	Cyr61	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0666	Rabbit REG3γ(Regenerating Islet Derived Protein 3 Gamma) ELISA Kit	REG3γ	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0667	Rabbit S100A8(S100 Calcium Binding Protein A8) ELISA Kit	S100A8	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0668	Rabbit S100A9(S100 Calcium Binding Protein A9) ELISA Kit	S100A9	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0669	Rabbit GPR37(G Protein Coupled Receptor 37) ELISA Kit	GPR37	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0670	Rabbit S100A11(S100 Calcium Binding Protein A11) ELISA Kit	S100A11	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0671	Rabbit S100A12(S100 Calcium Binding Protein A12) ELISA Kit	S100A12	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0672	Rabbit S100A13(S100 Calcium Binding Protein A13) ELISA Kit	S100A13	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0673	Rabbit S100A6(S100 Calcium Binding Protein A6) ELISA Kit	S100A6	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0674	Rabbit S100A7(S100 Calcium Binding Protein A7) ELISA Kit	S100A7	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0675	Rabbit S100B(S100 Calcium Binding Protein B) ELISA Kit	S100B	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0676	Rabbit SDC1(Syndecan 1) ELISA Kit	SDC1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0677	Rabbit SCF(Stem Cell Factor) ELISA Kit	SCF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0678	Rabbit SCFR(Stem Cell Growth Factor Receptor) ELISA Kit	SCFR	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0679	Rabbit GAL(Galanin) ELISA Kit	GAL	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-RB0680	Rabbit GALR1(Galanin Receptor 1) ELISA Kit	GALR1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0681	Rabbit PK(Pyruvate Kinase, Liver And RBC) ELISA Kit	PK	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0682	Rabbit PYGL(Glycogen Phosphorylase, Liver) ELISA Kit	PYGL	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0683	Rabbit SFRP2(Secreted Frizzled Related Protein 2) ELISA Kit	SFRP2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0684	Rabbit GP-73(Golgi Protein 73) ELISA Kit	GP-73	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0685	Rabbit SIAE(Sialic Acid Acetylerase) ELISA Kit	SIAE	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0686	Rabbit ST2(Syntenin 2) ELISA Kit	ST2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0687	Rabbit SIGLEC1(Sialic Acid Binding Immunoglobulin Like Lectin 1) ELISA Kit	SIGLEC1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0688	Rabbit SIGLEC10(Sialic Acid Binding Immunoglobulin Like Lectin 10) ELISA Kit	SIGLEC10	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0689	Rabbit SIGLEC2(Sialic Acid Binding Immunoglobulin Like Lectin 2) ELISA Kit	SIGLEC2	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0690	Rabbit PADI4(Peptidyl Arginine Deiminase Type IV) ELISA Kit	PADI4	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0691	Rabbit SIGLEC8(Sialic Acid Binding Immunoglobulin Like Lectin 8) ELISA Kit	SIGLEC8	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0693	Rabbit GCLM(Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit	GCLM	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0694	Rabbit GSTθ1/GSTt1(Glutathione S Transferase Theta 1) ELISA Kit	GSTθ1/GSTt1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0695	Rabbit GSTθ2/GSTt2(Glutathione S Transferase Theta 2) ELISA Kit	GSTθ2/GSTt2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0696	Rabbit GSTω1(Glutathione S Transferase Omega 1) ELISA Kit	GSTω1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0697	Rabbit GSTpi(Glutathione S Transferases Pi) ELISA Kit	GSTpi	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0698	Rabbit OPG(Osteoprotegerin) ELISA Kit	OPG	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0699	Rabbit OC/BGP(Osteocalcin) ELISA Kit	OC/BGP	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0700	Rabbit ELA3B(Elastase 3B) ELISA Kit	ELA3B	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0701	Rabbit MUSK(Muscle Skeletal Receptor Tyrosine Kinase) ELISA Kit	MUSK	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0702	Rabbit GSTα5(Glutathione S Transferase Alpha 5) ELISA Kit	GSTα5	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0703	Rabbit OPN(Osteopontin) ELISA Kit	OPN	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0704	Rabbit MPIF-1/CCL23(Myeloid Progenitor Inhibitory Factor 1) ELISA Kit	MPIF-1/CCL23	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0705	Rabbit SOX9(Sex Determining Region Y Box Protein 9) ELISA Kit	SOX9	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0706	Rabbit GREM1(Gremlin 1) ELISA Kit	GREM1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0707	Rabbit ON(Osteonectin) ELISA Kit	ON	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0708	Rabbit SF1(Steroidogenic Factor 1) ELISA Kit	SF1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0709	Rabbit PAPP2(Pappalysin 2) ELISA Kit	PAPP2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0710	Rabbit P-CK(Pan-Cytokeratin) ELISA Kit	P-CK	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0711	Rabbit ALDOC(Aldolase C, Fructose Bisphosphate) ELISA Kit	ALDOC	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0712	Rabbit PMP3(Peroxisomal Membrane Protein 3) ELISA Kit	PMP3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0713	Rabbit PPARγC1α(Peroxisome Proliferator Activated Receptor Gamma Coactivator 1) ELISA Kit	PPARγC1α	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0714	Rabbit PPAR-α(Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit	PPAR-α	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0715	Rabbit PPAR-γ(Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit	PPAR-γ	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0716	Rabbit NXF1(Nuclear RNA Export Factor 1) ELISA Kit	NXF1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0717	Rabbit NMP-22(Nuclear Matrix Protein 22) ELISA Kit	NMP-22	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0718	Rabbit NUP107(Nucleoporin 107kDa) ELISA Kit	NUP107	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0719	Rabbit NUP133(Nucleoporin 133kDa) ELISA Kit	NUP133	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0720	Rabbit NUP155(Nucleoporin 155kDa) ELISA Kit	NUP155	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0721	Rabbit NUP188(Nucleoporin 188kDa) ELISA Kit	NUP188	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL

FY-RB0722	Rabbit NUP205(Nucleoporin 205kDa) ELISA Kit	NUP205	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0723	Rabbit NUP214(Nucleoporin 214kDa) ELISA Kit	NUP214	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0724	Rabbit NUP35(Nucleoporin 35kDa) ELISA Kit	NUP35	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0725	Rabbit NUP37(Nucleoporin 37kDa) ELISA Kit	NUP37	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0726	Rabbit NUP50(Nucleoporin 50kDa) ELISA Kit	NUP50	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0727	Rabbit NUP85(Nucleoporin 85kDa) ELISA Kit	NUP85	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0728	Rabbit NUP98(Nucleoporin 98kDa) ELISA Kit	NUP98	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0729	Rabbit NRIP1(Nuclear Receptor Interacting Protein 1) ELISA Kit	NRIP1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0730	Rabbit RI(Ribonuclease Inhibitor) ELISA Kit	RI	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0731	Rabbit RNASE(Ribonuclease) ELISA Kit	RNASE	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0732	Rabbit RPL10(Ribosomal Protein L10) ELISA Kit	RPL10	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0733	Rabbit RPL26(Ribosomal Protein L26) ELISA Kit	RPL26	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0734	Rabbit MART/Melan-A(Melanoma Marker A) ELISA Kit	MART/Melan-A	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0735	Rabbit MAGE(Melanoma-associated Antigen) ELISA Kit	MAGE	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0736	Rabbit MLZE(Melanoma Derived Leucine Zipper Extra Nuclear Factor) ELISA Kit	MLZE	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0737	Rabbit THRSP(Thyroid Hormone Responsive) ELISA Kit	THRSP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0738	Rabbit EPOR(Erythropoietin Receptor) ELISA Kit	EPOR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0739	Rabbit EXT1(Exostoses 1) ELISA Kit	EXT1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0740	Rabbit GSTκ1(Glutathione S Transferase Kappa 1) ELISA Kit	GSTκ1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0741	Rabbit LTA4H(Leukotriene A4 Hydrolase) ELISA Kit	LTA4H	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0742	Rabbit GSTα4(Glutathione S Transferase Alpha 4) ELISA Kit	GSTα4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0743	Rabbit SNCγ(Synuclein Gamma) ELISA Kit	SNCγ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0744	Rabbit TXLNa(Taxilin Alpha) ELISA Kit	TXLNa	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0745	Rabbit TXLNb(Taxilin Beta) ELISA Kit	TXLNb	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0746	Rabbit RNF4(Ring Finger Protein 4) ELISA Kit	RNF4	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0747	Rabbit FABP6(Fatty Acid Binding Protein 6, Ileal) ELISA Kit	FABP6	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0748	Rabbit TMPRSS2(Transmembrane Protease, Serine 2) ELISA Kit	TMPRSS2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0749	Rabbit TNNI1(Troponin I Type 1, Slow Skeletal) ELISA Kit	TNNI1	Sandwich	120pg/mL	1.13 pg/mL	1.88-120pg/mL
FY-RB0750	Rabbit FSCN(Fascin) ELISA Kit	FSCN	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0751	Rabbit TNNC2(Troponin C Type 2) ELISA Kit	TNNC2	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0752	Rabbit TNC(Tenascin C) ELISA Kit	TNC	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0753	Rabbit TNR(Tenascin R) ELISA Kit	TNR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0754	Rabbit TNX(Tenascin X) ELISA Kit	TNX	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0755	Rabbit MYL1(Myosin Light Chain 1, Alkali, Fast Skeletal) ELISA Kit	MYL1	Sandwich	300pg/mL	2.81 pg/mL	4.69-300pg/mL
FY-RB0756	Rabbit MLCK/MYLK(Myosin Light Chain Kinase) ELISA Kit	MLCK/MYLK	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0757	Rabbit MYH1(Myosin Heavy Chain 1, Skeletal Muscle, Adult) ELISA Kit	MYH1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0758	Rabbit FABP3(Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit	FABP3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0759	Rabbit PYGM(Glycogen Phosphorylase, Muscle) ELISA Kit	PYGM	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0760	Rabbit MYOG(Myogenin) ELISA Kit	MYOG	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0761	Rabbit DMD(Dystrophin) ELISA Kit	DMD	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0762	Rabbit MSTN(Myostatin) ELISA Kit	MSTN	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-RB0763	Rabbit PUMA(P53 Upregulated Modulator of Apoptosis) ELISA Kit	PUMA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0764	Rabbit MGP(Matrix Gla Protein) ELISA Kit	MGP	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0765	Rabbit MMP-11(Matrix Metalloproteinase 11) ELISA Kit	MMP-11	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0766	Rabbit MGST1(Microsomal Glutathione S Transferase 1) ELISA Kit	MGST1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0767	Rabbit MMP-2(Matrix Metalloproteinase 2) ELISA Kit	MMP-2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0768	Rabbit MMP-3(Matrix Metalloproteinase 3) ELISA Kit	MMP-3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0769	Rabbit MMP-17(Matrix Metalloproteinase 17) ELISA Kit	MMP-17	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0770	Rabbit MMP-24(Matrix Metalloproteinase 24) ELISA Kit	MMP-24	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0771	Rabbit MMP-7(Matrix Metalloproteinase 7) ELISA Kit	MMP-7	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0772	Rabbit MMP-8(Matrix Metalloproteinase 8) ELISA Kit	MMP-8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0773	Rabbit ELA2B(Elastase 2B) ELISA Kit	ELA2B	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0774	Rabbit TIMP-2(Tissue Inhibitors of Metalloproteinase 2) ELISA Kit	TIMP-2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0775	Rabbit TIMP-3(Tissue Inhibitors of Metalloproteinase 3) ELISA Kit	TIMP-3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0776	Rabbit TIMP-4(Tissue Inhibitors of Metalloproteinase 4) ELISA Kit	TIMP-4	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB0777	Rabbit LGALS3BP(Lectin Galactoside Binding, Soluble 3 Binding Protein) ELISA Kit	LGALS3BP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0778	Rabbit TNNI2(Troponin I Type 2, Fast Skeletal) ELISA Kit	TNNI2	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0779	Rabbit TP/TXA2R(Thromboxane Receptor) ELISA Kit	TP/TXA2R	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0780	Rabbit SDF-2(Stromal Cell Derived Factor 2) ELISA Kit	SDF-2	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0781	Rabbit SDF2L1(Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit	SDF2L1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0782	Rabbit SDF-4(Stromal Cell Derived Factor 4) ELISA Kit	SDF-4	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0783	Rabbit TDGF1(Teratocarcinoma Derived Growth Factor 1) ELISA Kit	TDGF1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0784	Rabbit KLK8(Kallikrein 8) ELISA Kit	KLK8	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0785	Rabbit VLDLR(Very Low Density Lipoprotein Receptor) ELISA Kit	VLDLR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0786	Rabbit GAL3(Galectin 3) ELISA Kit	GAL3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0787	Rabbit FPR2(Formyl Peptide Receptor 2) ELISA Kit	FPR2	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0788	Rabbit PTHR2(Parathyroid Hormone Receptor 2) ELISA Kit	PTHrP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0789	Rabbit PTHrP(Parathyroid Hormone Related Protein) ELISA Kit	PTHrP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0790	Rabbit TPP1(Tripeptidyl Peptidase I) ELISA Kit	TPP1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0791	Rabbit TBG(Thyroxine Binding Globulin) ELISA Kit	TBG	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0792	Rabbit PAH(Phenylalanine Hydroxylase) ELISA Kit	PAH	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0793	Rabbit GJβ1(Gap Junction Protein Beta 1) ELISA Kit	GJβ1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0794	Rabbit PCT(Procalcitonin) ELISA Kit	PCT	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB0795	Rabbit PCPE1(Procollagen C Proteinase Enhancer 1) ELISA Kit	PCPE1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0796	Rabbit CDSN(Corneodesmosin) ELISA Kit	CDSN	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0797	Rabbit CNTN1(Contactin 1) ELISA Kit	CNTN1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0798	Rabbit CNTN2(Contactin 2) ELISA Kit	CNTN2	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0799	Rabbit CNTN3(Contactin 3) ELISA Kit	CNTN3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0800	Rabbit CNTNAP1(Contactin Associated Protein 1) ELISA Kit	CNTNAP1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0801	Rabbit CNTNAP4(Contactin Associated Protein 4) ELISA Kit	CNTNAP4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0802	Rabbit ELMO2(Engulfment And Cell Motility 2) ELISA Kit	ELMO2	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0803	Rabbit Des(Desmin) ELISA Kit	Des	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL

FY-RB0804	Rabbit SSRP1(Structure Specific Recognition Protein 1) ELISA Kit	SSRP1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0805	Rabbit HPR(Haptoglobin Related Protein) ELISA Kit	HPR	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0806	Rabbit TNNT1(Troponin T Type 1, Slow Skeletal) ELISA Kit	TNNT1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0807	Rabbit TNPO2(Transportin 2) ELISA Kit	TNPO2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0808	Rabbit UCP-2(Uncoupling Protein 2, Mitochondrial) ELISA Kit	UCP-2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0809	Rabbit HLTF(Helicase Like Transcription Factor) ELISA Kit	HLTF	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0810	Rabbit MT2(Metallothionein 2) ELISA Kit	MT2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0811	Rabbit TJP1(Tight Junction Protein ZO-1) ELISA Kit	TJP1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0812	Rabbit CLC(Charcot Leyden Crystal Protein) ELISA Kit	CLC	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0813	Rabbit CRY α B(Crystallin Alpha B) ELISA Kit	CRY α B	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0814	Rabbit Sp17(Sperm Protein 17) ELISA Kit	Sp17	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0815	Rabbit MSP(Macrophage Stimulating Protein) ELISA Kit	MSP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0816	Rabbit MCSFR(Macrophage Colony Stimulating Factor Receptor) ELISA Kit	MCSFR	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0817	Rabbit MIP-5(Macrophage Inflammatory Protein 5) ELISA Kit	MIP-5	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0818	Rabbit FZD8(Frizzled Homolog 8) ELISA Kit	FZD8	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0819	Rabbit CRHBP(Corticotropin Releasing Hormone Binding Protein) ELISA Kit	CRHBP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0820	Rabbit α HSP(Alpha-Hemoglobin Stabilizing Protein) ELISA Kit	α HSP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0821	Rabbit SPTAN1(Alpha-Fodrin) ELISA Kit	SPTAN1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0822	Rabbit MSLN(Mesothelin) ELISA Kit	MSLN	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0823	Rabbit PGF(Placental Growth Factor) ELISA Kit	PGF	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB0824	Rabbit OTR(Oxytocin Receptor) ELISA Kit	OTR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0825	Rabbit PDIA6(Protein Disulfide Isomerase A6) ELISA Kit	PDIA6	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0826	Rabbit Pro-MMP-10(Pro-Matrix Metalloproteinase 10) ELISA Kit	Pro-MMP-10	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0827	Rabbit Pro-MMP-13(Pro-Matrix Metalloproteinase 13) ELISA Kit	Pro-MMP-13	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0828	Rabbit PDGF-AB(Platelet Derived Growth Factor AB) ELISA Kit	PDGF-AB	Sandwich	100ng/ml	0.94 ng/ml	1.56-100ng/ml
FY-RB0829	Rabbit Surv(Survivin) ELISA Kit	Surv	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0830	Rabbit TSP-1(Thrombospondin-1) ELISA Kit	TSP-1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0831	Rabbit TSP-2(Thrombospondin-2) ELISA Kit	TSP-2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0832	Rabbit TRAIL/TNFSF10(Tumor Necrosis Factor Related Apoptosis Inducing Ligand) ELISA Kit	TRAIL/TNFSF10	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0833	Rabbit TREM-1(Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit	TREM-1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB0834	Rabbit TSLP(Thymic Stromal Lymphopoietin) ELISA Kit	TSLP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0835	Rabbit PLAUR/uPAR(Plasminogen Activator, Urokinase Receptor) ELISA Kit	PLAUR/uPAR	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0836	Rabbit VEGF-C(Vascular Endothelial Growth Factor C) ELISA Kit	VEGF-C	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0837	Rabbit VEGF-D(Vascular Endothelial Growth Factor D) ELISA Kit	VEGF-D	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0838	Rabbit VEGFR-2/KDR(Vascular Endothelial Growth Factor Receptor 2) ELISA Kit	VEGFR-2/KDR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0839	Rabbit GzmA(Granzyme A) ELISA Kit	GzmA	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0840	Rabbit GzmB(Granzyme B) ELISA Kit	GzmB	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0841	Rabbit GNLY(Granulysin) ELISA Kit	GNLY	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-RB0842	Rabbit CR2/CD21(Complement Receptor 2) ELISA Kit	CR2/CD21	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0843	Rabbit cADPRH/CD38(Cyclic ADP Ribose Hydrolase) ELISA Kit	cADPRH/CD38	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB0844	Rabbit LOX-1(Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit	LOX-1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL

FY-RB0845	Rabbit sAC(Soluble Adenylate Cyclase) ELISA Kit	sAC	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0846	Rabbit sFMC(Soluble Fibrin Monomer Complex) ELISA Kit	sFMC	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0847	Rabbit Tyk-2(Tyrosine Kinase 2) ELISA Kit	Tyk-2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0848	Rabbit TH(Tyrosine Hydroxylase) ELISA Kit	TH	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0849	Rabbit PSPN(Persephin) ELISA Kit	PSPN	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB0850	Rabbit TGFβI(Transforming Growth Factor Beta Induced Protein) ELISA Kit	TGFβI	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0851	Rabbit TSC22(Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA Kit	TSC22	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0852	Rabbit PKCη(Protein Kinase C Eta) ELISA Kit	PKCη	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0853	Rabbit TYR(Tyrosinase) ELISA Kit	TYR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0854	Rabbit FGL1(Fibrinogen Like Protein 1) ELISA Kit	FGL1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0855	Rabbit PRG4(Proteoglycan 4) ELISA Kit	PRG4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0856	Rabbit CX40(Connexin 40) ELISA Kit	CX40	Sandwich	25ng/mL	0.23 ng/mL	0.39-25ng/mL
FY-RB0857	Rabbit CX43(Connexin 43) ELISA Kit	CX43	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0858	Rabbit SDC4(Syndecan 4) ELISA Kit	SDC4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0859	Rabbit LAP(Leucine Aminopeptidase) ELISA Kit	LAP	Sandwich	12ng/mL	0.11 ng/mL	0.19-12ng/mL
FY-RB0860	Rabbit LTβ(Lymphotoxin Beta) ELISA Kit	LTβ	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0861	Rabbit LFA-1(Lymphocyte Function Associated Antigen 1) ELISA Kit	LFA-1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0862	Rabbit LFA-2/CD2(Lymphocyte Function Associated Antigen 2) ELISA Kit	LFA-2/CD2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0863	Rabbit Lptn/XCL1(Lymphotoactin) ELISA Kit	Lptn/XCL1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0864	Rabbit PGAM2(Phosphoglycerate Mutase 2) ELISA Kit	PGAM2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0865	Rabbit PFKM(Phosphofruktokinase, Muscle) ELISA Kit	PFKM	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0866	Rabbit pIκBα(Phospho-Inhibitory Subunit Of NF Kappa B Alpha) ELISA Kit	pIκBα	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0867	Rabbit PGM(Phosphoglucomutase) ELISA Kit	PGM	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0868	Rabbit PCK1(Phosphoenolpyruvate Carboxykinase 1) ELISA Kit	PCK1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0869	Rabbit PLSCR1(Phospholipid Scramblase 1) ELISA Kit	PLSCR1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0870	Rabbit PLSCR2(Phospholipid Scramblase 2) ELISA Kit	PLSCR2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0871	Rabbit PI3K2α(Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit	PI3K2α	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0872	Rabbit GPC1(Glypican 1) ELISA Kit	GPC1	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0873	Rabbit GPC2(Glypican 2) ELISA Kit	GPC2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0874	Rabbit GPC3(Glypican 3) ELISA Kit	GPC3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0875	Rabbit GPC4(Glypican 4) ELISA Kit	GPC4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0876	Rabbit PLTP(Phospholipid Transfer Protein) ELISA Kit	PLTP	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0877	Rabbit PLCγ1(Phospholipase C Gamma 1) ELISA Kit	PLCγ1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0878	Rabbit DS(Dermatan Sulfate) ELISA Kit	DS	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0879	Rabbit Trx(Thioredoxin) ELISA Kit	Trx	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0880	Rabbit TrxR1(Thioredoxin Reductase 1) ELISA Kit	TrxR1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0881	Rabbit TBP2(Thioredoxin Binding Protein 2) ELISA Kit	TBP2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0883	Rabbit FcαR/CD89(Receptor I for the Fc Fragment of IgA) ELISA Kit	FcαR/CD89	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0884	Rabbit FcγR II/CD32(Receptor II for the Fc Fragment of IgG) ELISA Kit	FcγR II/CD32	Sandwich	120pg/mL	1.13 pg/mL	1.88-120pg/mL
FY-RB0885	Rabbit FcγR3B(Fc Fragment of IgG Low Affinity IIIb Receptor) ELISA Kit	FcγR3B	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0886	Rabbit Fcγ(Fc Fragment of IgG) ELISA Kit	Fcγ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL

FY-RB0887	Rabbit MCP/CD46(Membrane Cofactor Protein) ELISA Kit	MCP/CD46	Sandwich	20pg/mL	0.19 pg/mL	0.31-20pg/mL
FY-RB0888	Rabbit TSP-4(Thrombospondin-4) ELISA Kit	TSP-4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0889	Rabbit TDT(Terminal Deoxynucleotidyl Transferase) ELISA Kit	TDT	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB0890	Rabbit SMN1(Survival of Motor Neuron 1, Telomeric) ELISA Kit	SMN1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0891	Rabbit ETRA(Endothelin Receptor A) ELISA Kit	ETRA	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0892	Rabbit ECE1(Endothelin Converting Enzyme 1) ELISA Kit	ECE1	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0893	Rabbit EL/LIPG(Endothelial Lipase) ELISA Kit	EL/LIPG	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0894	Rabbit VASPIN(Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit	VASPIN	Sandwich	60ng/mL	0.56 ng/mL	0.94-60ng/mL
FY-RB0895	Rabbit VF(Visfatin) ELISA Kit	VF	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0896	Rabbit ERN1(Endoplasmic Reticulum to Nucleus Signalling 1) ELISA Kit	ERN1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0897	Rabbit CFTR(Cystic Fibrosis Transmembrane Conductance Regulator) ELISA Kit	CFTR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0898	Rabbit NGB(Neuroglobin) ELISA Kit	NGB	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0899	Rabbit MADCAM1(Mucosal Addressin Cell Adhesion Molecule 1) ELISA Kit	MADCAM1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0900	Rabbit FAK(Focal Adhesion Kinase) ELISA Kit	FAK	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0901	Rabbit OTC(Ornithine Carbamoyl Transferase) ELISA Kit	OTC	Sandwich	1200pg/mL	11.25 pg/mL	18.75-1200pg/mL
FY-RB0902	Rabbit ODC(Ornithine Decarboxylase) ELISA Kit	ODC	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0903	Rabbit GBP4(Guanylate Binding Protein 4) ELISA Kit	GBP4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0904	Rabbit UGCG(UDP-glucose ceramide glucosyltransferase) ELISA Kit	UGCG	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0905	Rabbit PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit	PLAU/uPA	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0906	Rabbit UPP1(Uridine Phosphorylase 1) ELISA Kit	UPP1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0907	Rabbit ETRB(Endothelin Receptor B) ELISA Kit	ETRB	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0908	Rabbit GS(Glutamine synthetase) ELISA Kit	GS	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0909	Rabbit F1+2(Prothrombin Fragment 1+2) ELISA Kit	F1+2	Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-RB0910	Rabbit KERA(Keratocan) ELISA Kit	KERA	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0911	Rabbit VCL(Vinculin) ELISA Kit	VCL	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB0912	Rabbit EED(Embryonic Ectoderm Development) ELISA Kit	EED	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0913	Rabbit DPT(Dermatopontin) ELISA Kit	DPT	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0914	Rabbit SYK(Spleen Tyrosine Kinase) ELISA Kit	SYK	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0915	Rabbit PDPN(Podoplanin) ELISA Kit	PDPN	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0916	Rabbit OSF(Osteoclast Stimulating Factor 1) ELISA Kit	OSF	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0917	Rabbit OSCAR(Osteoclast Associated Receptor) ELISA Kit	OSCAR	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0918	Rabbit G6PD(Glucose 6 Phosphate Dehydrogenase) ELISA Kit	G6PD	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0919	Rabbit GPI(Glucose 6 Phosphate Isomerase) ELISA Kit	GPI	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0920	Rabbit GCK(Glucokinase) ELISA Kit	GCK	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB0921	Rabbit GKRP(Glucokinase Regulatory Protein) ELISA Kit	GKRP	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0922	Rabbit GLUT1(Glucose Transporter 1) ELISA Kit	GLUT1	Sandwich	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-RB0923	Rabbit GLUT2(Glucose Transporter 2) ELISA Kit	GLUT2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0924	Rabbit GLUT3(Glucose Transporter 3) ELISA Kit	GLUT3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0925	Rabbit GLUT4(Glucose Transporter 4) ELISA Kit	GLUT4	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB0926	Rabbit GTF2B(General Transcription Factor II B) ELISA Kit	GTF2B	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0927	Rabbit GTF2F1(General Transcription Factor II F, Polypeptide 1) ELISA Kit	GTF2F1	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL

FY-RB0928	Rabbit PK2(Prokineticin 2) ELISA Kit	PK2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0929	Rabbit PKR1(Prokineticin Receptor 1) ELISA Kit	PKR1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0930	Rabbit PKR2(Prokineticin Receptor 2) ELISA Kit	PKR2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0931	Rabbit EP3(Prostaglandin E Receptor 3) ELISA Kit	EP3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0933	Rabbit FP(Prostaglandin F Receptor) ELISA Kit	FP	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0934	Rabbit PTGIS(Prostaglandin I Synthase) ELISA Kit	PTGIS	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0935	Rabbit PTGES2(Prostaglandin E Synthase 2) ELISA Kit	PTGES2	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0936	Rabbit PTGS1(Prostaglandin Endoperoxide Synthase 1) ELISA Kit	PTGS1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0937	Rabbit PTGS2/COX-2(Prostaglandin Endoperoxide Synthase 2) ELISA Kit	PTGS2/COX-2	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB0939	Rabbit DSP(Desmoplakin) ELISA Kit	DSP	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0940	Rabbit ERCC1(Excision Repair Cross Complementing Rodent Repair Deficiency Co	ERCC1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0941	Rabbit HSPβ6/HSP-20(Heat Shock Protein Beta 6) ELISA Kit	HSPβ6/HSP-20	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0942	Rabbit HSP-27/HSPB1(Heat Shock Protein 27) ELISA Kit	HSP-27/HSPB1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0943	Rabbit HSP-40(Heat Shock Protein 40) ELISA Kit	HSP-40	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0944	Rabbit HSP-60/HSPD1(Heat Shock Protein 60) ELISA Kit	HSP-60/HSPD1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0945	Rabbit HSP-90(Heat Shock Protein 90) ELISA Kit	HSP-90	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB0946	Rabbit HSF2(Heat Shock Transcription Factor 2) ELISA Kit	HSF2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0947	Rabbit PSG1/SP1(Pregnancy Specific Beta-1-Glycoprotein 1) ELISA Kit	PSG1/SP1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0948	Rabbit LZM(Lysozyme) ELISA Kit	LZM	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0949	Rabbit TFF1(Trefoil Factor 1) ELISA Kit	TFF1	Sandwich	300pg/mL	2.81 pg/mL	4.69-300pg/mL
FY-RB0950	Rabbit TFF2(Trefoil Factor 2) ELISA Kit	TFF2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0951	Rabbit PEDF(Pigment Epithelium Derived Factor) ELISA Kit	PEDF	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0952	Rabbit SDH(Sorbitol Dehydrogenase) ELISA Kit	SDH	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0953	Rabbit Ep-CAM/CD326(Epithelial Cell Adhesion Molecule) ELISA Kit	Ep-CAM/CD326	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0956	Rabbit GMFB(Glia Maturation Factor, Beta) ELISA Kit	GMFB	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0957	Rabbit P2(Myelin Protein 2) ELISA Kit	P2	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0959	Rabbit NCAM/CD56(Neural Cell Adhesion Molecule) ELISA Kit	NCAM/CD56	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB0960	Rabbit L1CAM/CD171(L1-Cell Adhesion Molecule) ELISA Kit	L1CAM/CD171	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-RB0961	Rabbit NTRK2(Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit	NTRK2	Sandwich	60ng/mL	0.56 ng/mL	0.94-60ng/mL
FY-RB0962	Rabbit NTRK3(Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit	NTRK3	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0963	Rabbit NAIP(Neuronal Apoptosis Inhibitory Protein) ELISA Kit	NAIP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0964	Rabbit NPHN(Nephrin) ELISA Kit	NPHN	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0965	Rabbit GTEAGE(Glomerular Tissue Extract Advanced Glycosylation End Products) I	GTEAGE	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0966	Rabbit GROγ/CXCL3(Growth Regulated Oncogene Gamma) ELISA Kit	GROγ/CXCL3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0967	Rabbit GDF1(Growth Differentiation Factor 1) ELISA Kit	GDF1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0968	Rabbit GDF10(Growth Differentiation Factor 10) ELISA Kit	GDF10	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0969	Rabbit GDF11(Growth Differentiation Factor 11) ELISA Kit	GDF11	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0970	Rabbit GDF3(Growth Differentiation Factor 3) ELISA Kit	GDF3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0971	Rabbit GDF6(Growth Differentiation Factor 6) ELISA Kit	GDF6	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0972	Rabbit GDF7(Growth Differentiation Factor 7) ELISA Kit	GDF7	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0973	Rabbit GDF9(Growth Differentiation Factor 9) ELISA Kit	GDF9	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL

FY-RB0974	Rabbit GHITM(Growth Hormone Inducible Transmembrane Protein) ELISA Kit	GHITM	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0975	Rabbit GAP43(Growth Associated Protein 43) ELISA Kit	GAP43	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB0977	Rabbit SSTR1(Somatostatin Receptor 1) ELISA Kit	SSTR1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0978	Rabbit Grb10(Growth Factor Receptor Bound Protein 10) ELISA Kit	Grb10	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0979	Rabbit Grb14(Growth Factor Receptor Bound Protein 14) ELISA Kit	Grb14	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0980	Rabbit Grb2(Growth Factor Receptor Bound Protein 2) ELISA Kit	Grb2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0981	Rabbit Grb7(Growth Factor Receptor Bound Protein 7) ELISA Kit	Grb7	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0983	Rabbit RB1(Retinoblastoma Protein 1) ELISA Kit	RB1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0984	Rabbit CYPB(Cyclophilin B) ELISA Kit	CYPB	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB0985	Rabbit CYPD(Cyclophilin D) ELISA Kit	CYPD	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0986	Rabbit ECFR(Eosinophil Chemotactic Factor Receptor) ELISA Kit	ECFR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0987	Rabbit EAR1(Eosinophil Associated Ribonuclease A family Member 1) ELISA Kit	EAR1	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB0988	Rabbit EAR2(Eosinophil Associated Ribonuclease A family Member 2) ELISA Kit	EAR2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0989	Rabbit EAR3(Eosinophil Associated Ribonuclease A family Member 3) ELISA Kit	EAR3	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0990	Rabbit PLN(Phospholamban) ELISA Kit	PLN	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB0991	Rabbit RIPK2(Receptor TNFRSF Interacting Serine Threonine Kinase 2) ELISA Kit	RIPK2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0992	Rabbit TNK2(Tyrosine Kinase, Non Receptor 2) ELISA Kit	TNK2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB0993	Rabbit DAF(Decay Accelerating Factor) ELISA Kit	DAF	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0994	Rabbit MAP2K6/MKK6(Mitogen Activated Protein Kinase Kinase 6) ELISA Kit	MAP2K6/MKK6	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0995	Rabbit RXFP1(Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit	RXFP1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0996	Rabbit TACR1(Tachykinin Receptor 1) ELISA Kit	TACR1	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB0997	Rabbit TACR2(Tachykinin Receptor 2) ELISA Kit	TACR2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0998	Rabbit TACR3(Tachykinin Receptor 3) ELISA Kit	TACR3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0999	Rabbit PLP1(Proteolipid Protein 1, Myelin) ELISA Kit	PLP1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1000	Rabbit MYD88(Myeloid Differentiation Factor 88) ELISA Kit	MYD88	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1001	Rabbit HBF(Fetal Hemoglobin) ELISA Kit	HBF	Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-RB1002	Rabbit PR3(Proteinase 3) ELISA Kit	PR3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1003	Rabbit PL(Placental Lactogen) ELISA Kit	PL	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB1004	Rabbit PP13(Placental Protein13) ELISA Kit	PP13	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1005	Rabbit PLAP/ALPP(Placental Alkaline Phosphatase) ELISA Kit	PLAP/ALPP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1006	Rabbit PGRP-S(Peptidoglycan Recognition Protein-S) ELISA Kit	PGRP-S	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1007	Rabbit PPIL1(Peptidylprolyl Isomerase/Cyclophilin Like 1) ELISA Kit	PPIL1	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1008	Rabbit GP5(Glycoprotein V, Platelet) ELISA Kit	GP5	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1009	Rabbit GP6(Glycoprotein VI, Platelet) ELISA Kit	GP6	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1010	Rabbit GlyCAM1(Glycosylation Dependent Cell Adhesion Molecule 1) ELISA Kit	GlyCAM1	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1011	Rabbit GR α (Glucocorticoid Receptor Alpha) ELISA Kit	GR α	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1012	Rabbit GR β (Glucocorticoid Receptor Beta) ELISA Kit	GR β	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1013	Rabbit GSK3 β (Glycogen Synthase Kinase 3 Beta) ELISA Kit	GSK3 β	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1014	Rabbit GSK3 α (Glycogen Synthase Kinase 3 Alpha) ELISA Kit	GSK3 α	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1015	Rabbit GLTP(Glycolipid Transfer Protein) ELISA Kit	GLTP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1016	Rabbit FTL(Ferritin, Light Polypeptide) ELISA Kit	FTL	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL

FY-RB1017	Rabbit FTH(Ferritin, Heavy Polypeptide) ELISA Kit	FTH	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-RB1018	Rabbit PIR(Pirin) ELISA Kit	PIR	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1019	Rabbit MRP1(Multidrug Resistance Associated Protein 1) ELISA Kit	MRP1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1020	Rabbit SYNPO(Synaptopodin) ELISA Kit	SYNPO	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1021	Rabbit SYP(Synaptophysin) ELISA Kit	SYP	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB1022	Rabbit SAP25(Synaptosome Associated Protein 25) ELISA Kit	SAP25	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1026	Rabbit DNase-I(Deoxyribonuclease I) ELISA Kit	DNase-I	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1028	Rabbit SIGLEC12(Sialic Acid Binding Ig Like Lectin 12) ELISA Kit	SIGLEC12	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1029	Rabbit SIGLEC3(Sialic Acid Binding Ig Like Lectin 3) ELISA Kit	SIGLEC3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1030	Rabbit ITLN1(Intelectin 1/Omentin) ELISA Kit	ITLN1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1031	Rabbit TUBα3C(Tubulin Alpha 3c) ELISA Kit	TUBα3C	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1032	Rabbit STMN1(Stathmin 1) ELISA Kit	STMN1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1033	Rabbit DCX(Doublecortin) ELISA Kit	DCX	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1034	Rabbit VDR(Vitamin D Receptor) ELISA Kit	VDR	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1035	Rabbit PP(Pepsin) ELISA Kit	PP	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1036	Rabbit PGA(Pepsinogen A) ELISA Kit	PGA	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1037	Rabbit PGC(Pepsinogen C) ELISA Kit	PGC	Sandwich	3000ng/mL	28.13 ng/mL	46.88-3000ng/mL
FY-RB1038	Rabbit GKN1(Gastrokine 1) ELISA Kit	GKN1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1039	Rabbit GKN2(Gastrokine 2) ELISA Kit	GKN2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1041	Rabbit GT-34(Gastrin 34) ELISA Kit	GT-34	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1042	Rabbit ProGRP(Pro-Gastrin Releasing Peptide) ELISA Kit	ProGRP	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1043	Rabbit GIP(Gastric Inhibitory Polypeptide) ELISA Kit	GIP	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB1044	Rabbit MFG-E8(Milk Fat Globule EGF Factor 8) ELISA Kit	MFG-E8	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1045	Rabbit CK-13/KRT13(Cytokeratin 13) ELISA Kit	CK-13/KRT13	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1046	Rabbit CK-19/KRT19(Cytokeratin 19) ELISA Kit	CK-19/KRT19	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1047	Rabbit CK-20/KRT20(Cytokeratin 20) ELISA Kit	CK-20/KRT20	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1048	Rabbit CYFRA21-1(Cytokeratin Fragment Antigen 21-1) ELISA Kit	CYFRA21-1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1049	Rabbit CVL/EZR(Cytovillin/Ezrin) ELISA Kit	CVL/EZR	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1050	Rabbit CYB561(Cytochrome b-561) ELISA Kit	CYB561	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1051	Rabbit CYP21B(Cytochrome P450 21B) ELISA Kit	CYP21B	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1052	Rabbit ECM1(Extracellular Matrix Protein 1) ELISA Kit	ECM1	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1053	Rabbit MEPE(Matrix Extracellular Phosphoglycoprotein) ELISA Kit	MEPE	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1054	Rabbit ERK1(Extracellular Signal Regulated Kinase 1) ELISA Kit	ERK1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1055	Rabbit ERK2(Extracellular Signal Regulated Kinase 2) ELISA Kit	ERK2	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB1056	Rabbit ERK3(Extracellular Signal Regulated Kinase 3) ELISA Kit	ERK3	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1057	Rabbit SOCS3(Suppressors of Cytokine Signaling 3) ELISA Kit	SOCS3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1058	Rabbit CYTH1(Cytohesin 1) ELISA Kit	CYTH1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1059	Rabbit CYTIP(Cytohesin Interacting Protein) ELISA Kit	CYTIP	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1060	Rabbit FLNβ(Filamin B, Beta) ELISA Kit	FLNβ	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1061	Rabbit FLNC(Filamin C, Gamma) ELISA Kit	FLNC	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1062	Rabbit FMOD(Fibromodulin) ELISA Kit	FMOD	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL

FY-RB1063	Rabbit FLRT2(Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit	FLRT2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1064	Rabbit PAP(Plasmin-Antiplasmin Complex) ELISA Kit	PAP	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB1065	Rabbit PLG(Plasminogen) ELISA Kit	PLG	Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-RB1066	Rabbit PAI1(Plasminogen Activator Inhibitor 1) ELISA Kit	PAI1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1067	Rabbit PAI2(Plasminogen Activator Inhibitor 2) ELISA Kit	PAI2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1068	Rabbit tPA(Tissue-type Plasminogen Activator) ELISA Kit	tPA	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-RB1069	Rabbit TAFI(Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit	TAFI	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1070	Rabbit FBL(Fibrillarlin) ELISA Kit	FBL	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1071	Rabbit GHR(Growth Hormone Receptor) ELISA Kit	GHR	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1072	Rabbit FGα(Fibrinogen Alpha) ELISA Kit	FGα	Sandwich	16ng/mL	0.15 ng/mL	0.25-16ng/mL
FY-RB1073	Rabbit FGβ(Fibrinogen Beta) ELISA Kit	FGβ	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1074	Rabbit FGL2(Fibrinogen Like Protein 2) ELISA Kit	FGL2	Sandwich	120pg/mL	1.13 pg/mL	1.88-120pg/mL
FY-RB1075	Rabbit SNRPC(Small Nuclear Ribonucleoprotein Polypeptide C) ELISA Kit	SNRPC	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1076	Rabbit SBEM(Small Breast Epithelial Mucin) ELISA Kit	SBEM	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1077	Rabbit SRP9(Signal Recognition Particle 9kDa) ELISA Kit	SRP9	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1078	Rabbit STAT1(Signal Transducer and Activator of Transcription 1) ELISA Kit	STAT1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1079	Rabbit STAT2(Signal Transducer and Activator of Transcription 2) ELISA Kit	STAT2	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1080	Rabbit STAT4(Signal Transducer and Activator of Transcription 4) ELISA Kit	STAT4	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1081	Rabbit STAT6(Signal Transducer and Activator of Transcription 6) ELISA Kit	STAT6	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1082	Rabbit SOX1(Sex Determining Region Y Box Protein 1) ELISA Kit	SOX1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1083	Rabbit SOX18(Sex Determining Region Y Box Protein 18) ELISA Kit	SOX18	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1084	Rabbit SOX3(Sex Determining Region Y Box Protein 3) ELISA Kit	SOX3	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1085	Rabbit TS/TYMS(Thymidylate Synthetase) ELISA Kit	TS/TYMS	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1086	Rabbit TLA(Thymus Leukemia Antigen) ELISA Kit	TLA	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1087	Rabbit TP(Thymidine Phosphorylase) ELISA Kit	TP	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1091	Rabbit GSTμ2/GSTm2(Glutathione S Transferase Mu 2, Muscle) ELISA Kit	GSTμ2/GSTm2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1092	Rabbit SUMO2(Small Ubiquitin Related Modifier 2) ELISA Kit	SUMO2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1093	Rabbit VEGF-B(Vascular Endothelial Cell Growth Factor B) ELISA Kit	VEGF-B	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1094	Rabbit VEGFR3/Flt4(Vascular Endothelial Cell Growth Factor Receptor 3) ELISA Kit	VEGFR3/Flt4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1095	Rabbit VWF(Von Willebrand Factor) ELISA Kit	VWF	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB1096	Rabbit HbC(Hemoglobin C) ELISA Kit	HbC	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1097	Rabbit HBβ(Hemoglobin Beta) ELISA Kit	HBβ	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1098	Rabbit HO1(Heme Oxygenase 1) ELISA Kit	HO1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1099	Rabbit HO2(Heme Oxygenase 2) ELISA Kit	HO2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1100	Rabbit SGK2(Serum/Glucocorticoid Regulated Kinase 2) ELISA Kit	SGK2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1101	Rabbit SEMA4C(Semaphorin 4C) ELISA Kit	SEMA4C	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1102	Rabbit SERT(Serotonin Transporter) ELISA Kit	SERT	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1103	Rabbit TPP(Thrombus Precursor Protein) ELISA Kit	TPP	Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-RB1105	Rabbit FDP(Fibrinogen Degradation Product) ELISA Kit	FDP	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1106	Rabbit FPA(Fibrinopeptide A) ELISA Kit	FPA	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1107	Rabbit FPB(Fibrinopeptide B) ELISA Kit	FPB	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL

FY-RB1108	Rabbit PHLDA2(Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit	PHLDA2	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB1109	Rabbit PAFAHIB3(Platelet Activating Factor Acetylhydrolase Ib3) ELISA Kit	PAFAHIB3	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1110	Rabbit PKP1(Plakophilin 1) ELISA Kit	PKP1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1111	Rabbit PKP2(Plakophilin 2) ELISA Kit	PKP2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1112	Rabbit PFKP(Phosphofructokinase, Platelet) ELISA Kit	PFKP	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1113	Rabbit PDGFRα(Platelet Derived Growth Factor Receptor Alpha) ELISA Kit	PDGFRα	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1114	Rabbit PF4V1(Platelet Factor 4 Variant 1) ELISA Kit	PF4V1	Sandwich	4ng/mL	0.04 ng/mL	0.06-4ng/mL
FY-RB1115	Rabbit GYPE(Glycophorin E) ELISA Kit	GYPE	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1116	Rabbit NOX1(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit	NOX1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1117	Rabbit EloA(Elongin A) ELISA Kit	EloA	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1118	Rabbit PI(Proinsulin) ELISA Kit	PI	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1119	Rabbit AMY2(Amylase Alpha 2, Pancreatic) ELISA Kit	AMY2	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1120	Rabbit GC(Glycocalicin) ELISA Kit	GC	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1121	Rabbit GLP-2(Glucagon Like Peptide 2) ELISA Kit	GLP-2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1122	Rabbit PTF1α(Pancreas Specific Transcription Factor 1α) ELISA Kit	PTF1α	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1123	Rabbit REG1α(Regenerating Islet Derived Protein 1 Alpha) ELISA Kit	REG1α	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1124	Rabbit REG4(Regenerating Islet Derived Protein 4) ELISA Kit	REG4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1125	Rabbit OSM(Oncostatin M) ELISA Kit	OSM	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1126	Rabbit PHB(Prohibitin) ELISA Kit	PHB	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1127	Rabbit CRY1(Cryptochrome 1) ELISA Kit	CRY1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1128	Rabbit KNG1(Kininogen 1) ELISA Kit	KNG1	Sandwich	20pg/mL	0.19 pg/mL	0.31-20pg/mL
FY-RB1129	Rabbit PCDH1(Protocadherin 1) ELISA Kit	PCDH1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1130	Rabbit PCDHβ2(Protocadherin Beta 2) ELISA Kit	PCDHβ2	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1131	Rabbit FBN1(Fibrillin 1) ELISA Kit	FBN1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1132	Rabbit FBN2(Fibrillin 2) ELISA Kit	FBN2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1133	Rabbit FBN3(Fibrillin 3) ELISA Kit	FBN3	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1134	Rabbit PAQR7(Progesterin and AdipoQ Receptor Family Member VII) ELISA Kit	PAQR7	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1135	Rabbit PGR(Progesterone Receptor) ELISA Kit	PGR	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1136	Rabbit PS2(Presenilin 2) ELISA Kit	PS2	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1137	Rabbit EGR1(Early Growth Response Protein 1) ELISA Kit	EGR1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1138	Rabbit EGR2(Early Growth Response Protein 2) ELISA Kit	EGR2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1139	Rabbit EGR3(Early Growth Response Protein 3) ELISA Kit	EGR3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1140	Rabbit EGR4(Early Growth Response Protein 4) ELISA Kit	EGR4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1141	Rabbit NFYC(Nuclear Transcription Factor Y Gamma) ELISA Kit	NFYC	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1142	Rabbit GM(Gastric Mucin) ELISA Kit	GM	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1143	Rabbit MUC5B(Mucin 5 Subtype B) ELISA Kit	MUC5B	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1144	Rabbit MUC7(Mucin 7) ELISA Kit	MUC7	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1145	Rabbit eIF3a(Eukaryotic Translation Initiation Factor 3a) ELISA Kit	eIF3a	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1146	Rabbit PC4(Positive Cofactor 4) ELISA Kit	PC4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1147	Rabbit OSMR(Oncostatin M Receptor) ELISA Kit	OSMR	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB1148	Rabbit MK(Midkine) ELISA Kit	MK	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-RB1149	Rabbit PON1(Paraoxonase 1) ELISA Kit	PON1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1150	Rabbit TP63(Tumor Protein p63) ELISA Kit	TP63	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB1151	Rabbit TNF-β(Tumor Necrosis Factor Beta) ELISA Kit	TNF-β	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB1152	Rabbit TL1A(Tumor Necrosis Factor Related Ligand 1A) ELISA Kit	TL1A	Sandwich	200pg/ml	1.88 pg/ml	3.13-200pg/ml
FY-RB1153	Rabbit TSTA(Tumor Specific Transplantation Antigen) ELISA Kit	TSTA	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1154	Rabbit RAG-1(Recombination Activating Gene 1) ELISA Kit	RAG-1	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1155	Rabbit RAG-2(Recombination Activating Gene 2) ELISA Kit	RAG-2	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1156	Rabbit CDK1(Cyclin Dependent Kinase 1) ELISA Kit	CDK1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1157	Rabbit CDK2(Cyclin Dependent Kinase 2) ELISA Kit	CDK2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1158	Rabbit CDK4(Cyclin Dependent Kinase 4) ELISA Kit	CDK4	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1159	Rabbit CDK5(Cyclin Dependent Kinase 5) ELISA Kit	CDK5	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1160	Rabbit CDK7(Cyclin Dependent Kinase 7) ELISA Kit	CDK7	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1161	Rabbit CDK9(Cyclin Dependent Kinase 9) ELISA Kit	CDK9	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1162	Rabbit CDKN2A(Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit	CDKN2A	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1163	Rabbit Ntn1(Netrin 1) ELISA Kit	Ntn1	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1164	Rabbit Ntn4(Netrin 4) ELISA Kit	Ntn4	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1165	Rabbit TGM2(Transglutaminase 2, Tissue) ELISA Kit	TGM2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1166	Rabbit TGM3(Transglutaminase 3, Epidermal) ELISA Kit	TGM3	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1167	Rabbit TTR(Transthyretin) ELISA Kit	TTR	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB1168	Rabbit TCF20(Transcription Factor 20) ELISA Kit	TCF20	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1169	Rabbit TFR2(Transferrin Receptor 2) ELISA Kit	TFR2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1170	Rabbit TNPO1(Transportin 1) ELISA Kit	TNPO1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1171	Rabbit Pax(Paxillin) ELISA Kit	Pax	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1172	Rabbit PCX(Podocalyxin) ELISA Kit	PCX	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1173	Rabbit TPS(Tissue Polypeptide Specific Antigen) ELISA Kit	TPS	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1174	Rabbit TFPI2(Tissue Factor Pathway Inhibitor 2) ELISA Kit	TFPI2	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1175	Rabbit PON3(Paraoxonase 3) ELISA Kit	PON3	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-RB1176	Rabbit DDT(D-Dopachrome Tautomerase) ELISA Kit	DDT	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1177	Rabbit sC5b-9(Soluble Terminal Complement Complex C5b-9) ELISA Kit	sC5b-9	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB1178	Rabbit UCHL1(Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit	UCHL1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1179	Rabbit UBR4(Ubiquitin Protein Ligase E3 Component N-Recognin 4) ELISA Kit	UBR4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1180	Rabbit UBE2A(Ubiquitin Conjugating Enzyme E2A) ELISA Kit	UBE2A	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1182	Rabbit BCAR1(Breast Cancer Anti Estrogen Resistance 1) ELISA Kit	BCAR1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1183	Rabbit PLA2G2D(Phospholipase A2, Group II D) ELISA Kit	PLA2G2D	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1184	Rabbit UBAP2(Ubiquitin Associated Protein 2) ELISA Kit	UBAP2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1185	Rabbit MASP1(Mannan Associated Serine Protease 1) ELISA Kit	MASP1	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-RB1186	Rabbit MASP2(Mannan Associated Serine Protease 2) ELISA Kit	MASP2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1187	Rabbit USP8(Ubiquitin Specific Peptidase 8) ELISA Kit	USP8	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1188	Rabbit DPP6(Dipeptidyl Peptidase VI) ELISA Kit	DPP6	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB1189	Rabbit MRC2(Mannose Receptor C Type 2) ELISA Kit	MRC2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1190	Rabbit MRC1(Mannose Receptor C Type 1) ELISA Kit	MRC1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL

FY-RB1191	Rabbit GADL1(Glutamate Decarboxylase Like Protein 1) ELISA Kit	GADL1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1192	Rabbit GAD1(Glutamate Decarboxylase 1) ELISA Kit	GAD1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1193	Rabbit PPAR- δ (Peroxisome Proliferator Activated Receptor Delta) ELISA Kit	PPAR- δ	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1194	Rabbit NUP62(Nucleoporin 62kDa) ELISA Kit	NUP62	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1195	Rabbit PCNA(Proliferating Cell Nuclear Antigen) ELISA Kit	PCNA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1196	Rabbit RAN κ (Receptor Activator of Nuclear Factor Kappa B) ELISA Kit	RAN κ	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1197	Rabbit MCAM(Melanoma Cell Adhesion Molecule) ELISA Kit	MCAM	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1198	Rabbit MIA1(Melanoma Inhibitory Activity Protein 1) ELISA Kit	MIA1	Sandwich	20pg/mL	0.19 pg/mL	0.31-20pg/mL
FY-RB1199	Rabbit MFI2(Melanotransferrin) ELISA Kit	MFI2	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1200	Rabbit SPT β N4(Spectrin Beta, Non Erythrocytic 4) ELISA Kit	SPT β N4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1201	Rabbit I κ B α (Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit	I κ B α	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1202	Rabbit HNRPA2B1(Heterogeneous Nuclear Ribonucleoprotein A2/B1) ELISA Kit	HNRPA2B1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1203	Rabbit HNRPU(Heterogeneous Nuclear Ribonucleoprotein U) ELISA Kit	HNRPU	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1204	Rabbit RAR α (Retinoic Acid Receptor Alpha) ELISA Kit	RAR α	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1205	Rabbit RNASEN(Ribonuclease III, Nuclear) ELISA Kit	RNASEN	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1206	Rabbit MCSP(Melanoma Associated Chondroitin Sulfate Proteoglycan) ELISA Kit	MCSP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1208	Rabbit THR α (Thyroid Hormone Receptor Alpha) ELISA Kit	THR α	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1209	Rabbit NFR κ B(Nuclear Factor Related To Kappa B Binding Protein) ELISA Kit	NFR κ B	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1210	Rabbit TNKS1(Tankyrase 1) ELISA Kit	TNKS1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1211	Rabbit TNKS2(Tankyrase 2) ELISA Kit	TNKS2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1212	Rabbit PARP4(Poly ADP Ribose Polymerase 4) ELISA Kit	PARP4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1213	Rabbit ME20M(Melanoma Associated ME20) ELISA Kit	ME20M	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1214	Rabbit PXR(Pregnane X Receptor) ELISA Kit	PXR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1215	Rabbit RELB(V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit	RELB	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1216	Rabbit IL-24(Interleukin 24) ELISA Kit	IL-24	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1217	Rabbit AIM2(Absent In Melanoma 2) ELISA Kit	AIM2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1218	Rabbit TSPAN30/CD63(Tetraspanin 30) ELISA Kit	TSPAN30/CD63	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1219	Rabbit BANP(BTG3 Associated Nuclear Protein) ELISA Kit	BANP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1220	Rabbit PRKG1(Protein Kinase, cGMP Dependent Type I) ELISA Kit	PRKG1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1221	Rabbit PDE4D(Phosphodiesterase 4D, cAMP Specific) ELISA Kit	PDE4D	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1222	Rabbit CKB(Creatine Kinase, Brain) ELISA Kit	CKB	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1223	Rabbit LOXL4(Lysyl Oxidase Like Protein 4) ELISA Kit	LOXL4	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1224	Rabbit TNFRSF1B(Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELIS	TNFRSF1B	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1225	Rabbit T α 1(Thymosin- α 1) ELISA Kit	T α 1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1226	Rabbit PACAP-38(Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Ki	PACAP-38	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB1227	Rabbit PACAP-27(Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Ki	PACAP-27	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB1228	Rabbit PTK6(Protein Tyrosine Kinase 6) ELISA Kit	PTK6	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1229	Rabbit DAG1(Dystrophin Associated Glycoprotein 1) ELISA Kit	DAG1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1230	Rabbit GPX1(Glutathione Peroxidase 1) ELISA Kit	GPX1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1231	Rabbit ANXA1(Annexin A1) ELISA Kit	ANXA1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1232	Rabbit XOD(Xanthine Oxidase) ELISA Kit	XOD	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL

FY-RB1233	Rabbit VCAM-1/CD106(Vascular Cell Adhesion Molecule 1) ELISA Kit	VCAM-1/CD106	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1234	Rabbit CA125(Carbohydrate Antigen 125) ELISA Kit	CA125	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1235	Rabbit MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	MCP-1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1236	Rabbit CFD(Complement Factor D) ELISA Kit	CFD	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB1237	Rabbit IL-8(Interleukin 8) ELISA Kit	IL-8	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1238	Rabbit CXCL16(Chemokine C-X-C-Motif Ligand 16) ELISA Kit	CXCL16	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1239	Rabbit ES(Endostatin) ELISA Kit	ES	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1240	Rabbit Hcp(Hepcidin) ELISA Kit	Hcp	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1241	Rabbit Gas6(Growth Arrest Specific Protein 6) ELISA Kit	Gas6	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1242	Rabbit sgp130(soluble Glucoprotein 130) ELISA Kit	sgp130	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1243	Rabbit LEP(Leptin) ELISA Kit	LEP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1244	Rabbit MAU(Microalbuminuria) ELISA Kit	MAU	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1245	Rabbit LH(Luteinizing Hormone) ELISA Kit	LH	Sandwich	100mIU/mL	0.94 mIU/mL	1.56-100mIU/mL
FY-RB1246	Rabbit C-P(C-Peptide) ELISA Kit	C-P	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1247	Rabbit OxLDL(Oxidized Low Density Lipoprotein) ELISA Kit	OxLDL	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1248	Rabbit MMP-10(Matrix Metalloproteinase 10) ELISA Kit	MMP-10	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1249	Rabbit MMP-13(Matrix Metalloproteinase 13) ELISA Kit	MMP-13	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1250	Rabbit MYO(Myoglobin) ELISA Kit	MYO	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB1251	Rabbit GLP-1(Glucagon Like Peptide 1) ELISA Kit	GLP-1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1252	Rabbit C4(Complement Component 4) ELISA Kit	C4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1253	Rabbit TF(Transferrin) ELISA Kit	TF	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1254	Rabbit KIM-1(Kidney Injury Molecule 1) ELISA Kit	KIM-1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1255	Rabbit PICP(Procollagen I C-terminal Propeptide) ELISA Kit	PICP	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1256	Rabbit AChE(Acetylcholinesterase) ELISA Kit	AChE	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB1257	Rabbit ALCAM(Activated Leukocyte Cell Adhesion Molecule) ELISA Kit	ALCAM	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1258	Rabbit ALDOA(Aldolase A, Fructose Bisphosphate) ELISA Kit	ALDOA	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1259	Rabbit ANGPTL2(Angiopietin Like Protein 2) ELISA Kit	ANGPTL2	Sandwich	1200pg/mL	11.25 pg/mL	18.75-1200pg/mL
FY-RB1260	Rabbit FGFR4(Fibroblast Growth Factor Receptor 4) ELISA Kit	FGFR4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1261	Rabbit CYP1A1(Cytochrome P450, family 1, subfamily A, polypeptide 1) ELISA Kit	CYP1A1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1262	Rabbit IGF-2(Insulin Like Growth Factor 2) ELISA Kit	IGF-2	Sandwich	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-RB1263	Rabbit IGFBP-2(Insulin Like Growth Factor Binding Protein 2) ELISA Kit	IGFBP-2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1264	Rabbit ApoA2(Apolipoprotein A2) ELISA Kit	ApoA2	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-RB1265	Rabbit ApoF(Apolipoprotein F) ELISA Kit	ApoF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1266	Rabbit ApoA5(Apolipoprotein A5) ELISA Kit	ApoA5	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1267	Rabbit bFGF/FGF2(Basic Fibroblast Growth Factor) ELISA Kit	bFGF/FGF2	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB1268	Rabbit FGF4(Fibroblast Growth Factor 4) ELISA Kit	FGF4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1269	Rabbit BMP-1(Bone Morphogenetic Protein 1) ELISA Kit	BMP-1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1270	Rabbit BMP-15(Bone Morphogenetic Protein 15) ELISA Kit	BMP-15	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1271	Rabbit BMP-6(Bone Morphogenetic Protein 6) ELISA Kit	BMP-6	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1272	Rabbit CEA(Carcinoembryonic Antigen) ELISA Kit	CEA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1273	Rabbit CKLF(Chemokine Like Factor) ELISA Kit	CKLF	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL

FY-RB1274	Rabbit COL3(Collagen Type III) ELISA Kit	COL3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1275	Rabbit COL10(Collagen Type X) ELISA Kit	COL10	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1276	Rabbit COLEC11(Collectin Sub-Family Member 11) ELISA Kit	COLEC11	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1277	Rabbit C1INH(Complement 1 Inhibitor) ELISA Kit	C1INH	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1278	Rabbit C1q(Complement 1q) ELISA Kit	C1q	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1279	Rabbit C3(Complement Component 3) ELISA Kit	C3	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB1280	Rabbit C2(Complement Component 2) ELISA Kit	C2	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB1281	Rabbit CFB(Complement Factor B) ELISA Kit	CFB	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB1282	Rabbit CFH(Complement Factor H) ELISA Kit	CFH	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1283	Rabbit CFI(Complement Factor I) ELISA Kit	CFI	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1284	Rabbit C4b(Complement Component 4b) ELISA Kit	C4b	Sandwich	12ng/mL	0.11 ng/mL	0.19-12ng/mL
FY-RB1285	Rabbit α 2-M(Alpha-2 Macroglobulin) ELISA Kit	α 2-M	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB1286	Rabbit GL β (Galactosidase Beta) ELISA Kit	GL β	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1287	Rabbit M1 γ /CXCL9(Monocyte Interferon Gamma Inducing Factor) ELISA Kit	M1 γ /CXCL9	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1288	Rabbit PKM2(Pyruvate Kinase, Muscle) ELISA Kit	PKM2	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1289	Rabbit VTN(Vitronectin) ELISA Kit	VTN	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1290	Rabbit D β H(Dopamine- β -Hydroxylase) ELISA Kit	D β H	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1291	Rabbit HIF-1 α (Hypoxia Inducible Factor 1 Alpha) ELISA Kit	HIF-1 α	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1292	Rabbit LRG1(Leucine Rich Alpha-2-Glycoprotein 1) ELISA Kit	LRG1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1293	Rabbit GDH/GLDH(Glutamate dehydrogenase) ELISA Kit	GDH/GLDH	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1294	Rabbit GAD2(Glutamate Decarboxylase 2) ELISA Kit	GAD2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1295	Rabbit GST α 1(Glutathione S Transferase Alpha 1) ELISA Kit	GST α 1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1296	Rabbit IgA(Immunoglobulin A) ELISA Kit	IgA	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB1297	Rabbit TLN(Talin) ELISA Kit	TLN	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1298	Rabbit MMP-1(Matrix Metalloproteinase 1) ELISA Kit	MMP-1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1299	Rabbit MMP-9(Matrix Metalloproteinase 9) ELISA Kit	MMP-9	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1300	Rabbit I-PTH(intact Parathormone) ELISA Kit	I-PTH	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-RB1301	Rabbit HP(Haptoglobin) ELISA Kit	HP	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-RB1302	Rabbit MIP-3 β /ELC/CCL19(Macrophage Inflammatory Protein 3 β) ELISA Kit	MIP-3 β /ELC/CCL19	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1303	Rabbit EGFR2(Epidermal Growth Factor Receptor 2) ELISA Kit	EGFR2	Sandwich	300pg/mL	2.81 pg/mL	4.69-300pg/mL
FY-RB1304	Rabbit ApoB48(Apolipoprotein B48) ELISA Kit	ApoB48	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1305	Rabbit sTfR1(Soluble Transferrin Receptor 1) ELISA Kit	sTfR1	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB1306	Rabbit ORM2(Orosomucoid 2) ELISA Kit	ORM2	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1307	Rabbit NT-proANP(N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit	NT-proANP	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1308	Rabbit BGN(Biglycan) ELISA Kit	BGN	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1309	Rabbit NRG-1(Neuregulin 1) ELISA Kit	NRG-1	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB1310	Rabbit GFAP(Glial Fibrillary Acidic Protein) ELISA Kit	GFAP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1311	Rabbit NT-4(Neurotrophin 4) ELISA Kit	NT-4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1312	Rabbit CYPA(Cyclophilin A) ELISA Kit	CYPA	Sandwich	3000ng/mL	28.13 ng/mL	46.88-3000ng/mL
FY-RB1313	Rabbit FETUB(Fetuin B) ELISA Kit	FETUB	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-RB1314	Rabbit SHBG(Sex Hormone Binding Globulin) ELISA Kit	SHBG	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-RB1315	Rabbit TECK/CCL25(Thymus Expressed Chemokine) ELISA Kit	TECK/CCL25	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1316	Rabbit OLFM4(Olfactomedin 4) ELISA Kit	OLFM4	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1317	Rabbit VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit	VE-Cadherin	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1318	Rabbit IgE(Immunoglobulin E) ELISA Kit	IgE	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1319	Rabbit ALB(Albumin) ELISA Kit	ALB	Sandwich	200µg/mL	1.88 µg/mL	3.13-200µg/mL
FY-RB1320	Rabbit FG(Fibrinogen) ELISA Kit	FG	Sandwich	24ng/mL	0.23 ng/mL	0.38-24ng/mL
FY-RB1321	Rabbit LBP(Lipopolysaccharide Binding Protein) ELISA Kit	LBP	Sandwich	1200ng/mL	11.25 ng/mL	18.75-1200ng/mL
FY-RB1322	Rabbit SAA2(Serum Amyloid A2) ELISA Kit	SAA2	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1323	Rabbit sICAM-1/CD54(Soluble Intercellular Adhesion Molecule 1) ELISA Kit	sICAM-1/CD54	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1324	Rabbit F8(Coagulation Factor VIII) ELISA Kit	F8	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1325	Rabbit ANG1(Angiopoietin 1) ELISA Kit	ANG1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1326	Rabbit ADP/Acrp30(Adiponectin) ELISA Kit	ADP/Acrp30	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1327	Rabbit IFNα1(Interferon Alpha 1) ELISA Kit	IFNα1	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB1328	Rabbit NGAL(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	NGAL	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB1329	Rabbit DPP4(Dipeptidyl Peptidase IV) ELISA Kit	DPP4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1330	Rabbit sCD14(Soluble Cluster of Differentiation14) ELISA Kit	sCD14	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1331	Rabbit EMMPRIN/CD147(Extracellular Matrix Metalloproteinase Inducer) ELISA K	EMMPRIN/CD147	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1332	Rabbit CTSB(Cathepsin B) ELISA Kit	CTSB	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1333	Rabbit FABP1(Fatty Acid Binding Protein 1, Liver) ELISA Kit	FABP1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1334	Rabbit BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	BDNF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1335	Rabbit CASP1(Caspase 1) ELISA Kit	CASP1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1336	Rabbit CASP3(Caspase 3) ELISA Kit	CASP3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1337	Rabbit MIP-3α(Macrophage Inflammatory Protein 3 Alpha) ELISA Kit	MIP-3α	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1338	Rabbit MDC/CCL22(Macrophage-Derived Chemokine) ELISA Kit	MDC/CCL22	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1339	Rabbit CLU(Clusterin) ELISA Kit	CLU	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1340	Rabbit CRP(C-reactive protein) ELISA Kit	CRP	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1341	Rabbit ENA-78(Epithelial Neutrophil Activating Peptide 78) ELISA Kit	ENA-78	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1342	Rabbit SDF-1/CXCL12(Stromal Cell Derived Factor 1) ELISA Kit	SDF-1/CXCL12	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1343	Rabbit αFP(Alpha-Fetoprotein) ELISA Kit	αFP	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1344	Rabbit APN(Aminopeptidase N) ELISA Kit	APN	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1345	Rabbit GM-CSF(Granulocyte Macrophage Colony Stimulating Factor) ELISA Kit	GM-CSF	Sandwich	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-RB1346	Rabbit IGF-1(Insulin-like Growth Factor 1) ELISA Kit	IGF-1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1347	Rabbit FGF7/KGF(Fibroblast Growth Factor 7) ELISA Kit	FGF7/KGF	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1348	Rabbit IL-2(Interleukin 2) ELISA Kit	IL-2	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB1349	Rabbit IL-17(Interleukin 17) ELISA Kit	IL-17	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB1350	Rabbit ACAT2(Acetyl Coenzyme A Acetyltransferase 2) ELISA Kit	ACAT2	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1351	Rabbit TnTf/TNNT3(Troponin T, fast skeletal muscle) ELISA Kit	TnTf/TNNT3	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB1352	Rabbit IFABP/FABP2(Intestinal Fatty Acid Binding Protein) ELISA Kit	IFABP/FABP2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1353	Rabbit LP-a(Lipoprotein a) ELISA Kit	LP-a	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1354	Rabbit FE(Ferritin) ELISA Kit	FE	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB1356	Rabbit CG(Chorionic Gonadotropin) ELISA Kit	CG	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL

FY-RB1358	Rabbit ACD(Adrenocortical Dysplasia Homolog) ELISA Kit	ACD	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1359	Rabbit ACP1(Acid Phosphatase 1) ELISA Kit	ACP1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1360	Rabbit ADCY3(Adenylate Cyclase 3) ELISA Kit	ADCY3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1361	Rabbit ADCY4(Adenylate Cyclase 4) ELISA Kit	ADCY4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1362	Rabbit AIG1(Androgen Induced Protein 1) ELISA Kit	AIG1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1363	Rabbit AKAP1(A Kinase Anchor Protein 1) ELISA Kit	AKAP1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1364	Rabbit ALK(Anaplastic Lymphoma Kinase) ELISA Kit	ALK	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1367	Rabbit ANGPTL3(Angiopoietin Like Protein 3) ELISA Kit	ANGPTL3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1368	Rabbit RARS(Arginyl tRNA Synthetase) ELISA Kit	RARS	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1369	Rabbit iPLA2(Phospholipidase A2, Calcium Independent) ELISA Kit	iPLA2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1370	Rabbit GUK1(Guanylate Kinase 1) ELISA Kit	GUK1	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1371	Rabbit HB(Hemoglobin) ELISA Kit	HB	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1372	Rabbit HRG(Histidine Rich Glycoprotein) ELISA Kit	HRG	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1373	Rabbit AT-III(Antithrombin III) ELISA Kit	AT-III	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB1374	Rabbit IGFBP-1(Insulin Like Growth Factor Binding Protein 1) ELISA Kit	IGFBP-1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1375	Rabbit IGFBP-5(Insulin Like Growth Factor Binding Protein 5) ELISA Kit	IGFBP-5	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1376	Rabbit IGFBP-6(Insulin Like Growth Factor Binding Protein 6) ELISA Kit	IGFBP-6	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1379	Rabbit ApoC3(Apolipoprotein C3) ELISA Kit	ApoC3	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB1380	Rabbit ITIH4(Inter-alpha-trypsin Inhibitor Heavy Chain H4) ELISA Kit	ITIH4	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1381	Rabbit ARG1(Arginase I) ELISA Kit	ARG1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB1382	Rabbit ATXN2(Ataxin 2) ELISA Kit	ATXN2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1383	Rabbit ATF7(Activating Transcription Factor 7) ELISA Kit	ATF7	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1384	Rabbit ABCG1(ATP Binding Cassette Transporter G1) ELISA Kit	ABCG1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1385	Rabbit ABCA1(ATP-binding cassette transporter A1) ELISA Kit	ABCA1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1386	Rabbit ABCG2(ATP Binding Cassette Transporter G2) ELISA Kit	ABCG2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1387	Rabbit Aβ1-40(Amyloid Beta 1-40) ELISA Kit	Aβ1-40	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1388	Rabbit Aβ1-42(Amyloid Beta 1-42) ELISA Kit	Aβ1-42	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1389	Rabbit Bcl-3(B-cell Lymphoma/Leukemia 3) ELISA Kit	Bcl-3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1390	Rabbit BLNK(B-Cell Linker Protein) ELISA Kit	BLNK	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1391	Rabbit BTC(Betacellulin) ELISA Kit	BTC	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-RB1392	Rabbit BMPR1A(Bone Morphogenetic Protein Receptor 1A) ELISA Kit	BMPR1A	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB1393	Rabbit IBSP(Bone Sialoprotein 2) ELISA Kit	IBSP	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB1394	Rabbit BFP(Brain Finger Protein) ELISA Kit	BFP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1395	Rabbit C4BP(C4 Binding Protein) ELISA Kit	C4BP	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1396	Rabbit CGRP1(Calcitonin Gene Related Peptide 1) ELISA Kit	CGRP1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1397	Rabbit CaSR(Calcium Sensing Receptor) ELISA Kit	CaSR	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1398	Rabbit CAP1(Adenylate Cyclase Associated Protein 1) ELISA Kit	CAP1	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1399	Rabbit CMLC-1(Cardiac Myosin Light Chains 1) ELISA Kit	CMLC-1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1400	Rabbit cTn-I/TNNI3(cardiac Troponin I) ELISA Kit	cTn-I/TNNI3	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB1401	Rabbit CT1(Cardiotrophin 1) ELISA Kit	CT1	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1402	Rabbit CASP10(Caspase 10) ELISA Kit	CASP10	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL

FY-RB1403	Rabbit CASP12(Caspase 12) ELISA Kit	CASP12	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1404	Rabbit Cav-1(Caveolin-1) ELISA Kit	Cav-1	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-RB1405	Rabbit CCR1(CC-Chemokine Receptor 1) ELISA Kit	CCR1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1406	Rabbit CDC23(Cell Division Cycle Protein 23) ELISA Kit	CDC23	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1407	Rabbit CENPE(Centromere Protein E) ELISA Kit	CENPE	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1408	Rabbit CENPH(Centromere Protein H) ELISA Kit	CENPH	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1409	Rabbit CENPB(Centromere Protein B) ELISA Kit	CENPB	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1410	Rabbit CEP110(Centrosomal Protein 110kDa) ELISA Kit	CEP110	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1411	Rabbit CPN10(Chaperonin 10) ELISA Kit	CPN10	Sandwich	300pg/mL	2.81 pg/mL	4.69-300pg/mL
FY-RB1412	Rabbit CCL1(Chemokine C-C-Motif Ligand 1) ELISA Kit	CCL1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1413	Rabbit CCL4L1(Chemokine C-C-Motif Ligand 4 Like Protein 1) ELISA Kit	CCL4L1	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1414	Rabbit CCRL1(Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit	CCRL1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1415	Rabbit CX3CR1(Chemokine C-X3-C-Motif Receptor 1) ELISA Kit	CX3CR1	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1416	Rabbit CKLFSF7(Chemokine Like Factor Superfamily 7) ELISA Kit	CKLFSF7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1417	Rabbit CKLFSF8(Chemokine Like Factor Superfamily 8) ELISA Kit	CKLFSF8	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1418	Rabbit CCKBR(Cholecystokinin B Receptor) ELISA Kit	CCKBR	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1419	Rabbit CHAc(Choline acetylase) ELISA Kit	CHAc	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1420	Rabbit PC/CPG(Choline Phosphoglyceride) ELISA Kit	PC/CPG	Sandwich	16ng/mL	0.15 ng/mL	0.25-16ng/mL
FY-RB1421	Rabbit CHD3(Chromodomain Helicase DNA Binding Protein 3) ELISA Kit	CHD3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1422	Rabbit CHGA(Chromogranin A) ELISA Kit	CHGA	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1423	Rabbit CLDN14(Claudin 14) ELISA Kit	CLDN14	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1424	Rabbit CD109(Cluster of Differentiation 109) ELISA Kit	CD109	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1425	Rabbit GRN(Granulin) ELISA Kit	GRN	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1426	Rabbit CSF2R α (Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit	CSF2R α	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1427	Rabbit GCSFR(Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit	GCSFR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1428	Rabbit C1qR1(Complement Component 1q Receptor 1) ELISA Kit	C1qR1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1429	Rabbit C1r(Complement Component 1R) ELISA Kit	C1r	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1430	Rabbit C8b(Complement Component 8b) ELISA Kit	C8b	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1431	Rabbit C8g(Complement Component 8g) ELISA Kit	C8g	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1432	Rabbit CTAPIII(Connective Tissue Activating Peptide III) ELISA Kit	CTAPIII	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB1433	Rabbit CX31(Connexin 31) ELISA Kit	CX31	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB1434	Rabbit HELIX- II (Type II Collagen Helical Peptide) ELISA Kit	HELIX- II	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1435	Rabbit MTHFR(Methylenetetrahydrofolate Reductase) ELISA Kit	MTHFR	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1436	Rabbit PRKA α 1(Protein Kinase, AMP Activated Alpha 1) ELISA Kit	PRKA α 1	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1437	Rabbit I κ B β (Inhibitory Subunit of NF Kappa B Beta) ELISA Kit	I κ B β	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1438	Rabbit PTPRB(Protein Tyrosine Phosphatase Receptor Type B) ELISA Kit	PTPRB	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1439	Rabbit CREB1(Cyclic AMP Response Element Binding Protein 1) ELISA Kit	CREB1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1440	Rabbit CXCR3(CXC-Chemokine Receptor 3) ELISA Kit	CXCR3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1441	Rabbit CLEC2C(C-Type Lectin Domain Family 2 Member C) ELISA Kit	CLEC2C	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1442	Rabbit SELE(E-selectin) ELISA Kit	SELE	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1443	Rabbit FBXO2(F-Box Protein 2) ELISA Kit	FBXO2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-RB1444	Rabbit Flt3(FMS Like Tyrosine Kinase 3) ELISA Kit	Flt3	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1445	Rabbit GATA1(GATA Binding Protein 1) ELISA Kit	GATA1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1446	Rabbit GATA2(GATA Binding Protein 2) ELISA Kit	GATA2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1447	Rabbit GATA5(GATA Binding Protein 5) ELISA Kit	GATA5	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1448	Rabbit GDI1(GDP Dissociation Inhibitor 1) ELISA Kit	GDI1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1449	Rabbit GAP(GTPase Activating Protein) ELISA Kit	GAP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1450	Rabbit RGS4(Regulator of G Protein Signaling 4) ELISA Kit	RGS4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1451	Rabbit SHH(Sonic Hedgehog Protein) ELISA Kit	SHH	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1452	Rabbit TDP43(Tar DNA Binding Protein 43KDa) ELISA Kit	TDP43	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1453	Rabbit TAF1(Transcription initiation factor TFIID subunit 1) ELISA Kit	TAF1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1454	Rabbit TAF2(TATA Box Binding Protein Associated Factor 2) ELISA Kit	TAF2	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1455	Rabbit TAF13(TATA Box Binding Protein Associated Factor 13) ELISA Kit	TAF13	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1456	Rabbit TLR2(Toll-Like Receptor 2) ELISA Kit	TLR2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1457	Rabbit TU3A(TU3A Protein) ELISA Kit	TU3A	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1458	Rabbit LAT(Linker For Activation of T-cell) ELISA Kit	LAT	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1459	Rabbit TALLA-1(T-cell Acute Lymphoblastic Leukemia Antigen) ELISA Kit	TALLA-1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1460	Rabbit OCT1(Octamer-binding protein 1) ELISA Kit	OCT1	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1461	Rabbit IL1R1(Interleukin 1 Receptor Type I) ELISA Kit	IL1R1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1462	Rabbit LCA/CD45(Leukocyte Common Antigen) ELISA Kit	LCA/CD45	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1463	Rabbit NSE(Neuron-Specific Enolase) ELISA Kit	NSE	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1464	Rabbit cPLA2(Cytosolic Phospholipase A2) ELISA Kit	cPLA2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1465	Rabbit MIgA(Membrane Immunoglobulin A) ELISA Kit	MIgA	Sandwich	16ng/mL	0.15 ng/mL	0.25-16ng/mL
FY-RB1466	Rabbit MIgD(Membrane Immunoglobulin D) ELISA Kit	MIgD	Sandwich	16ng/mL	0.15 ng/mL	0.25-16ng/mL
FY-RB1467	Rabbit MIgM(Membrane Immunoglobulin M) ELISA Kit	MIgM	Sandwich	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-RB1468	Rabbit VIM(Vimentin) ELISA Kit	VIM	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1469	Rabbit FOXP1(Forkhead Box Protein P1) ELISA Kit	FOXP1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1470	Rabbit FGFBP2(Fibroblast Growth Factor Binding Protein 2) ELISA Kit	FGFBP2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1471	Rabbit FRS2(Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit	FRS2	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1472	Rabbit VLA4(Very Late Appearing Antigen 4) ELISA Kit	VLA4	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB1473	Rabbit PRF1(Perforin 1) ELISA Kit	PRF1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1474	Rabbit FSH(Follicle Stimulating Hormone) ELISA Kit	FSH	Sandwich	100mIU/mL	0.94 mIU/mL	1.56-100mIU/mL
FY-RB1476	Rabbit PYGB(Glycogen Phosphorylase, Brain) ELISA Kit	PYGB	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1477	Rabbit MAOA(Type A Monoamine Oxidase) ELISA Kit	MAOA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1478	Rabbit MCP-2/CCL8(Monocyte Chemotactic Protein 2) ELISA Kit	MCP-2/CCL8	Sandwich	16pg/mL	0.15 pg/mL	0.25-16pg/mL
FY-RB1479	Rabbit SREBP(Sterol Regulatory Element Binding Proteins) ELISA Kit	SREBP	Sandwich	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-RB1480	Rabbit EMLIN 1(Elastin Microfibril Interface Located Protein 1) ELISA Kit	EMLIN 1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1481	Rabbit PKB(Protein Kinase B) ELISA Kit	PKB	Sandwich	20ng/ml	0.19 ng/ml	0.31-20ng/ml
FY-RB1482	Rabbit PKD2(Protein Kinase D2) ELISA Kit	PKD2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1483	Rabbit LMP7/PSMB8(Low Molecular-weight Protein/proteasome Beta-type Subunit) ELISA Kit	LMP7/PSMB8	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1484	Rabbit VS(Versican/PG-M) ELISA Kit	VS	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1485	Rabbit PYY(Peptide YY) ELISA Kit	PYY	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-RB1486	Rabbit DPYSL3(Dihydropyrimidinase Like Protein 3) ELISA Kit	DPYSL3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1487	Rabbit FDFT1(Farnesyl Diphosphate Farnesyltransferase 1) ELISA Kit	FDFT1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1488	Rabbit SP-A1(Pulmonary surfactant-associated protein A1) ELISA Kit	SP-A1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1489	Rabbit SPB(Pulmonary Surfactant Associated Protein B) ELISA Kit	SPB	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1490	Rabbit SPD(Pulmonary Surfactant Associated Protein D) ELISA Kit	SPD	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1491	Rabbit PRCP(Prolylcarboxypeptidase) ELISA Kit	PRCP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1492	Rabbit RPA1(Replication protein A1) ELISA Kit	RPA1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1493	Rabbit MBP-C/MBL2(Mannose Binding Protein C/Mannose Binding Lectin) ELISA	MBP-C/MBL2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1494	Rabbit RNASEH1(Ribonuclease H1) ELISA Kit	RNASEH1	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1495	Rabbit RNASET2(Ribonuclease T2) ELISA Kit	RNASET2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1496	Rabbit NFKB-p105(Nuclear factor NF-kappa-B p105 subunit) ELISA Kit	NFKB-p105	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1497	Rabbit NF-κB p65(Nuclear Factor Kappa B p65) ELISA Kit	NF-κB p65	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1498	Rabbit βMSH(beta-Melanocyte Stimulating Hormone) ELISA Kit	βMSH	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1499	Rabbit GRIK2(Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit	GRIK2	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1500	Rabbit CKMB(Creatine Kinase MB Isoenzyme) ELISA Kit	CKMB	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1501	Rabbit CNTN4(Contactin 4) ELISA Kit	CNTN4	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1502	Rabbit UPP2(Uridine Phosphorylase 2) ELISA Kit	UPP2	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1503	Rabbit MT1(Metallothionein 1) ELISA Kit	MT1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1504	Rabbit SSFA2(Sperm Specific Antigen 2) ELISA Kit	SSFA2	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1505	Rabbit τPK1(Tau-Protein Kinase 1) ELISA Kit	τPK1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1506	Rabbit TRACP-5b(Tartrate Resistant Acid Phosphatase 5b)ELISA Kit	TRACP-5b	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1507	Rabbit HMGB-1(High Mobility Group Protein B1) ELISA Kit	HMGB-1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1508	Rabbit PDGF-BB(Platelet Derived Growth Factor-BB) ELISA Kit	PDGF-BB	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB1509	Rabbit RBP4(Retinol Binding Protein 4) ELISA Kit	RBP4	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1510	Rabbit TGF-β2(Transforming Growth Factor β2) ELISA Kit	TGF-β2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1511	Rabbit TAP1(Transporter associated with Antigen Processing 1) ELISA Kit	TAP1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1512	Rabbit CD28(Cluster of Differentiation 28) ELISA Kit	CD28	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1513	Rabbit CD30L(Cluster of Differentiation 30 Ligand) ELISA Kit	CD30L	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB1514	Rabbit CD86(T-Lymphocyte Activation Antigen CD86) ELISA Kit	CD86	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1515	Rabbit sTNF(Soluble Tumor Necrosis Factor) ELISA Kit	sTNF	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1516	Rabbit PDGFA(Platelet Derived Growth Factor Subunit A) ELISA Kit	PDGFA	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1517	Rabbit CX37(Connexin 37) ELISA Kit	CX37	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1518	Rabbit PLAA(Phospholipase A2 Activating Protein) ELISA Kit	PLAA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1519	Rabbit PLA1(Phospholipase A1) ELISA Kit	PLA1	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB1520	Rabbit PLCβ1/PI-PLC(Phospholipase C Beta 1, Phosphoinositide Specific) ELISA Kit	PLCβ1/PI-PLC	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1523	Rabbit PGDS2(Prostaglandin D2 Synthase)ELISA Kit	PGDS2	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1524	Rabbit DEFα5(Defensin Alpha 5, Paneth Cell Specific) ELISA Kit	DEFα5	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB1525	Rabbit IgM(Immunoglobulin M) ELISA Kit	IgM	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB1526	Rabbit EP2(Prostaglandin E Receptor 2) ELISA Kit	EP2	Sandwich	25ng/mL	0.23 ng/mL	0.39-25ng/mL
FY-RB1527	Rabbit ACP3/PAP(Prostatic Acid Phosphatase) ELISA Kit	ACP3/PAP	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1528	Rabbit HSP-70/HSPA9(Heat Shock Protein 70) ELISA Kit	HSP-70/HSPA9	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL

FY-RB1529	Rabbit HSF1(Heat Shock Factor 1) ELISA Kit	HSF1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1530	Rabbit WARS(Tryptophanyl tRNA Synthetase) ELISA Kit	WARS	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1532	Rabbit NT-3(Neurotrophin-3) ELISA Kit	NT-3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1533	Rabbit GROβ/CXCL2(Growth Regulated Oncogene Beta) ELISA Kit	GROβ/CXCL2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1534	Rabbit NE/ELA2(Neutrophil Elastase/Elastase-2) ELISA Kit	NE/ELA2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1535	Rabbit STK39(Serine/Threonine Kinase 39) ELISA Kit	STK39	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1536	Rabbit DAP6(Death Associated Protein 6) ELISA Kit	DAP6	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB1537	Rabbit TPPP(Tubulin Polymerization Promoting Protein) ELISA Kit	TPPP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1538	Rabbit VIPR1(Vasoactive Intestinal Peptide Receptor 1) ELISA Kit	VIPR1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1539	Rabbit CTLA4(Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit	CTLA4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1540	Rabbit CK-17/KRT17(Cytokeratin 17) ELISA Kit	CK-17/KRT17	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1541	Rabbit CK-18/KRT18(Cytokeratin 18) ELISA Kit	CK-18/KRT18	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1542	Rabbit CYP7A1(Cytochrome P450 7A1) ELISA Kit	CYP7A1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1543	Rabbit CCND1(Cyclin D1) ELISA Kit	CCND1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1544	Rabbit CCND2(Cyclin D2) ELISA Kit	CCND2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1545	Rabbit FLNα(Filamin A, Alpha) ELISA Kit	FLNα	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1546	Rabbit Fbn(Fibrinolysin) ELISA Kit	Fbn	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1547	Rabbit SNRPD1(Small Nuclear Ribonucleoprotein Polypeptide D1) ELISA Kit	SNRPD1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1548	Rabbit STAT5A(Signal Transducer and Activator of Transcription 5A) ELISA Kit	STAT5A	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1549	Rabbit SOX2(Sex Determining Region Y Box Protein 2) ELISA Kit	SOX2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1551	Rabbit HbH(Hemoglobin H) ELISA Kit	HbH	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1552	Rabbit SEPP1(Selenoprotein P) ELISA Kit	SEPP1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1553	Rabbit BMG/β2-MG(Beta-2-Microglobulin) ELISA Kit	BMG/β2-MG	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1554	Rabbit PDGFC(Platelet Derived Growth Factor C) ELISA Kit	PDGFC	Sandwich	160pg/mL	1.5 pg/mL	2.5-160pg/mL
FY-RB1555	Rabbit PDGFD(Platelet Derived Growth Factor D) ELISA Kit	PDGFD	Sandwich	30pg/mL	0.28 pg/mL	0.47-30pg/mL
FY-RB1556	Rabbit NOX4(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit	NOX4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1557	Rabbit NOX5(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 5) ELISA Kit	NOX5	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1559	Rabbit FNDC5(Fibronectin type III domain-containing protein 5) ELISA Kit	FNDC5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1560	Rabbit MEC/CCL28(Mucosae Associated Epithelia Chemokine) ELISA Kit	MEC/CCL28	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1561	Rabbit LpPLA2(Phospholipase A2, Lipoprotein Associated) ELISA Kit	LpPLA2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1562	Rabbit TACE/ADAM17(TNF α Converting Enzyme) ELISA Kit	TACE/ADAM17	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1563	Rabbit TSGF(Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit	TSGF	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1564	Rabbit TFR1(Transferrin Receptor 1) ELISA Kit	TFR1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1565	Rabbit CD8(Cluster of Differentiation 8) ELISA Kit	CD8	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1566	Rabbit SEPCR/PROCR(Soluble Endothelial Protein C Receptor) ELISA Kit	SEPCR/PROCR	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1567	Rabbit PLA2G2A(Phospholipase A2, Group II A) ELISA Kit	PLA2G2A	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1568	Rabbit PS1(Presenilin 1) ELISA Kit	PS1	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB1569	Rabbit MAP-2(Microtubule Associated Protein 2) ELISA Kit	MAP-2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1570	Rabbit OPRM1(Mu-type opioid receptor) ELISA Kit	OPRM1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1571	Rabbit INS(Insulin) ELISA Kit	INS	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1572	Rabbit TSH(Thyroid Stimulating Hormone) ELISA Kit	TSH	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL

FY-RB1573	Rabbit EPO(Erythropoietin) ELISA Kit	EPO	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1574	Rabbit TWEAK(Tumour Necrosis Factor Related Weak Inducer of Apoptosis) ELISA Kit	TWEAK	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1575	Rabbit APT(Abnormal Prothrombin) ELISA Kit	APT	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1576	Rabbit RANKL(Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit	RANKL	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1577	Rabbit CP(Ceruloplasmin) ELISA Kit	CP	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-RB1578	Rabbit NT-ProBNP(N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit	NT-ProBNP	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1579	Rabbit PTX 3/TSG-14(Pentraxin 3) ELISA Kit	PTX 3/TSG-14	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1580	Rabbit IL-6(Interleukin 6) ELISA Kit	IL-6	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB1581	Rabbit IL-10(Interleukin 10) ELISA Kit	IL-10	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB1582	Rabbit GDF2(Growth Differentiation Factor 2) ELISA Kit	GDF2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1583	Rabbit NRP1(Neuropilin 1) ELISA Kit	NRP1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB1585	Rabbit DBP(Vitamin D Binding Protein) ELISA Kit	DBP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1586	Rabbit α -Glu(Alpha Glucosidase) ELISA Kit	α -Glu	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1587	Rabbit MIF(Macrophage Migration Inhibitory Factor) ELISA Kit	MIF	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1588	Rabbit ENG(Endoglin) ELISA Kit	ENG	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1589	Rabbit LEPR(Leptin Receptor) ELISA Kit	LEPR	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1590	Rabbit VEGFR1/FLT1(Vasculoar Endothelial Growth Factor Receptor 1) ELISA Kit	VEGFR1/FLT1	Sandwich	300pg/mL	2.81 pg/mL	4.69-300pg/mL
FY-RB1591	Rabbit ACE(Angiotensin I Converting Enzyme) ELISA Kit	ACE	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1592	Rabbit RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	RANTES	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1593	Rabbit SELP(P-Selectin) ELISA Kit	SELP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1594	Rabbit PCSK9(Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit	PCSK9	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1595	Rabbit PDGF-AA(Platelet Derived Growth Factor AA) ELISA Kit	PDGF-AA	Sandwich	2000pg/ml	18.75 pg/ml	31.25-2000pg/ml
FY-RB1596	Rabbit PF4(Platelet Factor 4) ELISA Kit	PF4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1597	Rabbit ACV-A(Activin A) ELISA Kit	ACV-A	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1598	Rabbit FRS3(Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit	FRS3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1599	Rabbit GDF5(Growth Differentiation Factor 5) ELISA Kit	GDF5	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1600	Rabbit PAL(L-Phenylalanine Ammonia-Lyase) ELISA Kit	PAL	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-RB1601	Rabbit NUP153(Nucleoporin 153kDa) ELISA Kit	NUP153	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1602	Rabbit NUP160(Nucleoporin 160kDa) ELISA Kit	NUP160	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1603	Rabbit PSA(Prostate Specific Antigen) ELISA Kit	PSA	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1604	Rabbit TR(Thrombin Receptor) ELISA Kit	TR	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1605	Rabbit α 2-PI(α 2-Plasmin Inhibitor) ELISA Kit	α 2-PI	Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-RB1606	Rabbit ACT α 2(Actin Alpha 2, Smooth Muscle) ELISA Kit	ACT α 2	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB1607	Rabbit APC(Activated Protein C) ELISA Kit	APC	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1608	Rabbit Cys-C(Cystatin C) ELISA Kit	Cys-C	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1609	Rabbit F2(Coagulation Factor II) ELISA Kit	F2	Sandwich	24ng/mL	0.23 ng/mL	0.38-24ng/mL
FY-RB1610	Rabbit CCND3(Cyclin D3) ELISA Kit	CCND3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1611	Rabbit TGF- α (Transforming Growth Factor α) ELISA Kit	TGF- α	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1612	Rabbit TGF- β 3(Transforming Growth Factor Beta 3) ELISA Kit	TGF- β 3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1613	Rabbit Collagenase I ELISA Kit	Collagenase I	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1614	Rabbit MAT(Matrilysin) ELISA Kit	MAT	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL

FY-RB1615	Rabbit CD68(Cluster of Differentiation 68) ELISA Kit	CD68	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1616	Rabbit ACTα1(Actin Alpha 1, Skeletal Muscle) ELISA Kit	ACTα1	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB1617	Rabbit HSPG(Heparan Sulfate Proteoglycan) ELISA Kit	HSPG	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1618	Rabbit Gal(α-Galactosyl Epitopes) ELISA Kit	Gal	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1619	Rabbit TM(Thrombomodulin) ELISA Kit	TM	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1620	Rabbit ANG(Angiostatin) ELISA Kit	ANG	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1621	Rabbit MYC(V-Myc Myelocytomatosis Viral Oncogene Homolog) ELISA Kit	MYC	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1622	Rabbit PL(Pancreatic Lipase) ELISA Kit	PL	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB1623	Rabbit PP(Pancreatic Polypeptide) ELISA Kit	PP	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1624	Rabbit GS(Gelsolin) ELISA Kit	GS	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-RB1625	Rabbit GR(Glutathione Reductase) ELISA Kit	GR	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1626	Rabbit TPO(Thrombopoietin) ELISA Kit	TPO	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1627	Rabbit TG(Thyroglobulin) ELISA Kit	TG	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB1628	Rabbit PROS(Protein S) ELISA Kit	PROS	Sandwich	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-RB1629	Rabbit ucMGP(undercarboxylated Matrix Gla Protein) ELISA Kit	ucMGP	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1630	Rabbit PDH E1(Pyruvate Dehydrogenase E1) ELISA Kit	PDH E1	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-RB1631	Rabbit DBP(DNA Binding Protein) ELISA Kit	DBP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1632	Rabbit BTA(Bladder Tumor Antigen) ELISA Kit	BTA	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB1633	Rabbit CASP11(Caspase 11) ELISA Kit	CASP11	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1634	Rabbit GAD(Glutamic Acid Decarboxylase) ELISA Kit	GAD	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB1635	Rabbit PROZ(Protein Z) ELISA Kit	PROZ	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1636	Rabbit FSA(Fetal Sulfoslycoprotein Antigen) ELISA Kit	FSA	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1637	Rabbit MX1(Myxovirus Resistance 1) ELISA Kit	MX1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1638	Rabbit MFAP4(Microfibrillar Associated Protein 4) ELISA Kit	MFAP4	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-RB1639	Rabbit VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit	VEGF-A	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-RB0065	Rabbit ACTH(Adrenocorticotropic Hormone) ELISA Kit	ACTH	Competitive	800pg/mL	7.5 pg/mL	12.5-800pg/mL
FY-RB0071	Rabbit α-CTx(Alpha Crosslaps) ELISA Kit	α-CTx	Competitive	60ng/mL	0.56 ng/mL	0.94-60ng/mL
FY-RB0099	Rabbit 11-DH-TXB2(11-Dehydro-thromboxane B2) ELISA Kit	11-DH-TXB2	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0123	Rabbit ADH(Antidiuretic Hormone) ELISA Kit	ADH	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0150	Rabbit IAPP(Islet Amyloid Polypeptide) ELISA Kit	IAPP	Competitive	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0194	Rabbit GC(Glucagon) ELISA Kit	GC	Competitive	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0229	Rabbit APLN(Apelin) ELISA Kit	APLN	Competitive	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-RB0272	Rabbit ANP(Atrial Natriuretic Peptide) ELISA Kit	ANP	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0289	Rabbit β-EP(Beta-Endorphin) ELISA Kit	β-EP	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0298	Rabbit BK(Bradykinin) ELISA Kit	BK	Competitive	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-RB0304	Rabbit CNP(C-type Natriuretic Peptide) ELISA Kit	CNP	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0473	Rabbit Salusin α(Salusin Alpha) ELISA Kit	Salusin α	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0490	Rabbit β-CTx(Beta Crosslaps) ELISA Kit	β-CTx	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB0517	Rabbit OXA(Orexin A) ELISA Kit	OXA	Competitive	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-RB0545	Rabbit NP(Nandrolone Phenylpropionate) ELISA Kit	NP	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-RB0569	Rabbit SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit	SOD1	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-RB0570	Rabbit OGP(Osteogenic Growth Peptide) ELISA Kit	OGP	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB0583	Rabbit TRH(Thyrotropin Releasing Hormone) ELISA Kit	TRH	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0585	Rabbit δ SIP- α (Delta Sleep Inducing Peptide α) ELISA Kit	δ SIP- α	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB0587	Rabbit SCT(Secretin) ELISA Kit	SCT	Competitive	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB0692	Rabbit OFQ/N(Orphanin FQ/Nociceptin) ELISA Kit	OFQ/N	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0882	Rabbit MPA(Mycophenolic Acid) ELISA Kit	MPA	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB0932	Rabbit PGF2 α (Prostaglandin F2 α) ELISA Kit	PGF2 α	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB0938	Rabbit Dyn(Big Dynorphin) ELISA Kit	Dyn	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-RB0954	Rabbit NKA(Neurokinin A) ELISA Kit	NKA	Competitive	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-RB0955	Rabbit NKB(Neurokinin B) ELISA Kit	NKB	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-RB0958	Rabbit NPY(Neuropeptide Y) ELISA Kit	NPY	Competitive	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB0976	Rabbit SST(Somatostatin) ELISA Kit	SST	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB0982	Rabbit OXB(Orexin B) ELISA Kit	OXB	Competitive	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1023	Rabbit MT(Melatonin) ELISA Kit	MT	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1024	Rabbit DHEA-S(Dehydroepiandrosterone Sulfate) ELISA Kit	DHEA-S	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1025	Rabbit DPD(Deoxypyridinoline) ELISA Kit	DPD	Competitive	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-RB1027	Rabbit SA(Sialic Acid) ELISA Kit	SA	Competitive	500 μ g/mL	4.69 μ g/mL	7.81-500 μ g/mL
FY-RB1040	Rabbit MTL(Motilin) ELISA Kit	MTL	Competitive	60ng/mL	0.56 ng/mL	0.94-60ng/mL
FY-RB1088	Rabbit TMS β 10(Thymosin Beta 10) ELISA Kit	TMS β 10	Competitive	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-RB1089	Rabbit TP5(Thymopentin) ELISA Kit	TP5	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1090	Rabbit VIP(Vasoactive Intestinal Peptide) ELISA Kit	VIP	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-RB1104	Rabbit TXB2(Thromboxane B2) ELISA Kit	TXB2	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-RB1181	Rabbit SOD3(Superoxide Dismutase 3, Extracellular) ELISA Kit	SOD3	Competitive	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1207	Rabbit α MSH(Alpha-Melanocyte Stimulating Hormone) ELISA Kit	α MSH	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1355	Rabbit IgG(Immunoglobulin G) ELISA Kit	IgG	Competitive	40 μ g/mL	0.38 μ g/mL	0.63-40 μ g/mL
FY-RB1357	Rabbit AP17(Apelin 17) ELISA Kit	AP17	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-RB1365	Rabbit Ang- II (Angiotensin II) ELISA Kit	Ang- II	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1366	Rabbit Ang- I (Angiotensin I) ELISA Kit	Ang- I	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-RB1377	Rabbit AP13(Apelin 13) ELISA Kit	AP13	Competitive	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1378	Rabbit AP36(Apelin 36) ELISA Kit	AP36	Competitive	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-RB1475	Rabbit CRH(Corticotropin Releasing Hormone) ELISA Kit	CRH	Competitive	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-RB1521	Rabbit MS(Melatonin Sulfate) ELISA Kit	MS	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-RB1522	Rabbit EM2(Endomorphin 2) ELISA Kit	EM2	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-RB1531	Rabbit NT(Neurotensin) ELISA Kit	NT	Competitive	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-RB1550	Rabbit TMS β 4(Thymosin Beta 4) ELISA Kit	TMS β 4	Competitive	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-RB1558	Rabbit TAP(Trypsinogen Activation Peptide) ELISA Kit	TAP	Competitive	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-RB1584	Rabbit BNP(Brain Natriuretic Peptide) ELISA Kit	BNP	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0001	Chicken IgY(Immunoglobulin Y) ELISA Kit	IgY	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0002	Chicken ORM1(Orosomucoid 1) ELISA Kit	ORM1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0003	Chicken α -Glu(Alpha Glucosidase) ELISA Kit	α -Glu	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0004	Chicken AT-III(Antithrombin III) ELISA Kit	AT-III	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-Ch0005	Chicken bFGF/FGF2(Basic Fibroblast Growth Factor) ELISA Kit	bFGF/FGF2	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0006	Chicken EGF(Epidermal Growth Factor) ELISA Kit	EGF	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0007	Chicken CCL1(Chemokine C-C-Motif Ligand 1) ELISA Kit	CCL1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0008	Chicken MCSP(Melanoma Associated Chondroitin Sulfate Proteoglycan) ELISA Kit	MCSP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0009	Chicken PLSCR2(Phospholipid Scramblase 2) ELISA Kit	PLSCR2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0011	Chicken FLNβ(Filamin B, Beta) ELISA Kit	FLNβ	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0012	Chicken IL-4(Interleukin 4) ELISA Kit	IL-4	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0013	Chicken MYL9(Myosin Light Chain 9, Regulatory) ELISA Kit	MYL9	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0014	Chicken C/EBPε(CCAAT/Enhancer Binding Protein, Epsilon) ELISA Kit	C/EBPε	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0015	Chicken GUSβ(β-glucuronidase) ELISA Kit	GUSβ	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0016	Chicken HRG(Histidine Rich Glycoprotein) ELISA Kit	HRG	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0017	Chicken GPC2(Glypican 2) ELISA Kit	GPC2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0018	Chicken GzmA(Granzyme A) ELISA Kit	GzmA	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0019	Chicken IFN-γ(Interferon Gamma) ELISA Kit	IFN-γ	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0020	Chicken ADAM10(A Disintegrin And Metalloprotease 10) ELISA Kit	ADAM10	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0021	Chicken TMPRSS2(Transmembrane Protease, Serine 2) ELISA Kit	TMPRSS2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0022	Chicken CLCF1(Cardiotrophin Like Cytokine Factor 1) ELISA Kit	CLCF1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0023	Chicken PP(Pepsin) ELISA Kit	PP	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0024	Chicken DDT(D-Dopachrome Tautomerase) ELISA Kit	DDT	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0025	Chicken BAG3(Bcl2 Associated Athanogene 3) ELISA Kit	BAG3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0026	Chicken BMP-6(Bone Morphogenetic Protein 6) ELISA Kit	BMP-6	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0027	Chicken PCX(Podocalyxin) ELISA Kit	PCX	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0028	Chicken TXLNa(Taxilin Alpha) ELISA Kit	TXLNa	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0029	Chicken SDC4(Syndecan 4) ELISA Kit	SDC4	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0030	Chicken Grb14(Growth Factor Receptor Bound Protein 14) ELISA Kit	Grb14	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0031	Chicken FcγR3B(Fc Fragment of IgG Low Affinity IIIb Receptor) ELISA Kit	FcγR3B	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0032	Chicken CKLFSF8(Chemokine Like Factor Superfamily 8) ELISA Kit	CKLFSF8	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0033	Chicken HbH(Hemoglobin H) ELISA Kit	HbH	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0034	Chicken ASGR2(Asialoglycoprotein Receptor 2) ELISA Kit	ASGR2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0035	Chicken GUCA2A(Guanylate Cyclase Activator 2A) ELISA Kit	GUCA2A	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0036	Chicken TSP-4(Thrombospondin-4) ELISA Kit	TSP-4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0037	Chicken CYSLTR1(CysteinyL Leukotriene Receptor 1) ELISA Kit	CYSLTR1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0038	Chicken PCDHβ16(Protocadherin Beta 16) ELISA Kit	PCDHβ16	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0039	Chicken MLCK/MYLK(Myosin Light Chain Kinase) ELISA Kit	MLCK/MYLK	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0040	Chicken PLA2G4D(Phospholipase A2, Group IV D) ELISA Kit	PLA2G4D	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0041	Chicken CCND2(Cyclin D2) ELISA Kit	CCND2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0042	Chicken APT(Abnormal Prothrombin) ELISA Kit	APT	Sandwich	10ng/ml	0.1 ng/ml	0.16-10ng/ml
FY-Ch0043	Chicken PDGF-BB(Platelet Derived Growth Factor-BB) ELISA Kit	PDGF-BB	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0044	Chicken GPER1(G Protein Coupled Estrogen Receptor 1) ELISA Kit	GPER1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0045	Chicken PDGFD(Platelet Derived Growth Factor D) ELISA Kit	PDGFD	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-Ch0046	Chicken ADP/Acrp30(Adiponectin) ELISA Kit	ADP/Acrp30	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL

FY-Ch0047	Chicken GAL(Galanin) ELISA Kit	GAL	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0048	Chicken GAS(N-Acetylgalactosamine 6-Sulfatase) ELISA Kit	GAS	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0049	Chicken CKMB(Creatine Kinase MB Isoenzyme) ELISA Kit	CKMB	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0050	Chicken MEC/CCL28(Mucosae Associated Epithelia Chemokine) ELISA Kit	MEC/CCL28	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0051	Chicken HPR(Haptoglobin Related Protein) ELISA Kit	HPR	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0052	Chicken GREM1(Gremlin 1) ELISA Kit	GREM1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0053	Chicken CPN10(Chaperonin 10) ELISA Kit	CPN10	Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-Ch0054	Chicken ALCAM(Activated Leukocyte Cell Adhesion Molecule) ELISA Kit	ALCAM	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch0055	Chicken MAT(Matrilysin) ELISA Kit	MAT	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0056	Chicken SYK(Spleen Tyrosine Kinase) ELISA Kit	SYK	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0057	Chicken TP(Thymidine Phosphorylase) ELISA Kit	TP	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0058	Chicken GATA3(GATA Binding Protein 3) ELISA Kit	GATA3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0059	Chicken GPLD1(Glycosylphosphatidylinositol Specific Phospholipase D1) ELISA Ki	GPLD1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0060	Chicken COL9a3(Collagen Type IX Alpha 3) ELISA Kit	COL9a3	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0061	Chicken LTβR(Lymphotoxin Beta Receptor) ELISA Kit	LTβR	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0062	Chicken PDGFC(Platelet Derived Growth Factor C) ELISA Kit	PDGFC	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0063	Chicken SDH(Sorbitol Dehydrogenase) ELISA Kit	SDH	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0064	Chicken PLCγ2(Phospholipase C Gamma 2, Phosphatidylinositol Specific) ELISA Ki	PLCγ2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0065	Chicken MCP-4/CCL13(Monocyte Chemotactic Protein 4) ELISA Kit	MCP-4/CCL13	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0066	Chicken FXIIa(Activated Coagulation Factor XII) ELISA Kit	FXIIa	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0067	Chicken NUP35(Nucleoporin 35kDa) ELISA Kit	NUP35	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0068	Chicken CLDN3(Claudin 3) ELISA Kit	CLDN3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0069	Chicken CD40L(Cluster of Differentiation 40 Ligand) ELISA Kit	CD40L	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0070	Chicken PKCγ(Protein Kinase C Gamma) ELISA Kit	PKCγ	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0071	Chicken CAM(Calmodulin) ELISA Kit	CAM	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0072	Chicken GSTα1(Glutathione S Transferase Alpha 1) ELISA Kit	GSTα1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0073	Chicken GROα/CXCL1(Growth Regulated Oncogene Alpha) ELISA Kit	GROα/CXCL1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0074	Chicken AQP-1(Aquaporin 1) ELISA Kit	AQP-1	Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-Ch0075	Chicken HbC(Hemoglobin C) ELISA Kit	HbC	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0076	Chicken SREBP(Sterol Regulatory Element Binding Proteins) ELISA Kit	SREBP	Sandwich	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-Ch0077	Chicken PAH(Phenylalanine Hydroxylase) ELISA Kit	PAH	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0078	Chicken P-Cadherin(Cadherin, Placental) ELISA Kit	P-Cadherin	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0079	Chicken BMPR2(Bone Morphogenetic Protein Receptor II) ELISA Kit	BMPR2	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0080	Chicken NUP50(Nucleoporin 50kDa) ELISA Kit	NUP50	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0081	Chicken FGL2(Fibrinogen Like Protein 2) ELISA Kit	FGL2	Sandwich	150pg/mL	1.4 pg/mL	2.34-150pg/mL
FY-Ch0082	Chicken IGF-1(Insulin-like Growth Factor 1) ELISA Kit	IGF-1	Sandwich	1200pg/mL	11.25 pg/mL	18.75-1200pg/mL
FY-Ch0083	Chicken MMP-1(Matrix Metalloproteinase 1) ELISA Kit	MMP-1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0084	Chicken ApoA1(Apolipoprotein A1) ELISA Kit	ApoA1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0085	Chicken FGF19(Fibroblast Growth Factor 19) ELISA Kit	FGF19	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0086	Chicken IL-2(Interleukin 2) ELISA Kit	IL-2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0088	Chicken REG3γ(Regenerating Islet Derived Protein 3 Gamma) ELISA Kit	REG3γ	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL

FY-Ch0089	Chicken MPIF2(Myeloid Progenitor Inhibitory Factor 2) ELISA Kit	MPIF2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0090	Chicken GDI1(GDP Dissociation Inhibitor 1) ELISA Kit	GDI1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0091	Chicken CLDN1(Claudin 1) ELISA Kit	CLDN1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0092	Chicken Ank-B(Ankyrin B) ELISA Kit	Ank-B	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0093	Chicken IFN α 1(Interferon Alpha 1) ELISA Kit	IFN α 1	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0094	Chicken CK-13/KRT13(Cytokeratin 13) ELISA Kit	CK-13/KRT13	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0095	Chicken PKR2(Prokineticin Receptor 2) ELISA Kit	PKR2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0096	Chicken FAK(Focal Adhesion Kinase) ELISA Kit	FAK	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0097	Chicken BACE2(Beta-Site APP Cleaving Enzyme 2) ELISA Kit	BACE2	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch0098	Chicken GCSFR(Colony Stimulating Factor Receptor, Granulocyte) ELISA Kit	GCSFR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0099	Chicken ELMO2(Engulfment And Cell Motility 2) ELISA Kit	ELMO2	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0100	Chicken TFPI(Tissue Factor Pathway Inhibitor) ELISA Kit	TFPI	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-Ch0101	Chicken PEDF(Pigment Epithelium Derived Factor) ELISA Kit	PEDF	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0102	Chicken REG1 α (Regenerating Islet Derived Protein 1 Alpha) ELISA Kit	REG1 α	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0103	Chicken PP13(Placental Protein13) ELISA Kit	PP13	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0104	Chicken EPCR(Endothelial Protein C Receptor) ELISA Kit	EPCR	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0105	Chicken PHLDA2(Pleckstrin Homology Like Domain Family A, Member 2) ELISA Kit	PHLDA2	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0106	Chicken SNC γ (Synuclein Gamma) ELISA Kit	SNC γ	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0107	Chicken Ub(Ubiquitin) ELISA Kit	Ub	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0108	Chicken PMP3(Peroxisomal Membrane Protein 3) ELISA Kit	PMP3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0109	Chicken OIT3(Oncoprotein Induced Transcript 3) ELISA Kit	OIT3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0110	Chicken TIMP-3(Tissue Inhibitors of Metalloproteinase 3) ELISA Kit	TIMP-3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0111	Chicken RAPGEF1(Rap Guanine Nucleotide Exchange Factor 1) ELISA Kit	RAPGEF1	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0112	Chicken EPYC(Epiphycan) ELISA Kit	EPYC	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0113	Chicken PPAR γ C1 α (Peroxisome Proliferator Activated Receptor Gamma Coactivator 1) ELISA Kit	PPAR γ C1 α	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0114	Chicken THR α (Thyroid Hormone Receptor Alpha) ELISA Kit	THR α	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0115	Chicken ACTN3(Actinin Alpha 3) ELISA Kit	ACTN3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0117	Chicken EMMPRIN/CD147(Extracellular Matrix Metalloproteinase Inducer) ELISA Kit	EMMPRIN/CD147	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0118	Chicken GCP-2(Granulocyte Chemotactic Protein 2) ELISA Kit	GCP-2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0119	Chicken ACT α 2(Actin Alpha 2, Smooth Muscle) ELISA Kit	ACT α 2	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch0120	Chicken GasR(Gastrin Receptor) ELISA Kit	GasR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0121	Chicken DLAT(Dihydrolipoyl Transacetylase) ELISA Kit	DLAT	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0122	Chicken HSP-40(Heat Shock Protein 40) ELISA Kit	HSP-40	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0123	Chicken PICP(Procollagen I C-terminal Propeptide) ELISA Kit	PICP	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0124	Chicken STAT2(Signal Transducer and Activator of Transcription 2) ELISA Kit	STAT2	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0125	Chicken EMA(Epithelial Membrane Antigen) ELISA Kit	EMA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0126	Chicken CDK2(Cyclin Dependent Kinase 2) ELISA Kit	CDK2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0127	Chicken TNPO1(Transportin 1) ELISA Kit	TNPO1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0128	Chicken CDK9(Cyclin Dependent Kinase 9) ELISA Kit	CDK9	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0129	Chicken PRKG1(Protein Kinase, cGMP Dependent Type 1) ELISA Kit	PRKG1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0130	Chicken IGF-2(Insulin Like Growth Factor 2) ELISA Kit	IGF-2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-Ch0131	Chicken F13(Coagulation Factor X III) ELISA Kit	F13	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-Ch0133	Chicken NRP1(Neuropilin 1) ELISA Kit	NRP1	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0134	Chicken STAT5A(Signal Transducer and Activator of Transcription 5A) ELISA Kit	STAT5A	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0135	Chicken HYOU1(Hypoxia Up Regulated 1) ELISA Kit	HYOU1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0136	Chicken MLZE(Melanoma Derived Leucine Zipper Extra Nuclear Factor) ELISA Kit	MLZE	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0137	Chicken CPR(Cytochrome P450 Reductase) ELISA Kit	CPR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0138	Chicken CDKN3(Cyclin Dependent Kinase Inhibitor 3) ELISA Kit	CDKN3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0139	Chicken CXCR1(CXC-Chemokine Receptor 1) ELISA Kit	CXCR1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0140	Chicken Ntn1(Netrin 1) ELISA Kit	Ntn1	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0141	Chicken γ GT1(Gamma Glutamyltransferase 1) ELISA Kit	γ GT1	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0142	Chicken GAL12(Galectin 12) ELISA Kit	GAL12	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0143	Chicken BMF(Bcl-2 Modifying Factor) ELISA Kit	BMF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0144	Chicken GMFB(Glia Maturation Factor, Beta) ELISA Kit	GMFB	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-Ch0145	Chicken EPHA3(Ephrin Type A Receptor 3) ELISA Kit	EPHA3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0147	Chicken CCND3(Cyclin D3) ELISA Kit	CCND3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0148	Chicken DPT(Dermatopontin) ELISA Kit	DPT	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0149	Chicken ADAM8(A Disintegrin And Metalloprotease 8) ELISA Kit	ADAM8	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0150	Chicken NEF3(Neurofilament 3) ELISA Kit	NEF3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0151	Chicken LUM(Lumican) ELISA Kit	LUM	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0152	Chicken GP-73(Golgi Protein 73) ELISA Kit	GP-73	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0153	Chicken COLEC11(Collectin Sub-Family Member 11) ELISA Kit	COLEC11	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0154	Chicken PDXK(Pyridoxal Kinase) ELISA Kit	PDXK	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0155	Chicken TNF- α (Tumor Necrosis Factor Alpha) ELISA Kit	TNF- α	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0156	Chicken GST ω 1(Glutathione S Transferase Omega 1) ELISA Kit	GST ω 1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0157	Chicken GDF11(Growth Differentiation Factor 11) ELISA Kit	GDF11	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0158	Chicken CD30L(Cluster of Differentiation 30 Ligand) ELISA Kit	CD30L	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0159	Chicken TGM2(Transglutaminase 2, Tissue) ELISA Kit	TGM2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0160	Chicken FBL(Fibrillarlin) ELISA Kit	FBL	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0161	Chicken RGS4(Regulator of G Protein Signaling 4) ELISA Kit	RGS4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0162	Chicken C5b(Complement Fragment 5b) ELISA Kit	C5b	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0163	Chicken GRO γ /CXCL3(Growth Regulated Oncogene Gamma) ELISA Kit	GRO γ /CXCL3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0164	Chicken RUNX2(Runt-Related Transcription Factor 2) ELISA Kit	RUNX2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0165	Chicken IL-6(Interleukin 6) ELISA Kit	IL-6	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0166	Chicken FZD8(Frizzled Homolog 8) ELISA Kit	FZD8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0167	Chicken sgp130(soluble Glucoprotein 130) ELISA Kit	sgp130	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0168	Chicken PKN1(Protein Kinase N1) ELISA Kit	PKN1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0169	Chicken β -LPH(β -Lipotrophic Hormone) ELISA Kit	β -LPH	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0170	Chicken Smad7/MADH7(Mothers Against Decapentaplegic Homolog 7) ELISA Kit	Smad7/MADH7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0171	Chicken RHOA(Ras Homolog Gene Family, Member A) ELISA Kit	RHOA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0172	Chicken TU3A(TU3A Protein) ELISA Kit	TU3A	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0173	Chicken MAN1A1(Mannosidase Alpha Class 1A Member 1) ELISA Kit	MAN1A1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-Ch0174	Chicken EPHA1(Ephrin Type A Receptor 1) ELISA Kit	EPHA1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0175	Chicken CT1(Cardiotrophin 1) ELISA Kit	CT1	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0176	Chicken GP5(Glycoprotein V, Platelet) ELISA Kit	GP5	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0177	Chicken GP4/CD36(Platelet Membrane Glycoprotein IV) ELISA Kit	GP4/CD36	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0178	Chicken ON(Osteonectin) ELISA Kit	ON	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0179	Chicken PPAR-γ(Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit	PPAR-γ	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0180	Chicken Gal(α-Galactosyl Epitopes) ELISA Kit	Gal	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0182	Chicken POMC(Pro-Opiomelanocortin) ELISA Kit	POMC	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0183	Chicken ApoD(Apolipoprotein D) ELISA Kit	ApoD	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0184	Chicken IL-13(Interleukin 13) ELISA Kit	IL-13	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0185	Chicken NUP153(Nucleoporin 153kDa) ELISA Kit	NUP153	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0186	Chicken CKLFSF7(Chemokine Like Factor Superfamily 7) ELISA Kit	CKLFSF7	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0187	Chicken SRF/IF(Skin Reactive Factor/Inflammatory Factor) ELISA Kit	SRF/IF	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0188	Chicken LTβ(Lymphotoxin Beta) ELISA Kit	LTβ	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0189	Chicken CDK7(Cyclin Dependent Kinase 7) ELISA Kit	CDK7	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0190	Chicken TrxR1(Thioredoxin Reductase 1) ELISA Kit	TrxR1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0191	Chicken GKR(Protein Kinase C-δ) ELISA Kit	GKR	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0192	Chicken TACR1(Tachykinin Receptor 1) ELISA Kit	TACR1	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-Ch0193	Chicken C-P(C-Peptide) ELISA Kit	C-P	Sandwich	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-Ch0194	Chicken GTF2F1(General Transcription Factor II F, Polypeptide 1) ELISA Kit	GTF2F1	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0195	Chicken IDE(Insulin Degrading Enzyme) ELISA Kit	IDE	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0196	Chicken CYPA(Cyclophilin A) ELISA Kit	CYPA	Sandwich	3000ng/mL	28.13 ng/mL	46.88-3000ng/mL
FY-Ch0197	Chicken FABP1(Fatty Acid Binding Protein 1, Liver) ELISA Kit	FABP1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0198	Chicken ECFR(Eosinophil Chemotactic Factor Receptor) ELISA Kit	ECFR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0199	Chicken SHBG(Sex Hormone Binding Globulin) ELISA Kit	SHBG	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0200	Chicken D1R/DRD1(Dopamine Receptor D1) ELISA Kit	D1R/DRD1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0201	Chicken MIF(Macrophage Migration Inhibitory Factor) ELISA Kit	MIF	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0202	Chicken PLCγ1(Phospholipase C Gamma 1) ELISA Kit	PLCγ1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0203	Chicken MYC(V-Myc Myelocytomatosis Viral Oncogene Homolog) ELISA Kit	MYC	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0204	Chicken CX3CR1(Chemokine C-X3-C-Motif Receptor 1) ELISA Kit	CX3CR1	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0206	Chicken ARHGDIα(Rho GDP Dissociation Inhibitor Alpha) ELISA Kit	ARHGDIα	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0207	Chicken MSTN(Myostatin) ELISA Kit	MSTN	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0208	Chicken DBP(Vitamin D Binding Protein) ELISA Kit	DBP	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-Ch0209	Chicken CGRP1(Calcitonin Gene Related Peptide 1) ELISA Kit	CGRP1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0210	Chicken TG(Thyroglobulin) ELISA Kit	TG	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0211	Chicken TLN(Talin) ELISA Kit	TLN	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0213	Chicken SEMA4C(Semaphorin 4C) ELISA Kit	SEMA4C	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0214	Chicken LRP(Lung Resistance-Related Protein) ELISA Kit	LRP	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0215	Chicken PDGFRα(Platelet Derived Growth Factor Receptor Alpha) ELISA Kit	PDGFRα	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0216	Chicken TDT(Terminal Deoxynucleotidyl Transferase) ELISA Kit	TDT	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0217	Chicken PKBG(Protein Kinase B Gamma) ELISA Kit	PKBG	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL

FY-Ch0218	Chicken CA5B(Carbonic Anhydrase VB, Mitochondrial) ELISA Kit	CA5B	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0219	Chicken COL14(Collagen Type X IV) ELISA Kit	COL14	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0220	Chicken C1q(Complement 1q) ELISA Kit	C1q	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0221	Chicken C1r(Complement Component 1R) ELISA Kit	C1r	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0222	Chicken BMP-4(Bone Morphogenetic Protein 4) ELISA Kit	BMP-4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0223	Chicken GDH/GLDH(Glutamate dehydrogenase) ELISA Kit	GDH/GLDH	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0224	Chicken DAO(Diamine Oxidase) ELISA Kit	DAO	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch0225	Chicken MMP-25(Matrix Metalloproteinase 25) ELISA Kit	MMP-25	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0226	Chicken GDF5(Growth Differentiation Factor 5) ELISA Kit	GDF5	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0227	Chicken PACAPR(Pituitary Adenylate Cyclase Activating Polypeptide Receptor 1) E	PACAPR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0228	Chicken PAQR7(Progestin and AdipoQ Receptor Family Member VII) ELISA Kit	PAQR7	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0229	Chicken NUP214(Nucleoporin 214kDa) ELISA Kit	NUP214	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0230	Chicken XCR1(XC-Chemokine Receptor 1) ELISA Kit	XCR1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0231	Chicken AATF(Apoptosis Antagonizing Transcription Factor) ELISA Kit	AATF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0232	Chicken DAT(Dopamine Transporter) ELISA Kit	DAT	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0233	Chicken FGF6(Fibroblast Growth Factor 6) ELISA Kit	FGF6	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0234	Chicken CK-19/KRT19(Cytokeratin 19) ELISA Kit	CK-19/KRT19	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0235	Chicken NRG-3(Neuregulin 3) ELISA Kit	NRG-3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0236	Chicken ACD(Adrenocortical Dysplasia Homolog) ELISA Kit	ACD	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0237	Chicken ARO(Aromatase) ELISA Kit	ARO	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0238	Chicken Cyt-C(Cytochrome C) ELISA Kit	Cyt-C	Sandwich	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-Ch0239	Chicken CT(Calcitonin) ELISA Kit	CT	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0240	Chicken sICAM-1/CD54(Soluble Intercellular Adhesion Molecule 1) ELISA Kit	sICAM-1/CD54	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0241	Chicken Pro-MMP-1(pro-Matrix Metalloproteinase 1) ELISA Kit	Pro-MMP-1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0242	Chicken ANGPTL6(Angiopoietin Like Protein 6) ELISA Kit	ANGPTL6	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0243	Chicken GAP(GTPase Activating Protein) ELISA Kit	GAP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0244	Chicken TIMP-1(Tissue Inhibitors of Metalloproteinase 1) ELISA Kit	TIMP-1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0245	Chicken LDHC(Lactate Dehydrogenase C) ELISA Kit	LDHC	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch0246	Chicken D4R/DRD4(Dopamine Receptor D4) ELISA Kit	D4R/DRD4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0247	Chicken Cyr61(Cysteine Rich Protein, Angiogenic Inducer 61) ELISA Kit	Cyr61	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0248	Chicken CD28(Cluster of Differentiation 28) ELISA Kit	CD28	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0249	Chicken PTHR2(Parathyroid Hormone Receptor 2) ELISA Kit	PTH2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0250	Chicken RAG-1(Recombination Activating Gene 1) ELISA Kit	RAG-1	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0251	Chicken DCX(Doublecortin) ELISA Kit	DCX	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0252	Chicken CLC(Charcot Leyden Crystal Protein) ELISA Kit	CLC	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0253	Chicken PSA(Prostate Specific Antigen) ELISA Kit	PSA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0254	Chicken CHGA(Chromogranin A) ELISA Kit	CHGA	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0255	Chicken ADAM9(A Disintegrin And Metalloprotease 9) ELISA Kit	ADAM9	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0256	Chicken ECM1(Extracellular Matrix Protein 1) ELISA Kit	ECM1	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0257	Chicken MIgD(Membrane Immunoglobulin D) ELISA Kit	MIgD	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0258	Chicken COL6(Collagen Type VI) ELISA Kit	COL6	Sandwich	25ng/mL	0.23 ng/mL	0.39-25ng/mL

FY-Ch0259	Chicken TNNT1(Troponin T Type 1, Slow Skeletal) ELISA Kit	TNNT1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0260	Chicken CAD(Caspase Activated Deoxyribonuclease) ELISA Kit	CAD	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0261	Chicken S100A8(S100 Calcium Binding Protein A8) ELISA Kit	S100A8	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0262	Chicken LRP6(Low Density Lipoprotein Receptor Related Protein) ELISA Kit	LRP6	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0263	Chicken VASPIN(Visceral Adipose Specific Serine Protease Inhibitor) ELISA Kit	VASPIN	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0264	Chicken PK(Pyruvate Kinase, Liver And RBC) ELISA Kit	PK	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0265	Chicken MMP-11(Matrix Metalloproteinase 11) ELISA Kit	MMP-11	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0266	Chicken CG(Chorionic Gonadotropin) ELISA Kit	CG	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0267	Chicken IGF2R(Insulin Like Growth Factor 2 Receptor) ELISA Kit	IGF2R	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0268	Chicken BMPR1A(Bone Morphogenetic Protein Receptor 1A) ELISA Kit	BMPR1A	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0269	Chicken GLβ(Galactosidase Beta) ELISA Kit	GLβ	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0270	Chicken PTPRO(Protein Tyrosine Phosphatase Receptor Type O) ELISA Kit	PTPRO	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0271	Chicken E-Cad(E-Cadherin) ELISA Kit	E-Cad	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0272	Chicken Tgtp(T-cell Specific Gtpase) ELISA Kit	Tgtp	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0273	Chicken C/EBPβ(CCAAT/Enhancer Binding Protein Beta) ELISA Kit	C/EBPβ	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0274	Chicken LAP(Leucine Aminopeptidase) ELISA Kit	LAP	Sandwich	15ng/mL	0.14 ng/mL	0.23-15ng/mL
FY-Ch0275	Chicken GSK3α(Glycogen Synthase Kinase 3 Alpha) ELISA Kit	GSK3α	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0276	Chicken GRN(Granulin) ELISA Kit	GRN	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0277	Chicken ADCY3(Adenylate Cyclase 3) ELISA Kit	ADCY3	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0278	Chicken NOX4(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 4) ELISA Kit	NOX4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0279	Chicken GSTα5(Glutathione S Transferase Alpha 5) ELISA Kit	GSTα5	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0280	Chicken IL-2Rβ(Interleukin 2 Receptor Beta) ELISA Kit	IL-2Rβ	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0281	Chicken IL-16(Interleukin 16) ELISA Kit	IL-16	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0282	Chicken IgM(Immunoglobulin M) ELISA Kit	IgM	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0283	Chicken LFA-2/CD2(Lymphocyte Function Associated Antigen 2) ELISA Kit	LFA-2/CD2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0284	Chicken CHGB(Chromogranin B) ELISA Kit	CHGB	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0285	Chicken PLAP/ALPP(Placental Alkaline Phosphatase) ELISA Kit	PLAP/ALPP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0287	Chicken M-AChR M2(Muscarinic Acetylcholine Receptor M2) ELISA Kit	M-AChR M2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0288	Chicken GATA1(GATA Binding Protein 1) ELISA Kit	GATA1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0289	Chicken TFR1(Transferrin Receptor 1) ELISA Kit	TFR1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0290	Chicken KNG1(Kininogen 1) ELISA Kit	KNG1	Sandwich	25pg/mL	0.23 pg/mL	0.39-25pg/mL
FY-Ch0291	Chicken UBR4(Ubiquitin Protein Ligase E3 Component N-Recognin 4) ELISA Kit	UBR4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0292	Chicken COL6α1(Collagen Type VI Alpha 1) ELISA Kit	COL6α1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0293	Chicken BANP(BTG3 Associated Nuclear Protein) ELISA Kit	BANP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0294	Chicken ALK(Anaplastic Lymphoma Kinase) ELISA Kit	ALK	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0295	Chicken Aβ1-40(Amyloid Beta 1-40) ELISA Kit	Aβ1-40	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0296	Chicken GLTP(Glycolipid Transfer Protein) ELISA Kit	GLTP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0297	Chicken GSTθ2/GSTt2(Glutathione S Transferase Theta 2) ELISA Kit	GSTθ2/GSTt2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0298	Chicken PK2(Prokineticin 2) ELISA Kit	PK2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0299	Chicken NTRK2(Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit	NTRK2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0300	Chicken FLNC(Filamin C, Gamma) ELISA Kit	FLNC	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL

FY-Ch0301	Chicken EXT1(Exostoses 1) ELISA Kit	EXT1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0302	Chicken VIPR1(Vasoactive Intestinal Peptide Receptor 1) ELISA Kit	VIPR1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0303	Chicken CD99(Cluster of Differentiation 99) ELISA Kit	CD99	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0304	Chicken ARRβ2(Arrestin Beta 2) ELISA Kit	ARRβ2	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0305	Chicken DPP10(Dipeptidyl Peptidase X) ELISA Kit	DPP10	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0306	Chicken Smad4/MADH4(Mothers Against Decapentaplegic Homolog 4) ELISA Kit	Smad4/MADH4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0307	Chicken CFH(Complement Factor H) ELISA Kit	CFH	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-Ch0308	Chicken TWEAK(Tumour Necrosis Factor Related Weak Inducer of Apoptosis) ELISA Kit	TWEAK	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0309	Chicken PLP1(Proteolipid Protein 1, Myelin) ELISA Kit	PLP1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0310	Chicken CaSR(Calcium Sensing Receptor) ELISA Kit	CaSR	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0311	Chicken ERN1(Endoplasmic Reticulum to Nucleus Signalling 1) ELISA Kit	ERN1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0312	Chicken TTR(Transthyretin) ELISA Kit	TTR	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch0313	Chicken M-AChR M3(Muscarinic Acetylcholine Receptor M3) ELISA Kit	M-AChR M3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0315	Chicken PRF1(Perforin 1) ELISA Kit	PRF1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0316	Chicken NUP155(Nucleoporin 155kDa) ELISA Kit	NUP155	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0317	Chicken COL8α1(Collagen Type VIII Alpha 1) ELISA Kit	COL8α1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0318	Chicken PCPE1(Procollagen C Proteinase Enhancer 1) ELISA Kit	PCPE1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0319	Chicken TP/TXA2R(Thromboxane Receptor) ELISA Kit	TP/TXA2R	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0320	Chicken MlgM(Membrane Immunoglobulin M) ELISA Kit	MlgM	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0321	Chicken TS/TYMS(Thymidylate Synthetase) ELISA Kit	TS/TYMS	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0322	Chicken ATM(Ataxia Telangiectasia Mutated) ELISA Kit	ATM	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0323	Chicken SAA2(Serum Amyloid A2) ELISA Kit	SAA2	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0324	Chicken SDF-2(Stromal Cell Derived Factor 2) ELISA Kit	SDF-2	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0325	Chicken γ-ECS(γ-Glutamyl Systeine Synthetase) ELISA Kit	γ-ECS	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0326	Chicken ME20M(Melanoma Associated ME20) ELISA Kit	ME20M	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0327	Chicken I-PTH(intact Parathormone) ELISA Kit	I-PTH	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-Ch0328	Chicken BMP-15(Bone Morphogenetic Protein 15) ELISA Kit	BMP-15	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0329	Chicken ADA(Adenosine Deaminase) ELISA Kit	ADA	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0330	Chicken FGFR3(Fibroblast Growth Factor Receptor 3) ELISA Kit	FGFR3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0331	Chicken CCL3L1(Chemokine C-C-Motif Ligand 3 Like Protein 1) ELISA Kit	CCL3L1	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0333	Chicken CD97(Cluster of Differentiation 97) ELISA Kit	CD97	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0334	Chicken CKB(Creatine Kinase, Brain) ELISA Kit	CKB	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0335	Chicken GPR37(G Protein Coupled Receptor 37) ELISA Kit	GPR37	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0336	Chicken IL-9(Interleukin 9) ELISA Kit	IL-9	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0337	Chicken AGC(Aggregan) ELISA Kit	AGC	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0338	Chicken C2(Complement Component 2) ELISA Kit	C2	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch0339	Chicken GUK1(Guanylate Kinase 1) ELISA Kit	GUK1	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0340	Chicken SIGLEC1(Sialic Acid Binding Immunoglobulin Like Lectin 1) ELISA Kit	SIGLEC1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0342	Chicken LTA4H(Leukotriene A4 Hydrolase) ELISA Kit	LTA4H	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0343	Chicken PRCP(Prolylcarboxypeptidase) ELISA Kit	PRCP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0344	Chicken Bcl2L11(Bcl-2 Like Protein 11) ELISA Kit	Bcl2L11	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL

FY-Ch0345	Chicken SAP25(Synaptosome Associated Protein 25) ELISA Kit	SAP25	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0346	Chicken NOX5(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 5) ELISA Kit	NOX5	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0347	Chicken ADIPOR1(Adiponectin Receptor 1) ELISA Kit	ADIPOR1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0348	Chicken GHRH(Growth Hormone Releasing Hormone) ELISA Kit	GHRH	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0349	Chicken MRC2(Mannose Receptor C Type 2) ELISA Kit	MRC2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0350	Chicken BMP-7(Bone Morphogenetic Protein 7) ELISA Kit	BMP-7	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0351	Chicken PRB(Retinoblastoma Tumor Suppressor Protein) ELISA Kit	PRB	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0352	Chicken SIGLEC3(Sialic Acid Binding Ig Like Lectin 3) ELISA Kit	SIGLEC3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0353	Chicken GAL8(Galectin 8) ELISA Kit	GAL8	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0355	Chicken MT1E(Metallothionein 1E) ELISA Kit	MT1E	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0356	Chicken VLDLR(Very Low Density Lipoprotein Receptor) ELISA Kit	VLDLR	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0357	Chicken MUC1(Mucin 1) ELISA Kit	MUC1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0358	Chicken CTSC(Cathepsin C) ELISA Kit	CTSC	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0359	Chicken VF(Visfatin) ELISA Kit	VF	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0360	Chicken PRG4(Proteoglycan 4) ELISA Kit	PRG4	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch0361	Chicken GDF1(Growth Differentiation Factor 1) ELISA Kit	GDF1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0362	Chicken MIP-1 β (Macrophage Inflammatory Protein 1 Beta) ELISA Kit	MIP-1 β	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0363	Chicken TPS(Tissue Polypeptide Specific Antigen) ELISA Kit	TPS	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0364	Chicken GDF15(Growth Differentiation Factor 15) ELISA Kit	GDF15	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0365	Chicken RARS(Arginyl tRNA Synthetase) ELISA Kit	RARS	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0366	Chicken TAF2(TATA Box Binding Protein Associated Factor 2) ELISA Kit	TAF2	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0367	Chicken NSE(Neuron-Specific Enolase) ELISA Kit	NSE	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0368	Chicken CCRL1(Chemokine C-C-Motif Receptor Like Protein 1) ELISA Kit	CCRL1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0369	Chicken COL3(Collagen Type III) ELISA Kit	COL3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0370	Chicken GHITM(Growth Hormone Inducible Transmembrane Protein) ELISA Kit	GHITM	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0371	Chicken GZMK(Granzyme K) ELISA Kit	GZMK	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0372	Chicken ETRB(Endothelin Receptor B) ELISA Kit	ETRB	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0373	Chicken GPC3(Glypican 3) ELISA Kit	GPC3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0374	Chicken NUP85(Nucleoporin 85kDa) ELISA Kit	NUP85	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0375	Chicken PLD2(Phospholipase D2) ELISA Kit	PLD2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0376	Chicken MMP-17(Matrix Metalloproteinase 17) ELISA Kit	MMP-17	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0377	Chicken ProGRP(Pro-Gastrin Releasing Peptide) ELISA Kit	ProGRP	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0378	Chicken EPOR(Erythropoietin Receptor) ELISA Kit	EPOR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0379	Chicken ALDM(Aldehyde Dehydrogenase, Mitochondrial) ELISA Kit	ALDM	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0380	Chicken NT-4(Neurotrophin 4) ELISA Kit	NT-4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0381	Chicken FGFR4(Fibroblast Growth Factor Receptor 4) ELISA Kit	FGFR4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0382	Chicken LDL(Low Density Lipoprotein) ELISA Kit	LDL	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0383	Chicken TLA(Thymus Leukemia Antigen) ELISA Kit	TLA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0384	Chicken β MSH(beta-Melanocyte Stimulating Hormone) ELISA Kit	β MSH	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0385	Chicken SEPP1(Selenoprotein P) ELISA Kit	SEPP1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0386	Chicken ACT α 1(Actin Alpha 1, Skeletal Muscle) ELISA Kit	ACT α 1	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL

FY-Ch0387	Chicken ANG(Angiogenin) ELISA Kit	ANG	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0388	Chicken NCAM/CD56(Neural Cell Adhesion Molecule) ELISA Kit	NCAM/CD56	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0389	Chicken GAL14(Galectin 14) ELISA Kit	GAL14	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0390	Chicken Collagenase I ELISA Kit	Collagenase I	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0391	Chicken UBE1L(Ubiquitin Activating Enzyme E1 Like Protein) ELISA Kit	UBE1L	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0392	Chicken UBE2C(Ubiquitin Conjugating Enzyme E2C) ELISA Kit	UBE2C	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0393	Chicken MMP-24(Matrix Metalloproteinase 24) ELISA Kit	MMP-24	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0394	Chicken MCSFR(Macrophage Colony Stimulating Factor Receptor) ELISA Kit	MCSFR	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0395	Chicken ACV-A(Activin A) ELISA Kit	ACV-A	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0396	Chicken EP2(Prostaglandin E Receptor 2) ELISA Kit	EP2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0397	Chicken IL-1β(Interleukin 1 Beta) ELISA Kit	IL-1β	Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-Ch0398	Chicken ANG II R2/AGTR2(Angiotensin II Receptor 2) ELISA Kit	ANG II R2/AGTR2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0399	Chicken TUBα3C(Tubulin Alpha 3c) ELISA Kit	TUBα3C	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0400	Chicken MSLN(Mesothelin) ELISA Kit	MSLN	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0401	Chicken CSF2Rβ(Colony Stimulating Factor 2 Receptor Beta) ELISA Kit	CSF2Rβ	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0402	Chicken NUP188(Nucleoporin 188kDa) ELISA Kit	NUP188	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0403	Chicken CKLFSF4(Chemokine Like Factor Superfamily 4) ELISA Kit	CKLFSF4	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0404	Chicken ANG2(Angiopoietin 2) ELISA Kit	ANG2	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0405	Chicken C4BPβ(C4 Binding Protein Beta) ELISA Kit	C4BPβ	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0406	Chicken AQP-3(Aquaporin 3) ELISA Kit	AQP-3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0407	Chicken COL11α1(Collagen Type XI Alpha 1) ELISA Kit	COL11α1	Sandwich	120pg/mL	1.13 pg/mL	1.88-120pg/mL
FY-Ch0408	Chicken FPB(Fibrinopeptide B) ELISA Kit	FPB	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0409	Chicken CTNNα1(Catenin, Alpha 1) ELISA Kit	CTNNα1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0410	Chicken CALD(Caldesmon) ELISA Kit	CALD	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0411	Chicken HBF(Fetal Hemoglobin) ELISA Kit	HBF	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0412	Chicken PCIII(Procollagen III) ELISA Kit	PCIII	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0413	Chicken FOXO4(Forkhead Box Protein O4) ELISA Kit	FOXO4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0414	Chicken COL6α3(Collagen Type VI Alpha 3) ELISA Kit	COL6α3	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0415	Chicken Pro-MMP-13(Pro-Matrix Metalloproteinase 13) ELISA Kit	Pro-MMP-13	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0416	Chicken CACYBP(Calcyclin Binding Protein) ELISA Kit	CACYBP	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0417	Chicken TYR(Tyrosinase) ELISA Kit	TYR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0418	Chicken CENPB(Centromere Protein B) ELISA Kit	CENPB	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0419	Chicken HSP-60/HSPD1(Heat Shock Protein 60) ELISA Kit	HSP-60/HSPD1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0420	Chicken DCN(Decorin) ELISA Kit	DCN	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0421	Chicken MYH10(Myosin Heavy Chain 10, Non Muscle) ELISA Kit	MYH10	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0422	Chicken MVP(Major Vault Protein) ELISA Kit	MVP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0423	Chicken MCP-2/CCL8(Monocyte Chemotactic Protein 2) ELISA Kit	MCP-2/CCL8	Sandwich	20pg/mL	0.19 pg/mL	0.31-20pg/mL
FY-Ch0424	Chicken BFP(Brain Finger Protein) ELISA Kit	BFP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0425	Chicken AFGF/FGF1(Acidic Fibroblast Growth Factor 1) ELISA Kit	AFGF/FGF1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0426	Chicken MRP1(Multidrug Resistance Associated Protein 1) ELISA Kit	MRP1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0427	Chicken BMP-10(Bone Morphogenetic Protein 10) ELISA Kit	BMP-10	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-Ch0428	Chicken PAI2(Plasminogen Activator Inhibitor 2) ELISA Kit	PAI2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0429	Chicken GRIK2(Glutamate Receptor, Ionotropic, Kainate 2) ELISA Kit	GRIK2	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0430	Chicken Tyk-2(Tyrosine Kinase 2) ELISA Kit	Tyk-2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0431	Chicken ApoM(Apolipoprotein M) ELISA Kit	ApoM	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0432	Chicken UBE2A(Ubiquitin Conjugating Enzyme E2A) ELISA Kit	UBE2A	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0433	Chicken EAR1(Eosinophil Associated Ribonuclease A family Member 1) ELISA Kit	EAR1	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0434	Chicken RPL10(Ribosomal Protein L10) ELISA Kit	RPL10	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0435	Chicken Preptin ELISA Kit	Preptin	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0436	Chicken Bcl-10(B-cell Lymphoma/Leukemia 10) ELISA Kit	Bcl-10	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0437	Chicken BMP-2(Bone Morphogenetic Protein 2) ELISA Kit	BMP-2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0438	Chicken RXFP1(Relaxin/Insulin Like Family Peptide Receptor 1) ELISA Kit	RXFP1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0439	Chicken PSPN(Persephin) ELISA Kit	PSPN	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0440	Chicken CX40(Connexin 40) ELISA Kit	CX40	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0441	Chicken NAGase(N-Acetyl Beta-D-Glucosaminidase) ELISA Kit	NAGase	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0442	Chicken CNTF(Ciliary Neurotrophic Factor) ELISA Kit	CNTF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0443	Chicken APC(Activated Protein C) ELISA Kit	APC	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0444	Chicken EL(EuglobulinLysis) ELISA Kit	EL	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0445	Chicken AMBP(Alpha-1-Microglobulin/Bikunin Precursor) ELISA Kit	AMBP	Sandwich	15ng/mL	0.14 ng/mL	0.23-15ng/mL
FY-Ch0446	Chicken CRY α B(Crystallin Alpha B) ELISA Kit	CRY α B	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0447	Chicken PSAP(Prosaposin) ELISA Kit	PSAP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0448	Chicken IL-5(Interleukin 5) ELISA Kit	IL-5	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0449	Chicken SELE(E-selectin) ELISA Kit	SELE	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0450	Chicken PADI4(Peptidyl Arginine Deiminase Type IV) ELISA Kit	PADI4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0451	Chicken IP-10/CXCL10(Interferon Gamma Induced Protein 10kDa) ELISA Kit	IP-10/CXCL10	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0452	Chicken FADS2(Fatty Acid Desaturase 2) ELISA Kit	FADS2	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0453	Chicken IL-20(Interleukin 20) ELISA Kit	IL-20	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0454	Chicken GLUT2(Glucose Transporter 2) ELISA Kit	GLUT2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0455	Chicken PKI β (Protein Kinase Inhibitor Beta) ELISA Kit	PKI β	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0456	Chicken SUMO2(Small Ubiquitin Related Modifier 2) ELISA Kit	SUMO2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0458	Chicken IgA(Immunoglobulin A) ELISA Kit	IgA	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0459	Chicken PF4(Platelet Factor 4) ELISA Kit	PF4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0460	Chicken ARR β 1(Arrestin Beta 1) ELISA Kit	ARR β 1	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0461	Chicken MOG(Myelin Oligodendrocyte Glycoprotein) ELISA Kit	MOG	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0462	Chicken FETUB(Fetuin B) ELISA Kit	FETUB	Sandwich	300pg/mL	2.81 pg/mL	4.69-300pg/mL
FY-Ch0463	Chicken GTF2A1(General Transcription Factor II A, Polypeptide 1) ELISA Kit	GTF2A1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0464	Chicken PGRP-S(Peptidoglycan Recognition Protein-S) ELISA Kit	PGRP-S	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0465	Chicken MBP-C/MBL2(Mannose Binding Protein C/Mannose Binding Lectin) ELISA Kit	MBP-C/MBL2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0466	Chicken PDGFB(Platelet Derived Growth Factor Subunit B) ELISA Kit	PDGFB	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-Ch0467	Chicken RIPK2(Receptor TNFRSF Interacting Serine Threonine Kinase 2) ELISA Kit	RIPK2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0468	Chicken GAL2(Galectin 2) ELISA Kit	GAL2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0469	Chicken C-jun(Proto-oncogene c-Jun) ELISA Kit	C-jun	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL

FY-Ch0470	Chicken MAP2K6/MKK6(Mitogen Activated Protein Kinase Kinase 6) ELISA Kit	MAP2K6/MKK6	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0471	Chicken PIIIICP(Procollagen III C-Terminal ProPeptide) ELISA Kit	PIIIICP	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-Ch0472	Chicken REG4(Regenerating Islet Derived Protein 4) ELISA Kit	REG4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0474	Chicken C4BP α (C4 Binding Protein Alpha) ELISA Kit	C4BP α	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0475	Chicken OPN(Osteopontin) ELISA Kit	OPN	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0476	Chicken cTn-I/TNNI3(cardiac Troponin I) ELISA Kit	cTn-I/TNNI3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0477	Chicken MPIF-1/CCL23(Myeloid Progenitor Inhibitory Factor 1) ELISA Kit	MPIF-1/CCL23	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0478	Chicken SRP9(Signal Recognition Particle 9kDa) ELISA Kit	SRP9	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0479	Chicken SMAF(Specific Macrophage Arming Factor) ELISA Kit	SMAF	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0480	Chicken IGFBP-1(Insulin Like Growth Factor Binding Protein 1) ELISA Kit	IGFBP-1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0481	Chicken GSTs(Glutathione S-Transferases) ELISA Kit	GSTs	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0482	Chicken Fc ϵ R 1(Receptor 1 for the Fc Fragment of IgE) ELISA Kit	Fc ϵ R 1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0483	Chicken GAL9(Galectin 9) ELISA Kit	GAL9	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0484	Chicken ApoA4(Apolipoprotein A4) ELISA Kit	ApoA4	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0485	Chicken MAP2K1(Mitogen Activated Protein Kinase Kinase 1) ELISA Kit	MAP2K1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0486	Chicken HHIP(Hedgehog Interacting Protein) ELISA Kit	HHIP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0487	Chicken PON3(Paraoxonase 3) ELISA Kit	PON3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0488	Chicken NT-proANP(N-Terminal Pro Atrial Natriuretic Peptide) ELISA Kit	NT-proANP	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0489	Chicken SNC α (Synuclein Alpha) ELISA Kit	SNC α	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0490	Chicken CD59(Protectin) ELISA Kit	CD59	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0491	Chicken MIP-1 α (Macrophage Inflammatory Protein 1 Alpha) ELISA Kit	MIP-1 α	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0492	Chicken GTEAGE(Glomerular Tissue Extract Advanced Glycosylation End Products)	GTEAGE	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0493	Chicken PRSS2(Protease, Serine, 2) ELISA Kit	PRSS2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0494	Chicken FBN2(Fibrillin 2) ELISA Kit	FBN2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0495	Chicken Fc γ R3A(Fc Fragment of IgG Low Affinity IIIa Receptor) ELISA Kit	Fc γ R3A	Sandwich	25ng/mL	0.23 ng/mL	0.39-25ng/mL
FY-Ch0496	Chicken TNK2(Tyrosine Kinase, Non Receptor 2) ELISA Kit	TNK2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0497	Chicken Des(Desmin) ELISA Kit	Des	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0498	Chicken PCSK9(Proprotein Convertase Subtilisin/Kexin Type 9) ELISA Kit	PCSK9	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0499	Chicken CPE(Carboxypeptidase E) ELISA Kit	CPE	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0500	Chicken MCAM(Melanoma Cell Adhesion Molecule) ELISA Kit	MCAM	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0501	Chicken OPGL(Osteoprotegerin Ligand) ELISA Kit	OPGL	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0502	Chicken COL10(Collagen Type X) ELISA Kit	COL10	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0503	Chicken PGA(Pepsinogen A) ELISA Kit	PGA	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0504	Chicken F9(Coagulation Factor IX) ELISA Kit	F9	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-Ch0505	Chicken SCG2(Secretogranin II) ELISA Kit	SCG2	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-Ch0506	Chicken BPI(Bactericidal/Permeability Increasing Protein) ELISA Kit	BPI	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0507	Chicken RAR α (Retinoic Acid Receptor Alpha) ELISA Kit	RAR α	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0508	Chicken C3b(Complement Fragment 3b) ELISA Kit	C3b	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0509	Chicken iPLA2(Phospholipidase A2, Calcium Independent) ELISA Kit	iPLA2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0510	Chicken NMP-22(Nuclear Matrix Protein 22) ELISA Kit	NMP-22	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0511	Chicken CTGF(Connective Tissue Growth Factor) ELISA Kit	CTGF	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-Ch0512	Chicken THRSP(Thyroid Hormone Responsive) ELISA Kit	THRSP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0513	Chicken FABP3(Fatty Acid Binding Protein 3, Muscle and Heart) ELISA Kit	FABP3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0515	Chicken HSPβ6/HSP-20(Heat Shock Protein Beta 6) ELISA Kit	HSPβ6/HSP-20	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0516	Chicken TCF20(Transcription Factor 20) ELISA Kit	TCF20	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0517	Chicken CNX(Calnexin) ELISA Kit	CNX	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0518	Chicken PRB2(Basic Salivary Proline Rich Protein 2) ELISA Kit	PRB2	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0519	Chicken CCNB(Cyclin-B) ELISA Kit	CCNB	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0520	Chicken FBN1(Fibrillin 1) ELISA Kit	FBN1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0521	Chicken KLK6(Kallikrein 6) ELISA Kit	KLK6	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0522	Chicken APAF1(Apoptosis Protease Activating Factor 1) ELISA Kit	APAF1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0523	Chicken RB1(Retinoblastoma Protein 1) ELISA Kit	RB1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0524	Chicken H-Cadherin(Cadherin, Heart) ELISA Kit	H-Cadherin	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0525	Chicken CYP26A1(Cytochrome P450, family 26, subfamily A, polypeptide 1) ELISA Kit	CYP26A1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0526	Chicken FLOT2(Flotillin 2) ELISA Kit	FLOT2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0527	Chicken CYP24A1(Cytochrome P450, family 24, subfamily A, polypeptide 1) ELISA Kit	CYP24A1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0528	Chicken F11(Coagulation Factor XI) ELISA Kit	F11	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0529	Chicken MMP-2(Matrix Metalloproteinase 2) ELISA Kit	MMP-2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0530	Chicken FG(Fibrinogen) ELISA Kit	FG	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-Ch0531	Chicken HMGCS1(3-Hydroxy-3-Methylglutaryl-Coenzyme A Synthase 1) ELISA Kit	HMGCS1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0532	Chicken CASP7(Caspase 7) ELISA Kit	CASP7	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0533	Chicken Sp17(Sperm Protein 17) ELISA Kit	Sp17	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0534	Chicken CFB(Complement Factor B) ELISA Kit	CFB	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch0535	Chicken LDHA(Lactate Dehydrogenase A) ELISA Kit	LDHA	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0536	Chicken αHSP(Alpha-Hemoglobin Stabilizing Protein) ELISA Kit	αHSP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0537	Chicken DAP6(Death Associated Protein 6) ELISA Kit	DAP6	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0538	Chicken DSE(Dermatan Sulfate Epimerase) ELISA Kit	DSE	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0540	Chicken DT(Dopachrome Tautomerase) ELISA Kit	DT	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0541	Chicken SOX1(Sex Determining Region Y Box Protein 1) ELISA Kit	SOX1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0542	Chicken MK(Midkine) ELISA Kit	MK	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0543	Chicken sTNF-1(Soluble Tumor Necrosis Factor 1) ELISA Kit	sTNF-1	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0544	Chicken MAPKAPK3(Mitogen Activated Protein Kinase Activated Protein Kinase 3) ELISA Kit	MAPKAPK3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0545	Chicken IκBβ(Inhibitory Subunit of NF Kappa B Beta) ELISA Kit	IκBβ	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0546	Chicken TFPI2(Tissue Factor Pathway Inhibitor 2) ELISA Kit	TFPI2	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0547	Chicken PAFAHIB3(Platelet Activating Factor Acetylhydrolase Ib3) ELISA Kit	PAFAHIB3	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0548	Chicken CASP4(Caspase 4) ELISA Kit	CASP4	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch0549	Chicken SAG(S Antigen) ELISA Kit	SAG	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0550	Chicken Lptn/XCL1(Lymphotactin) ELISA Kit	Lptn/XCL1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0551	Chicken RI(Ribonuclease Inhibitor) ELISA Kit	RI	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-Ch0552	Chicken PKCθ(Protein Kinase C Theta) ELISA Kit	PKCθ	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0553	Chicken HBμ(Hemoglobin, Mu) ELISA Kit	HBμ	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0554	Chicken IL-6R(Interleukin 6 Receptor) ELISA Kit	IL-6R	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL

FY-Ch0555	Chicken IL-11(Interleukin 11) ELISA Kit	IL-11	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0556	Chicken SP-A1(Pulmonary surfactant-associated protein A1) ELISA Kit	SP-A1	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0557	Chicken PLAA(Phospholipase A2 Activating Protein) ELISA Kit	PLAA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0558	Chicken ARSB(Arylsulfatase B) ELISA Kit	ARSB	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0559	Chicken Try2(Trypsinogen II) ELISA Kit	Try2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0560	Chicken CETP(Cholesteryl Ester Transfer Protein) ELISA Kit	CETP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0561	Chicken RANTES(Regulated On Activation, Normal T-Cell Expressed and Secreted) ELISA Kit	RANTES	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0562	Chicken AIG1(Androgen Induced Protein 1) ELISA Kit	AIG1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0563	Chicken LIPF(Lipase, Gastric) ELISA Kit	LIPF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0564	Chicken GAL7(Galactin 7) ELISA Kit	GAL7	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0565	Chicken TDGF1(Teratocarcinoma Derived Growth Factor 1) ELISA Kit	TDGF1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0566	Chicken PKCζ(Protein Kinase C Zeta) ELISA Kit	PKCζ	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0567	Chicken ERK2(Extracellular Signal Regulated Kinase 2) ELISA Kit	ERK2	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0568	Chicken NUP205(Nucleoporin 205kDa) ELISA Kit	NUP205	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0569	Chicken TNPO2(Transportin 2) ELISA Kit	TNPO2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0570	Chicken ATXN1(Ataxin 1) ELISA Kit	ATXN1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0571	Chicken CA2(Carbonic Anhydrase II) ELISA Kit	CA2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0572	Chicken PGD2S(Prostaglandin D2 Synthase, Brain) ELISA Kit	PGD2S	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch0573	Chicken αFP(Alpha-Fetoprotein) ELISA Kit	αFP	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-Ch0574	Chicken MRP(Multidrug Resistance-Associated Protein) ELISA Kit	MRP	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0575	Chicken AGRN(Agrin) ELISA Kit	AGRN	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0576	Chicken SCTR(Secretin Receptor) ELISA Kit	SCTR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0577	Chicken TACE/ADAM17(TNF α Converting Enzyme) ELISA Kit	TACE/ADAM17	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0578	Chicken CLU(Clusterin) ELISA Kit	CLU	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0579	Chicken TREM-1(Triggering Receptor Expressed on Myeloid Cells-1) ELISA Kit	TREM-1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0580	Chicken ApoA5(Apolipoprotein A5) ELISA Kit	ApoA5	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0581	Chicken MYOG(Myogenin) ELISA Kit	MYOG	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0582	Chicken CKLFSF6(Chemokine Like Factor Superfamily 6) ELISA Kit	CKLFSF6	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0583	Chicken GIF(Gastric Intrinsic Factor) ELISA Kit	GIF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0584	Chicken BMG/β2-MG(Beta-2-Microglobulin) ELISA Kit	BMG/β2-MG	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0585	Chicken Bcl-3(B-cell Lymphoma/Leukemia 3) ELISA Kit	Bcl-3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0586	Chicken MUC6(Mucin 6) ELISA Kit	MUC6	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0587	Chicken IFABP/FABP2(Intestinal Fatty Acid Binding Protein) ELISA Kit	IFABP/FABP2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0588	Chicken IL-1ra/IL-1F3(Interleukin 1 Receptor Antagonist) ELISA Kit	IL-1ra/IL-1F3	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0589	Chicken ENA-78(Epithelial Neutrophil Activating Peptide 78) ELISA Kit	ENA-78	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0590	Chicken PRELP(Proline Arginine Rich End Leucine Rich Repeat Protein) ELISA Kit	PRELP	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0591	Chicken FcεR II/CD23(Receptor II for the Fc Region of Immunoglobulin E) ELISA Kit	FcεR II/CD23	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0592	Chicken MXD1(MAX Dimerization Protein 1) ELISA Kit	MXD1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0593	Chicken UBAP2(Ubiquitin Associated Protein 2) ELISA Kit	UBAP2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0594	Chicken GDF6(Growth Differentiation Factor 6) ELISA Kit	GDF6	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0595	Chicken TBG(Thyroxine Binding Globulin) ELISA Kit	TBG	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL

FY-Ch0596	Chicken ALDOA(Aldolase A, Fructose Bisphosphate) ELISA Kit	ALDOA	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0597	Chicken SGK2(Serum/Glucocorticoid Regulated Kinase 2) ELISA Kit	SGK2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0598	Chicken SOX3(Sex Determining Region Y Box Protein 3) ELISA Kit	SOX3	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0599	Chicken GTF2B(General Transcription Factor II B) ELISA Kit	GTF2B	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0600	Chicken cPLA2(Cytosolic Phospholipase A2) ELISA Kit	cPLA2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0601	Chicken ALDOC(Aldolase C, Fructose Bisphosphate) ELISA Kit	ALDOC	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0602	Chicken FGF9(Fibroblast Growth Factor 9) ELISA Kit	FGF9	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0603	Chicken USP8(Ubiquitin Specific Peptidase 8) ELISA Kit	USP8	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0604	Chicken sVCAM-1/CD106(soluble Vascular Cell Adhesion Molecule 1) ELISA Kit	sVCAM-1/CD106	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0605	Chicken MYH1(Myosin Heavy Chain 1, Skeletal Muscle, Adult) ELISA Kit	MYH1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0606	Chicken PC/CPG(Choline Phosphoglyceride) ELISA Kit	PC/CPG	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0607	Chicken Ntn4(Netrin 4) ELISA Kit	Ntn4	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0608	Chicken FABP4(Fatty Acid Binding Protein 4, Adipocyte) ELISA Kit	FABP4	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0609	Chicken OCT2A(Octamer-binding protein 2A) ELISA Kit	OCT2A	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0610	Chicken pIκBα(Phospho-Inhibitory Subunit Of NF Kappa B Alpha) ELISA Kit	pIκBα	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0611	Chicken CX43(Connexin 43) ELISA Kit	CX43	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0612	Chicken IL-23(Interleukin 23) ELISA Kit	IL-23	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0613	Chicken SYCN(Syncollin) ELISA Kit	SYCN	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0614	Chicken GzmB(Granzyme B) ELISA Kit	GzmB	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0615	Chicken LEP(Leptin) ELISA Kit	LEP	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0616	Chicken FUCA(Alpha-L-Fucosidase, Tissue) ELISA Kit	FUCA	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0617	Chicken IL-18(Interleukin 18) ELISA Kit	IL-18	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0618	Chicken SSRP1(Structure Specific Recognition Protein 1) ELISA Kit	SSRP1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0619	Chicken DKK2(Dickkopf Related Protein 2) ELISA Kit	DKK2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0620	Chicken VIM(Vimentin) ELISA Kit	VIM	Sandwich	25ng/mL	0.23 ng/mL	0.39-25ng/mL
FY-Ch0621	Chicken COL1α1(Collagen Type I Alpha 1) ELISA Kit	COL1α1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0622	Chicken IL-22(Interleukin 22) ELISA Kit	IL-22	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0623	Chicken SOX9(Sex Determining Region Y Box Protein 9) ELISA Kit	SOX9	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0624	Chicken LZM(Lysozyme) ELISA Kit	LZM	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0625	Chicken DKK4(Dickkopf Related Protein 4) ELISA Kit	DKK4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0626	Chicken TLR9(Toll-Like Receptor 9) ELISA Kit	TLR9	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0627	Chicken TPO(Thrombopoietin) ELISA Kit	TPO	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0628	Chicken SAP(Serum Amyloid P Component) ELISA Kit	SAP	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0629	Chicken CAP1(Adenylate Cyclase Associated Protein 1) ELISA Kit	CAP1	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0630	Chicken IGFBP-3(Insulin Like Growth Factor Binding Protein 3) ELISA Kit	IGFBP-3	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch0631	Chicken CX3CL1(Cheokine C-X3-C-Motif Ligand 1) ELISA Kit	CX3CL1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0632	Chicken sTfR1(Soluble Transferrin Receptor1) ELISA Kit	sTfR1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0633	Chicken DBP(DNA Binding Protein) ELISA Kit	DBP	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-Ch0634	Chicken CASP 9(Caspase 9) ELISA Kit	CASP 9	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0635	Chicken KLK8(Kallikrein 8) ELISA Kit	KLK8	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0636	Chicken PRSS33(Protease, Serine 33) ELISA Kit	PRSS33	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL

FY-Ch0637	Chicken ICAM4(Intercellular Adhesion Molecule 4) ELISA Kit	ICAM4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0638	Chicken PL(Pancreatic Lipase) ELISA Kit	PL	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0639	Chicken CYP21B(Cytochrome P450 21B) ELISA Kit	CYP21B	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0640	Chicken FcαR/CD89(Receptor I for the Fc Fragment of IgA) ELISA Kit	FcαR/CD89	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0641	Chicken GAL1(Galectin 1) ELISA Kit	GAL1	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch0642	Chicken NT-3(Neurotrophin-3) ELISA Kit	NT-3	Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-Ch0643	Chicken VEGFR1/FLT1(Vasuoar Endothelial Growth Factor Receptor 1) ELISA Kit	VEGFR1/FLT1	Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-Ch0644	Chicken β-Glu(Beta Glucosidase) ELISA Kit	β-Glu	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0645	Chicken AFP-L3(Alpha-Fetoprotein Lens Culinaris Agglutinin 3) ELISA Kit	AFP-L3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0646	Chicken TE(Telomerase) ELISA Kit	TE	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0647	Chicken AMY2(Amylase Alpha 2, Pancreatic) ELISA Kit	AMY2	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0648	Chicken MSR1(Macrophage Scavenger Receptor 1) ELISA Kit	MSR1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0649	Chicken PAR4(Protease Activated Receptor 4) ELISA Kit	PAR4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0650	Chicken PKD2(Protein Kinase D2) ELISA Kit	PKD2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0651	Chicken PHB(Prohibitin) ELISA Kit	PHB	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0652	Chicken S100A7(S100 Calcium Binding Protein A7) ELISA Kit	S100A7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0653	Chicken ACO1(Aconitase 1) ELISA Kit	ACO1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0654	Chicken BALP(Bone Alkaline Phosphatase) ELISA Kit	BALP	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0655	Chicken FLRT2(Fibronectin Leucine Rich Transmembrane Protein 2) ELISA Kit	FLRT2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0656	Chicken ApoF(Apolipoprotein F) ELISA Kit	ApoF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0657	Chicken TPO(Thyroid Peroxidase) ELISA Kit	TPO	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0658	Chicken FcγR II/CD32(Receptor II for the Fc Fragment of IgG) ELISA Kit	FcγR II/CD32	Sandwich	150pg/mL	1.4 pg/mL	2.34-150pg/mL
FY-Ch0660	Chicken sFMC(Soluble Fibrin Monomer Complex) ELISA Kit	sFMC	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0661	Chicken SOCS3(Suppressors of Cytokine Signaling 3) ELISA Kit	SOCS3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0662	Chicken EphA2(Tyrosine protein kinase receptor A2) ELISA Kit	EphA2	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0663	Chicken F7(Coagulation Factor VII) ELISA Kit	F7	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0664	Chicken MX1(Myxovirus Resistance 1) ELISA Kit	MX1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0665	Chicken ApoE(Apolipoprotein E) ELISA Kit	ApoE	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch0666	Chicken CDSN(Corneodesmosin) ELISA Kit	CDSN	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0667	Chicken ITIH4(Inter-alpha-trypsin Inhibitor Heavy Chain H4) ELISA Kit	ITIH4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0668	Chicken GLUT1(Glucose Transporter 1) ELISA Kit	GLUT1	Sandwich	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-Ch0669	Chicken DNase-I(Deoxyribonuclease I) ELISA Kit	DNase-I	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0670	Chicken FRS2(Fibroblast Growth Factor Receptor Substrate 2) ELISA Kit	FRS2	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0671	Chicken MUC17(Mucin 17) ELISA Kit	MUC17	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0672	Chicken EGR2(Early Growth Response Protein 2) ELISA Kit	EGR2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0673	Chicken Slit2(Slit Homolog 2) ELISA Kit	Slit2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0675	Chicken P27(P27 Protein) ELISA Kit	P27	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0676	Chicken GPX1(Glutathione Peroxidase 1) ELISA Kit	GPX1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0677	Chicken OSMR(Oncostatin M Receptor) ELISA Kit	OSMR	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0678	Chicken PDIA3(Protein Disulfide Isomerase A3) ELISA Kit	PDIA3	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch0679	Chicken SCFR(Stem Cell Growth Factor Receptor) ELISA Kit	SCFR	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL

FY-Ch0680	Chicken G-CSF(Granulocyte Colony Stimulating Factor 3) ELISA Kit	G-CSF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0682	Chicken UCRP(Ubiquitin Cross Reactive Protein) ELISA Kit	UCRP	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0683	Chicken GR α (Glucocorticoid Receptor Alpha) ELISA Kit	GR α	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0684	Chicken C8b(Complement Component 8b) ELISA Kit	C8b	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0685	Chicken AR(Androgen Receptor) ELISA Kit	AR	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0686	Chicken RAN κ L(Receptor Activator of Nuclear Factor Kappa B Ligand) ELISA Kit	RAN κ L	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0687	Chicken TPP(Thrombus Precursor Protein) ELISA Kit	TPP	Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-Ch0688	Chicken CASP12(Caspase 12) ELISA Kit	CASP12	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0689	Chicken OCT1(Octamer-binding protein 1) ELISA Kit	OCT1	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0690	Chicken TSPAN30/CD63(Tetraspanin 30) ELISA Kit	TSPAN30/CD63	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0691	Chicken ANKRD1(Ankyrin Repeat Domain 1) ELISA Kit	ANKRD1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0692	Chicken CTR(Calcitonin Receptor) ELISA Kit	CTR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0693	Chicken ETRA(Endothelin Receptor A) ELISA Kit	ETRA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0694	Chicken MIgA(Membrane Immunoglobulin A) ELISA Kit	MIgA	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0695	Chicken CYP1B1(Cytochrome P450, family 1, subfamily B, polypeptide 1) ELISA Kit	CYP1B1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0696	Chicken ACE2(Angiotensin I Converting Enzyme 2) ELISA Kit	ACE2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0697	Chicken RNASEH1(Ribonuclease H1) ELISA Kit	RNASEH1	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0698	Chicken DNM1(Dynamin 1) ELISA Kit	DNM1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0699	Chicken PP(Pancreatic Polypeptide) ELISA Kit	PP	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0700	Chicken VLA4(Very Late Appearing Antigen 4) ELISA Kit	VLA4	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0701	Chicken PPAR- δ (Peroxisome Proliferator Activated Receptor Delta) ELISA Kit	PPAR- δ	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0702	Chicken PTN(Pleiotrophin) ELISA Kit	PTN	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0703	Chicken OPRM1(Mu-type opioid receptor) ELISA Kit	OPRM1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0704	Chicken RBP4(Retinol Binding Protein 4) ELISA Kit	RBP4	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0705	Chicken EFNA5(Ephrin A5) ELISA Kit	EFNA5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0706	Chicken CTSL(Cathepsin L) ELISA Kit	CTSL	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch0708	Chicken DNM2(Dynamin 2) ELISA Kit	DNM2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0709	Chicken BAD(Bcl2 Associated Death Promoter) ELISA Kit	BAD	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0710	Chicken GLUT3(Glucose Transporter 3) ELISA Kit	GLUT3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0711	Chicken MCP/CD46(Membrane Cofactor Protein) ELISA Kit	MCP/CD46	Sandwich	25pg/mL	0.23 pg/mL	0.39-25pg/mL
FY-Ch0712	Chicken HNRPA2B1(Heterogeneous Nuclear Ribonucleoprotein A2/B1) ELISA Kit	HNRPA2B1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0713	Chicken CD30(Cluster of Differentiation 30) ELISA Kit	CD30	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0714	Chicken P-GP(Permeability Glycoprotein) ELISA Kit	P-GP	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0715	Chicken NFR κ B(Nuclear Factor Related To Kappa B Binding Protein) ELISA Kit	NFR κ B	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0716	Chicken S100B(S100 Calcium Binding Protein B) ELISA Kit	S100B	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0717	Chicken HLTf(Helicase Like Transcription Factor) ELISA Kit	HLTF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0718	Chicken ATA(Angiotensinase A) ELISA Kit	ATA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0719	Chicken PGM(Phosphoglucomutase) ELISA Kit	PGM	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0720	Chicken EloA(Elongin A) ELISA Kit	EloA	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0721	Chicken PTGS2/COX-2(Prostaglandin Endoperoxide Synthase 2) ELISA Kit	PTGS2/COX-2	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-Ch0722	Chicken CKLFSF5(Chemokine Like Factor Superfamily 5) ELISA Kit	CKLFSF5	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL

FY-Ch0723	Chicken CD109(Cluster of Differentiation 109) ELISA Kit	CD109	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0724	Chicken GP39(Glycoprotein 39, Cartilage) ELISA Kit	GP39	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0725	Chicken TSP-2(Thrombospondin-2) ELISA Kit	TSP-2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0726	Chicken AZU1(Azurocidin 1) ELISA Kit	AZU1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0727	Chicken NFIL3(Nuclear Factor, Interleukin 3 Regulated) ELISA Kit	NFIL3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0728	Chicken PGC(Pepsinogen C) ELISA Kit	PGC	Sandwich	3000ng/mL	28.13 ng/mL	46.88-3000ng/mL
FY-Ch0729	Chicken SOX2(Sex Determining Region Y Box Protein 2) ELISA Kit	SOX2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0730	Chicken GST01/GSTt1(Glutathione S Transferase Theta 1) ELISA Kit	GST01/GSTt1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0731	Chicken CFI(Complement Factor I) ELISA Kit	CFI	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0732	Chicken NPHN(Nephrin) ELISA Kit	NPHN	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0733	Chicken C4a(Complement Component 4a) ELISA Kit	C4a	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0734	Chicken SSFA2(Sperm Specific Antigen 2) ELISA Kit	SSFA2	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0735	Chicken EGR1(Early Growth Response Protein 1) ELISA Kit	EGR1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0736	Chicken F5(Coagulation Factor V) ELISA Kit	F5	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0737	Chicken CLDN11(Claudin 11) ELISA Kit	CLDN11	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0738	Chicken IκBα(Inhibitory Subunit of NF Kappa B Alpha) ELISA Kit	IκBα	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0739	Chicken HTRA1(Htra Serine, Peptidase 1) ELISA Kit	HTRA1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0741	Chicken ApoC3(Apolipoprotein C3) ELISA Kit	ApoC3	Sandwich	4000ng/ml	37.5 ng/ml	62.5-4000ng/ml
FY-Ch0742	Chicken TJPI(Tight Junction Protein ZO-1) ELISA Kit	TJPI	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0743	Chicken RAG-2(Recombination Activating Gene 2) ELISA Kit	RAG-2	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0744	Chicken RANκ(Receptor Activator of Nuclear Factor Kappa B) ELISA Kit	RANκ	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0745	Chicken CSTA(Cystatin A) ELISA Kit	CSTA	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0746	Chicken Bax(Bcl-2 Associated X Protein) ELISA Kit	Bax	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0747	Chicken EGFR(Epidermal Growth Factor Receptor) ELISA Kit	EGFR	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0748	Chicken PAR3(Protease Activated Receptor 3) ELISA Kit	PAR3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0749	Chicken LpPLA2(Phospholipase A2, Lipoprotein Associated) ELISA Kit	LpPLA2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0750	Chicken β2GP1-IgG(β2 Glycoprotein 1 IgG) ELISA Kit	β2GP1-IgG	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0751	Chicken CHAc(Choline acetylase) ELISA Kit	CHAc	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0752	Chicken CTX I (Cross Linked C-telopeptide of Type I Collagen) ELISA Kit	CTX I	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0753	Chicken NUP62(Nucleoporin 62kDa) ELISA Kit	NUP62	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0754	Chicken ADH1(Alcohol Dehydrogenase 1) ELISA Kit	ADH1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0755	Chicken RPS6Kβ1(Ribosomal Protein S6 Kinase Beta 1) ELISA Kit	RPS6Kβ1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0756	Chicken CRN(Corin) ELISA Kit	CRN	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0757	Chicken MMP-3(Matrix Metalloproteinase 3) ELISA Kit	MMP-3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0758	Chicken EFNA3(Ephrin A3) ELISA Kit	EFNA3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0759	Chicken BLNK(B-Cell Linker Protein) ELISA Kit	BLNK	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0760	Chicken PCI(Protein C Inhibitor) ELISA Kit	PCI	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0762	Chicken GPC4(Glypican 4) ELISA Kit	GPC4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0763	Chicken UBD(Ubiquitin D) ELISA Kit	UBD	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0764	Chicken MMP-13(Matrix Metalloproteinase 13) ELISA Kit	MMP-13	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0765	Chicken GAL4(Galectin 4) ELISA Kit	GAL4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL

FY-Ch0766	Chicken PLA2R1(Phospholipase A2 Receptor 1) ELISA Kit	PLA2R1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0767	Chicken GROβ/CXCL2(Growth Regulated Oncogene Beta) ELISA Kit	GROβ/CXCL2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0768	Chicken CNTN4(Contactin 4) ELISA Kit	CNTN4	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0770	Chicken PLA2G2A(Phospholipase A2, Group II A) ELISA Kit	PLA2G2A	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0771	Chicken S100A12(S100 Calcium Binding Protein A12) ELISA Kit	S100A12	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0772	Chicken GRβ(Glucocorticoid Receptor Beta) ELISA Kit	GRβ	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0773	Chicken Smad3/MADH3(Mothers Against Decapentaplegic Homolog 3) ELISA Kit	Smad3/MADH3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0774	Chicken NE/ELA2(Neutrophil Elastase/Elastase-2) ELISA Kit	NE/ELA2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0776	Chicken CYPD(Cyclophilin D) ELISA Kit	CYPD	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0777	Chicken PLA2G2D(Phospholipase A2, Group II D) ELISA Kit	PLA2G2D	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0778	Chicken SPB(Pulmonary Surfactant Associated Protein B) ELISA Kit	SPB	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0779	Chicken TNFRSF1A(Tumor Necrosis Factor Receptor Superfamily, Member 1A) ELISA Kit	TNFRSF1A	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0780	Chicken DPYSL2(Dihydropyrimidinase Like Protein 2) ELISA Kit	DPYSL2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0781	Chicken CASP1(Caspase 1) ELISA Kit	CASP1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0782	Chicken ACP5(Tartrate Resistant Acid Phosphatase 5) ELISA Kit	ACP5	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0783	Chicken MTHFR(Methylenetetrahydrofolate Reductase) ELISA Kit	MTHFR	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0784	Chicken CNTN3(Contactin 3) ELISA Kit	CNTN3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0785	Chicken MUC2(Mucin 2) ELISA Kit	MUC2	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0786	Chicken Ras-GAP1(Rac-GTPase Activating Protein 1) ELISA Kit	Ras-GAP1	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0787	Chicken COL15(Collagen Type X V) ELISA Kit	COL15	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0788	Chicken C3a(Complement Component 3a) ELISA Kit	C3a	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch0789	Chicken TGF-α(Transforming Growth Factor α) ELISA Kit	TGF-α	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0790	Chicken P II CP(Procollagen II C-Terminal ProPeptide) ELISA Kit	P II CP	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0791	Chicken DPP6(Dipeptidyl Peptidase VI) ELISA Kit	DPP6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0792	Chicken GATA5(GATA Binding Protein 5) ELISA Kit	GATA5	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0793	Chicken GYPE(Glycophorin E) ELISA Kit	GYPE	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0794	Chicken PSG1/SP1(Pregnancy Specific Beta-1-Glycoprotein 1) ELISA Kit	PSG1/SP1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0795	Chicken MYL1(Myosin Light Chain 1, Alkali, Fast Skeletal) ELISA Kit	MYL1	Sandwich	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-Ch0796	Chicken CA1(Carbonic Anhydrase I) ELISA Kit	CA1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0797	Chicken C3c(Complement Fragment 3c) ELISA Kit	C3c	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0798	Chicken BTC(Betacellulin) ELISA Kit	BTC	Sandwich	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-Ch0799	Chicken MEPE(Matrix Extracellular Phosphoglycoprotein) ELISA Kit	MEPE	Sandwich	60ng/mL	0.56 ng/mL	0.94-60ng/mL
FY-Ch0800	Chicken BDNF(Brain Derived Neurotrophic Factor) ELISA Kit	BDNF	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0801	Chicken Pro-MMP-10(Pro-Matrix Metalloproteinase 10) ELISA Kit	Pro-MMP-10	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0802	Chicken ELA2B(Elastase 2B) ELISA Kit	ELA2B	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0803	Chicken PIK3Cδ(Phosphoinositide-3-Kinase Catalytic Delta Polypeptide) ELISA Kit	PIK3Cδ	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0804	Chicken MANF(Mesencephalic Astrocyte Derived Neurotrophic Factor) ELISA Kit	MANF	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0805	Chicken PLP(Parathyroid Hormone-Like Protein) ELISA Kit	PLP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0806	Chicken PTGIS(Prostaglandin I Synthase) ELISA Kit	PTGIS	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0807	Chicken S100A9(S100 Calcium Binding Protein A9) ELISA Kit	S100A9	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0808	Chicken DAF(Decay Accelerating Factor) ELISA Kit	DAF	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL

FY-Ch0809	Chicken NRG-2(Neuregulin 2) ELISA Kit	NRG-2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0810	Chicken ARSA(Arylsulfatase A) ELISA Kit	ARSA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0811	Chicken INSR(Insulin Receptor) ELISA Kit	INSR	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0812	Chicken FMOD(Fibromodulin) ELISA Kit	FMOD	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0813	Chicken P2(Myelin Protein 2) ELISA Kit	P2	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0814	Chicken Sp100(Sp100 Nuclear Antigen) ELISA Kit	Sp100	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0815	Chicken Surv(Survivin) ELISA Kit	Surv	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0816	Chicken OC/BGP(Osteocalcin) ELISA Kit	OC/BGP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0817	Chicken COL13(Collagen Type X III) ELISA Kit	COL13	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0818	Chicken GPRC5A(G Protein Coupled Receptor, Family C, Group 5, Member A) ELISA Kit	GPRC5A	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0819	Chicken FPR2(Formyl Peptide Receptor 2) ELISA Kit	FPR2	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0821	Chicken CTLA4(Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit	CTLA4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0822	Chicken CVL/EZR(Cytovillin/Ezrin) ELISA Kit	CVL/EZR	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0823	Chicken NOS2/iNOS(Nitric Oxide Synthase 2, Inducible) ELISA Kit	NOS2/iNOS	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0824	Chicken GDF9(Growth Differentiation Factor 9) ELISA Kit	GDF9	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0825	Chicken GAP43(Growth Associated Protein 43) ELISA Kit	GAP43	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0826	Chicken GAD(Glutamic Acid Decarboxylase) ELISA Kit	GAD	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0827	Chicken SDF-1/CXCL12(Stromal Cell Derived Factor 1) ELISA Kit	SDF-1/CXCL12	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0829	Chicken LATS(Long Acting Thyroid Stimulator) ELISA Kit	LATS	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0830	Chicken BLC(B-Lymphocyte Chemoattractant) ELISA Kit	BLC	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0832	Chicken ASGR1(Asialoglycoprotein Receptor 1) ELISA Kit	ASGR1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch0833	Chicken SMN1(Survival of Motor Neuron 1, Telomeric) ELISA Kit	SMN1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0834	Chicken TNX(Tenascin X) ELISA Kit	TNX	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0835	Chicken F12(Coagulation Factor XII) ELISA Kit	F12	Sandwich	150pg/mL	1.4 pg/mL	2.34-150pg/mL
FY-Ch0836	Chicken IL-12(Interleukin 12) ELISA Kit	IL-12	Sandwich	800pg/mL	7.5 pg/mL	12.5-800pg/mL
FY-Ch0837	Chicken FSTL1(Follistatin Like Protein 1) ELISA Kit	FSTL1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0838	Chicken COL4 α 1(Collagen Type IV Alpha 1) ELISA Kit	COL4 α 1	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0839	Chicken AGR2(Anterior Gradient Protein 2) ELISA Kit	AGR2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0840	Chicken CKLF(Chemokine Like Factor) ELISA Kit	CKLF	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0841	Chicken D5R/DRD5(Dopamine Receptor D5) ELISA Kit	D5R/DRD5	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0842	Chicken PDCN(Podocin) ELISA Kit	PDCN	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0843	Chicken AhR(Aryl Hydrocarbon Receptor) ELISA Kit	AhR	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0844	Chicken MYD88(Myeloid Differentiation Factor 88) ELISA Kit	MYD88	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0845	Chicken α 1ACT(Alpha1 Antichymotrypsin) ELISA Kit	α 1ACT	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-Ch0846	Chicken HSF2(Heat Shock Transcription Factor 2) ELISA Kit	HSF2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0847	Chicken c-myc(c-myc Oncogene Product) ELISA Kit	c-myc	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0848	Chicken PYGB(Glycogen Phosphorylase, Brain) ELISA Kit	PYGB	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0849	Chicken M-AChR M1(Muscarinic Acetylcholine Receptor M1) ELISA Kit	M-AChR M1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0850	Chicken L1CAM/CD171(L1-Cell Adhesion Molecule) ELISA Kit	L1CAM/CD171	Sandwich	300pg/mL	2.81 pg/mL	4.69-300pg/mL
FY-Ch0851	Chicken RND1(Rho Family GTPase 1) ELISA Kit	RND1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0852	Chicken MUC3(Mucin 3) ELISA Kit	MUC3	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL

FY-Ch0853	Chicken TF(Transferrin) ELISA Kit	TF	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0854	Chicken CCND1(Cyclin D1) ELISA Kit	CCND1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0855	Chicken LOXL4(Lysyl Oxidase Like Protein 4) ELISA Kit	LOXL4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0856	Chicken D3R/DRD3(Dopamine Receptor D3) ELISA Kit	D3R/DRD3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0857	Chicken GLUT4(Glucose Transporter 4) ELISA Kit	GLUT4	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch0858	Chicken ApoC2(Apolipoprotein C2) ELISA Kit	ApoC2	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch0859	Chicken MIP/CXCL9(Monocyte Interferon Gamma Inducing Factor) ELISA Kit	MIP/CXCL9	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0860	Chicken CLDN9(Claudin 9) ELISA Kit	CLDN9	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0861	Chicken VEGFR3/Flt4(Vascular Endothelial Cell Growth Factor Receptor 3) ELISA Kit	VEGFR3/Flt4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0862	Chicken SPARCL1(SPARC Like Protein 1) ELISA Kit	SPARCL1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0863	Chicken VCL(Vinculin) ELISA Kit	VCL	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0864	Chicken SDF2L1(Stromal Cell Derived Factor 2 Like Protein 1) ELISA Kit	SDF2L1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0865	Chicken SDF-4(Stromal Cell Derived Factor 4) ELISA Kit	SDF-4	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0866	Chicken NAIP(Neuronal Apoptosis Inhibitory Protein) ELISA Kit	NAIP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0867	Chicken PLAUR/uPAR(Plasminogen Activator, Urokinase Receptor) ELISA Kit	PLAUR/uPAR	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0868	Chicken PTPRJ(Protein Tyrosine Phosphatase Receptor Type J) ELISA Kit	PTPRJ	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0869	Chicken HO2(Heme Oxygenase 2) ELISA Kit	HO2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0870	Chicken PMN Elastase(Polymorphonuclear Elastase) ELISA Kit	PMN Elastase	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0871	Chicken GDF7(Growth Differentiation Factor 7) ELISA Kit	GDF7	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0872	Chicken CRT(Calreticulin) ELISA Kit	CRT	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0873	Chicken MFI2(Melanotransferrin) ELISA Kit	MFI2	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0874	Chicken CDNF(Cerebral Dopamine Neurotrophic Factor) ELISA Kit	CDNF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0875	Chicken FcεR II (Receptor II for the Fc Fragment of IgE) ELISA Kit	FcεR II	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0876	Chicken DEFα1(Defensin Alpha 1, Neutrophil) ELISA Kit	DEFα1	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0878	Chicken CDX2(Caudal Type Homeobox Transcription Factor 2) ELISA Kit	CDX2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0879	Chicken PKP1(Plakophilin 1) ELISA Kit	PKP1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0880	Chicken SPOP(Speckle Type POZ Protein) ELISA Kit	SPOP	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0881	Chicken SAA(Serum Amyloid A) ELISA Kit	SAA	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0882	Chicken PFKM(Phosphofructokinase, Muscle) ELISA Kit	PFKM	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0883	Chicken GDF3(Growth Differentiation Factor 3) ELISA Kit	GDF3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0884	Chicken ORM2(Orosomucoid 2) ELISA Kit	ORM2	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0885	Chicken HSPG(Heparan Sulfate Proteoglycan) ELISA Kit	HSPG	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0886	Chicken AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit	AGER	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0887	Chicken P-CK(Pan-Cytokeratin) ELISA Kit	P-CK	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0888	Chicken LILRB4(Leukocyte Immunoglobulin Like Receptor Subfamily B Member 4) ELISA Kit	LILRB4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0889	Chicken CNR-1(Cadherin Related Neuronal Receptor1) ELISA Kit	CNR-1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0890	Chicken FLNα(Filamin A, Alpha) ELISA Kit	FLNα	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0891	Chicken PIK3Cβ(Phosphoinositide-3-Kinase Catalytic Beta Polypeptide) ELISA Kit	PIK3Cβ	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0893	Chicken DKK1(Dickkopf Related Protein 1) ELISA Kit	DKK1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0894	Chicken CFD(Complement Factor D) ELISA Kit	CFD	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch0895	Chicken GPI(Glucose 6 Phosphate Isomerase) ELISA Kit	GPI	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL

FY-Ch0896	Chicken TNF-β(Tumor Necrosis Factor Beta) ELISA Kit	TNF-β	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0897	Chicken PCDH1(Protocadherin 1) ELISA Kit	PCDH1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0898	Chicken NEFL(Neurofilament, Light Polypeptide) ELISA Kit	NEFL	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0899	Chicken MUC7(Mucin 7) ELISA Kit	MUC7	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0900	Chicken ANG-R-Tie2(Angiopoietin Receptor Tie2) ELISA Kit	ANG-R-Tie2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0901	Chicken SLPI(Secretory Leukocyte Protease Inhibitor) ELISA Kit	SLPI	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0902	Chicken CHD5(Chromodomain Helicase DNA Binding Protein 5) ELISA Kit	CHD5	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0903	Chicken IL-8(Interleukin 8) ELISA Kit	IL-8	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0904	Chicken C/EBPα(CCAAT/Enhancer Binding Protein Alpha) ELISA Kit	C/EBPα	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0905	Chicken C5(Complement Component 5) ELISA Kit	C5	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0906	Chicken VWF(Von Willebrand Factor) ELISA Kit	VWF	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0907	Chicken IRAK3(Interleukin 1 Receptor Associated Kinase 3) ELISA Kit	IRAK3	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0908	Chicken CEBPγ(CCAAT/Enhancer Binding Protein Gamma) ELISA Kit	CEBPγ	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0909	Chicken VEGF-C(Vascular Endothelial Growth Factor C) ELISA Kit	VEGF-C	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0910	Chicken OCT2B(Octamer-binding protein 2B) ELISA Kit	OCT2B	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch0911	Chicken CR2/CD21(Complement Receptor 2) ELISA Kit	CR2/CD21	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0912	Chicken FTL(Ferritin, Light Polypeptide) ELISA Kit	FTL	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0913	Chicken MASP1(Mannan Associated Serine Protease 1) ELISA Kit	MASP1	Sandwich	300pg/mL	2.81 pg/mL	4.69-300pg/mL
FY-Ch0914	Chicken CDK8(Cyclin Dependent Kinase 8) ELISA Kit	CDK8	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0915	Chicken UBC(Urinary Bladder Cancer Antigen) ELISA Kit	UBC	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0916	Chicken MAGE(Melanoma-associated Antigen) ELISA Kit	MAGE	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0917	Chicken PFKP(Phosphofructokinase, Platelet) ELISA Kit	PFKP	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0919	Chicken FGγ(Fibrinogen Gamma chain) ELISA Kit	FGγ	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0920	Chicken TDP43(Tar DNA Binding Protein 43KDa) ELISA Kit	TDP43	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0921	Chicken FSHR(Follicle Stimulating Hormone Receptor) ELISA Kit	FSHR	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0922	Chicken NOX1(Nicotinamide Adenine Dinucleotide Phosphate Oxidase 1) ELISA Kit	NOX1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0923	Chicken HELIX- II (Type II Collagen Helical Peptide) ELISA Kit	HELIX- II	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0924	Chicken MAPK14(Mitogen Activated Protein Kinase 14) ELISA Kit	MAPK14	Sandwich	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-Ch0925	Chicken EAR12(Eosinophil Associated Ribonuclease A family Member 12) ELISA Kit	EAR12	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0926	Chicken HNRPU(Heterogeneous Nuclear Ribonucleoprotein U) ELISA Kit	HNRPU	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0927	Chicken OTR(Oxytocin Receptor) ELISA Kit	OTR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0928	Chicken ADRA1A(Adrenergic Receptor alpha-1A) ELISA Kit	ADRA1A	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0929	Chicken ApoA2(Apolipoprotein A2) ELISA Kit	ApoA2	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-Ch0930	Chicken CTNNβ1(Catenin, Beta 1) ELISA Kit	CTNNβ1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0931	Chicken GlyCAM1(Glycosylation Dependent Cell Adhesion Molecule 1) ELISA Kit	GlyCAM1	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0932	Chicken CENPA(Centromere Protein A) ELISA Kit	CENPA	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0934	Chicken TBP(TATA Binding Protein) ELISA Kit	TBP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0935	Chicken TNFRSF1B(Tumor Necrosis Factor Receptor Superfamily, Member 1B) ELISA Kit	TNFRSF1B	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0936	Chicken CTX- II (Cross Linked C-telopeptide of Type II Collagen) ELISA Kit	CTX- II	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0937	Chicken TBP2(Thioredoxin Binding Protein 2) ELISA Kit	TBP2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0938	Chicken PAPP(A) (Pregnancy Associated Plasma Protein A) ELISA Kit	PAPP(A)	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL

FY-Ch0939	Chicken AVPR1A(Arginine Vasopressin Receptor 1A) ELISA Kit	AVPR1A	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch0940	Chicken PDGFRL(Platelet Derived Growth Factor Receptor Like Protein) ELISA Kit	PDGFRL	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0941	Chicken PGDS(Prostaglandin D Synthase) ELISA Kit	PGDS	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0942	Chicken PAR1(Protease Activated Receptor 1) ELISA Kit	PAR1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0943	Chicken MUC13(Mucin 13) ELISA Kit	MUC13	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0944	Chicken DAG1(Dystrophin Associated Glycoprotein 1) ELISA Kit	DAG1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0945	Chicken HCCS(Holocytochrome C Synthase) ELISA Kit	HCCS	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0946	Chicken CA3(Carbonic Anhydrase III, Muscle Specific) ELISA Kit	CA3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0947	Chicken MCP-3(Monocyte Chemotactic Protein 3) ELISA Kit	MCP-3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0948	Chicken CFP(Complement Factor P) ELISA Kit	CFP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0949	Chicken CDK6(Cyclin Dependent Kinase 6) ELISA Kit	CDK6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0950	Chicken ANG(Angiostatin) ELISA Kit	ANG	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0951	Chicken PLA2G12B(Phospholipase A2, Group XII B) ELISA Kit	PLA2G12B	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0952	Chicken DPYD(Dihydropyrimidine Dehydrogenase) ELISA Kit	DPYD	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0953	Chicken EED(Embryonic Ectoderm Development) ELISA Kit	EED	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0954	Chicken REN(Renin) ELISA Kit	REN	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0955	Chicken FSCN(Fascin) ELISA Kit	FSCN	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0956	Chicken MYH9(Myosin Heavy Chain 9, Non Muscle) ELISA Kit	MYH9	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0957	Chicken PCK1(Phosphoenolpyruvate Carboxykinase 1) ELISA Kit	PCK1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0958	Chicken CENPF(Centromere Protein F) ELISA Kit	CENPF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0959	Chicken ANG-R-Tie1(Angiopoietin Receptor Tie1) ELISA Kit	ANG-R-Tie1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0960	Chicken FASN(Fatty Acid Synthase) ELISA Kit	FASN	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0961	Chicken NOS3/eNOS(Nitric Oxide Synthase 3, Endothelial) ELISA Kit	NOS3/eNOS	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0962	Chicken α 2-PI(α 2-Plasmin Inhibitor) ELISA Kit	α 2-PI	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0963	Chicken CK-18/KRT18(Cytokeratin 18) ELISA Kit	CK-18/KRT18	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0964	Chicken PKM2(Pyruvate Kinase, Muscle) ELISA Kit	PKM2	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch0965	Chicken PTHrP(Parathyroid Hormone Related Protein) ELISA Kit	PTHrP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0966	Chicken PKC ϵ (Protein Kinase C Epsilon) ELISA Kit	PKC ϵ	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0968	Chicken GC(Glycocalicin) ELISA Kit	GC	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0969	Chicken NTX1(Cross Linked N-telopeptide of Type I Collagen) ELISA Kit	NTX1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0970	Chicken FN(Fibronectin) ELISA Kit	FN	Sandwich	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-Ch0971	Chicken HSP-27/HSPB1(Heat Shock Protein 27) ELISA Kit	HSP-27/HSPB1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0972	Chicken PDK4(Pyruvate Dehydrogenase Kinase Isozyme 4) ELISA Kit	PDK4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0973	Chicken PROC(Protein C) ELISA Kit	PROC	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0974	Chicken NRIP1(Nuclear Receptor Interacting Protein 1) ELISA Kit	NRIP1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0975	Chicken VEGF-D(Vascular Endothelial Growth Factor D) ELISA Kit	VEGF-D	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0976	Chicken PTGES2(Prostaglandin E Synthase 2) ELISA Kit	PTGES2	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch0977	Chicken UGT2B4(Uridine Diphosphate Glucuronosyltransferase 2 Family Polypeptide) ELISA Kit	UGT2B4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0978	Chicken MIP-3 β /ELC/CCL19(Macrophage Inflammatory Protein 3 β) ELISA Kit	MIP-3 β /ELC/CCL19	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0979	Chicken VTN(Vitronectin) ELISA Kit	VTN	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0980	Chicken GATA4(GATA Binding Protein 4) ELISA Kit	GATA4	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL

FY-Ch0981	Chicken FGFBP1(Fibroblast Growth Factor Binding Protein 1) ELISA Kit	FGFBP1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0982	Chicken CPP(Copeptin) ELISA Kit	CPP	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0983	Chicken IGFBP-5(Insulin Like Growth Factor Binding Protein 5) ELISA Kit	IGFBP-5	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0984	Chicken COL3a1(Collagen Type III Alpha 1) ELISA Kit	COL3a1	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch0985	Chicken CASP14(Caspase 14) ELISA Kit	CASP14	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0986	Chicken SHC/SLP-76(SH2-Containing Protein) ELISA Kit	SHC/SLP-76	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0987	Chicken PLN(Phospholamban) ELISA Kit	PLN	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch0988	Chicken SEPCR/PROCR(Soluble Endothelial Protein C Receptor) ELISA Kit	SEPCR/PROCR	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch0989	Chicken GCK(Glucokinase) ELISA Kit	GCK	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch0990	Chicken PDGFRβ(Platelet Derived Growth Factor Receptor Beta) ELISA Kit	PDGFRβ	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0991	Chicken EREG(Epiregulin) ELISA Kit	EREG	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0992	Chicken TALLA-1(T-cell Acute Lymphoblastic Leukemia Antigen) ELISA Kit	TALLA-1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch0993	Chicken FGF7/KGF(Fibroblast Growth Factor 7) ELISA Kit	FGF7/KGF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0994	Chicken ECF/CCL11(Eosinophil Chemotactic Factor) ELISA Kit	ECF/CCL11	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0995	Chicken OT(Ovotransferrin) ELISA Kit	OT	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0996	Chicken SIGLEC10(Sialic Acid Binding Immunoglobulin Like Lectin 10) ELISA Kit	SIGLEC10	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0997	Chicken FABP7(Fatty Acid Binding Protein 7, Brain) ELISA Kit	FABP7	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch0998	Chicken PGAM2(Phosphoglycerate Mutase 2) ELISA Kit	PGAM2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0999	Chicken PPAR-α(Peroxisome Proliferator Activated Receptor Alpha) ELISA Kit	PPAR-α	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1000	Chicken CSF2Rα(Colony Stimulating Factor 2 Receptor Alpha) ELISA Kit	CSF2Rα	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1001	Chicken TNNI2(Troponin I Type 2, Fast Skeletal) ELISA Kit	TNNI2	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1002	Chicken TAGLN(Transgelin) ELISA Kit	TAGLN	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1003	Chicken FSH(Follicle Stimulating Hormone) ELISA Kit	FSH	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1004	Chicken PRB1(Basic Salivary Proline Rich Protein 1) ELISA Kit	PRB1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1005	Chicken TACR2(Tachykinin Receptor 2) ELISA Kit	TACR2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1006	Chicken α Actin(Alpha-Actin) ELISA Kit	α Actin	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1007	Chicken PIAS1(Protein Inhibitor of Activated STAT 1) ELISA Kit	PIAS1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1008	Chicken TGM3(Transglutaminase 3, Epidermal) ELISA Kit	TGM3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1009	Chicken CRY1(Cryptochrome 1) ELISA Kit	CRY1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1010	Chicken CNTN1(Contactin 1) ELISA Kit	CNTN1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1011	Chicken TERT(Telomerase Reverse Transcriptase) ELISA Kit	TERT	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1012	Chicken CD3d(Cluster of Differentiation 3d) ELISA Kit	CD3d	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1013	Chicken STAT4(Signal Transducer and Activator of Transcription 4) ELISA Kit	STAT4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1014	Chicken IL-2Rα(Interleukin 2 Receptor Alpha) ELISA Kit	IL-2Rα	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1015	Chicken SERT(Serotonin Transporter) ELISA Kit	SERT	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1016	Chicken MIP-5(Macrophage Inflammatory Protein 5) ELISA Kit	MIP-5	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1017	Chicken CXCL16(Chemokine C-X-C-Motif Ligand 16) ELISA Kit	CXCL16	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1019	Chicken COL17(Collagen Type X VII) ELISA Kit	COL17	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1020	Chicken PRSS23(Protease, Serine 23) ELISA Kit	PRSS23	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1021	Chicken VDR(Vitamin D Receptor) ELISA Kit	VDR	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1022	Chicken MMSAM(Melanoma Metastasis Surface Adhesion Molecule) ELISA Kit	MMSAM	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL

FY-Ch1023	Chicken ERCC1(Excision Repair Cross Complementing Rodent Repair Deficiency C ERCC1		Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1024	Chicken FBLN4(Fibulin 4) ELISA Kit	FBLN4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1025	Chicken C/EBPδ(CCAAT/Enhancer Binding Protein, Delta) ELISA Kit	C/EBPδ	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1026	Chicken ADCY2(Adenylate Cyclase 2) ELISA Kit	ADCY2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1027	Chicken TNNC2(Troponin C Type 2) ELISA Kit	TNNC2	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1028	Chicken FGFR2(Fibroblast Growth Factor Receptor 2) ELISA Kit	FGFR2	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1029	Chicken NUP160(Nucleoporin 160kDa) ELISA Kit	NUP160	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1030	Chicken ERK1(Extracellular Signal Regulated Kinase 1) ELISA Kit	ERK1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1031	Chicken CYSLTR2(CysteinyL Leukotriene Receptor 2) ELISA Kit	CYSLTR2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1032	Chicken Ly6a(Lymphocyte Antigen 6 Complex Locus A) ELISA Kit	Ly6a	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1033	Chicken CENPE(Centromere Protein E) ELISA Kit	CENPE	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1034	Chicken MYBPC3(Myosin Binding Protein C, Cardiac) ELISA Kit	MYBPC3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1035	Chicken CEA(Carcinoembryonic Antigen) ELISA Kit	CEA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1036	Chicken CYP1A1(Cytochrome P450, family 1, subfamily A, polypeptide 1) ELISA K	CYP1A1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1037	Chicken tPA(Tissue-type Plasminogen Activator) ELISA Kit	tPA	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1038	Chicken FCN1(Ficolin 1) ELISA Kit	FCN1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1039	Chicken IL-15(Interleukin 15) ELISA Kit	IL-15	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1040	Chicken LEPR(Leptin Receptor) ELISA Kit	LEPR	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1041	Chicken PYGM(Glycogen Phosphorylase, Muscle) ELISA Kit	PYGM	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1042	Chicken AIF1(Allograft Inflammatory Factor 1) ELISA Kit	AIF1	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1043	Chicken GDF10(Growth Differentiation Factor 10) ELISA Kit	GDF10	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1044	Chicken PDAP1(Pdgfa Associated Protein 1) ELISA Kit	PDAP1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1045	Chicken Salusin α(Salusin Alpha) ELISA Kit	Salusin α	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1046	Chicken ABCG2(ATP Binding Cassette Transporter G2) ELISA Kit	ABCG2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1047	Chicken F-Box Protein 2(FBXO2) ELISA Kit	F-Box Protein 2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1048	Chicken BLVRB(Biliverdin Reductase B) ELISA Kit	BLVRB	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1049	Chicken TUBβ(β-Tubulin) ELISA Kit	TUBβ	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1050	Chicken EAR3(Eosinophil Associated Ribonuclease A family Member 3) ELISA Kit	EAR3	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch1051	Chicken PTF1α(Pancreas Specific Transcription Factor 1α) ELISA Kit	PTF1α	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1052	Chicken SF1(Steroidogenic Factor 1) ELISA Kit	SF1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1053	Chicken SDC3(Syndecan 3) ELISA Kit	SDC3	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1055	Chicken AQP-4(Aquaporin 4) ELISA Kit	AQP-4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1056	Chicken PAI1(Plasminogen Activator Inhibitor 1) ELISA Kit	PAI1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1057	Chicken TnTf/TNNT3(Troponin T, fast skeletal muscle) ELISA Kit	TnTf/TNNT3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1058	Chicken PLG(Plasminogen) ELISA Kit	PLG	Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-Ch1059	Chicken DBH(Dopamine-β-Hydroxylase) ELISA Kit	DBH	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1060	Chicken FE(Ferritin) ELISA Kit	FE	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch1061	Chicken PLP(Proteolipid Protein) ELISA Kit	PLP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1062	Chicken OxLDL(Oxidized Low Density Lipoprotein) ELISA Kit	OxLDL	Sandwich	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-Ch1063	Chicken CA724(Tumor Marker CA724) ELISA Kit	CA724	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1064	Chicken MYO(Myoglobin) ELISA Kit	MYO	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL

FY-Ch1065	Chicken TP53(Tumor Protein p53) ELISA Kit	TP53	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1066	Chicken C1INH(Complement 1 Inhibitor) ELISA Kit	C1INH	Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-Ch1067	Chicken CLEC4E(C-Type Lectin Domain Family 4 Member E) ELISA Kit	CLEC4E	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1068	Chicken CNTN2(Contactin 2) ELISA Kit	CNTN2	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1069	Chicken CCR8(Chemokine C-C-Motif Receptor 8) ELISA Kit	CCR8	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1070	Chicken ALB(Albumin) ELISA Kit	ALB	Sandwich	200µg/mL	1.88 µg/mL	3.13-200µg/mL
FY-Ch1071	Chicken OPG(Osteoprotegerin) ELISA Kit	OPG	Sandwich	8000pg/mL	75 pg/mL	125-8000pg/mL
FY-Ch1072	Chicken TRAIL/TNFSF10(Tumor Necrosis Factor Related Apoptosis Inducing Ligan TRAIL/TNFSF10	TRAIL/TNFSF10	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1073	Chicken GSK3β(Glycogen Synthase Kinase 3 Beta) ELISA Kit	GSK3β	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1074	Chicken Gas6(Growth Arrest Specific Protein 6) ELISA Kit	Gas6	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1075	Chicken COL7(Collagen Type VII) ELISA Kit	COL7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1076	Chicken IL-1R2(Interleukin 1 Receptor Type II) ELISA Kit	IL-1R2	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1077	Chicken FGα(Fibrinogen Alpha) ELISA Kit	FGα	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1078	Chicken TECK/CCL25(Thymus Expressed Chemokine) ELISA Kit	TECK/CCL25	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1079	Chicken PKB(Protein Kinase B) ELISA Kit	PKB	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1080	Chicken UCHL1(Ubiquitin Carboxyl Terminal Hydrolase L1) ELISA Kit	UCHL1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1081	Chicken LIFR(Leukemia Inhibitory Factor Receptor) ELISA Kit	LIFR	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1082	Chicken ANGPTL4(Angiopoietin Like Protein 4) ELISA Kit	ANGPTL4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1084	Chicken KIM-1(Kidney Injury Molecule 1) ELISA Kit	KIM-1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1085	Chicken C4BP(C4 Binding Protein) ELISA Kit	C4BP	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1086	Chicken MUC4(Mucin 4) ELISA Kit	MUC4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1087	Chicken TIMP-2(Tissue Inhibitors of Metalloproteinase 2) ELISA Kit	TIMP-2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1088	Chicken CREB1(Cyclic AMP Response Element Binding Protein 1) ELISA Kit	CREB1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1089	Chicken D2R/DRD2(Dopamine Receptor D2) ELISA Kit	D2R/DRD2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1090	Chicken ApoC1(Apolipoprotein C1) ELISA Kit	ApoC1	Sandwich	5000ng/mL	46.88 ng/mL	78.13-5000ng/mL
FY-Ch1091	Chicken DR-70TM(Tumor Marker DR-70 For Lung Cancer) ELISA Kit	DR-70TM	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1092	Chicken AIM2(Absent In Melanoma 2) ELISA Kit	AIM2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1093	Chicken CNR1(Cannabinoid Receptor 1, Brain) ELISA Kit	CNR1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1094	Chicken PDP(Pyruvate Dehydrogenase Phosphatase) ELISA Kit	PDP	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1095	Chicken P11NP(N-Terminal Procollagen III Propeptide) ELISA Kit	P11NP	Sandwich	80pg/mL	0.75 pg/mL	1.25-80pg/mL
FY-Ch1096	Chicken APP(Amyloid Precursor Protein) ELISA Kit	APP	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1097	Chicken F1+2(Prothrombin Fragment 1+2) ELISA Kit	F1+2	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-Ch1098	Chicken ApoH(Apolipoprotein H) ELISA Kit	ApoH	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1099	Chicken MYCBP(C-Myc Binding Protein) ELISA Kit	MYCBP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1100	Chicken NRG-4(Neuregulin 4) ELISA Kit	NRG-4	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch1101	Chicken PAC3(Platelet Associated Complement 3) ELISA Kit	PAC3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1102	Chicken UGCG(UDP-glucose ceramide glucosyltransferase) ELISA Kit	UGCG	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1103	Chicken EFNA4(Ephrin A4) ELISA Kit	EFNA4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1104	Chicken FDFT1(Farnesyl Diphosphate Farnesyltransferase 1) ELISA Kit	FDFT1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1105	Chicken ERβ(Estrogen Receptor Beta) ELISA Kit	ERβ	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1106	Chicken SIAE(Sialic Acid Acetyltransferase) ELISA Kit	SIAE	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-Ch1107	Chicken CFTR(Cystic Fibrosis Transmembrane Conductance Regulator) ELISA Kit	CFTR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1108	Chicken CLDN14(Claudin 14) ELISA Kit	CLDN14	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1109	Chicken BAFF/CD257(B-Cell Activating Factor) ELISA Kit	BAFF/CD257	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1111	Chicken PEX2(Peroxisomal Biogenesis Factor 2) ELISA Kit	PEX2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1112	Chicken TNNI1(Troponin I Type 1, Slow Skeletal) ELISA Kit	TNNI1	Sandwich	150pg/mL	1.4 pg/mL	2.34-150pg/mL
FY-Ch1113	Chicken IL-3(Interleukin 3) ELISA Kit	IL-3	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1114	Chicken PYGL(Glycogen Phosphorylase, Liver) ELISA Kit	PYGL	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1115	Chicken PLA2G10(Phospholipase A2, Group X) ELISA Kit	PLA2G10	Sandwich	25pg/mL	0.23 pg/mL	0.39-25pg/mL
FY-Ch1116	Chicken OTC(Ornithine Carbamoyl Transferase) ELISA Kit	OTC	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1117	Chicken CX26(Connexin 26) ELISA Kit	CX26	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1118	Chicken MT2(Metallothionein 2) ELISA Kit	MT2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1120	Chicken CYTIP(Cytohesin Interacting Protein) ELISA Kit	CYTIP	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1121	Chicken PDE4D(Phosphodiesterase 4D, cAMP Specific) ELISA Kit	PDE4D	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1122	Chicken 5-NT(5-Nucleotidase) ELISA Kit	5-NT	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1123	Chicken ACTg1(Actin Gamma 1) ELISA Kit	ACTg1	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1124	Chicken CDC23(Cell Division Cycle Protein 23) ELISA Kit	CDC23	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1125	Chicken GSTk1(Glutathione S Transferase Kappa 1) ELISA Kit	GSTk1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1126	Chicken ANG II R1(Angiotensin II Receptor 1) ELISA Kit	ANG II R1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1127	Chicken CHD3(Chromodomain Helicase DNA Binding Protein 3) ELISA Kit	CHD3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1128	Chicken LBP(Lipoplysaccharide Binding Protein) ELISA Kit	LBP	Sandwich	1200ng/mL	11.25 ng/mL	18.75-1200ng/mL
FY-Ch1129	Chicken SIGLEC8(Sialic Acid Binding Immunoglobulin Like Lectin 8) ELISA Kit	SIGLEC8	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1130	Chicken HGF(Hepatocyte Growth Factor) ELISA Kit	HGF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1132	Chicken sCD14(Soluble Cluster of Differentiation14) ELISA Kit	sCD14	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1133	Chicken EPO(Erythropoietin) ELISA Kit	EPO	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1134	Chicken MRC1(Mannose Receptor C Type 1) ELISA Kit	MRC1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1135	Chicken IGF2BP3(Insulin Like Growth Factor 2 mRNA Binding Protein 3) ELISA Kit	IGF2BP3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1136	Chicken LAT(Linker For Activation of T-cell) ELISA Kit	LAT	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1137	Chicken AKAP1(A Kinase Anchor Protein 1) ELISA Kit	AKAP1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1138	Chicken TSGF(Tumor Specific Growth Factor/Tumor Supplied Group of Factor) ELISA Kit	TSGF	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1139	Chicken GKN2(Gastrokine 2) ELISA Kit	GKN2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1140	Chicken PDH E1(Pyruvate Dehydrogenase E1) ELISA Kit	PDH E1	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1142	Chicken FSA(Fetal Sulfoslycoprotein Antigen) ELISA Kit	FSA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1143	Chicken TH(Tyrosine Hydroxylase) ELISA Kit	TH	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1144	Chicken SPTAN1(Alpha-Fodrin) ELISA Kit	SPTAN1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1145	Chicken GIP(Glucose Dependent Insulin Releasing Polypeptide) ELISA Kit	GIP	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch1146	Chicken FDP(Fibrinogen Degradation Product) ELISA Kit	FDP	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1147	Chicken CCKAR(Cholecystokinin A Receptor) ELISA Kit	CCKAR	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1148	Chicken OSF(Osteoclast Stimulating Factor 1) ELISA Kit	OSF	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1149	Chicken MCSF(Macrophage Colony Stimulating Factor 1) ELISA Kit	MCSF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1150	Chicken INSL3(Insulin Like Protein 3) ELISA Kit	INSL3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1151	Chicken S100A6(S100 Calcium Binding Protein A6) ELISA Kit	S100A6	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-Ch1152	Chicken GDNF(Glial Cell Line Derived Neurotrophic Factor) ELISA Kit	GDNF	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1153	Chicken FP(Prostaglandin F Receptor) ELISA Kit	FP	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1154	Chicken GKN1(Gastrokine 1) ELISA Kit	GKN1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1155	Chicken CYB561(Cytochrome b-561) ELISA Kit	CYB561	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1156	Chicken DPYSL3(Dihydropyrimidinase Like Protein 3) ELISA Kit	DPYSL3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1157	Chicken ETFβ(Electron Transfer Flavoprotein Beta Polypeptide) ELISA Kit	ETFβ	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1158	Chicken BMP-1(Bone Morphogenetic Protein 1) ELISA Kit	BMP-1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1159	Chicken LH(Luteinizing Hormone) ELISA Kit	LH	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1160	Chicken LIPH(Lipase H) ELISA Kit	LIPH	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1161	Chicken MMP-10(Matrix Metalloproteinase 10) ELISA Kit	MMP-10	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1162	Chicken PDGF-AB(Platelet Derived Growth Factor AB) ELISA Kit	PDGF-AB	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1163	Chicken RNASET2(Ribonuclease T2) ELISA Kit	RNASET2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1164	Chicken FOXO3(Forkhead Box Protein O3) ELISA Kit	FOXO3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1165	Chicken TLR2(Toll-Like Receptor 2) ELISA Kit	TLR2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1167	Chicken TTN(Titin) ELISA Kit	TTN	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1168	Chicken CEP110(Centrosomal Protein 110kDa) ELISA Kit	CEP110	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1169	Chicken GDF2(Growth Differentiation Factor 2) ELISA Kit	GDF2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1170	Chicken NCAD(Neural Cadherin) ELISA Kit	NCAD	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1171	Chicken DLK-1(Delta Like 1 Homolog) ELISA Kit	DLK-1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1172	Chicken SCF(Stem Cell Factor) ELISA Kit	SCF	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1173	Chicken TR(Thrombin Receptor) ELISA Kit	TR	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1174	Chicken CALU(Calumenin) ELISA Kit	CALU	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1175	Chicken DPP4(Dipeptidyl Peptidase IV) ELISA Kit	DPP4	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1177	Chicken α1BG(Alpha-1-B-Glycoprotein) ELISA Kit	α1BG	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1178	Chicken UPP2(Uridine Phosphorylase 2) ELISA Kit	UPP2	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1179	Chicken FETUA(Fetuin A) ELISA Kit	FETUA	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-Ch1180	Chicken PLA1(Phospholipase A1) ELISA Kit	PLA1	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch1181	Chicken TAP1(Transporter associated with Antigen Processing 1) ELISA Kit	TAP1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1182	Chicken NFYC(Nuclear Transcription Factor Y Gamma) ELISA Kit	NFYC	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1183	Chicken VE-Cadherin(Vascular Endothelial Cadherin) ELISA Kit	VE-Cadherin	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1184	Chicken DLD(Dihyrolipoyl Dehydrogenase) ELISA Kit	DLD	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1185	Chicken SPD(Pulmonary Surfactant Associated Protein D) ELISA Kit	SPD	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1186	Chicken NNE(Non-Neuronal Enolase) ELISA Kit	NNE	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1188	Chicken TXLNb(Taxilin Beta) ELISA Kit	TXLNb	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1189	Chicken ABCG1(ATP Binding Cassette Transporter G1) ELISA Kit	ABCG1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1190	Chicken TNKS1(Tankyrase 1) ELISA Kit	TNKS1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1191	Chicken CTAPIII(Connective Tissue Activating Peptide III) ELISA Kit	CTAPIII	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1192	Chicken CASP11(Caspase 11) ELISA Kit	CASP11	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1193	Chicken CD4(Cluster of Differentiation 4) ELISA Kit	CD4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1194	Chicken PGR(Progesterone Receptor) ELISA Kit	PGR	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1195	Chicken CCL4L1(Chemokine C-C-Motif Ligand 4 Like Protein 1) ELISA Kit	CCL4L1	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL

FY-Ch1196	Chicken SFRP2(Secreted Frizzled Related Protein 2) ELISA Kit	SFRP2	Sandwich	80ng/mL	0.75 ng/mL	1.25-80ng/mL
FY-Ch1197	Chicken GNLY(Granulysin) ELISA Kit	GNLY	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1198	Chicken PLSCR1(Phospholipid Scramblase 1) ELISA Kit	PLSCR1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1199	Chicken Grb7(Growth Factor Receptor Bound Protein 7) ELISA Kit	Grb7	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1200	Chicken IL-10(Interleukin 10) ELISA Kit	IL-10	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1201	Chicken PC I (Procollagen I) ELISA Kit	PC I	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1202	Chicken VEGF-B(Vascular Endothelial Cell Growth Factor B) ELISA Kit	VEGF-B	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1203	Chicken RPL26(Ribosomal Protein L26) ELISA Kit	RPL26	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1204	Chicken PCDH γ A2(Protocadherin Gamma A2) ELISA Kit	PCDH γ A2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1205	Chicken FOXO1(Forkhead Box Protein O1) ELISA Kit	FOXO1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1206	Chicken CK-17/KRT17(Cytokeratin 17) ELISA Kit	CK-17/KRT17	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1207	Chicken GARS(Glycyl tRNA Synthetase) ELISA Kit	GARS	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1208	Chicken Grb2(Growth Factor Receptor Bound Protein 2) ELISA Kit	Grb2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1209	Chicken FAIM3(Fas Apoptotic Inhibitory Molecule 3) ELISA Kit	FAIM3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1210	Chicken TNFSF14(Tumor Necrosis Factor Ligand Superfamily, Member 14) ELISA Kit	TNFSF14	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1211	Chicken BCAN(Brevican) ELISA Kit	BCAN	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1212	Chicken SPT β N4(Spectrin Beta, Non Erythrocytic 4) ELISA Kit	SPT β N4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1213	Chicken MGP(Matrix Gla Protein) ELISA Kit	MGP	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch1214	Chicken PLAU/uPA(Urokinase-Type Plasminogen Activator) ELISA Kit	PLAU/uPA	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1215	Chicken SPJ(Surfactant Associated Protein J) ELISA Kit	SPJ	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1216	Chicken CYPB(Cyclophilin B) ELISA Kit	CYPB	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch1217	Chicken COL1 α 2(Collagen Type I Alpha 2) ELISA Kit	COL1 α 2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1218	Chicken ET-1(Endothelin 1) ELISA Kit	ET-1	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1219	Chicken COL4(Collagen Type IV) ELISA Kit	COL4	Sandwich	300ng/mL	2.81 ng/mL	4.69-300ng/mL
FY-Ch1220	Chicken LIPC(Lipase, Hepatic) ELISA Kit	LIPC	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1221	Chicken Grb10(Growth Factor Receptor Bound Protein 10) ELISA Kit	Grb10	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1222	Chicken DCLK1(Doublecortin Like Kinase 1) ELISA Kit	DCLK1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1223	Chicken IL-12 p40(Interleukin 12 p40) ELISA Kit	IL-12 p40	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1224	Chicken Trx(Thioredoxin) ELISA Kit	Trx	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1225	Chicken AGT(Angiotensinogen) ELISA Kit	AGT	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-Ch1226	Chicken PPP1R1A(Protein Phosphatase 1 Regulatory Subunit A) ELISA Kit	PPP1R1A	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1227	Chicken PIAS3(Protein Inhibitor of Activated STAT3) ELISA Kit	PIAS3	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1228	Chicken PDIA6(Protein Disulfide Isomerase A6) ELISA Kit	PDIA6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1229	Chicken ABCA1(ATP-binding cassette transporter A1) ELISA Kit	ABCA1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1230	Chicken PLC β 1/PI-PLC(Phospholipase C Beta 1, Phosphoinositide Specific) ELISA Kit	PLC β 1/PI-PLC	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1231	Chicken CDK1(Cyclin Dependent Kinase 1) ELISA Kit	CDK1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1232	Chicken ApoC4(Apolipoprotein C4) ELISA Kit	ApoC4	Sandwich	25pg/mL	0.23 pg/mL	0.39-25pg/mL
FY-Ch1233	Chicken FST(Follistatin) ELISA Kit	FST	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1234	Chicken TSH(Thyroid Stimulating Hormone) ELISA Kit	TSH	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1235	Chicken EPHA10(Ephrin Type A Receptor 10) ELISA Kit	EPHA10	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1236	Chicken PKC δ 1(Protein Kinase C Delta) ELISA Kit	PKC δ 1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-Ch1237	Chicken ADAMTS4(ADAM with Thrombospondin Type 1 Motif 4) ELISA Kit	ADAMTS4	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1238	Chicken EGR3(Early Growth Response Protein 3) ELISA Kit	EGR3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1239	Chicken ACHAP(Acetylcholinesterase Associated Protein) ELISA Kit	ACHAP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1240	Chicken AChE(Acetylcholinesterase) ELISA Kit	AChE	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1241	Chicken STAT1(Signal Transducer and Activator of Transcription 1) ELISA Kit	STAT1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1243	Chicken STAT3(Signal Transducer and Activator of Transcription 3) ELISA Kit	STAT3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1245	Chicken NUP133(Nucleoporin 133kDa) ELISA Kit	NUP133	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1246	Chicken ANXA2(Annexin A2) ELISA Kit	ANXA2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1247	Chicken TL1A(Tumor Necrosis Factor Related Ligand 1A) ELISA Kit	TL1A	Sandwich	250pg/ml	2.35 pg/ml	3.91-250pg/ml
FY-Ch1248	Chicken GLRX3(Glutaredoxin 3) ELISA Kit	GLRX3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1249	Chicken TIEG1(TGF Beta Inducible Early Response Gene 1) ELISA Kit	TIEG1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1250	Chicken OCLN(Occludin) ELISA Kit	OCLN	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1251	Chicken DMD(Dystrophin) ELISA Kit	DMD	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1252	Chicken C1qR1(Complement Component 1q Receptor 1) ELISA Kit	C1qR1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1253	Chicken S100A13(S100 Calcium Binding Protein A13) ELISA Kit	S100A13	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1254	Chicken FOXP3(Forkhead Box Protein P3) ELISA Kit	FOXP3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1255	Chicken AR(Amphiregulin) ELISA Kit	AR	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1256	Chicken LCA/CD45(Leukocyte Common Antigen) ELISA Kit	LCA/CD45	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1257	Chicken INSL5(Insulin Like Protein 5) ELISA Kit	INSL5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1258	Chicken G6PD(Glucose 6 Phosphate Dehydrogenase) ELISA Kit	G6PD	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch1259	Chicken HSP-70/HSPA9(Heat Shock Protein 70) ELISA Kit	HSP-70/HSPA9	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch1260	Chicken NPYR5(Neuropeptide Y Receptor Y5) ELISA Kit	NPYR5	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1261	Chicken CNN1(Calponin 1, Basic, Smooth Muscle) ELISA Kit	CNN1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1262	Chicken ATF6(Activating Transcription Factor 6) ELISA Kit	ATF6	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch1263	Chicken FTH(Ferritin, Heavy Polypeptide) ELISA Kit	FTH	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-Ch1264	Chicken ECE1(Endothelin Converting Enzyme 1) ELISA Kit	ECE1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1265	Chicken IFN- β (Interferon Beta) ELISA Kit	IFN- β	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1266	Chicken PKI α (Protein Kinase Inhibitor Alpha) ELISA Kit	PKI α	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1267	Chicken C3(Complement Component 3) ELISA Kit	C3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1268	Chicken Ep-CAM/CD326(Epithelial Cell Adhesion Molecule) ELISA Kit	Ep-CAM/CD326	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1269	Chicken ANG1(Angiopoietin 1) ELISA Kit	ANG1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1270	Chicken F10(Coagulation Factor X) ELISA Kit	F10	Sandwich	300pg/mL	2.81 pg/mL	4.69-300pg/mL
FY-Ch1271	Chicken Alk-Smase(Alkaline Sphingomyelinase) ELISA Kit	Alk-Smase	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1272	Chicken ANGPTL1(Angiopoietin Like Protein 1) ELISA Kit	ANGPTL1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1273	Chicken DDC(Dopamine Decarboxylase) ELISA Kit	DDC	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1274	Chicken CMLC-1(Cardiac Myosin Light Chains 1) ELISA Kit	CMLC-1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1275	Chicken C4(Complement Component 4) ELISA Kit	C4	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch1276	Chicken LMP7/PSMB8(Low Molecular-weight Protein/proteasome Beta-type Subunit) ELISA Kit	LMP7/PSMB8	Sandwich	25ng/mL	0.23 ng/mL	0.39-25ng/mL
FY-Ch1277	Chicken GSTpi(Glutathione S Transferases Pi) ELISA Kit	GSTpi	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1278	Chicken ELA3B(Elastase 3B) ELISA Kit	ELA3B	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1279	Chicken NXF1(Nuclear RNA Export Factor 1) ELISA Kit	NXF1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-Ch1280	Chicken PDI(Protein Disulfide Isomerase) ELISA Kit	PDI	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1281	Chicken NPYR2(Neuropeptide Y Receptor Y2) ELISA Kit	NPYR2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1282	Chicken PARP4(Poly ADP Ribose Polymerase 4) ELISA Kit	PARP4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1283	Chicken PS2(Presenilin 2) ELISA Kit	PS2	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1284	Chicken STMN1(Stathmin 1) ELISA Kit	STMN1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1285	Chicken ANGPT4(Angiopoietin 4) ELISA Kit	ANGPT4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1286	Chicken FASL/TNFSF6(Factor Related Apoptosis Ligand) ELISA Kit	FASL/TNFSF6	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch1287	Chicken CLEC2C(C-Type Lectin Domain Family 2 Member C) ELISA Kit	CLEC2C	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1288	Chicken MUSK(Muscle Skeletal Receptor Tyrosine Kinase) ELISA Kit	MUSK	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1289	Chicken TP63(Tumor Protein p63) ELISA Kit	TP63	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1290	Chicken PRSS1(Protease, Serine, 1) ELISA Kit	PRSS1	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch1291	Chicken ASPN(Asporin) ELISA Kit	ASPN	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1292	Chicken AQP-2(Aquaporin 2) ELISA Kit	AQP-2	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1293	Chicken SELP(P-Selectin) ELISA Kit	SELP	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1294	Chicken PLB(Phospholipase B) ELISA Kit	PLB	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1295	Chicken SIGLEC12(Sialic Acid Binding Ig Like Lectin 12) ELISA Kit	SIGLEC12	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1296	Chicken NOS1/nNOS(Nitric Oxide Synthase 1, Neuronal) ELISA Kit	NOS1/nNOS	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1298	Chicken CTSK(Cathepsin K) ELISA Kit	CTSK	Sandwich	100fmol/mL	0.94 fmol/mL	1.56-100fmol/mL
FY-Ch1299	Chicken COL5a2(Collagen Type V Alpha 2) ELISA Kit	COL5a2	Sandwich	25ng/mL	0.23 ng/mL	0.39-25ng/mL
FY-Ch1300	Chicken sAC(Soluble Adenylate Cyclase) ELISA Kit	sAC	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1301	Chicken LOX-1(Lectin Like Oxidized Low Density Lipoprotein Receptor 1) ELISA Kit	LOX-1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1303	Chicken CD163(Cluster of Differentiation 163) ELISA Kit	CD163	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1304	Chicken PAR2(Protease Activated Receptor 2) ELISA Kit	PAR2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch1305	Chicken UBE3A(Ubiquitin Protein Ligase E3A) ELISA Kit	UBE3A	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1306	Chicken PKIγ(Protein Kinase Inhibitor Gamma) ELISA Kit	PKIγ	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1307	Chicken FBLN5(Fibulin 5) ELISA Kit	FBLN5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1309	Chicken ELA4(Elastase 4) ELISA Kit	ELA4	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1310	Chicken EAR2(Eosinophil Associated Ribonuclease A family Member 2) ELISA Kit	EAR2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1311	Chicken AGGF1(Angiogenic Factor with G Patch and FHA Domains 1) ELISA Kit	AGGF1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1312	Chicken NT-ProBNP(N-Terminal Pro-Brain Natriuretic Peptide) ELISA Kit	NT-ProBNP	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1313	Chicken GP6(Glycoprotein VI, Platelet) ELISA Kit	GP6	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1314	Chicken PICK1(Protein Interacting With Protein Kinase C, Alpha 1) ELISA Kit	PICK1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1315	Chicken PKP2(Plakophilin 2) ELISA Kit	PKP2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1316	Chicken CLDN4(Claudin 4) ELISA Kit	CLDN4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1317	Chicken N-MID-OT(N-MID Osteocalcin) ELISA Kit	N-MID-OT	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1318	Chicken ASK1(Apoptosis Signal Regulating Kinase 1) ELISA Kit	ASK1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1319	Chicken GATA2(GATA Binding Protein 2) ELISA Kit	GATA2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1320	Chicken Flt3(FMS Like Tyrosine Kinase 3) ELISA Kit	Flt3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1321	Chicken BID(BH3 Interacting Domain Death Agonist) ELISA Kit	BID	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1322	Chicken HSF1(Heat Shock Factor 1) ELISA Kit	HSF1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1323	Chicken Rock-2(Rho Associated Coiled Coil Containing Protein Kinase 2) ELISA Kit	Rock-2	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL

FY-Ch1324	Chicken C7(Complement Component 7) ELISA Kit	C7	Sandwich	40pg/mL	0.38 pg/mL	0.63-40pg/mL
FY-Ch1325	Chicken MMP-9(Matrix Metalloproteinase 9) ELISA Kit	MMP-9	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1326	Chicken BAFFR(B-Cell Activation Factor Receptor) ELISA Kit	BAFFR	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1327	Chicken CNTNAP4(Contactin Associated Protein 4) ELISA Kit	CNTNAP4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1328	Chicken MFAP4(Microfibrillar Associated Protein 4) ELISA Kit	MFAP4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1329	Chicken PKC η (Protein Kinase C Eta) ELISA Kit	PKC η	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1330	Chicken SEMA3A(Semaphorin 3A) ELISA Kit	SEMA3A	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1331	Chicken LGALS3BP(Lectin Galactoside Binding, Soluble 3 Binding Protein) ELISA	LGALS3BP	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1332	Chicken CG(Cryoglobulin) ELISA Kit	CG	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1333	Chicken CYTH1(Cytohesin 1) ELISA Kit	CYTH1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1334	Chicken Myelin-IgA(Myelin IgA) ELISA Kit	Myelin-IgA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1335	Chicken EIPS2(Estrogen Induced Protein pS2) ELISA Kit	EIPS2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1336	Chicken ACP3/PAP(Prostatic Acid Phosphatase) ELISA Kit	ACP3/PAP	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1337	Chicken WARS(Tryptophanyl tRNA Synthetase) ELISA Kit	WARS	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1338	Chicken RNF4(Ring Finger Protein 4) ELISA Kit	RNF4	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1339	Chicken PR3(Proteinase 3) ELISA Kit	PR3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1340	Chicken GTC(General Transcription Complex) ELISA Kit	GTC	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1341	Chicken IL-19(Interleukin 19) ELISA Kit	IL-19	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1342	Chicken TSP-1(Thrombospondin-1) ELISA Kit	TSP-1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1343	Chicken PAPP2(Pappalysin 2) ELISA Kit	PAPP2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1344	Chicken FOXC1(Forkhead Box Protein C1) ELISA Kit	FOXC1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1345	Chicken HSP-90(Heat Shock Protein 90) ELISA Kit	HSP-90	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1346	Chicken PIR(Pirin) ELISA Kit	PIR	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1347	Chicken GSP(Glycosylated Serum Protein) ELISA Kit	GSP	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1349	Chicken OCT2(Octamer-binding protein 2) ELISA Kit	OCT2	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch1350	Chicken A β 1-42(Amyloid Beta 1-42) ELISA Kit	A β 1-42	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1351	Chicken PTX 3/TSG-14(Pentraxin 3) ELISA Kit	PTX 3/TSG-14	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1352	Chicken MUC5AC(Mucin 5 Subtype AC) ELISA Kit	MUC5AC	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1353	Chicken ODC(Ornithine Decarboxylase) ELISA Kit	ODC	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1354	Chicken TARC(Thymus Activation Regulated Chemokine) ELISA Kit	TARC	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1355	Chicken BECN1(Beclin 1) ELISA Kit	BECN1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1356	Chicken Fbn(Fibrinolysin) ELISA Kit	Fbn	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1357	Chicken NES1(Nesfatin 1) ELISA Kit	NES1	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch1358	Chicken UPP1(Uridine Phosphorylase 1) ELISA Kit	UPP1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1359	Chicken PKC β 1(Protein Kinase C Beta 1) ELISA Kit	PKC β 1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1360	Chicken TEP1(Telomerase Associated Protein 1) ELISA Kit	TEP1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1361	Chicken ADIPOR2(Adiponectin Receptor 2) ELISA Kit	ADIPOR2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch1362	Chicken GL α (Galactosidase Alpha) ELISA Kit	GL α	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1363	Chicken SNRPD1(Small Nuclear Ribonucleoprotein Polypeptide D1) ELISA Kit	SNRPD1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1364	Chicken GT-34(Gastrin 34) ELISA Kit	GT-34	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1365	Chicken ACTN2(Actinin Alpha 2) ELISA Kit	ACTN2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL

FY-Ch1366	Chicken CA4(Carbonic Anhydrase IV) ELISA Kit	CA4	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch1367	Chicken CTACK(Cutaneous T-cell Attracting Chemokine) ELISA Kit	CTACK	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1369	Chicken PIGR/SC(Polymeric Immunoglobulin Receptor/Membrane Secretory Compo PIGR/SC		Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1370	Chicken PLVAP(Plasmalemma Vesicle Associated Protein) ELISA Kit	PLVAP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1371	Chicken AMH(Anti-Mullerian Hormone) ELISA Kit	AMH	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1372	Chicken PAIg(Platelet Associated Immunoglobulins) ELISA Kit	PAIg	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1373	Chicken STK39(Serine/Threonine Kinase 39) ELISA Kit	STK39	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch1374	Chicken GHR(Growth Hormone Receptor) ELISA Kit	GHR	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1375	Chicken FGB(Fibrinogen Beta) ELISA Kit	FGB	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1376	Chicken HBβ(Hemoglobin Beta) ELISA Kit	HBβ	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1377	Chicken ADM(Adrenomedullin) ELISA Kit	ADM	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1378	Chicken SEMA3B(Semaphorin 3B) ELISA Kit	SEMA3B	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1379	Chicken PTPRN(Protein Tyrosine Phosphatase Receptor Type N) ELISA Kit	PTPRN	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1380	Chicken PUMA(P53 Upregulated Modulator of Apoptosis) ELISA Kit	PUMA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1381	Chicken DSP(Desmoplakin) ELISA Kit	DSP	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1382	Chicken LFA-1(Lymphocyte Function Associated Antigen 1) ELISA Kit	LFA-1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1383	Chicken TSC22(Transforming Growth Factor Beta Stimulated Protein Clone 22) ELISA Kit	TSC22	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1384	Chicken EL/LIPG(Endothelial Lipase) ELISA Kit	EL/LIPG	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1385	Chicken α2-AP(α2-Antiplasmin) ELISA Kit	α2-AP	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch1386	Chicken BGN(Biglycan) ELISA Kit	BGN	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1387	Chicken STAT6(Signal Transducer and Activator of Transcription 6) ELISA Kit	STAT6	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch1388	Chicken TPPP(Tubulin Polymerization Promoting Protein) ELISA Kit	TPPP	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1389	Chicken MMP-7(Matrix Metalloproteinase 7) ELISA Kit	MMP-7	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1390	Chicken PCT(Procalcitonin) ELISA Kit	PCT	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1391	Chicken MMP-8(Matrix Metalloproteinase 8) ELISA Kit	MMP-8	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1392	Chicken EGFR2(Epidermal Growth Factor Receptor 2) ELISA Kit	EGFR2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1393	Chicken PON1(Paraoxonase 1) ELISA Kit	PON1	Sandwich	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-Ch1394	Chicken TEF(Thyrotroph Embryonic Factor) ELISA Kit	TEF	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1395	Chicken NAGLU(N-Acetyl Alpha-D-Glucosaminidase) ELISA Kit	NAGLU	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1396	Chicken HIF-1α(Hypoxia Inducible Factor 1 Alpha) ELISA Kit	HIF-1α	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1397	Chicken ConA(Concanavalin A) ELISA Kit	ConA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1398	Chicken PTPRB(Protein Tyrosine Phosphatase Receptor Type B) ELISA Kit	PTPRB	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1399	Chicken CYP7A1(Cytochrome P450 7A1) ELISA Kit	CYP7A1	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-Ch1400	Chicken β2GP1-IgM(β2 Glycoprotein 1 IgM) ELISA Kit	β2GP1-IgM	Sandwich	5ng/mL	0.05 ng/mL	0.08-5ng/mL
FY-Ch1401	Chicken Cav-1(Caveolin-1) ELISA Kit	Cav-1	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1402	Chicken FGL1(Fibrinogen Like Protein 1) ELISA Kit	FGL1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1403	Chicken PA2G4(Proliferation Associated Protein 2G4) ELISA Kit	PA2G4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1404	Chicken ANGPTL2(Angiopoietin Like Protein 2) ELISA Kit	ANGPTL2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1405	Chicken TACR3(Tachykinin Receptor 3) ELISA Kit	TACR3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1406	Chicken CRLF1(Cytokine Receptor Like Factor 1) ELISA Kit	CRLF1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1407	Chicken CCKBR(Cholecystokinin B Receptor) ELISA Kit	CCKBR	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL

FY-Ch1408	Chicken MMRN1(Multimerin 1) ELISA Kit	MMRN1	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1409	Chicken ADCY1(Adenylate Cyclase 1, Brain) ELISA Kit	ADCY1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch1410	Chicken LRG1(Leucine Rich Alpha-2-Glycoprotein 1) ELISA Kit	LRG1	Sandwich	25ng/mL	0.23 ng/mL	0.39-25ng/mL
FY-Ch1411	Chicken Smad1/MADH1(Mothers Against Decapentaplegic Homolog 1) ELISA Kit	Smad1/MADH1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1412	Chicken I-TAC/CXCL11(Interferon Inducible T-cell Alpha Chemoattractant) ELISA Kit	I-TAC/CXCL11	Sandwich	3000pg/mL	28.13 pg/mL	46.88-3000pg/mL
FY-Ch1413	Chicken TNKS2(Tankyrase 2) ELISA Kit	TNKS2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1414	Chicken GPC1(Glypican 1) ELISA Kit	GPC1	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1415	Chicken ANXA5(Annexin A5) ELISA Kit	ANXA5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1416	Chicken MASP2(Mannan Associated Serine Protease 2) ELISA Kit	MASP2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1417	Chicken F2(Coagulation Factor II) ELISA Kit	F2	Sandwich	30ng/mL	0.28 ng/mL	0.47-30ng/mL
FY-Ch1418	Chicken CASP3(Caspase 3) ELISA Kit	CASP3	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1419	Chicken BRAK/CXCL14(Breast and Kidney Expressed Chemokine) ELISA Kit	BRAK/CXCL14	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1420	Chicken PKR1(Prokineticin Receptor 1) ELISA Kit	PKR1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1421	Chicken ADRP(Adipose Differentiation Related Protein) ELISA Kit	ADRP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1422	Chicken GAD1(Glutamate Decarboxylase 1) ELISA Kit	GAD1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1423	Chicken eIF3a(Eukaryotic Translation Initiation Factor 3a) ELISA Kit	eIF3a	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1424	Chicken Pax(Paxillin) ELISA Kit	Pax	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1425	Chicken FABP5(Fatty Acid Binding Protein 5, Epidermal) ELISA Kit	FABP5	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1426	Chicken MT1(Metallothionein 1) ELISA Kit	MT1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1427	Chicken CHRM2(Cholinergic Receptor, Muscarinic 2) ELISA Kit	CHRM2	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1428	Chicken SBEM(Small Breast Epithelial Mucin) ELISA Kit	SBEM	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1429	Chicken CCR1(CC-Chemokine Receptor 1) ELISA Kit	CCR1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1430	Chicken SPC(Pulmonary Surfactant Associated Protein C) ELISA Kit	SPC	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch1431	Chicken MAU(Microalbuminuria) ELISA Kit	MAU	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1432	Chicken NGB(Neuroglobin) ELISA Kit	NGB	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1433	Chicken DS(Dermatan Sulfate) ELISA Kit	DS	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1434	Chicken TIMP-4(Tissue Inhibitors of Metalloproteinase 4) ELISA Kit	TIMP-4	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1435	Chicken TNR(Tenascin R) ELISA Kit	TNR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1436	Chicken PCDHβ2(Protocadherin Beta 2) ELISA Kit	PCDHβ2	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1437	Chicken CK-20/KRT20(Cytokeratin 20) ELISA Kit	CK-20/KRT20	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1438	Chicken FOXP1(Forkhead Box Protein P1) ELISA Kit	FOXP1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1439	Chicken IL-24(Interleukin 24) ELISA Kit	IL-24	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1440	Chicken GADL1(Glutamate Decarboxylase Like Protein 1) ELISA Kit	GADL1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1441	Chicken FPA(Fibrinopeptide A) ELISA Kit	FPA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1442	Chicken MUC20(Mucin 20) ELISA Kit	MUC20	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1443	Chicken IIGFALS(Insulin Like Growth Factor Binding Protein, Acid Labile Subunit) ELISA Kit	IIGFALS	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1444	Chicken DEFα5(Defensin Alpha 5, Paneth Cell Specific) ELISA Kit	DEFα5	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch1445	Chicken RHOB(Ras Homolog Gene Family, Member B) ELISA Kit	RHOB	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1446	Chicken REG1β(Regenerating Islet Derived Protein 1 Beta) ELISA Kit	REG1β	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1447	Chicken CBX3(Chromobox Homolog 3) ELISA Kit	CBX3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1448	Chicken PDGFA(Platelet Derived Growth Factor Subunit A) ELISA Kit	PDGFA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-Ch1449	Chicken NUP107(Nucleoporin 107kDa) ELISA Kit	NUP107	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1450	Chicken GCLM(Glutamate Cysteine Ligase, Modifier Subunit) ELISA Kit	GCLM	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1451	Chicken TFR2(Transferrin Receptor 2) ELISA Kit	TFR2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1452	Chicken CENPI(Centromere Protein I) ELISA Kit	CENPI	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1453	Chicken VS(Versican/PG-M) ELISA Kit	VS	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1454	Chicken RELB(V-Rel Reticuloendotheliosis Viral Oncogene Homolog B) ELISA Kit	RELB	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1455	Chicken PCDH15(Protocadherin 15) ELISA Kit	PCDH15	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1456	Chicken NF-κB p65(Nuclear Factor Kappa B p65) ELISA Kit	NF-κB p65	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1457	Chicken MAP-2(Microtubule Associated Protein 2) ELISA Kit	MAP-2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1458	Chicken OSCAR(Osteoclast Associated Receptor) ELISA Kit	OSCAR	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1459	Chicken GALR1(Galanin Receptor 1) ELISA Kit	GALR1	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1460	Chicken BACE1(Beta-site APP Cleaving Enzyme 1) ELISA Kit	BACE1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1461	Chicken TAT(Thrombin/Antithrombin Complex) ELISA Kit	TAT	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1462	Chicken NGF(Nerve Growth Factor) ELISA Kit	NGF	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1463	Chicken PDPN(Podoplanin) ELISA Kit	PDPN	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1464	Chicken TNC(Tenascin C) ELISA Kit	TNC	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1465	Chicken CDKN2A(Cyclin Dependent Kinase Inhibitor 2A) ELISA Kit	CDKN2A	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1466	Chicken NF-κB p105(Nuclear factor NF-kappa-B p105 subunit) ELISA Kit	NF-κB p105	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1467	Chicken GAL3(Galectin 3) ELISA Kit	GAL3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1468	Chicken RNASEN(Ribonuclease III, Nuclear) ELISA Kit	RNASEN	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1470	Chicken RCOR3(Rest Corepressor 3) ELISA Kit	RCOR3	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1471	Chicken MCRS1(Microspherule Protein 1) ELISA Kit	MCRS1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1472	Chicken EGR4(Early Growth Response Protein 4) ELISA Kit	EGR4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1473	Chicken CRHBP(Corticotropin Releasing Hormone Binding Protein) ELISA Kit	CRHBP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1474	Chicken GSTα4(Glutathione S Transferase Alpha 4) ELISA Kit	GSTα4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1475	Chicken EP3(Prostaglandin E Receptor 3) ELISA Kit	EP3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1476	Chicken PXR(Pregnane X Receptor) ELISA Kit	PXR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1477	Chicken CDK4(Cyclin Dependent Kinase 4) ELISA Kit	CDK4	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1478	Chicken RNASE(Ribonuclease) ELISA Kit	RNASE	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1479	Chicken NTRK3(Neurotrophic Tyrosine Kinase Receptor Type 3) ELISA Kit	NTRK3	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1480	Chicken GSTμ2/GSTm2(Glutathione S Transferase Mu 2, Muscle) ELISA Kit	GSTμ2/GSTm2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1481	Chicken PDGF-AA(Platelet Derived Growth Factor-AA) ELISA Kit	PDGF-AA	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch1482	Chicken SNRPC(Small Nuclear Ribonucleoprotein Polypeptide C) ELISA Kit	SNRPC	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1483	Chicken AIRE(Autoimmune Regulator) ELISA Kit	AIRE	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1484	Chicken PTGS1(Prostaglandin Endoperoxide Synthase 1) ELISA Kit	PTGS1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1485	Chicken PYY(Peptide YY) ELISA Kit	PYY	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1486	Chicken CTSB(Cathepsin B) ELISA Kit	CTSB	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1487	Chicken PCNA(Proliferating Cell Nuclear Antigen) ELISA Kit	PCNA	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1488	Chicken ELN(Elastin) ELISA Kit	ELN	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch1489	Chicken PDIA5(Protein Disulfide Isomerase A5) ELISA Kit	PDIA5	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1490	Chicken FGF4(Fibroblast Growth Factor 4) ELISA Kit	FGF4	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL

FY-Ch1491	Chicken EFNA1(Ephrin A1) ELISA Kit	EFNA1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1492	Chicken CP(Ceruloplasmin) ELISA Kit	CP	Sandwich	4000ng/mL	37.5 ng/mL	62.5-4000ng/mL
FY-Ch1493	Chicken GNβ1(G Protein Beta 1) ELISA Kit	GNβ1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1494	Chicken CABIN1(Calcineurin Binding Protein 1) ELISA Kit	CABIN1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1496	Chicken PKD1(Protein Kinase D1) ELISA Kit	PKD1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1497	Chicken GAD2(Glutamate Decarboxylase 2) ELISA Kit	GAD2	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch1498	Chicken CYP1A2(Cytochrome P450, family 1, subfamily A, polypeptide 2) ELISA Kit	CYP1A2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1499	Chicken GFAP(Glial Fibrillary Acidic Protein) ELISA Kit	GFAP	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1500	Chicken ACAT2(Acetyl Coenzyme A Acetyltransferase 2) ELISA Kit	ACAT2	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1501	Chicken IGFBP-2(Insulin Like Growth Factor Binding Protein 2) ELISA Kit	IGFBP-2	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1502	Chicken SYNPO(Synaptopodin) ELISA Kit	SYNPO	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1503	Chicken BCL2L1/BCL-X(Bcl-2 Like Protein 1) ELISA Kit	BCL2L1/BCL-X	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1504	Chicken GS(Gelsolin) ELISA Kit	GS	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1505	Chicken PPOX(Protoporphyrinogen Oxidase) ELISA Kit	PPOX	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1506	Chicken APN(Aminopeptidase N) ELISA Kit	APN	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1507	Chicken CAMK2α(Calcium/Calmodulin-Dependent Protein Kinase II Alpha) ELISA Kit	CAMK2α	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1508	Chicken TGF-β2(Transforming Growth Factor β2) ELISA Kit	TGF-β2	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1509	Chicken PTGES(Prostaglandin E Synthase, Microsomal) ELISA Kit	PTGES	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1510	Chicken CCL16(Chemokine C-C-Motif Ligand 16) ELISA Kit	CCL16	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1511	Chicken GM-CSF(Granulocyte Macrophage Colony Stimulating Factor) ELISA Kit	GM-CSF	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1512	Chicken MDC/CCL22(Macrophage-Derived Chemokine) ELISA Kit	MDC/CCL22	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1513	Chicken GR(Glutathione Reductase) ELISA Kit	GR	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1514	Chicken βTG/PBP/CXCL7/NAP2(Thromboglobulin, Beta) ELISA Kit	βTG/PBP/CXCL7/NAP2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1515	Chicken FGFR1(Fibroblast Growth Factor Receptor 1) ELISA Kit	FGFR1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1516	Chicken AMY1(Amylase Alpha 1, Salivary) ELISA Kit	AMY1	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1517	Chicken MIA1(Melanoma Inhibitory Activity Protein 1) ELISA Kit	MIA1	Sandwich	25pg/mL	0.23 pg/mL	0.39-25pg/mL
FY-Ch1518	Chicken sIgA(Secretory Immunoglobulin A) ELISA Kit	sIgA	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1519	Chicken CR(Calretinin) ELISA Kit	CR	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1520	Chicken PF4V1(Platelet Factor 4 Variant 1) ELISA Kit	PF4V1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1521	Chicken AST(Aspartate Aminotransferase) ELISA Kit	AST	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1522	Chicken F13A1(Coagulation Factor X III A1 Polypeptide) ELISA Kit	F13A1	Sandwich	60pg/mL	0.56 pg/mL	0.94-60pg/mL
FY-Ch1523	Chicken NUP37(Nucleoporin 37kDa) ELISA Kit	NUP37	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1524	Chicken MFG-E8(Milk Fat Globule EGF Factor 8) ELISA Kit	MFG-E8	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1525	Chicken MSMβ(Microsminoprotein Beta) ELISA Kit	MSMβ	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1526	Chicken TPP1(Tripeptidyl Peptidase I) ELISA Kit	TPP1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1527	Chicken BMP-3(Bone Morphogenetic Protein 3) ELISA Kit	BMP-3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1529	Chicken MAPτ(Microtubule Associated Protein Tau/Tau Protein) ELISA Kit	MAPτ	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1530	Chicken PRKAα1(Protein Kinase, AMP Activated Alpha 1) ELISA Kit	PRKAα1	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1531	Chicken NPYR1(Neuropeptide Y Receptor Y1) ELISA Kit	NPYR1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1532	Chicken Rac-GAP1(Rac-GTPase Activating Protein 1) ELISA Kit	Rac-GAP1	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1533	Chicken SEMA7A(Semaphorin 7A) ELISA Kit	SEMA7A	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL

FY-Ch1534	Chicken FGFBP2(Fibroblast Growth Factor Binding Protein 2) ELISA Kit	FGFBP2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1535	Chicken ESF(Erythropoiesis Stimulating Factor) ELISA Kit	ESF	Sandwich	100pg/mL	0.94 pg/mL	1.56-100pg/mL
FY-Ch1536	Chicken CLL1(Collectin Liver 1) ELISA Kit	CLL1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1537	Chicken TGM1(Transglutaminase 1, Keratinocyte) ELISA Kit	TGM1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1538	Chicken PRSS8(Protease, Serine 8) ELISA Kit	PRSS8	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1539	Chicken PDIA2(Protein Disulfide Isomerase A2) ELISA Kit	PDIA2	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1540	Chicken ERK3(Extracellular Signal Regulated Kinase 3) ELISA Kit	ERK3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1541	Chicken MT3(Metallothionein 3) ELISA Kit	MT3	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1542	Chicken C4b(Complement Component 4b) ELISA Kit	C4b	Sandwich	15ng/mL	0.14 ng/mL	0.23-15ng/mL
FY-Ch1543	Chicken LCN1(Lipocalin 1) ELISA Kit	LCN1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch1544	Chicken ATF7(Activating Transcription Factor 7) ELISA Kit	ATF7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1545	Chicken PI3K2 α (Phosphoinositide-3-Kinase Class-2-Alpha Polypeptide) ELISA Kit	PI3K2 α	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1546	Chicken GJ β 1(Gap Junction Protein Beta 1) ELISA Kit	GJ β 1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1547	Chicken CYFRA21-1(Cytokeratin Fragment Antigen 21-1) ELISA Kit	CYFRA21-1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1548	Chicken GH(Growth Hormone) ELISA Kit	GH	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1549	Chicken RETN(Resistin) ELISA Kit	RETN	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1550	Chicken SSTR1(Somatostatin Receptor 1) ELISA Kit	SSTR1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1551	Chicken IL-17A(Interleukin 17A) ELISA Kit	IL-17A	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1552	Chicken FGF21(Fibroblast Growth Factor 21) ELISA Kit	FGF21	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1553	Chicken MSP(Macrophage Stimulating Protein) ELISA Kit	MSP	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1554	Chicken CENPH(Centromere Protein H) ELISA Kit	CENPH	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1555	Chicken SOX18(Sex Determining Region Y Box Protein 18) ELISA Kit	SOX18	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1556	Chicken PC4(Positive Cofactor 4) ELISA Kit	PC4	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1557	Chicken SYP(Synaptophysin) ELISA Kit	SYP	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch1558	Chicken CA12(Carbonic Anhydrase XII) ELISA Kit	CA12	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1559	Chicken MMP-12(Matrix Metalloproteinase 12) ELISA Kit	MMP-12	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1560	Chicken TGF β 1(Transforming Growth Factor Beta Induced Protein) ELISA Kit	TGF β 1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1561	Chicken S100A11(S100 Calcium Binding Protein A11) ELISA Kit	S100A11	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1562	Chicken EMILIN 1(Elastin Microfibril Interface Located Protein 1) ELISA Kit	EMILIN 1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1563	Chicken PRSS12(Protease, Serine 12) ELISA Kit	PRSS12	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1565	Chicken Cys-C(Cystatin C) ELISA Kit	Cys-C	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1566	Chicken NUP98(Nucleoporin 98kDa) ELISA Kit	NUP98	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1567	Chicken CNTNAP2(Contactin Associated Protein Like Protein 2) ELISA Kit	CNTNAP2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1568	Chicken ACP1(Acid Phosphatase 1) ELISA Kit	ACP1	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1569	Chicken F13B(Coagulation Factor X III B Polypeptide) ELISA Kit	F13B	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1570	Chicken ASM(Acid Sphingomyelinase) ELISA Kit	ASM	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1571	Chicken IGFBP-4(Insulin Like Growth Factor Binding Protein 4) ELISA Kit	IGFBP-4	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1572	Chicken PP1L1(Peptidylprolyl Isomerase/Cyclophilin Like 1) ELISA Kit	PP1L1	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch1573	Chicken KLK7(Kallikrein 7) ELISA Kit	KLK7	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1574	Chicken CA6(Carbonic Anhydrase VI) ELISA Kit	CA6	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1575	Chicken UCP-2(Uncoupling Protein 2, Mitochondrial) ELISA Kit	UCP-2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL

FY-Ch1576	Chicken GMP-140(Alpha-granular Membrane Protein) ELISA Kit	GMP-140	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1577	Chicken MBP(Myelin Basic Protein) ELISA Kit	MBP	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1578	Chicken TM(Thrombin) ELISA Kit	TM	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1579	Chicken TAF13(TATA Box Binding Protein Associated Factor 13) ELISA Kit	TAF13	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1580	Chicken CASP10(Caspase 10) ELISA Kit	CASP10	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch1581	Chicken GIP(Gastric Inhibitory Polypeptide) ELISA Kit	GIP	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch1582	Chicken MART/Melan-A(Melanoma Marker A) ELISA Kit	MART/Melan-A	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1583	Chicken TGF-β3(Transforming Growth Factor β3) ELISA Kit	TGF-β3	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1584	Chicken CA9(Carbonic Anhydrase IX) ELISA Kit	CA9	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1585	Chicken ITLN1(Intelectin 1/Omentin) ELISA Kit	ITLN1	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1586	Chicken CX37(Connexin 37) ELISA Kit	CX37	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1587	Chicken GBP4(Guanylate Binding Protein 4) ELISA Kit	GBP4	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1588	Chicken TM(Thrombomodulin) ELISA Kit	TM	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1589	Chicken MAG(Myelin Associated Glycoprotein) ELISA Kit	MAG	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1590	Chicken CX31(Connexin 31) ELISA Kit	CX31	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch1591	Chicken Bcl-2(B-cell Leukemia/Lymphoma 2) ELISA Kit	Bcl-2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1592	Chicken GPR120(G Protein Coupled Receptor 120) ELISA Kit	GPR120	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1593	Chicken TSTA(Tumor Specific Transplantation Antigen) ELISA Kit	TSTA	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1594	Chicken MT1M(Metallothionein 1M) ELISA Kit	MT1M	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1595	Chicken ES(Endostatin) ELISA Kit	ES	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1596	Chicken TAFI(Thrombin Activatable Fibrinolysis Inhibitor) ELISA Kit	TAFI	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1597	Chicken HO1(Heme Oxygenase 1) ELISA Kit	HO1	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1598	Chicken LDHB(Lactate Dehydrogenase B) ELISA Kit	LDHB	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1599	Chicken IGF1R(Insulin Like Growth Factor 1 Receptor) ELISA Kit	IGF1R	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1600	Chicken MGST1(Microsomal Glutathione S Transferase 1) ELISA Kit	MGST1	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch1601	Chicken EG-VEGF(Endocrine Gland Derived Vascular Endothelial Growth Factor) ELISA Kit	EG-VEGF	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1602	Chicken FABP6(Fatty Acid Binding Protein 6, Ileal) ELISA Kit	FABP6	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1603	Chicken RPA1(Replication protein A1) ELISA Kit	RPA1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1605	Chicken CLS(Cardiolipin Synthase) ELISA Kit	CLS	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch1606	Chicken FBN3(Fibrillin 3) ELISA Kit	FBN3	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1607	Chicken CTSD(Cathepsin D) ELISA Kit	CTSD	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1610	Chicken ATXN2(Ataxin 2) ELISA Kit	ATXN2	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1611	Chicken CNTNAP1(Contactin Associated Protein 1) ELISA Kit	CNTNAP1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1612	Chicken FGF23(Fibroblast Growth Factor 23) ELISA Kit	FGF23	Sandwich	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1613	Chicken PAP(Plasmin-Antiplasmin Complex) ELISA Kit	PAP	Sandwich	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch1614	Chicken FRS3(Fibroblast Growth Factor Receptor Substrate 3) ELISA Kit	FRS3	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1615	Chicken SDC1(Syndecan 1) ELISA Kit	SDC1	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1616	Chicken CD42(Cluster of Differentiation 42) ELISA Kit	CD42	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1617	Chicken SELL(L-Selectin) ELISA Kit	SELL	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1618	Chicken MUC5B(Mucin 5 Subtype B) ELISA Kit	MUC5B	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1619	Chicken LP-a(Lipoprotein a) ELISA Kit	LP-a	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL

FY-Ch1620	Chicken FAS/CD95(Factor Related Apoptosis) ELISA Kit	FAS/CD95	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1621	Chicken OSM(Oncostatin M) ELISA Kit	OSM	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1622	Chicken PS(Plasma Anticoagulant Protein S) ELISA Kit	PS	Sandwich	50ng/mL	0.47 ng/mL	0.78-50ng/mL
FY-Ch1623	Chicken TARS(Threonyl tRNA Synthetase, cytoplasmic) ELISA Kit	TARS	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1624	Chicken NFKBIE(Nuclear Factor Of Kappa Light Polypeptide Gene Enhancer In B-cells) ELISA Kit	NFKBIE	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1625	Chicken C8g(Complement Component 8g) ELISA Kit	C8g	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1626	Chicken VEGFR-2/KDR(Vascular Endothelial Growth Factor Receptor 2) ELISA Kit	VEGFR-2/KDR	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1627	Chicken ANGPTL3(Angiopoietin Like Protein 3) ELISA Kit	ANGPTL3	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1628	Chicken CGα(Chorionic Gonadotropin Alpha Polypeptide) ELISA Kit	CGα	Sandwich	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1629	Chicken KERA(Keratocan) ELISA Kit	KERA	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1630	Chicken AIF(Apoptosis Inducing Factor) ELISA Kit	AIF	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1631	Chicken β Manase(β Mannosidase) ELISA Kit	β Manase	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1632	Chicken TSLP(Thymic Stromal Lymphopoietin) ELISA Kit	TSLP	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1633	Chicken ADCY4(Adenylate Cyclase 4) ELISA Kit	ADCY4	Sandwich	40ng/mL	0.38 ng/mL	0.63-40ng/mL
FY-Ch1634	Chicken PTPRK(Protein Tyrosine Phosphatase Receptor Type K) ELISA Kit	PTPRK	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1635	Chicken NRG-1(Neuregulin 1) ELISA Kit	NRG-1	Sandwich	400ng/mL	3.75 ng/mL	6.25-400ng/mL
FY-Ch1637	Chicken IDUα(Alpha-L-Iduronidase) ELISA Kit	IDUα	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1638	Chicken α2-M(Alpha-2 Macroglobulin) ELISA Kit	α2-M	Sandwich	6000ng/mL	56.25 ng/mL	93.75-6000ng/mL
FY-Ch1639	Chicken F8(Coagulation Factor VIII) ELISA Kit	F8	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1640	Chicken ST2(Syntenin 2) ELISA Kit	ST2	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1641	Chicken cTnT/TNNT2(Troponin T Type 2, Cardiac) ELISA Kit	cTnT/TNNT2	Sandwich	1500pg/mL	14.06 pg/mL	23.44-1500pg/mL
FY-Ch1642	Chicken CDK5(Cyclin Dependent Kinase 5) ELISA Kit	CDK5	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1643	Chicken MUC21(Mucin 21) ELISA Kit	MUC21	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1644	Chicken MAdCAM1(Mucosal Addressin Cell Adhesion Molecule 1) ELISA Kit	MAdCAM1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1645	Chicken AXIN2(Axis Inhibition Protein 2) ELISA Kit	AXIN2	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1646	Chicken IBSP(Bone Sialoprotein 2) ELISA Kit	IBSP	Sandwich	500ng/mL	4.69 ng/mL	7.81-500ng/mL
FY-Ch1647	Chicken PI(Proinsulin) ELISA Kit	PI	Sandwich	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch1648	Chicken COMP(Cartilage Oligomeric Matrix Protein) ELISA Kit	COMP	Sandwich	1000ng/mL	9.38 ng/mL	15.63-1000ng/mL
FY-Ch1649	Chicken GUCA2B(Guanylate Cyclase Activator 2B) ELISA Kit	GUCA2B	Sandwich	50pg/mL	0.47 pg/mL	0.78-50pg/mL
FY-Ch1650	Chicken MR(Mannose Receptor) ELISA Kit	MR	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1651	Chicken D-LDH(D-Lactate Dehydrogenase) ELISA Kit	D-LDH	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1652	Chicken IGFBP-7(Insulin Like Growth Factor Binding Protein 7) ELISA Kit	IGFBP-7	Sandwich	60ng/mL	0.56 ng/mL	0.94-60ng/mL
FY-Ch1654	Chicken OLFM4(Olfactomedin 4) ELISA Kit	OLFM4	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1655	Chicken C5a(Complement Component 5a) ELISA Kit	C5a	Sandwich	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1656	Chicken HB(Hemoglobin) ELISA Kit	HB	Sandwich	50μg/mL	0.47 μg/mL	0.78-50μg/mL
FY-Ch1657	Chicken ERα(Estrogen Receptor Alpha) ELISA Kit	ERα	Sandwich	200pg/mL	1.88 pg/mL	3.13-200pg/mL
FY-Ch1658	Chicken NLR(NOD-Like Receptor9) ELISA Kit	NLR	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1659	Chicken τPK1(Tau-Protein Kinase 1) ELISA Kit	τPK1	Sandwich	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1660	Chicken BTA(Bladder Tumor Antigen) ELISA Kit	BTA	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1661	Chicken Ta1(Thymosin-α1) ELISA Kit	Ta1	Sandwich	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1662	Chicken PS1(Presenilin 1) ELISA Kit	PS1	Sandwich	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL

FY-Ch1663	Chicken PACAP-38(Pituitary Adenylate Cyclase Activating Polypeptide 38) ELISA Kit	PACAP-38	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1664	Chicken PACAP-27(Pituitary Adenylate Cyclase Activating Polypeptide 27) ELISA Kit	PACAP-27	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1665	Chicken INS(Insulin) ELISA Kit	INS	Sandwich	20ng/mL	0.19 ng/mL	0.31-20ng/mL
FY-Ch1667	Chicken SI(Sucrase Isomaltase) ELISA Kit	SI	Sandwich	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0010	Chicken MT(Melatonin) ELISA Kit	MT	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0087	Chicken ANP(Atrial Natriuretic Peptide) ELISA Kit	ANP	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0116	Chicken GLP-1(Glucagon Like Peptide 1) ELISA Kit	GLP-1	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0132	Chicken δ SIP- α (Delta Sleep Inducing Peptide α) ELISA Kit	δ SIP- α	Competitive	400pg/mL	3.75 pg/mL	6.25-400pg/mL
FY-Ch0146	Chicken α -CTx(Alpha Crosslaps) ELISA Kit	α -CTx	Competitive	6000pg/mL	56.25 pg/mL	93.75-6000pg/mL
FY-Ch0181	Chicken MPA(Mycophenolic Acid) ELISA Kit	MPA	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0205	Chicken 6-OHDA(6-Hydroxydopamine) ELISA Kit	6-OHDA	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0212	Chicken Ang- I (Angiotensin I) ELISA Kit	Ang- I	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0286	Chicken AP36(Apelin 36) ELISA Kit	AP36	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0314	Chicken NKB(Neurokinin B) ELISA Kit	NKB	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-Ch0332	Chicken PGF2 α (Prostaglandin F2 α) ELISA Kit	PGF2 α	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch0341	Chicken TXB2(Thromboxane B2) ELISA Kit	TXB2	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0354	Chicken Dyn(Big Dynorphin) ELISA Kit	Dyn	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-Ch0457	Chicken MS(Melatonin Sulfate) ELISA Kit	MS	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0473	Chicken β -EP(Beta-Endorphin) ELISA Kit	β -EP	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0514	Chicken CNP(C-type Natriuretic Peptide) ELISA Kit	CNP	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0539	Chicken TMS β 10(Thymosin Beta 10) ELISA Kit	TMS β 10	Competitive	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0659	Chicken TAP(Trypsinogen Activation Peptide) ELISA Kit	TAP	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0674	Chicken TMS β 4(Thymosin Beta 4) ELISA Kit	TMS β 4	Competitive	250ng/mL	2.35 ng/mL	3.91-250ng/mL
FY-Ch0681	Chicken BNP(Brain Natriuretic Peptide) ELISA Kit	BNP	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0707	Chicken CRH(Corticotropin Releasing Hormone) ELISA Kit	CRH	Competitive	200ng/mL	1.88 ng/mL	3.13-200ng/mL
FY-Ch0740	Chicken SOD3(Superoxide Dismutase 3, Extracellular) ELISA Kit	SOD3	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0761	Chicken NPY(Neuropeptide Y) ELISA Kit	NPY	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0769	Chicken ACTH(Adrenocorticotropic Hormone) ELISA Kit	ACTH	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch0775	Chicken 8-epi-PGF2 α (8-Epi-Prostaglandin F2 Alpha) ELISA Kit	8-epi-PGF2 α	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0820	Chicken NADPH(Nicotinamide Adenine Dinucleotide Phosphate) ELISA Kit	NADPH	Competitive	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch0828	Chicken β -CTx(Beta Crosslaps) ELISA Kit	β -CTx	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0831	Chicken AP13(Apelin 13) ELISA Kit	AP13	Competitive	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch0877	Chicken OGP(Osteogenic Growth Peptide) ELISA Kit	OGP	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0892	Chicken Ang- II (Angiotensin II) ELISA Kit	Ang- II	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch0918	Chicken SST(Somatostatin) ELISA Kit	SST	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch0933	Chicken GLP-2(Glucagon Like Peptide 2) ELISA Kit	GLP-2	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch0967	Chicken SA(Sialic Acid) ELISA Kit	SA	Competitive	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1018	Chicken MTL(Motilin) ELISA Kit	MTL	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1054	Chicken DPD(Deoxyypyridinoline) ELISA Kit	DPD	Competitive	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch1083	Chicken VMA(Vanillylmandelic Acid) ELISA Kit	VMA	Competitive	2000ng/mL	18.75 ng/mL	31.25-2000ng/mL
FY-Ch1110	Chicken SCT(Secretin) ELISA Kit	SCT	Competitive	100ng/mL	0.94 ng/mL	1.56-100ng/mL

FY-Ch1119	Chicken NKA(Neurokinin A) ELISA Kit	NKA	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-Ch1131	Chicken TP5(Thymopentin) ELISA Kit	TP5	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1141	Chicken DHEA-S(Dehydroepiandrosterone Sulfate) ELISA Kit	DHEA-S	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1166	Chicken OXB(Orexin B) ELISA Kit	OXB	Competitive	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1176	Chicken SOD1(Superoxide Dismutase 1, Soluble) ELISA Kit	SOD1	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1187	Chicken EM2(Endomorphin 2) ELISA Kit	EM2	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-Ch1242	Chicken APLN(Apelin) ELISA Kit	APLN	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1244	Chicken OXA(Orexin A) ELISA Kit	OXA	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1297	Chicken GC(Glucagon) ELISA Kit	GC	Competitive	2000pg/mL	18.75 pg/mL	31.25-2000pg/mL
FY-Ch1302	Chicken 11-DH-TXB2(11-Dehydro-thromboxane B2) ELISA Kit	11-DH-TXB2	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1308	Chicken BK(Bradykinin) ELISA Kit	BK	Competitive	800ng/mL	7.5 ng/mL	12.5-800ng/mL
FY-Ch1348	Chicken VIP(Vasoactive Intestinal Peptide) ELISA Kit	VIP	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-Ch1368	Chicken SOD2(Superoxide Dismutase 2, Mitochondrial) ELISA Kit	SOD2	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1469	Chicken IAPP(Islet Amyloid Polypeptide) ELISA Kit	IAPP	Competitive	4000pg/mL	37.5 pg/mL	62.5-4000pg/mL
FY-Ch1495	Chicken AP17(Apelin 17) ELISA Kit	AP17	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1528	Chicken P II NP(Procollagen II N-Terminal ProPeptide) ELISA Kit	P II NP	Competitive	100ng/mL	0.94 ng/mL	1.56-100ng/mL
FY-Ch1564	Chicken α MSH(Alpha-Melanocyte Stimulating Hormone) ELISA Kit	α MSH	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1604	Chicken TRH(Thyrotropin Releasing Hormone) ELISA Kit	TRH	Competitive	5000pg/mL	46.88 pg/mL	78.13-5000pg/mL
FY-Ch1608	Chicken CTX-III(Cross Linked C-telopeptide of Type III Collagen) ELISA Kit	CTX-III	Competitive	10ng/mL	0.1 ng/mL	0.16-10ng/mL
FY-Ch1609	Chicken NT(Neurotensin) ELISA Kit	NT	Competitive	250pg/mL	2.35 pg/mL	3.91-250pg/mL
FY-Ch1636	Chicken OFQ/N(Orphanin FQ/Nociceptin) ELISA Kit	OFQ/N	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL
FY-Ch1653	Chicken NP(Nandrolone Phenylpropionate) ELISA Kit	NP	Competitive	500pg/mL	4.69 pg/mL	7.81-500pg/mL
FY-Ch1666	Chicken HPX(Hemopexin) ELISA Kit	HPX	Competitive	1000pg/mL	9.38 pg/mL	15.63-1000pg/mL